

NGS Single Cell RNA Analysis

NEBNext® SINGLE CELL/LOW INPUT LIBRARY PREP KIT

*How low
can you go?*



be **INSPIRED**
drive **DISCOVERY**
stay **GENUINE**


NEBNext® Single Cell/Low Input RNA Library Prep

High-quality transcript sequencing from a single cell or ultra-low input RNA

This unique workflow meets the demand for a highly sensitive, yet robust method that consistently generates high-quality, full-length transcript sequencing data from a single cell or ultra-low input RNA.

Optimized cDNA synthesis and amplification steps incorporate template switching, as well as utilize a unique protocol and suite of reagents. Even low-abundance transcripts are represented in the high yields of cDNA obtained. Subsequent library construction incorporates the Ultra™ II FS enzymatic DNA fragmentation/end repair/dA-tailing mix in a simple and efficient workflow.

NEBNext Single Cell/Low Input RNA Library Prep workflow



Advantages

- Generate the highest yields of high-quality full-length transcript sequencing libraries from single cells, or as little as 2 pg – 200 ng total RNA
- Experience unmatched detection of low abundance transcripts
- Rely on consistent transcript detection for a wide range of input amounts and sample types
- Obtain full length, uniform transcript coverage, regardless of input amount or sample type
- Use with cultured or primary cells, or total RNA
- Save time with a fast, streamlined workflow, minimal handling steps and hands-on time
 - Single-tube protocol from cell lysis to cDNA
 - Enzymatic DNA fragmentation, end repair and dA-tailing reagents in a single enzyme mix, with a single protocol, regardless of GC content




Information on automation can be found here:
www.neb.com/automation



To learn more about the technology, please visit:
www.neb.com/nebnextsinglecell
and download the extensive technical note.




The NEBNext Single Cell/Low Input RNA Library Prep Kit provides uniform coverage across the length of transcripts



Sequencing libraries were generated from HeLa, Jurkat and M1 single cells, or 10 pg of Universal Human Reference (UHR) RNA (Agilent® #740000) with recommended amounts of ERCC RNA Spike-In Mix I (Thermo Fisher Scientific® #4456740). The NEBNext Single Cell/Low Input RNA Library Prep Kit, or the SMART-Seq v4 Ultra® Low Input RNA Kit for Sequencing (Clontech® # 634891) plus the Nextera XT DNA Library Prep Kit (Illumina® #FC-131-1096) were used. Libraries were sequenced on an Illumina NextSeq® 500 using paired-end mode (2x76 bp). Gene body coverage shown is an average of four replicates and was calculated using Picard tools. The global view of the 5' to 3' coverage of the RefSeq transcripts reveals both consistency across different sample types and uniformity across the transcript length in the NEBNext libraries.


The NEBNext Single Cell/Low Input RNA Library Prep Kit increases transcript detection

Sequencing libraries were generated from Jurkat single cells (6 replicates) using the NEBNext Single Cell/Low Input RNA Library Prep Kit, or the SMART-Seq v4 Ultra® Low Input RNA Kit for Sequencing (Clontech® # 634891) plus the Nextera XT DNA Library Prep Kit (Illumina® #FC-131-1096). Libraries were sequenced on an Illumina NextSeq® 500 using paired-end mode (2x76 bp). TPM = Transcripts per Kilobase Million. Each dot represents the number of transcripts identified at the given TPM range, and each box represents the median, first and third quartiles per replicate. Salmon 0.6 was used for read mapping and quantification of all GENCODE v25 transcripts. Panels show the number of transcripts detected within the following TPM ranges: 1-5, 5-10, 10-50 and >50 TPM. Increased identification of low abundance transcripts is observed with the NEBNext libraries.




Generate higher library yields with the NEBNext Single Cell/Low Input RNA Library Prep Kit

Sequencing libraries were generated from HeLa, Jurkat and M1 single cells or 10 pg of Universal Human Reference (UHR) RNA (Agilent® #740000) with recommended amounts of ERCC RNA Spike-In Mix I (Thermo Fisher Scientific® #4456740). The NEBNext Single Cell/Low Input RNA Library Prep Kit, or the SMART-Seq v4 Ultra® Low Input RNA Kit for Sequencing (Clontech® #634891) plus the Nextera XT DNA Library Prep Kit (Illumina® #FC-131-1096) were used. For the NEBNext workflow ~80% of the cDNA was used as input into sequencing library preparation, and libraries were amplified with 8 PCR cycles. For the SMART-Seq v4/Nextera XT workflow, as recommended, 125 pg of cDNA was used as input into sequencing library preparation and 12 PCR cycles were used for amplification. Error bars indicate standard deviation for 6-11 replicates.



Fast NEBNext Single Cell/Low Input RNA Library Prep workflow with only 26 min hands-on time!



What users are saying:

“The implementation of this new single cell RNA-seq pipeline is an important addition to our core Scientific Operations and will enable our scientists to resolve at scale the transcriptional variation within a wide variety of single cells.

– Sarah Teichmann, PhD, head of cellular genetics at the Wellcome Sanger Institute

“NEB libraries show greater consistency between cells across all metrics studied (number of reads, mapping rate, 5'-3' coverage, genome features, and transcript detection).

– Technical Note Fluidigm C1™ <https://www.fluidigm.com/c1openapp/scripthub/script/2018-07/nebnext-mrna-sequencing-1530810776075-7>



To request a free sample, please visit www.neb.com/nebnextsinglecell

ORDERING INFORMATION

| PRODUCTS | NEB # | SIZE |
|--|----------|------------|
| NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina | E6420S/L | 24/96 rxns |
| NEBNext Single Cell/ Low Input cDNA Synthesis & Amplification Module | E6421S/L | 24/96 rxns |

| ALSO AVAILABLE | NEB # | SIZE |
|---|-------------------------------|-------------|
| NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs) | E6440S/L | 96/384 rxns |
| NEBNext Multiplex Oligos for Illumina Set 2 (96 Unique Dual Index Primer Pairs) | E6442S/L | 96/384 rxns |
| NEBNext Multiplex Oligos for Illumina Set 3 (96 Unique Dual Index Primer Pairs) | E6444S/L | 96/384 rxns |
| NEBNext Multiplex Oligos for Illumina Set 4 (96 Unique Dual Index Primer Pairs) | E6446S/L | 96/384 rxns |
| NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 1) | E7600S | 96 rxns |
| NEBNext Multiplex Oligos for Illumina (Dual Index Set 2) | E7780S | 96 rxns |
| NEBNext Multiplex Oligos for Illumina (96 Index Primers) | E6609S/L | 96/384 rxns |
| NEBNext Multiplex Oligos for Illumina (Index Primers Set 1, 2, 3, 4) | E7335, E7500, E7710, E7730S/L | 24/96 rxns |



NEBNext Selector

Use this tool to guide you through selection of NEBNext reagents for next generation sequencing sample preparation.

One or more of these products are covered by patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at gbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications. ILLUMINA®, NEXT-SEQ® and NEXTERA® are registered trademarks of Illumina, Inc. SMART-SEQ® and ULTRA® are registered trademarks of Takara Bio, Inc. Fluidigm C1™ is a trademark of Fluidigm Corporation.

© Copyright 2022, New England Biolabs, Inc.; all rights reserved.

Your local NEB distributor:

| | |
|---|--|
| <p>BeNeLux: BIOKE Tel: (+31) 71 720 0220 (BE: 0800 - 71640) info@bioke.com www.bioke.com</p> <p>BULGARIA ELTA 90M Ltd. Tel: (+359) 2 983 96 49 Fax: (+359) 2983 2211, 9274 232 office@elta90.eu www.elta90.com</p> <p>CROATIA: Diagnosticska Skopeli d.o.o. Tel: (+385) 1 377 8484 Fax: (+385) 1 377 8585 skopeli@skopeli.hr www.skopeli.hr</p> <p>CZECH REPUBLIC: BioTech a.s. Tel: (+420) 2727 017 39 info@ibiotech.cz www.ibiotech.cz</p> <p>DENMARK: BioNordika Denmark A/S Tel: (+45) 39 56 20 00 Fax: (+45) 39 56 19 42 info@bionordika.dk www.bionordika.dk</p> <p>FINLAND: BioNordika Oy Tel: (+358) 207 410 270 info@bionordika.fi www.bionordika.fi</p> <p>GREECE: BIOLINE SCIENTIFIC Tel: (+30) 210 522 6547 Fax: (+30) 210 524 4744 info@bioline.gr www.bioline.gr</p> <p>HUNGARY: Kvalitex Kft. Tel: (+361) 340 4700 Fax: (+361) 339 8274 info@kvalitex.hu www.kvalitex.hu</p> <p>ITALY: Euroclone S.p.A. Tel: (+39) 02 38 1951 info@euroclone.it www.euroclonengroup.it</p> <p>ISRAEL: Ornat BioChemicals & Lab Equipment Tel: (+972) 8 947 7077 Fax: (+972) 8 936 3034 ornatbio@ornat.co.il www.ornat.co.il</p> <p>LITHUANIA: Nanodiagnostika, Ltd. Tel: (+370) 525 052 44 lina@nanodiagnostika.lt www.nanodiagnostika.lt</p> <p>Norway: BioNordika AS Tel: (+47) 2303 5800 info@bionordika.no www.bionordika.no</p> | <p>POLAND: Lab-JOT Ltd. Sp. z o.o. Sp.k. Tel: (+48) 2233 598 84 Fax: (+48) 2233 598 19 biuro@labjot.com www.labjot.com</p> <p>PORTUGAL: Werfen Tel: (+351) 214 247 300 dac-pt@werfen.com www.pt.werfen.com</p> <p>ROMANIA: BioZyme Tel: (+40) 264 523 281 contact@biozyme.ro www.biozyme.ro</p> <p>RUSSIA and CIS: SkyGen LLC Tel: (+88) 0033 312 26 info@skygen.com www.skygen.com</p> <p>SERBIA: Alfa Genetics d.o.o. Tel: (+38) 1011 2084 252 office@alfagenetics.rs www.alfagenetics.rs</p> <p>SPAIN: Werfen Tel: (+34) 900 822 290 customerservice-es@werfen.com www.es.werfen.com</p> <p>SLOVAK REPUBLIC BioTech s.r.o. Tel: (+421) 254 7744 88 Fax: (+421) 232 2010 48 info@ibiotech.sk www.ibiotech.sk</p> <p>SLOVENIA: Mikro+Polo d.o.o. Tel: (+386) 2614 33 00 Fax: (+386) 2614 33 20 podpora@mikro-polo.si www.mikro-polo.si</p> <p>South Africa: Inqaba Biotechnical Industries (Pty) Ltd Tel: (+27) 1234 358 29 info@inqababiotec.co.za www.inqababiotec.co.za</p> <p>SWEDEN, ESTONIA, LATVIA, LITHUANIA: BioNordika Sweden AB Tel: (+46) 830 60 10 info@bionordika.se www.bionordika.se</p> <p>SWITZERLAND: BioConcept Ltd. Tel: (+41) 6148 680 80 Fax: (+41) 6148 680 00 info@bioconcept.ch www.bioconcept.ch</p> <p>TURKEY: EKA BIOLAB TEKNOLOJİ Tel: (+90) 312 248 13 24 Email: info@ekabiolab.com www.ekabiolab.com</p> |
|---|--|

Distributors

For a complete list of all NEB subsidiaries and distributors, please go to:

www.neb.com/international-support

