

# One *Taq* DNA Polymerase


ENZYMES & KITS FOR PCR



# OneTaq DNA Polymerase

The *One* you've been waiting for!

Choose OneTaq DNA Polymerase for all your amplification of standard, AT- and GC-rich templates, at a price-point that beats the competition. OneTaq DNA Polymerase is an optimized blend of Taq and Deep Vent DNA polymerases for use with routine and difficult PCR experiments. The 3'–5' exonuclease activity of Deep Vent DNA Polymerase increases the fidelity and robust amplification of Taq DNA Polymerase. The OneTaq Reaction Buffers and High GC Enhancer have been formulated for robust yields with minimal optimization, regardless of a template's GC content (see additional data on page 4).



## Achieve robust amplification for standard, AT- and GC-rich templates with OneTaq DNA Polymerase

Amplification of a selection of sequences with varying AT and GC content from human and *C. elegans* genomic DNA using OneTaq DNA Polymerase. GC content is indicated above gel. Marker M is the 1 kb DNA Ladder (NEB #N3232).

## It's that easy to choose the right buffer for maximum performance:

OneTaq DNA Polymerase is supplied with two 5× buffers (Standard and GC), as well as a High GC Enhancer solution to ensure maximum performance for routine, AT- or GC-rich amplicons.

AMPLICON % GC CONTENT	RECOMMENDED DEFAULT BUFFER	OPTIMIZATION NOTES
<50% GC	OneTaq Standard Reaction Buffer	Adjust annealing temperature, primer/template concentration, etc., if needed.
50–65% GC	OneTaq Standard Reaction Buffer	OneTaq GC Reaction Buffer can be used to enhance performance of difficult amplicons.
>65% GC	OneTaq GC Reaction Buffer	OneTaq GC Reaction Buffer with 10–20% OneTaq High GC Enhancer can be used to enhance performance of difficult amplicons.



## Scan the QR code

to read our guidelines for PCR optimization with OneTaq DNA Polymerases.

## ADVANTAGES

- Exceptional performance in endpoint PCR across a wide range of templates (Standard, AT- and GC-rich)
- Robust yields with minimal optimization
- Convenient product formats (stand-alone enzyme, master mixes, and Quick-Load formats)

## DETAILS & APPLICATIONS:

### Details

Extension Rate	1 kb/min
Amplicon Size	≤ 6 kb
Fidelity	2X Taq
Units/50 µl rxn	1.25 units
Resulting Ends	3' A/Blunt
3'→5' Exonuclease Activity	Yes
5'→3' Exonuclease Activity	Yes
Supplied Buffer	OneTaq Std Rxn Buffer, OneTaq GC Rxn Buffer
Supplied Enhancer	OneTaq High GC Enhancer

### Product Formats

Hot Start Available	Yes
- Activation Required	No
Master Mix Available	Yes
Direct Gel-loading Available	Yes
PCR Kit Available	No

### Applications


Routine PCR	Yes
SNP Detection	Yes
T/A, U/A Cloning	Yes
Colony PCR	Yes
High-Fidelity PCR	No

# OneTaq Quick-Load DNA Polymerase

Optimal for “standard” PCR and direct gel-loading

For direct and fast agarose gel-loading after “standard” PCRs such as genotyping or colony PCR etc., OneTaq DNA Polymerase is also available in a Quick-Load format. It is supplied with a density and tracking dye containing 5× OneTaq Quick-Load Reaction Buffer for direct gel loading in addition to the regular “colorless” 5× OneTaq Reaction Buffer.

## Convenient direct gel-loading feature of OneTaq Quick-Load DNA Polymerase is not compromising on performance



Amplification of a variety of DNA targets demonstrates strong performance of the OneTaq Quick-Load DNA Polymerase. Product yield (bars, left axis) and purity (circles, right axis) were calculated via microfluidic analysis from triplicate reactions after 30 cycles of PCR. Standard deviation is indicated by error bars (yield) or shaded bands (purity). GoTaq was cycled according to manufacturer's recommendations.

### TIPS

#### Use your OneTaq PCR products in Sanger sequencing?

The dye in OneTaq Quick-Load 2X Master Mix Buffer doesn't interfere with Sanger sequencing. Prepare your samples with the fast and easy Exo-CIP Rapid PCR Cleanup Kit (#E1050) and proceed directly to Sanger sequencing.

## Convenient Master Mix Formulations

Faster set-up, less pipetting, less errors & contaminations

OneTaq and OneTaq Hot Start DNA Polymerases are also available in convenient Master Mix and Quick-Load Master Mix formats. Master Mix formulations include dNTPs, MgCl<sub>2</sub> and other buffers and stabilizers. The Quick-Load Master Mix formats allow for direct gel-loading. OneTaq Master Mix formulations are optimally suited for faster reaction set-up with less pipetting steps – they increase the reliability of each PCR and reduce the risk of contamination.

ONE TAQ 2X MASTER MIXES	ART.-NR.	CONVENIENT	GEL-LOADING
with Standard Buffer	M0482S/L	•	
OneTaq Quick-Load 2× Master Mix with Standard Buffer	M0486S/L/X	•	•

ONE TAQ HOT START 2X MASTER MIXES	ART.-NR.	CONVENIENT	GEL-LOADING
with Standard Buffer	M0484S/L	•	
with GC Buffer	M0485S/L	•	
OneTaq Hot Start Quick-Load 2× Master Mix with Standard Buffer	M0488S/L	•	•
OneTaq Hot Start Quick-Load 2× Master Mix with GC Buffer	M0489S/L	•	•

# Robust PCR results you can trust!

Unparalleled robustness on all tested templates!

While OneTaq DNA Polymerase works perfectly well on AT-rich DNA (see page 2), it's the challenging templates that unveil the true performance and quality of any PCR polymerase. Therefore, OneTaq DNA Polymerase and its Hot Start counterpart have been stringently and systematically tested under selected demanding conditions. In the figure below, we demonstrate some examples of the convincing performance of OneTaq on GC-rich templates.

Choose OneTaq DNA Polymerase for all your templates regardless of GC-content to benefit from the unrivaled robustness and reliability, so you can always trust your results!

## RECOMMENDATION

- Use OneTaq DNA Polymerase for all your templates
- Benefit from unrivaled reliability and performance
- Obtain PCR results you can really trust
- Choose the OneTaq format (e.g. convenient Master Mix) that suits you best
- Benefit from our low and fair prices

## Better than the competition:

Polymerase	Additives	GC%					
		55	65	66	67	78	79
OneTaq <sup>®</sup> DNA Polymerase (NEB)	None*						
OneTaq HotStart DNA Polymerase (NEB)	None*						
AmpliTaq Gold <sup>™</sup> 360 DNA Polymerase (Thermo Fisher)	None 360 GC Enhancer						
DreamTaq <sup>™</sup> Hot Start DNA Polymerase (Thermo Fisher)	(Not provided)						
FastStart <sup>™</sup> Taq DNA Polymerase (Roche)	None GC-RICH solution						
GoTaq <sup>®</sup> G2 Hot Start Polymerase (Promega)	(Not provided)						
GoTaq Hot Start Polymerase (Promega)	(Not provided)						
HotStarTaq <sup>®</sup> DNA Polymerase (Qiagen)	Q-Solution None						
HotStarTaq Plus DNA Polymerase (Qiagen)	Q-Solution None						
iTaq <sup>™</sup> DNA polymerase (Bio-Rad)	(Not provided)						
JumpStart <sup>™</sup> Taq DNA Polymerase (Sigma)	None						
Platinum <sup>™</sup> II Taq Hot-Start DNA Polymerase (Thermo Fisher)	Platinum GC Enhancer None						
Platinum Taq DNA Polymerase High Fidelity (Thermo Fisher)	(Not provided)						
Platinum Taq DNA Polymerase (Thermo Fisher)	None KB Extender						
Ex Taq DNA Polymerase, hot-start version (TaKaRa)	(Not provided)						
Titanium <sup>®</sup> Taq DNA Polymerase (TaKaRa)	(Not provided)						

Yield (ng/ul)

· 0.0 · 1.0 · 2.0 · 3.0 · 4.0 · ≥5.0

% Purity

0 100

\* OneTaq products are supplied with OneTaq Standard Reaction Buffer and OneTaq GC Reaction Buffer. The GC reaction buffer was used to amplify the targets shown in the table. For other products, amplification reactions were conducted both with and without GC enhancers (if provided).

Amplification of a selection of high GC human genomic DNA targets demonstrates OneTaq performance. All polymerases were cycled according to manufacturer's recommendations, including the use of additives to enhance the amplification of targets with high GC content. Yield (dot size) and purity (color) of reaction product were quantified from triplicate reactions on a Perkin Elmer LabChip. A large, dark green dot represents the highest yield and purity.

# OneTaq Hot Start DNA Polymerase

## Room temperature reaction setup with no activation step

NEB's OneTaq Hot Start utilizes aptamer technology. This unique modified oligonucleotide binds to the polymerase through non-covalent interactions, blocking polymerase activity at temperatures below 45°C. The polymerase is activated during normal cycling conditions, allowing reactions to be set up at room temperature. OneTaq Hot Start DNA Polymerase does not require a separate high temperature incubation step to activate the enzyme. This ultimately shortens reaction times and increases ease of use.


### Recommended time for enzyme activation of commercially available Hot Start Taq products

MANUFACTURER	ENZYME	ACTIVATION STEP*	HOT START FORM
Thermo Fisher	Ampli Taq Gold 360	10', 95°C	Modified
Thermo Fisher	Platinum Taq	30"-2', 94°C	Ab
Promega	Go Taq Hot Start	2', 94-95°C	Ab
Qiagen	HotStar Taq	15', 95°C	Modified
Roche	FastStart Taq	4', 95°C	Modified
Sigma	JumpStart Taq	1', 94°C	Ab
Thermo Fisher	Thermo-Start Taq	15', 95°C	Modified
NEB	OneTaq	None	Aptamer

\* May include initial denaturation step

### ADVANTAGES

- Allows room temperature reaction setup
- Does not require a separate activation step
- Compatible with standard Taq protocols



### Extended room temperature incubation does not affect performance or specificity of OneTaq Hot Start DNA Polymerase

Amplification of a selection of sequences with varying GC content from human and *C.elegans* genomic DNA using OneTaq Hot Start DNA Polymerase. The presence or absence of an extended room temperature incubation does not affect performance. GC content is indicated under gel. Marker M is the 1 kb DNA Ladder (#N3232).



### PCR Selector

Use this tool to help select the right DNA polymerase for your PCR setup. Whether your amplicon is long, complex, GC-rich or present in a single copy, the PCR selection tool will identify the perfect DNA polymerase for your reaction.

[PCRselector.neb.com](http://PCRselector.neb.com)



### PCR Fidelity Estimator

Estimate the percentage of correct DNA copies (those without base substitution errors) per cycle of PCR for selected DNA polymerases.

[PCRFidelityEstimator.neb.com](http://PCRFidelityEstimator.neb.com)



GERMANY & AUSTRIA  
New England Biolabs GmbH  
Brüningstr. 50, Geb B852  
65926 Frankfurt/Main, Germany  
Tel: +49/(0)69/305-23140  
Fax: +49/(0)69/305-23149

Free Call: 0800/246 5227 (Germany)  
Free Call: 00800/246 52277 (Austria)  
info.de@neb.com

**www.neb-online.de**

#### HEADQUARTERS: USA

New England Biolabs, Inc.  
Telephone: (978) 927-5054

Toll Free (USA Orders): 1-800-632-5227  
Toll Free (USA Tech): 1-800-632-7799  
Fax: (978) 921-1350  
info@neb.com

**www.neb.com**

FRANCE  
New England Biolabs France SAS  
Genopole Campus 1, Bâtiment 6  
5 rue Henri Desbrères  
91030 Evry cedex, France  
Tel.: 0800 100 632 (Customer Service)

Tel.: 0800 100 633 (Technical Service)  
info.fr@neb.com

**www.neb-online.fr**

## Ordering Information

PRODUCTS	NEB #	SIZE
OneTaq DNA Polymerase	M0480S/L/X	200/1,000/ 5,000 units
OneTaq Quick-Load DNA Polymerase	M0509L/X	500/ 2,500 units
OneTaq 2X Master Mix with Standard Buffer	M0482S/L	100/500 rxns
OneTaq Quick-Load 2x Master Mix with Standard Buffer	M0486S/L/X	100/500/ 2,500 rxns
OneTaq Hot Start DNA Polymerase	M0481S/L/X	200/1,000/ 5,000 units
OneTaq Hot Start 2x Master Mix with Standard Buffer	M0484S/L	100/500 rxns
OneTaq Hot Start 2x Master Mix with GC Buffer	M0485S/L	100/500 rxns
OneTaq Hot Start Quick-Load 2x Master Mix with Standard Buffer	M0488S/L	100/500 rxns
OneTaq Hot Start Quick-Load 2x Master Mix with GC Buffer	M0489S/L	100/500 rxns
OneTaq RT-PCR Kit	E5310S	30 rxns
OneTaq One-Step RT-PCR Kit	E5315S	30 rxns

Please ask for larger packing sizes or quantities: info.de@neb.com. Purchase of this product provides the purchaser with a non-exclusive license to use OneTaq DNA Polymerase products for research purposes only.

One or more of these products are covered by patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at gbo@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications.

Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at [www.neb.com/support/terms-of-sale](http://www.neb.com/support/terms-of-sale). NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

AMPLITAQ GOLD, PLATINUM, THERMO-START, DREAMTAQ, are trademarks of Thermo Fisher Scientific. FASTSTART is a trademark of Roche. EX TAQ, TITANIUM are trademarks of TaKaRa Bio. GO TAQ is a registered trademark of Promega. HOTSTAR is a trademark of Qiagen. ITAQ is a trademark of BIO-RAD. JUMPSTART is a trademark of Sigma-Aldrich.

## Your local NEB distributor:

AZERBAIJAN &  
TAJIKISTAN:  
Bayramli Group LLC  
info@bayramligroup.com  
www.bayramligroup.com

BeNeLux:  
BIOKÉ  
order@bioke.com  
www.bioke.com

BULGARIA:  
ELTA 90M Ltd.  
office@elta90.eu  
www.elta90.com

CROATIA:  
Diagnostica Skalpeli d.o.o.  
skalpel@skalpel.hr  
www.skapel.hr

CZECH REPUBLIC:  
BioTech a.s.  
info@ibitech.cz  
www.ibitech.cz

DENMARK:  
BioNordika Denmark A/S  
info@bionordika.dk  
www.bionordika.dk

FINLAND:  
BioNordika Oy  
info@bionordika.fi  
www.bionordika.fi

GEORGIA:  
Oxgen-Import, LLC  
info@oxgenimport.ge  
www.oxgenimport.ge

GREECE:  
BIOLINE SCIENTIFIC  
info@bioline.gr  
www.bioline.gr

HUNGARY:  
Kvalitex Kft.  
info@kvalitex.hu  
www.kvalitex.hu

ITALY:  
EUROCLONE S.P.A.  
info@euroclone.it  
www.euroclonegroup.it

ISRAEL:  
Ornat BioChemicals &  
Lab Equipment  
ornatbio@ornat.co.il  
www.ornat.co.il

KAZAKHSTAN:  
INOS, LLC  
market@inos.kz  
www.inos.kz

LITHUANIA:  
Nanodiagnostika, Ltd.  
info@nanodiagnostika.lt  
www.nanodiagnostika.lt

NORWAY:  
BioNordika AS  
info@bionordika.no  
www.bionordika.no

POLAND:  
Lab-JOT Ltd. Sp.z o.o. Sp.k.  
biuro@labjot.com  
www.labjot.com

PORTUGAL:  
Werfen  
dac-pt@werfen.com  
www.pt.werfen.com

ROMANIA:  
BioZyme  
contact@biozyme.ro  
www.biozyme.ro

SERBIA:  
Alfa Genetics d.o.o.  
office@alfagenetics.rs  
www.alfagenetics.rs

SPAIN:  
Werfen  
customerservice-es@werfen.com  
www.es.werfen.com

SLOVAK REPUBLIC:  
BioTech s.r.o.  
info@ibitech.sk  
www.ibitech.sk

SLOVENIA:  
Mikro+Polo d.o.o.  
info@mikro-polo.si  
www.mikro-polo.si

SOUTH AFRICA:  
Inqaba Biotechnical Industries  
(Pty) Ltd  
info@inqababiotec.co.za  
www.inqababiotec.co.za

SWEDEN, ESTONIA, LATVIA,  
LITHUANIA:  
BIONORDIKA SWEDEN AB  
info@bionordika.se  
www.bionordika.se

SWITZERLAND:  
BioConcept Ltd.  
info@bioconcept.ch  
www.bioconcept.ch

TURKEY:  
EKA BIOLAB TEKNOLOJİ  
info@ekabiolab.com  
www.ekabiolab.com

## Distributors

For a list of NEB subsidiaries and distributors, please go to:

**www.neb.com/international-support**



## Tm Calculator

Use this tool when designing PCR reaction protocols to help determine the optimal annealing temperature for your amplicon. Simply input your DNA polymerase, primer concentration and your primer sequence and the Tm Calculator will guide you to successful reaction conditions.

**TmCalculator.neb.com**



Print product with financial  
**climate contribution**  
ClimatePartner.com/13845-2502-1005