


NEBNext UltraExpress[®] Library Prep Kits


STREAMLINED FOR SPEED



NEBNext UltraExpress Streamlined for speed

Sometimes speed is required to set you apart from the pack. The NEBNext UltraExpress DNA and RNA Library Prep Kits have been carefully optimized for speed, while providing the high yields and high quality that you've come to rely on. Each kit has a single-condition protocol, with fixed universal adaptor concentration and number of PCR cycles, for the ultimate in streamlining. With fewer workflow and cleanup steps and automation friendly transfer volumes, the kits were built for ease of use and automation compatibility. And, they do this all while reducing consumables waste, making your discoveries at the bench greener.

Savings* with NEBNext UltraExpress



* As compared to NEBNext Ultra II

“ Protocols designed to save time usually only save 20 - 30 minutes. The NEBNext UltraExpress RNA Library Prep protocol saved us just over an hour in processing time. This is quite significant. We were especially impressed with the new cleanup approach that resulted in squeaky clean libraries. ”

~ Director, Large Sequencing Core

5 Reasons to choose NEBNext UltraExpress

- 1 Fast workflow
- 2 High quality libraries from a wide input range
- 3 Single protocol for all input amounts
- 4 Fewer steps, consumables, and cleanups
- 5 Automation friendly

TOOLS & RESOURCES



Visit [NEBNext.com](https://www.neb.com) to find:

- The full list of products available
- Video protocols & Workflow animations
- Online tutorials to help with product selection, general handling tips and more
- Access to **NEBNext Selector Tool**, our online tool for help with selecting the right NEBNext product
- Over 30,000 NEBNext citations
- Protocols & FAQs
- Refer to NEBNext Automation Compatibility Selection Chart for details about successful kit automation

NEBNext UltraExpress DNA Library Prep Kit

The NEBNext UltraExpress DNA Library Prep Kit features a fast, streamlined workflow that enables the creation of high yield, high quality libraries from a broad input amount range in under 2 hours. The workflow incorporates master mixed reagents, reduced incubation times, and fewer cleanup steps, all conducted in a single tube, which minimizes plastic consumable waste. With simplicity at the forefront, the kit includes a single protocol for all input amounts, making it well-suited for all throughputs.

DNA Library Prep Workflow




The NEBNext UltraExpress DNA Library Prep Kit provides robust library yields over a wide input range

Libraries were prepared from 10, 50, 100 or 200 ng of Human NA19240 genomic DNA (Coriell Institute for Medical Research) using the same adaptor amount and 8 PCR cycles. Yields exceeded the minimum requirement (40 ng) for a single Illumina® NovaSeq® 6000 run to achieve whole genome sequencing with at least 30X coverage.

ADVANTAGES

- Fast workflow (< 2 hours)
- Fewer steps and consumables
- Fewer cleanups
- Single protocol for all input amounts
- High quality libraries from a wide input range (10 - 200 ng pre-sheared DNA)
- Single protocol for all input amounts
- Automation friendly



The NEBNext UltraExpress DNA Library Prep Kit produces libraries with uniform GC coverage and insert size from a range of input amounts

Libraries were prepared from 10, 50, 100 or 200 ng of Human NA19240 genomic DNA (Coriell Institute for Medical Research) using the same adaptor amount and 8 PCR cycles. Libraries were pooled and sequenced on an Illumina® MiSeq® (2 x 75 bases). Data showed consistent GC coverage (A) and insert size (B). 2 million paired-end reads from each library were sampled (seqtk v1.3), adapter-trimmed (seqprep v0.1) and mapped to the GRCh38 reference (bowtie2 v2.4.5), and GC coverage information was calculated using Picard's CollectGCBiasMetrics (v1.56.0). In (A), the horizontal grey line indicates the expected normalized coverage of 1.0, and the dots in shades of green represent read numbers at each GC%. The grey area plot is a histogram representing the distribution of GC content in 100 bp windows of the reference genome.



Need help with product selection?



If you need assistance navigating the ever-growing NEBNext product line, NEBNext Selector Tool can help guide your product selection.

Visit nebnextselector.neb.com

NEBNext UltraExpress FS DNA Library Prep Kit

The NEBNext UltraExpress FS DNA Library Prep Kit offers a fast, streamlined workflow, including enzymatic fragmentation, end prep, and dA-tailing with a single enzyme mix. In under 2 hours, the protocol enables the creation of high yield, high quality libraries, from a broad input amount range, while generating less plastic waste. Designed for simplicity and efficiency, the kit features a single protocol for all input amounts, making NEBNext UltraExpress FS DNA Library Prep well suited for all throughputs.

FS DNA Library Prep Workflow




The NEBNext UltraExpress FS DNA Library Prep Kit provides robust library yields over a wide input range

Libraries were prepared in triplicate from 10, 50, 100 and 200 ng of a 9:1 Human NA19240 genomic DNA (Coriell Institute for Medical Research) and *Escherichia coli* gDNA (Lostrand Labs Limited) mixed sample, using the NEBNext UltraExpress FS DNA single-protocol workflow (e.g., same adaptor amount and 6 PCR cycles for all input amounts). Yields exceeded the minimum requirement (40 ng) for a single Illumina® NovaSeq® 6000 run to achieve whole genome sequencing with at least 30X coverage.

ADVANTAGES

- Fast workflow (< 2 hours), including FS "Fragmentation System" enzymatic fragmentation
- Fewer steps, requiring only a single reaction tube
- Fewer cleanup steps required
- High quality libraries from a wide input range (10 - 200 ng intact DNA)
- Single protocol for all input amounts
- Automation friendly



The NEBNext UltraExpress FS DNA Library Prep Kit produces libraries with uniform GC coverage and insert size from a range of input amounts

Libraries were prepared from 10, 50, 100 and 200 ng of a 9:1 Human NA19240 genomic DNA (Coriell Institute for Medical Research) and *Escherichia coli* gDNA (Lostrand Labs Limited) mixed sample, using the NEBNext UltraExpress FS DNA single-protocol workflow (e.g., same adaptor amount and 6 PCR cycles for all input amounts). Libraries were pooled and sequenced on an Illumina NextSeq® 500/550 (2 x 75 bases). Data showed consistent (A) GC coverage and (B) insert size. 2 million paired-end reads from each library were sampled (seqtk v1.0), adaptor-trimmed (seqprep v0.1) and mapped to a composite reference containing GRCh38 and *E. coli* MG1655 contigs (bowtie2 v2.5.0). GC coverage and insert size distributions were calculated using Picard's CollectGCBiasMetrics and Picard CollectInsertSizeMetrics (v1.56.0). Picard CollectGCBiasMetrics (v1.56.0) was run on human autosomes only due to the even copy number assumption of the tool. In (A), the horizontal grey line indicates the expected normalized coverage of 1.0, and the dots in shades of green represent read numbers at each GC%. The grey area plot is a histogram representing the distribution of GC content in 100 bp windows of the reference genome.




Did you know that many NEB products can be purchased in large volumes and custom formats?

Learn more at www.neb.com/customizedsolutions


NEBNext UltraExpress RNA Library Prep Kit

The NEBNext UltraExpress RNA Library Prep Kit is the latest generation of NEBNext RNA library prep, with a fast, streamlined workflow. The kit is compatible with mRNA isolation and rRNA depletion workflows and a wide range of sample types. With a 3 hour library prep protocol, the kit enables creation of high quality RNA libraries in a single day, in conjunction with mRNA or rRNA depletion kits.

RNA Library Prep Workflow



The NEBNext UltraExpress RNA Library Prep Kit produces high library yields for a range of inputs, in poly(A) enrichment and rRNA depletion workflows.



Universal human reference RNA with the indicated input amounts was (A) enriched for poly(A) mRNA (NEB #E7490) or (B) depleted of ribosomal RNA (NEBNext rRNA Depletion Kit v2 (Human/Mouse/Rat) (NEB #E7400), followed by creation of strand-specific libraries using the NEBNext UltraExpress RNA Library Preparation Kit. Library yields were assessed using TapeStation® 4200 and values shown are for three replicates for each input amount.

ADVANTAGES

- Fast workflow (3 hours)
- Fewer steps and consumables
- Fewer cleanups
- Single protocol for all input amounts
- High quality directional libraries from a broad input range: 25 - 250 ng total RNA
- Automation friendly

NEBNext Poly(A) mRNA Magnetic Isolation Module

- Isolation of intact poly(A)+ RNA from 1 - 5 µg of previously isolated total RNA
- Protocol can be adapted for automated high-throughput applications
- Permits elution of intact mRNA in small volumes

Fast workflow: streamlined express protocol of only 30 minutes

Note that for input amounts >5 µg we recommend the NEBNext High Input Poly(A) mRNA Isolation Module.

NEBNext rRNA Depletion Kit v2 (Human/Mouse/Rat)

- Compatible with a broad range of input amounts of 10 ng - 1 µg
- Superior depletion of rRNA from human, mouse and rat RNA
- Depletes cytoplasmic (5S, 5.8S, 18S, 28S, human ITS, ETS) and mitochondrial (12S, 16S) rRNA
- Fast workflow: 2 hours, with less than 10 minutes hands-on time

GERMANY & AUSTRIA
New England Biolabs GmbH
Brüningstr. 50, Geb B852
65926 Frankfurt/Main, Germany
Tel: +49/(0)69/305-23140
Fax: +49/(0)69/305-23149

Free Call: 0800/246 5227 (Germany)
Free Call: 00800/246 52277 (Austria)
info.de@neb.com

www.neb-online.de

HEADQUARTERS: USA

New England Biolabs, Inc.
Telephone: (978) 927-5054

Toll Free (USA Orders): 1-800-632-5227
Toll Free (USA Tech): 1-800-632-7799
Fax: (978) 921-1350
info@neb.com

www.neb.com

FRANCE
New England Biolabs France SAS
Genopole Campus 1, Bâtiment 6
5 rue Henri Desbruères
91030 Evry cedex, France
Tel.: 0800 100 632 (Customer Service)

Tel.: 0800 100 633 (Technical Service)
info.fr@neb.com

www.neb-online.fr

Your local NEB distributor:

AZERBAIJAN &
TAJIKISTAN:
Bayramli Group LLC
info@bayramligroup.com
www.bayramligroup.com

BeNeLux:
BIOKÉ
order@bioke.com
www.bioke.com

BULGARIA:
ELTA 90M Ltd.
office@elta90.eu
www.elta90.com

CROATIA:
Diagnostica Skalpel d.o.o.
skalpel@skalpel.hr
www.skalpel.hr

CZECH REPUBLIC:
BioTech a.s.
info@ibiotech.cz
www.ibiotech.cz

DENMARK:
BioNordika Denmark A/S
info@bionordika.dk
www.bionordika.dk

FINLAND:
BioNordika Oy
info@bionordika.fi
www.bionordika.fi

GEORGIA:
Oxgen-Import, LLC
info@oxgenimport.ge
www.oxgenimport.ge

GREECE:
BIOLINE SCIENTIFIC
info@bioline.gr
www.bioline.gr

HUNGARY:
Kvalitex Kft.
info@kvalitex.hu
www.kvalitex.hu

ITALY:
EUROCLONE S.P.A.
info@euroclone.it
www.euroclonegroup.it

ISRAEL:
Ornat BioChemicals &
Lab Equipment
ornatbio@ornat.co.il
www.ornat.co.il

KAZAKHSTAN:
INOS, LLC
market@inos.kz
www.inos.kz

LITHUANIA:
Nanodiagnostika, Ltd.
info@nanodiagnostika.lt
www.nanodiagnostika.lt

NORWAY:
BioNordika AS
info@bionordika.no
www.bionordika.no

POLAND:
Lab-JOT Ltd. Sp.z o.o. Sp.k.
biuro@labjot.com
www.labjot.com

PORTUGAL:
Werfen
dac-pt@werfen.com
www.pt.werfen.com

ROMANIA:
BioZyme
contact@biozyme.ro
www.biozyme.ro

SERBIA:
Alfa Genetics d.o.o.
office@alfagenetics.rs
www.alfagenetics.rs

SPAIN:
Werfen
customerservice-es@werfen.com
www.es.werfen.com

SLOVAK REPUBLIC:
BioTech s.r.o.
info@ibiotech.sk
www.ibiotech.sk

SLOVENIA:
Mikro+Polo d.o.o.
info@mikro-polo.si
www.mikro-polo.si

SOUTH AFRICA:
Inqaba Biotechnical Industries
(Pty) Ltd
info@inqababiotec.co.za
www.inqababiotec.co.za

SWEDEN, ESTONIA, LATVIA,
LITHUANIA:
BIONORDIKA SWEDEN AB
info@bionordika.se
www.bionordika.se

SWITZERLAND:
BioConcept Ltd.
info@bioconcept.ch
www.bioconcept.ch

TURKEY:
EKA BIOLAB TEKNOLOJİ
info@ekabiolab.com
www.ekabiolab.com

Ordering Information

DNA	NEB #	SIZE
NEBNext UltraExpress DNA Library Prep Kit	E3325S/L	24/96 rxns
NEBNext UltraExpress FS DNA Library Prep Kit	E3340S/L	24/96 rxns
NEBNext Magnetic Separation Rack	S1515S	24 tubes

RNA	NEB #	SIZE
NEBNext UltraExpress RNA Library Prep Kit	E3330S/L	24/96 rxns
NEBNext Poly(A) mRNA Magnetic Isolation Module	E7490S/L	24/96 rxns
NEBNext rRNA Depletion Kit v2 (Human/Mouse/Rat)	E7400S/L/X	6/24/96 rxns
NEBNext High Input Poly(A) mRNA Magnetic Isolation Module	E3370S	24 rxns
NEBNext Magnetic Separation Rack	S1515S	24 tubes

Products and content are covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc (NEB). The use of trademark symbols does not necessarily indicate that the name is trademarked in the country where it is being read; it indicates where the content was originally developed. See www.neb.com/trademarks. The use of these products may require you to obtain additional third-party intellectual property rights for certain applications. For more information, please email busdev@neb.com.

Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at www.neb.com/support/terms-of-sale. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

© Copyright 2025, New England Biolabs, Inc.; all rights reserved



NEBNext Multiplex Oligos Selection Chart

The chart simplifies the process of selecting the appropriate indexing primers and adaptors for your sequencing applications. It provides a comprehensive overview of available options, including single, dual, and unique dual indices and addresses specific needs like reducing index hopping and incorporating unique molecular identifiers (UMIs) for enhanced accuracy.

Visit www.neb.com/oligos

Distributors

For a list of NEB subsidiaries and distributors, please go to:

www.neb.com/international-support

