

NGS Target Enrichment

NEBNEXT[®] DIRECT[™] CUSTOM READY PANELS

*Flexible
Precision.*




 NEW ENGLAND
BioLabs[®]

be **INSPIRED**
drive **DISCOVERY**
stay **GENUINE**

NEBNext Direct Custom Ready Panels for NGS target enrichment

NEBNext Direct combines a unique hybridization-based technology for highly specific target enrichment of genomic regions of interest with a superior library preparation workflow for Illumina® sequencing. This innovative approach to target enrichment balances the speed and precision of multiplexed PCR-based approaches with the content scalability typical of hybridization-based methods. The flexibility of NEBNext Direct allows a single workflow for assays ranging from single gene tests to comprehensive panels including several hundred genes. Regardless of sample type or assay content, NEBNext Direct allows you to enrich your targets with precision.

NEBNext Direct target enrichment workflow



Advantages of Technology

- Save time with a 1-day workflow that combines enrichment with library preparation
- Generate a higher percentage of your sequencing reads aligning to your targets
- Eliminate the need to over-sequence, reducing cost per sample
- Obtain uniform sequencing of all targets, regardless of gene content
- Generate high quality libraries with limited input amounts and degraded DNA samples, including FFPE and ctDNA
- Distinguish molecular duplicates, reducing false positive variants and improving sensitivity

Customer Feedback:

“The kit and its technology are easy to use and easy to automate, allowing us to get up and running quickly. The protocol itself is fast and efficient to obtain deep coverage of targets, giving homogeneous results for FFPE and frozen tumors, therefore opening doors for customized panels.”

Francis Rousseau, Ph.D., Director of Genomics for IntegraGen SA


“NEBNext Direct enrichment technology is by far the fastest and most automation friendly protocol available today. I can have samples on the sequencer in 6 hours starting from genomic DNA.(...)”

Eric C. Olivares, Founder, SEQanswers.com

Sequence only what you need!


Employing the unique NEBNext Direct hybridization-based enrichment method, NEBNext Direct Custom Ready Panels allow rapid customization of targeted gene panels. Select from a list of genes for which baits have been carefully designed and optimized to produce complete coverage of the full coding regions. High quality panels can be designed by you and rapidly delivered from any combination of genes. NEBNext Direct Custom Ready Panels provide the content you want with the performance you need.

**Unique coverage profile of NEBNext Direct:
Even coverage across target with less padding than other hybridization based approaches.**



IGV image showing coverage profile for exon 3 from the gene *KIT*, targeted using NEBNext Direct. Coverage profile demonstrates the ability to optimize individual baits to produce even coverage across longer exons. 100 ng of DNA was used as input for NEBNext Direct enrichment. Libraries were sequenced using Illumina 2 x 150 basepair sequencing

NEBNext Direct Custom Ready Panels demonstrate optimum performance across a wide range of panel sizes



Key target enrichment metrics demonstrate consistent performance across a range of panel sizes. 100 ng of DNA was tested against panels of 1, 10, 25, 50 and 100 genes, and sequenced using Illumina paired-end 150 bp sequencing. Larger panels included all genes present in smaller panels.

Advantages of Custom Ready Panels

- Built your own panel by choosing from single to hundreds of readily-available genes
- Panel content scalable from a single gene up to 1.5 MB
- No minimum order quantities: Pack sizes with 8, 24 and 96 reactions available
- Extremely fast delivery time
- Experience unmatched specificity and coverage uniformity
- Eliminate synthesis and optimization steps for faster turnaround
- Improve sensitivity with Unique Molecule Indexes (UMI)

Customer Feedback:

“NEB was fantastic while developing our panel or updating an existing one. The protocol is simple and fast and the results have been superb.”

Luca Magnani, Ph.D, CRUK Fellow, Imperial Centre for Translational and Experimental Medicine

“NEBNext Direct Custom Ready Panels have allowed my research to focus on the specific genes we need to explore. In addition to the convenience of easily selecting genes for focused panels, NEBNext Direct enrichment has provided the necessary reliability and depth of coverage to enable robust somatic variant calling.”

Guang Peng, M.D., PhD, MD Anderson Cancer Center

See more testimonials on www.NEBNextDirect.com

See more performance data on the back cover and at
www.NEBNextDirect.com

Currently available NEBNext Direct custom ready genes:

Select from a list of genes and design your panel of interest. As all baits for the listed genes have already been carefully designed and optimized to give complete coverage of the full coding regions, your custom panel is rapidly delivered.

NEBNext Direct Custom Ready Panels are scalable from a single gene up to 1.5 MB of content and can be ordered with no minimum ordering quantity in pack sizes with 8, 24 or 96 reactions!

| | | | | | | | | | |
|----------|----------|----------|---------|---------|--------|--------|---------|-----------|---------|
| AARS | ASXL1 | BSCL2 | CDK12 | CUL1 | EGFR | FANCA | FOXP1 | HDAC4 | IRF1 |
| ABCC9 | ATL1 | BTG1 | CDK4 | CUL3 | EGLN1 | FANCB | FRS2 | HDAC7 | IRF2 |
| ABL1 | ATM | BTG2 | CDK6 | CUX1 | EGR2 | FANCC | FUBP1 | HGF | IRF4 |
| ABL2 | ATP2A1 | BTK | CDK8 | CXCR4 | ELAC1 | FANCD2 | FZR1 | HINT1 | IRF8 |
| ABRAXAS1 | ATP7A | BTLA | CDKN1A | CYLD | ELAC2 | FANCE | GAA | HIST1H1C | IRS2 |
| ACD | ATP7B | BTNL2 | CDKN1B | DAG1 | ELP1 | FANCF | GABRA6 | HIST1H1D | ISPD |
| ACTA1 | ATR | BUB1B | CDKN1C | DAXX | ELP2 | FANCG | GADD45B | HIST1H1E | ITGA7 |
| ACTA2 | ATRX | CACNA1C | CDKN2A | DCLRE1C | EMD | FANCI | GALNT12 | HIST1H2AC | JAK1 |
| ACTB | AURKA | CACNA1S | CDKN2B | DDB2 | EMSY | FANCL | GAN | HIST1H2AG | JAK2 |
| ACTC1 | AURKB | CACNA2D1 | CDKN2C | DDR2 | ENG | FANCM | GARS | HIST1H2AL | JAK3 |
| ACTN2 | AXIN1 | CACNB2 | CEBPA | DDX3X | EP300 | FAS | GATA1 | HIST1H2AM | JARID2 |
| ACVR1B | AXIN2 | CAD | CEP112 | DDX41 | EP400 | FBN1 | GATA2 | HIST1H2BC | JPH2 |
| ADA | AXL | CALR | CEP57 | DES | EPCAM | FBXO11 | GATA3 | HIST1H2BJ | JUN |
| ADGRA2 | B2M | CALR3 | CFL2 | DHX29 | EPHA3 | FBXO31 | GATA4 | HIST1H2BK | JUP |
| ADGRB3 | B3GALNT2 | CAPN3 | CFTR ** | DICER1 | EPHA5 | FBXO7 | GATA6 | HIST1H2BO | KAT6A |
| AIP | B3GNT2 | CARD11 | CHD1 | DIS3L2 | EPHA7 | FBXW7 | GATAD1 | HIST1H3B | KBTBD13 |
| AK2 | BAG3 | CASQ2 | CHD2 | DKC1 | EPHB1 | FGD4 | GDAP1 | HNF1A | KCNE1 |
| AKAP9 | BAP1 | CASR | CHD4 | DMD | ERBB2 | FGF10 | GDNF | HNF1B | KCNE2 |
| AKT1 | BARF1 | CAV3 | CHEK1 | DNAJB2 | ERBB3 | FGF14 | GID4 | HOXB13 | KCNE3 |
| AKT2 | BCL10 | CBFB | CHEK2 | DNAJB6 | ERBB4 | FGF19 | GJB1 | HRAS | KCNH2 |
| AKT3 | BCL11B | CBL | CHKB | DNM2 | ERCC1 | FGF23 | GLA | HSD3B1 | KCNJ2 |
| ALK | BCL2 | CCND1 | CIC | DNMT1 | ERCC2 | FGF3 | GLI1 | HSP90AA1 | KCNJ5 |
| ALMS1 | BCL2L1 | CCND2 | CIITA | DNMT3A | ERCC3 | FGF4 | GMPPB | HSPB1 | KCNJ8 |
| AMER1 | BCL2L2 | CCND3 | CKS1B * | DNMT3B | ERCC4 | FGF6 | GNA11 | HSPB8 | KCNQ1 |
| ANK2 | BCL6 | CCNE1 | CLCN1 | DOCK8 | ERCC5 | FGFR1 | GNA12 | ICK | KDM2B |
| ANKRD1 | BCL7A | CCT6B | CNTN1 | DOT1L | ERG | FGFR2 | GNA13 | ID3 | KDM4B |
| ANO5 | BCOR | CD22 | COL3A1 | DPM1 | ERRF1 | FGFR3 | GNAQ | IDH1 | KDM4C |
| APC | BCORL1 | CD247 | COL6A1 | DPM2 | ESR1 | FGFR4 | GNAS | IDH2 | KDM5A |
| APH1A | BICD2 | CD274 | COL6A3 | DPM3 | ETNK1 | FH | GNB1 | IGF1R | KDM5C |
| APOA4 | BIN1 | CD36 | CORO1A | DPYD | ETS1 | FHIT | GNE | IGF2 | KDM6A |
| APOA5 | BIRC3 | CD3D | CPS1 | DSC2 | ETV1 | FHL1 | GPC3 | IGHMBP2 | |
| APOB | BIRC6 | CD3E | CREBBP | DSG2 | ETV4 | FIG4 | GPD1L | IKBKE | |
| APOC2 | BLM | CD58 | CRKL | DSP | ETV5 | FKRP | GREM1 | IKZF1 | |
| AR | BMPR1A | CD70 | CRLF2 | DTNA | ETV6 | FKTN | GRIN2A | IKZF2 | |
| ARAF | BRAF | CD79A | CRYAB | DTX1 | EWSR1 | FLCN | GRM3 | IKZF3 | |
| ARFRP1 | BRCA1 | CD79B | CSF1R | DUSP2 | EXO1 | FLNC | GRM8 | IL2RG | |
| ARHGAP26 | BRCA2 | CD82 | CSF3R | DUSP9 | EXOSC6 | FLT1 | GSK3B | IL7R | |
| ARID1A | BRD2 | CD83 | CSRFP3 | DVL3 | EXT1 | FLT3 | GTSE1 | ILK | |
| ARID1B | BRD3 | CDC73 | CTCF | DYNC1H1 | EXT2 | FLT4 | H3F3A * | INF2 | |
| ARID2 | BRD4 | CDH1 | CTNNA1 | EBF1 | EZH2 | FOXL2 | H3F3B | INHBA | |
| ASIP | BRIP1 | CDH2 | CTNNA1 | ECT2L | FAF1 | FOXP1 | HCN4 | INPP4B | |
| ASMTL | BRSK1 | CDH4 | CTRC | EED | FAM46C | FOXO1 | HDAC1 | INPP5D | |

* This gene has strong sequence homology to one or more other locations in the genome which can interfere with read mapping. Special care should be taken when interpreting sequencing results for this gene.

** CFTR targets include non-coding, CF-causing variants from the website CFTR2.org.

| | | | | | | | | |
|---------|--------|----------|---------|----------|---------|---------|-----------|--------|
| KDR | MDH2 | NEBL | PEX2 | PTPN1 | RPL35A | SMAD3 | TCF12 | UBA1 |
| KEAP1 | MDM2 | NEFL | PHF21A | PTPN11 | RPL5 | SMAD4 | TCF3 | UGT1A1 |
| KEL | MDM4 | NEXN | PHF6 | PTPN2 | RPS10 * | SMARCA1 | TCL1A | UTP6 |
| KIF1A | MED12 | NF1 | PHIP | PTPN6 | RPS19 | SMARCA4 | TERC | VCL |
| KIF1B | MEF2B | NF2 | PHKA1 | PTPRC | RPS24 | SMARCB1 | TERF2IP | VCP |
| KIF5A | MEF2C | NFE2L2 | PHOX2B | PTPRD | RPS26 * | SMARCE1 | TERT | VEGFA |
| KIT | MEGF10 | NFIB | PICK1 | PTPRO | RPS7 * | SMC1A | TET2 | VHL |
| KLHL40 | MEN1 | NFKBIA | PIGA | PYGM | RPTOR | SMC3 | TFG | VRK1 |
| KLHL6 | MET | NGF | PIK3C2B | QKI | RSPO2 | SMO | TGFBR1 | WAS |
| KMT2A | MFN2 | NHEJ1 | PIK3CA | RAB35 | RUNX1 | SMOX | TGFBR2 | WDR90 |
| KMT2B | MGA | NHP2 | PIK3CB | RAB7A | RUNX1T1 | SNCAIP | TINF2 | WISP3 |
| KMT2C | MIB1 | NKX2-1 | PIK3CG | RAC1 | RYR1 | SNTA1 | TLL2 | WNK1 |
| KMT2D | MITF | NOD1 | PIK3R1 | RAC2 | RYR2 | SNW1 | TMEM127 | WRN |
| KRAS | MKI67 | NOP10 | PIK3R2 | RAD21 | S1PR2 | SOCS1 | TMEM30A | WT1 |
| LAMA2 | MLH1 | NOTCH1 | PIM1 | RAD50 | SBDS | SOCS2 | TMEM43 | XBP1 |
| LAMA4 | MLH3 | NOTCH2 | PKHD1 | RAD51 | SBF2 | SOCS3 | TMEM5 | XPA |
| LAMP2 | MORC3 | NOTCH3 | PKP2 | RAD51C | SCN1B | SOS1 | TMPO | XPC |
| LARGE1 | MPL | NPM1 * | PLCG2 | RAD51D | SCN3B | SOX10 | TMPRSS2 | XPO1 |
| LDB1 | MPZ | NRAS | PLEC | RAF1 | SCN4A | SOX2 | TNFAIP3 | XRCC2 |
| LDB3 | MRE11 | NSD1 | PLEKHG5 | RAG1 | SCN4B | SOX9 | TNFRSF11A | XRCC3 |
| LDLR | MRPL36 | NSD2 | PLN | RAG2 | SCN5A | SPEN | TNFRSF14 | YARS |
| LEF1 | MSH2 | NT5C2 | PMS1 | RANGRF | SCN9A | SPINK1 | TNFRSF17 | YLPM1 |
| LIG4 | MSH3 | NTHL1 | PMS2 | RARA | SDHA | SPOP | TNNC1 | YTHDC1 |
| LITAF | MSH6 | NTRK1 | PNP | RARB | SDHAF2 | SPRED1 | TNNI3 | YY1AP1 |
| LMNA | MSR1 | NTRK2 | POLD1 | RASAL1 | SDHB | SPTA1 | TNNT1 | ZAP70 |
| LMO1 | MTM1 | NTRK3 | POLE | RASGEF1A | SDHC | SPTLC2 | TNNT2 | ZBTB2 |
| LRP1B | MTMR2 | NUP93 | POLH | RASSF1 | SDHD | SRC | TNPO3 | ZBTB33 |
| LRRFIP2 | MTOR | NUP98 | POMGNT1 | RB1 | SELENON | SRCAP | TOP1 | ZFHX3 |
| LRRK2 | MUTYH | ORAI1 | POMT1 | RBM10 | SERP2 | SRSF2 | TOP2A | ZMYM3 |
| LRSAM1 | MXI1 | OTC | POMT2 | RBM20 | SET | STAG2 | TP53 | ZNF217 |
| LTN1 | MYBPC3 | P2RY8 | POT1 | RECQL4 | SETBP1 | STAT3 | TP63 | ZNF24 |
| LYN | MYC | PAG1 | PPM1D | REEP1 | SETD2 | STAT4 | TPM1 | ZNF318 |
| LZTR1 | MYCL | PAK3 | PPP2R1A | RELN | SF1 | STAT5A | TPM2 | ZNF703 |
| MAF | MYCN | PALB2 | PRDM1 | RET | SF3A1 | STAT5B | TPM3 | ZNRF3 |
| MAFB | MYD88 | PALLD | PREX2 | RETREG1 | SF3B1 | STAT6 | TRAF2 | ZRSR2 |
| MAGED1 | MYH10 | PASK | PRF1 | RFFL | SGCA | STIM1 | TRAF3 | |
| MAGI1 | MYH11 | PAX5 | PRKAG2 | RFX7 | SGCB | STK11 | TRAF5 | |
| MAGI2 | MYH7 | PBRM1 | PRKAR1A | RHBDF2 | SGCD | STXBP5 | TRIM32 | |
| MALT1 | MYL2 | PBX1 | PRKCI | RHEB | SGCG | SUFU | TRPV4 | |
| MAP2K1 | MYL3 | PC | PRKDC | RHOA | SGK1 | SUZ12 | TRPV5 | |
| MAP2K2 | MYLK | PCBP1 | PRKN | RICTOR | SH2B3 | SYK | TSC1 | |
| MAP2K4 | MYLK2 | PCLO | PRPF40B | RINT1 | SH3TC2 | SYNE1 | TSC2 | |
| MAP3K1 | MYO18A | PCSK9 | PRPF8 | RIT1 | SHOC2 | TAF1 | TSHR | |
| MAP3K14 | MYOT | PDCD1 | PRPS1 | RNASSEL | SIL1 | TAF4 | TTR | |
| MAP3K6 | MYOZ2 | PDCD11 | PRSS1 | RNF2 | SLC12A6 | TANC2 | TUSC3 | |
| MAP3K7 | MYPN | PDCD1LG2 | PRSS8 | RNF38 | SLC52A2 | TAZ | TYK2 | |
| MAPK1 | NBN | PDGFRA | PRX | RNF43 | SLIT2 | TBL1XR1 | TYR | |
| MAX | NCOR2 | PDGFRB | PTCH1 | ROS1 | SLTM | TBX1 | TYRP1 | |
| MC1R | NCSTN | PDK1 | PTCH2 | RPL11 | SLX4 | TBX3 | U2AF1 * | |
| MCL1 | NDRG1 | PDLIM3 | PTEN | RPL26 | SMAD2 | TCAP | U2AF2 | |

New to NGS Target Enrichment?

Watch our tutorial videos on NEBNextDirect.com incl.:

- Webinar "Challenges and Opportunities for NGS target enrichment"
- NEBNext Direct Workflow overview
- NEB TV – Episode 11 about target enrichment in clinical applications

OVERVIEW

What are NEBNext Direct Custom Ready Panels?

NEBNext Direct Custom Ready Panels allow users to select from an extensive list of genes to create customized target enrichment panels for Illumina® sequencing. Modular bait sets for each gene have been designed, synthesized and optimized for exceptional specificity and target coverage uniformity, and bait sets of these genes can be mixed and matched into customized panels with rapid turnaround time.

What is the minimum reaction commitment for NEBNext Direct Custom Ready Panels?

NEBNext Direct Custom Ready Panels do not have a minimum reaction commitment. Kits are available in 8, 24 and 96 reaction sizes.

What sequencers are NEBNext Direct Custom Ready Panels compatible with?

NEBNext Direct Custom Ready Panels are compatible with the full range of Illumina® sequencing instrumentation.

GENE PANELS

What genes can be included in NEBNext Direct Custom Ready Panels?

The genes available through the NEBNext Direct Custom Ready offering will be continually updated, and currently include those associated with a variety of translational research areas, including cancer, neurological disorders, cardiological disease, autism, severe combined immunodeficiency, cystic fibrosis and the recommended genes for incidental findings by the American College of Medical Genetics. The full list of genes currently available can be found at www.neb.com/CustomReadyPanelForm.

How many genes can I include in a NEBNext Direct Custom Ready Panel, and are there any limitations as to how genes can be combined?

NEBNext Direct Custom Ready Panels can include anywhere from a single specific gene up to 1.5 megabases of total target territory. There are no limitations on genes that can be combined together in a Custom Ready Panel.

What if my gene or region of interest is not on the list?

For the inclusion of genes or other regions of interest that are not available through the NEBNext Direct Custom Ready menu, please contact NEBsolutions@neb.com. Please note that the addition of content beyond what is currently available through NEBNext Direct Custom Ready Panels will impact the pricing and turnaround time of the panels, and minimum reaction quantities may apply.

INPUT

What is the input DNA requirement?

NEBNext Direct Custom Ready Panels are compatible with DNA inputs ranging from 10 ng to 1 µg of genomic DNA. For panels where sequencing data will be used for somatic variant calling, we recommend using 100 ng of DNA or greater.

COVERAGE

What genomic regions are covered for each gene?

Each gene contains baits covering the full coding regions (all exons) for each gene selected, as defined by RefSeq. Baits are designed with variable padding of 0 to 60 bases into the intronic regions.

Can I select specific exons from each gene?

It is possible to select specific exons from each gene, however these requests will be treated as a custom panel, and will impact the turnaround time and pricing of the panel. Please contact nebnextdirect@neb.com for any inquiries.

What is the specificity and coverage uniformity of a NEBNext Direct Custom Ready Panel?

Bait sets for each gene included in the panel have undergone a rigorous development and optimization process to maximize specificity and target coverage uniformity. Because each individual panel contains a unique subset of genes, official specifications for specificity, coverage uniformity, and other performance metrics can not be provided. However, panels typically demonstrate specificity with >85% of the reads mapping to targets and coverage uniformity with >95% of the targeted bases having coverage >25% of the mean target coverage of the panel. Please see the data below for typical specificity and uniformity metrics that can be achieved using NEBNext Direct Custom Ready Panels.

Configure & order your custom panel easily online:

www.neb.com/CustomReadyPanelForm




How much coverage can I expect from my panel?


Target coverage is dependent on the input DNA amount, the total target territory of the panel, and the depth of sequencing. Data shown below demonstrates the typical performance of NEBNext Direct Custom Ready Panels. However, coverage for each individual panel may vary based on the specific genes that are requested.

Mean target depth achieved after PCR duplicate filtering using 10-500 ng of DNA and panels containing 1-100 genes

1.5-hour hybridization



16-hour hybridization



ANALYSIS

What are the recommended sequencing read lengths?

We recommend paired end 150 base sequencing to adequately cover the targets. In addition, 8 bases of i7 index sequencing, to read the sample index and 12 bases of i5 index sequencing, to read the unique molecule ID.

What are the recommendations for analyzing data from NEBNext Direct Custom Ready Panels?

We have developed and optimized a pipeline for internal processing of data directly from FASTQ files using open-source bioinformatics tools. Details on this pipeline can be found at: <https://github.com/DirectedGenomics/DemoPipeline>. Additionally, NEB has an agreement with Bluebee, who have hosted the demo pipeline on their platform, where NEB can provide activation codes for a limited number of samples to run through the pipeline. In order to request activation codes, please email nebnextdirect@neb.com. Please note that these pipelines are intended for evaluation use only, and that any production pipeline should be implemented and verified appropriately, with parameters optimized for the intended use.

PRODUCT DETAILS

Do NEBNext Direct Custom Ready Panels undergo verification prior to shipment?

Yes. For each NEBNext Direct Custom Ready Panel that is purchased, an NEBNext Direct library is created and sequenced internally to ensure that the expected target regions are included and that specificity and target coverage uniformity are optimal. A performance report will be emailed to you along with shipment of the panel.

How much does a NEBNext Direct Custom Ready Panel cost?

Tiered pricing is offered based on the total target territory included in the panel. Volume-based pricing is also available. Please contact your local New England Biolabs representative for official pricing, or visit www.neb.com/CustomReadyPanelForm.

How do I order an NEBNext Direct Custom Ready Panel?

Panels can be designed and ordered by visiting www.neb.com/CustomReadyPanelForm, launching the web tool, and selecting the genes of interest. Once a design has been submitted, you will receive target bed files for the coordinates included in the panel and a quotation for the panel that includes pricing and a unique part number for order placement.

How long will it take to receive my NEBNext Direct Custom Ready Panel?


Because the development and optimization of the baits has already been performed, panels are typically shipped within two weeks of order receipt.

NEBNext Direct is also available as
Genotyping-by-Sequencing Solution!

Please contact us to learn more:


info.de@neb.com

Sensitivity in detection of variants across panel size and DNA input amount



24 HapMap samples were blended to create a range of variant allele frequencies (VAF) down to 2%. 25, 50, 100, 200 and 500 ng of this blended DNA was enriched using NEBNext Direct Custom Ready Panels of 1, 10, 50, and 100 genes. Larger panels were inclusive of the genes in smaller panels. Resulting libraries were sequenced using 2 x 150 bp Illumina sequencing and variants were called using Mutect and Vardict variant calling algorithms.

NEBNext Direct Custom Ready Panels demonstrate retention of target behavior across panel sizes



IGV image of coverage profile for 4 BRAF exons included in panels of 1, 10, 25, 50 and 100 genes, demonstrate consistent target behavior with the addition of gene targets. 100 ng of DNA was used as input for NEBNext Direct enrichment using the 5 panels, including the BRAF gene. Libraries were sequenced using Illumina 2 x 150 basepair sequencing.

ORDERING INFORMATION

| PRODUCTS | NEB # | SIZE |
|-----------------------------------|------------|--------------|
| NEBNext Direct Custom Ready Panel | E6631S/L/X | 8/24/96 rxns |

| ALSO AVAILABLE | NEB # | SIZE |
|-------------------------------------|------------|--------------|
| NEBNext Direct Cancer HotSpot Panel | E7000S/L/X | 8/24/96 rxns |
| NEBNext Direct BRCA1/BRCA2 Panel | E6627S/L/X | 8/24/96 rxns |

For research use only.

One or more of these products are covered by patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at gbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications.

Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at www.neb.com/support/terms-of-sale.

NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

ILLUMINA® is a registered trademark of Illumina, Inc.

© Copyright 2019, New England Biolabs, Inc.; all rights reserved.

x/NEB131 – 05/19

GERMANY & AUSTRIA
 New England Biolabs GmbH
 Brünigstr. 50, Geb B852
 65926 Frankfurt/Main, Germany
 Tel: +49/(0)69/305-23140
 Fax: +49/(0)69/305-23149

Free Call: 0800/246 5227 (Germany)
 Free Call: 00800/246 52277 (Austria)
info.de@neb.com

www.neb-online.de

FRANCE
 New England Biolabs France SAS
 Genopole Campus 1, Bâtiment 6
 5 rue Henri Desbruères
 91030 Evry cedex, France
 Tel.: 0800 100 632 (Customer Service)
 Tel.: 0800 100 633 (Technical Service)
 FAX.: 0800 100 610
info.fr@neb.com

www.neb-online.fr

HEADQUARTERS:
USA

New England Biolabs, Inc.
 Telephone: (978) 927-5054

Toll Free (USA Orders): 1-800-632-5227
 Toll Free (USA Tech): 1-800-632-7799
 Fax: (978) 921-1350
info@neb.com
www.neb.com

www.neb.com

AUSTRALIA
 New England Biolabs, Inc.
 Toll Free: +61 1800 934 218
info.au@neb.com

CANADA
 New England Biolabs, Ltd.
 Toll Free: 1-800-387-1095
info.ca@neb.com

CHINA, PEOPLE'S REPUBLIC
 New England Biolabs (Beijing), Ltd.
 Telephone: 010-82378265/82378266
info@neb-china.com

JAPAN
 New England Biolabs Japan, Inc.
 Telephone: +81 (0)3 5669 6191
info.jp@neb.com

SINGAPORE
 New England Biolabs, PTE. Ltd.
 Telephone: +65 63859623
sales.sg@neb.com

UNITED KINGDOM
 New England Biolabs (UK), Ltd.
 Call Free: 0800 318486
info.uk@neb.com

