

Isothermal DNA Amplification

ROBUST TECHNOLOGIES FOR RAPID NUCLEIC ACID DETECTION

Update
2015/16




What is isothermal DNA amplification?

The Polymerase Chain Reaction (PCR) is a well-known approach for amplifying a specific DNA sequence. PCR involves the reiterative cycling of a reaction cocktail between different temperatures to achieve amplification. As routine as PCR is in the molecular biology and molecular diagnostic laboratory, there are other methods of sequence-specific DNA amplification.

These alternative approaches often do not require changing the reaction temperature and are, therefore, often referred to as isothermal amplification protocols. Isothermal amplification protocols are varied and have different advantages. In general, isothermal techniques are extremely fast and do not require thermocyclers, making them particularly well suited for field applications and point-of-care molecular diagnostics assays.

Advantages

- Fast
- Minimal equipment required
- Robust reactions in the presence of inhibitors
- Simplified optical detection



Interested in learning how NEB scientists are using isothermal amplification in their research?

Visit www.neb.com/IsothermalAmplification to find videos, protocols and recent publications, including a recent publication from NEB scientists, describing a pH-sensitive isothermal detection method.

Featured Products for Isothermal


Bst 3.0 DNA Polymerase

✓ Validated for LAMP

Bst 3.0 DNA Polymerase (NEB #M0374) is an *in silico* designed homolog of *Bacillus stearothermophilus* DNA Polymerase I, Large Fragment, (NEB #M0275) engineered for improved isothermal amplification performance and increased reverse transcription activity. *Bst* 3.0 contains 5' → 3' DNA polymerase activity with either DNA or

RNA templates and strong strand displacement activity, but lacks 5' → 3' and 3' → 5' exonuclease activity. *Bst* 3.0 demonstrates robust performance even in high concentrations of amplification inhibitors and features significantly increased reverse transcriptase activity compared to *Bst* DNA Polymerase.

Fast, single-enzyme RT-LAMP can be performed using *Bst* 3.0



RT-LAMP was performed using indicated DNA polymerase and Jurkat total RNA and primers for two genes (ACTB, left; HMBS2, right). Fastest results were observed with a 2-enzyme system, *Bst* 2.0 and WarmStart RTx, but robust amplification was also observed using *Bst* 3.0 without additional RT. *Bst* LF, *Bst* 2.0 and competitor enzymes showed highly variable performance, with slow threshold times or reaction failure on one of the two targets.

Not sure which product will work best for your experiment?

NEB offers a selection of *Bst* DNA Polymerase-based products for isothermal DNA amplification. Use this chart to determine which product will work best for your needs.

	5' → 3' EXO ACTIVITY	AMPLIFICATION SPEED	ROOM TEMPERATURE SETUP	REVERSE TRANSCRIPTASE ACTIVITY	INHIBITOR TOLERANCE	APPLICATIONS
<i>Bst</i> DNA Polymerase, Full Length	★★	N/A	N/A	N/A	★	Nick translation reactions at elevated temperatures
<i>Bst</i> DNA Polymerase, Large Fragment	N/A	★	N/A	★	★	General strand-displacement reactions, original polymerase for LAMP and other diagnostic amplifications
<i>Bst</i> 2.0 DNA Polymerase	N/A	★★	N/A	★★	★★	Improved LAMP, SDA, and other amplification reactions
<i>Bst</i> 2.0 WarmStart DNA Polymerase	N/A	★★	★★★	★★	★★	Consistent, room-temperature, and high-throughput amplification assays
<i>Bst</i> 3.0 DNA Polymerase	N/A	★★★	★★	★★★	★★★	Fastest, most robust LAMP and RT-LAMP reactions. High reverse transcriptase activity up to 72°C

- ★★★ Optimal, recommended product for selected application
- ★★ Works well for selected application
- ★ Will perform selected application, but is not recommended
- N/A Not applicable to this application



Did you know that many of these products can be purchased in larger volumes? Contact bulks.de@neb.com to find out more.


DNA Amplification from NEB

WarmStart RTx Reverse Transcriptase

✓ Validated for RT-LAMP


WarmStart RTx Reverse Transcriptase (NEB #M0380) is a unique *in silico* designed RNA-directed DNA polymerase coupled with a reversibly-bound aptamer that inhibits RTx activity below 40°C. This enzyme can synthesize a complementary DNA strand initiating from a primer using RNA (cDNA synthesis) or single-stranded DNA as a template. RTx is a robust enzyme for RNA detection in amplification reactions and is particularly well-suited for use in loop-mediated isothermal amplification (LAMP). The WarmStart property enables high throughput applications, room temperature setup, and increases the consistency and specificity of amplification reactions. RTx contains intact RNase H activity.

WarmStart control of WarmStart RTx



cDNA synthesis was performed for 10 minutes, followed by qPCR analysis. Resulting Ct's were normalized to a "no RT" control for 0% activity and fastest Ct for 100% activity. WarmStart RTx is inhibited by a reversibly bound aptamer at temperatures below 40°C, and is fully active at temperatures 42°C and higher.

WarmStart improves speed and sensitivity in RT-LAMP



RT-LAMP reactions with *Bst* 2.0 WarmStart DNA Polymerase and the indicated reverse transcriptase were incubated at 65°C with 1 pg – 100 ng of Jurkat total RNA. Reactions were monitored with real-time fluorescence, and the resulting curves are shown on left, with the corresponding threshold times on right. WarmStart RTx provides faster reaction threshold times for improved consistency and sensitivity with lower input RNA amounts. RT-LAMP reactions performed with AMV Reverse Transcriptase resulted in inconsistent detection, as indicated by wide variation at lower RNA input concentrations (blue curves).

What is WarmStart?

"WarmStart" is the term we use to describe a mesophilic enzyme that is inactive at room temperature, and activated when the reaction is warmed above approximately 40°C. NEB currently supplies two enzymes, *Bst* 2.0 DNA Polymerase (NEB #M0537) and WarmStart RTx Reverse Transcriptase (NEB #M0380), that have this property.

For room temperature reaction setup, such as in high-throughput workflows, we recommend pairing *Bst* 2.0 WarmStart with WarmStart RTx for RT-LAMP. The double WarmStart mixture ensures reproducibility.

Optimization tips for LAMP

- Use LAMP primer design software (e.g., Primer Explorer – primerexplorer.jp/e/). Select 2–3 sets for each target and compare performance in a LAMP assay.
- Include loop primers for faster reactions
- Use high magnesium (6–8 mM) and dNTP (1–1.4 mM) concentrations for best reactions
- Omit betaine, unless it has a demonstrated benefit
- Optimize the reaction temperature (60–65°C for *Bst* LF and 63–70°C for *Bst* 2.0/3.0)
- To prevent contamination, use *Bst* 3.0 or Antarctic Thermo-labile UDG (NEB #M0372), which denatures rapidly


Examples of isothermal technologies

Loop-mediated Isothermal Amplification (LAMP & RT-LAMP)


LAMP is designed to detect a target nucleic acid without sophisticated equipment. LAMP uses 4-6 primers recognizing 6-8 distinct regions of the target DNA. A strand-displacing DNA polymerase initiates synthesis and two of the primers form loop structures to facilitate subsequent rounds of amplification. LAMP provides high sensitivity (fg or <10 copies of target), and reactions can be performed in as little as 5–10 minutes. Additionally, reactions can be performed with limited resources (e.g., using a water bath for incubation and detection of results by eye), or with real-time measurement and high-throughput instruments.

Detection of RNA targets is accomplished by simple addition of a reverse transcriptase to the LAMP reaction (e.g., WarmStart® RTx Reverse Transcriptase), with RT-LAMP performed as a true one-step, isothermal workflow.

Overview of LAMP



Overview of SDA



* Target amplification, shown above for SDA_F, will also occur simultaneously with SDA_R.

Strand Displacement Amplification (SDA)

SDA relies on a strand-displacing DNA polymerase, typically *Bst* DNA Polymerase, Large Fragment (NEB #M0275) or Klenow Fragment (3'→5' exo⁻) (NEB # M0212), to initiate at nicks created by a strand-limited restriction endonuclease or nicking enzyme (e.g., Nt.BstNBI, NEB # R0607) at a site contained in a primer. The nicking site is regenerated with each polymerase displacement step, resulting in exponential amplification. SDA is typically used in clinical diagnostics.

Helicase-dependent Amplification (HDA)

HDA employs the double-stranded DNA unwinding activity of a helicase to separate strands, enabling primer annealing and extension by a strand-displacing DNA polymerase. Like PCR, this system requires only two primers. HDA has been employed in several diagnostic devices and FDA-approved tests.

Nicking Enzyme Amplification Reaction (NEAR)

NEAR employs a strand-displacing DNA polymerase initiating at a nick created by a nicking enzyme, rapidly producing many short nucleic acids from the target sequence. This process is extremely rapid and sensitive, enabling detection of small target amounts in minutes. NEAR is commonly used for pathogen detection in clinical and biosafety applications.

Choose from our selection of products

For your isothermal DNA application.

PRODUCT	NEB #	SIZE
<i>Bst</i> 3.0 DNA Polymerase	M0374S/L/M	1,600/8,000/8,000 units
<i>Bst</i> 2.0 WarmStart DNA Polymerase	M0538S/M/L	1,600/8,000 units
<i>Bst</i> 2.0 DNA Polymerase	M0537S/M/L	1,600/8,000 units
<i>Bst</i> DNA Polymerase, Large Fragment	M0275S/M/L	1,600/8,000 units
<i>Bst</i> DNA Polymease, Full Length	M0328S/L	500/2,500 units
WarmStart RTx Reverse Transcriptase	M0380S/L	50/250 reactions
Nt.BstNBI	R0607S/L	1,000/5,000 units
IsoAmp® II Universal tHDA Kit	H0110S	50 reactions
AMV Reverse Transcriptase	M0277S/T/L	200/500/1,000 units
Antarctic Thermolabile UDG	M0372S/L	100/500 units
Deoxynucleotide (dNTP) Solution Mix	N0447S/L	8/40 µmol of each
Deoxynucleotide (dNTP) Solution Set	N0446S	25 µmol of each

NEW ENGLAND BIOLABS®, NEB® and WARMSTART® are registered trademarks of New England Biolabs, Inc.

OMNIAMP® is a registered trademark of Lucigen, Inc.

ISOAMP® is a registered trademark of BioHelix Corporation. The IsoAmp® II Universal tHDA Kit was developed and produced by BioHelix Corporation, now a wholly owned subsidiary of Quidel Corporation.

The purchase of NEB RTx products conveys to the purchaser the limited, nontransferable right to use the purchased products to perform reverse transcription loop-mediated isothermal amplification ("RT-LAMP") for research use only. LAMP is a patented technology belonging to Eiken Chemical Co., Ltd., and any use other than research may require a license from Eiken Chemical Co., Ltd. A patent is pending for NEB's RTx product.

Germany & Austria

New England Biolabs GmbH
Brüningstr. 50, Geb B852
65926 Frankfurt/Main, Germany
Tel: +49/(0)69/305-23140
Fax: +49/(0)69/305-23149

Free Call: 0800/246 5227 (Germany)
Free Call: 00800/246 52277 (Austria)
info.de@neb.com
www.neb-online.de

USA

New England Biolabs, Inc.
Telephone (978) 927-5054
Toll Free (USA Orders) 1-800-632-5227
Toll Free (USA Tech) 1-800-632-7799
Fax (978) 921-1350
info@neb.com
www.neb.com

Canada

New England Biolabs, Ltd.
Toll Free: 1-800-387-1095
info.ca@neb.com

China, People's Republic

New England Biolabs (Beijing), Ltd.
Telephone: 010-82378265/82378266
info@neb-china.com

France

New England Biolabs France
Telephone : 0800 100 632
info.fr@neb.com

Japan

New England Biolabs Japan, Inc.
Telephone: +81 (0)3 5669 6191
info@neb-japan.com

Singapore

New England Biolabs, PTE. Ltd.
Telephone +65 6776 0903
sales.sg@neb.com

United Kingdom

New England Biolabs (UK), Ltd.
Call Free: 0800 318486
info.uk@neb.com



www.neb-online.de