

Robust Colony PCR from Multiple *E. coli* Strains using OneTaq[®] Quick-Load[®] Master Mixes

Yan Xu, Ph.D., New England Biolabs, Inc.

Introduction

Colony PCR is a commonly used method to quickly screen for plasmids containing a desired insert directly from bacterial colonies. This method eliminates the need to culture individual colonies and prepare plasmid DNA before analysis. However, the presence of bacterial cell contents and culture media in colony PCR reactions often results in polymerase inhibition. A robust polymerase is required to perform colony PCR with high efficiency in many different bacterial strains.

OneTaq DNA Polymerase, an optimized blend of Taq and Deep Vent[™]_R DNA polymerases, has been formulated for robust yields with minimal optimization. This robustness makes OneTaq ideal for use in demanding applications, such as colony PCR.

Furthermore, the OneTaq Quick-Load Master Mix product format increases the ease-of-use for colony PCR. The master mix formulation contains dNTPs, MgCl₂, buffer components and stabilizers, as well as two commonly used tracking dyes for DNA gels. On a 1% agarose gel in 1X TBE, Xylene Cyanol FF migrates at ~4 kb and Tartrazine migrates at ~10 bp. Both dyes are present in concentrations that do not mask any co-migrating DNA bands.

General Protocol

1. Transform ligation mix or other plasmid-containing reaction mixture into the desired bacterial strain, and incubate agar plates overnight at the appropriate temperature.
2. Set up 50 µl reactions as follows:

| | |
|-------------------|----------|
| OneTaq Master Mix | 25 µl |
| PCR primer | 200 nM |
| H ₂ O | to 50 µl |

Note: If OneTaq Hot Start Quick-Load 2X Master Mix is used, reactions can be set up at room temperature. If OneTaq Quick-Load 2X Master Mix is used, reactions should be set up on ice.

3. Use a sterile toothpick to pick up individual colonies and dip into each reaction tube.
4. As soon as the solution looks cloudy, remove the toothpick. To create a stock of each individual colony either:
 - a.) Dip the toothpick into 3 ml growth media with appropriate antibiotics and culture overnight.
 - or
 - b.) Streak the toothpick onto another agar plate containing the appropriate antibiotics and grow overnight.
5. Transfer reactions to a PCR cyclor, and perform PCR following the guidelines on the right for cycling conditions:

| | |
|------------------------------|-----------------|
| Initial denaturation: | |
| 94°C | 2 minutes |
| 30 cycles: | |
| 94°C | 15–30 seconds |
| 45–68°C | 15–60 seconds |
| 68°C | 1 minute per kb |
| Final hold: | |
| 68°C | 5–10 minutes |
| 10°C | hold |

6. Load 4–6 µl of each PCR reaction directly onto an agarose gel, alongside an appropriate DNA ladder.

Materials

- Well-isolated bacterial colonies, ideally 1–2 mm in diameter
- Sterile toothpicks or pipette tips
- Additional agar plate, or culture tubes with growth media for retention of original colonies.
- OneTaq Quick-Load 2X Master Mix with Standard Buffer (M0486) or OneTaq Hot Start Quick-Load 2X Master Mix with Standard Buffer (M0488)*
- Sterile H₂O
- PCR primers


**For amplicons with a GC content over 65% GC, OneTaq Quick-Load 2X Master Mix with GC Buffer or OneTaq Hot Start Quick-Load 2X Master Mix with GC Buffer may be used.*

Results


Colony PCR was performed in 2 separate experiments using the protocol described above, with the following colonies:

1. Colonies obtained from transformation of a plasmid with a 4.5 kb insert into 18 different *E. coli* strains. Amplification of the plasmid insert was achieved in each case. *OneTaq* Quick-Load 2X Master Mix with Standard Buffer and *OneTaq* Hot Start Quick-Load 2X Master Mix with Standard Buffer were used.
2. Colonies from *E. coli* library clones with inserts ranging from 0.8 kb to 10 kb. *OneTaq* Quick-Load 2X Master Mix with Standard Buffer was used, and results illustrate the robustness of the *OneTaq* Quick-Load 2X Master Mix in this application.

Colony PCR of library clones with inserts of 0.8 kb – 10 kb, using *OneTaq* Quick-Load 2X Master Mix



Colony PCR of a 4.5 kb insert using *OneTaq* and *OneTaq* Hot Start Quick-Load 2X Master Mixes with Standard Buffer and 18 different *E. coli* strains



Reactions were set up according to the protocol and analyzed by agarose electrophoresis. Marker M is the 1 kb DNA Ladder (NEB #N3232)

| Lane | Name | NEB # | Lane | Name | NEB # | Lane | Name | NEB # |
|------|-----------------|-------|------|--------------------------------|-------|------|------------------------------|-------|
| 1 | NEB 10-beta | C3019 | 7 | Lemo21(DE3) | C2528 | 13 | T7 Express lysY | C3010 |
| 2 | NEB 5-alpha | C2987 | 8 | NIc021(DE3) | C2529 | 14 | T7 Express | C2566 |
| 3 | NEB 5-alpha F1q | C2992 | 9 | NEB Express F ⁺ | C3037 | 15 | T7 Express Crystal | C3022 |
| 4 | dam-/dcm- | C2925 | 10 | NEB Express | C2523 | 16 | SHuffle [®] Express | C3028 |
| 5 | NEB Turbo | C2984 | 11 | T7 Express F ⁺ | C3016 | 17 | SHuffle T7 Express lysY | C3030 |
| 6 | BL21(DE3) | C2527 | 12 | T7 Express lysY/F ⁺ | C3013 | 18 | SHuffle T7 Express | C3029 |

Summary

OneTaq and *OneTaq* Hot Start Quick-Load Master Mixes provide reliable performance in colony PCR, and are compatible with multiple *E. coli* strains. Reliable performance has been seen with amplicons up to 10 kb. The Quick-Load format offers additional convenience by enabling direct loading of the PCR reaction onto an agarose gel for analysis. Lastly, the Hot Start formulation provides additional functionality by reducing interference from primer-dimers and secondary amplification products.

Ordering Information:

| PRODUCT | NEB # | SIZE | PRICE | |
|---|------------|-----------------------|------------------------|----------------------------|
| <i>OneTaq</i> [®] DNA Polymerase | M0480S/L/X | 200/1,000/5,000 units | 38 € / 152 € / 600 € | |
| <i>OneTaq</i> [®] 2X Master Mix with Standard Buffer | M0482S/L | 100/500 rxns | 39 € / 156 € | |
| <i>OneTaq</i> [®] 2X Master Mix with GC Buffer | M0483S/L | 100 | 39 € | |
| <i>OneTaq</i> [®] Quick-Load [®] 2X Master Mix with Standard Buffer | M0486S/L/X | 100/500/5x500 rxns | 40 € / 160 € / 708 € | Featured in this App Note. |
| <i>OneTaq</i> [®] Hot Start DNA Polymerase | M0481S/L/X | 200/1,000/5,000 units | 78 € / 312 € / 1,248 € | |
| <i>OneTaq</i> [®] Hot Start 2X Master Mix with Standard Buffer | M0484S/L | 100/500 rxns | 66 € / 272 € | |
| <i>OneTaq</i> [®] Hot Start 2X Master Mix with GC Buffer | M0485S/L | 100/500 rxns | 66 € / 272 € | |
| <i>OneTaq</i> [®] Hot Start Quick-Load [®] 2X Master Mix with Standard Buffer | M0488S/L | 100/500 rxns | 69 € / 282 € | Featured in this App Note. |
| <i>OneTaq</i> [®] Hot Start Quick-Load [®] 2X Master Mix with GC Buffer | M0489S/L | 100/500 rxns | 69 € / 282 € | |
| <i>OneTaq</i> [®] RT-PCR Kit | E5310S | 30 rxns | 142 € | |

NEW ENGLAND BIOLABS[®], ONETAQ[®], QUICK-LOAD[®] and SHUFFLE[®] are registered trademarks of New England Biolabs, Inc. DEEP VENT[™] is a trademark of New England Biolabs, Inc.



www.neb-online.de



be INSPIRED
drive DISCOVERY
stay GENUINE