



Stay genuine.

Since our establishment in 1974, New England Biolabs has been different. From our founding principles — placing the advancement of science and stewardship of the environment as our highest priorities — to our unique corporate culture, NEB's philosophy can be distilled down to three core values: passion, humility and being genuine.

As part of our ongoing commitment to these values, in 2014, we established the Passion in Science Awards® to recognize the unsung heroes in the scientific community who contribute to making the world a better place through their inspirational and innovative work in artistic expression, humanitarian service, environmental stewardship and scientific mentorship. We have always believed that science is more than just a vocation — it embodies an ethos that inspires acts of compassion, brilliance and originality.

Thank you for your ongoing trust and support. If there is anything you believe we should be doing differently, please share your thoughts with us. We wish you continued success in your research.

Want to experience the NEB catalog in augmented reality?



Download

the NEB Augmented Reality (AR) app for iPhone® or iPad® at the Apple® App Store or for Android™ on Google Play™.



Explore

the NEB Catalog & Technical Reference and keep an eye out for the augmented reality icon.



Discover

videos, tutorials and immersive experiences that can be accessed by scanning the icon using the AR app.



be INSPIRED

drive **DISCOVERY**

stay **GENUINE**

At New England Biolabs, we are motivated by a set of core values that are still as true today as they were when the company was founded 45 years ago. These principles continue to guide us both as a company and as individuals.

Advancement of Science

We believe that basic research and the cultivation of scientific knowledge is critical for us to stay connected with our customers and to drive scientific breakthroughs. At NEB, over 30 labs participate in research projects, which are aided by post-doctoral fellows and students in Masters and Ph.D. programs. NEB researchers have authored or co-authored over 1,200 publications to date. most of which are in peer-reviewed journals.

Environmental Stewardship

We continuously strive to promote ecologically sound practices and environmental sustainability in order to protect our natural resources, both locally and globally. Further, it is our goal to continuously improve our business processes to minimize our impact on the environment.

Humanitarianism

We see opportunities where science can be used to improve lives, and are inspired by researchers who help push the social barometer to a kinder, more compassionate, and healthier future. This philosophy lies behind NEB's longstanding commitment to its parasitology research program, which contributes to the understanding and treatment of poorly-funded and understudied tropical diseases. NEB has also helped to establish several foundations devoted to humanitarian efforts.

Delivering the Highest Quality Product

In order to accelerate your research, it is our goal to deliver best-in-class product quality and technical support. With a reliance on recombinant technologies, our products are designed and manufactured in our ISO 13485:2016 and ISO 9001:2015 certified facility in Ipswich, MA, USA. We are constantly improving the stringency and range of our quality controls to ensure that our products will perform to your expectations, every

time. In addition, we recently expanded our manufacturing footprint by opening a facility in Rowley, MA, for the production of GMP-grade materials for customers requiring an enhanced level of quality documentation and support.









A Unique Approach to Wastewater Treatment

Our state-of-the-art Solar Aquatics System™ utilizes and accelerates the process found in streams and wetlands to treat the campus' wastewater, making it clean enough for groundwater recharge.

Solar Aquatics System™ is a trademark of Ecological Engineering Associates.





The NEB Facility

Our research and production facility, located in Ipswich, MA, USA, is LEED® certified, which is awarded based on environmentally-focused standards that include:

- · Water efficiency
- Energy conservation
- · Atmospheric protection
- Sustainable building materials and resources
- · Indoor environmental quality
- Innovation and building design

 $\textit{LEED}{}^{\otimes} \textit{ is a registered trademark of the U.S. Green Building Council Corporation}.$

Recycling at NEB

NEB established the first shipping box recycling program over 40 years ago, in order to divert polystyrene from landfills.

Our extensive in-house recycling program includes:

- Paper
- Aluminum
- Plastics
- Batteries
- Glass
- Electronics

NEB also composts its waste, diverting as much as 75% of the waste from landfills.



Partnering with NEB

With over 40 years of leadership in the life science industry, our scientific expertise, global reach, and proven track record of turning innovative ideas into successful products positions NEB as a compelling partner. With experience in fields as diverse as next generation sequencing, RNA biology, qPCR and protein engineering, NEB is ready to work with you to develop custom solutions specific to your needs, and to help bring your technologies to market. Further, our global distribution network can help to ensure that your products will have worldwide reach.

Custom Solutions

From development to commercialization, NEB provides the technical expertise, consistent scalable manufacturing, quality systems and a global distribution network to enable a successful long-term partnership. Our dedicated team is ready to work with you to develop novel, high performance enzymes tailored to your application, optimize these enzymes in your workflow, enable small to large scale production, generate quality controls and customize packaging. With our ISO 13485 and ISO 9001 certified manufacturing processes, as well as the ability to manufacture GMP-grade products, you can be confident in our robust process, documentation and risk mitigation for the product you need. For more information, contact custom@neb.com.

Global Business Development

The business development team at NEB operates on a global basis to enable innovation in molecular technologies through strategic partnerships, licensing and new ventures. We do so by leveraging the talents and assets of NEB, including our scientific, commercial and international resources. Co-development collaborations benefit from access to our expertise and proprietary technologies, ensuring commercial outcomes with shorter innovation cycles. Further success is then achieved due to our privileged market position as products enter a well-managed global distribution network. For more information, contact busdev@neb.com.

International Network

NEB has extensive worldwide distribution capabilities through a network of exclusive distributors and agents, plus wholly-owned subsidiaries located in Australia, Canada, China, France, Germany, Japan, Singapore and the United Kingdom. For more information, please see the back inside cover.

NEBnow Freezer Program Network

Your research doesn't have to be put on hold waiting for that critical reagent to be delivered. With NEBnow® on-site freezers, enjoy convenient and affordable access to NEB's high quality reagents, anytime. Our NEBnow Freezer Program Team works closely with your institution to customize inventory best suited to your research program. Save time and avoid shipping fees with consolidated shipments. For more information, contact freezers@neb.com.

Enzymes for Innovation

The NEB catalog highlights a wide variety of enzyme functionalities found in nature or engineered for specific purposes. However, in molecular biology, new tools can often lead to new discoveries. Taking advantage of the enzymology expertise at NEB, we now offer a growing selection of novel enzymes with interesting and unique activities for manipulating DNA, RNA, proteins and glycans, even if specific applications for them have yet to be discovered. Our hope is that by engaging researchers' imaginations, our "Enzymes for Innovation" initiative will enable the discovery of new molecular tools and workflows. If you are looking for an enzyme functionality that it is not currently available, visit www.enzymesforinnovation.com or contact enzymesforinnovation@neb.com.



We see ourselves as an extension of your manufacturing and operations team, dedicated to enabling you to develop new technologies, products and services, and delivering them with exceptional service and support.

- Director, OEM & Custom Solutions, New England Biolabs





Trust, transparency and timely communications are key to achieving a mutually beneficial outcome.
Our goal is to exceed your expectations as we work together to carry your innovations forward.

- Executive Director, Global Business Development, New England Biolabs

The NEBnow freezer program is an amazing addition; having 24/7 access to the NGS reagents has been extremely helpful. NEB allows you to customize the freezer to hold the items you use most, saving us valuable time by cutting out the ordering process!

– Assistant Director, Genomic Sequencing and Analysis Facility, University of Texas, Austin



Non-Profits and Research Foundations

New England Biolabs has played a role in the establishment of several organizations that are advancing humanitarian efforts and environmental stewardship.







Creative Action Institute transforms the way people communicate and collaborate to catalyze community-driven solutions that advance gender equality and build a sustainable planet. Through our experiential trainings, convenings and coaching, we develop leadership, build networks and support grassroots advocacy. Visit CreativeActionInstitute.org to learn more.



The Ocean Genome Legacy Center of New England Biolabs is a non-profit research center dedicated to the conservation of marine genome diversity and maintaining a repository of genomic DNA from marine organisms around the world. Learn more at northeastern.edu/ogl.

Bestellinformationen (Deutschland & Österreich):

Kontakt

New England Biolabs GmbH Brüningstr. 50, Geb. B852 65926 Frankfurt/Main Deutschland

Telefon: +49/(0)69/305-23140 Fax: +49/(0)69/305-23149

Email: info.de@neb.com

Online-Shop und technische Information:

www.neb-online.de

Kostenfreie Servicenummern

Deutschland

Telefon: 0800/246-5227 Fax: 0800/246-5229

Österreich

Telefon: 00800/246-52277

Öffnungszeiten

Wir stehen Ihnen montags - freitags von 8:30 bis 17:00 Uhr persönlich zur Verfügung.

Kostenfreier Technischer Support

Nutzen Sie bei technischen Fragen gerne kostenfrei die Expertise unseres wissenschaftlichen Beratungsteams.

Telefon: 0800/246-5227 (in D) bzw. 00800/246-52277 (in A) oder: Email: techsupport.de@neb.com

24 Stunden Lieferservice

Bei Bestelleingang werktags (Mo-Do) bis 16:00 Uhr erhalten Sie Ihre Ware am nächsten Tag! Bestellungen, die vor einem Wochenende/ Feiertag eingehen, werden am kommenden Montag/ Werktag versendet.

Mehrweg-Versandboxen

Die Versandboxen inkl. Kühlakkus sind Mehrwegverpackungen und werden in Deutschland von unseren Transportunternehmen regelmäßig bei Ihnen abgeholt und wiederverwendet.

Verpackung- & Transportpauschale

Deutschland

Brief/Paketlieferungen frei. Kurier-/Kühlsendungen 15 €; frachtkostenfrei ab einem Nettowert von 250 €.

Österreich

Brief/Paketlieferungen frei. Kurier-/Kühlsendungen 20 €; frachtkostenfrei ab einem Nettowert von 350 €.

NEB erhebt keine zusätzlichen Mindermengen- oder Trockeneiszuschläge!

Garantie und Gewährleistung

NEB GmbH garantiert für alle Produkte die Spezifikation entsprechend des jeweiligen Produktdatenblatts. Bei Nichterfüllung liefern wir kostenfreien Ersatz. Die Reklamation von fehlerhaften Produkten muss innerhalb von 60 Tagen nach Erhalt der Ware erfolgen.

Preise

Die in diesem Katalog niedergelegten Angaben bezüglich der Preise, Mengen und Liefermöglichkeiten sind freibleibend. Preise können ohne vorherige Ankündigung geändert werden. Es gelten die jeweils zum Zeitpunkt der Bestellung gültigen Preise. Die unverbindlichen Preise sind Netto-Preise zuzüglich der gültigen gesetzlichen Mehrwertsteuer.

Zahlung

Die Rechnungsbeträge sind innerhalb von 30 Tagen nach Erhalt der Ware rein netto zur Zahlung fällig.

Erstausstattungsrabatt

Wir sind Ihnen bei der Erstausstattung Ihres neuen Labors gerne behilflich und entwickeln mit Ihnen gemeinsam eine optimale Start-Strategie.

Konsignationslager/ Freezer Programm

Wir stellen Ihnen gerne ein individuell zugeschnittenes Depot der von Ihnen häufig benötigten Produkte vor Ort zur Verfügung. Unmittelbarer Zugriff auf diese Produkte und minimaler Verwaltungsaufwand sind die Leitgedanken dieses NEB Freezer Programms. Bitte fordern Sie unsere detaillierten Unterlagen an!

Sonderformulierungen/ Mengenlieferungen/OEM

Wir bieten Ihnen unsere Katalogprodukte auch auf Ihre Bedürfnisse zugeschnitten in unterschiedlichen Mengen, Formulierungen oder Verpackungen an. Bitte richten Sie Ihre Anfragen an unsere kostenfreie Servicenummer oder an Email: custom de@neb.com.

Beschränkungen und Haftung

Die Produkte und Angebote von New England Biolabs GmbH richten sich ausschließlich an gewerbliche Nutzer, Universitäten, Forschungsund Bildungseinrichtungen etc., nicht jedoch an Verbraucher/Privatpersonen. NEB Produkte sind Forschungsreagenzien und ausschließlich für in vitro Anwendungen bestimmt. Die Verwendung durch ausgebildetes Laborpersonal wird vorausgesetzt. Weder NEB GmbH noch NEB Inc. haften für Schäden, die sich aus dem Gebrauch von NEB Produkten, direkt oder indirekt, ergeben oder ergeben könnten.

Markenrechte/Patente etc.

Die Produkte in diesem Katalog sind durch ein oder mehrere Patente, Markenrechte bzw. das Urheberrecht geschützt. Diese Rechte bzw. deren Rechtekontrolle obliegen NEB. Zwar entwickelt und validiert NEB seine Produkte für diverse Anwendungen, jedoch kann für den Einsatz in besonderen Anwendungen eventuell die Zustimmung Dritter (Patentrechtinhaber) notwendig sein.

Für weitere Informationen kontaktieren Sie bitte NEBs Global Business Development Team unter Email: gbd@neb.com.

Unsere vollständigen Allgemeinen Geschäftsbedingungen finden Sie im Internet unter www.neb-online.de.

Für Irrtümer oder Druckfehler im Katalog übernehmen wir keine Haftung.

Featured Products from NEB

NEBNext® Reagents for NGS Sample Prep

- Streamline workflows, minimize inputs, and improve library yields & quality with our growing portfolio.
- Avoid DNA damage generated with bisulfite sequencing! Our NEBNext Enzymatic Methyl-seq Kit (EM-seq™)generates high quality libraries for superior detection of 5-mC and 5-hmC from fewer sequencing reads.
- Generate high quality, full length transcript sequencing libraries from single cells or as little as 2 pg of total RNA with our NEBNext Single Cell/Low Input RNA Library Prep Kit.
- The NEBNext Ultra[™] II FS DNA Library Prep Kit incorporates our novel fragmentation reagent and offers a fast and reliable solution for library construction.

See page 134 for details





Monarch® Nucleic Acid Purification Kits

- Maximize performance and minimize your environmental impact.
- Choose from kits for genomic DNA & total RNA extraction, DNA & RNA cleanup, plasmid miniprep and gel extraction. Buffers and columns are also available separately.
- Obtain highly-pure DNA & RNA from a wide variety of sample types.
- Elute in small volumes and prevent buffer carryover with our unique column designs.
- Save time with fast, user-friendly protocols.

See page 122 for details

Cloning & Synthetic Biology

- Assemble multiple DNA fragments and transform in under two hours, regardless of fragment length or compatibility, with the NEBuilder HiFi DNA Assembly Kit. This versatile kit can be used for a variety of DNA assembly methods.
- Achieve 20+ fragment assembly with high efficiency and accuracy with our NEB Golden Gate Assembly Kit (Bsal-HF[®]v2). Couple with our Ligase Fidelity Viewer online tool to ensure the highest fidelity ligation.

See pages 88–90 for details



Table of Contents

Corporate Philosop Partnering with NE Restriction	В	2 4 ases		n-Profits & Research I dering Information	Foundations	6 6		Luna qPCR & RT-PCR Luna Universal qPCR Master Mix Luna Universal Probe qPCR Master Mix LunaScript RT SuperMix Kit Luna Universal One-Step RT-qPCR Kit	69 70 70 70 71
Aatll	BceAl		Dral	Hpy188III	Notl-HF	Sbfl-HF		Luna Universal Probe One-Step RT-qPCR Kit	71
AbaSI AccI Acc65I	Bcgl BciVl Bcll	BspDI	DrallI-HF Drdl Eael	HpyAV HpyCH4III HpyCH4IV	Nrul Nrul-HF Nsil	Scal-HF ScrFI SexAl		RT-PCR One Tag One-Step RT-PCR Kit	72
Accibil Acili Acili Acili Acili Acili Acili Afel Afili Afilii Afilii Agel Agel-HF Ahdi Alel-v2 Alui Alwi Alwi Alwi Alwi Apal Apal Apal Apal Apol-HF Asci AsiSi Aval Avali Bael	BCII BCII-HF BCODI BfaI BfuAI BgII BgIII BgIII BmgBI BmrI BmtI-HF BpmI Bpu10I BpuEI BsaI BsaI-HFv2 BsaAI BsaBI BsaHI BsaJI BsaWI BsaXI BsaXI BsaXI BsaXI BsaXI BsaXI BsaXI BsaSI BsaXI	BspHI BspMI BspQI BsrI BsrBI BsrBI BsrFI-v2 BsrGI BsrGI-HF BssHII BssSI-v2 BstAPI BstEII BstEII-HF BstVI BstVI BstVI BstVI BstVI BstZ171-HF Bsu361 BtgI BtgZI BtsI-v2	Eael Eagl Eagl Eagl-HF Earl Ecil Eco53kl EcoNI EcoO109I EcoP15I EcoRI-HF EcoRV EcoRV-HF Esp3I Fatl Faul Fnu4HI Fokl Fsel Fspl FspEl Haell HaellI Hgal Hincll	HyyCH4IV HpyCH4V KasI KpnI KpnI-HF LpnPI MboI MfeI MfeI-HF MluI MluI-HF MluCI MlyI MmeI MnII MscI MseI MspI MspA1I MspA1I MspJI NarI Ncil	Nsil-HF Nspil-HF Nspil-HF Paccl PaeR71 Pcil PfilFI PfilMI PleI PluTI Pmel PmlI PpuMI PshAI Psil PspGI PspOMI PsptI-HF PvuI PvuII-HF PvuII PvuII-HF Rsal RsrII	Sexal Stanii Stol Stol SgrAl Smal Smil SnaBi Spel Spel-HF Sphl-HF Stfl Sspl-HF Sttl Styl-HF Styl Styl-HF Styl Styl-HF StyD4l Swal Taq ^e l Tfil Tsel TspMI		One Taq RT-PCR Kit Isothermal Amplification and Strand Displacemen Bst DNA Polymerase, (Large Frag. and Full-Length) Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase Bst 3.0 DNA Polymerase WarmStart Colorimetric LAMP Master Mix (DNA & RNA) IsoAmp II Universal tHDA Kit phi29 DNA Polymerase Polymerases for DNA Manipulation (labeling, blunting, etc.) PreCR Repair Mix Sulfolobus DNA Polymerase IV Therminator DNA Polymerase DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Klenow Fragment (3'→5' exo-) T4 DNA Polymerase, Large Fragment Terminal Transferase Nucleotide Solutions & Buffers Reaction Buffers and Diluents	73 73 73 74 74 75 75 76 76 76 77 77 77 77 78 78 78
BaeGI	BsiEl	BtsCI	HindIII	Ncol	Sacl	TspRI		Acyclonucleotide Set	79 79
BamHI BamHI-HF BanI BanII BbsI BbsI-HF	BsiHKAI BsiWI BsiWI-HF BsII BsmI BsmAI	Clal CspCl CviAll CviKl-1	HindIII-HF HinfI HinP1I HpaI HpaII HphI	Ncol-HF Ndel NgoMIV Nhel Nhel-HF NlallI	SacI-HF SacII SaII SaII-HF SapI Sau3AI	Tth111I Xbal Xcml XhoI Xmal XmnI		Deoxynucleotide (dNTP) Solution Mix/Set Ribonucleotide Solution Mix/Set 7-deaza-dGTP Adenosine-5´ Triphosphate (ATP) 5-methyl-dCTP dATP Solution	79 79 79 79 79 79
Bbvl BbvCl Bccl High-Fidelity Re Methylation Ser AbaSl, Dpnl, Dpnll Nicking Endonu	BsmBl BsmFl BsoBl estriction Enzym nsitive Restriction , FspEl, Hpall, Lpr cleases	Ddel Dpnl Dpnll es on Enzymes for Ep iPl, McrBC, Mspl, M	Hpy99I Hpy166II Hpy188I igenetics Str	NIaIV NmeAIII NotI	Sau96l Sbfl	Zral 301 264–265		cDNA Synthesis cDNA Synthesis Selection Chart LunaScript RT SuperMix Kit ProtoScript II First Strand cDNA Synthesis Kit ProtoScript First Strand cDNA Synthesis Kit Template Switching RT Mix ProtoScript II Reverse Transcriptase M-MuLV Reverse Transcriptase AMV Reverse Transcriptase WarmStart RTx Reverse Transcriptase	80 70 194 194 193 192 192 192
Homing Endonu I-Ceul, I-Scel, PI-F	cleases	ssSI, Nd.BīsI, Nī.AIWI	, Nt.BbvCI, Nt.	BsmAI, Nt.BspQI, Nt.E	istnbi, nt.cvipii	53–55 55–56		PCR Cleanup Monarch PCR & DNA Cleanup Kit	81
NEBuffers, Diluents, Gel Loading Dyes, BSA & Rec Molecular Biology Grade, NEB Tube Opener			ecombinant /	Albumin,		56–57		Exo-CIP Rapid PCR Kit DNA Modifying Enzymes & Cloning Technologies	81
Polymerases & Amplification Technologies						66	Cloning & Synthetic Biology Cloning Workflow Comparison NEBuilder HiFi DNA Assembly Master Mix & Cloning Kit	86 88	
Q5U Hot Start High Q5 High-Fidelity P Q5 Site-Directed N (with or without co	INA Polymerase & Fidelity DNA Polymerase & Fidelity DNA Polymerase & Fidelity DNA Polymerase & Master Management Cells) Stilly DNA Polymerase & Master Management Cells Polymerase & Master Manageme	merase & Master Mixymerase ase & Master Mixes ase & Master Mix Mixes Master Mixes or Std. Buffer	62 63 63 63 63 63 92 64 64 64 65 65 66 66	Taq Master Mixes (2) Taq PCR Kit Multiplex PCR 5X Ma Hot Start Taq DNA PC Specialty PCR LongAmp Taq DNA PC LongAmp Taq PCR K Hemo KlenTaq EpiMark Hot Start Ta Other PCR Polyme Vent & Vent (exo") DI Deep Vent & Deep Vet qPCR & RT-qPCR What is Difference Be and Two-step RT-PCf	olymerase & A2X Nolymerase & Manag DNA Polymerasit 7 DNA Polymerases NA Polymerases nt (exo^) DNA Potween One-step	Master Mix ster Mix ase & 2X Master Mix ee	66 66 66 66 67 67 67 67 68 68	NEBuilder HiFi DNA Assembly Bundle for Large Fragments Gibson Assembly Master Mix & Cloning Kit NEB Golden Gate Assembly Kit (Bsal-HFv2) BioBrick Assembly Kit NEB PCR Cloning Kit (with and without competent cells) USER & Thermolabile USER II Enzyme Q5 Site Directed Mutagenesis Kit (with and without competent cells) Quick Blunting Kit DNA Ligases DNA Ligases Selection Chart T4 DNA Ligase Quick Ligase Selection Chart T4 DNA Ligase Quick Ligation Kit Instant Sticky-end Ligase Master Mix Blunt T/A Ligase Master Mix	93 94 95 96

ElectroLigase	96	EcoGII Methyltransferase	115	Methylation Analysis	
T3 DNA Ligase	96	EcoRI Methyltransferase	115	NEBNext Enzymatic Methyl-seq Kit	141
7 DNA Ligase	96	HaeIII Methyltransferase	115	NEBNext Enzymatic Methyl-seq Conversion Module	141
HiFi <i>Taq</i> DNA Ligase	97	Hhal Methyltransferase	115	NEBNext Q5U Master Mix	141
Taq DNA Ligase	97	Hpall Methyltransferase	116	NEBNext Multiplex Oligos for Enzymatic Methyl-seq	
E. coli DNA Ligase	98	Mspl Methyltransferase	116	(Unique Dual Index Primer Pairs)	141
9°N DNA Ligase	98	Taql Methyltransferase	116	Ultra II RNA	
SplintR Ligase	199	Human Dnmt1 Methyltransferase	270	NEBNext Ultra II Directional RNA Library Prep Kit	
DNA Kinases		Single-stranded DNA Binding Proteins		with or without Sample Purification Beads	142
74 Polynucleotide Kinase	98	RecA, RecA _f	117	NEBNext Ultra II RNA Library Prep Kit with or without Sample Purification Beads	1/10
4 Polynucleotide Kinase (3´ phosphatase minus)	98	T4 Gene 32 Protein	117	with or without Sample Furnication beaus	142
5-hydroxymethyluridine Kinase	99	ET SSB (Extreme Thermostable ssBinding Protein)	117	Single Cell RNA	
ONA Phosphatases & Sulfurylases		Cloning Plasmids & DNAs		NEBNext Single/Cell Low Input RNA Library	
Phosphatase Selection Chart	99	Cloning Plasmids (pBR322, pUC19)	118	Prep Kit	144
Quick Dephosphorylation Kit	99	DNAs(M13mp18 DNA, Lambda, Lambda		NEBNext Single/Cell Low Input cDNA Synthesis & Amplification Module	144
Shrimp Alkaline Phosphatase (rSAP)	100	(N ⁶ -Methyladenine-free), φX174)	118	NEBNext Single Cell Lysis Module	144
Antarctic Phosphatase	100	Other		,	
Alkaline Phosphatase (CIP)	100	M13K07 Helper Phage	118	Small RNA	4.45
Pyrophosphatase, Inorganic (<i>E. coli</i> , yeast)	189	Competent Cell Selection Chart for Cloning	119	NEBNext Multiplex Small RNA Library Prep Set (Set 1, 2)	145
Thermostable Inorganic Pyrophosphatase	189 101	Monarch Nucleic Acid Purification Kits	119	NEBNext Multiplex Small RNA Library Prep Kit (Index Primers 1-48)	145
Apyrase <i>Tte</i> UvrD Helicase	101	Nuclaia Asid Durification		NEBNext Small RNA Library Prep Set (Multiplex Compatible	
	101	Nucleic Acid Purification			, 1 10
Nucleases		Plasmid Purification		Modules & Enzymes	
Properties of Exonucleases and Endonucleases	102	Monarch Plasmid Miniprep Kit	125	NEBNext Ultra II FS DNA Module	140
Common Apps for Exonucleases & Endonucleases	103	Companion Products		3	, 143
DNase I (RNase-Free) Exonuclease I	202 103	Monarch Plasmid Miniprep Columns	125	NEBNext Ultra II Ligation Module 140 Ultra II RNA First Strand Synthesis Module	, 143 143
Exonuclease i Thermolabile Exonuclease I	103	Monarch Plasmid Buffers (Resuspension (B1), Lysis (B2),		Ultra II Directional RNA Second Strand Synthesis Module	
Exonuclease III	103	Neutralization (B3), Wash 1, Wash 2)	125	Ultra II Non-directional RNA Second Strand	170
Exonuclease V (RecBCD)	104	Monarch DNA Elution Buffer	125	Synthesis Module	143
Exonuclease VII	104	DNA Cleanup		NEBNext Ultra II Q5 Master Mix	155
Exonuclease VIII, truncated	104	Monarch DNA Gel Extraction Kit	126	Adaptors & Primers	
Exonuclease T	105	Monarch PCR & DNA Cleanup Kit (5 μg)	127	NEBNext Multiplex Oligos	
Thermostable FEN1	105	Companion Products		(96 Unique Dual Index Primer Pairs)	146
ambda Exonuclease	105	•	, 127	NEBNext Multiplex Oligos (Dual Index Primers Set 1, 2)	
Micrococcal Nuclease	105	Monarch Gel Dissolving Buffer	126	NEBNext Multiplex Oligos (Index Primers Set 1, 2, 3, 4)	
Mung Bean Nuclease	106	Monarch DNA Buffers		NEBNext Multiplex Oligos (96 Index Primers)	146
Nuclease BAL-31	106	(Cleanup Binding, Wash, Elution) 126	, 127	NEBNext Multiplex Oligos	
Nuclease P1	106 106	Genomic DNA Purification		(Methylated Adaptor, Index Primers Set 1)	146
RecJ _f 75 Exonuclease	106	Monarch Genomic DNA Purification Kit	128	NEBNext Adaptor Dilution Buffer	146
7 Exonuclease	107		.20	Target Enrichment	
Nucleoside Digestion Mix	202	Companion Products	400	NEBNext Direct Cancer HotSpot Panel	150
·		Monarch gDNA Purification Columns	128 128	NEBNext Direct BRCA1/BRCA2	150
ONA Repair Proteins	107	Monarch Collection Tubes II Monarch gDNA Buffers (Tissue Lysis, Cell Lysis,	120	NEBNext Direct Custom Ready Panels	151
Properties of DNA Repair Enzymes DNA Repair Glycosylases on Various	107	Blood Lysis, Binding, Wash, Elution)	128	•	
Damaged Bases	108	Monarch RNase A	128	Library Quantification NEBNext Library Quant Kit	156
APE 1	109	Proteinase K, Molecular Biology	128	NEBNext Library Dilution Buffer	156
Endonuclease IV	109	RNA Purification		•	100
Tth Endonuclease IV	109	Monarch Total RNA Miniprep Kit	130	Reagents for Ion Torrent	
Endonuclease III (Nth)	109	Monarch RNA Cleanup Kit (10 μg, 50 μg, 500 μg)	131	NEBNext Fast DNA Library Prep Set	157
Tma Endonuclease III	110	1 (13, 13, 13,		NEBNext Fast DNA Fragmentation & Library Prep Set	157
Endonuclease V	110	Companion Products	100	Reagents for rRNA Depletion	
Endonuclease VIII	110	Monarch RNA Purification Columns	130	NEBNext rRNA Depletion Kit (Human/Mouse/Rat)	
Fpg DAAC	110 111	Monarch gDNA Removal Columns Monarch DNA/RNA Protection Reagent	130 130	with or without RNA Sample Purification Beads	147
nAAG F4 PDG (T4 Endonuclease V)	111	Monarch RNA Buffers (Lysis, Priming, Wash,	130	NEBNext Poly(A) mRNA Magnetic Isolation Module	147
Jracil-DNA Glycosylase (UDG)	111		, 131	NEBNext Globin RNA & rRNA Depletion Kit (Human/Mouse/Rat) with or without RNA	
Afu Uracil-DNA Glycosylase (UDG)	111	Monarch Total RNA Miniprep Enzyme Pack	130	Sample Purification Beads	148
Antarctic Thermolabile UDG	112	Monarch RNA Cleanup Columns (10 μg, 50 μg, 500 μg)	131	NEBNext Magnetic Separation Rack	148
nSMUG1	112	Monarch Collection Tubes II 130	, 131	Descents Suitable for any Blotform	
PreCR Repair Mix	76	Nuclease-free Water 130	, 131	Reagents Suitable for any Platform DNA Separation	
Other Enzymes		Accessories		NEBNext Microbiome DNA Enrichment Kit	152
3-Agarase I	112	Monarch Microfuge Tube EcoRack	127		
Cre Recombinase	113	•		DNA Repair	150
7 Endonuclease I	113	NEBNext Reagents for		NEBNext FFPE DNA Repair Mix	153
TelN Protelomerase	113	Next Generation Sequencing		Modules & Enzymes	
「opoisomerase I (<i>E. coli</i>)	113	Reagents for Illumina Sequencing		NEBNext Magnesium RNA Fragmentation Module	143
DNA Methyltransferases			, 137	NEBNext dsDNA Fragmentase	154
CpG Methyltransferase (M.SssI)	114			NEBNext dsDNA Fragmentase Reaction Buffer v2	154
GpC Methyltransferase (M.CviPI)	114	Ultra II DNA NEPNovt Ultra II DNA Library Prop Kit		NEBNext Ultra II Q5 Master Mix	155
Alul Methyltransferase	114	NEBNext Ultra II DNA Library Prep Kit with or without Sample Purification Beads	139		
BamHI Methyltransferase	114	NEBNext Ultra II FS DNA Library Prep Kit	100		
dam Methyltransferase	115	with or without Sample Purification Beads	139		

Table of Contents (Cont.)

`	,				
Markers and Ladders		Pyrophosphatase, Inorganic (<i>E. coli</i> , yeast)	189	Protein Expression &	
(DNA, RNA & Protein)		Thermostable Inorganic Pyrophosphatase	189	Purification Technologies	
		Vaccinia Capping System	190	Protein Expression Selection Chart	212
DNA Markers & Ladders DNA Ladders		mRNA Cap 2´-O-Methyltransferase	190	·	212
Size Ranges of DNA Ladders	166	RNA Cap Analogs	190	Expression Systems	
Gel Loading Dye, Purple (6X) (with or without SDS)	166	3´-Desthiobiotin-GTP, 3´-Biotin-GTP	191	E. coli pMAL Protein Fusion & Purification System	213
1 kb DNA Ladder	167	yDcpS	191	pivial Frotein Fusion & Furnication System	213
100 bp DNA Ladder	167	cDNA Synthesis, RT-PCR & RT-qPCR		pMAL Companion Products:	
1 kb Plus DNA Ladder	167	cDNA Synthesis Selection Chart	191	Amylose Resin, Amylose Resin High Flow	213
50 bp DNA Ladder	167	ProtoScript II Reverse Transcriptase	192	Amylose Magnetic Beads	213
Fast DNA Ladder	167	LunaScript RT SuperMix Kit	70	Anti-MBP Magnetic Beads	213
Low Molecular Weight DNA Ladder	167	M-MuLV Reverse Transcriptase	192	Anti-MBP Monoclonal Antibody	213
PCR Marker	167	AMV Reverse Transcriptase	192	Factor Xa Protease	213
Quick-Load Purple 1 kb Plus DNA Ladder	168	WarmStart RTx Reverse Transcriptase	193	pMAL Vectors	213
Quick-Load Purple 1 kb DNA Ladder	168	Template Switching RT Enzyme Mix	193	IMPACT Kit	214
Quick-Load Purple 100 bp DNA Ladder	168	Primers for cDNA Synthesis	193	IMPACT Vector Comparison	214
Quick-Load Purple 50 bp DNA Ladder	168	ProtoScript II First Strand cDNA Synthesis Kit	194	IMPACT Companion Products:	
Quick-Load Purple Low Molecular Weight DNA Ladder	168	ProtoScript First Strand cDNA Synthesis Kit	194	Anti-CBD Monoclonal Antibody	214
Quick-Load DNA Ladders	168	Luna Universal One-Step RT-qPCR Kit	71	Chitin Resin	214
Quick-Load 1 kb Extend DNA Ladder	168	Luna Universal Probe One-Step RT-qPCR Kit	71	IMPACT Vectors	214
TriDve DNA Ladders	168	One Taq One-Step RT-PCR Kit	72		
PFG Ladders	168	One Taq RT-PCR Kit	72	Yeast	045
Supercoiled DNA Ladder	170	RNA Ligases & Modifying Enzymes		K. lactis Protein Expression Kit	215
·		Reported Activities for T4 RNA Ligases	195	K. lactis Companion Products:	
Conventional DNA Markers	1) 100	RNA Ligase Selection Chart	196	BstXI, SacI	32, 47
Lambda DNA (Mono Cut Mix, HindIII Digest, BstEII Digest	,	T4 RNA Ligase 1	197	Yeast Medium Pack	215
φX174 DNA–HaelII Digest	169	T4 RNA Ligase 2	197	Cell-free Expression	
pBR322 DNA (BstNI Digest, MspI Digest)	169	T4 RNA Ligase 2, truncated	197	PURExpress <i>In Vitro</i> Protein Synthesis Kit	216
RNA Markers & Ladders		T4 RNA Ligase 2, truncated K227Q and KQ	198	PURExpress Δ Ribosome Kit	216
dsRNA Ladder	170	RtcB Ligase	198	PURExpress Δ (aa, tRNA) Kit	216
microRNA Marker	170	Thermostable 5´ AppDNA/RNA Ligase	198	PURExpress \(\Delta \) (aa, thive) Kit	216
ssRNA Ladder	170	5' DNA Adenylation Kit	199	,	210
Low Range ssRNA Ladder	170	SplintR Ligase	199	PURExpress Companion Products:	
Dustain Maulaus 9 Laddaus		RNA 5´ Pyrophosphohydrolase (RppH)	199	PURExpress Disulfide Bond Enhancer	216
Protein Markers & Ladders	171	5´ Deadenylase	199	E. coli Ribosome	216
Unstained Protein Standard, Broad Range (10–200 kDa)	171	RNase I,	200	Competent Cells	
Color Prestained Protein Standard, Broad Range (10–250 kDa)	171	RNase H, Thermostable RNase H	200	Competent Cells for Protein Expression	229-233
Blue Prestained Protein Standard, Broad Range	17.1	RNase HII	200	Competent Cell Expression Strain Selection Chart	217
(11–250 kDa)	171	Quick Dephosphorylation Kit	99	Magnetic Metrices and Dooks	
(11 200 Kbu)		Antarctic Phosphatase	100	Magnetic Matrices and Racks	218
Genome Editing		Alkaline Phosphatase Calf Intestinal (CIP)	100	Protein A Magnetic Beads	218
Featured Products Supporting CRISPR Workflows		Shrimp Alkaline Phosphatase (rSAP)	100	Protein G Magnetic Beads Streptavidin Magnetic Beads	218
EnGen Spy Cas9 NLS	• 177	T4 Polynucleotide Kinase	98	Hydrophilic Streptavidin Magnetic Beads	218
EnGen Mutation Detection Kit	178	ShortCut RNase III	201	Goat Anti-Rabbit IgG Magnetic Beads	219
EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i>	179	XRN-1	201	Goat Anti-Mabbit 19d Magnetic Beads	219
EnGen Spy Cas9 Nickase	177	RNase Control		Goat Anti-Rat IgG Magnetic Beads	219
EnGen Spy dCas9 (SNAP-tag)	177	RNase Inhibitor, Murine/ Human Placenta	202, 203	Chitin Magnetic Beads	219
EnGen Lba Cas12a (Cpf1)	177	Ribonucleoside Vanadyl Complex	202, 203	Amylose Magnetic Beads	219
EnGen Sau Cas9	177	Tilbondolcoside Valladyi Complex	200	Ni-NTA Magnetic Beads	219
Cas9 Nuclease, S. pyogenes	177	RNA Library Prep for Next Gen Squencing		Anti-MBP Magnetic Beads	219
Monarch Total RNA Miniprep Kit	130	NEBNext Ultra II Directional RNA Library Prep Kit		Magnetic Separation Racks (6-tube, 12-tube, 50 ml)	219
Q5 Site-directed Mutagenesis Kit		with or without Sample Purification Beads	142	96-Well Microtiter Plate Magnetic Separation Rack	219
(with or without competent cells)	92	NEBNext Ultra II RNA Library Prep Kit	4.40	SNAP-Capture Magnetic Beads	283
Q5 High-Fidelity DNA Polymerases	63	with or without Sample Purification Beads	142		
NEBuilder HiFi DNA Assembly Master Mix & Cloning Kit	88	NEBNext Single Cell/Low Input RNA Library Prep Kit	144	Competent Cells	
HiScribe T7 ARCA mRNA Kit (with or without tailing)	186	NEBNext Single Cell/Low Input cDNA	144	Strain Properties Chart	224
HiScribe T7 High Yield RNA Synthesis Kit	185	Synthesis & Amplification Module Exonuclease T	202	Convenient Formats	225
HiScribe T7 Quick High Yield RNA Synthesis Kit	186	Nucleoside Digestion Mix	202		LLO
HiScribe SP6 RNA Synthesis Kit	185	DNase I (RNase-free)	202	Cloning Strains	000
T7 Endonuclease I	113	Divase i (ilivase-liee)	202	NEB Cloning Competent <i>E. coli</i> Sampler	226
DNA D		RNA Cleanup & Isolation		NEB Turbo Competent E. coli	226
RNA Reagents		EpiMark N6-Methyladenosine Enrichment Kit	204	(High Efficiency) Electrocompetent	226
Avoiding Ribonuclease Contamination	184	p19 siRNA Binding Protein	204	'	220
RNA Synthesis		Magnetic mRNA Isolation Kit	204	NEB 10-beta Competent <i>E. coli</i>	227
HiScribe T7 Quick High Yield RNA Synthesis Kit	185	Oligo d(T) ₂₅ Magnetic Beads	205	(High Efficiency)	227
HiScribe T7 High Yield RNA Synthesis Kit	185	Oligo d(T) ₂₅ Cellulose Beads	205	Electrocompetent	221
HiScribe SP6 RNA Synthesis Kit	185	polyA Spin mRNA Isolation Kit	205	NEB 5-alpha Competent <i>E. coli</i>	227
HiScribe T7 ARCA mRNA Kit (with or without tailing)	186	Streptavidin Magnetic Beads	205	(High Efficiency)	227
EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i>	179	Monarch Total RNA Miniprep Kit	130	(Subcloning Efficiency)	227
Recommended HiScribe Kits by Application	187	Monarch RNA Cleanup Kit (10 μg, 50 μg, 500 μg)	131	Electrocompetent	221
T3, T7, SP6 RNA Polymerase	188	RNA Markers & Ladders		NEB 5-alpha F [*] / ₁ Competent <i>E. coli</i>	227
Hi-T7 RNA Polymerase	188	dsRNA Ladder	207	(High Efficiency)	221
E. coli Poly(A) Polymerase	188	microRNA Marker	207	NEB Stable Competent <i>E. coli</i>	228
Poly(U) Polymerase	188	ssRNA Ladder	207	(High Efficiency)	
E. coli RNA Polymerase, (Core and Holo Enzyme)	189	Low Range ssRNA Ladder	207	dam ⁻ /dcm ⁻ Competent <i>E. coli</i>	228
Ribonucleotide Solution Mix/Set	189	RNA Loading Dye (2X)	207		
	100	Universal miRNA Cloning Linker	207		

Universal miRNA Cloning Linker

Protein Expression Strains		Furin	255	Biotin Labels, Vista Label	282
BL21, BL21(DE3) Competent E. coli	229	TEV Protease	255	SNAP-Capture, Building Blocks	283
Lemo21(DE3) Competent E. coli	229	Protein Phosphatases		Deference Annandiy	
NiCo21(DE3) Competent E. coli	230	Properties of Protein Phosphatases from NEB	256	Reference Appendix	
NEB Express Competent E. coli	230	Lambda Protein Phosphatase (Lambda PP)	256	Online Interactive Tools	288, 289
(High Efficiency)	230	<i>p</i> -Nitrophenylphosphate (PNPP)	256	Databases, Mobile Apps & Freezer Program Locator	289
NEB Express ¹⁹ Competent <i>E. coli</i> (High Efficiency)	230	Sodium Orthovanadate	256	Restriction Endonucleases	
T7 Express Competent <i>E. coli</i>	200	Protein Kinases		Optimizing Restriction Enzyme Reactions	290
(High Efficiency)	231	cAMP-dependent Protein Kinase (PKA)	257	Double Digestion	291
T7 Express /ysY Competent E. coli		Casein Kinase II (CK2)	257	Types of Restriction Enzymes	291
(High Efficiency)	231	Adenosine 5´ Triphosphate	257	Troubleshooting Guide	292
T7 Express <i>lysY/l^q</i> Competent <i>E. coli</i>				Performance	
(High Efficiency)	231	Epigenetics		Performance Chart for Restriction Enzymes	293-298
SHuffle Express Competent <i>E. coli</i>	232	Epigenetics		Activity of Enzymes at 37°C	299
SHuffle T7 Express Competent <i>E. coli</i>	232	EpiMark 5-hmC and 5-mC Analysis Kit	262	Activity of DNA Modifying Enzymes in CutSmart Buffe	er 299
SHuffle T7 Competent <i>E. coli</i>	233	T4 Phage β-glucosyltransferase	262	Tips for Avoiding Star Activity	300
SHuffle T7 Express <i>lysY</i> Competent <i>E. coli</i>	233	EpiMark Methylated DNA Enrichment Kit	263	High-Fidelity (HF) Restriction Enzymes	301
How do I express my protein in SHuffle?	233	EpiMark N6-Methyladenosine Enrichment Kit	204	Reduced Star Activities of HF Enzymes	301
Glycobiology & Protein Tools		EpiMark Bisulfite Conversion Kit	263	Time-Saver Qualified Restriction Enzymes	302–303
		EpiMark Hot Start <i>Taq</i> DNA Polymerase	68	Isoschizomers/Recognition Sequences	
Endoglycosidases	000	Methylation-Dependent Restriction Enzymes		Cross Index of Recognition Sequences	304-306
Endo F2	238 238	AbaSI, FspEI, LpnPI, MspJI	264	Palindromic Penta-nucleotide Recognition Sequences	305
Endo F3 Endo D	239	McrBC	265	Specificities Greater Than 6 Bases	305
Endo H, Endo H,	239	Additional Restriction Enzymes for Epigenetic Ana	lvsis	Interrupted Palindromes	305-306
Endo S	239	Dpnl, Dpnll	34	Multiple & Nonpalindromic Recognition Sequences	306-309
PNGase A	240	Hpall	38	Alphabetized List of NEB Recognition Sequences	309-310
PNGase F and PNGase F, recombinant	240	Mspl	41	Isoschizomers	311–327
RNase B	240	·		Properties of Restriction Enzymes	
Rapid PNGase F and Rapid PNGase F (non-reducing format)		Histones	067	Survival in a Reaction	328-329
Rapid PNGase F Antibody Standard	241	EpiMark Nucleosome Assembly Kit Histone H1 ^o Human, Recombinant	267 267	Cleavage of Supercoiled DNA	329
Remove-iT PNGase F	241	Histone H3.1/H4 Tetramer Human, Recombinant	267	Generating New Cleavage Sites	
O-Glycosidase	242	Histone H2A/H2B Dimer Human, Recombinant	267	Recleavable Filled-in 5´ Overhangs	330
Protein Deglycosylation Mix II	242	Nucleosome Control DNA	267	Recleavable Blunt Ends	331
Fetuin	242	Histone H2A, H2B Human, Recombinant	268	Compatible Cohesive Ends	332–333
Exoglycosidases		Histone H3.1, H3.2, H3.3 Human, Recombinant	268	·	
N-Glycan Sequencing Kit	243	Histone H4 Human, Recombinant	269	Methylation Sensitivity	004 000
α-N-Acetylgalactosaminidase	243	DNA Mothyltronoforago		Dam, Dcm and CpG Methylation	334–336
B-N-Acetylglucosaminidase S	243	DNA Methyltransferases CpG Methyltransferase (M.SssI)	114	PCR, qPCR & RT-qPCR	
B- <i>N</i> -Acetylhexosaminidase,	244	GpC Methyltransferase (M.CviPI)	114	General Guidelines for PCR Optimization	337
α1-2 Fucosidase	244	Alul Methyltransferase	114	PCR Troubleshooting Guide	338
∝1-3,4 Fucosidase	244	BamHI Methyltransferase	114	Optimization Tips for Luna qPCR	339
α1-2,3,4,6 Fucosidase	244	dam Methyltransferase	115	Optimization Tips for Luna One-Step RT-qPCR	340
α1-2,4,6 Fucosidase 0	245	EcoGII Methyltransferase	115	Luna qPCR Troubleshooting Guide	341
α1-3,6 Galactosidase	245	EcoRI Methyltransferase	115	Luna One-Step RT-qPCR Troubleshooting Guide	342
cx1-3,4,6 Galactosidase	245	Haelli Methyltransferase	115	Cleavage Close to the End of DNA Fragments Activity of Restriction Enzymes in PCR Buffers	343 344–345
B1-3 Galactosidase B1-3.4 Galactosidase	245 246	Hhal Methyltransferase	115	Activity of nestriction Enzymes in Fon Duners	344-343
B1-4 Galactosidase S	246	Hpall Methyltransferase	116	Cloning and Transformation	
α1-2,3 Mannosidase	246	Mspl Methyltransferase	116	Getting Started with Molecular Cloning	346
ca1-2,3,6 Mannosidase	246	Human DNA (cytosine-5) Methyltransferase (Dnmt1)	270	Traditional Cloning Quick Guide	347–348
∞1-6 Mannosidase	247	Control DNA		Troubleshooting Guide for Cloning	349–351
cc2-3,6,8 Neuraminidase	247	HeLa Genomic DNA	271	Optimization Tips for Your Cloning Reaction	352–353
α2-3,6,8,9 Neuraminidase A	247	NEDNovt Boogonto for ChID Son Library Bron		Troubleshooting Guide for DNA Cleanup & Plasmid Purification Using Monarch Kits	354
cx2-3 Neuraminidase S	247	NEBNext Reagents for ChIP-Seq Library Prep 5-methyl-dCTP	271	Choosing Sample Input Amounts	001
Endoglycoceramidase I (EGCase I)	257	NEBNext Ultra II DNA Library Prep Kit for Illumina	211	When Using the Monarch gDNA Purification Kit	355
Heparin Lyases		with or without Sample Purification Beads	139	Troubleshooting Guide for Genomic DNA	
Bacteroides Heparinase I, II, III	249	NEBNext Multiplex Oligos	146	Purification Using Monarch Kits	356
•		NEBNext Ultra II End Repair/dA-Tailing Module	140	Choosing Sample Input Amounts	
Phage Display	050	NEBNext Ultra II Ligation Module	140	When Using the Monarch Total RNA Miniprep Kit	357
Ph.D7 Phage Display Peptide Library Kit Ph.D12 Phage Display Peptide Library Kit	250	NEBNext Ultra II Q5 Master Mix	155	Troubleshooting Guide for Total RNA Extraction	250
Ph.DC7C Phage Display Peptide Library Kit	250 250	NEBNext Enzymatic Methyl-seq Kit	273	& Purification Using Monarch Kits Genetic Markers	358 359–360
Ph.D12 Library	250	NEBNext Enzymatic Methyl-seq Conversion Module	273	Enhancing Transformation Efficiency	361
Ph.D. Peptide Display Cloning System	251	NEBNext Q5U Master Mix	273	Protein Expression with T7 Express Strains	362
	201	NEBNext Multiplex Oligos for Enzymatic	116	·	
Proteases & Proteome Analysis	0.40	Taql Methyltransferase Methyl-seq (Unique Dual Index Primer Pairs)	116 273	NGS Library Preparation	000
IdeZ Protease	248	woonyr-oog (onique duai illuex i illilet Falis)	210	DNA/RNA Input Guidelines for Library Prep	363
Trypsin-digested BSA MS Standards (CAM Modified)	252	Cellular Analysis		Protein Expression	
Trypsin-ultra, Mass Spectrometry Grade Endoproteinase GluC	252 253	Cellular Imaging & Analysis		Labeling with SNAP-tag Technology	
Endoproteinase Gido Endoproteinase AspN	253	Overview of SNAP- and CLIP-tag Technologies	278	Troubleshooting Guide	364
Endoproteinase Aspiv Endoproteinase LysC	253	Comparison of SNAP-tag/CLIP-tag Technology to GFP	279	Frequencies of Restriction Sites in Sequenced DNAs	365–368
Proteinase K, Molecular Biology Grade	254	SNAP-tag Starter Kit	279	DNA Maps	369–381
Thermolabile Proteinase K	254	Fluorescent Substrates	280	Odds and Ends	
Factor Xa Protease	254	Blocking Agents	281	General Molecular Biology Information	382-385
Enterokinase, light chain	255	SNAP-tag Antibodies, Purified Protein	281	Index	386
		Cloning Vectors	282	III WA	500

Application Reference Guide

DNA Cloning		SNP Analysis	76, 100, 103	Genome Editing	
•	78	Synthetic Biology	88–91	Cas9-based Genome Editing	
A-tailing Alkaline Elution	78 109–111	T/A Cloning	86, 95, 96	Mutation Detection) 176–177 178
Alkaline Unwinding	109–111	Transcriptional Mapping	106	sgRNA Synthesis	179
Chromosomal DNA Removal	103–111	Transformation Tips	361–362	Synna Synnesis	119
Cleavage Close to the Ends	343	Transformation/Competent cells	119, 222–233		
Cloning Optimization	352–353	USER Cloning	91	Epigenetics	
Cloning Plasmids	118	Vectors	118	5-mC & 5-hmC Detection, E	nrichment
Cloning Quick Guide	347–348	Vector Maps	369-381	& Quantitation	262–263, 273
Cloning Troubleshooting Guide	349–350			5-hmC Labeling	262
Cloning Workflow	86			Bisulfite Sequencing	263, 273
DNA Assembly	88–91	DNA Amplification	& PCR	ChIP Assay	267
DNA Cleanup	81, 119, 127	cDNA Synthesis	80, 191-194	ChIP-Seq Sample Preparation	
DNA Digestion	16–57	Chain Termination	76, 79		146, 155, 272
DNA Extraction from Agarose	112, 119, 126	Colorimetric LAMP	74	Combined Bisulfide Restrict Analysis (COBRA)	ion 271
DNA Ligase Selection Chart	93, 94	DNA Polymerases	60–81	Control DNA Substrates	271
DNA Maps	369-381	Damaged Templates	76, 153	Glucosylation of 5-hmC	262
DNA Markers and Ladders	164-170	DNA Removal	100, 103–107	HAT Assay	267
DNA Repair		DNA Repair		Histone Modification Studie	
Enzymes	107-113	Polymerases	76	Immunodetection of 5-hmC	262
Polymerases	76	Repair Mix	76, 153	m6A Enrichment	204
Removal	103, 202	DNA Sequencing	134–161	Methyl-Seq	273
Repair Mix	76, 153	Genome Amplification	74	Methylation Dependent Clea	
Electrophoresis		Helicase Dependent Amplification	74	Methylation Sensitive Restri	
Loading Dye	57, 166	Incorporation of Modified Nucleotic		Enzyme Cleavage	17, 264, 334–336
Markers & Ladders	164–171	Isothermal Amplification	73, 74	Methylation Sensitive Single	
End Treatments		LAMP	73, 74, 101	Primer Extension (Ms-SNuF	PE) 271
Blunting	77–78, 92, 105, 106	Lesion Bypass	76 134–161	Methylation Specific PCR (M	,
dA-Tailing	78	Library Construction Next Generation Sequencing	134–161	Methyltransferases	270
Dephosphorylation	99–100	Nick Translation	73, 77	Nucleosome Assembly	266–269
End Labeling	98	Nucleic Acid Purification	122–131	SNP Analysis	271
Phosphorylation	98	Nucleotides	79	Southern Blotting	271
Genomic DNA Purification	128–129	One-Step RT-PCR	69, 72		
Gibson Assembly	88–89	PCR	03,72	DNA Analysis	
Golden Gate Assembly	90	AT-Rich Targets	62-65, 68	RNA Analysis	100 100
Library Construction	134–161, 203, 273	Cleanup	81, 119, 127	Adaptor Ligation	196–199
Ligase Chain Reaction	97, 98	Colony PCR	62, 65–66	Blunting	202 190–191
Ligase Detection Reaction	97, 98 93–98	Extraction-free PCR	62, 67	Capping cDNA Synthesis	80, 191–194
Ligation Selection Chart	93–96 93, 94	Fast	62–64	Dephosphorylation	99–100
Ligation Independent Cloning	93, 94	GC-Rich Targets	62, 63, 65, 67, 68	m6A Enrichment	204
Methyltransferases	114–116, 270	High-Fidelity	62-64, 92, 155	Gene Silencing	201
Nick Repair	95, 96, 97, 98	High Throughput	62-64, 68	In vitro Translation	185–189
Nicking Endonucleases	17, 53–55	Long Amplicon	62-64, 67	Isolation and Enrichment	204–205
Nucleic Acid Purification	81, 118, 122–131	Multiplex PCR	62-66	Labeling	188, 190, 197
Padlock Probes	97	Optimization	337	Library Construction	137, 142–144, 146, 147,
PCR Cloning	86, 91	Routine	65-66		160–161, 203
Plasmid Purification	119, 125	Selection Chart	62	Ligation	195-198
Polymerases	62–68	Troubleshooting	338	Activity Chart	195
Recombination	113, 117	Probe Labeling	77	Selection Chart	196
Recombinational Cloning	86	qPCR	69, 70	Loading Dye	207
Restriction Digests		Optimization	339	mRNA Isolation	204-205
Activity in PCR Buffers	344-345	Troubleshooting	341	Phosphorylation	98
Double Digestion	291	RT-LAMP	73, 193	Poly(A) Tailing	188
High Fidelity Restriction Enzy	mes 17, 301	RT-PCR	69, 72	Poly(U) Tailing	188
Methylation Dependent	264-265	RT-qPCR	69, 70–71	Probe Generation	188
Methylation Sensitivity	334–336	Optimization	340	RACE	193
Optimization	290	Troubleshooting	342	Reverse Transcription	80, 191–194
Restriction Enzymes	16–57	Second Strand cDNA Synthesis	77	RNAi	201
Time Saver Qualified Restrict		SNP Analysis	76, 271	RNA Cleanup	131, 206
Enzymes	302–303	ssDNA Binding Proteins	117	RNA Library Propagation	188
Troubleshooting	292	Strand Displacement Amplification		RNA Library Preparation	137, 142–144, 146, 147, 160–161, 203
Seamless Cloning	86, 88, 90	Nicking Enzymes Site Directed Mutagenesis	53–55 88, 92	RNA Markers and Ladders	170, 207
Single Cell Gel Electrophoresis	110–111	Two Step RT-PCR	69, 72	RNA Miniprep	130, 206
Site Directed Mutagenesis	88, 92, 104	Whole Genome Amplification	69, 72 75	RNA Removal	200, 202
Site Mapping	106	Tribio denome Ampinication	13	RNA Synthesis	185–189
					100 100

RNase Inhibit	ors 202–203
RT-PCR	72
RT-qPCR	71
Sequencing	137, 142–144, 146, 147, 160–161, 203
siRNA	201, 204
Tailing	188

Sample Prep for Next Gen Sequencing

Work don boquonong					
Adenylation of ssDNA Linkers 1					
Barcodes 14	46				
Enzymatic Fragmentation 13	39				
FFPE-treated DNA 15	53				
Library Construction					
Amplification 15	55				
ChIP-Seq 136, 139, 140, 146, 155, 2	72				
DNA 136, 138–141, 146, 152–15	59				
dsDNA Fragmentation 139, 140, 15	54				
mRNA 137, 142–144, 146–147, 160–161, 20	03				
Single Cell RNA 14	44				
Small RNA 145, 160–16	61				

Library Quantitation	156
Magnetic Separation	148
Microbiome Enrichment	152
Sample Preparation	134-161
Target Enrichment	149-151

Protein Expression & Analysis

Expression	210-219
Bacterial	213-214
Cell-free	216
Difficult to Express	
Proteins	213, 216, 217, 229–233
Selection Chart	212
Yeast	215
Glycosylation	236-257
Endoglycosidases	238-243
Exoglycosidases	243-247
MS Standards	252
Peptide Ligation	214
Phage Display	250-251
Proteases	248, 252-255

Protein Labeling	214, 276-283
Protein Standards	171
Protein Phosphorylation	256-257
Protein Purification	
Expression & Purification Kits	213-216
Magnetic Matrices and Racks	218-219
Selection Chart	212, 217
Pull-Down Assays	288
Transformation	
Competent Cells	217, 229–233
T7 Expression Strains	231
Troubleshooting & Tips	361-362
Yeast	215
Strain Properties	224

Cellular Analysis

Protein Labeling (e.g., SNAP-tag)	276-283
Protein Localization	276-283
Pulse Chase	276-283
Receptor Internalization	276-283
Protein Pull-Down Assays	283
Protein Detection (SDS-PAGE)	282



Conservation of Biodiversity

Each edition of the New England Biolabs Catalog contains a collection of mini-reviews that addresses various scientific, environmental and/or humanitarian topics. The theme of the 2019-20 Catalog is "Conservation of Biodiversity".

14 What is Biodiversity and Why Conserve it	162 The Threat to Pollinators	Taking Molecular Tools to the Jungle
Habitat Fragmentation	172 The Ripple Effect	The Importance of the Top Predator
82 Whisky and Oysters	180 The Downstream Effects of Overfishing	274 The Importance and Fragility of the Amazon Biome
120 The Race to Save the Reefs	208 Microbes Cleaning up our Mess	284 Backyard Biodiversity
132 Choking the Oceans with Plastic	220 The Effects of Dams on River Systems	

Mini-reviews were authored by Joanne Gibson, a Science Writer at New England Biolabs, Inc. Joanne received her Ph.D. in Molecular Biology from the University of Sydney, Faculty of Medicine, Sydney, Australia. She then became a postdoctoral associate at the Massachusetts Institute of Technology, Cambridge, MA. Prior to joining NEB, Joanne served as an Editor for American Journal Experts. When she is not working, Joanne enjoys spending time with her family, traveling & camping, and volunteering for local educational outreach opportunities.





What is Biodiversity and Why Conserve it?

Biodiversity is a broad concept that refers to the "web of life". It describes all of the organisms on Earth, the variety of organisms in an ecosystem and the genetic make-up of those organisms. Biodiversity is responsible for sustaining life as we know it — from producing oxygen, to pollinating our food sources and cleaning our water supplies.

Relationships between organisms have evolved over millions of years into a complex network that functions in harmony with the Earth. All species within an ecosystem have a role to play, and the productivity and stability of any given ecosystem depends on the ability of all of the organisms to work together to maintain balance and ensure each other's survival.

Since humans have inhabited the Earth, there has been a reduction in biodiversity and genetic diversity, yet our health and well-being are intricately tied to its survival. For example, our food sources require a vast range of native pollinators. Diverse ecosystems are responsible for purifying water and recycling nutrients. Our economic well-being relies on the availability of natural resources, such as timber and crude oil. Rainforests help manage the reduction of CO_2 to regulate our climate. Fifty percent of modern pharmaceutical products used in developed countries are derived from plants, animals and microorganisms, and this number increases to 80% if you take into consideration traditional medicines used across the globe. Biodiversity enriches our cultural and recreational experiences: it pleases the senses and gives us the opportunity to exercise, whether it be hiking, bird watching or boating on a lake.

Ecosystems are vulnerable to collapse, and a loss of biodiversity makes them more susceptible to disease and sudden environmental changes. It renders the ecosystem less adaptable.

Humans have caused rapid changes to ecosystems and a decline in biodiversity. Habitat destruction, particularly in biodiverse-rich tropical regions, is the primary cause of an accelerated level of extinction of many species. Upstream of habitat destruction is overpopulation, which has led to an unsustainable level of resource consumption. Downstream of habitat destruction is climate change — to which the burning of fossil fuels consequently warming the atmosphere, the destruction of forests that play a critical role in regulating climate, ocean acidification and changes in plant morphology all contribute. Other factors, such as the accidental and deliberate introduction of non-native species to an ecosystem and overfishing of our lakes and oceans, are also significant threats to biodiversity.

One key challenge we face is to better understand the complexity of all biodiversity and the network of interactions that occur within ecosystems — this knowledge will help us to better evaluate the threats. Some of the most significant disruptions to ecosystems have been the result of a lack of awareness of the complexity within that ecosystem. Our attempts to control an ecosystem for human benefit has, in many cases, been disastrous.

We also need to expand protected areas of land and ocean, with a desire to both preserve ecosystems and benefit from what they have to offer. Precious resources must be carefully utilized, while fully considering the implications.

Biodiversity is declining, and human population is growing exponentially; however, we are now gaining tools and the necessary knowledge to protect our planet for future generations. Scientific and community involvement is continually expanding, which has the potential to lead to policy changes designed to conserve biodiversity.

Restriction Endonucleases



The leader in the discovery & production of restriction enzymes.

Having supplied restriction enzymes to the research community for over 40 years, NEB has earned the reputation of being a leader in enzyme technologies. Working continuously to be worthy of that distinction, NEB strives to develop enzymes of the highest purity and unparalleled quality.

NEB scientists continue to improve our existing portfolio, as well as explore the utility of NEB reagents in new technologies. As a result, NEB scientists continue to publish scientific papers and be awarded grants in this area. With the industry's largest research and development group dedicated to restriction enzymes, we are proud to have been there first: the first to commercialize a recombinant enzyme, the first to introduce a nicking enzyme, and the first to supply a true restriction enzyme master mix. In addition, NEB has a continuing history of innovation by engineering restriction enzymes with altered specificities and improved performance. Through ongoing research in these areas, we are committed to driving the innovations that allow us to offer maximum performance and convenience.

Featured Tools and Resources

Performance/Activity Chart for Restriction Enzymes

Tips for Restriction Enzyme Optimization

Restriction Enzyme
Troubleshooting Guide

302 Time-Saver™ List



Visit NEBrestrictionenzymes.com to find additional online tools, video tech tips and tutorials to help you in setting up restriction enzyme reactions.

Icon Descriptions

The gene encoding this enzyme was cloned at NEB.

This enzyme is purified from a recombinant source.

e This enzyme has been engineered for maximum performance.

Time-Saver qualified enzymes will digest 1 µg of substrate DNA in 5–15 minutes using 1 µl of enzyme under recommended reaction conditions. These enzymes can also be used overnight with no loss of sample. For more information, see pages 302–303.

2*site Indicates that the restriction enzyme requires two or more sites for cleavage.

NEB 1.1 NEB 2.1 NEB 3.1 CutSmart NEB U

Indicates which reaction buffer is supplied with the enzyme for optimal activity. Enzymes with buffer requirements not met by one of the four standard NEBuffers (1.1, 2.1, 3.1 or CutSmart®) are supplied with their own unique NEBuffer (NEB U). NEBuffers are color-coded (NEB 1.1-yellow, NEB 2.1-blue, NEB 3.1-red, CutSmart-green) and supplied as 10X stocks with each enzyme. For more information, consult the Performance Chart on pages 293–298.

Epi This enzyme is EpiMark validated for epigenetics studies.

This enzyme is supplied with a separate tube of S-adenosylmethionine (SAM). To obtain 100% activity, SAM should be added to the 1X reaction mix as indicated. When required, a concentrated stock of SAM is supplied with the enzyme.

dam dcm CpG This restriction enzyme is sensitive to dam, dcm, or CpG methylation. (Note that CpG methylation is applicable to eukaryotic genomic DNA only.) For more information, see pages 334–336.

186° 180° 1 No. Indicates whether or not the enzyme can be heat inactivated. Enzymes are first tested by incubation at 65°C for 20 minutes; any enzyme not inactivated at 65°C is then tested by incubation at 80°C for 20 minutes. If an enzyme can be heat inactivated, the temperature is indicated in the icon.

25° 37° 50° 55° 60° 65° 75°

Indicates the enzyme's optimal incubation temperature.

dil A dil B dil C Indicates which diluent buffer (A, B or C) is recommended for making dilutions of restriction enzymes. For more information see pages 293–298.



Enzymes for Innovation.



```
RR Aatli
                 ₽ BfuAl
                                   R BssSI-v2
                                                     RR Fokl
                                                                         RR Nael
                                                                                           RR Sall
R AbaSi
                 B BgII
                                   R BstAPI
                                                     RR Fsel
                                                                         RR Nari
                                                                                           RR Sall-HF
                 B BgIII
                                                                         RR Ncil
RR Acci
                                   RII BstBl
                                                      RR Fspl
                                                                                            RR Sapi
Acc65I
                                   BstEII
                                                     RR FspEl
                                                                         RR Ncol
                                                                                           ■ Sau3Al
                 R Blol
RR Acil
                 I BmgBl
                                   RR BstEII-HF
                                                     RR Haell
                                                                         RN Ncol-HF
                                                                                           RN Sau961
RR AcII
                 RR Bmrl
                                   R BstNI
                                                     RR Haelli
                                                                         RR Ndel
                                                                                           RR Sbfl
                                                                                           RR Sbfl-HF
RR Acul
                                                                         RR NgoMIV
                 RM Bmtl
                                      BstUI
                                                     ₽ Hgal
                 RR Bmtl-HF

    BstXI

RR Afel
                                                     RR Hhal
                                                                         RR Nhel
                                                                                           RN Scal-HF
RR AfIII
                 R Bpml
                                   R BstYl
                                                                         RR Nhel-HF
                                                                                           RR ScrFI
                                                     ₩ HincII
RR AfIIII
                 Bpu10I
                                   RR BstZ17I-HF
                                                     RR HindIII
                                                                         RR NIalli
                                                                                           RR SexAl
R Agel
                 RR BouEl
                                   RN Bsu361
                                                     RR HindIII-HF
                                                                         RR NIaIV
                                                                                           RR SfaNI
RR Agel-HF
                 RR Bsal
                                   RN Btgl
                                                     R# Hinfl
                                                                         RR NmeAlli
                                                                                           RR SfcI
RR Ahdi
                 RR Bsal-HFv2
                                   RR BtgZI
                                                     RR HinP11
                                                                         R Notl
                                                                                           RR Sfil
RR Alel-v2
                                                                         ™ Notl-HF
                 RR BsaAl
                                   RR BtsI-v2
                                                     RR Hpal
                                                                                           RR Sfol
RR Alui
                 RN BsaBl
                                   RR BtslMutl
                                                                         RR Nrul
                                                                                           RR SgrAl
                                                     RR Hpall
R Alwi
                 RR BsaHl
                                   RR BtsCI
                                                     RR Hphl
                                                                         RR Nrul-HF
                                                                                           RR Smal
RR AlwNi
                 RN BsaJI
                                      Cac8I
                                                     RR Hpy991
                                                                         RR Nsil
                                                                                           RR Smll
RR Apal
                 RN BsaWl
                                   RR Clai
                                                     RR Hpy166II
                                                                         RR Nsil-HF
                                                                                           RR SnaBl
RR ApaLl
                     BsaXI
                                   RR CspCI
                                                     RR Hpy1881
                                                                         RR Nspl
                                                                                           RR Spel
■ ApeKI
                 RR BseRI
                                   RR CviAII
                                                     Hpy188III
                                                                         RR Paci
                                                                                           Spel-HF
ı Apol
                 I BseYI
                                   RR CviKI-1
                                                     RR HpyAV
                                                                         RR PaeR7I
                                                                                           ■ Sphl
RR Apol-HF
                 RR Bsgl
                                   RR CviQ1
                                                     RR HpyCH4III
                                                                         RR Pcil
                                                                                           RR Sphl-HF
RR Ascl
                 RR BsiEl
                                   RR Ddel
                                                     RR HpyCH4IV
                                                                         RR PfIFI
                                                                                           RR Srfl
                                                                         PfIMI
                                                                                           III Sspl
R Asel
                 RR BsiHKAI
                                                     RR HpyCH4V
                                   RR Dpnl
                 ■ BsiWI
                                                                                           Sspl-HF
RR AsiSI
                                   ■ DpnII
                                                     RR Kasi
                                                                         RR Plei
R Aval
                 RR BsiWI-HF
                                   RR Drai
                                                     R Kpnl
                                                                         RR PluTi
                                                                                           RR Stul
RR Avall
                 RR BsII
                                   RR Dralli-HF
                                                     RR Kpnl-HF
                                                                         RR Pmel
                                                                                           R Styl
RR Avrii
                 RR Bsml
                                                     RR LpnPl
                                                                         RR PmII
                                                                                           RR Styl-HF
                                     Drdl
                                   RR Eael
RN Bael
                 RN BsmAl
                                                     RR Mbol
                                                                         RR PpuMI
                                                                                           RR StyD4I
                                                                                            RR Swal
■ BaeGI
                 I BsmBl
                                   RR Eagl
                                                     RR Mboll
                                                                         RR PshAl
                                                                         RR Psil
R BamHI
                 RR BsmFl
                                   RR Eagl-HF
                                                     RR Mfel
                                                                                           R Taq<sup>α</sup>l
RR BamHI-HF
                 RR BsoBI
                                   R# Earl
                                                     RR Mfel-HF
                                                                         R PspGI
                                                                                            RR Tfil
                 RR Bsp1286I
                                   RR Ecil
                                                     RR Mlul
RR Banl
                                                                         RR PspOMI
                                                                                               Tsel
                                                     RR MIul-HF
                                                                         RR PspXI
RN Banli
                 RR BspCNI
                                   RR Eco53kl
                                                                                               Tsp45I
RR Bbsl
                 RR BspDI
                                   RR EcoNI
                                                     RR MluCl
                                                                         R Pstl
                                                                                               TspMI
BbsI-HF
                 ■ BspEl
                                   RR Eco0109I
                                                     RR Miyi
                                                                         RR PstI-HF
                                                                                           RR TspRI
RR Bbvl
                 RR BspHI
                                   EcoP15I
                                                     RR Mmel
                                                                         Pvul
                                                                                           RR Tth1111
RN BbvCl
                 BspMI
                                                     RR MnII
                                                                         RR Pvul-HF
                                                                                           RR Xbal
                                   R EcoRI
RR Bccl
                 ■ BspQI
                                   RR EcoRI-HF
                                                     R# MscI
                                                                         ₽ Pvull
                                                                                           ™ Xcml
BceAl
                    Bsrl
                                   ™ EcoRV
                                                     RR Msel
                                                                         RR Pvull-HF
                                                                                           RR Xhol
Bcgl
                 ₩ BsrBl
                                   RR EcoRV-HF
                                                     RR MsII
                                                                         RR Rsal
                                                                                           R Xmal
                 ■ BsrDI
                                   RR Esp31
                                                     RR Mspl
                                                                         RR RsrII
                                                                                           RR Xmnl
RR BciVI
                                   Fatl
                                                     ₩ MspA1I
R Bell
                 RR RsrFI-v2
                                                                         R Sacl
                                                                                           RR 7rai
RR BcII-HF
                 ■ BsrGI
                                   RR Faul
                                                     RM MspJi
                                                                         RR SacI-HF
RN BcoDI
                 RR BsrGI-HF
                                   RR Fnu4HI
                                                      RR Mwol
                                                                         RR Sacil
```

High-Fidelity Restriction Enzymes

RR Bfal

RR BssHII

```
Agel-HF, EM Apol-HF, EM BamHi-HF, EM Bbsi-HF, EM Bcil-HF, EM Bmti-HF, EM Bsai-HFv2, EM BsiWi-HF, EM BsrGi-HF, EM BstEil-HF, EM BstZ171-HF, EM Draili-HF, EM Eagl-HF, EM EcoRi-HF, EM EcoRv-HF, EM Hindili-HF, EM Kpni-HF, EM Mfel-HF, EM Ncol-HF, EM Nhel-HF, EM Noti-HF, EM Nsil-HF, EM PstI-HF, EM Pvul-HF, EM Pvul-HF, EM Pvul-HF, EM Pvul-HF, EM Pvul-HF, EM Sail-HF, EM Sbfi-HF, EM Sphi-HF, EM Sphi-HF, EM Styl-HF 19, 301
```

Methylation Sensitive Restriction Enzymes for Epigenetics Studies

```
am Abasi, am Dpni, am Dpnii, am Fspei, am Hpali, am LpnPi, am McrBC, am Mspi, am Mspji
264–265
```

Nicking Endonucleases

```
■ Nb.BbvCl, ■ Nb.Bsml, ■ Nb.BsrDl, ■ Nb.BssSl, ■ Nb.Btsl, ■ Nt.Alwi, ■ Nt.BbvCl, ■ Nt.BsmAl, ■ Nt.BspQl, ■ Nt.BstNBl, ■ Nt.CviPll 53-55
```

Homing Endonucleases

```
Im I-Ceul, Im I-Scel, Im PI-Pspl, Im PI-Scel 55-56
```

NEBuffers, Diluents, Gel Loading Dyes, BSA & Recombinant Albumin, Molecular Biology Grade, NEB Tube Opener 56–57

Enzymes in green are all 100% active in CutSmart Buffer (see page 293-298).

One or more of these products are covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at gbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications. Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at https://www.neb.com/support/terms-of-sale. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

Looking to bring convenience to your workflow?

Speed up digestions with Time-Saver™ Qualified Restriction Enzymes


There are 195 NEB restriction enzymes that can digest DNA in 5–15 minutes; many of which are CutSmart or High-Fidelity (HF°) Restriction Enzymes. If you prefer, you can also digest overnight with no unwanted star activity. All of our enzymes are rigorously tested for nuclease contamination. Only NEB can offer enzymes with the power to digest in 5–15 minutes, and the flexibility to withstand overnight digestions with no loss of sample (see page 302–303).

For more information, visit www.neb.com/timesaver



Improve your analysis with our Purple Gel Loading Dye

Our Gel Loading Dye, Purple (6X), which is supplied with most restriction enzymes and all HF enzymes, sharpens bands and eliminates the UV shadow seen with other dyes. This solution contains SDS, which often results in sharper bands, as some restriction enzymes are known to remain bound to DNA following cleavage.



Bring flexibility to your workflow

NEB offers the largest selection of restriction enzymes commercially available. With an evergrowing list to choose from, currently at 286 restriction enzymes – including traditional restriction enzymes, nicking endonucleases, homing endonucleases and methylation-sensitive enzymes for epigenetics studies – there is no need to look anywhere else.



Simplify reaction setup and double digestion with CutSmart® Buffer

Over > 215 enzymes are 100% active in a single buffer, CutSmart Buffer, making it significantly easier to set up double digest reactions. Since CutSmart Buffer includes BSA, there are fewer tubes and pipetting steps to worry about. Additionally, many DNA modifying enzymes are 100% functionally active in CutSmart Buffer, eliminating the need for subsequent purification (see page 299).

For more information, visit www.NEBCutSmart.com


Looking to optimize performance in your reaction?


Choose a High-Fidelity (HF) Restriction Enzyme

NEB High-Fidelity (HF) restriction enzymes have the same specificity as native enzymes, with the added benefits of reduced star activity, rapid digestion (5–15 minutes), and 100% activity in CutSmart Buffer. Enjoy the improved performance of our engineered enzymes at the same price as the native enzymes!

For more information, visit www.neb.com/HF

EcoRI-HF (NEB #R3101) shows no star activity in overnight digests, even when used at higher concentrations. 50 µl rxns were set up using 1 µg of Lambda DNA, the indicated amount of enzyme and the recommended reaction buffer. Rxns were incubated overnight at 37°C. Marker M is the 1 kb DNA Ladder (NEB #N3232).






Benefit from industry-leading quality controls

NEB's reputation as a leader in enzyme technologies stems from the quality and reliability of our restriction enzymes. All of our restriction enzymes undergo stringent quality control testing, ensuring the highest levels of purity and lot-to-lot consistency.

For more information on quality at NEB, visit www.neb.com/quality

Restriction Enzyme Competitor Study: Nuclease Contamination



EcoRI, NotI, and BamHI from multiple suppliers were tested in reactions containing a fluorescent labeled single stranded, double stranded blunt, 3' overhang or 5' overhang containing oligonucleotides. The percent degradation is determined by capillary electrophoresis and peak analysis. The resolution is at the single nucleotide level.



AatII

#R0117L

#R0117S 500 units61 € 2,500 units244 €

5′... G A C G T C ... 3′ 3'... C_T G C A G ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 20,000 units/ml

CutSmart Ri O dil B 37° W CpG NEBuffer 1.1 2.1 3.1 CutSmart

% Activity <10 50* 50

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Note: *May exhibit star activity in this

3′... T T G C A A ... 5′

5′... A A^TC G T T ... 3′

Reaction Conditions: CutSmart

Buffer, 37°C

AclI

#R0598S

#R0598L

CutSmart RR dil B 37° WW CpG

NEBuffer 1.1 2.1 3.1 CutSmart 300 units 69 € 1,500 units276 € % Activity <10 <10 <10 100

Concentration: 5,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

AbaSI

#R0665S 1,000 units113 € See page 264 for more information.

AccI

#R0161S 1,000 units74 € #R0161L 5,000 units296 €

5′... G T[▼]M K A C ... 3′ 3′... C A K M_AT G ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

CutSmart Ri Epi dil C 25° 65

CutSmart RX dilA 37° 60 CpG

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 50 50 10 100

Methylation Sensitivity: Cleavage of mammalian genomic

methylation (see p. 334).

Concentration: 10,000 units/ml

DNA is blocked by overlapping CpG

AcuI

#R0641S 300 units66 € #R0641L 1,500 units264 €

5′... C T G A A G (N)₁₆▼... 3′ 3′... G A C T T C (N) 14... 5′

Reaction Conditions: CutSmart Buffer + SAM, 37°C, Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 50 100

Concentration: 5,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

CutSmart RR SAM dii B 37° 165

Note: Star activity may result from a glycerol concentration of > 5%.

Acc65I

#R0599S #R0599L

5′... G^TG T A C C ... 3′ 3'... C C A T G G ... 5'

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10.000 units/ml

2,000 units69 € NEBuffer 1.1 2.1 3.1 CutSmart 10,000 units276 € % Activity 10 75* 100 25

> Methylation Sensitivity: Blocked by some combinations of overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

Note: *May exhibit star activity in this buffer.

AfeI

#R0652S 200 units 72 € #R0652L 1,000 units 288 €

5'... AGC GCT ... 3' 3'... T C G C G A ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 25 100

Concentration: 10,000 units/ml

CutSmart Rill dil B 37° 65 CpG

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

AciI

#R0551S 200 units69 € #R0551L 1,000 units276 €

5′... C[▼]C G C ... 3′ 3'... G G C₄G ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RR OdiA 37° CpG NEBuffer 1.1 2.1 3.1 CutSmart

Concentration: 10,000 units/ml

% Activity <10 25 100 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

AfIII

#R0520S #R0520L

5′... C[▼]T T A A G ... 3′ 3′... G A A T T_AC ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RX dilA 37° 65

2,000 units66 € NEBuffer 1.1 2.1 3.1 CutSmart 10,000 units264 € % Activity 50 100 10 100

Concentration: 20,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.



























AfIII

1,250 units272 €

NEBuffer 1.1 250 units68 €

5′... A^TCRYGT...3′ 3′... T G Y R C A ... 5′

Reaction Conditions: NEBuffer 3.1. 37°C. Heat inactivation: 80°C for 20 minutes.

RX NEB 3.1 dil B 37° 180°

2.1 3.1 CutSmart % Activity 10 50 100

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

NEW AleI-v2

#R0685S 500 units69 € #R0685L 2,500 units276 €

5'...CACNNNNGTG...3' 3'...GTGNNNNCAC...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 10 10 10 100

CutSmart Rill e dil B 37° CpG

Concentration: 10,000 units/ml

Methylation Sensitivity: Impaired by overlapping CpG methylation (see p. 334).

AgeI

#R0541S

#R0541L

#R0552S 300 units72 € #R0552L 1,500 units 288 €

5′... A C C G G T ... 3′ 3′... T G G C CA ... 5′

Reaction Conditions: NEBuffer 1.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

RX NEB 1.1 dil C 37° 655 CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity **100** 75 25

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Note: Star activity may result from extended digestion.

AluI

#R0137S 1,000 units68 € #R0137L 5,000 units272 €

5′... A G[™]C T ... 3′ 3'... T C_AG A ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

CutSmart Ri G dil B 37° W NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 50 100

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

AgeI-HF®

#R3552S 300 units72 € #R3552L 1,500 units288 €

5′... A C C G G T ... 3′ 3'... T G G C C A ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Rill e dilA 37° 655 CpG

High-Fidelity

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 50 10 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Concentration: 20,000 units/ml

AlwI

#R0513S 500 units70 € #R0513L 2,500 units280 €

5′... G G A T C (N)₄▼... 3′ 3′... C C T A G (N)₅... 5′

Reaction Conditions: CutSmart Buffer, 37°C

Concentration: 10.000 units/ml

CutSmart Ri dilA 37° dil dam

NEBuffer 1.1 2.1 3.1 CutSmar % Activity 50 50 10 100

Methylation Sensitivity: Blocked by dam methylation (see p. 334).

Note: Star activity may result from extended digestion, high enzyme concentration or glycerol concentration of > 5%.

AhdI

#R0584S 1,000 units76 € #R0584L 5,000 units304 €

5'...GACNNNNGTC...3' $3^\prime...\,C\,TG\,NNNN\,N\,C\,AG\,...\,5^\prime$

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

CutSmart RR dilA 37° 655 CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 25 10

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by some combinations of overlapping CpG methylation (see p. 334).

AlwNI

#R0514S #R0514L

5...CAGNNNCTG...3 3′... G T CNNNG A C ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 10,000 units/ml

CutSmart RR diA 37° dim

500 units69 € NEBuffer 1.1 2.1 3.1 CutSmar 2,500 units276 € % Activity 10 100 50 100

> Methylation Sensitivity: Blocked by overlapping dcm methylation (see p. 334).

Note: Conditions of low ionic strength, high enzyme concentration, glycerol concentration > 5%, or pH > 8.0 may result in star activity (see p. 300).

ApaI CutSmart Ri dilA 25° dil dcm CpG

#R0114S 5,000 units72 € #R0114L 25,000 units288 €

5′... G G G C C^{*}C ... 3′ 3′... C_AC C G G G ... 5′

Reaction Conditions: CutSmart Buffer, 25°C. Heat inactivation: 65°C for 20 minutes

Concentration: 50,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 25 <10 100

Activity at 37°C: 100% However, the half-life of Apal at 37°C is only 30 minutes.

Methylation Sensitivity: Blocked by overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

ApoI-HF®

#R3566S 1,000 units72 € #R3566L 5,000 units288 €

5′... R^TA A T T Y ... 3′ 3′... Y T T A A R ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 100 10 100

Concentration: 20,000 units/ml

CutSmart RR e dil B 37° Wh

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

ApaLI

#R0507M

#R0507S 2,500 units66 € #R0507L 12,500 units264 € for high (5X) concentration

12,500 units264 €

5'... GTGCAC...3' $3^\prime...~C~A~C~G~T_{\blacktriangle}G~\dots~5^\prime$

Reaction Conditions: CutSmart

Buffer, 37°C

CutSmart RR dilA 37° Mb CpG NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 100 100 10 100 Concentration: 10,000 and 50.000 units/ml

Methylation Sensitivity:

Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

AscI

#R0558S 500 units68 € #R0558L 2,500 units272 €

5′... G G C G C C ... 3′ $3'\dots$ C C G C G C $_{\mathbf{A}}$ G G \dots 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes

CutSmart RR OdiA 37° CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 10 10 100

Concentration: 10,000 units/ml Methylation Sensitivity: Cleavage of mammalian genomic DNA is

blocked by CpG methylation (see p. 334).

ApeKI

#R0643S 250 units72 € #R0643L 1,250 units288 €

5′... G[▼]C W G C ...3′ 3'... C G W C₄G ... 5'

Reaction Conditions: NEBuffer 3.1,

Concentration: 5,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 50 100 10

Activity at 37°C: 10%

Methylation Sensitivity:

Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

AseI

#R0526S 2,000 units66 € #R0526L 10.000 units264 € for high (5X) concentration

#R0526M 10,000 units264 €

5′... AT T A A T ... 3′ 3′... T A A T_AT A ... 5′

Reaction Conditions: NEBuffer 3.1. 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 50 100

Concentration: 10.000 and 50,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%.

ApoI

#R0566S 1,000 units72 € 5,000 units288 € #R0566L

5′... R^VA A T T Y ... 3′ 3′... Y T T A A R ... 5′

Reaction Conditions: NEBuffer 3.1, 50°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 75 100

Activity at 37°C: 50%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

AsiSI

#R0630S 500 units68 € #R0630L 2,500 units272 €

5'... G C G A T C G C ... 3' 3'... C G C,T A G C G ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 10,000 units/ml

CutSmart Rill dil B 37° 80 CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 25 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Note: Star activity may result from extended digestion.





























AvaI CutSmart 💓 RR 🕜 dii A 37° 🚧 CpG

#R0152S 2,000 units ...65 € #R0152L 10,000 units ...260 € for high (5X) concentration #R0152T 2,000 units ...65 € #R0152M 10,000 units ...260 €

5′... C^TY C G R G ... 3′ 3′... G R G C Y_AC ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer	1.1	2.1	3.1	CutSmart
% Activity	<10	100	25	100

Concentration: 10,000 and 50,000

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

BaeGI

#R0708S

500 units65 €

5′... G K G C M C ... 3′ 3′... C M C G K G ... 5′

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 **3.1** CutSmart % Activity 75 75 100 25

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

AvaII

#R0153S 2,000 units65 € #R0153L 10,000 units260 €

for high (5X) concentration #R0153M 10,000 units 260 €

 $5'...G^{\P}G W C C ...3'$ $3'...C C W G_{\P}G ...5'$

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes. NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 50 75 10 100

CutSmart RR O dil A 37° W dcm CpG

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Blocked by overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

BamHI

#R0136S 10,000 units50 € #R0136L 50,000 units200 € for high (5X) concentration #R0136T 10,000 units50 € #R0136M 50,000 units200 €

5′... G G A T C C ... 3′ 3′... C C T A G G ... 5′

Reaction Conditions: NEBuffer 3.1, 37°C

 NEBuffer
 1.1
 2.1
 3.1
 CutSmart

 % Activity
 75
 100
 100
 100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%.

AvrII

#R0174S 100 units72 € #R0174L 500 units288 €

5′... C^TC T A G G ... 3′ 3′... G G A T C_AC ... 5′

Reaction Conditions: CutSmart

Buffer, 37°C

CutSmart RR dilB 37° Mb

 NEBuffer
 1.1
 2.1
 3.1
 CutSmart

 % Activity
 100
 50
 50
 100

Concentration: 5,000 units/ml

Methylation Sensitivity: Not sensitive to dam. dcm or mammalian

CpG methylation.

BamHI-HF®

#R3136S 10,000 units50 € #R3136L 50,000 units200 € for high (5X) concentration #R3136T 10,000 units50 € #R3136M 50,000 units200 €

 $5'\dots$ $G^{\mathsf{T}}G$ A T C C \dots 3' $3'\dots$ C C T A $G_{\bullet}G$ \dots 5'

Reaction Conditions: CutSmart

Buffer, 37°C

High-Fidelity

CutSmart RR e e dilA 37° Mb

 NEBuffer
 1.1
 2.1
 3.1
 CutSmart

 % Activity
 100
 50
 10
 100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

BaeI

#R0613S 250 units 66 € NEBuffer 1.1 2.1 3.1 Cutsma
#R0613L 1,250 units 264 € % Activity 50 100 50 100

 $5^{'}.\overset{\blacktriangledown}{\underset{10}{\downarrow}}(N) \ A \ C \ (N)_{4} \ G \ T \ A \ Y \ C \ (N)_{12}\overset{\blacktriangledown}{\underset{12}{\downarrow}}.3^{'} \ 3^{'}.\underset{15}{\underset{15}{\downarrow}}(N) \ T \ G \ (N)_{4} \ C \ A \ T \ R \ G \ (N)_{7}\overset{\blacktriangledown}{\underset{1}{\downarrow}}.5^{'}$

Reaction Conditions: CutSmart Buffer + SAM, 25°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

CutSmart Ril O SAM diiA 25° 666 CpG

Activity at 37°C: 20%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

BanI

#R0118S 5,000 units68 €

5′...G^TG Y R C C ...3′ 3′... C C R Y G_AG ...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 20,000 units/ml

CutSmart RR dil A 37° 65 dcm CpG

 NEBuffer
 1.1
 2.1
 3.1
 CutSmart

 % Activity
 10
 25
 <10</td>
 100

Methylation Sensitivity: Blocked by some combinations of overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

500 units264 €

BanII

CutSmart Ril dil A 37° 180°

#R0119S 2,000 units68 €

5′... G R G C Y C ... 3′ 3'... C_AY C G R G ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 10.000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 50

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from extended digestion.

BbvCI

#R0601S

#R0601L

Buffer, 37°C.

100 units66 €



NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 100 50 100

Methylation Sensitivity: Cleavage 5′... C C[▼]T C A G C ... 3′ of mammalian genomic DNA is 3′... G G A G T₄C G ... 5′ impaired by overlapping CpG methylation (see p. 334). Reaction Conditions: CutSmart

> Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

BbsI

#R0539S 300 units68 € #R0539L 1,500 units272 €

5'... G A A G A C $(N)_2^{\blacktriangledown}$... 3' 3'... C T T C T G $(N)_{6_{\blacktriangle}}$... 5'

Reaction Conditions: NEBuffer 2.1, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 25

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Store at -80°C.

#R0704S 1,000 units69 € #R0704L 5,000 units276 €

Concentration: 2,000 units/ml

3'... G G T A G (N)₅... 5'

Concentration: 10,000 units/ml

BccI



5′... C C A T C (N)₄...3′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.



NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 50 10 100

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%.

High-Fidelity

BbsI-HF®

#R3539S 300 units 68 € #R3539L 1,500 units272 €

5'... G A A G A C $(N)_2^{\blacktriangledown}...3'$ 3'... CTTCTG (N)₆...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RR e dil B 37° 65

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 10 <10 100

Concentration: 20,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BceAI

#R0623S 50 units66 € #R0623L 250 units264 €

5'... A C G G C $(N)_{12}^{\blacktriangledown}...3'$ 3'... T G C C G $(N)_{14}^{\blacktriangledown}...5'$

Reaction Conditions: NEBuffer 3.1. 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 2,000 units/ml



NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 100 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

BbvI

#R0173S 300 units70 €

5′... G C A G C (N)₈ ·... 3′ 3′... C G T C G (N)₁₂... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes

CutSmart Rill 2*site dil B 37° 165

% Activity 100 100 25 100 Concentration: 2.000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%.

BcgI



Reaction Conditions: NEBuffer 3.1 + SAM, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 2,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 75* 100 50*

RX 2+site NEB 3.1 SAM dil A 37° 1664 dam CpG

Methylation Sensitivity: Impaired by overlapping dam methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

Note: *May exhibit star activity in this buffer.

NEB 1.1





























BciVI CutSmart Ril O dil C 37° 1864

#R0596S 200 units66 € #R0596L 1,000 units264 €

5′... G T A T C C (N), ♥...3′ 3′... C A T A G G (N)₅ ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 25 <10 100

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

#R0568L 2.500 units304 € 5′... C[▼]T A G ... 3′ 3′... G A T_AC ... 5′

BfaI

#R0568S

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

500 units76 €

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 10 <10

CutSmart Ri dil B 37° 180°

Concentration: 10.000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from extended digestion.

BclI

#R0160S 3,000 units 62 € 15,000 units248 € #R0160L

5′... T G A T C A ... 3′ 3′... A C T A G₄T ... 5′

Reaction Conditions: NEBuffer 3.1,

50°C

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 100 75

Activity at 37°C: 50%

Methylation Sensitivity: Blocked by dam methylation (see p. 334).

BfuAI

#R0701S 250 units66 € 1,250 units264 € #R0701L

5′... A C C T G C (N)₄ ... 3′ 3′... T G G A C G (N)₈...5′

Reaction Conditions: NEBuffer 3.1, 50°C. Heat inactivation: 65°C for

20 minutes.

Concentration: 5,000 units/ml

% Activity <10 25 100 10

NEBuffer 1.1 2.1 3.1 CutSmart

Activity at 37°C: 50%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by overlapping CpG methylation (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.

NEW

BclI-HF®

#R3160S 3,000 units62 € #R3160L 15,000 units248 €

5'... TGATCA...3' 3′... A C T A G_AT ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

High-Fidelity

CutSmart Rill e dil B 37° 655 dam NEBuffer 1.1 2.1 3.1 CutSmart

Concentration: 20,000 units/ml

% Activity 100 100 10 100

Methylation Sensitivity: Blocked by dam methylation (see p. 334).

BfuCI

BfuCl has been replaced by its isoschizomer Sau3Al.

BglI

#R0143S 2,000 units60 € #R0143L 10,000 units240 €

5'...GCCNNNNNGGC...3' 3′... CGGNNNNNCCG...5′

Reaction Conditions: NEBuffer 3.1. 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart

Concentration: 10,000 units/ml

% Activity 10 25 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

BcoDI

CutSmart RR O dil B 37° W CpG 1,000 units68 € #R0542S

5′... G T C T C (N), [▼]... 3′ 3′... C A G A G (N)₅...5′

Reaction Conditions: CutSmart

Buffer, 37°C

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 75 75 100

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by some combinations of overlapping CpG methylation (see p. 334).

BglII

#R0144S 2.000 units61 € #R0144L 10,000 units244 € for high (5X) concentration #R0144M 10,000 units244 €

5'... A G A T C T ... 3' 3′... T C T A GA... 5′

Reaction Conditions: NEBuffer 3.1, 37°C.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 10 100

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BlpI CutSmart RX dilA 37° Vib

#R0585S 500 units76 € 2,500 units304 € #R0585L

5′... G C^TT N A G C ... 3′ 3′... C G A N T C G ... 5′

Reaction Conditions: CutSmart Buffer, 37°C

NEBuffer 1.1 2.1 3.1 CutSmart

Concentration: 10,000 units/ml

% Activity 50 100 10 100

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BmtI-HF®

#R3658S 300 units 68 € #R3658L 1,500 units 272 €

5′... G C T A G^TC ... 3′ 3′... C_AG A T C G ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 10 100

CutSmart RR C dil B 37° 65

High-Fidelity

Concentration: 20,000 units/ml.

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BmgBI

#R0628S 500 units70 € #R0628L 2,500 units280 €

5′... C A C G T C ... 3′ 3′... G T G C A G ... 5′

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 10 100 10

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.

BpmI

#R0565S #R0565L

5′... C T G G A G (N)₁₆ ... 3′ 3′... G A C C T C (N)₁₄... 5′

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 2,000 units/ml.

RX 2+site NEB 3.1 dil B 37° 1654

100 units68 € NEBuffer 1.1 2.1 3.1 CutSmart 500 units272 € % Activity 75 100 100 100

> Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from extended digestion.

BmrI

#R0600S 100 units68 €

5′... A C T G G G (N)₅ ... 3′ 3′... T G A C C C (N)₄ ... 5′

Reaction Conditions: NEBuffer 2.1. 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5.000 units/ml

RR NEB 2.1 dil B 37° 166

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 75 100 75

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: *May exhibit star activity in this buffer.

Bpu10I

#R0649S 200 units 69 € #R0649L 1.000 units 276 €

5'... C C T N A G C ... 3' 3'... GGANT CG...5'

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 5,000 units/ml

RN MEB 3.1 dil B 37° 186

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 25 100 25

Methylation Sensitivity: Not sensitive to dam. dcm or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%.

BmtI

#R0658S 300 units 68 €

5′... G C T A G[▼]C ... 3′ 3'... C G A T C G ... 5

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 65°C for 20 minutes

Concentration: 10,000 units/ml

RX NEB 3.1 dil B 37° 655

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 100 100

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from extended digestion.

BpuEI

#R0633S 500 units69 €

5′... C T T G A G (N)₁₆ ... 3′ 3′... G A A C T C (N)₁₄ ... 5′

Reaction Conditions: CutSmart Buffer + SAM, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

CutSmart RI SAM dil B 37° 165

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50* 100 50* 100

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: *May exhibit star activity in this buffer.





























BsaI CutSmart Rill dil B 37° 655 dcm CpG

#R0535S 1,000 units69 € #R0535L 5,000 units276 €

5′... G G T C T C (N), ▼...3′ 3′... C C A G A G (N)₅... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for

20 minutes.

Concentration: 10.000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 75 75 100 100

Methylation Sensitivity: Impaired by some combinations of overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.

High-Fidelity

BsaHI

CutSmart Ril O dilC 37° 66 dem CpG

#R0556S 2,000 units66 €

5′... G R C G Y C ... 3′ 3′... C Y G C R G ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 100 100

Concentration: 10,000 units/ml

Methylation Sensitivity: Blocked by some combinations of overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

BsaI-HF®v2

#R3733S 1,000 units69 € #R3733L 5,000 units276 €

5′... G G T C T C (N), ▼... 3′ $3'\dots$ C C A G A G (N)₅ ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 100 100

CutSmart Rill e dil B 37° til dcm CpG

Concentration: 20,000 units/ml

Methylation Sensitivity: Impaired by some combinations of overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

BsaJI

#R0536S 1,000 units69 € #R0536L 5,000 units276 €

5′... C^{*}C N N G G ... 3′ 3′... G G N N C C ... 5′

Reaction Conditions: CutSmart Buffer, 60°C. Heat inactivation: 80°C for 20 minutes.

CutSmart Rik dil A 60° km²

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 100

Concentration: 10,000 units/ml

Activity at 37°C: 20%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BsaAI

#R0531S 500 units69 € #R0531L 2,500 units276 €

5′... Y A C[▼]G T R ... 3′ 3'... RTG_CAY ... 5'

Reaction Conditions: CutSmart

Buffer, 37°C

CutSmart RR dilC 37° Wb CpG NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 100 100 100 100 Concentration: 5.000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

BsaWI

250 units74 € #R0567S

5′... W C C G G W ... 3′ 3′... W G G C C₄W ... 5′

Reaction Conditions: CutSmart Buffer, 60°C. Heat inactivation: 80°C for 20 minutes.

CutSmart RR 60° km

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 100 50

Concentration: 10,000 units/ml

Activity at 37°C: 20%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BsaBI

#R0537S 2,000 units 69 € #R0537L 10,000 units276 €

5...GATNNNATC...3 3′... C T A N N N N T A G ... 5′

Reaction Conditions: CutSmart Buffer, 60°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 10,000 units/ml

Activity at 37°C: 20%

CutSmart Rill dil B 60° 60° dam CpG NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 50 100 75

Methylation Sensitivity: Blocked by overlapping dam methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

Note: Star activity may result from extended digestion.

BsaXI

#R0609S 100 units69 € #R0609L 500 units276 €

 $\begin{array}{l} 5^{\prime}...\overset{\blacktriangledown}{{}_{g}}(N) \ A \ C \ (N)_{_{5}} \ C \ T \ C \ C \ (N)_{_{10}}\overset{\blacktriangledown}{{}_{\cdots}} \ldots 3 \\ 3^{\prime}...\underset{12}{\overset{\longleftarrow}{}_{12}}(N) \ T \ G \ (N)_{_{5}} \ G \ A \ G \ G \ (N)_{_{7}} \end{array}$

Reaction Conditions: CutSmart Buffer, 37°C

Concentration: 2,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50* 100* 10

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

CutSmart dilC 37° \dil

Note: *May exhibit star activity in this

buffer.

BseRI

#R0581S 200 units74 € #R0581L 1,000 units296 €

5′...G A G G A G (N)₁₀...3′ 3′... C T C C T C (N)₈... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes

CutSmart RR dilA 37° W

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 75

Concentration: 5,000 units/ml.

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BsiHKAI

#R0570S 1,000 units67 €

5′... G W G C W C ... 3′ 3′...C_w C G W G ... 5′

Reaction Conditions: CutSmart

Buffer, 65°C.

Concentration: 10.000 and 50,000 units/ml.

Methylation Sensitivity: Not

Activity at 37°C: 5%

sensitive to dam, dcm or mammalian CpG methylation.

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 25 100 100 100

CutSmart RR dilA 65° W6

BseYI

#R0635S 100 units69 € #R0635L 500 units276 €

5'...C"CCAGC...3' 3'...GGGTC₄G...5'

Reaction Conditions: NEBuffer 3 1, 37°C. Heat inactivation: 80°C for 20 minutes.

RX NEB 3.1 dil B 37° VM CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 50 100

Concentration: 5,000 units/ml.

Methylation Sensitivity:

Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

BsiWI

#R0553S 300 units69 € #R0553L 1,500 units276 €

5'... C'G T A C G ... 3' 3'... G C A T G₁C ... 5'

Reaction Conditions: NEBuffer 3.1, 55°C. Heat inactivation: 65°C for

20 minutes.

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 25 50* 100

Activity at 37°C: 50%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Note: *May exhibit star activity in this buffer.

High-Fidelity

BsgI

50 units69 € #R0559S #R0559L 250 units276 €

 $5'\dots GTGCAG(N)_{16}^{\blacktriangledown}\dots 3'$ 3′... C A C G T C (N)₁₄... 5′

Reaction Conditions: CutSmart Buffer + SAM, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Ril 2+site SAM dil B 37° 1654 NEBuffer 1.1 2.1 3.1 CutSmart

Concentration: 5,000 units/ml

% Activity 25 50 25

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BsiWI-HF®

#R3553S 300 units69 € #R3553L 1,500 units276 €

5′... C[▼]G T A C G ... 3′ 3′... G C A T G_AC ... 5′

Reaction Conditions: CutSmart

Buffer, 37°C.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 10 100

Concentration: 20,000 units/ml

CutSmart RR C G dil B 37° Mb CpG

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

BsiEI

#R0554S 1,000 units70 €

5′...CGRY CG...3′ 3′... G C_AY R G C ... 5′

Reaction Conditions: CutSmart Buffer, 60°C.

Concentration: 10,000 units/ml

CutSmart Ril O dil A 60° Wb CpG NEBuffer 1.1 2.1 3.1 CutSmart

Activity at 37°C: 30%

% Activity 25 50 <10 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

BslI

#R0555S 1.000 units67 € #R0555L 5,000 units268 €

5...CCNNNNNNNGG...3 3′... G G N N N N N N C C ... 5′

Reaction Conditions: CutSmart Buffer, 55°C

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 75 100 100

CutSmart RR dilA 55° th dcm CpG

Activity at 37°C: 30%

Methylation Sensitivity: Blocked by some combinations of overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).





























CutSmart RR dilA 65° 866 **BsmI**

#R0134S 500 units70 € 2,500 units280 € #R0134L

5′...G A A T G C N ... 3′ 3′... C T T A C₄G N ... 5′

Reaction Conditions: CutSmart Buffer, 65°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 <10 100

Concentration: 10,000 units/ml

Activity at 37°C: 20%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BsoBI

#R0586S 10,000 units67 €

5′...C[▼]Y C G R G ...3′ 3′...GRGCY_C...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 100

Concentration: 10,000 units/ml

CutSmart RR 6 dil A 37° 866

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BsmAI

#R0529S 1,000 units 69 € #R0529L 5,000 units276 €

5′...G T C T C (N), ▼...3′ 3′... C A G A G (N)₅...5′

Reaction Conditions: CutSmart

Buffer, 55°C

Concentration: 5,000 units/ml

CutSmart RR Odil B 55° WW CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 100 100

Activity at 37°C: 50%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

Bsp1286I

#R0120S 500 units66 €

5′... G D G C H^TC ... 3′ 3′... C_AH C G D G ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

CutSmart Ril G dilA 37° Y654 NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 25 25 25

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

BsmBI

#R0580S 200 units74 € #R0580L 1,000 units296 €

5′... C G T C T C (N), √... 3′ $3' \dots G C A G A G (N)_{5_{A}} \dots 5'$

Reaction Conditions: NEBuffer 3.1. 55°C. Heat inactivation: 80°C for 20 minutes

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 50* 100

Activity at 37°C: 20%

Methylation Sensitivity:

Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Note: *May exhibit star activity in this buffer.

BspCNI

#R0624S 100 units67 €

5′... C T C A G (N)₁₀ ... 3′ 3′... G A G T C (N)₈ ... 5′ and

5′... C T C A G (N)₉ ... 3′ 3′...GAGTC(N)₇...5′

Note: The cleavage site of BspCNI varies. Two equally represented species of fragments are produced from BspCNI cleavage.

CutSmart RR SAM dilA 25° 866

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 75 10

Reaction Conditions: CutSmart Buffer + SAM, 25°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 2.000 units/ml

Activity at 37°C: 75%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BsmFI

#R0572S 100 units\$72 € #R0572L 500 units ... \$288 €

5′... G G G A C (N)₁₀ ▼... 3′ 3′... C C C T G (N)₁₄... 5′

Reaction Conditions: CutSmart Buffer, 65°C. Heat inactivation: 80°C for 20 minutes

Concentration: 2,000 units/ml Activity at 37°C: 50%

CutSmart Ri dil A 65° dcm CpG

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 25 50 50 100

Methylation Sensitivity: Blocked by overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG

Note: Star activity may result from extended digestion, high enzyme concentration or glycerol concentration of > 5%

methylation (see p. 334).

BspDI

#R0557S 2,000 units66 €

5′... A T[▼]C G A T ... 3′ 3′... T A G C T A ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 75 50 100

CutSmart Ri dii A 37° dam CpG

Concentration: 10,000 units/ml Methylation Sensitivity: Blocked

by overlapping dam methylation. Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

BspEI

#R0540S 1,000 units69 € #R0540L 5,000 units276 €

5'... TCCGGA...3' 3′... A G G C C_AT ... 5′

Reaction Conditions: NEBuffer 3.1. 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 10 100 <10

Concentration: 10,000 units/ml

Methylation Sensitivity: Blocked by overlapping dam methylation. Cleavage of mammalian genomic DNA is impaired by CpG methylation (see p. 334).

BsrI

#R0527S 1,000 units66 € #R0527L 5,000 units264 €

5′... A C T G G N[▼]...3′ 3′... T G A C₄C N ... 5′

Reaction Conditions: NEBuffer 3.1, 65°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart

● NEB3.1 dil B 65° 📆

Concentration: 10.000 units/ml

Activity at 37°C: 20%

% Activity <10 50 100

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BspHI

#R0517S 500 units68 € #R0517L 2,500 units272 €

5'... TCATGA...3' 3′... A G T A C_AT ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

CutSmart Ri O dil A 37° W dam

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 50 25

Concentration: 10,000 units/ml

Methylation Sensitivity: Impaired by overlapping dam methylation (see p. 334).

BsrBI

#R0102S 1,000 units68 € #R0102L 5,000 units272 €

5′...CCG^TCTC...3′ 3′...GG C₄G A G ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 100 100

CutSmart RR dilA 37° 66 CpG

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

BspMI

#R0502S 100 units69 €

5′... A C C T G C (N)₄ ▼... 3′ 3′... T G G A C G (N)₈... 5′

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 50* 100

Concentration: 2,000 units/ml

Rii 2+site NEB 3.1 dil B 37° (65)

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian

Note: *May exhibit star activity in this buffer.

CpG methylation.

BsrDI

#R0574S 200 units66 € #R0574L 1,000 units264 €

5'... G C A A T G N N ... 3' 3′... CGTTACNN...5′

Reaction Conditions: NEBuffer 2.1. 65°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 5.000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 100 75

Activity at 37°C: 30%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%.

BspQI

#R0712S 500 units72 € #R0712L 2,500 units288 €

5′...GCTCTTC(N), ▼...3′ 3′... C G A G A A G (N), 5′

Reaction Conditions: NEBuffer 3.1, 50°C. Heat inactivation: 80°C for 20 minutes

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100* 100* 100 100*

Activity at 37°C: 10%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: *May exhibit star activity in this buffer.

5′... R[▼]C C G G Y ... 3′ 3'... Y G G C C₄R ... 5'

BsrFI-v2

#R0682S

Reaction Conditions: CutSmart

Buffer 37°C

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 25 0 100

CutSmart RR e dilC 37° W CpG

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see

p. 334).



































1,000 units66 €

BsrGI

#R0575S 1,000 units68 € #R0575L 5,000 units272 €

5′... T[▼]G T A C A ...3′ 3'... A C A T G_T ... 5'

Reaction Conditions: NEBuffer 2.1, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart

25

High-Fidelity

Concentration: 10,000 units/ml

% Activity 25 100 100

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

> 20 minutes. Concentration: 5,000 units/ml

Reaction Conditions: CutSmart

Buffer, 60°C. Heat inactivation: 80°C for

200 units 96 € NEBuffer 1.1 2.1 3.1 CutSmart 1,000 units 384 € % Activity 50 100 25

5′... G C A N N N N T G C ... 3′ 3′... C G T N N N N N A C G ... 5′ Activity at 37°C: 10%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

CutSmart Ri dilA 60° 60° CpG

BsrGI-HF®

1,000 units #R3575S ..68€ #R3575L 5,000 units272 €

5′... T[▼]G T A C A ... 3′ 3'... A C A T G T ... 5'

Reaction Conditions: CutSmart Buffer, 37°C, Heat inactivation: 80°C for 20 minutes.

CutSmart RR C dilA 37° km NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 10 100 100 100

Concentration: 20.000 units/ml.

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BstBI

#R0519S 2,500 units62 €

5'... T T C G A A ... 3' 3′... A A G C_AT T ... 5′

Reaction Conditions: CutSmart

Buffer, 65°C

BstAPI

#R0654S

#R0654L

#R0519L 12,500 units248 €

Concentration: 20,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 75 100 10

CutSmart RR dilA 65° WW CpG

Activity at 37°C: 10%

Methylation Sensitivity:

Cleavage of mammalian genomic DNA is blocked by CpG methylation

(see p. 334).

BssHII

#R0199S 500 units 68 € #R0199L 2,500 units272 € for high (5X) concentration #R0199M 2,500 units272 €

5′...GCGCGC...3′ 3′...CGCGC**,**G...5′

Reaction Conditions: CutSmart Buffer, 50°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 100 100

CutSmart RR O dil B 50° 655 CpG

Concentration: 5,000 and 25,000 units/ml

Activity at 37°C: 75%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

BstEII

#R0162S #R0162L for high (5X) concentration

#R0162M 10,000 units248 €

5'...G"GTNACC...3'

Reaction Conditions: NEBuffer 3.1,

Concentration: 10,000 and 50.000 units/ml

2,000 units62 € 10,000 units248 €

3'...CCANTG.G...5'

60°C

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 75 100

Activity at 37°C: 50%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%.

BssSI-v2

#R0680S 200 units74 € #R0680L 1,000 units296 €

5'... C'A C G A G ... 3' 3′...G T G C T C ... 5′

Reaction Conditions: CutSmart Buffer, 37°C

CutSmart RR e dil B 37° VIII

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 25 <10 100

Concentration: 10.000 units/ml Methylation Sensitivity: Not sensitive to dam, dcm or mammalian

CpG methylation.

BstEII-HF®

#R3162S 2,000 units62 € #R3162L 10,000 units248 €

for high (5X) concentration

#R3162M 10,000 units248 €

5'...GGTNACC...3' 3′...CCANTG.G...5′

Reaction Conditions: CutSmart

Buffer, 37°C

High-Fidelity

CutSmart RR C dilA 37° Mb

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 10 <10 100

Concentration: 20.000 and 100,000 units/ml.

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian

CpG methylation.

BstNI



#R0168S 3,000 units61 € #R0168L 15,000 units244 €

5'... C C WGG ... 3' 3′... G G W C C ... 5′

Reaction Conditions: NEBuffer 3.1,

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 100 100

Activity at 37°C: 30%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BstZ17I-HF®

CutSmart RR C O dil A 37° Mb CpG



% Activity 100 100 10 100

#R3594S 1,000 units70 € #R3594L 5,000 units280 €

5′...GTA*TAC...3′ 3'... C A TAT G ... 5'

Reaction Conditions: CutSmart Buffer, 37°C.

Concentration: 20,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

BstUI

#R0518S 1,000 units63 € #R0518L 5,000 units252 €

5′... C G[▼]C G ... 3′ 3′...G C₄G C ... 5′

Reaction Conditions: CutSmart Buffer, 60°C

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 50 100 25 Activity at 37°C: 20%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Bsu36I

#R0524S 1,000 units66 € #R0524L 5,000 units264 €

5'... C C"T N A G G ... 3' 3′... G G A N T C C ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

CutSmart Ril O dil C 37° 1804

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 100 100

Concentration: 10,000 units/ml Methylation Sensitivity: Not

sensitive to dam, dcm or mammalian CpG methylation.

BstXI

#R0113S 1,000 units 68 € #R0113L 5,000 units272 €

5′... C C A N N N N N N T G G ... 3′ 3′... G G T N N N N N N A C C ... 5′

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 50 100

Concentration: 10,000 units/ml

Methylation Sensitivity: Blocked by some combinations of overlapping dcm methylation (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.

BtgI

#R0608S 1,000 units66 €

5′...C^TCRYGG...3 3'... G G Y R C₄C ...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

CutSmart RR 6 dil B 37° W

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 100 100

CpG methylation.

Concentration: 10.000 units/ml Methylation Sensitivity: Not sensitive to dam, dcm or mammalian

BstYI

#R0523S 2,000 units66 € #R0523L 10,000 units264 €

5'...R GATCY...3' 3′... Y C T A G R ... 5′

Reaction Conditions: NEBuffer 2.1, 60°C

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 75

Activity at 37°C: 30%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

BtgZI

#R0703S 100 units69 € #R0703L 500 units276 €

5′...G C G A T G (N)₁₀ ... 3′ 3′...CGCTAC(N)₁₄...5′

Reaction Conditions: CutSmart Buffer, 60°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 5,000 units/ml

CutSmart Ri dil A 60° 80 CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 25 <10 100

Activity at 37°C: 75%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by CpG methylation (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.





























BtsI-v2 Cutsmart 💥 Ril e o dnA 55° Wb

#R0667S 500 units71 € #R0667L 2,500 units284 €

5'...GCAGTGNN...3' 3'...CGTCACNN...5'

Reaction Conditions: CutSmart

Buffer, 55°C.

Concentration: 10,000 units/ml

Activity at 37°C: 75%


NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 25 100

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

ClaI #R0197S

#R0197L



5′...AT CGAT...3′ 3′...TAGC TA...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.
 NEBuffer
 1.1
 2.1
 3.1
 CutSmart

 % Activity
 10
 50
 50
 100

CutSmart RX dilA 37° dil dam CpG

Concentration: 10,000 units/ml

Methylation Sensitivity: Blocked

by overlapping *dam* methylation. Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

BtsIMutI

#R0664S 100 units104 €

5′... C AGTG N N ... 3′ 3′... G T C A C N N ... 5′

Reaction Conditions: CutSmart Buffer, 55°C. Heat inactivation: 80°C for 20 minutes. CutSmart Ril e dilA 55° 💖

 NEBuffer
 1.1
 2.1
 3.1
 CutSmart

 % Activity
 100
 50
 10
 100

Concentration: 1.000 units/ml

Activity at 37°C: 50%

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

CspCI

#R0645S 500 units70 €

 $5'..._{10\cdot11}^{\mathbf{Y}}$ (N) C A A (N)₅ G T G G (N)_{12·13}...3' $3'..._{12\cdot13}$ (N) G T T (N)₅ C A C C (N)_{10·11}...5'

Note: The cleavage point may shift one base pair depending on the DNA sequence context before and after the recognition site. For a given sequence, one site will predominate. For details, see www.neb.com.

 NEBuffer
 1.1
 2.1
 3.1
 CutSmart

 % Activity
 10
 100
 10
 100

CutSmart R 2*site SAM dil A 37° 65

Reaction Conditions: CutSmart Buffer + SAM, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

Methylation Sensitivity: Not
sensitive to dam, dcm or mammalian
CpG methylation.

BtsCI

#R0647S 2,000 units72 €

5′...G G AT G N N^V...3′ 3′...C C T A C₁N N ...5′

Reaction Conditions: CutSmart Buffer, 50°C. Heat inactivation: 80°C for 20 minutes.

CutSmart RX dil B 50° 1864

 NEBuffer
 1.1
 2.1
 3.1
 CutSmart

 % Activity
 10
 100
 25
 100

Concentration: 20,000 units/ml

Activity at 37°C: 50%

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

CviAII

#R0640S 200 units66 € #R0640L 1,000 units264 €

5′...C*ATG...3′ 3′...GTA₄C...5′

Reaction Conditions: CutSmart Buffer, 25°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RR O dii C 25° 654

 NEBuffer
 1.1
 2.1
 3.1
 CutSmart

 % Activity
 50
 50
 10
 100

Concentration: 10,000 units/ml

Activity at 37°C: 20%

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Cac8I

#R0579S 100 units74 € #R0579L 500 units296 €

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes. CutSmart dil B 37° S CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 75 100 100

Concentration: 5,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

CviKI-1

#R0710S 250 units72 €

5′...RG^{*}CY...3′ 3′...YC₄GR...5′

Reaction Conditions: CutSmart

Buffer, 37°C

Concentration: 5,000 units/ml

CutSmart RR dilA 37° M6 NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 25 100 100 100

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

CviQI

#R0639S 2.000 units72 € #R0639L 10,000 units288 €

5'... GTAC ... 3' 3'...CAT_G...5'

Reaction Conditions: NEBuffer 3.1,

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 75 100* 100 75* Activity at 37°C: 10%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: *May exhibit star activity in

this buffer.

DraI

2,000 units66 € #R0129S #R0129L 10,000 units264 €

5′... T T T T A A A ... 3′ 3′... A A A A T T T ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for

20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 75 75 50 100

CutSmart RR dilA 37° Kith

Concentration: 20,000 units/ml.

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian

CpG methylation.

DdeI

#R0175S 1,000 units68 € #R0175L 5,000 units272 €

5′...CTNAG...3′ 3′...G A N T_AC ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes

CutSmart RR 6 dil B 37° 65

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 75 100 100 100

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

DraIII-HF®

#R3510S 1,000 units69 € #R3510L 5,000 units276 €

5'... C A C N N N G T G ... 3' 3'... G T G N N N C A C ... 5'

Reaction Conditions: CutSmart

Buffer, 37°C

CutSmart RR e dil B 37° th CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 50 10 100

High-Fidelity

Concentration: 20.000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by some combinations of overlapping CpG methylation (see p. 334).

DpnI

#R0176S 1,000 units65 € 5,000 units260 € #R0176L

CH, 5′... G A T C ... 3′ 3'... C T_A G ... 5' CH,

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 75 100

CutSmart Ri Epi 🔮 dii B 37° 🕍 CpG

Concentration: 20,000 units/ml

Methylation Sensitivity: Dpnl cleaves only when its recognition site is methylated. DNA purified from a dam+ strain will be a substrate for Dpnl. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

DrdI

#R0530S 300 units69 € 1,500 units276 € #R0530L

5'...GACNNNNNNGTC...3' 3'...CTGNNNNNNCAG...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

CutSmart Odil A 37° CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 50 10 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.

DpnII

#R0543S 1,000 units72 € #R0543L 5,000 units288 € for high (5X) concentration #R0543T 1,000 units72 € #R0543M 5,000 units288 €

5′...*G AT C ... 3′ 3′...CTAG...5′

Reaction Conditions: NEBuffer DpnII, 37°C. Heat inactivation: 65°C for 20 minutes.

RX Epi 🚱 NEBU dil B 37° 🙀 dam

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 25 100*

Concentration: 10,000 and 50.000 units/ml

Methylation Sensitivity: Blocked by dam methylation (see p. 334).

Note: *May exhibit star activity in this buffer.

EaeI

#R0508S 200 units66 € #R0508L 1,000 units264 €

5′...Y[™]G G C C R ... 3′ 3'... R C C G G Y ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Rill dil A 37° 65 dcm CpG NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 10 50 <10 100 Concentration: 5,000 units/ml

Methylation Sensitivity: Blocked by overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

NEB 1.1































EagI

#R0505S 500 units 66 € 2,500 units264 € #R05051 for high (5X) concentration #R0505M 2,500 units264 €

5′... C[▼]G G C C G ... 3′ 3′...G C C G G C ... 5′

Reaction Conditions: NEBuffer 3.1. 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 25 100

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Eco53kI

#R0116S 1,000 units60 €

5′...GAG[▼]CTC...3′ 3'...CTC_GAG...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 <10

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

CutSmart RR dilA 37° CpG

Note: Star activity may result from a glycerol concentration of > 5%.

EagI-HF® CutSmart RR e dil B 37° 655 CpG

#R3505S 500 units66 € #R3505L 2,500 units264 €

for high (5X) concentration #R3505M 2,500 units264 €

5′...C^{*}G G C C G ...3′ 3'... G C C G G C ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes

NEBuffer 1.1 2.1 3.1 CutS % Activity 25 100 100 100

High-Fidelity

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

EcoNI

#R0521S 1,000 units69 € #R0521L 5,000 units276 €

5′... C C T N N N N A G G ... 3′ 3′... G G A N N N N T C C ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes

CutSmart Ril O dil A 37° 165 NEBuffer 1.1 2.1 3.1 CutSmart

Concentration: 10,000 units/ml

% Activity 50 100 75

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

EarI

500 units66 € #R0528S 2,500 units264 € #R0528L

5′... C T C T T C (N), ▼... 3′ 3′...G A G A A G (N)₄....5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 10 <10 100

Concentration: 20,000 units/ml

CutSmart Ril O dil B 37° K CpG

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by overlapping CpG methylation (see p. 334).

EcoO109I

#R0503S 2,000 units66 €

5'...R G N C C Y ... 3' 3′... Y C C N G₄G R ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 20,000 units/ml

CutSmart RR dilA 37° dim NEBuffer 1.1 2.1 3.1 CutSmart

Methylation Sensitivity: Blocked by overlapping dcm methylation (see

p. 334).

% Activity 50 100 50

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

EciI

100 units66 € #R0590S #R0590L 500 units264 €

5′...GGCGGA(N),, ▼...3′ 3'...CCGCCT(N), ...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes

Concentration: 2,000 units/ml

CutSmart Rik dil A 37° CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 50 50

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see

Note: Star activity may result from extended digestion.

EcoP15I

#R0646S 500 units69 € #R0646L 2,500 units276 €

5′... C A G C A G (N)₂₅... 3′ 3′...GTCGTC(N)₂₇...5′

Reaction Conditions: NEBuffer 3.1 + ATP, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 75 100 100 100

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian

CpG methylation.

EcoRI

#R0101S	10,000 units	50€
#R0101L	50,000 units	200€
for high (5X) cond	entration	

#R0101T 10,000 units50 € #R0101M 50,000 units200 €

5′...G*A A T T C ...3′ 3′... C T T A A G ... 5′

Reaction Conditions: NEBuffer EcoRI, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart

Concentration: 20,000 and 100.000 units/ml

% Activity 25 100* 50

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combination of overlapping CpG methylation (see p. 334).

Note: *May exhibit star activity in this buffer.

Esp3I



 $5^{\prime}...$ CGTCTC $(N)_1$ $^{\blacktriangledown}...$ 3 $^{\prime}$ 3′... G C A G A G (N)₅ ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 100 100 <10 100

CutSmart RR dil B 37° 65° CpG

Concentration: 10,000 units/ml Methylation Sensitivity: Cleavage

of mammalian genomic DNA is

blocked by CpG methylation (see p. 334).

High-Fidelity EcoRI-HF®

10,000 units50 €
50,000 units200 €
ntration
10,000 units50 €
50,000 units200 €

5′...G*A A T T C ... 3′ 3′... C T T A A_AG ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 100 <10 100

CutSmart Rill e dill 37° 655 CpG

Concentration: 20,000 and 100,00 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

FatI

#R0650S 50 units 94 € #R0650L 250 units 376 €

5'... CATG... 3' 3′...GTAC₄...5′

Reaction Conditions: NEBuffer 2.1, 55°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 100 50

RX NEB 2.1 dil A 55° 866

Concentration: 2,000 units/ml

Activity at 37°C: 20%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

EcoRV

#R0195S 4,000 units\$61 € #R0195L 20,000 units244 € for high (5X) concentration #R0195T 4,000 units61 € #R0195M 20,000 units244 €

5′... G A T[™]A T C ... 3′ 3′... C T A T A G ... 5′

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 80°C for 20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 50 100

Concentration: 20,000 and 100,000 units/ml.

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by some combinations of overlapping CpG methylation (see p. 334).

Faul

#R0651S 200 units 72 € #R0651L 1,000 units 288 €

5′... C C C G C (N)₄ ... 3′ 3′...GGGCG(N)₆...5′

Reaction Conditions: CutSmart Buffer, 55°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

CutSmart Rill dil A 55° CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 50 10 100

Activity at 37°C: 20%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation

(see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%

High-Fidelity

EcoRV-HF®

#R3195S 4,000 units61 € #R3195L 20,000 units244 € for high (5X) concentration 4,000 units61 € #R3195T #R3195M 20,000 units244 €

5'... G A T*A T C ... 3' 3′... C T A T A G ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RR e dil B 37° Km CpG NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 100 100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by some combinations of overlapping CpG methylation (see p. 334).

Fnu4HI

#R0178S #R0178L

5′...G C[▼]N G C ... 3′ 3′... C G N_AC G ... 5′

Reaction Conditions: CutSmart

Buffer, 37°C

CutSmart RR OdiA 37° W CpG

NEBuffer 1.1 2.1 3.1 CutSmart 200 units68 € % Activity <10 <10 <10 100 1,000 units272 €

Concentration: 10,000 units/ml

Methylation Sensitivity:

Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).































FokI CutSmart Ri 2+site dil A 37° dcm CpG

#R0109S 1,000 units68 € #R0109L 5,000 units272 €

5′...G G A T G (N)₉ ▼...3′ 3′... C C T A C (N)₁₃... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 75 100

Methylation Sensitivity: Impaired by overlapping dcm methylation. Cleavage of mammalian genomic DNA is impaired by overlapping CpG methylation (see p. 334).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

HaeIII

#R0108S 3,000 units62 € #R0108L 15,000 units248 € for high (5X) concentration

3,000 units62 € #R0108T #R0108M 15,000 units248 €

5′...GG^TCC...3′ 3′...CC₄GG...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C

for 20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 25 100

CutSmart RR OdiA 37° W

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

FseI

#R0588S 100 units74 € #R05881 500 units296 €

5'...GGCCGG CC...3' 3′...CC,GGCCGG...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Rill O dil B 37° 655 dcm CpG NEBuffer 1.1 2.1 3.1 CutSmart

Concentration: 2,000 units/ml

% Activity 100 75 <10

Methylation Sensitivity: Impaired by some combinations of overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

5′... G A C G C (N), ▼... 3′ 3′... C T G C G (N)₁₀... 5′

Reaction Conditions: NEBuffer 1.1, 37°C. Heat inactivation: 65°C for

HgaI

#R0154S 100 units71 € #R0154L 500 units284 €

20 minutes

Concentration: 2,000 units/ml

RX NEB 1.1 dil A 37° 655 CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 25

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

FspI

#R0135S 500 units68 € #R0135L 2,500 units272 €

5′...TGC*****GCA...3′ 3′...ACGCGT...5′

Reaction Conditions: CutSmart Buffer, 37°C

CutSmart RR dilC 37° Mb CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 100 10 100

Concentration: 10,000 units/ml. Methylation Sensitivity: Cleavage

of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

HhaI

#R0139S 2,000 units60 € #R0139L 10,000 units240 €

5′...G C G^TC ...3′ 3′...C₄G C G ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RR dilA 37° CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 100 100

Concentration: 20,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

FspEI

CutSmart Ri Epi dil B 37° 186

#R0662S 200 units 106 € See page 264 for more information.

HaeII

#R0107S 2,000 units62 € #R0107L 10,000 units248 €

5'... RGCGC Y ... 3' 3'... Y C G C G R ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

CutSmart RR dilA 37° 66 CpG

% Activity 25 100 10 100 Concentration: 20.000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

HincII

#R0103S 1,000 units62 € #R0103L 5,000 units248 €

5′... G T Y R A C ... 3′ 3′... C A R Y T G ... 5′

Reaction Conditions: NEBuffer 3.1. 37°C. Heat inactivation: 65°C for 20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 25 100 100 100 Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see

p. 334).

500 units61 €

2,500 units244 €

HindIII



#R0104S	10,000 units	50€
#R0104L	50,000 units	200€
for high (5X) cond	entration	
#R0104T	10,000 units	50€
"Do 10 11 1		

#R0104M 50,000 units200 €

5´... A*A G C T T ... 3´ 3'... T T C G A A ... 5'

Reaction Conditions: NEBuffer 2.1, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 50 50

Concentration: 20,000 and 100.000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from extended digestion.

HpaI #R0105S

#R0105L

Buffer, 37°C

Concentration: 5,000 units/ml







NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 75 25 100

Methylation Sensitivity: Cleavage 5′... G T T*A A C ... 3′ of mammalian genomic DNA is 3′... C A A T T G ... 5′ blocked by some combinations of overlapping CpG methylation (see Reaction Conditions: CutSmart

p. 334).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

High-Fidelity

HindIII-HF®

10,000 units50 € #R3104S 50,000 units200 € #R3104L for high (5X) concentration #R3104T 10,000 units50 €

50,000 units200 €

5′... A*A G C T T ... 3′ 3′... T T C G A A ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

CutSmart Ril e dil B 37° Will

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 100 10 100

Concentration: 20.000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

HpaII

#R0171S 2,000 units65 € #R0171L 10,000 units260 € for high (5X) concentration

#R0171M 10,000 units260 €

5′...C^{*}CGG...3′ 3'... G G C C ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 50 <10

CutSmart Ri Epi dilA 37° til CpG

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Hinfl

#R3104M

#R0155S 5,000 units63 € #R0155L 25,000 units252 € for high (5X) concentration #R0155T 5,000 units63 € #R0155M 25,000 units252 €

5′... G[™]A N T C ... 3′ 3′... C T N A₄G ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

CutSmart Ril O dil A 37° 186 CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 100 100 Concentration: 10,000 and

50,000 units/ml Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see

p. 334).

HphI

#R0158S 1,000 units65 € #R0158L 5,000 units260 €

5′...G G T G A (N)₈ ... 3′ 3′... C C A C T (N)₇... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes

Concentration: 5.000 units/ml

CutSmart Ri dil B 37° dil B 37° dam dem

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 50 <10

Methylation Sensitivity: This enzyme is blocked by overlapping dam and dcm methylation (see p. 334).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

HinP1I

2,000 units68 € #R0124S

5′...G^{*}CGC...3′ 3'... C G C G ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Ril O dil A 37° 655 CpG NEBuffer 1.1 2.1 3.1 CutSmart

Concentration: 10,000 units/ml

% Activity 100 100 100 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Hpy99I

#R0615S 100 units 69 € #R0615L 500 units276 €

5′... C G W C G[▼]...3′ 3'... GCWGC ... 5'

Reaction Conditions: CutSmart Buffer Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Ri dil A 37° CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 10 <10 100

Concentration: 2,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

























Hpy166II

#R0616S 1,000 units72 €

5′... G T N N A C ... 3′ 3′... C A N N T G ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RR diiC 37° CpG

 NEBuffer
 1.1
 2.1
 3.1
 CutSmart

 % Activity
 100
 100
 50
 100

Concentration: 10,000 units/ml

Methylation Sensitivity:

Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

HpyCH4III

#R0618S 250 units70 € #R0618L 1,250 units280 €

5'...ACNGT...3' 3'...TGNCA...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RRI driA 37° 66° NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 25 <10</td> 100

Concentration: 5,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Hpy188I

#R0617S 1,000 units69 € #R0617L 5,000 units276 €

5′...TCN GA...3′ 3′...AGNCT...5′

Reaction Conditions: CutSmart Buffer, 37°C.Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

CutSmart Rik dil A 37° dil dam

 NEBuffer
 1.1
 2.1
 3.1
 CutSmart

 % Activity
 25
 100
 50
 100

Methylation Sensitivity: Blocked by overlapping *dam* methylation (see p. 334).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

HpyCH4IV

#R0619S 500 units70 € #R0619L 2,500 units280 €

5'... A C G T ... 3' 3'... T G C A ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

 NEBuffer
 1.1
 2.1
 3.1
 cutSmart

 % Activity
 100
 50
 25
 100

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Hpy188III

#R0622S 500 units69 € #R0622L 2,500 units276 €

5′...T C^TN N G A ...3′ 3′...A G N N_AC T ...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

CutSmart Ril dii B 37° 1654 dam CpG

% Activity 100 100 10 100

Methylation Sensitivity: Blocked by overlapping dam methylation.
Cleavage of mammalian genomic

NEBuffer 1.1 2.1 3.1 CutSmart

Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.

HpyCH4V

#R0620S 100 units69 € #R0620L 500 units276 €

5′...T G^VC A ...3′ 3′...A C_sG T ...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart

CutSmart Ril O dil A 37° 165

Concentration: 5,000 units/ml

% Activity 50 50 25 100

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

HpyAV

#R0621S 100 units72 € #R0621L 500 units288 €

 $5' \dots G G T T G (N)_6^{\blacktriangledown} \dots 3'$ $3' \dots G G A A G (N)_{5_{\blacktriangle}} \dots 5'$

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes

Concentration: 2,000 units/ml

CutSmart RX O 37° 654 CpG

 NEBuffer
 1.1
 2.1
 3.1
 CutSmart

 % Activity
 100
 100
 25
 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by overlapping CpG methylation. (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.

KasI

#R0544S 250 units69 € #R0544L 1,250 units276 €

5′... G^TG C G C C ... 3′ 3′... C C G C G₁G ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

 $\textbf{Concentration:}\ 5{,}000\ units/ml$

CutSmart RR dil B 37° 654 CpG

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 50 100 50 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.

glycerol concentration of > 5%.

KpnI RN NEB 1.1 dil A 37° 116 #R0142S 4,000 units62 €

#R0142L 20,000 units248 € for high (5X) concentration 20,000 units248 € #R0142M

5'...GGTAC C...3' 3'...C_CATGG...5'

Reaction Conditions: NEBuffer 1.1,

Concentration: 10,000 and 50.000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 75 <10 50

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

MfeI

#R0589S 500 units74 € NEBuffer 1.1 2.1 3.1 CutSmart #R0589L 2,500 units296 € % Activity 75 50 10 100

5'...C"A A T T G ... 3' 3′...GTTAAC...5′

Reaction Conditions: CutSmart

Buffer, 37°C

Concentration: 10.000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian

CutSmart Ri dil A 37° Vib

Note: Star activity may result from extended digestion.

CpG methylation.

High-Fidelity

KpnI-HF®

#R3142S 4.000 units62 € #R3142L 20,000 units248 € for high (5X) concentration

20,000 units248 €

5'...GGTAC C...3' 3'...CCATGG...5'

Reaction Conditions: CutSmart

Buffer, 37°C

#R3142M

CutSmart RR e e dil A 37° Wh

% Activity 100 25 <10 100 Concentration: 20,000 and 100.000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

LpnPI #R0663S 200 units108 € CutSmart Rill Epi dil B 37° 165

See page 264 for more information.

MboI CutSmart RR dilA 37° dilA CpG

#R0147S 500 units74 € #R0147L 2,500 units296 € for high (5X) concentration

#R0147M 2,500 units296 €

5″...**™**G A T C ... 3″ 3'... C T A G₄... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 75 100 100 100

Concentration: 5,000 and 25,000 units/ml

Methylation Sensitivity: Blocked by dam methylation. Its isoschizomer Sau3Al is not. Cleavage of mammalian genomic DNA is impaired by overlapping CpG methylation (see p. 334).

MfeI-HF®

#R3589S #R3589L

5′...C*A A T T G ... 3′ 3′...GTTAA,C...5′

Reaction Conditions: CutSmart

Buffer, 37°C

High-Fidelity

CutSmart RR e dii A 37° Wh

500 units74 € NEBuffer 1.1 2.1 3.1 CutSmart 2,500 units296 € % Activity 75 25 <10 100

Concentration: 20,000 units/ml Methylation Sensitivity:

Not sensitive to dam, dcm or mammalian CpG methylation.

MluI

#R0198S 1,000 units66 € #R0198L 5,000 units264 €

5'... A CGCGT...3' 3′...TGCGC₄A...5′

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 80°C for 20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 50 100

Concentration: 10,000 units/ml Methylation Sensitivity: Cleavage

of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

MboII

#R0148S 300 units68 € #R0148L 1,500 units272 €

5′...G A A G A (N), ♥...3′ 3′...CTTCT(N)₇...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Rik 2*site dil C 37° 655 dam

% Activity 100* 100 50 100 Concentration: 5,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart

Methylation Sensitivity: Blocked by overlapping dam methylation (see p. 334).

Note: *May exhibit star activity in this buffer.

MluI-HF®

5′... A CGCGT...3′

1,000 units 66 €

#R3198S #R3198L 5,000 units264 €

3′...TGCGC₄A...5′

Reaction Conditions: CutSmart Buffer, 37°C

% Activity 25 100 100 100 Concentration: 20,000 units/ml

CutSmart Rill e dilA 37° th CpG

NEBuffer 1.1 2.1 3.1 CutSmart

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see

p. 334).





























High-Fidelity





MluCI CutSmart RR dilA 37° Wb

#R0538S 1,000 units64 € #R0538L 5,000 units256 €

5′...*A A T T ...3′ 3′... T T A A₄... 5′

Reaction Conditions: CutSmart

Buffer, 37°C

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 10 10 100

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian

CpG methylation.

MscI

#R0534S 250 units72 € #R0534L 1,250 units288 € for high (5X) concentration

#R0534M 1,250 units288 €

5′...TGG^{*}CCA...3′ 3'... A C C,G G T ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 100 100

CutSmart Ril dilC 37° km dcm

Concentration: 5,000 and 25.000 units/ml

Methylation Sensitivity: Blocked by overlapping *dcm* methylation (see p. 334). The single Mscl site in pBR322 overlaps a dcm methylation site; consequently, pBR322 which has been grown in a dcm host should be used for cloning.

MlyI

#R0610S 1,000 units67 € #R0610L 5,000 units268 €

5′... G A G T C (N), [▼]... 3′ 3′... C T C A G (N)₅....5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.



Concentration: 10.000 units/ml

% Activity 50 50 10

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

MseI

#R0525S 500 units69 € #R0525L 2,500 units276 € for high (5X) concentration #R0525M 2,500 units276 €

5′... T[▼]T A A ... 3′ 3′... A A T_AT ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 75 100 75 100

CutSmart Ril O dil A 37° (65)

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

MmeI

#R0637S 100 units69 € #R0637L 500 units276 €

5′...T C C R A C (N)₂₀ ... 3′ 3′...AGGYTG(N)₁₈...5′

Reaction Conditions: CutSmart Buffer + SAM, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 50 100

CutSmart Ri 2*site SAM dii B 37° 655 CpG

Concentration: 2,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

MslI

#R0571S #R0571L

5′...CAYNN[¶]NNRTG...3′ 3′...GTRNNNYAC...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

500 units72 € 2,500 units288 €

Concentration: 10,000 units/ml Methylation Sensitivity: Not sensitive to dam, dcm or mammalian

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 50 50 <10 100

CutSmart Ril O dil A 37° km

CpG methylation.

MnlI

#R0163S 500 units71 € #R0163L 2,500 units284 €

5′... C C T C (N)₇ ... 3′ 3′... G G A G (N)₆,... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes

CutSmart RX dil B 37° 65 NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 75 100 50 100 Concentration: 5,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

MspI

#R0106S 5,000 units 64 € #R0106L for high (5X) concentration #R0106T #R0106M

5'...C'CGG...3' 3'...GGC,C...5'

CutSmart Rik Epi 😻 dil A 37° 🚻

25,000 units256 € 5,000 units64 € 25,000 units256 €

Concentration: 20,000 and 100,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmar

% Activity 75 100 50 100

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Reaction Conditions: CutSmart

Buffer, 37°C

MspA1I

#R0577S

#R0577L

500 units72 € 2,500 units288 €

5′... C M G^{*}C K G ... 3′ 3′...G K C₄G M C ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Rill Odil B 37° 65 CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 50 10

Concentration: 10,000 units/ml

Methylation Sensitivity:

Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

NciI

#R0196S 2,000 units66 € #R0196L 10,000 units264 €

5'... C C S G G ... 3' 3'... G G S₄C C ... 5'

Reaction Conditions: CutSmart

Buffer, 37°C

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 25 10 100

CutSmart RR dilA 37° Mb CpG

Concentration: 20,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by overlapping CpG methylation (see p. 334).

MspJI

CutSmart Rill Epi dil B 37° (65)

#R0661S 200 units108 €

See page 264 for more information. #R0661L 1,000 units432 €

MwoI

#R0573S 500 units72 € #R0573L 2,500 units288 €

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 100 100 100

CutSmart RR Odil B 60° CpG

5'...GCNNNNNNNGC...3' $3^{\prime}\dots$ CGNNNNNNNCG...5 $^{\prime}$

Reaction Conditions: CutSmart

Buffer, 60°C

Concentration: 5,000 units/ml

Activity at 37°C: 10%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

NcoI

#R0193M

#R0193S 1,000 units65 € #R0193L 5,000 units260 € for high (5X) concentration #R0193T 1,000 units65 €

5,000 units260 €

5'...C"CATGG...3' 3'... G G T A C_AC ... 5'

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

50,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 100 100 100 100

Concentration: 10.000 and

NaeI

CutSmart Rill 2*site dil A 37° 16 CpG

#R0190S 500 units63 € #R0190L 2,500 units252 €

% Activity 25 25 <10 100

5′...GCC GGC ...3′ 3′... C G G C C G ... 5′

Reaction Conditions: CutSmart Buffer 37°C

NEBuffer 1.1 2.1 3.1 CutSmart

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

NcoI-HF®

#R3193S 1.000 units65 € #R3193L 5,000 units260 €

for high (5X) concentration #R3193M 5,000 units260 €

5'...C"CATGG...3' 3'...GGTACC...5'

Reaction Conditions: CutSmart Buffer, 37°C, Heat inactivation: 80°C for 20 minutes.

High-Fidelity

CutSmart RR C dil B 37° km

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 10

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

NarI

#R0191S 500 units68 € #R0191L 2,500 units272 €

5′... G G C G C C ... 3′ 3′...CCGC**_**GG...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Rik 2*site dil A 37° CpG

NEBuffer 1.1 2.1 3.1 CutSmar % Activity 100 100 10

Concentration: 5,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

NdeI

#R0111S 4,000 units 66 € #R0111L 20,000 units 264 €

5'...CA"TATG...3' 3'...GTATAC...5

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C

for 20 minutes.

CutSmart RR dilA 37° 165

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 75 100 100 100

Concentration: 20,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

































NgoMIV

#R0564S 1,000 units66 € #R0564L 5,000 units264 €

5'...G CGGC...3' 3'... C G G C C₄G ... 5'

Reaction Conditions: Cutsmart

Buffer, 37°C

Concentration: 10,000 units/ml

CutSmart RX 2+site dilA 37° CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 50 10 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

NlaIV #R0126S

200 units69 €

3′...CCNNGG...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 10 10 100

CutSmart Rik dil B 37° dim CpG

Concentration: 2,000 units/ml

Methylation Sensitivity: Blocked by overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

NheI

#R0131S 1,000 units68 € 5,000 units272 € #R01311 for high (5X) concentration

#R0131M 5,000 units272 €

5'...G'CTAGC...3' 3'... C G A T C₄G ... 5'

Reaction Conditions: NEBuffer 2.1, 37°C Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 10

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

High-Fidelity

NmeAIII

#R0711S 250 units72 €

5′...GCCGAG(N)₂₀₋₂₁...3′ 3'... CGGCTC(N)₁₈₋₁₉...5

Note: The cleavage point may shift one base pair depending on the DNA sequence context between the recognition site and the position of cleavage. For a given sequence, generally one site will predominate. For details, see www.neb.com

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 10 <10 100

CutSmart R 2+site SAM dil B 37° (65)

Reaction Conditions: CutSmart Buffer + SAM, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 2,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

NheI-HF®

#R3131S 1,000 units68 € #R3131L 5,000 units272 € for high (5X) concentration #R3131M 5,000 units272 €

5'...G"C TAGC...3' 3′...CGATC<u>.</u>G...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 25 <10 100

CutSmart RR e dilC 37° 60 CpG

Concentration: 20,000 and 100.000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

NotI

#R0189M

500 units72 € #R0189S #R0189L 2,500 units288 € for high (5X) concentration

2,500 units288 €

5′...GC*****GGCCGC...3′

3′...CGCCGGCG...5′

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 50 100 25 Concentration: 10,000 and

50 000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see

p. 334).

NlaIII

#R0125S 500 units69 € #R0125L 2,500 units276 €

5′... C A T G[▼]... 3′ 3′...<mark>G</mark> T A C ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RR G dil B 37° 1654

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 <10 <10 100

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

NotI-HF®

#R3189S 500 units72 € #R3189L 2,500 units288 €

for high (5X) concentration

#R3189M 2,500 units288 €

5'...GC GGCCGC...3' 3'...CGCCGGCG...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Rill e dilA 37° tt CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 25 100

High-Fidelity

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

NruI

#R0192S 1,000 units62 € #R0192L 5,000 units248 € for high (5X) concentration

1,000 units62 € #R0192T #R0192M 5,000 units248 €

5'... T C G C G A ... 3' 3'... A G C,G C T ... 5'

Reaction Conditions: NEBuffer 3.1,

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 10 100 10

Concentration: 10,000 and

50.000 units/ml

Methylation Sensitivity: Blocked by overlapping dam methylation. Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

NspI #R0602S

#R0602L

CutSmart RR OdiA 37°

250 units70 € NEBuffer 1.1 2.1 3.1 CutSmart 1,250 units280 € % Activity 100 100 <10

Concentration: 10,000 units/ml 5'...RCATG Y ...3' 3′...Y_GTACR...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C

for 20 minutes.

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

High-Fidelity

NruI-HF®

#R3192S 1,000 units62 € #R3192L 5,000 units248 €

5'... T C G C G A ... 3' 3′... A G C₄G C T ... 5′

Reaction Conditions: CutSmart Buffer, 37°C.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 0 25 50 100

CutSmart Rill e dilA 37° Mil dam CpG

Concentration: 20,000 units/ml

Methylation Sensitivity: Blocked by overlapping dam methylation. Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

PacI

#R0547S 250 units 69 € #R0547L 1,250 units276 €

5′... T T A A T[▼]T A A ... 3′ 3′... A A T_AT A A T T ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Ril O dil A 37° 165

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 75 10 100

Concentration: 10.000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

NsiI

#R0127S 1,000 units66 € #R0127L 5,000 units264 €

5′... A T G C A T ... 3′ 3′...T,ACGTA...5′

Reaction Conditions: NEBuffer 3.1. 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 75 100 25

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian

Concentration: 10,000 units/ml

CpG methylation.

PaeR7I

#R0177S 2,000 units62 €

5'...C"T C G A G ... 3' 3′...GAGCT₄C...5′

Reaction Conditions: CutSmart

Buffer, 37°C

CutSmart RR dilA 37° Mb CpG NEBuffer 1.1 2.1 3.1 CutSmart

Concentration: 20,000 units/ml

% Activity 25 100 10

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see

p. 334).

High-Fidelity

NsiI-HF®

..66€ #R3127S 1,000 units #R3127L 5,000 units264 €

5′... A T G C A[▼]T ... 3′ 3'... T_A C G T A ... 5

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

CutSmart Rik e e dil B 37° km

% Activity <10 20 <10 Concentration: 20,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

PciI

#R0655S 200 units 72 € #R0655L 1,000 units 288 €

5′... A CATGT...3′ 3'... T G T A CA... 5'

Reaction Conditions: NEBuffer 3.1. 37°C. Heat inactivation: 80°C for

20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 75 100

Rii NEB 3.1 dil B 37° 💖

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: *May exhibit star activity in this buffer.





























PflFI CutSmart RR dilA 37° 655

#R0595S 2,000 units66 €

5'...GACN NNGTC...3' 3'...CTGNNNCAG...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart

Methylation Sensitivity: Not

% Activity 25 100 25 100

sensitive to dam, dcm or mammalian CpG methylation.

PmeI

#R0560S 500 units76 € #R0560L 2,500 units304 €

5'...GTTT"AAAC...3' 3′...CAAATTTG...5′

Reaction Conditions: CutSmart Buffer, 37°C, Heat inactivation: 65°C for 20 minutes.

Concentration: 10.000 units/ml

CutSmart RR dilA 37° 655 CpG NEBuffer 1.1 2.1 3.1 CutSmart

% Activity <10 50 10 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

PflMI

#R0509S 1.000 units66 € #R0509L 5,000 units264 €

5'... C C A N N N N T G G ... 3' 3′... G G T N N N N N A C C ... 5′

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml



NEBuffer 1.1 2.1 3.1 CutSmart % Activity 0 100 100

Methylation Sensitivity: Blocked by overlapping dcm methylation (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.

#R0532S 2,000 units #R0532L 10,000 units272 €

5′...CAC GTG...3′

Reaction Conditions: CutSmart

PmlI

.68€

3'...GTGCAC...5'

Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 20,000 units/ml



% Activity 100 50 <10 100 Methylation Sensitivity: Cleavage

of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

PleI

#R0515S 1,000 units 66 €

5′... G A G T C (N)₄ ... 3′ 3′... C T C A G (N)₅... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

CutSmart Ri 2+site dil A 37° 655 CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 50 25 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

PpuMI

#R0506S 500 units69 € #R0506L 2,500 units276 €

5'...RG WCCY...3' 3′...YCCWG₄GR...5′

Reaction Conditions: CutSmart

Buffer, 37°C

Concentration: 10,000 units/ml

CutSmart RR dil B 37° Wb dcm

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 <10 <10 100

Methylation Sensitivity: Blocked

by overlapping dcm methylation

(see p. 334).

PluTI

500 units68 € #R0713S

5′...GGCGC^{*}C...3′ 3'...C,CGCGG...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

CutSmart Ri 2+site dil A 37° 655 CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 25 <10 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

PshAI

#R0593S 1.000 units70 € #R0593L 5,000 units280 €

5'... G A C N N N N G T C ... 3 3'...CTGNNNNCAG...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

CutSmart RR dilA 37° 655 CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 50 10 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see

p. 334).

10,000 units66 €

PsiI CutSmart Rill dil B 37° 65 200 units 106 € #R0657S NEBuffer 1.1 2.1 3.1 CutSmart #R0657L 1,000 units 424 €

5′...T T A[▼]T A A ...3′ 3′...A A T_A T T ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

% Activity 10 100 10


Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%.

PstI #R0140S

#R0140L

#R0140T



CpG methylation.

buffer.

Methylation Sensitivity: Not

sensitive to dam, dcm or mammalian

Note: *May exhibit star activity in this

50,000 units264 € for high (5X) concentration 10,000 units66 € Concentration: 20,000 and 100 000 units/ml #R0140M 50,000 units264 €

5′...CTGCA^{*}G...3′ 3'... GACGTC... 5'

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

PspGI

#R0611S 1,000 units 66 €

5′...*****CCWGG...3′ 3'... G G W C C ... 5'

Reaction Conditions: CutSmart

Buffer, 75°C

Concentration: 10,000 units/ml

CutSmart Ri dil A 75° Wid dcm

100

% Activity 25 100 50 Activity at 37°C: 10%

NEBuffer 1.1 2.1 3.1 CutSmart

Methylation Sensitivity: Blocked by dcm methylation (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.

PstI-HF®

#R3140S .66€ 10,000 units 50,000 units264 € #R3140L for high (5X) concentration

#R3140T 10,000 units66 € #R3140M 50,000 units264 €

5′...CTGCA G...3′ 3′...G,A C G T C ... 5′

Reaction Conditions: CutSmart

Buffer, 37°C

High-Fidelity

CutSmart RR e dilC 37° Mb

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 75 50 100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

PspOMI

#R0653S 1,500 units 62 € #R0653L 7,500 units 248 €

5′...G*****GGCCC...3′ 3′...CCCGG,G...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Rill dil B 37° 65 dcm CpG NEBuffer 1.1 2.1 3.1 CutSmart

Concentration: 20,000 units/ml

% Activity 10 10 <10

Methylation Sensitivity: Impaired by some combinations of overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

PvuI

#R0150S 500 units71 € #R0150L 2,500 units284 €

5'...CGATCG...3' 3′...GC_TAGC...5′

Reaction Conditions: NEBuffer 3.1.

37°C

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 25 100 <10

Methylation Sensitivity: Cleavage

of mammalian genomic DNA is blocked by CpG methylation (see

Concentration: 10,000 units/ml

p. 334).

PspXI

#R0656S 200 units 72 € #R0656L 1,000 units 288 €

5'... V CTC G A G B ... 3' 3′...BGAGCT_CV...5′

Reaction Conditions: CutSmart Buffer, 37°C

CutSmart Rik dil B 37° Wil CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 100 25 100 Concentration: 5,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by CpG methylation (see page 334).

PvuI-HF®





#R3150S 500 units71 € #R3150L 2,500 units284 €

5'...CGATCG...3' 3′...GC_TAGC...5′

Reaction Conditions: CutSmart

Buffer, 37°C

% Activity 25 100 100 100 Concentration: 20,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see

p. 334).

























High-Fidelity

High-Fidelity

PvuII		
#R0151S	5,000 units	60€
#R0151L	25,000 units	240 €
for high (5X) con	centration	
#R0151T	5,000 units	60€
#R0151M	25,000 units	240 €
5′ C A G	CTG3′	

Dynati

Reaction Conditions: NEBuffer 3.1, 37°C

54	R	•	NEB.	3.1 dil B	37° 166
NEBuffer	1.1	2.1	3.1	CutSma	rt

Concentration: 10,000 and 50,000 units/ml

% Activity 50 100 100 100*

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

*May exhibit star activity in this buffer.

SacI #R0156S

#R0156L



for high (5X) concentration

#R0156M 10,000 units244 €

5′...GAGCT^{*}C...3′ 3′... C_AT C G A G ... 5′

Reaction Conditions: NEBuffer 1.1, 37°C. Heat inactivation: 65°C for 20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 50 10

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

High-Fidelity

PvuII-HF®

3′...GTC₄GAC...5′

#R3151S 5.000 units 60 € 25,000 units240 € #R3151L for high (5X) concentration #R3151M 25,000 units240 €

5'... C A G C T G ... 3' 3'... G T C₄G A C ... 5'

Reaction Conditions: CutSmart

Buffer, 37°C.

CutSmart Rill e dil B 37° Wb

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 <10 <10 100

Concentration: 20,000 and 100 000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

SacI-HF®

#R3156S 2,000 units61 € 10,000 units244 € #R3156L for high (5X) concentration

#R3156M 10,000 units244 €

5′...GAGCT C...3′ 3′...C_AT C G A G ... 5′

Reaction Conditions: Cutsmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 50 <10 100

CutSmart Rill e dilA 37° 655 CpG

Concentration: 20,000 and 100 000 units/ml

Methylation Sensitivity: Not sensitive to dam or dcm methylation. Blocked by some combinations of overlapping CpG methylation.

RsaI

#R0167S 1.000 units60 € #R0167L 5,000 units240 €

5′...G T[™]A C ... 3′ 3′...C A T G ... 5′

Reaction Conditions: CutSmart

Buffer, 37°C

CutSmart RR O dil A 37° WW CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 50 <10 100

Concentration: 10,000 units/ml Methylation Sensitivity: Cleavage

of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

SacII

#R0157S 2.000 units63 € #R0157L 10,000 units252 €

5′...CCGC^{*}GG...3′ 3'...GG_CGCC...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes

CutSmart Ri 2*site dilA 37° CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 100 10 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Concentration: 20.000 units/ml

RsrII

#R0501S 500 units74 € #R0501L 2,500 units296 €

5´...CG G W C C G...3´ 3'...G C C W G,G C ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RR 2+site dil C 37° CpG NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 25 75 10 100 Concentration: 5,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

SalI

#R0138S 2,000 units59 € #R0138L 10,000 units236 € for high (5X) concentration 2,000 units59 € #R0138T #R0138M 10,000 units236 €

5'...G"T CGAC...3' 3′... C A G C T₄G ... 5′

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 65°C for 20 minutes.





NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 <10 100 <10

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

High-Fidelity Sall-HF® CutSmart Rill e dilA 37° 655 CpG

#R3138S 2.000 units 59 € #R3138L 10,000 units 236 € for high (5X) concentration

#R3138T 2,000 units 59 € #R3138M 10,000 units 236 €

5′...G^TT C G A C ... 3′ 3′...CAGCT_G...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 100 100 100

Concentration: 20,000 and 100.000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

SbfI

CutSmart Ril dil A 37° 180°

#R0642S 500 units 74 € #R0642L 2,500 units 296 €

5'...CCTGCA GG...3' 3′...GG,ACGTCC...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 25 <10

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%.

SapI

#R0569S 250 units 68 € #R0569L 1.250 units 272 €

5′...GCTCTTC(N), ▼...3′ 3′... C G A G A A G (N)₄... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Ril O dil B 37° 165

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 75 50 >10

Concentration: 10,000 units/ml Methylation Sensitivity: Not

sensitive to dam, dcm or mammalian CpG methylation.

SbfI-HF®

#R3642S 500 units 74 € 2,500 units 296 € #R3642L

5′...CCTGCA^VGG...3′ 3′...GG_AACGTCC...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

High-Fidelity CutSmart RR e dilB 37° Will

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 25 <10

Concentration: 20,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Sau3AI

#R0169S 200 units62 € #R0169L 1,000 units248 €

5′...[™]G A T C ... 3′ 3′... C T A G₄... 5′

Reaction Conditions: NEBuffer 1.1. 37°C. Heat inactivation: 65°C for 20 minutes

Concentration: 5,000 units/ml

RR MEB 1.1 dil A 37° 655 CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 50 10 100

Methylation Sensitivity: Unlike DpnII and Mbol, Sau3Al is not blocked by dam methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

ScaI-HF®

#R3122S 1,000 units 66€ #R3122L 5,000 units 264 €

for high (5X) concentration

#R3122M 5,000 units 264 €

5′... A G T*A C T ... 3′ 3′... T C A T G A ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

High-Fidelity CutSmart RR e dil B 37° W

NEBuffer 1.1 2.1 3.1 CutSmart

Concentration: 20,000 and 100,000 units/ml

% Activity 100 100 10 100

Methylation Sensitivity: Not sensitive to dam. dcm or mammalian CpG methylation.

Sau96I

#R0165S 1,000 units 62 €

5′...G G N C C ... 3′ 3′... C C N G₄G ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

CutSmart Ri dii A 37° dim CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 100 100

Methylation Sensitivity: Blocked by overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

ScrFI

#R0110S 1,000 units 68 €

5′...CC[™]NGG...3′ 3'...GGNCC...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart

CutSmart Ri dil C 37° dim CpG

% Activity 100 100 100 100 Methylation Sensitivity: Blocked

by overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 334).

Note: Star activity may result from extended digestion.



























SexAI

#R0605S 200 units 68 € #R0605L 1,000 units 272 €

5′...A CCWGGT...3′ 3′...TGGWCC_A...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Rik dil A 37° dim

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 75 50 100

Concentration: 5,000 units/ml

Methylation Sensitivity: Blocked by dcm methylation (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.

SfoI



#R0606S 500 units 66 € #R0606L 2,500 units 264 €

5′...GGC^{*}GCC...3′ 3′...CCG_{*}CGG...5′

Reaction Conditions: CutSmart

Buffer, 37°C

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 100 100

Methylation Sensitivity: Blocked by some combinations of overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

SfaNI

#R0172S 300 units 104 € #R0172L 1,500 units 416 €

5′...G C A T C (N)₅...3′ 3′... C G T A G (N)₉... 5′

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 2.000 units/ml

RX NEB 3.1 dil B 37° 165 CpG

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 75 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by some combinations of overlapping CpG methylation (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.

SgrAI

#R0603S 1,000 units 69 € #R0603L 5,000 units 276 €

5'... CR CGGYG...3' 3′...GYGGCC_RC...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

CutSmart RX 2+site dil A 37° CpG

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 100 100 10 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

SfcI

#R0561S 200 units 66 € #R0561L 1,000 units 264 €

5′...CTRYAG...3′ 3′...GAYRT_C...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RR dil B 37° 65

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 75 50 25

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%.

Smal

#R0141S 2,000 units 58 € #R0141L 10,000 units 232 €

5′...CCC**~**GGG...3′ 3′...GGG**_**CCC...5′

Reaction Conditions: CutSmart Buffer, 25°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RR dil B 25° 655 CpG NEBuffer 1.1 2.1 3.1 CutSmart

% Activity <10 <10 <10 100 Concentration: 20,000 units/ml

Activity at 37°C: 50% (15 minute

half-life)

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

SfiI

CutSmart Ri 2*site dil C 50° th dcm CpG

#R0123S 3,000 units 70 € #R0123L 15,000 units 280 €

5'...GGCCNNNN NGGCC...3' 3'...C CGGNNNNNCCGG...5'

Reaction Conditions: CutSmart Buffer, 50°C

Concentration: 20,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 50 100

Activity at 37°C: 10%

Methylation Sensitivity: Impaired by overlapping dcm methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

SmlI

#R0597S 500 units 69 € #R0597L 2,500 units 276 €

5′...C^TTYRAG...3′ 3′...GARYT₄C...5′

Reaction Conditions: CutSmart Buffer, 55°C

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 75 25 100

CutSmart Rill dil A 55° Will

Activity at 37°C: 10%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

SnaBI CutSmart Ri dil A 37° W CpG

#R0130S 500 units 62 € #R0130L 2,500 units 248 €

for high (5X) concentration 2,500 units 248 € #R0130M

5′... TAC[▼]G TA...3′ 3′... AT G_C AT ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 5,000 and 25,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 50 10 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

SphI-HF®

#R3182S 500 units 66 € #R3182L 2,500 units 264 € for high (5X) concentration #R3182M 2,500 units 264 €

5'...GCATG C...3' 3'...CGTACG...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RR C dil B 37° NEBuffer 1.1 2.1 3.1 CutSmart

100.000 units/ml

% Activity 50 25 10 100 Concentration: 20,000 and

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

SpeI

#R0133M

#R0133S 500 units 66 € #R0133L 2,500 units 264 € for high (5X) concentration

2,500 units 264 €

5'... A CTAGT...3' 3′... T G A T CA ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

CutSmart RR OdinC 37° W

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 75 100 25

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

SrfI

#R0629S 500 units 81 € #R0629L 2,500 units 360 €

3′...CGGG,CCCG...5′

for 20 minutes.

CutSmart Rill e G dil B 37° CpG

5'...GCCC GGGC...3'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C

Concentration: 20,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 50 <10

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p 334).

SpeI-HF®

#R3133S 500 units ... 66€ #R3133L 2,500 units 264 € for high (5X) concentration #R3133M 2,500 units 264 €

5′... A C T A G T ... 3′ 3′... T G A T CA ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 20,000 and 100,000 units/ml

High-Fidelity

CutSmart RR e dilC 37° Will

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 50 <10

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

SspI

#R0132S 1,000 units 71 € #R0132L 5.000 units 284 € for high (5X) concentration

#R0132M 5,000 units 284 €

5′... A A T^{*}A T T ... 3′ 3′... T T A_AT A A ... 5′

Reaction Conditions: NEBuffer Sspl. 37°C. Heat inactivation: 65°C for 20 minutes.

% Activity 50 100 50 Concentration: 5,000 and 25.000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

SphI

#R0182S 500 units 66 € #R0182L 2,500 units 264 € for high (8X) concentration 2,500 units 264 € #R0182M

5'...GCATG C ... 3' 3'... CAGTACG... 5'

Reaction Conditions: NEBuffer 2.1, 37°C. Heat inactivation: 65°C for 20 minutes.

RR NEB 2.1 dil B 37° KM

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 100 50 Concentration: 10,000 and

80,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from extended digestion.

SspI-HF®

#R3132S 1,000 units 71 € #R3132L 5,000 units 284 € for high (5X) concentration #R3132M 5,000 units 284 €

5′... A A T^TA T T ... 3′ 3′... T T A_AT A A ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

High-Fidelity

CutSmart Rill e dil B 37° 165

% Activity 25 100 <10 100 Concentration: 20.000 and 100,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian

CpG methylation.

































StuI CutSmart Ri diA 37° th dcm

#R0187S 1,000 units 62 € #R0187L 5,000 units 248 € for high (10X) concentration #R0187M 5,000 units 248 €

5'... A G G C C T ... 3' 3'... T C C₄G G A ... 5'

Reaction Conditions: CutSmart

Buffer, 37°C

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 50


Concentration: 10,000 and 100,000 units/ml

Methylation Sensitivity: Blocked by overlapping *dcm* methylation (see p. 334).

SwaI #R0604S

#R0604L







Activity at 37°C: 50%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Reaction Conditions: NEBuffer 3.1, 25°C. Heat inactivation: 65°C for

5′... A T T T[™]A A A T ... 3′

3′... T A A A A T T T A ... 5′

20 minutes.

Concentration: 10,000 units/ml

StyI

#R0500S 3,000 units 72 €

5'...C V W W G G ... 3' 3′...G G W W C₄C ...5′

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 25 100

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

$Taq^{\alpha}I$

#R0149S 4,000 units #R0149I for high (5X) concentration

4,000 units 62 € #R0149T #R0149M

5′... T[▼]C G A ... 3′ 3′... A G C_AT ... 5′

Reaction Conditions: CutSmart Buffer, 65°C. Heat inactivation: 80°C for 20 minutes

62€ NEBuffer 1.1 2.1 3.1 CutSmart 20,000 units 248 € % Activity 50 75 100 100

Concentration: 20,000 and 100,000 units/ml 20,000 units 248 €

Activity at 37°C: 10%

Methylation Sensitivity: Blocked by overlapping dam methylation (see p. 334).

CutSmart RR dil B 65° dil dam

StyI-HF®

3,000 units 72 € #R3500S #R3500L 15,000 units 288 €

5′...C^{*}C W W G G ...3′ 3′...G G W W C_AC ...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

High-Fidelity CutSmart Rill e dilA 37° 165

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 25 100

Concentration: 20,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

TfiI

#R0546S 500 units 70 €

5'... GAWTC...3' 3′... C T W A₄G ... 5′

Reaction Conditions: CutSmart

Buffer, 65°C

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 100 100

CutSmart RR OdilC 65° WM CpG

Activity at 37°C: 10%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

StvD4I

#R0638S 200 units 66 €

5′...*CCNGG...3′ 3'... G G N C C ... 5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

CutSmart Ri dil B 37° di dcm CpG NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 10 100 100 100

methylation. (see p. 334).

Methylation Sensitivity: Blocked by overlapping dcm methylation. Cleavage of mammalian genomic DNA is impaired by overlapping CpG

TseI

#R0591S 75 units 66 € #R0591L 375 units 264 €

5′...G WGC...3′ 3'... C G W C G ... 5'

Reaction Conditions: CutSmart Buffer, 65°C

Concentration: 5,000 units/ml. Activity at 37°C: 20%

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 75 100 100 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.

15,000 units 284 €

3,000 units 71 €

15,000 units 284 €

Tsp45I

#R0583S

#R0583L

200 units 72 € 1,000 units 288 €

5′... **™**G T S A C ... 3′ 3′... C A S T G₄...5′

Reaction Conditions: CutSmart

Concentration: 5,000 units/ml

CutSmart dil A 65°

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 100 50 <10

Activity at 37°C: 10%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

5′...T^{*}CTAGA...3′ 3′... A G A T C_AT ... 5′

for high (5X) concentration

XbaI

#R0145S

#R0145L

#R0145T

#R0145M

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RR O dii A 37° K dam 3,000 units 71 €

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 100 75 100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Blocked by overlapping dam methylation (see p. 334).

TspMI

#R0709S 200 units 72 €

5′... C'C C G G G ... 3′ 3′... GGGCC<u>C</u>C ... 5′

Reaction Conditions: CutSmart Buffer, 75°C

Concentration: 5,000 units/ml

Activity at 37°C: 20%

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50* 75* 50* 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Note: *May exhibit star activity in

this buffer.

XcmI

#R0533S 1,000 units 69 € #R0533L 5,000 units 276 €

37°C. Heat inactivation: 65°C for 20 minutes.

5'...C CANNNN NNNTGG...3' 3'...GGTNNNNNNNNACC...5'

Reaction Conditions: NEBuffer 2.1,

Concentration: 5,000 units/ml

RX NEB 2.1 dil C 37° 165

1.1 **2.1** 3.1 CutSmart NEBuffer % Activity 10 100 25

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Note: Star activity may result from extended digestion.

TspRI

#R0582S 1,000 units ... #R0582L

5′... N N C A S T G N N ▼... 3′ 3'...NNGTSACNN ...5'

Reaction Conditions: CutSmart Buffer, 65°C

Concentration: 10,000 units/ml

CutSmart RR 65° Mb

. 75€ NEBuffer 1.1 2.1 3.1 CutSmart 5,000 units 300 € % Activity 25 50 25

Activity at 37°C: 10%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

XhoI

#R0146S 5,000 units 71 € #R0146I 25,000 units 284 €

for high (5X) concentration

#R0146M 25,000 units 284 €

5′...C^TT C G A G ...3′ 3′...GAGCT₄C...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.



% Activity 75 100 100 100

Concentration: 20,000 and 100.000 units/ml Methylation Sensitivity: Cleavage

of mammalian genomic DNA is impaired by CpG methylation (see p. 334).

Tth 1111

#R0185S 400 units 60 €

5′...GACN<mark>'</mark>NNGTC...3′ 3′...C T G N N N C A G ... 5′

Reaction Conditions: CutSmart

Buffer, 65°C

CutSmart Ril O dil B 65° Mil

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 25 100

Concentration: 5,000 units/ml

Activity at 37°C: 10%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

XmaI

#R0180S 500 units 71 € #R0180L 2,500 units 284 €

for high (5X) concentration

#R0180M 2,500 units 284 €

5′...C^TCCGGG...3′ 3′...GGGCC₄C...5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 50 <10 100

CutSmart RR dilA 37° CpG

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by CpG methylation (see p. 334).

Note: Star activity may result from a glycerol concentration of > 5%.































CutSmart Rill dil B 37° 86 CpG

#R0659L 1,000 units 264

Concentration: 20,000 units/ml

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

#R0659L 1,000 units 264

5′... G A C G T C ... 3′
3′... C T G C A G ... 5′

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C

ZraI

for 20 minutes.

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Nicking Endonucleases

Reaction Conditions: CutSmart

for 20 minutes.

Buffer, 37°C. Heat inactivation: 65°C

As a rule, when restriction endonucleases bind to their recognition sequences in DNA, they cleave both strands of the duplex at the same time. Two independent hydrolytic reactions proceed in parallel, most often driven by the presence of two catalytic sites within each enzyme, one for hydrolyzing each strand. One of the focuses of our research program is to engineer restriction enzymes so that they hydrolyze only one strand of the duplex, generating DNA molecules that are "nicked", rather than cleaved. These conventional nicks (3´-hydroxyl, 5´-phosphate) can serve as initiation points for a variety of further enzymatic reactions such as replacement DNA synthesis, strand-displacement amplification (1), exonucleolytic degradation or the creation of small gaps (2).

Nicking endonucleases (NEases) are as simple to use as restriction endonucleases. Since the nicks generated by 6- or 7-base nicking endonucleases do not fragment DNA, their activities are monitored by conversion of supercoiled plasmids to open circles. Alternatively, substrates with nicking sites close enough on opposite strands to create a double-stranded cut can be used instead

The uses of nicking endonucleases are still being explored. NEases can generate nicked or gapped duplex DNA for DNA mismatch repair studies and for diagnostic applications. The long overhangs that nicking enzymes make can be used in DNA fragment assembly. Nt.BbvCl has been used to generate long and non-complementary overhangs when used with Xbal in the USER® cloning protocol from NEB. Nicking endonucleases are also useful for isothermal DNA amplification, which relies on the production of site-specific nicks. For example, isothermal DNA amplification using Nt.BstNBI in concert with Vent® (exo-) DNA Polymerase (NEB #M0257) (EXPAR) has been reported for detection of a specific DNA sequence in a sample (3). Another isothermal DNA amplification technique [Nicking Endonuclease Mediated- DNA Amplification (NEMDA)] (4) has been described using the 3-base cutter Nt.CviPII and Bst DNA Polymerase I. Frequent cutting NEases can generate short partial duplex DNA fragments from genomic DNA. These fragments can be used for cloning or used as probes for hybridization-based applications. Nicking enzymes have also been used for genome mapping.


NEB continues to engineer more nicking enzymes, particularly in response to specific customer needs and applications.

References:

for 20 minutes.

- (1) Walker, G.T. et al. (1992) *Proc. Natl. Acad. Sci. USA*, 89, 392–396.
- (2) Wang, H. and Hays, J.B. (2000) Mol. Biotechnol., 15, 97-104.
- (3) Van Ness, J. et al. (2003) Proc. Natl. Acad. Sci. USA, 89, 4504-4509.
- (2) Chan, S.H. et al. (2004) Nucl. Acids Res., 32, 6187-6199.

Nb.Bb	vCI	CutSmart Rill dii A 37° Killy					
#R0631S #R0631L	1,000 units 71 € 5,000 units 284 €	NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 100 100					
5′ C C T C		Concentration: 10,000 units/ml					
	onditions: CutSmart Heat inactivation: 80°C	Methylation Sensitivity: Not sensitive to <i>dam, dcm</i> or mammalian CpG methylation.					



Nb.BsrDI

5,000 units 284 €

CutSmart Rid dil A 65° 800 1,000 units 71 €

5′... G C A A T G N N ... 3′ 3'... CGTTACNN...5'

#R0648S

#R0648L

Reaction Conditions: CutSmart Buffer, 65°C. Heat inactivation: 80°C for 20 minutes

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 25 100 100

Concentration: 10,000 units/ml

Activity at 37°C: 75%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Nt.BbvCI

#R0632S 1.000 units 71 € #R0632L 5,000 units 284 €

5'... C CTCAGC...3' 3'... GGAGTCG...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 50 100 10 100

CutSmart Ri dil A 37° 80 CpG

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 334).

Nb.BssSI

#R0681T

#R0681S 1,000 units 72 € for high (5X) concentration

5,000 units 288 €

5'... CACGAG...3' 3'... GTGCT_C...5'

Reaction Conditions: NEBuffer 3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

RX NEB 3.1 dil B 37° W/6

Concentration: 20,000 and 100,000 units/ml.

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity <10 100 100

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Nt.BsmAI

#R0121S 500 units 72 €

5'... GTCTCN N ... 3' 3'... CAGAGNN...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart Ri dil A 37° CpG NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 100 50 10 100 Concentration: 5,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Nb.BtsI

#R0707S 1,000 units 72 €

5'... G C A G T G N N ... 3' 3'... CGTCAC,NN...5'

Reaction Conditions: CutSmart Buffer, 37°C

CutSmart Ri dil A 37° Kill

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 75 100 75 100

Concentration: 10,000 units/ml Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Nt.BspQI

#R0644S 1,000 units 72 €

5′...GCTCTTCN...3′ 3'...CGAGAAGN...5'

Reaction Conditions: NEBuffer 3.1. 50°C. Heat inactivation: 80°C for 20 minutes

RN NEB 3.1 dil B 50° WW NEBuffer 1.1 2.1 3.1 CutSmart

Concentration: 10,000 units/ml

Activity at 37°C: 80%

% Activity <10 25 100

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.

Nt.AlwI

#R0627S 500 units 71 €

5'...GGATCNNNNNN...3' 3'...CCTAGNNNNN...5'

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes

CutSmart Ri dil A 37° dam

Concentration: 10,000 units/ml

NEBuffer 1.1 2.1 3.1 CutSmart

% Activity 10 100 100 100

Methylation Sensitivity: Cleavage is blocked by dam methylation (see p. 334).

Nt.BstNBI

#R0607S 1,000 units 71 € #R0607L 5,000 units 284 €

5'...GAGTCNNNNN...3' 3'...CTCAGNNNNN...5'

Reaction Conditions: NEBuffer 3.1, 55°C. Heat inactivation: 80°C for 20 minutes.

RX NEB3.1 dil A 55° 1866

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 0 10 **100**

Concentration: 10,000 units/ml Activity at 37°C: 10%

Methylation Sensitivity: Not sensitive to dam, dcm or mammalian CpG methylation.



























Nt.CviPII

CutSmart Ril dil A 37° 65 CpG

#R0626S 40 units 130 €

5′...*CCD... 3′ 3'... GGH ... 5'

D = A or G or T (not C)H = A or C or T (not G)

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 100 25 100

Concentration: 2,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).

Homing Endonucleases

Homing endonucleases are double stranded DNases that have large, asymmetric recognition sites (12-40 base pairs) and coding sequences that are usually embedded in either introns or inteins (1). Introns are spliced out of precursor RNAs, while inteins are spliced out of precursor proteins (2,3). Homing endonucleases are named using conventions similar to those of restriction endonucleases with intron-encoded endonucleases containing the prefix, "I-" and intein endonucleases containing the prefix, "PI-"(1,7).

Homing endonuclease recognition sites are extremely rare. For example, an 18 base pair recognition sequence will occur only once in every 7 x 10¹⁰ base pairs of random sequence. This is equivalent to only one site in 20 mammalian-sized genomes (4). However, unlike standard restriction endonucleases, homing endonucleases tolerate some sequence degeneracy within their recognition sequence (5,6). As a result, their observed sequence specificity is typically in the range of 10–12 base pairs.

Homing endonucleases do not have stringently-defined recognition sequences in the way that restriction enzymes do. That is, single base changes do not abolish cleavage but reduce its efficiency to variable extents. The precise boundary of required bases is generally not known. The recognition sequence listed is one site that is known to be recognized and cleaved.

References:

- (1) Belfort, M. and Roberts, R.J. (1997) Nucleic Acids Res., 25, 3379–3388.
- (2) Dujon, B. et al. (1989) Gene, 82, 115-118.
- (3) Perler, F.B. et al. (1994) Nucleic Acids Res., 22, 1125-1127.
- (4) Jasin, M. (1996) Trends in Genetics, 12, 224-228
- (5) Gimble, F.S. and Wang, J. (1996) J. Mol. Biol., 263, 163-180.
- (6) Argast, M.G. et al. (1998) J. Mol. Biol., 280, 345-353.
- (7) Roberts, R.J. et al. (2003) Nucleic Acids Res., 31, 1805-1812.

I-CeuI

CutSmart RN dil B 37° (55)

#R0699S 500 units 74 € #R0699L 2,500 units 296 €

(Supplied with 5 µg of plasmid DNA)

Reaction Conditions: CutSmart Buffer, 37°C, Heat inactivation: 65°C for 20 minutes.

Specificity: The homing or recognition site for this endonuclease is shown below:

5′...TAACTATAACGGTCCTAAGGTAGCGAA...3 3′...ATTGATATTGCCAGGATTCCATCGCTT...5

NEBuffer 1.1 2.1 3.1 CutSmart % Activity 10 10 <10 100

Concentration: 5,000 units/ml

Note: Homing endonucleases do not have stringently-defined recognition sequences.

I-SceI

#R0694S 500 units 71 € #R0694L 2,500 units 284 €

(Supplied with 5 µg of plasmid DNA)

Reaction Conditions: CutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Specificity: The homing or recognition site for this endonuclease is shown below:

5′...TAGGGATAACAGGGTAATA...3′ 3'...ATCCCTATTGTCCCATTAT...5

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 50 25 100

Concentration: 5,000 units/ml

CutSmart RX dil B 37° 655

Note: Homing endonucleases do not have stringently-defined recognition sequences.

PI-PspI

RX NEBU BSA dil B 65° Wh

#R0695S 500 units 70 € (Supplied with 5 µg of plasmid DNA)

Reaction Conditions: NEBuffer PI-PspI + BSA, 65°C

Specificity: The homing or recognition site for this endonuclease is shown

5'...TGGCAAACAGCTATTATGGGTATTATGGGT...3
3'...ACCGTTTGTCGATAATACCCATAATACCCA...5

NEBuffer 1.1 2.1 3.1 CutSmart 10 10 % Activity <10

Concentration: 5,000 units/ml

Activity at 37°C: 5%

Note: Homing endonucleases do not have stringently-defined recognition

sequences.

PI-SceI

#R0696S 250 units 71 € #R0696L 1,250 units 284 €

(Supplied with 5 µg of plasmid DNA)

Reaction Conditions: NEBuffer PI-Scel + BSA, 37°C. Heat inactivation: 65°C for 20 minutes.

Specificity: The homing or recognition site for this endonuclease is shown

below:

5'...ATCTATGTCGGGTGCGGAGAAAGAGGTAATGAAATGG...3
3'...TAGATACAGCCCACGCCTCTTTCTCCATTACTTTACC...5

NEBuffer 1.1 2.1 3.1 CutSmart % Activity <10 <10 <10 <10

Concentration: 5,000 units/ml

Note: Homing endonucleases do not have stringently-defined recognition

RN NEBU BSA dil B 37° Kith

sequences.

Recombinant Albumin, Molecular Biology Grade

#B9200S

12 mg 42 €

Companion Product:


BSA, Molecular Biology Grade #R9000S 12 mg 30 € Reaction Conditions: Recombinant Albumin, Molecular Biology Grade, is a non-bovine derived albumin that can serve as an alternative to Bovine Serum Albumin (BSA). Like BSA, it has been shown to prevent adhesion of enzymes to reaction tubes and pipette surfaces. It also stabilizes some proteins during incubation. Choose Recombinant Albumin, when there is a need to avoid BSA sequences.

NEB Tube Opener

#C1008S

2 units 20 €

Description: Use to open a variety of microcentrifuge tubes. Can be used for snap-on caps or screw-on caps.



Reaction Buffers

NEBuffer™ 1.1 (10X)

#B7201S 5.0 ml 23 €

NEBuffer 2.1 (10X)

#B7202S 5.0 ml 23 €

NEBuffer 3.1 (10X)

5.0 ml 23 € #B7203S

CutSmart Buffer (10X)

#B7204S 5.0 ml 23 €

NEBuffer Set (EcoRI, DpnII) #B7006S 5.0 ml 23 € S-adenosylmethionine (SAM) #B9003S (32 mM) 0.5 ml 38 €

Nuclease-free Water

25 ml 25 € #B1500S 100 ml 61 € #B1500L

Description: New England Biolabs provides a color-coded 10X NEBuffer with each restriction endonuclease to ensure optimal (100%) activity. Most of our enzymes are supplied with one of four standard NEBuffers. Occasionally, an enzyme has specific buffer requirements not met by one of the four standard NEBuffers, in which case the enzyme is supplied with its own unique NEBuffer.

Nuclease-free Water is ideal for the preparation of reagents and for use in enzymatic reactions. No toxic agents, such as DEPC, are used in the production of this water, so as to avoid inhibition in enzymatic reactions.

For more information about NEBuffer compatibility, please turn to page 293.

Reaction Buffer Compositions: Visit www.neb.com for details.























Diluent Buffers

Diluent A #B8001S	5.0 ml 32 €
Diluent B #B8002S	5.0 ml 32 €
Diluent C #B8003S	5.0 ml 32 €

Description: Diluent Buffers are recommended for making dilutions of restriction endonucleases. When necessary, we recommend diluting enzymes just prior to use and suggest that the final concentration of diluted enzymes be at least 1,000 units/ml. Diluent preference for each restriction endonuclease is listed with its catalog entry; a complete listing of all restriction endonucleases and their diluents can be found on page 293.

Storage Conditions:

Store at -20°C.

Diluent Buffer Compositions: Visit www.neb.com for details.

Gel Loading Dyes

Gel Loading Dye, Bli	ue (6X)
#B7021S	4 ml 48 €
Gel Loading Dye, Ora	ange (6X)
#B7022S	4 ml 48 €
Gel Loading Dye, Pu	ırple (6X)
#B7024S	4 ml 43 €
Gel Loading Dye, Pu	rple (6X), no SDS
#B7025S	4 ml 43 €

Description: Pre-mixed loading dye solutions are available with a choice of tracking dyes, for agarose and non-denaturing polyacrylamide gel electrophoresis. Three of the solutions contain SDS, which often results in sharper bands, as some restriction enzymes remain bound to their DNA substrates following cleavage. Some interference may be observed when using SYBR® or GelRed® as precast dyes in the presence of increased concentrations of SDS. When using these dyes as precast dyes, NEB recommends using our Gel Loading Dye, Purple (6X), no SDS (NEB #B7025). Gel Loading Dyes contain EDTA to chelate up to 10 mM magnesium, thereby stopping the reaction. Bromophenol blue is the standard tracking dye for electrophoresis. It migrates at approximately 300 bp on a standard 1% TBE agarose gel. Orange G will not appear in gel photographs, and runs ahead of all but the smallest restriction fragments, migrating at approximately 50 bp on a standard 1% TBE agarose gel. NEB also offers a unique purple dye which migrates similarly to Bromophenol blue. Specifically chosen, this dye does not leave a shadow under UV light.

Gel Loading Dye Compositions: Visit www.neb.com for details.

Note: Use 5 μ I of Gel Loading Dye per 25 μ I reaction, or 10 μ I per 50 μ I reaction. Mix well before loading gel. Store at room temperature.

SYBR® is a registered trademark of Molecular Probes, Inc GELRED® is a registered trademark of Biotium.

Twice a year, we bring together NEB team members from around the world to learn about the latest NEB innovations. U.S. Sales Representatives Corina, Faraz, Lynn, and the General Manager of NEB France, Eric Beguec, join NEB's CEO Jim Ellard (center), as he welcomes them to NEB's Ipswich, Massachusetts headquarters.







Habitat Fragmentation

Habitat fragmentation occurs when a continuous habitat is disrupted and divided into smaller habitats that are physically isolated from each other. It carves the environment into patches and divides plant and animal populations. Fragmentation can occur when a forest is disrupted for urban development or agricultural purposes, following a forest fire, or by the addition of roadways through a previously continuous habitat. The result is that the number of plants and animals in the fragments decreases, and the further apart the fragments are, the greater the loss in species diversity and genetic diversity. This increases the likelihood of species becoming extinct, as they are less able to deal with environmental changes or disturbances.

Furthermore, "edge effects" are observed in fragmented habitats. Edge effects refer to the contrast between adjacent habitats that can hinder the safe travel of species. When the contrast is great, for example, with suburban development, abiotic factors, such as an increase in sunlight, wind and temperature are observed. These factors influence which plants and animals live near the edge. A sharp contrast between habitats affects the safety of animals as they move around to forage for food and water. It also makes seasonal migration more dangerous. If an animal has to leave its core habitat and enter a more hostile environment, they become more vulnerable to attack by a predator.

In a fragmented forest, such as the Middle Magdalena Valley in Colombia, the Brown Spider Monkey no longer has a continuous canopy to move around — 80% of its habitat has been cleared for farmland and cattle ranches. Spider monkeys have a prehensile tail that acts like a fifth limb and gives them the appearance of a spider; they require large tracts of undisturbed forest due to their large body size.

While the monkeys rely on the forest, the forest also relies on the monkeys to regenerate. Spider monkeys disperse the seeds of approximately 100 tropical plants through ingestion of the fruit in one location, and then elimination of the seeds in a different part of the forest. Spider monkeys defecate 13 to 17 times a day, and this represents the dispersal of hundreds to thousands of seeds. Many of the seeds that they disperse belong to the dense, hardwood trees that are the most effective in eliminating carbon dioxide from the atmosphere.

Habitat fragmentation limits the distance that can be traveled safely by the spider monkey, and consequently, the area of seed dispersal. Fragmentation is causing a decline in the spider monkey population, which causes concern regarding the future of the forest. In an effort to prevent this, researchers observe the monkeys and are directed to their feces by dung beetles. They collect, wash and categorize the seeds in order to predict which species of plants will be dispersed, the quantity of each species, and to see if there is a pattern between the type of seeds dispersed and the monkey's gender or age. They then plant the seeds in degraded areas that could potentially provide continuity, a nature corridor, between fragmented areas.

It is hoped that efforts to restore the continuity of the Middle Magdalena Valley will allow spider monkeys to again move freely around the canopy and carry out their crucial role in preserving the biodiversity of the forest. The tale of the spider monkey demonstrates how habitat fragmentation effects not just a single species, but an entire ecosystem.

DNA Polymerases & Amplification Technologies



NEB has pursued the discovery & development of DNA polymerases for over 25 years.

As the first company to sell *Taq* DNA polymerase and the first to discover a PCR-stable high-fidelity DNA polymerase, NEB established a foundation in developing innovative, high quality tools for the research and diagnostic community.

PCR & aPCR

NEB's product portfolio features a large selection of polymerases for PCR. Q5 High-Fidelity DNA Polymerase offers fidelity 280 times higher than *Taq*, along with superior performance with minimal optimization. One *Taq* is ideal for robust amplification in routine PCR applications. Both are available in hot start formulations, utilizing a novel aptamer-based approach that does not require an activation step.

Fluorescence-based quantitative real-time PCR (qPCR) is the gold standard for detection and quantitation of nucleic acids due to its specificity and sensitivity. Luna qPCR & RT-qPCR products feature an inert blue tracking dye for easy reaction setup and are available for intercalating dye or fluorescent probe-based detection methods.

Isothermal Amplification

Sequence specific isothermal amplification approaches eliminate the need for temperature cycling, providing advantages for certain point-of-care diagnostic needs. NEB's broad suite of reagents continue to enable advancement in isothermal amplification. Nicking enzymes, WarmStart enzymes, strand-displacing DNA polymerases, and RNA polymerases offer flexibility to assemble and design an isothermal amplification platform.

Featured Products

Q5® High-Fidelity DNA Polymerase

One *Taq*® DNA Polymerase

69 Luna® qPCR & RT-qPCR Products

74 WarmStart® LAMP Kit (DNA & RNA)

Featured Tools & Resources

62 PCR Polymerase Selection Chart

337 Guidelines for PCR Optimization

PCR Troubleshooting



Visit www.neb.com/PCR to find additional online tools, video tech tips and tutorials to help you in setting up your PCR experiments.





	DOD Deliverage Calestian Chart			DT DOD	
	PCR Polymerase Selection Chart	62		RT-PCR	
	High Fidelity PCR			One <i>Taq</i> One-Step RT-PCR Kit One <i>Taq</i> RT-PCR Kit	72 72
RR	Q5 High-Fidelity DNA Polymerase	63	m	Olie lay NI-FON KIL	12
RR	Q5 High-Fidelity 2X Master Mix	63		Isothermal Amplification	
RR	Q5 Hot Start High-Fidelity DNA Polymerase	63		and Strand Displacement	
RR	Q5 Hot Start High-Fidelity 2X Master Mix	63	RR	Bst DNA Polymerase, (Large Frag. and Full-Length)	73
RR	Q5U Hot Start High-Fidelity DNA Polymerase	63	RR	Bst 2.0 DNA Polymerase	73
RR	Q5 High-Fidelity PCR Kit	63	RR	Bst 2.0 WarmStart DNA Polymerase	73
	Q5 Site-Directed Mutagenesis Kit	92	RR	Bst 3.0 DNA Polymerase	73
RR	Q5 Site-Directed Mutagenesis Kit	00		WarmStart Colorimetric LAMP	
	(without competent cells)	92		Master Mix (DNA & RNA)	74
	Phusion High-Fidelity DNA Polymerase	64		WarmStart LAMP Kit (DNA & RNA)	74
	Phusion High-Fidelity PCR Master Mix w/HF Buffer	64		IsoAmp II Universal tHDA Kit	75
	Phusion High-Fidelity PCR Master Mix w/GC Buffer Phusion Hot Start Flex DNA Polymerase	64 64	RR	phi29 DNA Polymerase	75
	Phusion Hot Start Flex 2X Master Mix	64		Polymerases for DNA Manipulation	
	Phusion High-Fidelity PCR Kit	64		(labeling, blunting, etc.)	
11111	Triusion riigh riddity r on tat	U -1		PreCR Repair Mix	76
	Routine PCR		RR	Sulfolobus DNA Polymerase IV	76
RR	One <i>Taq</i> DNA Polymerase (Std., Quick-Load)		RR	Therminator DNA Polymerase	76
	One Taq Master Mixes (Std., Quick-Load)	65	RR	DNA Polymerase I (E. coli)	77
	One <i>Taq</i> Hot Start DNA Polymerase	65		DNA Polymerase I, Large (Klenow) Fragment	77
	One Taq Hot Start Master Mixes (Std., Quick-Load)		RR	Klenow Fragment (3´→5´ exo-)	77
	Taq DNA Polymerase w/ThermoPol Buffer	66	RR	T4 DNA Polymerase	77
	Taq DNA Polymerase w/Standard Taq Buffer	66		T7 DNA Polymerase (unmodified)	78
	Taq DNA Polymerase w/Standard Taq (Mg-free) Buffer Taq Master Mixes (2X, 5X, Quick-Load)	66		Bsu DNA Polymerase, Large Fragment	78
	Tag PCR Kit	66	RR	Terminal Transferase	78
	Multiplex PCR 5X Master Mix	66		Nucleotide Solutions & Buffers	
	Hot Start <i>Tag</i> DNA Polymerase	66		Reaction Buffers and Diluents	79
	Hot Start <i>Tag</i> 2X Master Mix	66		Acyclonucleotide Set	79
				Deoxynucleotide (dNTP) Solution Mix/Set	79
	Specialty PCR			Ribonucleotide Solution Mix/Set	79
	LongAmp Taq DNA Polymerase	67		7-deaza-dGTP	79
	LongAmp Hot Start <i>Taq</i> DNA Polymerase	67		Adenosine-5´ Triphosphate (ATP)	79
	LongAmp Taq 2X Master Mix	67		5-methyl-dCTP	79
	LongAmp Hot Start <i>Taq</i> 2X Master Mix	67		dATP Solution	79
	LongAmp <i>Taq</i> PCR Kit	67		cDNA Synthesis	
	Hemo KlenTaq	67		cDNA Synthesis Selection Chart	80
Kii	EpiMark Hot Start <i>Taq</i> DNA Polymerase	68	RR	LunaScript RT SuperMix Kit	70
	Other PCR Polymerases		RR	ProtoScript II First Strand cDNA Synthesis Kit	194
RR	Vent DNA Polymerase	68	RR	ProtoScript First Strand cDNA Synthesis Kit	194
RR	Vent (exo-) DNA Polymerase	68	RR	Template Switching RT Enzyme Mix	193
RR	Deep Vent DNA Polymerase	68	RR	ProtoScript II Reverse Transcriptase	192
RR	Deep Vent (exo-) DNA Polymerase	68	RR	M-MuLV Reverse Transcriptase	192
	aPCR & RT-aPCR				192
WŁ	nat is Difference Between One-step		RR	WarmStart RTx Reverse Transcriptase	193
	d Two-step RT-PCR?	69		PCR Cleanup	
	Luna qPCR & RT-qPCR	69		Monarch PCR & DNA Cleanup Kit	81
	Luna Universal qPCR Master Mix	70		Exo-CIP Rapid PCR Kit	81
RR	Luna Universal Probe qPCR Master Mix	70			
RR	LunaScript RT SuperMix Kit	70		Recombinant Enzy	/me
	Luna Universal One-Step RT-qPCR Kit	71		- In recombinant Enzy	HIC
RR	Luna Universal Probe One-Step RT-qPCR Kit	71			

One or more of these products are covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at gbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications. Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at https://www.neb.com/support/terms-of-sale. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

PCR Polymerase Selection Chart

For over 40 years, New England Biolabs, Inc. has been a world leader in the discovery and production of reagents for the life science industry. NEB offers a wide range of DNA polymerases, and through our commitment to research, ensures the development of innovative and high quality tools for PCR and related applications. The following table simplifies the selection of a polymerase that best suits your experiment.

High-fidelity polymerases benefit from a Tm+3 annealing temp. Use the **NEB Tm Calculator** to ensure successful PCR at TmCalculator.neb.com

★ indicates recommended choice for application

A LIE	STAND	ARD PCR	HIC	GH-FIDELITY	PCR	SPECIALTY PCR				
			HIGHEST Fidelity		MODERATE FIDELITY	LONG Amplicons		U Rance	BLOOD DIRECT PCR	
	One <i>Taq®/</i> One <i>Taq</i> Hot Start	<i>Taq </i> Hot Start <i>Taq</i>	Q5 [®] /Q5 Hot Start	Phusion ^{®(4)} / Phusion ⁽⁴⁾ Flex	Vent [®] / Deep Vent™	LongAmp [®] / LongAmp Hot Start <i>Taq</i>	Epimark [®] Hot Start <i>Taq</i>	Q5U™	Hemo KlenTaq®	
PROPERTIES										
Fidelity vs. <i>Taq</i>	2X	1X	~280X ⁽²⁾	> 39X	5-6X	2X	1X	ND	ND	
Amplicon Size	< 6 kb	≤ 5 kb	$\leq 20 \text{ kb}$	≤ 20 kb	≤ 6 kb	≤ 30 kb	≤ 1 kb	app-specific	≤ 2 kb	
Extension Time	1 kb/min	1 kb/min	6 kb/min	4 kb/min	1 kb/min	1.2 kb/min	1 kb/min	2 kb/min	0.5 kb/min	
Resulting Ends	3´ A/Blunt	3´ A	Blunt	Blunt	Blunt	3´ A/Blunt	3´ A	Blunt	3´ A	
3'→5' exo	Yes	No	Yes	Yes	Yes	Yes	No	Yes	No	
5'→3' exo	Yes	Yes	No	No	No	Yes	Yes	No	No	
Units/50 µl Reaction	1.25	1.25	1.0	1.0	0.5-1.0	5.0	1.25	1.0	N/A	
Annealing Temperature	Tm-5	Tm-5	Tm+3	Tm+3	Tm-5	Tm-5	Tm-5	Tm-3	Tm-5	
APPLICATIONS										

APPLICATIONS	APPLICATIONS								
Routine PCR	*	•	•	•	•	•			
Colony PCR	*	•							
Enhanced Fidelity	•		*	•	•	•			
High Fidelity			*	•					
High Yield	*	•	*	•					
Fast			*	•					
Long Amplicon			*	•		*			
GC-rich Targets	*		*		•	•			
AT-rich Targets	*	•	*	•		•	•	*	
High Throughput	•	•	•	•			•	*	
Multiplex PCR	•	★ (1)	•	•					
Extraction-free PCR									*
DNA Labeling		*							
Site-directed Mutagenesis			*	•					
Carryover Prevention							•	*	
USER® Cloning							•	*	

NGS APPLICATIONS										
	NGS Library Amplification			★(3)	•			★(5)		

FORMATS							
Hot Start Available	•	•	•	•	•	•	
Kit		•	•	•	•		
Master Mix Available	•	•	•	•	•		
Direct Gel Loading	•	•					

- (1) Use Multiplex PCR 5X Master Mix.
- (2) Due to the very low frequency of misincorporation events being measured, the error rate of high-fidelity enzymes like Q5 is challenging to measure in a statistically significant manner. We continue to investigate improved assays to characterize Q5's very low error rate to ensure that we present the most robust accurate fidelity data possible (Popatov, V. and Ong, J.L. (2017) PLoS One, 12(1):e0169774. doi 10.1371/journal. pone. 0169774).
- (3) Use NEBNext High-Fidelity 2X PCR Master Mix.
- (4) Phusion DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific.
- (5) Use NEBNext Enzymatic Methyl-seg Kit (EM-seg™) for Illumina.



Why choose Q5 for your PCR?















Q5® High-Fidelity DNA Polymerase

RX NEBU HIFI PCR 100 Tm+3

Q5 Hot Start High-Fidelity DNA Polymerase

Q5 High-Fidelity DNA	Polymerase
#M0491S	100 units 96 €
#M0491L	500 units 424 €
	Master Mix s (50 µl vol) 157 € s (50 µl vol) 628 €
Q5 Hot Start High-Fid	lelity DNA Polymerase
#M0493S	100 units 112 €
#M0493L	500 units 508 €
#M0494L 500 rxns (1	lelity 2X Master Mix 2 x 1.25 ml) 180 € 0 x 1.25 ml) 720 € 1 x 12.5 ml) 712 €

NOW AVAILABLE: NEBNext Ultra II Q5 Master Mix – see page 155.

Q5 POLYMERASE DETAILS	
Extension Rate	6 kb/min
Amplicon Size	≤ 20 kb
Fidelity	~ 280X <i>Taq</i>
Units / 50 µl rxn	1 unit
Resulting Ends	Blunt
3´→5´ Exonuclease Activity	Yes
5´→3´ Exonuclease Activity	No
Supplied Buffer	Q5 Reaction Buffer
Supplied Enhancer	Q5 High GC Enhancer
PRODUCT FORMATS	
Hot Start Available	Yes
Activation Required	No
Master Mix Available	Yes
PCR Kit Available	Yes
NGS Version Available	Yes
APPLICATIONS	
High-Fidelity PCR	Yes
Difficult PCR	Yes
High GC PCR	Yes
T/A, U/A Cloning	No
Colony PCR	No
Blunt Cloning	Yes
Multiplex PCR	Yes
USER Cloning	Yes (Q5U)
Carryover Prevention	Yes (Q5U)

Tm Calculator

High-Fidelity polymerases benefit from a Tm*3 annealing temp. Use the **NEB Tm Calculator** to ensure successful PCR at **TmCalculator.neb.com**

NEW

Q5U[™] Hot Start High-Fidelity DNA Polymerase #M0515S 100 units 152 € #M0515L 500 units 640 €

Q5 High-Fidelity PCR Kit

#E0555S 50 rxns (50 µl vol)92 € #E0555L 200 rxns (50 µl vol)312 €

Q5 Site-Directed Mutagenesis Kit

#E0554S 10 rxns 192€

Q5 Site-Directed Mutagenesis Kit (without competent cells)

#E0552S 10 rxns 132 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

For a complete listing of Reaction Buffers, see pages 79.

Description: Q5 High-Fidelity DNA Polymerase sets a new standard for both fidelity and performance. With the highest fidelity amplification available (~280 times higher than *Taq*), Q5 DNA Polymerase results in ultra-low error rates. Q5 DNA Polymerase is composed of a novel polymerase that is fused to the processivity-enhancing Sso7d DNA binding domain, improving speed, fidelity and reliability of performance.

The Q5 buffer system is designed to provide superior performance with minimal optimization across a broad range of amplicons, regardless of GC content. For routine or complex amplicons up to ~65% GC content, Q5 Reaction Buffer provides reliable and robust amplification. For amplicons with high GC content (>65% GC), addition of the Q5 High GC Enhancer ensures continued maximum performance.



Q5 Hot Start DNA Polymerase: In contrast to chemically-modified or antibody-based hot start polymerases, NEB's Q5 Hot Start utilizes a unique synthetic aptamer. This structure binds to the polymerase through non-covalent interactions, blocking activity during the reaction setup. The polymerase is activated during normal cycling conditions, allowing reactions to be set up at room temperature. Q5 Hot Start does not require a separate high temperature activation step, reducing the potential for sample degredation, shortening reaction times and increasing ease-of-use. Q5 Hot Start is an ideal choice for high specificity amplification and provides robust amplification of a wide variety of amplicons, regardless of GC content.


Q5U Hot Start High-Fidelity DNA Polymerase: A modified version of Q5 High-Fidelity DNA Polymerase capable of incorporating dUTP for carryover prevention. Q5U is also compatible with USER cloning methods and enables the amplification of bisulfite treated/deaminated DNA.

Additional Formats: For added convenience, Q5 DNA Polymerase is available in master mix format or as a kit. Master mix formulations include dNTPs, Mg+ and all necessary buffer components. The Q5 High-Fidelity PCR Kit contains the Q5 High-Fidelity 2X Master Mix, nuclease-free water and the Quick-Load Purple 1 kb Plus DNA Ladder. For information on the Q5 Site-Directed Mutagenesis Kit, with or without competent cells, see page 92.

Concentration: 2.000 units/ml


Visit www.Q5PCR.com for more information.

LABCHIP® is a registered trademark of Caliper Life Sciences, Inc.



Robust amplification with Q5 (A) and Q5 Hot Start (B) High-Fidelity DNA Polymerases, regardless of GC content.

Amplification of a variety of human genomic amplicons from low-to-high GC content using either Q5 or Q5 Hot Start High-Fidelity DNA Polymerase. Reactions using Q5 Hot Start were set up at room temperature. All reactions were conducted using 30 cycles of amplification. and visualized by microfluidic LabChip® analysis.



The stand-alone enzyme comes with a reaction buffer that supports robust amplification of high AT to routine targets. Addition of the High GC Enhancer allows amplification of GC rich and difficult targets. For added convenience, the master mix formats allow robust amplification of a broad range of targets with a single formulation.

Phusion® High-Fidelity DNA Polymerase

RX HiFi PCR WWw Tm+3

Phusion Hot Start Flex DNA Polymerase

Phusion High-Fidelity DNA Polymerase #M0530S 100 units98 € #M0530L 500 units 432 €

Phusion High-Fidelity PCR Master Mix with HF Buffer

#M0531S 100 rxns (50 µl vol) 164 € #M0531L 500 rxns (50 µl vol) 650 €

Phusion High-Fidelity PCR Master Mix with GC Buffer

#M0532S 100 rxns (50 µl vol) 164 € #M0532L 500 rxns (50 μl vol) 650 €

Phusion Hot Start Flex DNA Polymerase

#M0535S 100 units 114 € #M0535L 500 units 518 €

Phusion Hot Start Flex 2X Master Mix #M0536S 100 rxns (50 µl vol) 185 € #M0536L 500 rxns (50 μl vol) 735 €

Phusion High-Fidelity PCR Kit

#E0553S 50 rxns (50 μl vol) 68 € #E0553L 200 rxns (50 µl vol) 238 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

For a complete listing of Reaction Buffers, see pages 79.

PHUSION POLYMERASE D	ETAILS
Extension Rate	4 kb/min
Amplicon Size	≤ 20 kb
Fidelity	> 50X <i>Taq</i>
Units / 50 µl rxn	1 unit
Resulting Ends	Blunt
3´→5´ Exonuclease Activity	Yes
5´→3´ Exonuclease Activity	No
Supplied Buffers	- 5X Phusion HF Buffer - 5X Phusion GC Buffer
Supplied Enhancer	100% DMS0
PRODUCT FORMATS	
Hot Start Available Activation Required	Yes No
Master Mix Available	Yes
PCR Kit Available	Yes
APPLICATIONS	
High-Fidelity PCR	Yes
T/A, U/A Cloning	No
Colony PCR	No
Blunt Cloning	Yes

Description: DNA polymerases with high fidelity are important for applications in which the DNA sequence needs to be correct after amplification. Manufactured and quality controlled at New England Biolabs, Thermo Scientific® Phusion High-Fidelity DNA Polymerase offers both high fidelity and robust performance, and thus can be used for all PCR applications. Its unique structure, a novel Pvrococcus-like enzyme fused with a processivity-enhancing domain, increases fidelity and speed. Phusion Hot Start Flex DNA Polymerase is available as a standalone enzyme or in a master mix format, and enables high specificity amplification. Phusion DNA Polymerase is an ideal choice for cloning and can be used for long amplicons.

Additional Formats: Phusion and Phusion Hot Start Flex DNA Polymerases are also available in master mix format. Phusion master mixes are available with HF or GC Buffer.

RN HiFi PCR 100 MM Tm+3


The Phusion PCR Kit contains Phusion Polymerase, Phusion HF and GC buffers, deoxynucleotides, MgCl₂, DMSO and DNA size standards.

Concentration: 2.000 units/ml

* Phusion DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific. Phusion® is a registered trademark and property of Thermo Fisher Scientific.


Гm Calculator

High-Fidelity polymerases benefit from a Tm+3 annealing temp. Use the NEB Tm Calculator to ensure successful PCR at TmCalculator.neb.com



Phusion DNA Polymerase generates robust amplification even with short extension times.

A 1.2 kb C. elegans genomic amplicon was analyzed using different polymerases. Reactions were set up in duplicate and visualized by microfluidic LabChip analysis. All reactions were conducted according to manufacturer's instructions using 30 cycles and extension times of either 30 seconds or 2 minutes as indicated.



Phusion Hot Start Flex DNA Polymerase delivers robust amplification. All amplicons are from human Jurkat template except for the 1.3 kb C. elegans amplicon. Amplicon sizes are indicated above gel. All reactions were conducted in duplicate according to manufacturer's instructions using 30 cycles and visualized after microfluidic LabChip analysis.







One Tag® DNA Polymerase

RX NEBU PCR WW Tm-5

One Taq Hot Start DNA Polymerase

One?	Tan I	Prn	du	cts
Olic	ıay ı	10	uu	υιs

	One Tag DNA	Polymerase	
	#M0480S	200 units	38€
	#M0480L	1,000 units	152€
	#M0480X	5,000 units	600 €
	One Tag 2X M	laster Mix with Stan	dard Buffer
	#M0482S	100 rxns (50 µl vol)	39€
	#M0482L	500 rxns (50 µl vol)	156€
١	EW		
	One Taq Quick	k-Load® DNA Polym	nerase
	#M0509S		20€
	#M0509L		80 €
	#M0509X	2,500 units	340 €
	One Tag Quick	k-Load 2X Master M	lix with
	Standard Buf	fer	
	#M0486S	100 rxns (50 µl vol)	40€
	#M0486L	$500~\text{rxns}$ (50 μI vol)	160€
	One Tag RT-P	CR Kit	
	#E5310S		143€
	0 = 0	I DT DOD IVI	
	,	step RT-PCR Kit	4540
	#E5315S	30 rxns	154 €

ONE <i>Taq</i> POLYMERASE DE Extension Rate	1 kb/min
Amplicon Size	≤ 6 kb
Fidelity	> 2X Taq
Units / 50 µl rxn	1.25 units
Resulting Ends	3´ A/Blunt
3´→5´ Exonuclease Activity	Yes
5´→3´ Exonuclease Activity	
Supplied Buffers	- One <i>Taq</i> Std Rxn Buffer - One <i>Taq</i> GC Rxn Buffer
Supplied Enhancer	One <i>Taq</i> High GC Enhancer
PRODUCT FORMATS	
Hot Start Available Activation Required	Yes No
Master Mix Available	Yes
Direct Gel-loading Available	Yes
PCR Kit Available	Yes
APPLICATIONS	
Routine PCR	Yes
T/A, U/A Cloning	Yes
Colony PCR	Yes

One Taq Hot Start Products

•	
One <i>Taq</i> Hot #M0481S #M0481L #M0481X	Start DNA Polymerase 200 units78 € 1,000 units312 € 5,000 units1248 €
Standard Bu #M0484S	Start 2X Master Mix with uffer 100 rxns (50 µl vol)
One <i>Taq</i> Hot Buffer #M0485S #M0485L	Start 2X Master Mix with GC 100 rxns (50 µl vol)
with Standar #M0488S	Start Quick-Load 2X Master Mix rd Buffer 100 rxns (50 µl vol) 69 € 500 rxns (50 µl vol) 282 €
with GC Buf #M0489S	Start Quick-Load 2X Master Mix ffer 100 rxns (50 µl vol)
For a compl	ete listing of Deoxynucleotide

For a complete listing of Deoxynucleotide Solutions, see page 79.

For a complete listing of Reaction Buffers, see pages 79.

RR NEBU PCR WY WW Tm-5

Description: One *Taq* DNA Polymerase is an optimized blend of *Taq* and Deep Vent® DNA polymerases for use with routine and difficult PCR experiments. The 3′→ 5′ exonuclease activity of Deep Vent DNA Polymerase increases the fidelity and robust amplification of *Taq* DNA Polymerase. The One *Taq* reaction buffers and High GC Enhancer have been formulated for robust yields with minimal optimization, regardless of a template's GC content.

One *Taq* **Hot Start DNA Polymerase:** One *Taq* **Hot Start DNA** Polymerase utilizes an aptamer-based inhibitor. The hot start formulation combines convenience with decreased interference from primer-dimers and secondary products.

The aptamer-based inhibitor binds reversibly, blocking polymerase activity at temperatures below 45°C, allowing reactions to be set up at room temperature. One *Taq* Hot Start DNA Polymerase is activated during normal cycling conditions, eliminating the need for a separate high temperature incubation step to activate the enzyme. One *Taq* Hot Start DNA Polymerase can therefore be substituted into typical or existing *Taq*-based protocols.

One *Taq* and One *Taq* Hot Start are provided with Standard Reaction Buffer, GC Reaction Buffer and High GC Enhancer. Recommendations for buffer selection, based on % GC content, are shown below.


Quick-Load formats offer direct loading of PCR products onto gels, eliminating the need for a PCR clean-up step.

Additional Formats: For added convenience, One *Taq* and One *Taq* Hot Start DNA Polymerases are available in master mix format. Master mixes are available with Standard or GC Buffer. High GC Enhancer is also provided with master mixes containing GC Buffer. For direct gel loading, Quick-Load versions of master mixes are also available. For information on the One *Taq* RT-PCR Kit or One *Taq* One-step RT-PCR Kit, see page 72.

Concentration: 5,000 units/ml

One Taq Buffer Recommendations

AMPLICON % GC CONTENT	RECOMMENDED DEFAULT BUFFER	OPTIMIZATION NOTES
< 50% GC	One Taq Standard Reaction Buffer	Adjust annealing temperature, primer/ template concentration, etc., if needed.
50–65% GC	One Taq Standard Reaction Buffer	One Taq GC Reaction Buffer can be used to enhance performance of difficult amplicons.
> 65% GC	One Tag GC Reaction Buffer	One Taq GC Reaction Buffer with 10-20% One Taq High GC Enhancer can be used to enhance performance of difficult amplicons.



Amplification of a selection of sequences with varying GC content from human and C. elegans genomic DNA using OneTaq DNA Polymerase. GC content is indicated above gel. Marker M is the 1 kb DNA Ladder (NEB #N3232).

Taq DNA Polymerase

RX NEBU PCR WW Tm-5

Hot Start Tag DNA Polymerase

Tag Polymerases

Taq DNA Polymerase with				
ThermoPol®	Buffer			
#M0267S	400 units60 €			
#M0267L	2,000 units 240 €			
#M0267X	4,000 units 440 €			
#M0267E	20,000 units1200 €			
<i>Taq</i> DNA Pol	ymerase with Std <i>Taq</i> Buffer			
#M0273S	400 units60 €			
#M0273L	2,000 units 240 €			
#M0273X	4,000 units 440 €			
#M0273E	20,000 units1200 €			
Taq DNA Pol	ymerase with Std <i>Taq</i>			
(Mg-free) Bu	iffer			
#M0320S	400 units60 €			
#M0320L	2,000 units 240 €			
Tag PCR Kit				
#E5000S	200 rxns99 €			
Taq 2X Mast	er Mix			
#M0270L	500 rxns (50 μl vol) 141 €			

Taq DNA POLYMERASE DE	TAILS
Extension Rate	1 kb/min
Amplicon Size	≤ 5 kb
Units / 50 µl rxn	1.25 units
Resulting Ends	3´ A
3´→5´ Exonuclease Activity	No
5´→3´ Exonuclease Activity	Yes
Supplied Buffer (product dependent)	- Std <i>Taq</i> Rxn Buffer or - ThermoPol Rxn Buffer
PRODUCT FORMATS	
PRODUCT FORMATS Hot Start Available Activation Required	Yes No
Hot Start Available	1.00
Hot Start Available Activation Required	No
Hot Start Available Activation Required Master Mix Available	No Yes
Hot Start Available Activation Required Master Mix Available Direct Gel-loading Available	No Yes Yes
Hot Start Available Activation Required Master Mix Available Direct Gel-loading Available PCR Kit Available	No Yes Yes
Hot Start Available Activation Required Master Mix Available Direct Gel-loading Available PCR Kit Available APPLICATIONS	No Yes Yes Yes
Hot Start Available Activation Required Master Mix Available Direct Gel-loading Available PCR Kit Available APPLICATIONS Routine PCR	No Yes Yes Yes

	<i>Taq</i> 2X Master Mix 500 rxns (50 μl vol) 216 €	
<i>Taq</i> 5X Mast #M0285L	er Mix 500 rxns (50 µl vol) 141 €	
	CR 5X Master Mix 100 rxns (50 µl vol) 228 €	
Hot Start 7	aq Products	
Hot Start Tac	7 DNA Polymerase	
#M0495S	200 units66 €	
#M0495L	1,000 units 264 €	
Hot Start <i>Taq</i> 2X Master Mix		

For a complete listing of Deoxynucleotide Solutions, see page 79.

100 rxns (50 μl vol) 58 €

500 rxns (50 μl vol) 232 €

#M0496S

#M0496L

For a complete listing of Reaction Buffers, see pages 79.

RR NEBU PCR 165 YM6 Tm-5

Description: Taq DNA Polymerase is a thermostable DNA polymerase that possesses a $5 \rightarrow 3$ polymerase activity and a 5 flap endonuclease activity. It is the most widely used enzyme for PCR. To accomodate a variety of PCR applications, Tag is available with different reaction buffers. Standard Tag Reaction Buffer is designed to support existing PCR platforms, and is an ideal choice for DHPLC and high-throughput applications. ThermoPol Reaction Buffer was designed at NEB, and is formulated to promote high product yields, even under demanding conditions. For additional convenience, Taq DNA Polymerase is also available in kit and master mix formats. For direct gel loading, a Quick-Load version of the Taq 2X Master Mix is also available.

Hot Start Tag DNA Polymerase: With value pricing and attractive commercial terms, Hot Start Tag from NEB is an ideal choice for molecular diagnostics and other applications. In contrast to chemically modified or antibody-based hot start polymerases, NEB's Hot Start Taq utilizes an aptamer-based technology. The unique aptamer binds to the polymerase through non-covalent interactions, inhibiting polymerization at non-permissive temperatures. This novel method eliminates the need for an activation step, reducing the potential for sample degradation and decreasing overall reaction time.

Additional Formats: For added convenience, *Tag* and Hot Start Tag DNA Polymerase are available in master mix format. For direct gel loading, a Quick-Load version of the Tag 2X Master Mix is also available. The Tag PCR Kit contains Tag DNA Polymerase, dNTP Mix, Buffer, MgCl, and the Quick-Load Purple 1 kb Plus DNA Ladder.

The Multiplex PCR 5X Master Mix formulation has been specifically optimized for enhanced performance in multiplex PCR reactions.

Concentration: 5,000 units/ml

Taq Buffer Selection Chart

CHOICE OF BUFFER	AVAILABLE PRODUCTS	NEB#
ThermoPol Reaction Buffer: Designed for optimal yield and specificity	Taq DNA Polymerase with ThermoPol Buffer	M0267
Standard <i>Tag</i> Reaction Buffer:	Taq DNA Polymerase with Standard Taq Buffer	M0273
Detergent-free and designed to be compatible with existing assay systems	Taq DNA Polymerase with Standard Taq (Mg- free) Buffer	M0320











LongAmp® Taq DNA Polymerase

LongAmp Hot Start Taq DNA Polymerase

LongAmp Taq DNA Polymerase #M0323S 500 units79 € #M0323L 2,500 units316 € LongAmp Hot Start Taq DNA Polymerase #M0534S 500 units104 € #M0534L 2,500 units416 € LongAmp Taq 2X Master Mix

#M0287S 100 rxns (50 µl vol) 132 € #M0287L 500 rxns (50 µl vol) 528 €

Extension Rate	1.2 kb/min
Amplicon Size	≤ 30 kb
Fidelity	2X <i>Taq</i> DNA Polymerase
Units / 50 µl rxn	5 units
Resulting Ends	3´ A/Blunt
3´→5´ Exonuclease Activity	Yes
5´→3´ Exonuclease Activity	Yes
Supplied Buffer (product dependent)	LongAmp <i>Taq</i> Rxn Buffer
PRODUCT FORMATS	
Hot Start Available	Yes
Activation Required	No
Master Mix Available	Yes
Direct Gel-loading Available	No
PCR Kit Available	Yes
APPLICATIONS	
THE ELECTRICATION	
	Yes
Long Amplicons	Yes
Long Amplicons Routine PCR T/A, U/A Cloning	


LongAmp Hot Start *Taq* 2X Master Mix #M0533S 100 rxns (50 µl vol) 165 € #M0533L 500 rxns (50 µl vol) 660 €

LongAmp Tag PCR Kit

#E5200S 100 rxns (50 µl vol) 113 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

For a complete listing of Reaction Buffers, see pages 79.



Amplification of longer templates with LongAmpTaq. Amplification of specific sequences from human genomic DNA using LongAmp Taq DNA Polymerase. Amplicon sizes are indicated below gel. Ladder (M) is NEB 1 kb DNA Ladder (NEB #N3232).

RX NEBU PCR 1115 Tm-5

RR NEBU PCR 100 11m-5

Description: An optimized blend of *Taq* and Deep Vent DNA Polymerases, LongAmp *Taq* DNA Polymerase enables amplification of longer PCR products with higher fidelity than *Taq* DNA Polymerase alone.

LongAmp Hot Start Tag DNA Polymerase:

LongAmp Hot Start Taq DNA Polymerase utilizes a unique synthetic aptamer. This structure binds reversibly to the enzyme, inhibiting polymerase activity at temperatures below 45°C, but releases the enzyme during normal PCR cycling conditions.

Additional Formats: For added convenience, LongAmp Taq and LongAmp Hot Start Taq are available in master mix format. The LongAmp Taq PCR Kit includes LongAmp Taq DNA Polymerase (2,500 units/ml), dNTP Mix (10 mM), LongAmp Taq Reaction Buffer Pack (5X), MgSO, (100 mM) and nuclease-free water.

Concentration: 2,500 units/ml


Hemo KlenTaq®

#M0332S 200 rxns (25 µl vol) 107 € #M0332L 1,000 rxns (25 µl vol) 428 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

Extension Rate	0.5 kb/min
Amplicon Size	≤ 2 kb
Units / 50 µl rxn	4 units
Resulting Ends	3´ A
3´→5´ Exonuclease Activity	No
5´→3´ Exonuclease Activity	No
Supplied Buffer	Hemo KlenTaq Rxn Buffer
APPLICATIONS	
Extraction-free PCR	Yes
T/A, U/A Cloning	Yes
Colony PCR	Yes

Description: Hemo KlenTaq is a truncated version of *Taq* DNA Polymerase, lacking the first 280 amino acids. Hemo KlenTaq also contains mutations that make it resistant to inhibitors present in whole blood. It can amplify DNA from whole blood samples from humans and mice, without the need for DNA extraction. Hemo KlenTaq tolerates up to 20% whole blood in a 25 μ l reaction (30% in a 50 μ l reaction). Hemo KlenTaq works well with most common anticoagulants, including heparin, citrate and EDTA.



RN NEBU → MM PCR Tm⁻5

Source: An *E. coli* strain that carries a mutant *Taq* DNA polymerase gene. The protein lacks the N-terminal 5´ flap endonuclease domain and the gene has three internal point mutations.

Reaction Buffer: 1X Hemo KlenTaq Reaction Buffer

KLENTAQ® is a registered trademark of Wayne M.Barnes.

Amplification of human whole blood with Hemo KlenTaq. Lane 1: 5% blood + Na-EDTA; Lane 2: 5% blood + K-EDTA; Lane 3: 5% blood + Na-Heparin; Lane 4: 5% blood + Na-Citrate; Lane 5: 1.2 mm² FTA Guthrie Card containing dried human blood + Na-Heparin; Lane 6: 1.2 mm² PTA Guthrie Card containing dried human blood + Na-Heparin (washed with 50 µl H₂O at 50°C for 5 min.). Ladder (M) is the 1 kb Plus DNA Ladder (NEB #N3200).

EpiMark® Hot Start Taq DNA Polymerase

#M0490S 100 rxns (50 µl vol) 61 € #M0490L 500 rxns (50 μl vol) 244 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

EPIMARK POLYMERASE D	ETAILS
Extension Rate	1 kb/min
Amplicon Size	≤ 1 kb
Units / 50 µl rxn	1.25 units
Resulting Ends	3´A
3´→5´ Exonuclease Activity	No
5´→3´ Exonuclease Activity	Yes
Supplied Buffer	EpiMark Hot Start Taq Rxn Buffer
APPLICATIONS	
AT-rich Targets	Yes
Bisulfite-converted DNA	Yes
Routine PCR	Yes
T/A, U/A Cloning	Yes

Description: EpiMark Hot Start *Tag* DNA Polymerase is an excellent choice for use on bisulfite-converted DNA. It is a mixture of Taq DNA Polymerase and a temperaturesensitive, aptamer-based inhibitor. The inhibitor binds reversibly to the enzyme, inhibiting polymerase activity at temperatures below 45°C, but releases the enzyme during normal PCR cycling conditions. This permits room temperature reaction assembly with no separate hightemperature incubation step to activate the enzyme.

Companion Product:

EpiMark Bisulfite Conversion Kit #E3318S 48 reactions\$162

RR NEBU Epi PCR 165 1m-5

Source: An E. coli strain that carries the Tag DNA Polymerase gene from Thermus aquaticus YT-1.

Concentration: 5,000 units/ml

Vent[®] & Deep Vent[®] DNA Polymerases

#M0254S 200 units72 € 1,000 units 288 € #M0254L Vent (exo-) DNA Polymerase #M0257S 200 units72 € 1,000 units 288 € #M0257L Deep Vent DNA Polymerase #M0258S 200 units 102 € #M0258L 1,000 units 408 € Deep Vent (exo-) DNA Polymerase

Vent DNA Polymerase

#M0259S 200 units 102 € #M0259L 1,000 units 408 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

VENT/DEEP VENT POLYMER	RASES DETAILS
Extension Rate	1 kb/min
Amplicon Size	≤ 6 kb
Fidelity	5-6X <i>Taq</i>
Resulting Ends	Blunt
3´→5´ Exonuclease Activity	Yes (M0254, M0258)
5´→3´ Exonuclease Activity	No
Supplied Buffer	ThermoPol Rxn Buffer

Description: Vent DNA Polymerase was the first highfidelity thermophilic DNA polymerase available for PCR. The fidelity is 5-fold higher than that observed for Tag DNA Polymerase, and is derived in part from an integral $3 \rightarrow 5$ proofreading exonuclease activity. Greater than 90% activity remains following a 1 hour incubation at 95°C.

Deep Vent DNA Polymerase is a modified version of Vent DNA Polymerase. Deep Vent has similar fidelity with even more stability than Vent at temperatures of 95 to 100°C, with a half-life of 8 hours at 100°C.


Vent (exo-) DNA Polymerase has been genetically engineered to eliminate the $3 \rightarrow 5$ proofreading exonuclease activity associated with Vent DNA Polymerase. The fidelity of polymerization by this form is reduced to a level about 2-fold higher than that of Tag DNA Polymerase. Likewise, Deep Vent (exo-) DNA Polymerase has been genetically engineered to eliminate the 3´→ 5´ proofreading exonuclease activity associated with Deep Vent DNA Polymerase.



Source: Vent DNA Polymerase is purified from a strain of E. coli that carries the Vent DNA Polymerase gene from the archaea Thermococcus litoralis. Vent (exo-) is purified from an E. coli strain that carries the Vent (D141A/E143A) DNA Polymerase gene, a genetically engineered form of the native DNA polymerase.

Deep Vent DNA Polymerase is purified from a strain of E. coli that carries the Deep Vent DNA Polymerase gene from Pyrococcus species GB-D. Deep Vent (exo-) is purified from an E. coli strain that carries the Deep Vent (D141A/E143A) DNA Polymerase gene, a genetically engineered form of the native polymerase.

Concentration: 2,000 units/ml



Amplification of Jurkat genomic DNA with Vent DNA Polymerase. Amplicon Sizes are indicated below gel. Marker (M) is the 1 kb DNA Ladder (NEB #N3232).











Luna qPCR and RT-qPCR


- Convenient master mix and supermix formats with user-friendly protocols simplify reaction setup
- Non-interfering, visible tracking dye helps to eliminate pipetting errors
- Novel, thermostable reverse transcriptase (RT) improves performance
- One product per application streamlines selection process

Fluorescence-based quantitative real-time PCR (qPCR) is the gold standard for the detection and quantification of nucleic acids due to its sensitivity and specificity. Luna products from NEB are optimized for qPCR or RT-qPCR, and are available for either intercalating dye or probe-based detection methods. Each Hot Start *Taq*-based Luna qPCR master mix has been formulated with a unique passive reference dye that is compatible across a wide variety of thermal cyclers, regardless of ROX requirements. No additional components are required to ensure compatibility. For two-step RT-qPCR, the LunaScript® RT SuperMix Kit offers a fast (less than 15 minute), robust, and sensitive option for cDNA synthesis upstream of your Luna qPCR experiment.

What is the Difference Between One-Step and Two-Step RT-PCR?

In one-step RT-PCR, cDNA synthesis and PCR amplification are performed in a single reaction. This offers the benefit of a streamlined workflow with fewer chances of contamination. In two-step RT-PCR, cDNA synthesis and PCR amplification are done separately to allow for more flexibility and customization.

Find the right Luna product for your application



Avoiding pipetting errors with Luna



Luna® Universal qPCR & Probe qPCR Master Mixes

Luna Universai qPCR Master Mix			
#M3003S	200 reactions 105 €		
#M3003L	500 reactions 237 €		
#M3003X	1,000 reactions 424 €		
#M3003E	2.500 reactions 948 €		

Luna Universal Probe gPCR Master Mix

Luna Omvorsar	i iobe di ottivias	LUI IVIIA
#M3004S	200 reactions	88 €
#M3004L	500 reactions	198€
#M3004X	1,000 reactions	364 €
#M3004E	2.500 reactions	812 €

Description: The NEB Luna Universal gPCR Master Mix is an optimized 2X reaction mix for real-time qPCR detection and quantitation of target DNA sequences using the SYBR®/FAM channel of most real-time qPCR instruments.

The NEB Luna Universal Probe gPCR Master Mix is a 2X reaction mix optimized for real-time qPCR detection and quantitation of target DNA sequences using hydrolysis probes.

Each Hot Start Tag-based Luna qPCR master mix has been formulated with a unique passive reference dye that

RX PCR WW WW


is compatible across a wide variety of instrument platforms, including those that require a ROX reference signal. This means that no additional components are required to ensure machine compatibility. The mixes also contain dUTP, enabling carryover prevention when reactions are treated with NEB's Antarctic Thermolabile UDG. A blue visible dye assists in tracking the reagents when pipetting into clear, multi-welled PCR plates. These features, combined with rapid, sensitive and precise real-time aPCR performance, make Luna the universal choice for all your qPCR experiments.

Companion Product:

Antarctic Thermolabile	UDG	
#M0372S	100 units	76 €
#M0372L	500 units	304 €

- One product per application
- Convenient master mix formats and userfriendly protocols simplify reaction setup
- Non-interfering, visible tracking dye helps to eliminate pipetting errors
- Rigorously tested to optimize specificity, sensitivity, accuracy and reproducibility
- Unique passive reference dye for compatibility across wide range of instruments

Find the right Luna product for your application



LunaScript® RT SuperMix Kit

#E3010S	25 reactions 125 €
#E3010L	100 reactions 420 €

Companion Products:

Luna Universal qPCR Master Mix

#IVI30035	200 reactions 105 €		
#M3003L	500 reactions 237 €		
#M3003X	1,000 reactions 424 €		
#M3003E	2,500 reactions 948 €		
Luna Universal Probe qPCR Master Mix			
#M3004S	200 reactions 88 €		
#M3004L	500 reactions 198 €		
#M3004X	1,000 reactions 364 €		
#M3004E	2,500 reactions 812 €		

- Less than 15 minute first-strand cDNA synthesis protocol
- Combine with Luna qPCR master mixes for robust RT-qPCR results

Description: LunaScript RT SuperMix Kit is an optimized master mix containing all the necessary components for first strand cDNA synthesis in the context of a two-step RT-qPCR workflow. It features the thermostable Luna Reverse Transcriptase, which supports cDNA synthesis at elevated temperatures. Murine RNase Inhibitor is also included to protect template RNA from degradation. The LunaScript RT SuperMix Kit contains random hexamer and poly-dT primers, allowing for even coverage across the length of the RNA targets.

RX PCR W WW

In addition, the LunaScript RT SuperMix Kit contains an inert blue tracking dye, providing a visual indicator that can be followed throughout the two-step RT-qPCR process. The LunaScript RT SuperMix offers robust, linear, and sensitive detection using total RNA inputs as high as 1 µg and as low as single copies of RNA.

The Luna Universal One-Step RT-qPCR Kit Includes:

- LunaScript RT SuperMix
- No-RT Control Mix
- Nuclease-free Water



How can we ensure best in class performance with Luna?











Heat Inactivation



Luna Universal One-Step RT-qPCR Kit #E3005S 200 reactions 218 € #E3005L 500 reactions 498 € #E3005X 1,000 reactions 869 € #E3005E 2,500 reactions 1924 € Luna Universal Probe One-Step RT-qPCR Kit #E3006S 200 reactions 198 € 500 reactions 448 € #E3006L #E3006X 1.000 reactions 788 € #E3006E 2,500 reactions1738 €

Companion Product:

Antarctic Thermolabile I	JDG	
#M0372S	100 units	76 €
#M0372L	500 units	304 €

- Novel, thermostable reverse transcriptase (RT) improves performance
- WarmStart RT paired with Hot Start Taq increases reaction specificity and robustness
- Non-interfering, visible tracking dye helps to eliminate pipetting errors
- Products perform consistently across a wide variety of sample sources
- Unique passive reference dye for compatibility across wide range of instruments

Description: One-Step RT-qPCR provides a convenient and powerful method for RNA detection and quantitation. In a single tube, RNA is first converted to cDNA by a reverse transcriptase, and then a DNA-dependent DNA polymerase amplifies the cDNA, enabling quantitation via qPCR.

The Luna RT-qPCR kits contain a novel, *in silico*-designed reverse transcriptase (RT) engineered for improved performance. Both the Luna WarmStart Reverse Transcriptase and Hot Start *Taq* DNA Polymerase, included in these kits, utilize a temperature-sensitive, reversible aptamer, which inhibits activity below 45°C. This enables room temperature reaction setup and prevents undesired non-specific activity. Furthermore, the WarmStart RT has increased thermostability, improving performance at higher reaction temperatures.

The NEB Luna Universal One-Step RT-qPCR Kit is optimized for dye-based real-time quantitation of target RNA sequences via the SYBR/FAM fluorescence channel of most real-time instruments.

The NEB Luna Universal Probe One-Step RT-qPCR Kit is optimized for real-time quantitation of target RNA sequences using hydrolysis probes.


Each Hot Start *Taq*-based Luna RT-qPCR kit has been formulated with a unique passive reference dye that is compatible across a wide variety of instrument platforms, including those that require a ROX reference signal. This means that no additional components are required to ensure machine compatibility. The mixes also contain dUTP, enabling carryover prevention when reactions are treated with NEB's Antarctic Thermolabile UDG. A blue visible dye assists in tracking the reagents when pipeting into clear, multi-welled PCR plates. These features, combined with rapid, sensitive and precise real-time qPCR performance, make Luna the universal choice for all your qPCR experiments.

The Luna Universal One-Step RT-qPCR Kit Includes:

- Luna Universal One-Step Reaction Mix
- Luna WarmStart RT Enzyme Mix
- Nuclease-free Water

The Luna Universal Probe One-Step RT-qPCR Kit Includes:

- Luna WarmStart® RT Enzyme Mix
- Luna Universal Probe One-Step Reaction Mix
- Nuclease-free Water



Luna RT-qPCR products offer exceptional sensitivity, reproducibility and performance. RT-qPCR targeting human GAPDH was performed using the Luna Universal One-Step RT-qPCR Kit over an 8-log range of input template concentrations (1 µp - 0.1 pg Jurkat total RNA) with 8 replicates at each concentration. Reaction setup and cycling conditions followed recommended protocols, including a 10-minute RT step at 55°C for the thermostable Luna WarmStart Reverse Transcriptase. NTC = non-template control.



NEB's OEM & Custom Solutions team works closely with our customers to ensure they receive the innovative solutions they need to accelerate their own research. Denisa, Julie and Beth support the operations of this growing branch within the company.

One Tag One-Step RT-PCR Kit

3,000 units71 €

15,000 units 284 €

#E5315S 30 reactions 154 €

Companion Products:

#M0314S

ProtoScript II First Strand cDNA Synthesis Kit				
#E6560S	30 reactions 162 €			
#E6560L	150 reactions 648 €			
ProtoScript II Reven	se Transcriptase			
#M0368S	4,000 units 80 €			
#M0368L	10,000 units 160 €			
#M0368X	40,000 units 574 €			
RNase Inhibitor, Mu	ırine			

For a complete listing of One Tag products, see page 65.

- Combine cDNA synthesis and PCR in a single reaction
- Detect at little as 0.1 pg of a GAPH target
- Robust amplification from 100 bp to 9 kb
- Faster protocols with less hands-on time
- Quick-Load Reaction Mix allows instant gel loading

Description: The One Tag One-Step RT-PCR Kit offers sensitive and robust end-point detection of RNA templates. cDNA synthesis and PCR amplification steps are performed in a single reaction using gene-specific primers, resulting in a streamlined RT-PCR protocol.

The kit combines three optimized mixes: One Tag One-Step Enzyme Mix, One Taq One-Step Reaction Mix and One Tag One-Step Quick-Load Reaction Mix. The Enzyme Mix combines ProtoScript II Reverse Transcriptase, Murine RNase Inhibitor and One Tag Hot Start DNA Polymerase. ProtoScript II Reverse Transcriptase is a mutant M-MuLV reverse transcriptase with reduced RNase H activity and increased thermostability. One Tag Hot Start DNA Polymerase is mixture of a Hot Start Tag DNA Polymerase combined with a proof-reading DNA polymerase, resulting in high-yield amplification with minimal optimization. The One Tag One-Step RT-PCR Kit is capable of amplifying long transcripts up to 9 kb in length.


Two optimized reaction mixes are included, One Tag One-Step Reaction Mix and Quick-Load One Tag One-Step Reaction Mix. The reaction mixes offer robust conditions for both cDNA synthesis and PCR amplification. The unique Quick-Load One Tag One-Step Reaction Mix contains additional dyes, offering color indication for reaction setup as well as direct gel loading.

Both total RNA and mRNA can be used as template. The kit can detect a GAPDH target as low as 0.1 pg per reaction. It can routinely detect RNA targets up to 9 kb. The One Tag One-Step RT-PCR Kit is capable of multiplex detection of two or three targets.

RX PCR WW

The One Tag One-Step RT-PCR Kit Includes:

- One Tag One-Step Enzyme Mix
- One Tag One-Step Reaction Mix
- Quick-Load One Tag One-Step Reaction Mix
- Nuclease-free Water



Detection of RNA templates of different length. About 100 ng of Jurkat total RNA was used in 50 µl reactions following the standard protocol. The target sizes were Lane 1: 0.7 kb, Lane 2: 1.1 kb, Lane 3: 1.9 kb, Lane 4: 2.3 kb, Lane 5: 2.5 kb, Lane 6: 5.5 kb, Lane 7: 7.6 kb and Lane 8: 9.3 kb. The marker lane (M) contains Quick-Load 1 kb Plus DNA Ladder (NEB #N0469).

One Tag RT-PCR Kit

30 reactions 143 € #E5310S

Companion Products:

RNase Inhibitor,	Murine
#M0314S	3,000 units 71 €
#M0314L	15,000 units 284 €
M-MuLV Reverse	e Transcriptase
#M0253S	10,000 units 70 €
#M0253L	50,000 units 280 €
One Taq Hot Start #M0484S #M0484L	2X Master Mix with Standard Buffer 100 rxns (50 µl vol)
PolyA Spin mRN #S1560S	A Isolation Kit 8 isolations 224 €

For a complete listing of One Tag products, see page 65.

- Robust RT-PCR reactions
- Convenient master mix formats
- OneTag DNA Polymerase offers robust amplification for a wide range of templates


Description: One Taq RT-PCR Kit combines two powerful mixes. M-MuLV Enzyme Mix and One Taa Hot Start 2X Master Mix with Standard Buffer for 2-step RT-PCR applications. The two mixes require minimal handling during reaction setup and yet offer consistent and robust RT-PCR reactions.

The first strand cDNA synthesis is achieved by using two optimized mixes, M-MuLV Enzyme Mix and M-MuLV Reaction Mix. M-MuLV Enzyme Mix combines M-MuLV Reverse Transcriptase and RNase Inhibitor. Murine while M-MuLV Reaction Mix contains dNTPs and an optimized buffer. The kit also contains two optimized primers for reverse transcription and nuclease-free water. An anchored Oligo-dT primer [d(T)₂₂VN] forces the primer to anneal to the beginning of the polyA tail. The optimized Random Primer Mix provides random and consistent priming sites covering the entire RNA templates including both mRNAs and non-polyadenylated RNAs.

The One Tag RT-PCR Kit Includes:

- 10X M-MuLV Enzyme Mix
- 2X M-MuLV Reaction Mix
- One Tag Hot Start 2X Master Mix with Standard Buffer
- Random Primer Mix (60 μM)
- Oligo d(T)₂₃VN Primer (50 μM)**
- Nuclease-free Water

RX PCR YM6



First strand cDNA synthesis. Synthesis was carried out in the presence of 1X M-MuLV Enzyme Mix at 42°C using 0.5 μg of human spleen total RNA in the presence of dT₂₃VN (lanes 1, 4 and 7) or Random Hexamer Mix (lanes 2, 5 and 8). No-RT controls were lanes 3, 6 and 9. One Tag Hot Start Master Mix was used to amplify a 1.5 kb fragment of betaactin gene, a 0.6 kb fragment of GAPDH gene, and a 5.5 kb fragment from p532 gene in 35 cycles. The marker lane (M) contains 1 kb Plus DNA Ladder (NEB #N3200).







PCR PCR Enzyme

Hot Start/ WarmStar

BSA Requires BSA

65° Heat Inactivation

Tm⁻⁵ Annealing Temperature

Epi EpiMark

Extraction Free

^{**} Oligo $d(T)_{23}$ VN and Random Primer Mix contain 1 mM dNTP

Bst DNA Polymerase-based Products for Isothermal DNA Amplification

	5'→ 3' EXO ACTIVITY	AMPLIFICATION SPEED	ROOM TEMPERATURE SETUP	REVERSE Transcriptase Activity	INHIBITOR Tolerance	APPLICATIONS
Bst DNA Polymerase, Full Length	**	N/A	N/A	N/A	*	Nick translation reactions at elevated temperatures Primer extension
Bst DNA Polymerase, Large Fragment	N/A	*	N/A	*	*	General strand-displacement reactions
Bst 2.0 DNA Polymerase	N/A	**	N/A	**	*	Improved LAMP, SDA, and other amplification reactions Minimal effect of substitution of dTTP with dUTP
Bst 2.0 WarmStart DNA Polymerase	N/A	**	***	**	**	Consistent, room-temperature, and high-throughput amplification assays Minimal effect of substitution of dTTP with dUTP
Bst 3.0 DNA Polymerase	N/A	***	**	***	***	Fastest, most robust LAMP and RT-LAMP reactions High reverse transcriptase activity up to 72°C Strand displacement DNA synthesis

Optimal, recommended product for selected application

Bst DNA Polymerases

Bst DNA Polymerase	, Large Fragment		
#M0275S	1,600 units69 €		
#M0275L	8,000 units 276 €		
for high (15X) concentration #M0275M	8,000 units 276 €		
Bst DNA Polymerase, Full Length			
#M0328S	500 units69 €		

Bst 2.0 DNA Polymerase

#M0537S	1,600 units	69 €
#M0537L	8,000 units	276€
f 11 1 (45)() 1 11	_	

for high (15X) concentration

#M0537M 8,000 units 276 €

Bst 2.0 WarmStart DNA Polymerase #M0538S 1 600 units

#M0538S	1,600 units	77€
#M0538L	8,000 units	308 €
for high (15X) concentration	1	

#M0538M

8,000 units 308 €

Bst 3.0 DNA Polymerase

#IVIU3745	1,000 utilis	
#M0374L	8,000 units	280€
for high (15X) concentration	n	

// AOO 7 40

#M0374M 8,000 units 280 €

Companion Products:

WarmStart RTx Revers	e Transcriptase	
#M0380S	50 reactions	65€
#M0380L	250 reactions	260€

Tte UvrD Helicase

#M1202S 0.5 μg71 €


For a complete listing of Deoxynucleotide Solutions, see page 79.

Description: Bst DNA Polymerase, Large Fragment, is the portion of the Bacillus stearothermophilus DNA Polymerase protein that contains the $5 \rightarrow 3$ polymerase activity, but lacks $5 \rightarrow 3$ exonuclease activity.

Bst DNA Polymerase, Full Length is the full length polymerase from Bacillus stearothermophilus. It has 5 → 3 1 polymerase and double-strand specific 5 → 3 exonuclease activities, but lacks 3 → 5 exonuclease activity.

Bst 2.0 DNA Polymerase is an in silico designed homologue of Bst DNA Polymerase I. Large Fragment. It contains 5 → 3 DNA polymerase activity and strong strand displacement activity but lacks 5' -> 3' exonuclease activity. It has improved amplification speed, yield, salt tolerance and thermostability compared to wild-type Bst DNA Polymerase, Large Fragment.

Bst 2.0 WarmStart DNA Polymerase utilizes aptamer technology to inhibit activity at non-permissive temperatures (< 50°C). Like "Hot Start" PCR polymerases, the WarmStart feature enables room temperature set up and prevents non-templated addition of nucleotides, increasing reaction efficiencies. Additionally, no separate activation step is required to release the aptamer from the enzyme. Bst 2.0 WarmStart DNA Polymerase permits reaction temperatures from 60-72°C




RX NEBU 65° W

Bst 3.0 DNA Polymerase is a similarly designed in silico homologue engineered for improved performance and increased reverse transcriptase activity. Bst 3.0 DNA Polymerase contains $5 \rightarrow 3$ DNA polymerase activity with either DNA or RNA templates but lacks $5 \rightarrow 3$ and 3'→ 5' exonuclease activity. It demonstrates robust performance in the presence of inhibitors and significantly increase reverse transcriptase activity compared to Bst DNA Polymerase.

Concentration: Bst DNA Polymerase, Full Length: 5,000 units/ml. All others: 8,000 and 120,000 units/ml

Heat Inactivation: 80°C for 20 minutes

Usage Notes: No Bst DNA Polymerase-based products can be used for thermal cycle sequencing or PCR. Bst 2.0 WarmStart DNA Polymerase permits reaction temperatures from 60–72°C. Generally, reaction temperatures above 72°C are not recommended for any Bst DNA Polymerase-based product.



Benefits of Bst 2.0 WarmStart: Identical LAMP reactions were run either immediately after setup (solid line) or after a 2 hour incubation at 25°C. Without the protection from Bst 2.0 WarmStart, this room temperature incubation results in variable LAMP performance. Bst 2.0 WarmStart provides more consistent amplification reaction and enables room-temperature and high-throughput setup.

Works well for selected application

Will perform selected application, but is not recommended

Not applicable to this application

WarmStart® Colorimetric LAMP Master Mix (DNA & RNA)

#M1800S 100 reactions 212 € 500 reactions 848 € #M1800I

- LAMP amplification of RNA and DNA
- One-step field and point-of-need LAMP reactions
- Set up reactions at room temperature with our unique dual WarmStart formulation
- Simple visual indication of amplification for easy detection and analysis

Description: The WarmStart Colorimetric LAMP 2X Master Mix is an optimized formulation of Bst 2.0 Warm-Start DNA Polymerase and WarmStart RTx in a special low-buffer reaction solution containing a visible pH indicator for rapid and easy detection of Loop-Mediated Isothermal Amplification (LAMP) and RT-LAMP reactions. This system is designed to provide a fast, clear visual detection of amplification based on the production of protons and subsequent drop in pH that occurs from the extensive DNA polymerase activity in a LAMP reaction. producing a change in solution color from pink to yellow (an overview of LAMP and primer design can be found in our video library, https://www.neb.com/tools-andresources/video-library). This mix can be used for any LAMP or RT-LAMP reaction and requires only a heated chamber and samples, with readout of positive amplification judged by eye in 15-40 minutes.

Positive Negative Start 30 min. 65°C

RX 65° 86 65°

WarmStart LAMP Kit (DNA & RNA)

#F1700S 100 reactions 212 € #E1700L 500 reactions 848 €

Companion Product:

Tte UvrD Helicase #M1202S

0.5 μg\$71


- LAMP amplification of RNA and DNA targets
- Improve LAMP specificity and sensitivity with optimized master mixes
- Set up reactions at room temperature with our unique dual WarmStart formulation
- Use with a variety of detection methods including fluorescence, turbidity, visual detection and electrophoresis

Description: The WarmStart LAMP Kit (DNA & RNA) is designed to provide a simple, one-step solution for Loop-Mediated Isothermal Amplification (LAMP) of DNA or RNA (RT-LAMP) targets. LAMP and RT-LAMP are commonly used isothermal amplification techniques that provides rapid detection of a target nucleic acid using LAMP-specific primers (supplied by the user) and a strand-displacing DNA polymerase. This kit is supplied with the WarmStart LAMP 2X Master Mix, which contains a blend of Bst 2.0 WarmStart DNA Polymerase and WarmStart RTx Reverse Transcriptase in an optimized LAMP buffer solution. Both Bst 2.0 Warm-Start DNA Polymerase and WarmStart RTx Reverse Transcriptase have been engineered for improved performance in LAMP and RT-LAMP reactions. A fluorescent dye is also supplied to enable real-time fluorescence measurement of LAMP. The WarmStart LAMP Kit is compatible with multiple detection methods, including turbidity detection, real-time fluorescence detection (when used with LAMP fluorescent dye) and end-point visualization.

RN 65° W

The WarmStart LAMP Kit Includes:

- WarmStart LAMP 2X Master Mix
- LAMP Fluorescent Dye (50X)



















IsoAmp® II Universal tHDA Kit

#H0110S 50 reactions 340 €

- Easy-to-use for assay development
- Helicase eliminates need for thermocycler
- Reactions performed at constant temp
- Amplify & detect short DNA sequences (70-120 bp)
- Use with a variety of templates (microbial genomic DNA, viral DNA, plasmid DNA and cDNA)
- Amplify a single copy of target DNA by tHDA when optimized primers and buffer are used

Description: Thermophilic Helicase-Dependent Amplification (tHDA) is a novel method for isothermal amplification of nucleic acids. Like PCR, the tHDA reaction selectively amplifies a target sequence defined by two primers. However, unlike PCR, tHDA uses an enzyme called a helicase to separate DNA, rather than heat. This allows DNA amplification without the need for thermocycling. The tHDA reaction can also be coupled with reverse transcription for RNA analysis.

IsoAmp II Universal tHDA Kit is based on a secondgeneration thermophilic Helicase-Dependent Amplification platform. The reactions supported by IsoAmp II Universal tHDA Kit include tHDA, reverse transcription HDA (RT-HDA), real-time quantitative HDA (qHDA) and real-time quantitative RT-HDA (qRT-HDA), from a single reaction buffer.

The IsoAmp II Universal tHDA Kit Includes:

- IsoAmp dNTP solution and IsoAmp Enzyme Mix
- 10X Annealing Buffer II, 100 mM MgSO₄ and 500 mM NaCl
- Control template and amplification primers

RX

...........

HDA technology. Helicase Dependent Amplification: Step 1: Helicase unwinding and primer binding. Step 2: DNA polymerization. Step 3: DNA amplification.

Developed by BioHelix Corporation a NEB-affiliated company, now part of Quidel Corporation. ISOAMP® is a registered trademark of BioHelix Corporation

phi29 DNA Polymerase

#M0269S 250 units58 € #M0269L 1,250 units 232 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

- Extreme processivity
- Extreme strand displacement
- Replication requiring a high degree of strand displacement and/or processive synthesis
- High-fidelity replication at moderate temperatures

Description: phi29 DNA Polymerase is the replicative polymerase from the Bacillus subtilis phage phi29 (φ29). This polymerase has exceptional strand displacement and processive synthesis properties. The polymerase has an inherent 3´→ 5´ proofreading exonuclease activity.



Concentration: 10,000 units/ml Heat Inactivation: 65°C for 10 minutes



Steve is a Product Marketing Manager and has been with NEB for 2 years. In his role, Steve works with our scientists specializing in amplification technologies to bring the most useful products to market.



PreCR® Repair Mix

#M0309S 30 reactions 157 € #M03091 150 reactions 628 €

Companion Product:

beta-Nicotinamide adenine dinucleotide (NAD+) 0.2 ml 32 €

- Repair DNA prior to its use in DNArelated technologies
- Easy-to-use protocols
- Does not harm template

Need to repair FFPE-treated DNA prior to next gen sequencing? See page 153 for more information. **Description:** The PreCR Repair Mix is an enzyme cocktail formulated to repair damaged template DNA prior to its use in the polymerase chain reaction (PCR), microarrays, or other DNA technologies. The PreCR Repair Mix is active on a broad range of DNA damage, including those that block PCR (e.g., apurinic/apyrimidinic sites, thymidine dimers, nicks and gaps) and those that are mutagenic (e.g., deaminated cytosine and 8-oxo-guanine). In addition, it will remove a variety of moieties from the 3´ end of DNA leaving a hydroxyl group. The PreCR Repair Mix will not repair all damage that inhibits/interferes with PCR. It can be used in conjunction with any thermophilic polymerase.

Applications:

 Repair DNA prior to its use as a template in PCR or other DNA technologies

Reagents Supplied:


1X PreCR Repair Mix 10X ThermoPol Reaction Buffer 100X NAD+ solution Control Template (UV damaged λ DNA) PCR primers for control template Purified BSA

Repair of different types of DNA damage with the PreCR Repair Mix. The gel shows amplification of damaged DNA that was either not treated (-) or treated (+) with the PreCR Repair Mix. Type of DNA damage is shown. Note: heat treated DNA is incubated at 99°C for 3 minutes. Marker (M) is the 1 kb Plus DNA Ladder (NEB #N3200).

RR Wh

Types of DNA Damage

DNA DAMAGE	CAUSE	REPAIRED BY PRECR REPAIR MIX?
abasic sites	hydrolysis	yes
nicks	hydrolysis nucleases shearing	yes
thymidine dimers	UV radiation	yes
blocked 3´-ends	multiple	yes
oxidized guanine	oxidation	yes
oxidized pyrimidines	oxidation	yes
deaminated cytosine	hydrolysis	yes
fragmentation	hydrolysis nucleases shearing	no
protein-DNA crosslinks	formaldehyde	no



Sulfolobus DNA Polymerase IV

#M0327S

100 units 119 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

- Synthesis of DNA through DNA lesions (lesion bypass)
- DNA Repair

Description: Sulfolobus DNA Polymerase IV is a thermostable Y-family lesion-bypass DNA Polymerase that efficiently synthesizes DNA across a variety of DNA template lesions.

RX NEBU 55° WW

Concentration: 2.000 units/ml

Therminator™ DNA Polymerase

#M0261S 200 units 101 € #M0261L 1,000 units 404 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

- Incorporation of modified nucleotides
- DNA sequencing by partial ribosubstitution
- DNA sequencing or SNP analysis using dideoxy or acyclo chain terminators

Description: Therminator DNA Polymerase is a 9°N™ DNA Polymerase variant with an enhanced ability to incorporate modified substrates such as dideoxynucleotides, ribonucleotides and acyclonucleotides.

Source: An E. coli strain that carries the 9°N (D141A / E143A / A485L) DNA Polymerase gene, a genetically engineered form of the native DNA polymerase from Thermococcus species 9°N-7.



Concentration: 2.000 units/ml

Usage Notes: Amplification of extended regions may require optimization of reaction conditions.





















DNA Polymerase I (E. coli)

#M0209S 500 units68 € #M0209L 2,500 units272 €

- Nick translation of DNA
- Second strand cDNA synthesis

Description: DNA Polymerase I (*E. coli*) is a DNA-dependent DNA polymerase with inherent $3' \rightarrow 5'$ and $5' \rightarrow 3'$ exonuclease activities. The $5' \rightarrow 3'$ exonuclease activity removes nucleotides ahead of the growing DNA chain, allowing nick translation.

RX NEB 2 37° 1/54

Concentration: 10,000 units/ml

Heat Inactivation: 75°C for 20 minutes

Usage Notes: DNase I is not included with this enzyme and must be added for nick translation reactions.

DNA Polymerase I, Large (Klenow) Fragment

#M0210S 200 units63 € #M0210L 1,000 units252 €

for high (10X) concentration

#M0210M 1,000 units 252 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

- Generates probes using random primers
- Removal of 3´ overhangs or fill-in of 5´ overhangs to form blunt ends
- Second strand cDNA synthesis

Description: DNA Polymerase I, Large (Klenow) Fragment was originally derived as a proteolytic product of *E. coli* DNA Polymerase I that retains polymerase and $3' \rightarrow 5'$ exonuclease activity, but lacks $5' \rightarrow 3'$ exonuclease activity. Klenow retains the polymerization fidelity of the holoenzyme without degrading 5' termini.

Source: An *E. coli* strain that contains the *E. coli* polA gene that has had its $5 \rightarrow 3$ exonuclease domain removed.

RX NEB 2 25° 166

Concentration: 5,000 and 50,000 units/ml Heat Inactivation: 75°C for 20 minutes

Usage Notes: Elevated temperatures, excessive amounts of enzyme, failure to supplement with dNTPs or long reaction times will result in recessed ends due to the $3' \rightarrow 5'$ exonuclease activity of the enzyme.

Klenow Fragment $(3' \rightarrow 5' \text{ exo}^-)$

#M0212S 200 units63 € #M0212L 1,000 units252 €

for high (10X) concentration #M0212M 1,000 units 252 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

- Generates probes using random primers
- Random priming labeling
- Second strand cDNA synthesis

Description: Klenow Fragment ($3' \rightarrow 5'$ exo-) is an N-terminal truncation of DNA Polymerase I that retains polymerase activity, but has lost the $5' \rightarrow 3'$ exonuclease activity, and has mutations (D355A, E357A) that abolish the $3' \rightarrow 5'$ exonuclease activity.

Source: An *E. coli* strain containing a plasmid with a fragment of the *E. coli polA* (D355A, E357A) gene starting at codon 324.

RX NEB 2 37° 1/64

Concentration: 5,000 and 50,000 units/ml Heat Inactivation: 75°C for 20 minutes

Usage Notes: Klenow Fragment ($3' \rightarrow 5' \text{ exo}^-$) is not suitable for generating blunt ends because it lacks the $3' \rightarrow 5'$ exonuclease activity necessary to remove non-templated 3' additions.

T4 DNA Polymerase

#M0203S 150 units66 € #M0203L 750 units264 €

Companion Product:

Quick Blunting™ Kit

#E1201S 20 reactions 85 € #E1201L 100 reactions 340 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

- Gap filling (no strand displacement activity)
- Removal of 3´ overhangs or fill-in of 5´ overhangs to form blunt ends
- Probe labeling using replacement synthesis
- Single-strand deletion subcloning

Description: T4 DNA Polymerase catalyzes the synthesis of DNA in the $5' \rightarrow 3'$ direction and requires the presence of template and primer. This enzyme has a $3' \rightarrow 5'$ exonuclease activity which is much more active than that found in *E. coli* DNA Polymerase I. Unlike DNA Polymerase I, T4 DNA Polymerase does not have a $5' \rightarrow 3'$ exonuclease function.

Concentration: 3,000 units/ml

RX NEB 2.1

Heat Inactivation: 75°C for 20 minutes

Usage Notes: Elevated temperatures, excessive amounts of enzyme, failure to supplement with dNTPs or long reaction times will result in recessed ends due to the $3' \rightarrow 5'$ exonuclease activity of the enzyme.

T7 DNA Polymerase (unmodified)

#M0274S 300 units70 € #M0274L 1,500 units 280 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

 Gap-filling reaction (no strand displacement)

Description: T7 DNA Polymerase catalyzes the replication of T7 phage DNA during infection. The protein dimer has two catalytic activities: DNA polymerase activity and strong $3 \rightarrow 5$ exonuclease. The high fidelity and rapid extension rate of the enzyme make it particularly useful in copying long stretches of DNA template.

Source: T7 DNA Polymerase consists of two subunits: T7 gene 5 protein (84 kDa) and E. coli thioredoxin (12 kDa). Each protein is cloned and overexpressed in a T7 expression system in E. coli.

Reaction Conditions: 1X T7 DNA Polymerase Reaction Buffer. Supplement with BSA and dNTPs (not included). Incubate at 37°C. Heat inactivation: 75°C for 20 minutes.

RX NEBU BSA 37° 📆

Unit Definition: One unit is defined as the amount of enzyme that will incorporate 10 nmol of dNTP into acid insoluble material in 30 minutes at 37°C.

Concentration: 10,000 units/ml

Usage Notes: The high polymerization rate of the enzyme makes long incubations unnecessary. T7 DNA Polymerase is not suitable for DNA sequencing.

Bsu DNA Polymerase, Large Fragment

#M0330S 200 units66 € #M0330L 1,000 units 264 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

- Random primer labeling
- Second strand cDNA synthesis
- Single dA tailing
- Strand displacement DNA synthesis

Description: Bsu DNA Polymerase I, Large Fragment retains the $5 \rightarrow 3$ polymerase activity of the *Bacillus* subtilis DNA polymerase I, but lacks the 5 → 3 exonuclease domain. This large fragment naturally lacks 3'→ 5' exonuclease activity.

Source: An *E. coli* strain that contains the *Bacillus* subtilis DNA polymerase I gene (starting from codon 297 thus lacking the $5' \rightarrow 3'$ exonuclease domain).

Concentration: 5,000 units/ml

Heat Inactivation: 75°C for 20 minutes

RX NEB 2 37° 164

Usage Notes: Bsu DNA Polymerase, Large Fragment is not suitable for generating blunt ends because it lacks the 3'→ 5' exonuclease activity necessary to remove nontemplated 3' additions.

Bsu DNA Polymerase, Large Fragment retains 50% activity at 25°C and is twice as active as Klenow Fragment $(3 \rightarrow 5 \text{ exo})$ at this temperature.

Terminal Transferase

#M0315S 500 units70 € #M0315L 2,500 units 280 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

- Addition of homopolymer tails to the 3' ends of DNA
- Labeling the 3´ ends of DNA with modified nucleotides (e.g., ddNTP, DIG-dUTP)
- TUNEL assay (in situ localization of apoptosis)
- TdT dependent PCR

Description: Terminal Transferase (TdT) is a template independent polymerase that catalyzes the addition of deoxynucleotides to the 3' hydroxyl terminus of DNA molecules. Protruding, recessed or blunt-ended double or single-stranded DNA molecules serve as a substrate for TdT. The 58.3 kDa enzyme does not have 5' or 3' exonuclease activity. The addition of Co2+ in the reaction makes tailing more efficient.

RX NEBU 37° 😘

Concentration: 20,000 units/ml Heat Inactivation: 75°C for 20 minutes



78





PCR PCR Enzyme

Hot Start/ WarmStart



65° Heat Inactivation







Polymerase Reaction Buffers

Q5

Q5 Reaction Buffer Pack [Q5 Reaction Buffer (5X), Q5 High GC Enhancer (5X)] #B9027S 6.0 ml27 €

Phusion*

Phusion HF Buffer Pack [Phusion HF Reaction Buffer (5X), MgCl₂ (50 mM), DMSO]

#B0518S 6.0 ml23 €

Phusion GC Buffer Pack [Phusion GC Reaction Buffer (5X), MgCl₂ (50 mM), DMSO]

6.0 ml23 €

Taq

#B0519S

Standard $\overline{\it Taq}$ Reaction Buffer [Standard $\overline{\it Taq}$ Reaction Buffer (10X), MgCl $_2$ (25 mM)]

#B9014S 6.0 ml22 €

Standard *Taq* (Mg-free) Reaction Buffer Pack [Standard *Taq* (Mg-free) Reaction Buffer (10X), MgCl₂ (25 mM)]

#B9015S 6.0 ml22 €

Other

ThermoPol Reaction Buffer Pack [ThermoPol Reaction Buffer (10X), MgSO₄ (100 mM)]

#B9004S 6.0 ml22 €

Isothermal Amplification Buffer Pack [Isothermal Amplification Buffer (10X)] #B0537S 6.0 ml27 €

Isothermal Amplification Buffer II Pack [Isothermal Amplification Buffer II (10X)] #B0374S 6.0 ml26 € **Description:** Q5 Reaction Buffer and High GC Enhancer are provided with both Q5 and Q5 Hot Start High-Fidelity DNA Polymerases.

Phusion High-Fidelity DNA Polymerase is supplied with 5X Phusion HF Buffer, 5X Phusion GC Buffer, DMSO, and 50 mM MgCl₂.

Standard *Taq* Reaction Buffer is provided with *Taq* DNA Polymerase as an alternative to the ThermoPol Reaction Buffer.

ThermoPol Reaction Buffer is provided with *Taq*, Vent, Deep Vent, *Bst* Full Length and *Bst* Large Fragment, Sulfolobus IV and Therminator DNA Polymerases; this buffer contains 2 mM MgSO₄ when the buffer is diluted to its final 1X concentration.

Isothermal Amplification Buffer is supplied with *Bst* 2.0 and *Bst* 2.0 WarmStart DNA Polymerases.

Isothermal Amplification Buffer II is supplied with ${\it Bst}\,3.0$ DNA Polymerase.

Nucleotides

Acyclonucleotide Set

#N0460S 0.5 μmol of each70 €

Deoxynucleotide (dNTP) Solution Set

#N0446S 25 μmol of each 164 €

Deoxynucleotide (dNTP) Solution Mix

#N0447S 8 μ mol of each61 \in #N0447L 40 μ mol of each244 \in

Ribonucleotide Solution Set

#N0450S 10 μmol of each69 € #N0450L 50 μmol of each277 €

Ribonucleotide Solution Mix

#N0466S 10 µmol of each72 € #N0466L 50 µmol of each 288 €

7-deaza-dGTP

#N0445S 0.3 µmol of each67 € #N0445L 1.5 µmol of each268 €

Adenosine 5´-Triphosphate (ATP)

#P0756S 1.0 ml35 € #P0756S 5.0 ml 138 € 5-methyl-dCTP

#N0356S 1 μmol70 €

dATP Solution

#N0440S 25 μmol53 €

Description:

Deoxynucleotide Solution Set:

Four separate solutions of ultrapure deoxynucleotide (dATP, dCTP, dGTP and dTTP). Each deoxynucleotide is supplied at a concentration of 100 mM.

Deoxynucleotide Solution Mix:

An equimolar solution of ultrapure dATP, dCTP, dGTP and dTTP. Each deoxynucleotide is supplied at a concentration of 10 mM.

Ribonucleotide Solution Set:

Four separate solutions of ATP, CTP, GTP and UTP, pH 7.5, as sodium salts.

Ribonucleotide Solution Mix:

A buffered equimolar solution of ribonucleotide triphosphates: rATP, rCTP, rGTP and rUTP. Each nucleotide is supplied at a concentration of 25 mM (total rNTP concentration equals 100 mM).

7-deaza-dGTP:

7-deaza contains a 5 mM solution of 7-deaza-dGTP as a dilithium salt.

5-methyl-dCTP:

dm5CTP supplied as 10 mM solution at pH 7.0.

dATP Solution:

dATP Solution contains a 100 mM solution of dATP as a sodium salt at pH 7.4.

Acyclonucleotide Set:

Four separate tubes of acyNTPs (acyATP, acyCTP, acyGTP and acyTTP). Acyclonucleotides are supplied as a dry powder. Addition of 50 μ l of distilled or de-ionized (Milli-Q®) water will result in a final concentration of 10 mM acyNTP.

Acyclonucleotides (acyNTPs) act as chain terminators and are thus useful in applications that normally employ dideoxynucleotides such as DNA sequencing and SNP detection. AcyNTPs are especially useful in applications with archaeal DNA Polymerases, more specifically with Therminator DNA Polymerase. Therminator DNA Polymerase is an engineered enzyme with an increased capacity to incorporate analogs with altered sugars, such as ribonucleotides, dideoxynucleotides, 2´ deoxynucleotides and especially acyclo-base analogs.

For RNA Cap Analogs, see page 190.

MILLI-Q® is a registered trademark of Millipore, Inc.

^{*} Phusion DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific. Phusion[®] is a registered trademark and property of Thermo Fisher Scientific.

Products for cDNA Synthesis

CDNA SYNTHESIS	FEATURES	SIZE	PRICE
KITS	•		
NEW LunaScript® RT SuperMix Kit (NEB #E3010)	Ideal for cDNA synthesis in a two-step RT-qPCR workflow Single tube supermix contains random hexamer and oligo-dT primers, dNTPs, Murine RNase Inhibitor, and Luna Reverse Transcriptase Visible blue tracking dye for easy reaction setup Fast 13-minute protocol	25/100 rxns	125 €/420 €
ProtoScript® II First Strand cDNA Synthesis Kit (NEB #E6560)	Generates cDNA at least 10 kb in length Contains ProtoScript II Reverse Transcriptase, an enzyme with increased thermostability and reduced RNase H activity Convenient 2-tube kit includes dNTPs, Oligo-dT primer and Random Primer Mix	30/150 rxns	162 €/648 €
ProtoScript First Strand cDNA Synthesis Kit (NEB #E6300)	Generates cDNA at least 5 kb in length Contains M-MuLV Reverse Transcriptase Convenient 2-tube kit includes dNTPs, Oligo-dT primer and Random Primer Mix	30/150 rxns	155 €/620 €
NEW Template Switching RT Enzyme Mix (NEB #M0466)	Incorporates a universal adaptor sequence at the 3´ end of cDNA during the RT reaction Enzyme mix and buffer are optimzed for efficient template switching RT enzyme mix includes RNase Inhibitor High sensitivity for cDNA amplification — enables transcriptome analysis by RNA-seq from single cells or as low as 2 pg of human total RNA Robust and simple workflow for 5´ Rapid Amplification of cDNA Ends (RACE) Retains the complete 5´ end of transcripts for 2nd Strand cDNA Synthesis	20/100 rxns	88 €/352 €
STANDALONE REAGENTS			
ProtoScript II Reverse Transcriptase (NEB #M0368) An alternative to SuperScript® II	RNase H ⁻ mutant of M-MuLV Reverse Transcriptase with increased thermostability and reduced RNase H activity Increased reaction temperatures (37–50°C)	4,000/10,000/40,000 units	80 €/160 €/ 574 €
M-MuLV Reverse Transcriptase (NEB #M0253)	Robust reverse transcriptase for a variety of templates Standard reaction temperatures (37–45°C)	10,000/50,000 units	70 €/288 €
AMV Reverse Transcriptase (NEB #M0277)	Robust reverse transcriptase for a broad temperature range (37–52°C) Can be used for templates requiring higher reaction temperatures	200/1,000 units	71 €/286 €
WarmStart RTx Reverse Transcription (NEB #M0380)	Permits room temperature reaction setup Increased reaction temperatures (50–65°C) Optimized for RT-LAMP isothermal detection	50/250 rxns	65 €/260 €

More information about our reverse transcriptases and cDNA synthesis kits can be found in the RNA analysis chapter, pages 192–194.

Ellen has been with NEB for 29 years, and will retire this spring with much fanfare. NEB's Research and Development teams have been lucky to have a Principal Development Scientist with Ellen's dedication and expertise. She is also an advocate for science education and has brought joy to many children at the holiday party as Mrs. Claus. Ellen is looking forward to spending more time gardening, knitting and traveling in the coming years.











Monarch PCR & DNA Cleanup Kit (5 µg)

#T1030S 50 preps88 € #T1030L 250 preps398 €

Companion Products:

Monarch DNA Cleanup Columns (5 µg) #T1034L 100 columns 128 € Monarch DNA Cleanup Binding Buffer 235 ml97 € Monarch DNA Wash Buffer #T1032L 25 ml32 € Monarch DNA Elution Buffer #T1016L 25 ml32 € Monarch Plasmid Miniprep Kit #T1010S 50 preps 75 € #T1010L 250 preps 325 € Monarch DNA Gel Extraction Kit88€ #T1020S 50 preps #T1020L 250 preps 398 €

- Elute in as little as 6 μl
- Prevent buffer retention and salt carry-over with optimized column design
- Purify small DNA and oligos with a slight protocol modification
- Save time with fast, user-friendly protocol
- Purchase optimized kit formats or buffers & columns separately for your convenience

Description: The Monarch PCR & DNA Cleanup Kit (5 μg) is a rapid and reliable method for the purification and concentration of up to 5 µg of high-quality, doublestranded DNA from enzymatic reactions such as PCR, restriction digestion, ligation and reverse transcription. This method employs a bind/wash/elute workflow with minimal incubation and spin times, resulting in purification in less than 5 minutes. DNA Cleanup Binding Buffer is used to dilute the samples and ensure they are compatible for loading onto the proprietary silica matrix under high salt conditions. The DNA Wash Buffer ensures enzymes, short primers (≤ 40 nt), detergents and other low-molecular weight reaction components (e.g., nucleotides, DMSO, betaine) are removed, thereby allowing low-volume elution of concentrated, highpurity DNA. Eluted DNA is ready for use in restriction digests, DNA sequencing, ligation and other enzymatic manipulations. The unique column design ensures no buffer retention and no carryover of contaminants, allowing elution of sample in volumes as low as 6 µl. A slight protocol modification enables purification of small DNA and oligonucleotides.

Applications:

- PCR cleanup
- · Enzymatic reaction cleanup
- cDNA cleanup
- · Labeling cleanup
- Plasmid cleanup
- · Oligo cleanup

The Monarch PCR & DNA Cleanup Kit Includes:

- Monarch DNA Cleanup Columns (5 μg)
- Monarch DNA Cleanup Binding Buffer
- Monarch DNA Wash Buffer
- Monarch DNA Elution Buffer
- Monarch Collection Tubes

With Monarch PCR & DNA Cleanup Kit, you can purify your DNA in as little as 5 minutes.

NEW

Exo-CIP™ Rapid PCR Cleanup Kit

#E1050S 100 reactions82 € #E1050L 400 reactions270 €


- 5 minute protocol for enzymatic cleanup of primers and dNTPs
- Improves sequencing results, allowing longer reads

Description: The Exo-CIP Rapid PCR Cleanup Kit contains optimized formulations of thermolabile Exonuclease I and thermolabile Calf Intestinal Phosphatase and is used to rapidly degrade residual PCR primers and dephosphorylate excess dNTPs after amplification. Degradation occurs in only 4 minutes at 37°C, and is immediately followed by rapid inactivation of the enzymes by heating for 1 minute at at 80°C. In just 5 minutes, the PCR product is ready for downstream analysis such as Sanger sequencing, SNP detection, or library preparation for NGS. The Exo-CIP Rapid PCR Cleanup Kit is compatible with all commonly-used reaction buffers.

The Exo-CIP Rapid PCR Cleanup Kit Includes:

- Exo-CIP Tube A (thermolabile Exo I)
- Exo-CIP Tube B (thermolabile CIP)

RX



Exo-CIP Rapid PCR Cleanup Kit workflow. 1 µl of Exo-CIP Tube A (thermolabile Exo I) and 1 µl of Exo-CIP Tube B (thermolabile CIP) are added to the PCR product to degrade excess primers and dNTPs. The mixture is incubated at 37°C for 4 minutes, followed by a 1 minute incubation at 80°C to irreversibly inactivate both enzymes. The cleaned PCR product is ready for downstream applications or analysis.





Whisky and Oysters

Oysters are essential members of a diverse marine ecosystem because of their specific roles in contributing to biodiversity. They purify water by pumping it through their gills — obtaining food for themselves and filtering plankton and chemical contaminants from the water. Additionally, generations of oysters settle on top of each other, forming a reef of sorts, that provides habitats for other marine organisms, such as crabs.

While oyster reefs all over the world are threatened due to overfishing, in one protected location in northeast Scotland, there is a university, a marine conservation society, and a whisky distillery all striving to restore a long-lost oyster habitat.

The Dornoch Firth is a large and complex estuary that has been designated both a Special Protection Area (SPA) and Special Area of Conservation (SAC), because it is one of the northern-most estuaries for migrating and wintering birds. White sandy beaches at the mouth of the estuary lead to salt marshes, sandflats and mudflats that support an incredibly diverse plant and animal environment, including breeding osprey, waders and wildfowl, seal, dolphins, otters and mussel reefs.

The Native European Oyster thrived in the waters of the Dornoch Firth for over 10,000 years, until it was overfished and disappeared completely just over 100 years ago.

Glenmorangie Distillery is picturesquely situated on the banks of Dornoch Firth, producing single malt whisky for the past 175 years. In 2014, Glenmorangie pioneered a project to reduce their environmental footprint. The whisky distillery teamed up with Heriot-Watt University and the Marine Conservation Society in a partnership called DEEP (Dornoch Environmental Enhancement Project) to make some drastic changes to the water quality, re-introduce the Native European Oyster, and set a precedent for the re-establishment of disappearing reefs worldwide.

The DEEP project had two major initiatives. First was the introduction of an anaerobic digestion plant that purifies 95% of the by-products of whisky distillation that were previously released into the Firth. The remaining 5% consists of mostly organic compounds, such as barley, which is used as food by the oysters that then go on to further improve the water quality in the firth.

Next was the re-establishment of the Native European Oyster reef. Following an initial small scale exploratory step, twenty tons of waste mussel and scallop shells were laid down to stabilize the sediment and give the oysters a surface on which to grow. Divers then laid down 20,000 Native European Oysters in an area that covers 40 hectares of the firth. Their hope is that the oyster population will grow to 200,000 in three years, and then four million after five years, generating a fully sustainable oyster reef.

Ultimately, the introduction of an anaerobic digestion plant and the re-introduction of the Native European Oyster reef will drastically improve the water quality in Dornoch Firth. In turn, this supports all of the other organisms that form the diverse marine ecosystem there. A picturesque landscape for a renowned distillery on the edge of an ecosystem teaming with bird, plant and marine life, is a victory for both whisky connoisseurs and for the conservation of biodiversity.

DNA Modifying Enzymes & Cloning Technologies



The trusted source for DNA-modifying enzymes & cloning technologies.

Molecular biology, which also includes synthetic biology, is a fundamental area of research and development. Central to these advances has been the use of DNA modifying enzymes and novel cloning technologies. Some common examples of DNA-modifying enzymes include kinases, ligases and methylases, while newer cloning technologies include NEBuilder® HiFi DNA Assembly and Golden Gate Assembly.

NEB's 40+ year history as a leader in enzyme technologies gives you confidence in the products and support you'll receive. NEB continues to serve the scientific community by providing the tools to carry out the most innovative research, from start to finish. All NEB products pass stringent quality control assays to ensure the highest level of functionality and purity.

NEB offers several online tools to aid in your cloning experiments, including:

- NEBcloner® find the right products and protocols for each step of your traditional cloning experiment, including double digests
- NEBioCalculator® use this tool for your scientific calculations and conversions
- NEBuilder Assembly Tool use this tool for help with your DNA assembly primer design
- Thermostable Ligase Reaction Temperature Calculator estimate incubation temperature when using thermostable ligases
- NEB Golden Gate Assembly Tool use this tool for help with construct design for Golden Gate Assembly
- Ligase Fidelity Viewer visualize overhang ligation preferences for Golden Gate Assembly design

To view the full list of online tools available, see page 288.

Featured Products

- NEBuilder® HiFi DNA Assembly
 Master Mix & Cloning Kit
- 91 NEB® PCR Cloning Kit (with & without competent cells)
- 90 NEB Golden Gate Assembly Kit (Bsal-HF°v2)
- 97 HiFi Taq DNA Ligase
- **99** Quick Dephosphorylation Kit

Featured Tools & Resources

- 86 Cloning Overview
- 349 Troubleshooting Guide for Cloning
- Visit ClonewithNEB.com to view our online tutorials explaining each of the steps in the cloning workflow.




	Cloning & Synthetic Biology			Nuclease BAL-31	106
	Cloning Workflow Comparison	86		Nuclease P1	106
RR	NEBuilder HiFi DNA Assembly Master Mix	88		RecJ _f	106
RR	NEBuilder HiFi DNA Assembly Cloning Kit	88		T5 Exonuclease	107
RR	NEBuilder HiFi DNA Assembly Bundle	00	Rii	T7 Exonuclease	107
	for Large Fragments	88		Nucleoside Digestion Mix	202
	Gibson Assembly Cloping Kit	89 89		DNA Repair Proteins	
	Gibson Assembly Cloning Kit NEB Golden Gate Assembly Kit (Bsal-HFv2)	90		Properties of DNA Repair Enzymes	107
	BioBrick Assembly Kit	91		DNA Repair Glycosylases on Various	
	NEB PCR Cloning Kit	31		Damaged Bases	108
	(with and without competent cells)	91		APE 1 Endonuclease IV	109
RR	USER & Thermolabile USER II Enzymes	91		Tth Endonuclease IV	109 109
RR	Q5 Site Directed Mutagenesis Kit			Endonuclease III (Nth)	109
	(with and without competent cells)	92		Tma Endonuclease III	110
RR	Quick Blunting Kit	92		Endonuclease V	110
	DNA Ligases			Endonuclease VIII	110
	DNA Ligase Selection Chart	93	Rit	Fpg	110
	Substrate-based Ligase Selection Chart	94	Rit	hAAG	111
RR	T4 DNA Ligase	95	Rii	T4 PDG (T4 Endonuclease V)	111
RR	Quick Ligation Kit	95	Rit	Uracil-DNA Glycosylase (UDG)	111
RR	Instant Sticky-end Ligase Master Mix	95	Rii	Afu Uracil-DNA Glycosylase (UDG)	111
	Blunt T/A Ligase Master Mix	96		Antarctic Thermolabile UDG	112
	ElectroLigase	96		hSMUG1	112
	T3 DNA Ligase	96	Rii	PreCR Repair Mix	76
	T7 DNA Ligase	96		Other Enzymes	
	HiFi Taq DNA Ligase Tag DNA Ligase	97 97	RR	β-Agarase I	112
m	E. coli DNA Ligase	98	Rit	Cre Recombinase	113
RR	9°N DNA Ligase	98	RR	T7 Endonuclease I	113
	SplintR Ligase	199	Rii	TelN Protelomerase	113
			Rii	Topoisomerase I (E. coli)	113
	DNA Kinases			DNA Methyltransferases	
	T4 Polynucleotide Kinase	98	RR	CpG Methyltransferase (M.SssI)	114
	T4 Polynucleotide Kinase (3' phosphatase minus)	98 99		GpC Methyltransferase (M.CviPI)	114
m	5-hydroxymethyluridine Kinase	99		Alul Methyltransferase	114
	DNA Phosphatases & Sulfurylases		Rit	BamHI Methyltransferase	114
	Phosphatase Selection Chart	99		dam Methyltransferase	115
RR	Quick Dephosphorylation Kit	99	Rit	EcoGII Methyltransferase	115
RR	Shrimp Alkaline Phosphatase (rSAP)	100	Rit	EcoRI Methyltransferase	115
RR	Antarctic Phosphatase	100		HaeIII Methyltransferase	115
_	Alkaline Phosphatase (CIP)	100		Hhal Methyltransferase	115
	Pyrophosphatase, Inorganic (<i>E. coli</i> , yeast)	189 189		Hpall Methyltransferase	116
	Thermostable Inorganic Pyrophosphatase Apyrase	101		Mspl Methyltransferase	116
	Tte UvrD Helicase	101		Taql Methyltransferase Human Dnmt1 Methyltransferase	116 270
	710 OVID Hollouse	101		Tiutilali Dililiti Wettiyitialisielase	210
	Nucleases			Single-stranded DNA Binding Proteins	
	Properties of Exonucleases and Endonucleases	102	Rit	RecA, RecA _f	117
		103		T4 Gene 32 Protein	117
	DNase I (RNase-Free) Exonuclease I	202 103	Rii	ET SSB (Extreme Thermostable ssBinding Protein	າ) 117
	Thermolabile Exonuclease I	103		Cloning Plasmids & DNAs	
	Exonuclease III	103	R#	Cloning Plasmids	
	Exonuclease V (RecBCD)	104		(pBR322, pUC19)	118
	Exonuclease VII	104	RR	DNAs	
	Exonuclease VIII, truncated	104		(M13mp18 DNA, Lambda, Lambda	140
RR	Exonuclease T	105		(N ⁶ -Methyladenine-free), φX174)	118
RR	Thermostable FEN1	105		Other	
	Lambda Exonuclease	105		M13K07 Helper Phage	118
RR	Micrococcal Nuclease	105		Competent Cell Selection Chart for Cloning	119
	Mung Bean Nuclease	106		Monarch Nucleic Acid Purification Kits	119

Recombinant Enzyme

One or more of these products are covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at gbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications. Your purchase, acceptance, and/or payment of and for NBB's products is pursuant to NBB's Terms of Sale at https://www.neb.com/support/terms-of-sale. NBB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NBB.


Cloning Workflow Comparison

The figure below compares the steps for the various cloning methodologies, along with the time needed for each step in the workflows.





For help with choosing the right product for each step in the cloning workflow, try NEBcloner at **NEBcloner.neb.com**



RE Digest

Colony PCR

Sequencing

Site-Directed

Mutagenesis

RX

NEBuilder® HiFi DNA Assembly Master Mix & Cloning Kit

NEBuilder HiFi DNA Assembly Master Mix #E2621S 10 reactions162 € #E2621L 50 reactions648 € 250 reactions 2589 € #E2621X

NEBuilder HiFi DNA Assembly

Cloning Kit

#E5520S 10 reactions 190 €

NEBuilder HiFi DNA Assembly Bundle for Large Fragments

#E2623S 20 reactions 498 €

- Simple and fast seamless cloning
- Increased number of successful assembly products, particularly for longer or greater numbers of fragments
- Flexible sequence design, with no need to engineer cloning site
- Complex assembly achieved in an hour
- Less screening/re-sequencing of constructs, virtually error-free, highfidelity assembly
- Use in successive rounds of assembly; removes 5´ and 3´ restriction enzyme mismatches
- Bridge two ds-fragments with a synthetic ss-DNA oligo
- DNA can be used immediately for transformation or as template for PCR or RCA
- Switch from other systems easily; compatible with Gibson Assembly-. In-Fusion- designed fragments
- Adapts for multiple DNA manipulations, including site-directed mutagenesis.
- Assemble multiple DNA fragments and transform in just under 2 hours
- Clone into any vector with no additional sequence added (scarless)

To learn how simple NEBuilder HiFi is, visit NEBuilderHiFi.com **Description:** NEBuilder HiFi DNA Assembly Master Mix was developed to improve the efficiency and accuracy of DNA assembly. This method allows for seamless assembly of multiple DNA fragments, regardless of fragment length or end compatibility. This method has been used to assemble either single-stranded oligonucleotides or different sizes of DNA fragments with varied overlaps (15-80 bp). It has utility for the synthetic biology community, as well as those interested in one-step cloning of multiple fragments due to its ease of use, flexibility and simple master-mix format. The reaction features different enzymes that perform in the same buffer:

- Exonuclease creates single-stranded 3´ overhangs that facilitate the annealing of fragments that share complementarity at one end (the overlap region)
- The polymerase fills in gaps within each annealed fragment
- . The DNA ligase seals nicks in the assembled DNA

The end result is a double-stranded, fully sealed DNA molecule that can serve as template for PCR, RCA or a variety of other molecular biology applications, including direct transformation of E. coli.

The NEBuilder HiFi DNA Assembly Cloning Kit combines the power of the NEBuilder HiFi DNA Assembly Master Mix with NEB 5-alpha Competent E. coli, enabling fragment assembly and transformation in just under 2 hours.

NEBuilder HiFi Kits can be purchased with NEB 5-alpha Competent E. coli (Cloning Kit, NEB #E5520) or as a bundle with NEB 10-beta Competent E. coli (Bundle for Large Fragments, NEB #E2623). NEB 5-alpha competent cells are excellent for routine assemblies of 15 kb or less. NEB recommends NEB 10-beta competent cells for assemblies larger than 15 kb.

The NEBuilder HiFi DNA Assembly **Master Mix Includes:**

- NEBuilder HiFi DNA Assembly Master Mix
- NEBuilder Positive Control

The NEBuilder HiFi DNA Assembly Cloning Kit Includes:


- NEBuilder HiFi DNA Assembly Master Mix
- NEBuilder Positive Control
- NEB 5-alpha Competent E. coli (High Efficiency)
- SOC Outgrowth Medium
- pUC19 Control DNA

The NEBuilder HiFi DNA Assembly Bundle for Large Fragments Includes:


- NEBuilder HiFi DNA Assembly Master Mix
- NEBuilder Positive Control
- NEB 10-beta Competent E. coli (High Efficiency)
- NEB 10-beta/Stable Outgrowth Medium
- pUC19 Control DNA



Speed up your experimental design with our primer design tool at NEBuilder.neb.com



Overview of the NEBuilder HiFi DNA Assembly Cloning Method.













Gibson Assembly® Master Mix & Cloning Kit

Gibson Assembly Master Mix

#E2611S 10 reactions 164 € #E2611L 50 reactions 656 €

Gibson Assembly Cloning Kit

#E5510S 10 reactions 191 €

- Increased number of successful assembly products, particularly for longer or greater numbers of fragments
- Flexible sequence design with no need to engineer cloning sites
- Assemble multiple DNA fragments and transform in just under 2 hours
- Clone into any vector with no additional sequence added (scarless)
- No PCR clean-up step required

Description: Gibson Assembly Master Mix was developed by Dr. Daniel Gibson and his colleagues at the J. Craig Venter Institute and licensed to NEB by Synthetic Genomics, Inc. It allows for successful assembly of multiple DNA fragments, regardless of fragment length or end compatibility.

Gibson Assembly efficiently joins multiple overlapping DNA fragments in a single-tube isothermal reaction. The Gibson Assembly Master Mix includes three different enzymatic activities that perform in a single buffer:

- The exonuclease creates a single-stranded 3 overhang that facilitates the annealing of fragments that share complementarity at one end
- The polymerase fills in gaps within each annealed fragment
- . The DNA ligase seals nicks in the assembled DNA

The Gibson Assembly Cloning Kit combines the power of the Gibson Assembly Master Mix with NEB 5-alpha Competent E. coli, enabling fragment assembly and transformation in just under 2 hours.

The Gibson Assembly Cloning Kit has been optimized for the assembly and cloning of up to 6 fragments.

R₩

The Gibson Assembly Master Mix Includes:

- Gibson Assembly Master Mix
- NEBuilder Positive Control

The Gibson Assembly Cloning Kit Includes:

- Gibson Assembly Master Mix
- NEBuilder Positive Control
- NEB 5-alpha Competent E. coli (High Efficiency)
- SOC Outgrowth Medium
- pUC19 Control DNA

GIBSON ASSEMBLY® is a registered trademark of Synthetic Genomics Inc



Speed up your experimental design with our primer design tool at NEBuilder.neb.com

To learn how simple Gibson Assembly is, view our online tutorials at NEBGibson.com

NEB® Golden Gate

Synthetic Biology/DNA Assembly Selection Chart

	DNA Assembly	Assembly	Assembly Kit	USER® Enzyme
	(NEB #E2621) (NEB #E5520) (NEB #E2623)	(NEB #E5510) (NEB #E2611)	(Bsal-HFv2) (NEB #E1601)	(NEB #M5505)
PROPERTIES				
Removes 5´ or 3´ End Mismatches	***	*	N/A	N/A
Assembles with High Fidelity at Junctions	***	**	***	***
Tolerates Repetitive Sequences at Ends	*	*	***	***
Generates Fully Ligated Product	***	***	***	NR
Joins dsDNA with Single-stranded Oligo	***	**	NR	NR
Assembles with High Efficiency with Low Amounts of DNA	***	**	**	**
Accommodates Flexible Overlap Lengths	***	***	*	**
APPLICATIONS				
Simple Cloning (1-2 Fragments)	***	***	***	***
4-6 Fragment Assembly (one pot)	***	***	***	***
7-11 Fragment Assembly (one pot)	***	**	***	***
12-24 Fragment Assembly (one pot) (1)	*	*	***	NR
Template Construction for In vitro Transcription	***	***	***	***
Synthetic Whole Genome Assembly	***	*	*	*
Multiple Site-directed Mutagenesis	***	**	**	**
Library Generation	***	***	***	**
Metabolic Pathway Engineering	***	**	***	***
TALENS	**	**	***	**
Short Hairpin RNA Cloning (shRNA)	***	**	*	*
gRNA Library Generation	***	**	*	*
Large Fragment (> 10 kb) Assembly	***	***	***	**
Small Fragment (< 100 bp) Assembly	***	*	***	***
Use in Successive Rounds of Restriction Enzyme Assembly	***	*	NR	*

NEBuilder® HiFi NEB Gibson®

- Optimal, recommended product for selected application
 - Works well for selected application Will perform for selected application, but is not recommended
 - N/A Not applicable to this application
 - NR Not recommended
- (1) Please visit www.neb.com/GoldenGate for more information.







NEB Golden Gate Assembly Kit (BsaI-HF[®]v2)

#F1601S 20 reactions 158 € #E1601L 100 reactions 421 €

Companion Products:

NEB 5-	alpha	Competent	Е.	coli
(Hinh F	fficier	ncv)		

20 x 0.05 ml197 € #C2987H #C29871 6 x 0.2 ml 153 € #C2987P 1 x 96 well plate 493 € #C2987R 1 x 384 well plate1096 € #C2987U 96 x 50 µl/tube742 €

NFB 10-beta Competent F. coli.

(High Efficiency)

#C3019H 20 x 0.05 ml231 € #C3019I 6 x 0.2 ml179 €

NEB Cloning Competent E. coli Sampler #C1010S 8 tubes 120 €

Q5 Hot Start High-Fidelity 2X Master Mix #M0494S 100 rxns (50 µl vol) 180 € 500 rxns (50 µl vol)720 €

- Updated to include Bsal-HFv2 (optimized for Golden Gate)
- Seamless cloning no scar remains following assembly
- Includes destination plasmid with T7/ SP6 promoters
- Ordered assembly of multiple fragments (2-20+) in a single reaction
- Can also be used for cloning of single inserts and library preparations
- Efficient with regions of high GC content and areas of repeats
- Compatible with a broad range of fragment sizes (< 100 bps to > 15 kb)

Type IIS Enzymes used in Golden Gate:

- Bsal (NEB #R0535)
- Bsal-HFv2 (NEB #R3733)
- Bbsl (NEB #R0539)
- BbsI-HF (NEB #R3539)
- BsmBI (NEB #R0580)
- Esp3I (NEB #R0734)

Description: The NEB Golden Gate Assembly Kit (Bsal-HFv2) contains an optimized mix of Bsal-HFv2 and T4 DNA Ligase. Together these enzymes can direct the assembly of multiple inserts/modules using the Golden Gate approach. Also included is the pGGA destination plasmid, which provides a backbone for your assembly, features convenient restriction enzyme sites for subcloning, and has T7/SP6 promoter sequences to enable in vitro transcription.

The efficient and seamless assembly of DNA fragments, commonly referred to as Golden Gate assembly, has its origins in 1996, when for the first time it was shown that multiple inserts could be assembled into a vector backbone using only the sequential or simultaneous activities of a single Type IIS restriction enzyme and T4 DNA Ligase.

Type IIS restriction enzymes bind to their recognition sites but cut the DNA downstream from that site at a positional, not sequence-specific, cut site. Thus, a single Type IIS restriction enzyme can be used to generate DNA fragments with unique overhangs. As an example, Bsal has a recognition site of GGTCTC(N1/N5), where the GGTCTC represents the recognition/binding site, and the N1/N5 indicates the cut site is one base downstream on the top strand, and five bases downstream on the bottom strand. Assembly of digested fragments proceeds through annealing of complementary four base overhangs on adjacent fragments. The digested fragments and the final assembly no longer contain Type IIS restriction enzyme recognition sites, so no further cutting is possible. The assembly product accumulates with time.

RX

While particularly useful for multi-fragment assemblies such as Transcription Activator Like Effectors (TALEs) and TALEs fused to a Fokl nuclease catalytic domain (TALENs), the Golden Gate method can also be used for cloning of single inserts and inserts from diverse populations that enable library creation.

Advances in Ligase Fidelity: Research at NEB has led to increased understanding of ligase fidelity, including the development of a comprehensive method for profiling end-joining ligation fidelity in order to predict which overhangs have improved fidelity. This research allows careful choice of overhang sets, which is especially important for achieving complex Golden Gate Assemblies. To learn more, visit www.neb.com/goldengate.

The NEB Golden Gate Assembly Kit (Bsal-HFv2) Includes:


- NEB Golden Gate Assembly Mix
- T4 DNA Ligase Reaction Buffer (10X)
- pGGA Destination Plasmid

NEB Golden Gate Assembly Tool

Speed up your experimental design with our assembly tool at GoldenGate.neb.com

Ligase Fidelity Viewer*

Visualize overhang ligation preferences for Golden Gate Assembly design. *This is a NEBeta™ Tool which is not yet optimized for design and usability. Please share your feedback with us. Access this tool at www.neb.com/research/nebeta-tools.



In its simplest form, Golden Gate Assembly requires a Type IIS recognition site, in this case, Bsal (GGTCTC), added to both ends of a dsDNA fragment. After digestion, these sites are left behind, with each fragment bearing the designed 4-base overhangs that direct the assembly.













BioBrick® Assembly Kit

#E0546S 50 reactions 276 €

The BioBrick Assembly Kit was developed in partnership with Ginkgo BioWorks. For more details and for technical questions, please see: ginkgobioworks.com/support

Description: The BioBrick Assembly Kit provides a streamlined method for assembly of BioBrick parts into multi-component genetic systems. BioBrick parts are DNA sequences that encode a defined biological function and can be readily assembled with any other BioBrick part. The process for assembling any two BioBrick parts is identical and results in a new composite BioBrick part.

Please refer to the individual datacards for each reagent's recommended use and storage conditions.

BIOBRICK® is a registered trademark of The BioBricks Foundation (http://www.biobricks.org)

RX

The BioBrick Assembly Kit Includes:

- EcoRI-HF
- Xbal
- Spel
- Pstl
- 10X NEBuffer 2.1
- T4 DNA Ligase
- 10X T4 DNA Ligase Buffer

NEB® PCR Cloning Kit (with or without competent cells)

NEB PCR Cloning Kit

#E1202S 20 reactions 295 €

NEB PCR Cloning Kit (without competent cells)

#E1203S 20 reactions 145 €

- In vitro transcription with both SP6 & T7 promoters
- Easy cloning of all PCR products, including blunt and TA ends
- Fast cloning with 5-minute ligation step
- Simplified screening with low/no colony background and no blue/white selection
- Save time by eliminating purification steps
- More flanking restriction sites available for easy subcloning, including choice of two single digest options
- Bsal site removed to allow cloning of Golden Gate modules

Description: The NEB PCR Cloning Kit contains an optimized 2X Cloning Master Mix with a proprietary ligation enhancer and a linearized vector that uses a novel mechanism for background colony suppression to give a low background. It allows simple and quick cloning of any PCR amplicon, whether the amplification reactions are performed with proofreading DNA polymerases, such as Q5® or Phusion® which produce blunt ends; or nonproofreading DNA polymerases, such as Tag or Tag mixes (One Taq, Long Amp Taq) which produce singlebase overhangs. This is possible due to "invisible" end polishing components in the master mix that are active during the ligation step. The kit also allows direct cloning from amplification reactions without purification, and works well whether or not the primers used in the PCR possess 5´-phosphate groups.

- Provided analysis primers allow for downstream colony PCR screening or sequencing
- Ready-to-use kit components include 1 kb control amplicon, linearized cloning vector and single-use competent E. coli (NEB #E1202 only)
- Longer shelf life (12 months), as compared to some commercially available products

R₩

The PCR Cloning Kit Includes:

- Linearized pMiniT™ 2.0 Vector
- Cloning Mix 1
- Cloning Mix 2
- Amplicon Cloning Control (1 kb)
- Cloning Analysis Forward Primer
- Cloning Analysis Reverse Primer
- NEB 10-beta Competent E. coli (Cloning Efficiency) (NEB #E1202 only)
- NEB 10-beta/Stable Outgrowth Medium (NEB #E1202 only)
- pUC19 Control DNA

PHUSION® is a registered trademark of Thermo Fisher Scientific

USER® Enzyme

Thermolabile USER II Enzyme

USER Enzyme

#M5505S 50 units73 € #M5505L 250 units292 €

NEW

Thermolabile USER II Enzyme

#M5508S 50 units98 € #M5508L 250 units392 €

- USER Cloning
- Directional RNA-Seq
- NEBNext adaptor cleavage

Description: USER (Uracil-Specific Excision Reagent) Enzyme generates a single nucleotide gap at the location of a uracil. USER Enzyme is a mixture of Uracil DNA glycosylase (UDG) and the DNA glycosylase-lyase Endonuclease VIII. UDG catalyses the excision of a uracil base, forming an abasic (apyrimidinic) site while leaving the phosphodiester backbone intact. The lyase activity of Endonuclease VIII breaks the phosphodiester backbone at the 3' and 5' sides of the abasic site so that base-free deoxyribose is released.

Thermolabile USER (Uracil-Specific Excision Reagent) II Enzyme generates a single nucleotide gap at the location of a uracil residue. It can be 100% inactivated at temperatures > 65°C

CutSmart Ril 37° Kill

CutSmart RN 37° Wh

Reaction Conditions: CutSmart Reaction Buffer. Incubate at 37°C. Heat Inactivation of Thermolabile USER II Enzyme: 65°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to nick 10 pmol of a 34-mer oligonucleotide duplex containing a single uracil base, in 15 minutes at 37°C in a total reaction volume of 10 μl. Unit assay conditions can be found at www.neb.com.

Concentration: 1,000 units/ml









Q5[®] Site Directed Mutagenesis Kit (with or without competent cells)

RX

Q5 Site Directed Mutagenesis Kit

#E0554S 10 reactions 192 €

Q5 Site Directed Mutagenesis Kit (without competent cells)

#E0552S 10 reactions 132 €

Companion Products:

NEB PCR Cloning Kit #E1202S 20 reactions 295 €

NEB PCR Cloning Kit (without competent cells)

#E1203S 20 reactions 145 €

KLD Enzyme Mix

#M0554S 25 reactions 302 €

- Robust exponential amplification generates high yields of desired mutations from a wide range of templates.
- Low error rate of Q5 High-Fidelity DNA Polymerase reduces screening time.
- Room temperature reaction setup
- Use of standard primers eliminates need for phosphorylated or purified oligos
- Easy-to-use master mix format

For help with primer design, try NEBaseChanger at NEBaseChanger.neb.com


Description: The Q5 Site-Directed Mutagenesis Kit allows rapid site-specific mutagenesis of double-stranded plasmid DNA in less than 2 hours. The kit utilizes the robust Q5 Hot Start High-Fidelity DNA Polymerase along with custom mutagenic primers to create substitutions, deletions and insertions in a wide variety of plasmids. Transformation into high-efficiency NEB 5-alpha Competent E. coli, provided with (NEB #E0554), ensures robust results with plasmids up to 14.3 kb in length.

Applications:

· Generation of mutations, insertions or deletions in plasmid DNA

The Q5 Site Directed Mutagenesis Kit Includes:

- Q5 Hot Start High-Fidelity Master Mix (2X)
- KLD Enzyme Mix (10X) and Reaction Buffer (2X)
- Control Primer Mix (10 µM each) and Template DNA (5 ng/µl)
- NEB 5-alpha Competent E. coli (High Efficiency) (NEB #E0554 only)
- pUC19 Transformation Control Plasmid (5 pg/μl)
- SOC Outgrowth Medium



Q5 Site-Directed Mutagenesis Kit Overview. The first step is an exponential amplification using Q5 Hot Start High-Fidelity DNA Polymerase. The second step is a unique enzyme mix containing a kinase, ligase and Dpnl. Together, these enzymes allow for rapid circularization of the PCR product and removal of the template DNA. The last step is a high-efficiency transformation into chemically competent cells.

Quick Blunting™ Kit

#F1201S 20 reactions85 € #E1201L 100 reactions 340 €

Special Offer:

Quick Blunting and Quick Ligation Kits #E0542S 20 reactions..... 161 € #E0542L 100 reactions 644 €

See page 95 for details on the Quick Ligation Kit.

- Restriction enzyme-digested DNA is blunted in less than 30 minutes
- Reactions are performed at room temperature in a ready-to-use mix
- Suitable for restriction enzyme-digested DNA, sheared or nebulized DNA or PCR product

Description: The Quick Blunting Kit is used to convert DNA with incompatible 5 or 3 overhangs to 5 phosphorylated, blunt-ended DNA for efficient blunt-end ligation into DNA cloning vectors. DNA is blunted using T4 DNA Polymerase (NEB #M0203) which has both $3 \rightarrow 5$ exonuclease activity and 5´→ 3´ polymerase activity. T4 Polynucleotide Kinase (NEB #M0201) is included in the enzyme mix for phosphorylation of the 5'ends of bluntended DNA for subsequent ligation into a cloning vector. This kit is optimized for blunting up to 5 ug of DNA in a single reaction.

Applications:

- · Prepare sheared, nebulized or restriction enzymedigested DNA for blunt-end ligation into a plasmid, cosmid. fosmid or BAC vector
- · Prepare PCR products for efficient blunt-end cloning


The Quick Blunting Kit Includes:

- Blunting Enzyme Mix
- 10X Blunting Buffer
- 1 mM Deoxynucleotide (dNTP) Solution Mix


RX W

Heat Inactivation: 70°C for 10 minutes.

Notes: PCR generated DNA must be purified before blunting by using a commercial purification kit, such as Monarch PCR & DNA Cleanup Kit (NEB #T1030) phenol extraction/ethanol precipitation or gel electrophoresis. Restriction enzyme digested DNA can be blunted directly without purification.













DNA Ligase Selection Chart

NEB offers a variety of ligases for DNA research. Many of these enzymes are recombinant, and all offer the quality and value you have come to expect from our products. While more than one ligase may work for your application, the following selection chart presents our recommendations for optimal performance.

Visit NEBStickTogether.com for more information on DNA Ligases.

DNA APPLICATIONS	Blunt/TA Ligase Master Mix	Instant Sticky-end Ligase Master Mix	Electro Ligase™	T4 DNA Ligase	Quick Ligation™ Kit	T3 DNA Ligase	T7 DNA Ligase	HiFi <i>Taq</i> DNA Ligase	<i>E. coli</i> DNA Ligase	<i>Taq</i> DNA Ligase	9°N™ DNA Ligase	NEBNext Ultra II Ligation Module	SplintR® Ligase
Ligation of sticky ends	**	***	**	**	***	**	**	*	*	*	*		
Ligation of blunt ends	***	*	**	**	***	**							
T/A cloning	***	*	**	**	**	*	*						
Electroporation			***	**									
Ligation of sticky ends only							***						
Repair of nicks in dsDNA	**	**	**	***	**	**	**	**	**	**	**		**
High complexity library cloning	**	**	**	***	**								
Adaptor Ligation	***	**	**	*	**	*						A	
Ligation-Dependent DNA Sequence & SNP Detection (LCR, LDR & related methods)								***		**	**		
Ligation-Dependent RNA Sequence & SNP Detection				*									***
Ligation of adjacent ssDNAs on an RNA Splint													***
NGS APPLICATIONS		'											
NGS Library Prep dsDNA-dsDNA (ligation)	A			A		A						A	
FEATURES													
Salt tolerance (> 2X that of T4 DNA Ligase)						~							
Ligation in 15 min. or less	~	~		V	~	~	V	~		V	~	~	~
Master Mix Formulation	~	~										~	
Thermostable								~		V	~		

- $\star\star\star$ Optimal, recommended ligase for selected application
 - $\bigstar \bigstar$ Works well for selected application
 - \star Will perform selected application, but is not recommended
 - ▲ Please consult the specific NGS protocol to determine the optimal enzyme for your needs

Helpful Online Tools:

Recombinant

For help with choosing the best ligase for your experiment, try **NEBcloner** at **NEBcloner.neb.com**

For help with estimating incubation temperature when using thermostable ligases, try the Thermostable Ligase Reaction Temperature Calculator at LigaseCalc.neb.com

For help with scientific calculations and conversions, try **NEBioCalculator** at **NEBioCalculator.neb.com**

Find an overview of ligation.







Substrate-based Ligase Selection Chart

DNA ∼ This chart provides our recommendation for a choice of ligase to use in a reaction, based upon the type of substrate present. DNA, RNA and hybrid substrates are represented and require specific enzymes to achieve the highest-efficiency ligation.

	RECOMMENDED LIGASE	COMMENTS
NICKED DNA/RNA		
5' OH P 3'	T4 RNA Ligase 2	
5' ~~~~OH p~~~~ 3'	T4 RNA Ligase 2	
5' OH p 3'	T4 RNA Ligase 2	
5' OH p 3'	T4 DNA Ligase	
5' ——OH p 3'	N/A	No ligase optimized for this activity
5' ——OH p	N/A	No ligase optimized for this activity
2, OH b 3,	SplintR Ligase	100 – 1,000-fold higher efficiency than T4 DNA Ligase
5' — OH p — 3'	T4 DNA Ligase	For high temperatures, we recommend <i>Taq</i> DNA Ligase. For highest fidelity, we recommend HiFi <i>Taq</i> DNA Ligase.
ssDNA/RNA		
5° ——— OH p——— 3°	N/A	See CircLigase™
5' — OH p 3'	N/A	No ligase optimized for this activity
5' ~~~~ OH p——— 3'	T4 RNA Ligase 1	Supplement with ATP
5' OH p OM 3'	T4 RNA Ligase 1	
5' ~~~ OH App ——— 3'	T4 RNA Ligase 2 Truncated KQ	
5' OH App 3'	T4 RNA Ligase 2 Truncated KQ	
5' ——— OH App ——— 3'	Thermostable 5´ App DNA/RNA Ligase	We recommend a Proteinase K cleanup
5' — OH App ~ ~ 3'	Thermostable 5´ App DNA/RNA Ligase	We recommend a Proteinase K cleanup
5' — p OH ~ 3'	RtcB Ligase	Supplement with GTP and Mn ²⁺
2, OH 3,	RtcB Ligase	Supplement with GTP and Mn ²⁺
5' ~~~~~ OH pNp 3'	T4 RNA Ligase 1	
5' ~~~~~ OH pdNp 3'	T4 RNA Ligase 1	Reported to work, but ligates inefficiently. Consider pdCp.
5' ——— OH pNp 3'	T4 RNA Ligase 1	
5' ———— OH pdNp 3'	T4 RNA Ligase 1	
dsDNA/RNA		
5' ————————————————————————————————————	Blunt T/A Ligase Master Mix	
5' ————————————————————————————————————	Blunt T/A Ligase Master Mix	
3, ————————————————————————————————————	Quick Ligation Kit or Instant Sticky-end Ligase Master Mix	For ligating ends under high salt conditions, we recommend T3 DNA Ligase. For ligation of cohesive ends ONLY, we recommend T7 DNA Ligase.
2,b OH 2, 2,OH b 3,	Quick Ligation Kit or Instant Sticky-end Ligase Master Mix	For ligating ends under high salt conditions, we recommend T3 DNA Ligase. For ligation of cohesive ends ONLY, we recommend T7 DNA Ligase.
2, b OH b 2, 2, 2, 3, 4, 5, 6, 6, 6, 6, 6, 6, 6, 6, 6, 6, 6, 6, 6,	T4 RNA Ligase 2	
5' OH POHO 3'	T4 RNA Ligase 2	

CIRCLIGASE™ is a trademark of EpiCentre Technologies Corp.

















T4 DNA Ligase

Regular Concentration (400,000 cohesive end units/ml)

#M0202S 20,000 units 68 € #M0202L 100,000 units 272 €

High Concentration

(2.000.000 cohesive end units/ml)

#M0202T 20,000 units68 € #M0202M 100,000 units272 €

- Efficient ligation of sticky or blunt ends
- T/A Cloning
- Repair of nicks in dsDNA
- High complexity library cloning
- Ligation of RNA to DNA

For help with molar ratio calculations, try **NEBioCalculator** at **NEBioCalculator.neb.com**

Competitive Nuclease Contamination Study.

T4 DNA Ligase from multiple suppliers was tested in reactions containing a fluorescent labeled single stranded, double stranded blunt, 3 overhang or 5 overhang containing oligonucleotides. The percent degradation by contaminating nucleases is determined by capillary electrophoresis and peak analysis. The resolution is at the single nucleotide level.

RX NEBU 16° K

Description: Catalyzes the formation of a phosphodiester bond between juxtaposed 5 ´phosphate and 3´ hydroxyl termini in duplex DNA or RNA. This enzyme will join blunt end and cohesive end termini, as well as repair single stranded nicks in duplex DNA, RNA or DNA/RNA hybrids.

Reaction Conditions: 1X T4 DNA Ligase Reaction Buffer. Recommended DNA concentration (0.1 to 1 μ M of 5' termini). Incubate at 16°C. Heat inactivation: 65°C for 10 minutes.

Unit Definition (Cohesive End Unit):


One unit is defined as the amount of enzyme required to give 50% ligation of HindIII fragments of λ DNA (5´ DNA

termini concentration of 0.12 μM [300 $\mu g/ml$]) in 20 μl of 1X T4 DNA Ligase Reaction Buffer in 30 minutes at 16°C.

Concentration: 400,000 and 2,000,000 cohesive end

Usage Notes: ATP is an essential cofactor for the reaction. This contrasts with *E. coli* DNA Ligase which requires NAD*.

Ligation can also be performed in any of the four restriction endonuclease NEBuffers or in T4 Polynucleotide Kinase Buffer if they are supplemented with 1 mM ATP.



Quick Ligation™ Kit

#M2200S 30 reactions99 € #M2200L 150 reactions396 €

Special Offer:

Quick Blunting and Quick Ligation Kits
#E0542S 20 reactions.....161 €
#E0542L 100 reactions.....644 €

See page 106 for details on the Quick Blunting Kit.

- 5 minute ligation reactions for sticky or blunt ends
- Reactions performed at room temperature
- T/A Cloning
- Repair of nicks in dsDNA

Description: The Quick Ligation Kit enables ligation of cohesive or blunt end DNA fragments in 5 minutes at room temperature (25°C).

The Quick Ligation Kit Includes:

- Quick T4 DNA Ligase (recombinant)
- 2X Quick Ligation Reaction Buffer

Reaction Conditions: 1X Quick Ligation Reaction Buffer. Incubate at room temperature (25°C).

R% 25° ₩₩

Notes:

Most ligations performed using the Quick Ligation Kit reach an end point at 5 minutes or less at 25°C. Incubation beyond this time provides no additional benefit. Overnight incubations can result in lower transformation efficiencies.

For electroporation applications we recommend the use of ElectroLigase (NEB #M0369).

Instant Sticky-end Ligase Master Mix

#M0370S 50 reactions 100 € #M0370L 250 reactions 400 €

- Instant ligation of sticky ends
- Repair of nicks in dsDNA
- High-complexity library cloning

Description: Instant Sticky-end Ligase Master Mix is a ready-to-use 2X solution of T4 DNA ligase and a proprietary ligation enhancer in an optimized reaction buffer. It is specifically formulated to rapidly ligate cohesive-end (2–4 bp) substrates. No thawing of the master mix is required, as it maintains a liquid state during storage at –20°C, and no incubation time is necessary to achieve ligation efficiencies sufficient for successful cloning of sticky-end substrates.

RR 186

Reaction Conditions: 1X Instant Sticky-end Ligase Master Mix with DNA substrates in a 10 μ l reaction volume. A 10 μ l reaction contains 1,800 cohesive end units of T4 DNA Ligase.

Usage Note: Product maintains a liquid state at -20° C. Freeze-thaw testing at -70° C has confirmed that performance after 20 freeze/thaw cycles is close to that of the original mix.

Blunt/TA Ligase Master Mix

#M0367S 50 reactions 100 € #M0367L 250 reactions 400 €

- Ligation of blunt ends
- T/A Cloning
- Repair of nicks in dsDNA
- High-complexity library cloning
- Ligation of sticky ends

Description: Blunt/TA Ligase Master Mix is a ready-touse 2X solution of T4 DNA Ligase, proprietary ligation enhancer and optimized reaction buffer. This master mix is specifically formulated to improve ligation and transformation of both blunt-end and single-base overhang substrates. No thawing is necessary as it remains liquid during storage at -20°C.

RX 25° W

Reaction Conditions: 1X Blunt/TA Ligase Master Mix with DNA substrates in a 10 ul reaction volume incubated at 25°C. A 10 µl reaction contains 1,800 cohesive end units of T4 DNA Ligase.

Usage Note: Product maintains a liquid state at -20°C. Freeze-thaw testing at -70°C has confirmed that the performance is unchanged after 20 freeze/thaw cycles.

ElectroLigase®

#M0369S 50 reactions 121 €

- Ideal for transformation by electroporation
- Vector or library construction
- Linker ligation
- Fragment assembly
- T/A Cloning

Description: ElectroLigase combines T4 DNA ligase and an optimized, ready-to-use 2X reaction buffer containing a proprietary ligation enhancer and no PEG. This combination is specifically formulated to promote robust ligation of all types of DNA ends (blunt, sticky, TA). It is directly compatible, without desalting or purification, with electrocompetent cells used for transformation by electroporation.

RX NEBU 25° V654

Reaction Conditions: 1X ElectroLigase Reaction Buffer with DNA substrates and 1 µl ElectroLigase in an 11 µl reaction volume incubated at 25°C. Heat inactivation: 65°C for 15 minutes

Usage Note: Product maintains a liquid state at -20°C.

T3 DNA Ligase

#M0317S 100,000 units68 € #M0317L 750,000 units 272 €

- Ligation of sticky or blunt ends
- Increased salt tolerance
- Repair of nicks in dsDNA

Description: T3 DNA Ligase is an ATP-dependent ds-DNA ligase from bacteriophage T3. Cohesive ends, blunt ends, and nick sealing can all be efficiently catalyzed by T3 DNA Ligase. Blunt end ligation is enhanced by the addition of PEG 6000 to the reaction. T3 DNA Ligase exhibits a higher tolerance (2-fold) for NaCl in the reaction compared to T4 DNA Ligase, making the enzyme a versatile choice for in vitro molecular biology protocols requiring DNA ligase activity.

Reaction Conditions: 1X T3 DNA Ligase Reaction Buffer. Incubate at 25°C.

RX NEBU 25° 166

Unit Definition: One unit is defined as the amount of enzyme required to give 50% ligation of 100 ng HindIII fragments of λ DNA in a total reaction volume of 20 μ l in 1 minute at 25°C in 1X T3 DNA Ligase Reaction Buffer.

Concentration: 3,000,000 units/ml

Notes: ATP is an essential cofactor for the reaction.

T3 DNA Ligase is also active in buffers without PEG 6000, including T4 DNA Ligase Buffer for applications in which PEG 6000 is detrimental. Supplement with 1 mM ATP (final concentration). In these buffers T3 DNA Ligase exhibits ~10-fold reduction in activity. In applications where a high concentration of NaCl needs to be maintained, we suggest using a reaction buffer without PEG 6000.

T7 DNA Ligase

#M0318S 100,000 units68 € #M0318L 750,000 units 272 €

- Ligation of sticky ends only
- Repair of nicks in dsDNA

Description: T7 DNA Ligase is an ATP-dependent ds-DNA ligase from bacteriophage T7. It will catalyze the formation of a phosphodiester bond between adjacent 5' phosphate and 3' hydroxyl groups of duplex DNA. Cohesive end ligation and nick sealing can be efficiently catalyzed by T7 DNA Ligase. However, unlike T4 and T3 DNA Ligases, blunt end ligation is not efficiently catalyzed by T7 DNA Ligase, making it a good choice for applications in which blunt and sticky ends of DNA are present but only the sticky ends are to be joined.

Reaction Conditions: 1X T7 DNA Ligase Reaction Buffer. Incubate at 25°C.



Unit Definition: One unit is defined as the amount of enzyme required to give 50% ligation of 100 ng HindIII fragments of λ DNA in a total reaction volume of 20 μl in 30 minutes at 25°C in 1X T7 DNA Ligase Reaction Buffer.

Concentration: 3,000,000 units/ml

Notes: ATP is an essential cofactor for the reaction.

T7 DNA Ligase is also active in buffers without PEG 6000, including T4 DNA Ligase Buffer for applications in which PEG 6000 is detrimental. Supplement the reaction with 1 mM ATP (final concentration). Using these buffers, the activity of T7 DNA Ligase is reduced ~10-fold.















HiFi Taq DNA Ligase

#M0647S 50 reactions 158 €

- High fidelity, thermostable
- Repair of nicks in dsDNA
- Allele-specific gene detection using ligase-dependent methods, including the Ligase Chain Reaction (LCR) and Ligase Detection Reaction (LDR)
- Ligation of padlock probes


RX NEBU WW

Description: An optimized blend of a thermostable DNA Ligase and a proprietary additive, HiFi *Taq* DNA Ligase efficiently seals nicks in DNA with unmatched high fidelity. The formation of a phosphodiester bond between juxtaposed 5′ phosphate and 3′ hydroxyl termini of two adjacent oligonucleotides that are hybridized to a complementary target DNA is enhanced in the improved reaction buffer and mismatch ligation is dramatically reduced. The improved formulation allows higher

resolution discrimination between ligation donors and acceptors, enabling precise detection of SNPs and other allele variants. HiFi *Taq* DNA Ligase is active at elevated temperatures (37–75°C).

Reaction Conditions: 1X HiFi *Taq* DNA Ligase Reaction Buffer. Incubate at 25°C.

For help with calculating ligation temp, try our Thermostable Ligase Reaction Temperature Calculator at LigaseCalc.neb.com



HiFi Taq DNA Ligase displays increased fidelity. (A) Schematic of multiplexed substrate pools. Each substrate pool contained a single splint with a defined NN at the ligation junction (e.g., AA, AC, AG...) along with all four upstream probes and all four FAM-labeled downstream probes. Each probe that encodes the base at the ligation junction is of unique length allowing for separation and analysis by capillary electrophoresis. A total of 16 substrate pools were prepared, one for each unique splint. (B) Comparison of the ligation fidelity of Ampligase (Epicentre), Taq DNA Ligase and HiFi Taq DNA Ligase. Fidelity measurements were performed using 1 µl of ligase in a 50 µl reaction mixture in the supplied buffers at 1X concentration. Reactions were incubated 30 min at 55°C, using multiplexed substrate pools as outlined in (A). Rows represent a single template sequence, while columns indicate a particular ligation product resulting from a specific pair of probes ligating with the indicated bases at the ligation junction. A dot indicates detection of a product (see legend above). The diagonal from the top left to the bottom right represents Watson-Crick ligation products; all other spaces indicate mismatch ligation products. While Taq DNA Ligase and Ampligase perform similarly under these conditions, with a range of mismatch products detectable, HiFi Taq DNA Ligase shows dramatically fewer mismatch products while maintaining high yields.

Thermus aquaticus (Taq) DNA Ligase

#M0208S 2,000 units80 € #M0208L 10,000 units320 €

- Thermostable
- Repair of nicks in dsDNA
- Used in Gibson Assembly method
- Allele-specific gene detection using Ligase Detection Reaction and Ligase Chain Reaction

Description: *Taq* DNA Ligase is a thermostable ligase that catalyzes the formation of a phosphodiester bond between the 5´-phosphate and 3´-hydroxyl termini of two adjacent DNA strands. The strands to be ligated need to be hybridized and accurately paired, with no gap, to a complementary DNA strand; allowing resolution of single nucleotide variants. *Taq* DNA Ligase uses NAD as a cofactor and is active at elevated temperatures (37–75°C).

Reaction Conditions: 1X Taq DNA Ligase Reaction Buffer. Incubate at 45°C.

RX NEBU 45° 166

Requires NAD+ as a cofactor. NAD+ is supplied in the 10X Taq DNA Ligase Reaction Buffer; the buffer should be stored at –70°C to extend the half life of the NAD+ cofactor.

Unit Definition: (Cohesive End Unit)

One unit is defined as the amount of enzyme required to give 50% ligation of the 12-base pair cohesive ends of 1 μg of BstEll-digested λ DNA in a total reaction volume of 50 μl in 15 minutes at 45°C. Unit assay conditions can be found at www.neb.com.

Concentration: 40,000 units/ml

Notes: Will not ligate short 4-base overlaps (typical of restriction enzyme digests), while it efficiently ligates 12 base pair overlaps.

E. coli DNA Ligase

#M0205S 200 units62 € #M0205L 1,000 units 248 €

- Repair of nicks in dsDNA
- Okayama and Berg cDNA cloning

Description: E. coli DNA Ligase catalyzes the formation of a phosphodiester bond between the 5´-phosphate and 3´-hydroxyl of two adjacent DNA strands in duplex DNA with cohesive ends. It is not appreciably active on bluntended substrates. E. coli DNA Ligase uses NAD as a cofactor and can be heat-inactivated. E. coli DNA Ligase is active at a range of temperatures (4-37°C).

Reaction Conditions: 1X E. coli DNA Ligase Reaction Buffer. Optimal ligation occurs at 16°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to give 50% ligation of HindIII

RX NEBU 16° 165

fragments of λ DNA (5´ DNA termini concentration of 0.12 µM, 300 µg/ml) in a total reaction volume of 20 µl in 30 minutes at 16°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10.000 units/ml

Usage Notes:

Requires NAD+ (nicotinamide adenine dinucleotide) as a cofactor.

Ligation of blunt-ended fragments is extremely inefficient. For ligation of blunt-ended fragments we recommend the Blunt/TA Ligase Master Mix (NEB #M0367) or the Quick Ligation Kit (NEB #M2200).

9°N™ DNA Ligase

#M0238S 2,500 units80 €

- Repair of nicks in DNA while incubating at high temperatures
- Thermostable
- Allele-specific gene detection using Ligase Detection Reaction and Ligase Chain Reaction

Description: 9°N DNA Ligase is a thermostable ligase that catalyzes the formation of a phosphodiester bond between the 5´-phosphate and 3´-hydroxyl of two adjacent DNA strands that are hybridized and accurately paired, with no gap, to a complementary DNA strand. 9°N DNA Ligase uses ATP as a cofactor and it is active at elevated temperatures (45-70°C).

Reaction Conditions: 1X 9°N DNA Ligase Reaction Buffer, Incubate at 45°C.

RX NEBU 45° WW

Unit Definition: (Cohesive End Unit)

One unit is defined as the amount of enzyme required to give 50% ligation of the 12-base pair cohesive ends of 1 μg of BstEII-digested λ DNA in a total reaction volume of 50 µl in 15 minutes at 45°C. Unit assay conditions can be found at www.neb.com.

Concentration: 40,000 units/ml

Notes: Will not ligate short 4-base overlaps (typical of restriction enzyme digests), while it efficiently ligates 12 base pair overlaps.

SplintR® Ligase

#M0375S 1,250 units95 € #M0375L 6,250 units 428 €

The SplintR Ligase, also known as PBCV-1 DNA Ligase or Chlorella virus DNA Ligase, efficiently

catalyzes the ligation of adjacent, single-stranded DNA oligonucleotides splinted by a complementary

See page 199 for more information.

RNA strand.

RX 25° 1664

T4 Polynucleotide Kinase & T4 Polynucleotide Kinase (3' phosphatase minus) RX NEBU 37° KK

T4 Polynucleotide Kinase

#M0201S 500 units58 € #M0201L 2,500 units 232 €

T4 Polynucleotide Kinase (3' phosphatase minus)

#M0236S 200 units 111 € #M0236L 1,000 units 444 €

- 5´ phosphorylation of DNA/RNA for subsequent ligation
- End labeling DNA or RNA for probes and DNA sequencing
- Removal of 3 phosphoryl groups with T4 Polynucleotide Kinase (NEB #M0201)
- T4 PNK (3´ phosphatase minus) (NEB #M0236) can be used for the 5 phosphorylation of 3´ phosphorylated mononucleotide to generate a substrate (pNp) that can be added to the 3´ end of DNA or RNA
- 5´end labeling of 3´phosphorylated oligos with T4 PNK (3´ phosphatase minus) (NEB #M0236)

Description: T4 Polynucleotide Kinase catalyzes the transfer and exchange of P, from the y position of ATP to the 5´ hydroxyl terminus of polynucleotides (doubleand single-stranded DNA and RNA), as well as nucleoside 3´ monophosphates. The T4 Polynucleotide Kinase (NEB #M0201) also catalyzes the removal of 3' phosphoryl groups from 3' phosphoryl polynucleotides, deoxynucleoside 3´ monophosphates and deoxynucleoside 3' diphosphates. The modified version (NEB #M0236) exhibits full kinase activity with no 3' phosphatase activity.

Reaction Conditions: 1X T4 Polynucleotide Kinase Reaction Buffer. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Usage Notes:

Fresh buffer is required for optimal activity (in older buffers, loss of DTT due to oxidation lowers activity).

CTP. GTP. TTP. UTP. dATP or dTTP can be substituted for ATP as a phosphate donor.

Protocols for phosphorylation (radioactive and non-radioactive) of DNA & RNA can be found at www.neb.com.

The efficiencies of blunt and recessed 5' end phosphorylation can be improved by heating to 70°C for 5 minutes, then chilling on ice prior to kinase addition and by adding PEG-8,000 to 5% (w/v).

T4 Polynucleotide Kinase requires ATP for activity, but the supplied reaction buffer does not contain ATP to allow for high specific activity radiolabeling reactions.

Often, a kinase reaction is followed by a ligation reaction. In such cases, the T4 PNK reaction is performed in ligase buffer at 37°C for 30 minutes. The product of this reaction can be used directly in the ligation reaction without a buffer change or heat inactivation UNLESS there is a need to keep other DNA fragments dephosphorylated during ligation. When this is desirable, PNK should be heat inactivated prior to ligation.

Unit Definition: One Richardson unit is defined as the amount of enzyme catalyzing the incorporation of 1 nmol of acid-insoluble [32P] in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10.000 units/ml





















NEW

5-hydroxymethyluridine DNA Kinase

#M0659S

1,000 units73 €

This is an **Enzyme for Innovation** (EFI). EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit **www.neb.com/EnzymesforInnovation** to view the full list.

Description: 5-hydroxymethyluridine DNA Kinase (5-HMUDK) transfers the gamma phosphate from ATP to the hydroxymethyl moiety of 5-hydroxymethyluridine in polymeric DNA.

Reaction Conditions: 1X T4 DNA Ligase Reaction Buffer. Incubate at 37°C. Heat inactivation: 80°C for 10 minutes.

RR NEBU 🔅 37° 🚮

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μ g of Bacillus subtilis bacteriophage SP8 genomic DNA in 30 minutes at 37°C in a total reaction volume of 20 μ l against subsequent cleavage by Ncol-HF restriction endonuclease.

Concentration: 20,000 units/ml

Phosphatase Selection Chart

	Quick Dephosphorylation Kit	Recombinant Shrimp Alkaline Phosphatase (rSAP)	Antarctic Phosphatase	Alkaline Phosphatase Calf Intestinal (CIP)
FEATURES				
100% heat inactivation	2 minutes/80°C	5 minutes/65°C	2 minutes/80°C	No
High specific activity	•	•		•
Improved stability	•	•		
Works directly in NEBuffers	•	•	•	•
Requires additive			● (Zn ²⁺)	
Quick Protocol	•			

Quick Dephosphorylation Kit

#M0508S 100 reactions68 € #M0508L 500 reactions272 €

- Dephosphorylation of DNA and RNA
- Dephosphorylation of cloning vector DNA to prevent recircularization during ligation
- Removal of dNTPs and pyrophosphate from PCR reactions prior to sequencing or SNP analysis
- Dephosphorylation of DNA prior to 5´ endlabeling using T4 Polynucleotide Kinase

Description: The Quick Dephosphorylation Kit contains Quick CIP, a heat-labile recombinant version of calf intestinal alkaline phosphatase (CIP). Quick CIP nonspecifically catalyzes the dephosphorylation of 5′ and 3′ ends of DNA and RNA phosphomonoesters. It also hydrolyses ribo- and deoxyribonucleoside triphosphates (NTPs and dNTPs). Quick CIP is useful in many applications that require the dephosphorylation of DNA or RNA ends. In cloning, dephosphorylation prevents re-ligation of linearized plasmid DNA. The enzyme can quickly dephosphorylate 5′ protruding, 5′ recessed, and blunt ends in just 10 minutes. Quick CIP may also

CutSmart RN 37° 1866

be used to degrade unincorporated dNTPs in PCR reactions to prepare templates for DNA sequencing or SNP analysis.

Quick CIP is completely and irreversibly inactivated by heating it at 80°C for 2 minutes, unlike wild type CIP, which is not inactivated by heat. This makes removal of Quick CIP prior to ligation or end-labeling unnecessary.

The Quick Dephosphorylation Kit Includes:

- CutSmart Buffer
- Quick CIP



Shrimp Alkaline Phosphatase (rSAP)

500 units60 € #M0371S #M0371L 2,500 units 240 €

- Dephosphorylation of DNA and RNA
- Dephosphorylation of cloning vector DNA to prevent recircularization during ligation
- Removal of dNTPs and pyrophosphate from PCR reactions prior to sequencing or SNP analysis
- Dephosphorylation of DNA prior to 5´ endlabeling using T4 Polynucleotide Kinase

Description: Shrimp Alkaline Phosphatase (rSAP) is a heat labile alkaline phosphatase purified from a recombinant source. rSAP is identical to the native enzyme, and contains no affinity tags or other modifications. rSAP nonspecifically catalyzes the dephosphorylation of 5' and 3' ends of DNA and RNA phosphomonoesters. Also, rSAP hydrolyzes ribo-, as well as deoxyribonucleoside triphosphates (NTPs and dNTPs), rSAP is useful in many applications that require the dephosphorylation of DNA or RNA ends. In cloning, dephosphorylation prevents religation of linearized plasmid DNA. rSAP may also be used to degrade unincorporated dNTPs in PCR reactions to prepare templates for DNA sequencing

CutSmart RN 37° 65

or SNP analysis. rSAP is completely and irreversibly inactivated by heating at 65°C for 5 minutes, thereby making removal of rSAP prior to ligation or end-labeling unnecessary.

Reaction Conditions: 1X CutSmart Reaction Buffer. Incubate at 37°C. Heat inactivation: 65°C for 5 minutes.

Unit Definition: One unit is defined as the amount of enzyme that hydrolyzes 1 µmol of p-Nitrophenyl Phosphate, PNPP (NEB #P0757) in a total reaction volume of 1 ml in 1 minute at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 1,000 units/ml

Antarctic Phosphatase

#M0289S 1,000 units68 € #M0289L 5,000 units 272 €

- Dephosphorylation of DNA and RNA
- Dephosphorylation of cloning vector DNA to prevent recircularization during ligation
- Removal of dNTPs and pyrophosphate from PCR reactions prior to sequencing or SNP analysis
- Dephosphorylation of DNA prior to 5´ endlabeling using T4 Polynucleotide Kinase

Description: Antarctic Phosphatase catalyzes the dephosphorylation of 5' and 3' ends of DNA and RNA phosphomonoesters. Antarctic Phosphatase also hydrolyzes ribo-, as well as deoxyribonucleoside triphosphates (NTPs and dNTPs). Antarctic Phosphatase is useful in many applications that require the dephosphorylation of DNA or RNA ends. In cloning, dephosphorylation prevents religation of linearized plasmid DNA. The enzyme acts on 5' protruding, 5' recessed and blunt ends. Antarctic Phosphatase may also be used to degrade unincorporated dNTPs in PCR reactions to prepare templates for DNA sequencing. The enzyme is completely and irreversibly inactivated by heating at 70°C for 5 minutes, thereby making removal of Antarctic Phosphatase prior to ligation or end-labeling unnecessary.

RX NEBU 37° VSS

Reaction Conditions: 1X Antarctic Phosphatase Reaction Buffer. Incubate at 37°C. Heat inactivation: 80°C for 2 minutes.

Unit Definition: One unit is defined as the amount of enzyme that will dephosphorylate 1 µg of pUC19 vector DNA cut with a restriction enzyme generating 5' recessed ends in 30 minutes at 37°C. Dephosphorylation is defined as > 95% inhibition of recirculation in a self-ligation reaction and is measured by transformation into E. coli. Unit assay conditions can be found at www.neb.com.

Concentration: 5.000 units/ml

Alkaline Phosphatase, Calf Intestinal (CIP)

#M0290S 1,000 units72 € #M0290L 5,000 units 288 €

- Dephosphorylation of DNA and RNA
- Dephosphorylation of cloning vector DNA to prevent recircularization during ligation
- Removal of dNTPs and pyrophosphate from PCR reactions prior to sequencing or SNP analysis
- Dephosphorylation of DNA prior to 5´ endlabeling using T4 Polynucleotide Kinase

Description: Alkaline Phosphatase, Calf Intestinal (CIP) nonspecifically catalyzes the dephosphorylation of 5' and 3' ends of DNA and RNA phosphomonoesters. CIP also hydrolyses ribo-, as well as deoxyribonucleoside triphosphates (NTPs and dNTPs). CIP is useful in many applications that require the dephosphorylation of DNA or RNA ends. In cloning, dephosphorylation prevents religation of linearized plasmid DNA. In cloning, dephosphorylation prevents religation of linearized plasmid DNA. The enzyme acts on 5' protruding, 5' recessed and blunt ends. CIP may also be used to degrade unincorporated dNTPs in PCR reactions to prepare templates for DNA sequencing or SNP analysis.

CutSmart 37° Mb

Source: Calf intestinal mucosa

Reaction Conditions: 1X CutSmart Reaction Buffer. Incubate at 37°C.

Unit Definition: One unit is defined as the amount of enzyme that hydrolyzes 1 µmol of *p*-Nitrophenyl Phosphate, PNPP (NEB #P0757) in a total reaction volume of 1 ml in 1 minute at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10.000 units/ml





















Inorganic & Thermostable Inorganic Pyrophosphatases

R{{

See page 189 for more information.

Pyrophosphatase, Inorganic (*E. coli*) #M0361S 10 units 65 € #M0361L 50 units 260 €

Pyrophosphatase Inorganic (Yeast)

#M2403S 10 units67 € #M2403L 50 units268 € Thermostable Inorganic Pyrophosphatase #M0296S 250 units72 € #M0296L 1,250 units 288 €

Apyrase

#M0398S 10 units72 € #M0398L 50 units288 €

- Highly efficient degradation of ATP to AMP
- Removal of deoxynucleotides in DNA pyrosequencing between cycles
- Conversion of 5´ triphosphorylated RNA to ligatable monophosphorylated
- Conversion of 5' triphosphorylated RNA to 5' exonuclease XRN-1 (NEB #M0338) sensitive monophosphorylated RNA
- Supplied at 10-fold higher concentration

Description: Apyrase (recombinant, *E. coli*) is a highly active ATP-diphosphohydrolase that catalyses the sequential hydrolysis of ATP to ADP and ADP to AMP releasing inorganic phosphate. It is a recombinant version of one of several isoforms of apyrase. It can also hydrolyse 5′ tri- and diphosphate ribonucleosides and deoxyribonuclesides to their respective 5′ monophosphates. Apyrase can catalyse the conversion of 5′ triphosphorylated RNA to 5′ monophosphorylated RNA by sequential removal of γ and β phosphates.

Reaction Conditions: 1X Apyrase Reaction Buffer. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme that catalyses the release of 1 µmol of inorganic phosphate from ATP (1 mM, NEB #P0756) in 1X Apyrase Reaction Buffer in 1 minute at 30°C in a total reaction volume of 50 µl. Unit assay conditions can be found at www neb com

Concentration: 500 units/ml

Notes: Apyrase has a higher ratio of activity for ATP:ADP (14:1).

Apyrase is a calcium-activated enzyme. It is approximately 50% active when Mg²⁺ substitutes Ca²⁺ in Apyrase Reaction Buffer.

As a metal-dependent enzyme, Apyrase can be inhibited by EGTA and EDTA.

The activity of Apyrase is approximately 30% in NEBuffers 1.1, 2.1, 3.1 and CutSmart Buffer.

Apyrase does not remove 5´ caps from eukaryotic mRNA

NEW

Tte UvrD Helicase

#M1202S 0.5 μg71 €

This is an **Enzyme for Innovation** (EFI). EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit **www.neb.com/EnzymesforInnovation** to view the full list.


- Unwinds double-stranded DNA
- Thermostable to 65°C
- Reduces non-specific product formation in isothermal amplification (e.g. LAMP)

Description: *Tite* UvrD Helicase is a repair helicase from the thermophilic organism *Thermoanaerobacter teng-congensis*. It is capable of unwinding double-stranded DNA without a requirement for a specific flap or overhang structure. *Tite* UvrD Helicase is active on a wide range of DNA substrates and, along with its thermostability (active to 70°C), the *Tite* UvrD Helicase has been demonstrated to be a useful additive for improving specificity of isothermal amplification reactions.



Reaction Conditions: 1X Isothermal Amplification Buffer. Supplement with 1 mM ATP. Heat inactivation: 80°C for 20 minutes.

Concentration: 20 µg/ml



Tte UvrD Helicase activity.

Properties of Exonucleases and Endonucleases

			NA Trate	ACTIVITY	OR DIGITATION .		PARTIAL			END- LABELED PHOSPHO-	UNITS FOR		
ENZYME	POLARITY	SS	ds	WITHOUT 5' PHOSPHATE	5'ext	3´ext	BLUNT	NICK	DIGESTION SS EXTENSION ²	PRODUCTS PRODUCED ³	FAM CLEAVAGE ⁴	ROTHIOATE CLEAVAGE ⁵	90% DIGESTION OF 2 M OLIGO ⁶
Lambda Exonuclease	5´ → 3´	+/-	+	+/7	+/-	+	+	-	3´	ssDNA, dNMP	+	-	2
RecJ _f	5´ → 3´	+	-	+	+/-9	-	+/-10	-	NA	dNMP	+/-	-	> 1500
Exonuclease III	3´→5´	+/-	+	+	+	+/-	+	+	5´	ssDNA, dNMP	+	_	10
Exonuclease I	3´→5´	+	-	+	-	+/-9	+/-10	NR	NA	dNMP, dinucleotide ¹¹	+	-	20
Thermolabile Exonuclease I	3´→5´	+	_	+	-	+/-9	+/-10	NR	NA	dNMP, dinucleotide ¹¹	+	_	20
Exonuclease T	3´→5´	+	-	+	-	+12	+/-10	NR	NA	dNMP	+	-	> 100
Exonuclease V	both	+	+	+	+	+	+	-	Both 5´, 3´	short oligos	+	NA	NA
Exonuclease VIII, truncated	5´→3´	-	+	+	+	+	+	_	3′	ssDNA, dNMP	NA	NA	NA
Exonuclease VII	both	+	_	+	+/-	+/-	-	-	5´	short oligos	NA	+	> 100
BAL-31 Nuclease	3´→5´, endo¹³	+	+	+	+	+	+	+	NA	dsDNA, dNMP	NA	NA	NR
Mung Bean Nuclease	endo	+	_	+	+	+	_	-	NA	ssDNA, dsDNA	NA	NA	10
DNase I	endo	+	+	NA	NA	NA	NA	NA	NA	ssDNA, dsDNA oligonucleotides, di- and trinucleotides	NA	NA	0.1
Micrococcal Nuclease	endo	+	+	NA	NA	NA	NA	NA	NA	ssDNA, dsDNA 3´ monophospho- and diphosphonucleotides ¹⁴	NA	NA	NR
Nuclease P1	endo	+	-	+	+	+	-	-	NA	5´ mononucleotides	NA	NA	NR
T5 Exonuclease	5´→3´, ss endo	+	+	+	+	+	+	+	NA	dNMP to 6-mer	+	NA	2
T7 Exonuclease	5´→3´	+/-	+	+	+/-	+	+	+	3′	ssDNA, dNMP, dinucleotide ⁸	+	-	10
T7 Endonuclease I ¹⁵	endo	_	+	NA	NA	NA	NA	+/-	NA	dsDNA	NA	NA	NA
Thermostable FEN1	endo	NA	NA	+	-	-	-	-	NA	cleaved flap	?	?	?

This table is intended to be used as a guide. Not all reported activities and properties for each exonuclease or endonuclease are listed. The amount of enzyme, substrate and time of incubation can have a dramatic effect upon the desired outcome of the experiment.

Table Legend:

- + activity; preferred substrate
- no significant activity
- +/activity greatly reduced relative to preferred substrate
- NR not reported
- not applicable NA
- single-stranded SS
- ds double-stranded
- ext extension

dNMP deoxyribonucleoside monophosphate

Footnotes:

- The ability to act on short extensions and blunt ends distinguishes these enzymes; such ends are conveniently generated by restriction digestion. The 5' and 3' extensions tested were 4 nt in length.
- Partial digestion of dsDNA by Lambda Exonuclease, T7 Exonuclease and Exonuclease III will produce dsDNA products with ss extensions. Complete digestion produces ssDNA as a product.
- Complete hydrolysis of the preferred substrate will generate the listed products.
- The ability of an exonuclease to initiate on the end of the preferred DNA substrate (ss or ds) containing a fluorescein group linked to either the 5° or 3' end. Phosphoramidite chemistry was used to synthesize oligonucleotides with FAM groups.

The 5´FAM was added to the oligonucleotide as a [(3´,6´-dipivaloylfluoresceinyl)-6-carboxamidohexyl]-I-O-(2-cyanoethyl)-(N,N-diisopropyl)-phosphoramidite while the 3' FAM oligonucleotides were synthesized using the 1-dimethoxytrityloxy-3-[O-(N-carboxy-(di-O-pivaloyl-fluorescein)-3-aminopropyl)]-propyl-2-0succinoyl-long chain alkylamino-CPG support.

- The S_a stereoisomer of the phosphorothioate linkage is known to greatly inhibit cleavage of many nucleases while the R_n stereoisomer shows either no or less inhibition. Synthetic oligonucleotides containing approximately an equal ratio of both isomers necessitates the use of multiple phosphorothipates to block cleavage The presence of 6 consecutive phosphorothicates on oligonucleotides of the preferred substrate blocked all exonucleases effectively (5% or less degradation) except T7 exonuclease which had approximately 10% cleavage. Analysis of exonucleases with oligonucleotides having either one or two consecutive phosphorothicates revealed significant degradation (not shown). T5 Exonuclease is not recommened for this application since it also has ssDNA endonuclease activity.
- The amount of enzyme in units to cleave greater than 90% of 2 µM 35-mer oligonucleotide(s) of the preferred substrate (ss or dsDNA) in a 10 µl volume using the unit digest reaction conditions. All enzymes tested, except for Exo T and RecJ., could effectively cleave the substrate to completion
- Lambda exonuclease has a strong preference for initiating degradation on dsDNA containing a 5' phosphate. Thus, if linear dsDNA has a 5' phosphate at one end and lacks a 5' phosphate on the other end, then lambda exonuclease will preferentially degrade the DNA from the phosphorylated end.

- It has been reported that the initial first product hydrolyzed from dsDNA by T7 Exonuclease is a dinucleotide. Subsequent hydrolytic cleavage releases dNMP.
- RecJ, is not suitable for making 5' extensions blunt. Exo I is not suitable for making 3' extensions blunt. Both RecJ, and Exo I require longer length ssDNA extensions to initiate than those generated by restriction enzymes.
- 10. Depending upon the DNA sequence and amount of exonuclease, RecJ., Exo I and Exo T may remove a few nucleotides from flush termini.
- 11. Exo I releases dNMP from ssDNA, except at the last hydrolytic step where a dinucleotide is produced.
- 12. Exo T can be used to make 3' extensions blunt, however, it yields 2-4 fold fewer ligatable blunt ends when compared to Klenow polymerase plus dNTP on a test substrate.
- 13. BAL31 has been reported as having both ss endonuclease activity as well as 3'→5' ds exonuclease activity. Thus, any linear DNA is a substrate for BAL31
- 14. Products of Micrococcal Nuclease degradation have 3 phosphates.
- 15. T7 endonuclease recognizes and cleaves non-perfectly matched DNA, cruciform DNA, Holliday structures or junctions. It will act more slowly on nicked dsDNA.





















Common Applications for Exonucleases and Endonucleases

APPLICATION	RECOMMENDED ENZYME(s)
Removal of 3´ overhangs	Quick Blunting Kit T4 DNA Polymerase* + dNTPs, Klenow + dNTPs
5' overhang treatment: Fill in Cleavage	Quick Blunting Kit T4 DNA Polymerase* + dNTPs, Klenow + dNTPs Mung Bean Nuclease
Removal of oligonucleotides post PCR	Exonuclease I, Thermolabile Exonuclease I, Exonuclease VII
Removal of chromosomal DNA in plasmid preparations	Lambda Exonuclease (Exonuclease I can be added to remove ssDNA generated by Lambda Exonuclease)
	T5 Exonuclease (Degrades linear ss + dsDNA, nicked DNA)
	Exonuclease V (RecBCD) (Degrades linear ss + dsDNA)
Removal of DNA in RNA preparations	DNase I
Chromatin immunoprecipitation (ChIP) analysis	Micrococcal Nuclease
Generating ssDNA from linear dsDNA: If $5' \rightarrow 3'$ polarity required If $3' \rightarrow 5'$ polarity required Best general choice	Lambda Exonuclease Exonuclease III Lambda Exonuclease

^{*}T4 DNA Polymerase has a strong 3´→ 5´ exo activity.

DNase I (RNase-Free)

1,000 units71 € 5,000 units 284 €

DNase I (RNase-free) is an endonuclease that nonspecifically cleaves DNA to release di-, tri- and

oligonucleotide products with 5´-phosphorylated and 3´-hydroxylated ends.

See page 202 for more information.

Exonuclease I (E. coli)

#M0293S 3,000 units70 € #M0293L 15,000 units 280 €

- 3´→5´ single strand exonuclease
- Degradation of post-PCR primers

RX NEBU 37° 1864

RX

Description: Catalyzes the removal of nucleotides from single-stranded DNA in the 3° to 5° direction.

Reaction Conditions: 1X Exonuclease I Reaction Buffer. Incubate at 37°C. Heat inactivation: 80°C for 20 minutes. **Unit Definition:** One unit is defined as the amount of enzyme that will catalyze the release of 10 nmol of acid-soluble nucleotide in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 20,000 units/ml

NEW

#M0303S

#M0303L

Thermolabile Exonuclease I

#M0568S 3,000 units78 € #M0568L 15,000 units312 €

- Removal of single-stranded primers in PCR reactions prior to DNA sequencing or SNP analysis
- Removal of single-stranded primers for nested PCR reactions
- Removal of single-stranded DNA from dsDNA

Description: Catalyzes the removal of nucleotides from single-stranded DNA in the 3´ to 5´ direction.

Reaction Conditions: 1X NEBuffer 3.1. Incubate at 37°C. Heat inactivation: 80°C for 1 minute.

RX NEB 3.1 37° 1866

Unit Definition: One unit is defined as the amount of enzyme that will catalyze the release of 2 nmol of acid-soluble nucleotide in a total reaction volume of 100 μ l in 6 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 20,000 units/ml

Exonuclease III (E. coli)

#M0206S 5,000 units62 € 25,000 units 248 € #M0206L

- 3´→5´ exonuclease
- Production of unidirectional nested deletions
- Site-directed mutagenesis
- Preparation of strand-specific probes
- Preparation of single-stranded substrates for dideoxy sequencing

Description: Catalyzes the stepwise removal of mononucleotides from 3´ hydroxyl temini of duplex DNA. A limited number of nucleotides are removed during each binding event, resulting in coordinated progressive deletions within the population of DNA molecules.

The preferred substrates are blunt or recessed 3' termini, although the enzyme also acts at nicks in duplex DNA to produce single-strand gaps. The enzyme is not active on ssDNA, and thus 3' protruding termini are resistant to cleavage. The degree of resistance depends on the length of the extension, with 4 bases or longer being essentially resistant to cleavage. This can be exploited to produce unidirectional deletions from a linear molecule with one resistant (3 overhang) and one susceptible (blunt or 5' overhang) terminus.

Exonuclease III has also been reported to have RNase H, 3´ phosphatase and AP-endonuclease activities.

RX NEB1 37° 166

Reaction Conditions: 1X NEBuffer 1. Incubate at 37°C. Heat inactivation: 70°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to produce 1 nmol of acid-soluble total nucleotide in a total reaction volume of 50 ul in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 100,000 units/ml

Usage Notes: Phosphorothioate linkages are not cleaved by Exonuclease III. Unidirectional deletions can also be created by protecting one terminus by incorporation of an α -phosphorothioate-containing nucleotide.

Exonuclease V (RecBCD)

#M0345S 1,000 units74 € #M0345L 5,000 units 296 €

 Degradation of linear ssDNA and dsDNA, while preserving nicked and supercoiled plasmid DNA

Description: Exonuclease V, a RecBCD complex from E. coli, has several different enzyme activities, including an ATP-dependent single-stranded DNA endonuclease activity, ss- and ds- DNA exonuclease activity. The hydrolysis in each case is bi-directional (from both the 3' and 5' ends) and processive, producing oligonucleotides. All Exonuclease V activities have divalent cation requirements. Mg2+ is required for the exonuclease activity, while Ca2+ inhibits the exonuclease activity and allows double-stranded DNA unwinding (helicase activity) without hydrolysis.

RN NEBU 37° 166

Reaction Conditions: 1X NEBuffer 4. Supplement with 1 mM ATP. Incubate at 37°C. Heat inactivation: 70°C for 30 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to produce 1 nmol of acid-soluble deoxyribonucleotide from double-stranded DNA in 30 minutes at 37°C in a total reaction volume of 50 µl. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

Exonuclease VII

200 units 159 € #M0379L 1,000 units 636 €

- Removal of ss-oligonucleotide primers
- Removal of terminal phosphorothioated ss-oligonucleotide primers
- Mapping positions of introns in genomic DNA
- Removal of single-stranded DNA from dsDNA

Description: Exonuclease VII, (Exo VII) derived from E. coli, cleaves single-stranded DNA (ssDNA) from both $5' \rightarrow 3'$ and $3' \rightarrow 5'$ direction. This enzyme is not active on linear or circular dsDNA. It is useful for removal of ss-oligonucleotide primers from a completed PCR reaction when different primers are required for subsequent PCR reactions. Digestion of ssDNA by Exonuclease VII is metal-independent.

RX NEBU 37° 166

Reaction Conditions: 1X Exonuclease VII Reaction Buffer. Incubate at 37°C. Heat inactivation: 95°C for 10 minutes

Unit Definition: One unit is defined as the amount of enzyme that will catalyze the release of 1 nmol of acidsoluble nucleotide in a total reaction volume of 50 µl in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

Exonuclease VIII, truncated

#M0545S 1,000 units74 € #M0545L 5,000 units 296 €

 Degradation of linear dsDNA while maintaining double-stranded circular DNA

Description: Exonuclease VIII, truncated, is a genetically engineered active domain of exonuclease VIII from E. coli. Exonuclease VIII, truncated is able to initiate nucleotide removal from the 5' termini of linear double-stranded DNA in the 5' to 3' direction. The enzyme does not degrade supercoiled dsDNA and circular ssDNA.

Reaction Conditions: 1X NEBuffer 4. Incubate at 37°C. Heat inactivation: 70°C for 15 minutes.

RX NEB 4 37° 166

Unit Definition: One unit is defined as the amount of enzyme required to produce 1 nmol of acid-soluble deoxyribonucleotide from double-stranded DNA in a total reaction volume of 50 µl in 30 minutes at 37°C in 1X NEBuffer 4 with 0.15 mM sonicated duplex [3H] DNA. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

















CutSmart or NEB 1.1 Recommended Buffer



Exonuclease T

#M0265S 250 units72 € #M0265L 1,250 units 288 €

- Specific for single-stranded DNA or RNA
- Generation of blunt ends from RNA or DNA with a 3´ extension

Description: Exonuclease T (Exo T), also known as RNase T, is a single-stranded RNA or DNA specific nuclease that requires a free 3′ terminus and removes nucleotides in the 3′→5′ direction. Exonuclease T can be used to generate blunt ends from RNA or DNA

molecules that have 3' extensions.

Source: Exonuclease T is overexpressed and purified as a C-terminal fusion to maltose-binding protein (MBP). Following affinity chromatography, Exo T is cleaved from MBP leaving an additional amino acid on the N-terminus and a Phe instead of a Met (i.e., Glu-Phe-Exo T instead of Met-Exo T).

RX NEBU 25° KS

Reaction Conditions: 1X NEBuffer 4. Incubate at 25°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to produce 0.1 nmol of TCA soluble nucleotides from 1 nmol of [³H]-labeled polythymidine in a total reaction volume of 100 µl in 30 minutes at 25°C in 1X NEBuffer 4 with 1 nmol [³H]-labeled polythymidine DNA.

Concentration: 5,000 units/ml

Usage Note: Exo T is has different activity on RNA vs. DNA. For RNA, 1 unit of Exo T is required to completely digest 1.0 pmol of rA20 under standard reaction conditions.

Thermostable FEN1

#M0645S 1.600 units74 €

Cleavage of flap DNA structure

This is an **Enzyme for Innovation** (EFI). EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit **www.neb.com/EnzymesforInnovation** to view the full list.

Description: Thermostable Flap Endonuclease 1, FEN1, catalyzes the cleavage of 5´ DNA flaps from branched double stranded DNA substrates, creating a 5´ phosphate terminus. FEN1 products can be ligated by DNA ligase to create double stranded DNA. *In vivo*, FEN1 is an essential component of the Okazaki fragment maturation pathway, and also plays a role in base excision repair.

Reaction Conditions: 1X ThermoPol Reaction Buffer. Incubate at 65°C.

№ RX NEBU 🔅 65° ₩6

Unit Definition: One unit is defined as the amount of enzyme required to cleave 10 pmol of 5´ flap containing oligonucleotide substrate in a total reaction volume of 10 µl in 10 minutes at 65°C. Unit assay conditions can be found at www.neb.com.

Concentration: 32.000 units/ml

Lambda Exonuclease

#M0262S 1,000 units70 € #M0262L 5,000 units280 €

- Highly processive 5´→3´ exonuclease
- Removal of 5´ mononucleotides from duplex DNA

Description: A highly processive enzyme that catalyzes the removal of mononucleotides from duplex DNA in a 5′→ 3′ direction. The preferred substrate is 5′ phosphorylated double-stranded DNA, although it will also degrade single-stranded and non-phosphorylated substrates at a reduced rate. Lambda Exonuclease is unable to initiate digestion at nicks or gaps.

Source: A genetic fusion of the *E. coli* Lambda Exonuclease gene with the gene encoding maltose binding protein (MBP). Following affinity chromatography, Lambda Exonuclease is cleaved from the fusion construct and purified away from MBP.

RN NEBU 37° KB

Reaction Conditions: 1X Lambda Exonuclease Reaction Buffer. Incubate at 37°C. Heat inactivation: 75°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to produce 10 nmol of acid-soluble deoxyribonucleotide from double-stranded substrate in a total reaction volume of 50 µl in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 5,000 units/ml

Note: 5´-OH ends are digested 20X slower than 5´-PO₄ ends. ssDNA is digested 100X slower than dsDNA.

Micrococcal Nuclease

#M0247S 320,000 gel units74 €

- Degradation of nucleic acids present in protein preparations
- In vitro translation
- Reduction of the viscosity of cell lysates during non-mechanical cell lysis preparation
- Chromatin structure analysis
- Rapid RNA sequencing
- ChIP analysis

Description: Micrococcal Nuclease is derived from *Staphylococcus aureus* and is a relatively non-specific endo-exonuclease. The enzyme digests double-stranded, single-stranded, circular and linear nucleic acids. Cleavage preferences have been observed at sites rich in adenylate, deoxyadenylate or thymidylate. Both DNA and RNA are degraded to 3° phosphomononucleotides and dinucleotides.

Source: An *E. coli* strain containing a genetic fusion of the micrococcal nuclease gene (Gene ID: 3238436) and the gene coding for maltose binding protein (MBP). The micrococcal nuclease is cleaved from the fusion protein and purified away from MBP.

Reaction Conditions: 1X Micrococcal Nuclease Reaction Buffer. Supplement with 100 μ g/ml BSA. Incubate at 37°C.

RX NEBU BSA 37° 166

Unit Definition: (Agarose Gel Unit) One gel unit is defined as the amount of enzyme required to digest 1 μ g of lambda genomic DNA in 15 minutes at 37°C into molecular DNA fragments (100–400 base pairs) on a 1.2% agarose gel.

Note: 10,000 Gel Units are approximately equal to 1,000 Kunitz Units.

Concentration: 2 x 106 gel units/ml

Notes: 1–5 mM Ca²⁺ is required for activity. The enzyme is active in the pH range 7–10, with optimal activity at pH 9.2, as long as salt concentration is less than 100 mM. Enzyme can be inactivated by addition of excess EGTA.

Mung Bean Nuclease

#M0250S 1,500 units66 € #M0250L 7,500 units 264 €

- Removal of single-stranded extensions (3´ and 5´) to leave ligatable blunt ends
- Transcriptional mapping
- Cleavage of hairpin loops
- Excision of gene coding sequences from genomic DNA

Description: A single-strand specific DNA and RNA endonuclease which will degrade single-stranded extensions from the ends of DNA and RNA molecules, leaving blunt, ligatable ends.

Source: Mung bean sprouts

Reaction Conditions: 1X Mung Bean Nuclease Reaction Buffer. Incubate at 30°C.

Unit Definition: One unit is defined as the amount of enzyme required to produce 1 µg of acid-soluble total

NEBU 30° APP

nucleotide in a total reaction volume of 50 µl in 1 minute at 37°C in 1X Mung Bean Nuclease Reaction Buffer with 0.5 mg/ml denatured Calf Thymus DNA. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

Usage Note: Do not attempt to heat inactivate, DNA will "breathe" before enzyme inactivates, causing undesirable degradation.

Nuclease BAL-31

#M0213S 50 units67 €

- Progressive shortening of duplex DNA
- Restriction site mapping

Description: Nuclease BAL-31 exonuclease degrades both 3' and 5' termini of duplex DNA without generating internal scissions. The enzyme is also a highly-specific single-stranded endonuclease which cleaves at nicks, gaps and single-stranded regions of duplex DNA and RNA.

Source: Purified from the culture medium of *Alteromonas* espejiana BAL-31. Contains a mixture of "fast" and "slow" species of the enzyme.

Reaction Conditions: 1X Nuclease BAL-31 Reaction Buffer. Incubate at 30°C. Heat inactivation: 65°C for 10 minutes.

NEBU 30° V65

Unit Definition: One unit is defined as the amount of enzyme required to remove 200 base pairs from each end of linearized double-stranded φX174 DNA (40 μg/ ml) in 50 µl of 1X Nuclease BAL-31 Reaction Buffer in 10 minutes at 30°C

Concentration: 1,000 units/ml

Usage Notes: Duplex products of the exonuclease are a mixture of blunt and staggered ends. This mixture can be cloned directly, although maximal ligation efficiency requires repairing the staggered ends with a suitable DNA polymerase.

If necessary, the enzyme may be diluted in reaction buffer just prior to use.

Nuclease P1

#M0660S 10,000 units52 €

- Removal of single-stranded tails from DNA molecules to create blunt ends
- Cleavage of hairpin loops
- DNA or RNA base compositional analysis
- Removal of nucleic acids during protein purification

Description: Nuclease P1 (from P. citrinum) is a zinc-dependent single-strand specific nuclease which hydrolyzes 3'→ 5' phosphodiester bonds in RNA and ssDNA with no base specificity. Nuclease P1 also exhibits 3´-phosphomonoesterase activity.

Although a single-strand specific nuclease, Nuclease P1 does display some activity toward double-stranded DNA (dsDNA) in Nuclease P1 Reaction Buffer. If preferentially degrading single-stranded nucleic acids (ssDNA or RNA) in the presence of dsDNA, we recommend using Nuclease P1 in 1X NEBuffer 1.1, to limit activity on dsDNA while maintaining single-strand nuclease activity.

NEBU 37° 🐇

Reaction Conditions: 1X Nuclease P1 Reaction Buffer. Incubate at 37°C. Heat inactivation: 75°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to liberate 1.0 ug of acid soluble nucleotides from Torula Yeast total RNA per minute at 37°C in 1X Nuclease P1 Reaction Buffer.

Concentration: 100,000 units/ml

Usage Notes: Substrate specificity for Nuclease P1 is as follows: 3' AMP > RNA > ssDNA >> dsDNA.

The rate of hydrolysis of 2´ AMP is 3,000-fold less than that of 3' AMP.

Rec].

#M0264S 1,000 units70 € #M0264L 5,000 units 280 €

- Single-stranded DNA specific 5´→3´ exonuclease
- Removal of deoxynucleotide monophosphates from DNA

Description: RecJ, is a single-stranded DNA-specific exonuclease that catalyzes the removal of deoxy-dNMPs from DNA in the $5' \rightarrow 3'$ direction.

DNA substrate containing a 22 base 5´ extension results in products that are a mixture of DNA fragments that have blunt-ends, 5' extensions (1–5 nucleotides) and recessed 5' ends (1-8 nucleotides). RecJ, does not require a 5' phosphate.

Source: RecJ, is overexpressed and purified as a C-terminal fusion to MBP. MBP does not affect the catalytic activity of RecJ., but does enhance RecJ. solubility.

RX NEB 2 37° KK

Reaction Conditions: 1X NEBuffer 2. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to produce 0.05 nmol TCA soluble deoxyribonucleotide in a total reaction volume of 50 µl in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 30.000 units/ml





















106

BSA Requires BSA

37° Incubation Temperature

T5 Exonuclease

#M0363S 1,000 units66 € #M0363L 5,000 units264 €

- Degradation of linear ssDNA, linear dsDNA and nicked plasmid DNA while preserving supercoiled plasmid DNA
- Used in Gibson Assembly method

RX NEBU 37° WW

Description: T5 Exonuclease degrades DNA in the 5´ to 3´ direction. T5 Exonuclease is able to initiate nucleotide removal from the 5´ termini or at gaps and nicks of linear or circular dsDNA. However, the enzyme does not degrade supercoiled dsDNA. T5 Exonuclease also has ssDNA assay condition endonuclease activity.

Reaction Conditions: 1X NEBuffer 4. Incubate at 37°C.

Unit Definition: One unit is defined as the amount of enzyme required to produce 1 nmol of acid-soluble deoxyribonucleotide from double-stranded DNA in a total reaction volume of 50 μ l in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

T7 Exonuclease

#M0263S 1,000 units64 € #M0263L 5,000 units256 €

- 5´→3´ exonuclease
- Removal of 5´ mononucleotides from DNA

Description: T7 Exonuclease acts in the 5′ to 3′ direction, catalyzing the removal of 5′ mononucleotides from duplex DNA. T7 Exonuclease is able to initiate nucleotide removal from the 5′ termini or at gaps and nicks of double-stranded DNA. It will degrade both 5′ phosphorylated or 5′ dephosphorylated DNA. It has also been reported to degrade RNA and DNA from RNA/ DNA hybrids in the 5′ to 3′ direction, but it is unable to degrade either double-stranded or single-stranded RNA.

RX NEBU 25° YM

Reaction Conditions: 1X NEBuffer 4. Incubate at 25°C.

Unit Definition: One unit is defined as the amount of enzyme required to produce 1 nmol of acid-soluble deoxyribonucleotide from double-stranded DNA in a total reaction volume of $50\,\mu l$ in 30 minutes at $25^{\circ}C$. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

37°

Nucleoside Digestion Mix

#M0649S 50 reactions 119 €

The Nucleoside Digestion Mix is a mixture of enzymes that provides a convenient one-step method to generate single nucleosides from DNA or RNA. Optimized for

quantitative analysis by liquid chromatography-mass spectrometry (LC-MS), this reagent eliminates the need for sequential multi-step, time-consuming digestion protocols

See page 202 for more information.

Properties of DNA Repair Enzymes

Substrates and Cleavage Products of DNA Repair Glycosylases

				CREATED Leavage
ENZYME	SUBSTRATE	CLEAVAGE SITE	5´-terminus	3´-terminus
APE 1	AP sites	1st phosphodiester bond 5´ to the lesion	dR5P	OH
Endo III	Oxidized Pyrimidines, AP sites	1st phosphodiester bond 3´ to the lesion	Р	PA
Endo IV	AP sites	1st phosphodiester bond 5´ to the lesion	dR5P	OH
Endo V	Deoxyinosine, mismatches, hairpin/unpaired loop, flaps, pseudo Y structures	2nd phosphodiester bond 3´ to the lesion	Р	ОН
Endo VIII	Oxidized Pyrimidines, AP sites	1st phosphodiester bond both 5' and 3' to the lesion	Р	Р
Fpg	Oxidized Purines, AP sites	1st phosphodiester bond both 5' and 3' to the lesion	Р	Р
hOGG1	Oxidized Purines, AP sites	1st phosphodiester bond 3´ to the lesion	Р	PA
T7 Endo I	Cruciform, Holliday junctions, mismatches, herteroduplexes	1st, 2nd or 3rd phosphodiester bond 5´ to the mismatches	Р	OH
T4 PDG	Pyrimidine Dimers	N-glycosidic bond of the 5′ Thymine of the dimer and the 1st phosphodiester bond 3′ to the AP site	P*	PA
UDG				
Afu UDG			45 %	
Antarctic Thermolabile UDG	Deoxyuridine	N-glycosidic bond	AP site (i	no break)
hSMUG1	Deoxyuridine, 5-hyroxyuracil, 5-hydroxymethyluracil 5-formyluracil	N-glycosidic bond	AP site (no break)
hAAG	Deoxyinosine, Alkylated purines	N-glycosidic bond	AP site (no break)

*A pyrimidine dimer still covalently attached.

Table	Le	ge	nd:
AF)	an	uri

AP apurinic/apyrimidinic sites

P phosphate

OH hydroxyl

dR5P deoxyribose-5´-phosphate

PA 3´-phospho-α, β-unsaturated aldehyde

DNA Repair Glycosylases on Various Damaged Bases

		DOUBL	E-STRANDE	D DNA OLIGO	OS (34-MERS	3)						
	ENZYME	AP:A	DHT:A	5-hmU:A	5-hmU:G	l:T	6-MeA:T	8-0G:C	8-0G:G	U:A	U:G	THYMINE GLYCOL:A
	APE 1	++++	+	-	-	-	-	-	-	-	-	-
	Endo III	++++	+	-	-	-	-	-	-	-	-	+
	Tma Endo III	++++	++	-	-	-	-	+	+	-	-	++
	Endo IV	++++	+	-	-	-	-	-	-	-	-	-
	Tth Endo IV	++++	+	-	-	-	-	-	+	-	-	-
	Endo V*	+++	+	+	+	++++	+	++	+	+	++++	++
	Endo VIII	++++	++	-	-	-	-	-	-	-	-	+++
	Fpg	+	+	-	-	-	-	++++	++++	-	-	+
	hAAG	_	_	-	-	++++	-	-	-	-	-	-
	hNEIL1	++++	++	-	-	-	-	+	+	-	-	++
60	hOGG1	++	_	-	-	-	-	++++	+	-	-	-
Ü	T4 PDG	++++	-	-	-	-	-	-	-	-	-	-
REPAIR ENZYMES	UDG	N/A	_	-	-	-	-	-	-	++++	+	-
AIR	Afu UDG	N/A	-	-	-	-	-	-	-	++++	+	-
RE	hSMUG1	N/A	_	+++	+++	_	-	-	-	++++	++++	-

Standard reaction conditions were used to titer the enzymes with the alternate base.

		SINGLE	SINGLE-STRANDED DNA OLIGOS (34-MERS)						
	ENZYME	AP	DHT	5-hmU	ı	6-MeA	8-0G	U	THYMINE GLYCOL:A
	APE 1	++	-	-	-	-	-	-	-
	Endo III	++	-	-	-	-	-	-	-
	Tma Endo III	++	+	-	_	-	-	-	-
	Endo IV	-	-	-	-	-	-	-	-
	Tth Endo IV	-	-	-	-	-	-	-	-
	Endo V	+	-	-	++++	-	+	-	-
	Endo VIII	+++	=	-	-	-	-	=	-
	Fpg	+	+	-	-	-	+	-	+
	hAAG	-	_	-	+	_	_	_	-
	hNEIL1	+	+	-	-	-	-	-	+
s	hOGG1	++	_	-	_	_	+	_	-
YME	T4 PDG	-	-	-	-	-	-	-	-
ENZ	UDG	N/A	-	-	_	-	_	++++	_
REPAIR ENZYMES	Afu UDG	N/A	-	-	-	-	_	++++	-
RE	hSMUG1	N/A	-	++	_	_	-	+++	_

Table Legend:

AP apurinic/apyrimidinic site. The AP site is created by treating a uracil containing oligo with UDG.

DHT 5,6-dihydrothymine

5-hmU 5-hydroxymethyluracil

I inosine

6-MeA 6 methyladenine

8-0G 8 oxoguanine

II uridino

AP:A apurinic/apyrimidinic site base paired with adenine

DHT:A 5,6 dihydrothymine base paired with an adenine

5-hmU:A 5-hydroxymethyluracil base paired with an adenine

5-hmU:G 5-hydroxymethyluracil base paired with a guanine

I:T inosine base paired with a thymine

6-MeA:T 6-methyladenine base paired with a thymine

8-OG:C 8-oxoguanine base paired with a cytosine

8-OG:G 8-oxoguanine base paired with a guanine

U:A uridine paired with an adenine

U:G uridine paired with a guanine

Activity Level:

++++ 100%

+++ 50%

++ 10% - 25%

+ < 10%

no detectable enzyme activity (< 0.7%)

N/A not applicable

















^{*}Nicks only, does not remove damage

APE 1

#M0282S 1,000 units72 € #M0282L 5,000 units288 €

- Single-cell gel electrophoresis (Comet assay)
- Alkaline elution
- Alkaline unwinding
- Modified nick translation

Description: Human apurinic/apyrimidinic (AP) endonuclease, APE 1, also known as HAP 1 or Ref-1, shares homology with *E. coli* Exonuclease III. APE 1 cleaves the phosphodiester backbone immediately 5' to an AP site via hydrolytic mechanism, generating a ssDNA break that leaves a 3'-hydroxyl and 5'-deoxyribose phosphate terminus. Besides AP endonuclease activity, APE 1 has also been reported to have weak DNA 3'-diesterase, 3' to 5' exonuclease and RNase H activities.

Reaction Conditions: 1X NEBuffer 4. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

RX NEB 4 37° 664

Unit Definition: One unit is defined as the amount of enzyme required to cleave 20 pmol of a 34-mer oligonucleotide duplex containing a single AP site in a total reaction volume of 10 µl in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Recommended Dilution for the Comet Assay: 1:10³. A detailed protocol can be found at www.neb.com.

Concentration: 10,000 units/ml

Endonuclease IV

#M0304S 1,000 units76 € #M0304L 5,000 units304 €

- Single-cell gel electrophoresis (Comet assay)
- Alkaline elution
- Alkaline unwinding

Description: Endonuclease IV can act on several types of oxidative damage in DNA. The enzyme is an apurinic/apyrimidinic (AP) endonuclease that will hydrolyze intact AP sites in DNA. AP sites are cleaved at the first phosphodiester bond that is 5´ to the lesion leaving a hydroxyl group at the 3´ terminus and a deoxyribose 5´-phosphate at the 5´ terminus. The enzyme also has a 3´-diesterase activity and can release phosphoglycolal-dehyde, intact deoxyribose 5´-phosphate and phosphate from the 3´ end of DNA.

Reaction Conditions: 1X NEBuffer 3. Incubate at 37°C. Heat inactivation: 85°C for 20 minutes.

RN NEB3 37° 166

Unit Definition: One unit is defined as the amount of enzyme required to cleave 1 pmol of a 34-mer oligonucleotide duplex containing a single AP site in a total reaction volume of 10 µl in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Recommended Dilution for the Comet Assay: $1:10^4$ to $1:10^5$. A detailed protocol can be found at www.neb.com.

Concentration: 10,000 units/ml

Tth Endonuclease IV

#M0294S 500 units76 €

- Thermostable
- Alkaline elution
- Alkaline unwinding

Description: *Tth* Endonuclease IV is a thermostable apurinic/apyrimidinic (AP) endonuclease. *Tth* Endo IV will hydrolyze an AP site at the first phosphodiester bond 5′ to the lesion leaving a 3′ hydroxyl and a deoxyribose 5′-phosphate. The enzyme also has a 3′-diesterase activity.

Reaction Conditions: 1X ThermoPol Reaction Buffer. Incubate at 65°C.

RX NEBU 65° WW

Unit Definition: One unit is defined as the amount of enzyme required to cleave 1 pmol of a 60-mer oligonucleotide duplex containing a single AP site in a total reaction volume of $10 \, \mu l$ in 1 hour at 65°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

Endonuclease III (Nth)

#M0268S 1,000 units72 €

- Single-cell gel electrophoresis (Comet assay)
- Alkaline elution
- Alkaline unwinding

Description: Endonuclease III (Nth) protein from *E. coli* acts both as a *N*-glycosylase and an AP-lyase. The *N*-glycosylase activity releases damaged pyrimidines from dsDNA generating a basic (AP) site. The AP-lyase activity of the enzyme cleaves 3´ to the AP site leaving a 5´ phosphate and a 3´-phospho $-\alpha$, β unsaturated aldehyde.

Some of the damaged bases recognized and removed by Endonuclease III (Nth) include urea, 5, 6 dihydroxythymine, thymine glycol, 5-hydroxy-5 methylhydanton, uracil glycol, 6-hydroxy-5, 6-dihdrothimine and methyltartronylurea.

Source: An *E. coli* strain which carries the cloned *nth* gene.

RX NEBU 37° KK

Reaction Conditions: 1X Endonuclease III (Nth) Reaction Buffer. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to cleave 1 pmol of a 34-mer oligonucleotide duplex containing a single AP site in a total reaction volume of $10~\mu l$ in 1 hour at $37^{\circ}C$. Unit assay conditions can be found at www.neb.com.

Recommended Dilution for the Comet Assay: 1:10⁴ to 1:10⁵. A detailed protocol can be found at www.neb.com.

Concentration: 10,000 units/ml

Tma Endonuclease III

500 units76 € #M0291S

- Alkaline elution
- Alkaline unwinding

Description: A thermostable homolog of the *E. coli* Endonuclease III (Nth). It acts as an N-glycosylase and an AP-lyase. The N-glycosylase activity releases damaged pyrimidines from double-stranded DNA, generating an apurinic (AP) site. The AP-lysase activity then cleaves the resulting abasic site.

Tma Endonuclease III recognizes abasic sites, 5,6 dihydroxythymine and thymine glycol in DNA.

RX NEBU 65 VIS

Reaction Conditions: 1X ThermoPol Reaction Buffer. Incubate at 65°C.

Unit Definition: One unit is defined as the amount of enzyme required to cleave 1 pmol of a 60-mer oligonucleotide duplex containing a single AP site in a total reaction volume of 10 µl in 1 hour at 65°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10.000 units/ml

Endonuclease V

#M0305S 250 units76 €

- Cleavage of oligonucleotides containing deoxyinosines
- Mismatch cleavage

Description: Endonuclease V is a repair enzyme found in E. coli that recognizes deoxyinosine, a deamination product of deoxyadenosine in DNA. Endonuclease V. often called Deoxyinosine 3´ Endonuclease, recognizes DNA containing deoxyinosines (paired or not) on double-stranded DNA, single-stranded DNA with deoxyinosines and, to a lesser degree, DNA containing abasic sites (AP) or urea, base mismatches, insertion/ deletion mismatches, hairpin or unpaired loops, flaps and pseudo-Y structures.

Endonuclease V cleaves the second phosphodiester bond 3´ to the mismatch of deoxyinosine, leaving a nick with 3´-hydroxyl and 5´-phosphate.

RX NEB 4 37° 1654

Source: An E. coli strain containing a gene fusion of the Endo V gene and the gene coding for the maltose binding protein (MBP). The fusion protein is purified to near homogeneity and is active as a fusion. The protein contains 223 amino acids and has a molecular weight of 24 9 kDa

Reaction Conditions: 1X NEBuffer 4. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to cleave 1 pmol of a 34-mer oligonucleotide duplex containing a single deoxyinosine site in a total reaction volume of 10 µl in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

Endonuclease VIII

#M0299S 1,000 units76 € #M0299L 5.000 units 304 €

- Single-cell gel electrophoresis (Comet assay)
- Alkaline elution
- Alkaline unwinding

Description: Endonuclease VIII acts as both an Nglycosylase and an AP-lyase. The N-glycosylase activity releases damaged pyrimidines from double-stranded DNA, generating an apurinic (AP) site. The AP-Ivase activity cleaves 3´ and 5´ to the AP site leaving a 5´ phosphate and a 3´ phosphate. Damaged bases recognized and removed by Endonuclease VIII include urea, 5, 6- dihydroxythymine, thymine glycol, 5-hydroxy-5- methylhydanton, uracil glycol, 6-hydroxy-5, 6-dihydrothymine and methyltartronylurea. While Endonuclease VIII is similar to Endonuclease III, Endonuclease VIII has β and δ lyase activity while Endonuclease III has β lyase activity.

RX NEBU 37° 16

Reaction Conditions: 1X Endonuclease VIII Reaction Buffer. Incubate at 37°C. Heat inactivation: 75°C for 10 minutes

Unit Definition: One unit is defined as the amount of enzyme required to cleave 1 pmol of a 34 mer oligonucleotide duplex containing a single AP site in a total reaction volume of 10 µl in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Recommended Dilution for the Comet Assay: 1:104 to 1:105. A detailed protocol can be found at

Concentration: 10.000 units/ml

Fpg

#M0240S 500 units74 € #M0240L 2,500 units 296 €

- Single-cell gel electrophoresis (Comet assay)
- Alkaline elution
- Alkaline unwinding

Description: Fpg (formamidopyrimidine [fapy]-DNA glycosylase), also known as 8-oxoguanine DNA glycosylase, acts both as an N-glycosylase and an AP-lyase. The N-glycosylase activity releases damaged purines from double-stranded DNA, generating an apurinic (AP) site. The AP-lyase activity cleaves both 3' and 5' to the AP site thereby removing the AP site and leaving a 1-base gap with a 5' and 3' phosphate.

Some of the damaged bases recognized and removed by Fpg include 7, 8-dihydro-8-oxoguanine (8-oxoguanine), 8-oxoadenine, fapy-guanine, methy-fapy-guanine, fapyadenine, aflatoxin B1-fapy-guanine, 5-hydroxy-cytosine and 5-hydroxy-uracil.

RN NEB1 BSA 37° WW

www neb com

Reaction Conditions: 1X NEBuffer 1. Supplement with 100 μg BSA. Incubate at 37°C. Heat inactivation: 60°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to cleave 10 pmol of a 34-mer oligonucleotide duplex containing a single 8-oxoguanine base paired with a cytosine in a total reaction volume of 10 µl in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Recommended Dilution for the Comet Assay: 1:103 to 1:104. A detailed protocol can be found at www neb com

Concentration: 8,000 units/ml





















110

Recombinant Enzyme

hAAG

#M0313S 500 units76 €

- Single-cell gel electrophoresis (Comet Assay)
- Alkaline elution
- Alkaline unwinding

Description: Human Alkyladenine DNA Glycosylase (hAAG) excises alkylated and oxidative DNA damaged sites, including 3-methyladenine, 7-methylguanine, 1,N6-ethenoadenine and hypoxanthine. hAAG catalyzes the hydrolysis of the *N*-glycosidic bond to release the damaged base. hAAG is also known as methylpurine DNA glycosylase (MPG) or 3-methyladenine-DNA glycosylase (ANPG).

Reaction Conditions: 1X ThermoPol Reaction Buffer. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

RX NEBU 37° KM

Unit Definition: One unit is defined as the amount of enzyme required to create an AP site from 1 pmol of a 34-mer oligonucleotide duplex containing a single deoxyinosine site in a total reaction volume of 10 µl in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

T4 PDG (T4 Endonuclease V)

#M0308S

2,000 units76 €

- DNA damage studies
- Single-cell gel electrophoresis (Comet assay)

Description: T4 PDG (pyrimidine dimer glycosylase) has both DNA glycosylase and AP-lyase activity. The 16 kd protein recognizes cis-syn-cyclobutane pyrimidine dimers caused by UV irradiation. The enzyme cleaves the glycosyl bond of the 5' end of the pyrimidine dimer and the endonucleolytic activity cleaves the phosphodiester bond at the apurinic/apyrimidinic (AP) site.

Reaction Conditions: 1X T4 PDG Reaction Buffer. Supplement with BSA. Incubate at 37°C.

Unit Definition: One unit is defined as the amount of enzyme that catalyzes the conversion of $0.5 \mu g$

RX NEBU BSA 37° Mb

of UV irradiated supercoiled pUC19 DNA to > 95% nicked plasmid in a total reaction volume of 20 μl in 30 minutes at 37°C. Nicking is assessed by agarose gel electrophoresis. Irradiated plasmid contains an average of 3–5 pyrimidine dimers. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

Note: Incubation time should be ≤ 30 minutes for best results.

Uracil-DNA Glycosylase (UDG)

#M0280S #M0280L 1,000 units74 € 5,000 units 296 €

Companion Product:

Uracil Glycosylase Inhibitor (UGI)

#M0281S 200 units72 € #M0281L 1,000 units288 €

- Eliminates PCR carry-over contamination
- Release of uracil from ss- or ds- DNA

Description: *E. coli* Uracil-DNA Glycosylase (UDG) catalyzes the release of uracil from uracil-containing DNA. UDG efficiently hydrolyzes uracil from single-stranded or double-stranded DNA, but not from oligomers (6 or fewer bases).

Reaction Conditions: 1X UDG Reaction Buffer Incubate at 37°C.

Unit Definition: One unit is defined as the amount of enzyme that catalyzes the release of 60 pmol of uracil per minute from double-stranded, uracil-containing DNA.

RX NEBU 37° KM

Activity is measured by release of [3 H]-uracil in a 50 µl reaction containing 0.2 µg DNA (10^{4} - 10^{5} cpm/µg) in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 5,000 units/ml

Usage Notes: UDG is active over a broad pH range with an optimum at pH 8.0, does not require divalent cation, and is inhibited by high ionic strength (> 200 mM).

Afu Uracil-DNA Glycosylase (UDG)

#M0279S #M0279L 200 units76 € 1,000 units 304 €

- Thermostable
- Release of uracil from ss- or ds- DNA

Description: A thermostable homolog of the *E. coli* Uracil-DNA Glycosylase (UDG) from *Archaeoglobus fulgidus*. *Afu* UDG catalyses the release of uracil from uracil-containing DNA. *Afu* UDG efficiently hydrolyzes uracil from single-stranded or double-stranded DNA.

Reaction Conditions: 1X ThermoPol II (Mg-free) Reaction Buffer. Incubate at 65°C.

Unit Definition: One unit is defined as the amount of enzyme that catalyzes the release of 60 pmol of uracil per minute from double-stranded uracil-containing

RX NEBU 65° W6

DNA. Activity is measured by release of [3 H]-uracil in a 50 μ l reaction containing 0.2 μ g DNA (4 - 10 cpm/ 4) in 30 minutes at 65 $^{\circ}$ C. Unit assay conditions can be found at www.neb.com.

Concentration: 2,000 units/ml.

Usage Notes: Afu UDG retains 50% activity in the presence of 150 mM NaCl. Afu UDG retains less than 1% activity after boiling for 30 minutes in standard reaction conditions. Under standard reaction conditions, uracil glycosylase inhibitor (UGI) does not inhibit Afu UDG.

Antarctic Thermolabile UDG

#M0372S 100 units76 € #M0372L 500 units 304 €

- Eliminates PCR carry-over contamination
- Release of uracil from ss- or ds- DNA

Description: Antarctic Thermolabile UDG (Uracil-DNA Glycosylase) catalyzes the release of free uracil from uracil-containing single-stranded or double-stranded DNA. The resulting abasic sites are susceptible to the hydrolytic cleavage at the elevated temperature and high pH. This enzyme is sensitive to heat and can be rapidly and completely inactivated at temperatures above 50°C.

Reaction Conditions: 1X Standard Tag Reaction Buffer. Incubate at 37°C. Heat inactivation: 50°C for 10 minutes

RX NEBU 37° this

Unit Definition: One unit is defined as the amount of enzyme that catalyzes the release of 60 pmol of uracil per minute from double-stranded, uracil-containing DNA. Activity is measured by release of [3H]-uracil in а 50 µl Standard *Tag* reaction buffer containing 0.2 µg DNA (104-105 cpm/ug) in 30 minutes at 37°C.

Concentration: 1,000 units/ml

hSMUG1

#M0336S 500 units76 €

- Oxidative DNA damage studies
- Single-cell gel electrophoresis (Comet assay)

Description: Human single-strand-selective monofunctional uracil-DNA Glycosylase (hSMUG1) excises deoxyuracil and deoxyuracil-derivatives bearing an oxidized group at C5, such as 5-hydroxyuracil, 5-hydroxymethyluracil and 5-formyluracil in ssDNA and dsDNA

Reaction Conditions: 1X NEBuffer 1. Supplement with 100 μg BSA. Incubate at 37°C.Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to excise 1 pmol of deoxyuracil from a 34 mer oligonucleotide duplex containing a single dU site in a total reaction volume of 10 µl in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 5,000 units/ml.

RX NEB1 BSA 37° 1654

Usage Notes: hSMUG1 has 50% activity on 5-hydroxymethyluracil when compared to uracil. hSMUG1 has 50% activity on ssDNA compared to dsDNA.

PreCR® Repair Mix

#M0309S 30 reactions 157 € #M0309L 150 reactions 628 €

The PreCR Repair Mix is an enzyme cocktail formulated to repair damaged template DNA prior to its use in the

polymerase chain reaction (PCR), microarrays or other

See page 76 for more information.

β-Agarase I

#M0392S 100 units79 € #M0392L 500 units 316 €

- Extraction of DNA from agarose gels
- DNase and RNase free

Description: β-Agarase cleaves the agarose subunit, unsubstituted neoagarobiose [3,6-anhydro- α -Lgalactopyranosyl-1-3-D-galactose] to neoagarooligosaccharides. β-Agarase I digests agarose, releasing trapped DNA and producing carbohydrate molecules which can no longer gel. The remaining carbohydrate molecules and β-Agarase I will not, in general, interfere with subsequent DNA manipulations such as restriction endonuclease digestion, ligation and transformation.

Reaction Conditions: 1X β-Agarase I Reaction Buffer. Incubate at 42°C.

RX NEBU 42° 1654

RR W

DNA technologies.

Unit Definition: One unit is defined as the amount of enzyme required to digest 200 µl of molten low melting point or NuSieve agarose to nonprecipitable neoagarooligosaccharides in 1 hour at 42°C.

Concentration: 1,000 units/ml

Heat Inactivation: Incubation at 95°C for 2 minutes or incubation at 65°C for 15 minutes inactivates 50 units of $\beta\textsc{-Agarase}$ I. $\beta\textsc{-Agarase}$ I retains activity for several hours at 40-45°C and is stabilized by the presence of agarose in the reaction.



112



















Cre Recombinase

#M0298S 50 units70 € #M0298L 250 units280 €

for high (15X) concentration

#M0298M 250 units 280 €

- Excision of DNA between two loxP sites
- Fusion of DNA molecules containing loxP sites
- Inversion of DNA between loxP sites

Description: Cre Recombinase is a Type I topoisomerase from bacteriophage P1 that catalyzes the site-specific recombination of DNA between *loxP* sites. The enzyme requires no energy cofactors, and Cre-mediated recombination quickly reaches equilibrium between substrate and reaction products. The *loxP* recognition element is a 34 base pair (bp) sequence comprised of two 13 bp inverted repeats flanking an 8 bp spacer region which confers directionality. Recombination products depend on the location and relative orientation of the *loxP* sites. Two DNA species containing single *loxP* sites will be fused. DNA between repeated *loxP* sites will be excised in circular form while DNA between opposing *loxP* sites will be inverted with respect to external sequences.

RX NEBU 37° 📆

Reaction Conditions: 1X Cre Recombinase Reaction Buffer. Incubate at 37°C. Heat inactivation: 70°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme necessary to produce maximal site-specific recombination of 0.25 μ g pLox2+ control DNA in a total reaction volume of 50 μ l in 30 minutes at 37°C. Maximal recombination is determined by agarose gel analysis and by transformation of reactions followed by selection on ampicillin plates.

Concentration: 1,000 units/ml and 15,000 units/ml.

T7 Endonuclease I

#M0302S 250 units69 € #M0302L 1,250 units 276 €

- Recognition of mismatched DNA
- Resolve four-way junction or branched DNA
- Detection or cleavage of heteroduplex and nicked DNA
- Random cleavage of linear DNA for shotgun cloning

Description: T7 Endonuclease I recognizes and cleaves non-perfectly matched DNA, cruciform DNA structures, Holliday structures or junctions, heteroduplex DNA and more slowly, nicked ds- DNA. It cleaves at the first, second or third phosphodiester bond 5´ to the mismatch.

Source: An *E. coli* strain that carries a fusion of maltose binding protein (MBP) and T7 Endo I.

Reaction Conditions: 1X NEBuffer 2. Incubate at 37°C.

RX NEB 2 37° W/6

Unit Definition: One unit is defined as the amount of enzyme required to convert > 90% of 1 μ g of supercoiled cruciform pUC(AT)* to > 90% linear form in a total reaction volume of 50 μ l in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

*pUC(AT) is derived from pUC19 with a modification of the polylinker between the EcoRI and Pstl sites.

Concentration: 10,000 units/ml

Usage Note: It is important to control the amount of enzyme and the reaction time used for cleavage of a particular substrate. Temps. > 42°C increase nonspecific nuclease activity.

NEW

TelN Protelomerase

#M0651S 250 units74 €

- Isolated from phage N15
- Cuts dsDNA at a TeIN recognition sequence (56 bp)
- Leaves covalently closed ends at the site of cleavage

This is an **Enzyme for Innovation** (EFI). EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit **www.neb.com/EnzymesforInnovation** to view the full list.

Description: TelN Protelomerase, from phage N15, cuts dsDNA at a TelN recognition sequence (56 bp) and leaves covalently closed ends at the site of cleavage.

Reaction Conditions: 1X ThermoPol Reaction Buffer. Incubate at 30°C. Heat inactivation: 75°C for 5 minutes.

💥 RX MBU 🛊 30° 🐇

Unit Definition: One unit is defined as the amount of enzyme required to cleave 0.5 µg of pMiniT-TelN Bsallinearized control plasmid (313 fmol of TelN recognition sites) in a total reaction volume of 50 µl in 30 minutes at 30°C in 1X ThermoPol Reaction Buffer.

Concentration: 5,000 units/ml

Topoisomerase I (E. coli)

#M0301S 100 units72 € #M0301L 500 units288 €

- Recognition of mismatched DNA
- Catalyzes relaxation of negativelysupercoiled DNA

Description: Topoisomerase I (*E. coli*) catalyzes the relaxation of negatively-supercoiled DNA. Topoisomerase I has also been implicated in knotting and unknotting DNA, and in linking complementary rings of single-stranded DNA into double-stranded rings. The intact holoenzyme is a 97 kDa protein.

Reaction Conditions: 1X CutSmart Buffer. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

CutSmart RX 37°

Unit Definition: One unit is defined as the amount of enzyme that catalyzes the relaxation of > 95% of 0.5 μ g of negatively supercoiled pUC19 RF I DNA in a total reaction volume of 25 μ l in 15 minutes at 37°C. DNA supercoiling is assessed by agarose gel electrophoresis in the absence of ethidium bromide.

Concentration: 5,000 units/ml.

CpG Methyltransferase (M.SssI)

#M0226S	100 units	71€
#M0226L	500 units	284 €

for high (5X) concentration

#M0226M 500 units 284 €

- Blocking restriction enzyme cleavage
- Studying of CpG methylation-dependent gene expression
- Probing sequence-specific contacts within the major groove of DNA
- Altering the physical properties of DNA
- Uniform [3H]-labeling of DNA
- Decreasing the number of RE cut sites, yielding an apparent increase in specificity

Description: The CpG Methyltransferase (M.Sssl)

methylates all cytosine residues (C5) within the doublestranded dinucleotide recognition sequence 5´...CG...3´.

Reaction Conditions: 1X NEBuffer 2 + SAM. Supplement with 160 µM S-adenosylmethionine (supplied). Incubate at 37°C. Heat inactivation: 65° for 20 minutes.

RN NEB 2 SAM 37° 1656

Note: MaCl, is not required as a cofactor. In the presence of Mg²⁺, methylation becomes distributive rather than processive and exhibits topoisomerase activity.

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg of λ DNA in a total reaction volume of 20 µl in 1 hour at 37°C against cleavage by BstUI restriction endonuclease.

Concentration: 4,000 units/ml and 20,000 units/ml. Assayed on λ DNA.

See pages 334-336 for a complete list of restriction endonucleases blocked by CpG methylation.

GpC Methyltransferase (M.CviPI)

#M0227S 200 units74 € #M0227L 1,000 units 296 €

CH₃ 5'... G C ... 3' 3'... C G ... 5 CH,

- Blocking restriction enzyme cleavage
- Altering the physical properties of DNA
- Uniform [³H]-labeling of DNA

Description: The GpC Methyltransferase (M.CviPI) methylates all cytosine residues (C5) within the doublestranded dinucleotide recognition sequence 5'...GC...3'.

Reaction Conditions: 1X GC Reaction Buffer. Supplement with 160 µM SAM (supplied). Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Note: MgCl₂ is not required as a cofactor.

RN NEBU SAIM 37° 1654

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg of λ DNA in a total reaction volume of 20 μl in 1 hour at 37°C against cleavage by HaellI restriction endonuclease.

Concentration: 4,000 units/ml. Assayed on λ DNA.

For more information on products to study DNA methylation, see pages 260-273.

AluI Methyltransferase

#M0220S

100 units73 €

CH, 5'... A G C T ... 3' 3'... T C G A ... 5' CH₃

Description: Alul Methyltransferase modifies the cytosine residue (C5) in the sequence to the left.

Reaction Conditions: 1X Alul Methyltransferase Reaction Buffer + SAM. Supplement with 80 µM S-adenosylmethionine (supplied). Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

RR NEBU SAM 37° 1654

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μ g λ DNA in 1 hour at 37°C in a total reaction volume of 10 μl against cleavage by Alul restriction endonuclease.

Concentration: 5,000 units/ml

BamHI Methyltransferase

#M0223S #M0223L

100 units74 € 500 units 296 €

CH₃ 5'... GGATCC...3' 3'... CCTAGG...5' CH3

Description: BamHI Methyltransferase modifies the internal cytosine residue (N4) in the sequence to the left.

Reaction Conditions: 1X BamHI Methyltransferase Reaction Buffer + SAM. Supplement with 80 μM S-adenosylmethionine (supplied). Incubate at 37°C.

RX NEBU SAM 37° WW

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μ g λ DNA in 1 hour at 37°C in a total reaction volume of 10 µl against cleavage by BamHI restriction endonuclease.

Concentration: 4,000 units/ml















dam Methyltransferase

#M0222S 500 units75 € #M0222L 2,500 units300 €

5'... G Å T C ... 3' 3'... C T A G ... 5' **Description:** *dam* Methyltransferase modifies the adenine residue (N⁶) in the sequence to the left.

Reaction Conditions: 1X *dam* Methyltransferase Reaction Buffer + SAM. Supplement with $80 \mu M$ S-adenosylmethionine (supplied). Incubate at 37° C. Heat inactivation: 65° C for 20 minutes.

RN NEBU SAM 37° 155

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μ g λ DNA in 1 hour at 37°C in a total reaction volume of 10 μ l against cleavage by Mbol restriction endonuclease.

Concentration: 8,000 units/ml

EcoGII Methyltransferase

#M0603S 200 units73 €

 $\begin{array}{c} \text{CH}_{_{3}} \\ \text{5'} \dots \text{ A} \dots \text{ 3'} \\ \text{3'} \dots \text{ T} \dots \text{ 5'} \end{array}$

This is an Enzyme for Innovation (EFI).

EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit www.neb.com/EnzymesforInnovation to view the full list.

Description: EcoGII Methyltransferase is a non-specific methyltransferase that modifies adenine residues (N^6) in any sequence context.

Reaction Conditions: 1X CutSmart Buffer + SAM. Supplement with 160 μ M S-adenosylmethionine (supplied). Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to protect 100 ng FAM-labeled dsDNA in 30 minutes at 37°C in a total reaction volume of 20 μ l against cleavage by Mbol restriction endonuclease.

CutSmart RR SAM ★ 37° 1654

Concentration: 5,000 units/ml

Note: SAM is unstable at 37°C (pH 7.5) and should be replenished in reactions incubated > 4 hours. For use of methylation reaction SAM should be diluted 1:200 to a final concentration of 160 μ M. EcoGII Methyltransferase is sensitive to salt. Make sure DNA solution is low in salt concentration or that it makes up only a small % of the final reaction volume.

EcoRI Methyltransferase

#M0211S 10,000 units65 €

5´... G A A T T C ... 3´ 3´... C T T A A G ... 5´ CH₃ **Description:** EcoRI Methyltransferase modifies the internal adenine residue (N°) in the sequence to the left.

Reaction Conditions: 1X EcoRI Methyltransferase Reaction Buffer + SAM. Supplement with 80 µM S-adenosylmethionine (supplied). Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

RX NEBU SAM 37° 1654

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μ g λ DNA in 1 hour at 37°C in a total reaction volume of 10 μ l against cleavage by EcoRl restriction endonuclease.

Concentration: 40.000 units/ml

Note: EcoRI Methyltransferase is inhibited by MgCl₂. Only 50% activity is retained at a concentration of 4 mM MgCl₂.

HaeIII Methyltransferase

#M0224S 500 units73 €

CH₃ 5'... G G C C ... 3' 3'... C C G G ... 5' **Description:** HaelII Methyltransferase modifies the internal cytosine residue (C⁵) in the sequence to the left.

Reaction Conditions: 1X HaellI Methyltransferase Reaction Buffer + SAM. Supplement with 80 μM S-adenosylmethionine (supplied). Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

RX NEBU SAM 37° 1654

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μ g λ DNA in 1 hour at 37°C in a total reaction volume of 10 μ l against cleavage by Haelll restriction endonuclease.

Concentration: 10,000 units/ml

Note: HaellI Methyltransferase protects DNA against cleavage by Notl.

HhaI Methyltransferase

#M0217S 1,000 units73 €

5'... G C G C ... 3' 3'... C G C G ... 5' **Description:** Hhal Methyltransferase modifies the internal cytosine residue (C⁵) in the sequence to the left.

Reaction Conditions: 1X Hhal Methyltransferase Reaction Buffer + SAM. Supplement with $80 \mu M$ S-adenosylmethionine (supplied). Incubate at $37^{\circ}C$.

RR NEBU SAM 37° WAS

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μ g λ DNA in 1 hour at 37°C in a total reaction volume of 10 μ l against cleavage by Hhal restriction endonuclease.

Concentration: 25,000 units/ml

HpaII Methyltransferase

#M0214S

100 units69 €

CH. 5'... C C G G ... 3' 3'... GGCC...5' CH,

Description: Hpall Methyltransferase recognizes the same sequence as the Mspl Methyltransferase, but modifies the internal cytosine residue (C5) in the sequence to the left.

Reaction Conditions: 1X Hpall Methyltransferase Reaction Buffer + SAM. Supplement with 80 µM S-adenosylmethionine (supplied). Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

RX NEBU SAM 37° 1654

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μ g λ DNA in 1 hour at 37°C in a total reaction volume of 10 µl against cleavage by Hpall restriction endonuclease.

Concentration: 4,000 units/ml

MspI Methyltransferase

#M0215S

100 units73 €

CH₃ 5'... C'CGG...3' 3'... GGCC...5' CH₃

Description: Mspl Methyltransferase recognizes the same sequence as the Hpall Methyltransferase, but modifies the external cytosine residue (C5) in the sequence to the left

Reaction Conditions: 1X Mspl Methyltransferase Reaction Buffer + SAM. Supplement with 80 µM S-adenosylmethionine (supplied). Incubate at 37°C.

RX NEBU SAM 37° WW

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg λ DNA in 1 hour at 37°C in a total reaction volume of 10 µl against cleavage by Mspl restriction endonuclease.

Concentration: 5,000 units/ml

TaqI Methyltransferase

#M0219S

1,000 units73 €

5'... T C G A ... 3' 3′... A G C T ... 5′ CH₃

Description: *Taq* I Methyltransferase modifies the adenine residue (N6) in the sequence to the left.

Reaction Conditions: 1X CutSmart Buffer + SAM. Supplement with 80 µM S-adenosylmethionine (supplied). Incubate at 65°C.

Note: Activity at 37°C is 25%.

CutSmart RR SAM 65° Mb

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg λ DNA in 1 hour at 65°C in a total reaction volume of 20 µl against cleavage by Taq l restriction endonuclease

Concentration: 10.000 units/ml

Human DNA (cytosine-5) Methyltransferase (Dnmt1)

RX

See page 270 for more information.

#M0230S #M0230L

50 units 123 € 250 units 492 €

Dnmt1 methylates cytosine residues in hemimethylated DNA at 5'...CG...3'. Mammalian Dnmt1 is believed to be involved in carcinogenesis, embryonic development

and several other biological functions. The bulk of the methylation takes place during DNA replication in the S-phase of the cell cycle.

Colby and Andrew are the newest members of our Global Business Development group. Andrew specializes in next generation sequencing, while Colby brings expertise in RNA structural biology, metabolic engineering and microbiology.





















RecA & RecA

RecA #M0249S 200 μg72 € #M0249L 1,000 μg288 €

RecA_f

#M0355S 200 µg72 € #M0355L 1,000 µg 288 €

- Visualization of DNA structures with electron microscopy
- D-loop mutagenesis
- Screening libraries using RecA-coated probes
- Cleavage of DNA at a single predetermined site
- RecA mediated affinity capture for full length cDNA cloning

RX NEBU 37° 🛞

Reaction Conditions: 1X RecA Reaction Buffer. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Molecular Weight: RecA: 37,973 daltons

RecA_t: 38,907 daltons. **Concentration:** 2 mg/ml

duplex DNA and searches for a homologous region, (iii) the strands are exchanged.

RecA, is a N-terminal 6X His tagged recombinant protein.

Description: E. coli RecA is necessary for genetic

recombination, reactions involving DNA repair and UV-

induced mutagenesis. RecA promotes the autodiges-

tion of the lexA repressor, umuD protein and lambda

repressor. Cleavage of LexA derepresses more than 20

genes. In vitro studies indicate that in the presence of ATP,

RecA promotes the strand exchange of single-strand DNA fragments with homologous duplex DNA. The reaction has three distinct steps: (i) RecA polymerizes on the single-strand DNA, (ii) the nucleoprotein filament binds the

T4 Gene 32 Protein

#M0300S 100 μg78 € #M0300L 500 μg312 €

- Increase yield and processivity of reverse transcriptase during RT-PCR
- Increase yield and specificity of PCR products from soil samples
- Stabilization and marking of ssDNA structures

Description: T4 Gene 32 Protein is a single-stranded DNA (ssDNA) binding protein required for bacteriophage T4 replication and repair. It cooperatively binds to and stabilizes transiently formed regions of ssDNA and plays an important structural role during T4 phage replication. It also has been used extensively to stabilize and mark regions of ssDNA for electron microscopic examination of intracellular DNA structures. Recently, it has been shown to improve restriction enzyme digestion, improve the yield and efficiency of reverse transcription (RT) reactions during RT-PCR, enhance T4 DNA polymerase activity, as well as increase the yield of PCR products.

RX NEB 4 37° 655

Reaction Conditions: 1X NEBuffer 4. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: Sold by mass of pure protein as determined by OD_{280} (A_{280} =1.184 at 1 mg/ml, 1 cm)

Molecular Weight: 33,485 daltons

Concentration: 10 mg/ml

ET SSB (Extreme Thermostable Single-Stranded Binding Protein)

#M2401S 50 μg 161 €

- Improve the processivity of DNA polymerase
- Stabilization and marking of ssDNA structure
- Increase the yield and specificity of PCR
- Increase the yield and processivity of RT during RT-PCR
- Improve DNA sequencing through regions with strong secondary structure

Description: ET SSB (Extreme Thermostable Single-Stranded DNA Binding Protein) is a single-stranded DNA binding protein isolated from a hyperthermophilic microorganism, which remains fully active after incubation at 95°C for 60 minutes. Due to the extreme thermostability, ET SSB can be used in applications that require extremely high temperature conditions, such as nucleic acid amplification and sequencing.

Unit Definition: Sold by mass of pure protein as determined by OD_{280} ($A_{280} = 0.774$ at 1 mg/ml, 1 cm).

RX

Concentration: 500 µg/ml Molecular Weight: 16 kDa

Usage Note: ET SSB is active in any polymerase buffer. Add 200 ng of ET SSB per 50 µl reaction.

Cloning Plasmids and DNAs

CLONING PLASMID/DNA	NEB #	FEATURES	CONCENTRATION	MW/SIZE	SIZE	PRICE
pBR322 Vector	N3033S/L	Commonly used cloning vectors	1,000 μg/ml	2.83 x 10 ⁶ Da/4,361 bp	50/250 μg	73 €/292 €
pUC19 Vector	N3041S/L	Tet, Amp resistance	1,000 μg/ml	1.75 x 10 ⁶ Da/2,686 bp	50/250 μg	74 €/296 €
M13mp18 RF I DNA	N4018S	 Phage vectors derived from bacteriophage M13 DNA, covalently closed circular 13 Unique RE sites with β-gal gene Blue/white selection 	100 μg/ml	7,249 bp	10 µg	86€
M13mp18 Single-stranded DNA	N4040S		250 μg/ml	7,249 bp	10 µg	39€
Lambda DNA	N3011S/L		500 μg/ml	31.5 x 10 ⁶ Da/48,502 bp	250/1,250 μg	68 €/272 €
Lambda DNA (<i>N®</i> -methyladenine-free)	N3013S/L	Commonly used DNA substrate	500 μg/ml	31.5 x 10 ⁶ Da/48,502 bp	250/1,250 µg	68 €/272 €
φX174 RF I DNA	N3021S/L	Commonly used DNA substrate Covalently closed circular form of \$\phi\$X174	1,000 µg/ml	3.5 x 10 ⁶ Da/5,386 bp	30/150 µg	71 €/284 €
φX174 RF II DNA	N3022L	Commonly used DNA substrate Double-stranded nicked circular form of	1,000 µg/ml	3.5 x 10 ⁶ Da/5,386 bp	150 μg	284 €
φX174 Virion DNA	N3023S/L	Single-stranded viral DNA	1,000 μg/ml	1.7 x 10 ⁶ Da/5,386 bp	50/250 μg	78 €/312 €

NEB offers a selection of common cloning plasmids and DNAs for use as substrates.

Additional information for many of these DNAs can be found in the technical reference section of this catalog or at www.neb.com.

M13KO7 Helper Phage

#N0315S

1.8 ml 108 €

 Production of single-stranded phagemid DNA for sequencing and mutagenesis

Description: M13K07 is a derivative of M13 phage with the origin of replication from P15A and the kanamycin resistance gene from Tn903 both inserted within the M13 origin of replication, which is able to replicate in the absence of phagemid DNA. In the presence of phagemid bearing a wild-type M13 or f1 origin, single-stranded phagemid is packaged preferentially and secreted into the culture medium. This allows easy production of single-stranded phagemid DNA for mutagenesis or sequencing. M13K07 carries the kanamycin resistance marker.

Source: M13K07 phage supernatant is isolated from infected E. coli ER2738 by a standard procedure.

Concentration: 1.0 x 10¹¹ pfu/ml

Note: NEB does not recommend the use of M13K07 as a cloning vector. For cloning peptides displayed on M13 phage, we recommend the Ph.D.™ Peptide Display Cloning System (see page 251).

Companion Product:

Anti-M13 plll Monoclonal Antibody

Anti-M13 plll Monoclonal Antibody (mouse isotype IgG2a) is derived from the A23 hybridoma produced by the fusion of mouse myeloma cells and splenocytes from BALB/c mice immunized with the C-terminal half of M13 coat protein III (residues 259-406).

Note: ELISA against whole phage using this antibody is not recommended since the epitope lies in a region that is not accessible on intact M13 phage virions.

#E8033S

0.1 ml 278 €





















Competent Cell Selection Chart for Cloning

	NEB 5-alpha Competent E. coli (#C2987)	NEB Turbo Competent <i>E. coli</i> (#C2984)	NEB 5-alpha F´ / ^q Competent E. coli (#C2992)	NEB 10-beta Competent E. coli (#C3019)	dam ⁻ /dcm ⁻ Competent E. coli (#C2925)	NEB Stable Competent E. coli (#C3040)
FEATURES						
Versatile	•					•
Fast growth (< 8 hours)		•				
Toxic gene cloning		•	•			•
Large plasmid/BAC cloning				•		
Dam/Dcm-free plasmid growth					•	
Retroviral/lentiviral vector cloning						•
recA-	•		•	•		•
endA-	•	•	•	•	•	•
FORMATS						
Chemically competent	•	•	•	•	•	•
Electrocompetent	•	•		•		
Subcloning	•					
96-well format*	•					
384-well format*	•					
12 x 8-tube strips*	•					

^{*} Other strains are available upon request. For more information, contact custom@neb.com.

Monarch® Nucleic Acid DNA Purification Kits

Monarch kits are also available for Genomic DNA & RNA extraction, see pages 128–130.

Monarch kits provide fast and reliable purification of high quality DNA from bacterial cultures, agarose gels, and enzymatic reactions using best-in-class technology. Our unique column design offers elution in lower volumes than standard purification kits, providing concentrated, high quality DNA suitable for use in downstream applications such as DNA sequencing, PCR, restriction enzyme digests and other enzymatic manipulations.

Our column design also eliminates buffer retention and the risk of carryover contamination, providing fast, worry-free DNA purification. Designed with sustainability in mind, Monarch kits use significantly less plastic and responsibly-sourced packaging. Learn more about Monarch Nucleic Acid Purification Products on pages 122-131.

PRODUCT	NEB#	SIZE	PRICE
Monarch Plasmid Miniprep Kit	T1010S/L	50/250 preps	75 €/325 €
Monarch DNA Gel Extraction Kit	T1020S/L	50/250 preps	88 €/398 €
Monarch PCR & DNA Cleanup Kit (5 μg)	T1030S/L	50/250 preps	88 €/398 €
COLUMNS AVAILABLE SEPARATELY			
Monarch DNA Cleanup Columns (5 μg)	T1034L	100 columns	128 €
Monarch Plasmid Miniprep Columns	T1017L	100 columns	84 €
BUFFERS AVAILABLE SEPARATELY			
Monarch DNA Cleanup Binding Buffer	T1031L	235 ml	97 €
Monarch DNA Wash Buffer	T1032L	25 ml	32 €
Monarch DNA Elution Buffer	T1016L	25 ml	32 €
Monarch Gel Dissolving Buffer	T1021L	235 ml	99 €
Monarch Plasmid Lysis Buffer (B2)	T1012L	2 x 27 ml	32 €
Monarch Plasmid Neutralization Buffer (B3)	T1013L	110 ml	52 €
Monarch Plasmid Resuspension Buffer (B1)	T1011L	55 ml	32 €
Monarch Plasmid Wash Buffer 1	T1014L	2 x 27 ml	32 €
Monarch Plasmid Wash Buffer 2	T1015L	30 ml	32 €





The Race to Save the Reefs

Coral reefs occupy less than 1% of Earth's ocean floor; however, they are home to 25% of all marine wildlife. Coral is the central infrastructure of reefs, and the trees of an underwater forest, which provides microhabitats, shelter and breeding grounds for thousands of species.

There is an incredible network of collaborative relationships between thousands of reef species that underpin all life in this rich habitat, and each of the residents have a role. First and foremost, coral is made up of millions of polyps that house microalgae (*Zooxanthellae*), which photosynthesize and provide 90% of the coral's food. There is also an abundance of unrelated animals hunting together to share a meal, cleaning and protecting each other, recycling waste and defending the reef in a reciprocally altruistic fashion.

Coral is a keystone species, meaning it has a crucial role that no other species in its ecosystem can perform. This role is essential for the survival of the ecosystem, and therefore, if it becomes threatened, it jeopardizes the entire ecosystem.

The threats to the world's reefs include pollution, infectious disease, overfishing and climate change. Rising sea temperatures cause the polyps to eject the microalgae. Without its food source, the coral becomes photobleached and subsequently leaves all of the reef species without their habitat.

In the past 30 years, more than half the world's corals have been affected by bleaching, and the intervals between bleaching events have become shorter, leaving the coral without time to recover. Scientists have concluded that if water temperatures continue to rise at the current rate, all coral reefs will die by the turn of the century.

Can these fragile ecosystems be restored? A glimmer of hope comes from the pioneering work of marine biologists who are manually growing and planting corals that can tolerate higher temperatures and ocean acidification.

Researchers break slow-growing, massive reef-building corals into small pieces, which stimulates rapid healing and growth. The corals are then outplanted back onto the damaged reefs. However, this clonal method of reproduction is not enough; genetic diversity comes from sexual reproduction, and luckily, corals reproduce both asexually and sexually.

Biologists are also collecting eggs and sperm from colonies of Brain coral during spawning events that occur typically one week after a full moon. Swarms of butterflyfish direct the divers to coral that is about to spawn. They collect the sperm and eggs in tents and fertilize them in the lab. The reproductive success rate is approximately 0.2% in nature, but in the lab, it is upwards of 90%.

Other marine biologists are breeding "super coral" by crossing the most robust species of coral in the lab and then transplanting them back onto the reef, where they can hopefully withstand the current stresses that are leading to their decline.

It remains unclear at this point whether the scale of assisted coral transplantation can match that of coral loss, and ocean warming continues to pose a threat to these recovery efforts. Regardless, the passion of these coral enthusiasts and marine biologists to save the Earth's coral reefs is genuinely inspiring.

Nucleic Acid Purification



Time for change.

Nucleic acid purification is an important step in molecular biology workflows. Further, there are many commercially-available solutions from which to choose. Our Research and Development team spent time with customers to better understand what could be done to improve upon current nucleic acid purification kits. This feedback helped us develop our new line of Monarch® Nucleic Acid Purification kits, which have been optimized for maximum performance and minimal environmental impact.

Monarch kits are available for DNA and RNA extraction and cleanup, plasmid purification, and gel extraction. They utilize unique column designs which enable the isolation of highly-pure nucleic acids, free from contaminants and often in low volumes. Monarch kits are supported by a variety of validated, user-friendly protocols to support multiple workflows and applications.

We know that it can be difficult to be environmentally friendly in the lab, where sterility and convenience are of utmost importance. At times, it may seem that sustainability and benchwork are at odds with each other. Although we can't completely solve this problem, we can make changes to our product design to help move toward the goal of a greener lab, and that's exactly what we did with the design of our Monarch kits. These kits use less plastic, as well as responsibly-sourced and recyclable packaging. The columns have thinner walls, reducing total plastic usage without affecting performance. All bottles were carefully chosen to minimize plastic usage, and the kit boxes provide a detailed explanation of how to recycle each component.

Let's work together to clean up the world of nucleic acid purification, one prep at a time.

Featured Products

Monarch Plasmid Miniprep Kit

Monarch Genomic DNA
Purification Kit

130 Monarch Total RNA Miniprep Kit

Featured Tools & Resources



Visit NEBMonarch.com to view our online protocols and tips for optimization of nucleic acid purification



PLASMID PURIFICATION		Monarch gDNA Tissue Lysis Buffer	128
		Monarch gDNA Cell Lysis Buffer	128
Products		Monarch gDNA Blood Lysis Buffer	128
Monarch Plasmid Miniprep Kit	125	Monarch gDNA Binding Buffer	128
Componion Draduata		Monarch gDNA Wash Buffer	128
Companion Products		Monarch gDNA Elution Buffer	128
Monarch Plasmid Miniprep Columns	125	Monarch RNase A	128
Monarch Plasmid Resuspension Buffer (B1)	125	Proteinase K, Molecular Biology	128
Monarch Plasmid Lysis Buffer (B2)	125	Trotomaso K, Woloodiai Biology	
Monarch Plasmid Neutralization Buffer (B3)	125	RNA PURIFICATION	
Monarch Plasmid Wash Buffer 1	125	Products	
Monarch Plasmid Wash Buffer 2	125		
Monarch DNA Elution Buffer	125	Monarch Total RNA Miniprep Kit	130
		Monarch RNA Cleanup Kit (10 μg)	131
MA CLEANUD		Monarch RNA Cleanup Kit (50 μg)	131
DNA CLEANUP		Monarch RNA Cleanup Kit (500 μg)	131
Products		Companion Products	
Monarch DNA Gel Extraction Kit	126	Monarch RNA Purification Columns	130
Monarch PCR & DNA Cleanup Kit (5 μg)	127	Monarch gDNA Removal Columns	130
		Monarch DNA/RNA Protection Reagent	130
Companion Products		Monarch RNA Lysis Buffer	130
Monarch DNA Cleanup Columns (5 μg)	126, 127	Monarch Total RNA Miniprep Enzyme Pack	130
Monarch Gel Dissolving Buffer	126	Monarch RNA Priming Buffer	130
Monarch DNA Cleanup Binding Buffer	127	Monarch RNA Wash Buffer	130
Monarch DNA Wash Buffer	126, 127	Monarch RNA Cleanup Columns (10 µg)	131
Monarch DNA Elution Buffer	126, 127	Monarch RNA Cleanup Columns (50 µg)	131
		Monarch RNA Cleanup Columns (500 µg)	131
NEW COLUMN THE PROPERTY OF THE		Monarch Collection Tubes II	130, 131
GENOMIC DNA PURIFICATION		Monarch RNA Cleanup Binding Buffer	131
Products		Monarch RNA Cleanup Wash Buffer	131
Monarch Genomic DNA Purification Kit	128	Nuclease-free Water	130, 131
Companion Products			
Monarch gDNA Purification Columns	128	ACCESSORIES	
Assessed Astronomy Tolera II		Manarah Migrafuga Tuba Faa Baak	107

Make the right choice and migrate to Monarch®

Purification of DNA and RNA is an essential step in many molecular biology workflows, including enzyme digests, transformation, electrophoresis, PCR, qPCR, RT-PCR, RT-qPCR and library preparation for next gen sequencing. Monarch kits enable quick and easy purification of high-quality DNA and RNA, suitable for use in a variety of downstream applications. Recover pure, intact DNA and RNA in minutes with fast, user-friendly protocols and optimized buffer systems, and focus your time on the experiments that will drive your research forward. Monarch kits are all designed with sustainability in mind; whenever possible, kits and components are made with significantly less plastic and are packaged with responsibly-sourced, recyclable packaging.

Experience exceptional performance and streamlined workflows

- Efficient extraction of high quality DNA and RNA from a variety of samples
- Simplified DNA and RNA cleanup in low elution volumes
- Enhanced column designs for improved performance
- Fast, user-friendly protocols
- Optimized buffer systems

Choose Monarch Kits for pure value

- Buffers and columns available separately
- No additional shipping or handling charges**
- No hazardous materials fees**
- Competitive pricing

Reduce your impact on the environment

- Less plastic used in product design
- Responsibly sourced and recyclable packaging
- Packaging and protocol cards are printed with water and soy-based inks
- Reusable kit boxes made from post-consumer content

To learn more, visit NEBMonarch.com



* Visit www.NEBMonarchPackaging.com for details.



^{**} In the US and select subsidiary locations. Contact your local distributor for shipping policies.

Monarch Plasmid Miniprep Kit

#T1010S 50 preps75 € #T1010L 250 preps325 €

Companion Products:

Monarch Plasmid Miniprep #T1017L 100 columns + collection	
Monarch Plasmid Resuspen	sion Buffer (B1)
#T1011L	55 ml32 €
Monarch Plasmid Lysis Buff	er (B2)
#T1012L	2 x 27 ml32 €
Monarch Plasmid Neutraliza	tion Buffer (B3)
#T1013L	110 ml52 €
Monarch Plasmid Wash Buff	fer 1
#T1014L	2 x 27 ml32 €
Monarch Plasmid Wash Buff	fer 2
#T1015L	30 ml32 €
Monarch DNA Elution Buffer #T1016L	25 ml32 €
	n Kit 50 preps 88 € 50 preps 398 €
	up Kit (5 μg) 50 preps 88 € 50 preps 398 €


- Prevent buffer retention and salt carry-over with optimized column design
- Monitor completeness of certain steps using colored buffer system
- No need to add RNase before starting
- Elute in low volumes
- Purchase optimized kit formats or buffers & columns separately for your convenience
- Easily label columns using tab and frosted surfaces

Description: The Monarch Plasmid Miniprep Kit is a rapid and reliable method for the purification of high-quality plasmid DNA. This method employs standard cell resuspension, alkaline lysis and neutralization steps, with the additional benefit of color indicators at certain steps to easily monitor completion. After clarification of the lysate by centrifugation, the DNA is bound to the proprietary silica matrix under high salt conditions. Unique wash buffers ensure salts, proteins, RNA and other cellular components (endotoxins) are removed, allowing low-volume elution of concentrated, highly pure DNA, ready for use in restriction digests, DNA sequencing, PCR and other enzymatic manipulations.

The Monarch Plasmid Miniprep Kit Includes:

- Monarch Plasmid Miniprep Columns
- Monarch Plasmid Resuspension Buffer (B1)
- Monarch Plasmid Lysis Buffer (B2)
- Monarch Plasmid Neutralization Buffer (B3)
- Monarch Plasmid Wash Buffer 1
- Monarch Plasmid Wash Buffer 2 (5X)
- Monarch DNA Elution Buffer
- Monarch Collection Tubes

For your convenience, Monarch kit components, including columns and buffers, are available separately.



Monarch Plasmid Miniprep Kits consistently yield more concentrated plasmid DNA with equivalent purity and functionality as the leading supplier. Preps were performed according to recommended protocols using 1.5 ml aliquots of the same overnight culture. One microliter of each prep was digested with HindIII-HF (NEB #R3104) to linearize the vector and the digests were resolved on a 1% w/v agarose gel.



Monarch DNA Gel Extraction Kit

#T1020S	50 preps	88 €
#T1020L	250 preps	398€

Companion Products:

#T1034L 100 columns + co	(10)
Monarch Gel Dissolving E #T1021L	35 ml 99 €
Monarch DNA Wash Buffe	er
#T1032L	25 ml32 €
Monarch DNA Elution Buf	fer
#T1016L	25 ml32 €
Monarch Plasmid Minipre	ep Kit
#T1010S	50 preps 75 €
#T1010L	250 preps 325 €
Monarch PCR & DNA Cle #T1030S #T1030L	anup Kit (5 µg) 50 preps88 € 250 preps398 €

- Elute in as little as 6 μl
- Prevent buffer retention and salt carry-over with optimized column design
- Save time with fast, user-friendly protocol
- Purchase optimized kit formats or buffers & columns separately for your convenience


With Monarch DNA Cleanup Columns, DNA can be eluted in as little as 6 μ l.

Description: The Monarch DNA Gel Extraction Kit rapidly and reliably purifies up to 5 µg of concentrated high-quality, double-stranded DNA from agarose gels. This method employs a bind/wash/elute workflow with minimal incubation and spin times, resulting in purification in less than 15 minutes. The Monarch Gel Dissolving Buffer is used to digest the agarose gel slice and ensure the sample is compatible for loading the DNA onto the proprietary silica matrix under high salt conditions. The wash buffer ensures trace amounts of DNA binding dyes, electrophoresis buffer salts and gel loading buffer components are removed. Low-volume elution produces concentrated, highly pure DNA ready for use in restriction digests. DNA sequencing, ligation. and other enzymatic manipulations. The unique column design ensures no buffer retention and no carryover of contaminants, allowing elution of sample in volumes as low as 6 µl.


The Monarch DNA Gel Extraction Kit Includes:

- Monarch DNA Cleanup Columns (5 μg)
- Monarch Gel Dissolving Buffer
- Monarch DNA Wash Buffer
- Monarch DNA Elution Buffer
- Monarch Collection Tubes

SPECIFICATIONS	
Binding Capacity	5 μg
DNA Size Range	50 bp-25 kb
Elution Volume	≥ 6 µl
Typical Recovery	DNA 50 bp-10 kb 70-90%
	DNA 11-25 kb 50-70%



Monarch DNA Gel Extraction Kit reproducibly recovers DNA over a broad range of molecular weights. A mixture of 7 DNA fragments ranging from 10 kb down to 0.5 kb was prepared and one-half of the mixture was resolved on a 1% gel. Each fragment was manually excised from the agarose gel and processed using the Monarch DNA Gel Extraction Kit. The entire elution of each fragment was resolved on a new gel with the remainder of the original mixture for comparison.



Our optimized column design supplied with the Monarch Gel Extraction and PCR & DNA Cleanup Kits enables elution in as little as 6 μ l, and eliminates buffer retention

Monarch PCR & DNA Cleanup Kit (5 µg)

#T1030S 50 preps88 € #T1030L 250 preps398 €

Companion Products:

Monarch DNA Cleanup C #T1034L 100 columns + co	olumns (5 µg) ollection tubes 128 €
Monarch DNA Cleanup B #T1031L	inding Buffer 235 ml97 €
Monarch DNA Wash Buffe #T1032L	er 25 ml32 €
Monarch DNA Elution But #T1016L	ffer 25 ml32 €
Monarch Plasmid Minipre #T1010S #T1010L	ep Kit 50 preps 75 € 250 preps 325 €
Monarch DNA Gel Extract #T1020S	tion Kit 50 preps 88 €

Elute in as little as 6 μl

#T1020L

 Prevent buffer retention and salt carry-over with optimized column design

250 preps398 €

- Purify small DNA and oligos with a slight protocol modification
- Save time with fast, user-friendly protocol
- Purchase optimized kit formats or buffers & columns separately for your convenience

With the Monarch PCR & DNA Cleanup Kit, you can purify your DNA in as little as 5 minutes.


Description: The Monarch PCR & DNA Cleanup Kit (5 μg) is a rapid and reliable method for the purification and concentration of up to 5 µg of high-quality, doublestranded DNA from enzymatic reactions such as PCR, restriction digestion, ligation and reverse transcription. This method employs a bind/wash/elute workflow with minimal incubation and spin times, resulting in purification in less than 5 minutes. DNA Cleanup Binding Buffer is used to dilute the samples and ensure they are compatible for loading onto the proprietary silica matrix under high salt conditions. The DNA Wash Buffer ensures enzymes, short primers (≤ 40 nt), detergents and other low-molecular weight reaction components (e.g., nucleotides, DMSO, betaine) are removed, thereby allowing low-volume elution of concentrated, highpurity DNA. Eluted DNA is ready for use in restriction digests, DNA sequencing, ligation and other enzymatic manipulations. The unique column design ensures no buffer retention and no carryover of contaminants, allowing elution of sample in volumes as low as 6 µl. A slight protocol modification enables purification of small DNA and oligonucleotides.

Applications:

- PCR cleanup
- · Enzymatic reaction cleanup
- cDNA cleanup
- Labeling cleanup
- Plasmid cleanup
- Oligo cleanup

The Monarch PCR & DNA Cleanup Kit Includes:

- Monarch DNA Cleanup Columns (5 μg)
- Monarch DNA Cleanup Binding Buffer
- Monarch DNA Wash Buffer
- Monarch DNA Elution Buffer
- Monarch Collection Tubes



Monarch PCR & DNA Cleanup Kit (5 μg) removes low molecular weight primers from dsDNA samples. Three independent amplicons (267 bp, 520 bp, 1003 bp) were spiked with two oligonucleotides (16-mer, 24-mer) to a final concentration of 1 μM. Half of each mix was purified with the Monarch PCR & DNA Cleanup Kit (5 μg) following the included protocol. Equivalent fractions of the original mixture and the eluted material were resolved on a 20% TBE acrylamide gel at 100V for one hour and stained with SYBR® Green II.

Monarch Microfuge Tube EcoRack

#T5020S 2 racks40 €

Description: The Monarch Microfuge Tube EcoRack is a bench-top tube rack made from plastic recovered during the manufacture of Monarch Nucleic Acid Purification Columns. Plastic that would otherwise be discarded during the injection molding process is recovered and re-molded into this useful lab accessory that can hold up to 48 tubes each side. One side can accommodate tubes 1.5-2 mls and the other can accommodate 0.5 ml tubes.





Monarch Genomic DNA Purification Kit

#T3010S	50 preps	153 €
#T3010L	150 preps	395 €

Companion Products:

Monarch gDNA Purification	n Columns
#T3017L 100	D columns130 €
Monarch Collection Tubes #T2018L	II 100 tubes27 €
Monarch gDNA Tissue Lysi	is Buffer
#T3011L	34 ml40 €
Monarch gDNA Cell Lysis I	Buffer
#T3012L	20 ml 40 €
Monarch gDNA Blood Lysis	s Buffer
#T3013L	20 ml 40 €
Monarch gDNA Binding Bu	ıffer
#T3014L	65 ml110 €
Monarch gDNA Wash Buffe	er
#T3015L	60 ml32 €
Monarch gDNA Elution Buf	ffer
#T3016L	34 ml32 €
Monarch RNase A #T3018L	1 ml61 €
Proteinase K, Molecular Bio	ology
#P8107S	2 ml78 €


- Purify high quality gDNA from a wide variety of sample types (cells, blood, tissues, and more)
- Experience extremely low residual RNA contamination (typically < 1%)
- Isolate high molecular weight gDNA (peak size typically ≥ 50 kb)
- Take advantage of user-friendly protocols with fast and efficient lysis steps
- Additionally, use the kit to clean up genomic DNA
- Enjoy the flexibility to purchase kit components separately

The Monarch Genomic DNA Purification Kit is an excellant complement to the NEBNext Library Preparation products for NGS. The large peak sizes of purified DNA makes it an exceptional choice for purification upstream of long read sequencing platforms.


Description: The Monarch Genomic DNA Purification Kit is a comprehensive solution for cell lysis, RNA removal, and purification of intact genomic DNA (gDNA) from a wide variety of biological samples, including cultured cells, blood, and mammalian tissues. Additionally, bacteria and yeast can be processed with extra steps to enhance lysis in these tough-to-lyse samples. Protocols are also included to enable purification from clinically-relevant samples, such as saliva and cheek swabs, as well as rapid cleanup of previously extracted gDNA. Purified gDNA has high quality metrics, including $A_{260/280} > 1.8$ and A $_{260/230}$ > 2.0, high DIN scores and minimal residual RNA. The purified gDNA is suitable for downstream applications, such as endpoint PCR, qPCR and library prep for next generation sequencing (NGS). Typical peak size is 50-70 kb, making this kit an excellent choice upstream of long-read sequencing platforms.

The Monarch gDNA Purification Kit Includes:


- Monarch gDNA Purification Columns
- Monarch Collection Tubes II
- Monarch gDNA Tissue Lysis Buffer
- Monarch gDNA Cell Lysis Buffer
- Monarch gDNA Blood Lysis Buffer
- Monarch gDNA Binding Buffer
- Monarch gDNA Wash Buffer
- Monarch gDNA Elution Buffer
- Monarch RNase A
- Proteinase K, Molecular Biology



SPECIFICATIONS	
Input	Cultured mammalian cells: up to 5 x 10 ⁶ cells Mammalian whole blood: 100 µl Tissue: up to 25 mg, depending on tissue type Bacteria: up to 2 x 10 ⁶ Yeast: up to 5 x 10 ⁷ Saliva: up to 500 µl Buccal swabs Genomic DNA requiring cleanup
Binding Capacity	30 μg genomic DNA
Yield	Varies depending on sample type, see "Guidelines for Choosing Sample Input Amounts", (page 355)
Genomic DNA Size	Peak size > 50 kb for most sample types; may be lower for saliva and buccal swabs
RNA Content	< 1% (with included RNase A treatment)
Purity	$A_{260/280} \ge 1.8, A_{260/230} \ge 2.0$



The Monarch Genomic DNA Purification Kit efficiently purifies high-quality, high molecular weight gDNA from a variety of sample types. 100 ng of genomic DNA from each sample was loaded on a 0.75% agarose gel. gDNA was isolated following the standard protocols for blood, cultured cells and tissue, and the supplemental protocols for buccal swabs, saliva, Gram- and Gram+ bacteria. Starting material used: 1 x 10° HeLa cells, 100 µl human blood, 10 µl bird blood, 10 mg frozen tissue powder, 1 buccal swab, 500 µl saliva and ~1 x 10° bacterial cells. Lambda DNA-Hind III digest (NEB #N3012) was used as a marker in the last lane (M). Purified gDNA samples were analyzed using a Genomic DNA ScreenTape® on an Agilent Technologies® 4200 TapeStation®. Samples typically yield peak sizes 50—70 kb and DINs of -9. The cell fractions processed in the buccal swab and saliva preps contain dead cells, as expected, causing a smear like pattern with typical low molecular weight apoptotic bands.



The Monarch Genomic DNA Purification Kit generates high quality genomic DNA suitable for sensitive applications like long range PCR and gPCR.

A: Amplification reactions were set up with primer pairs specific for 6, 8, 10, 12, 16, 20 kb amplicons from human DNA. LongAmp® Hot Start Taq 2X Master Mix (NEB #M0533) was used and 25 ng template DNA was added to each sample. PCR reactions were carried out on an Applied Biosystems® 2720 Thermal Cycler. Monarch-purified gDNA isolated from HeLa cells and human blood were compared to commercially available reference DNA from the human cell line NA19240 F11. 10 µl was loaded on a 1.5% agarose gel, using the 1 kb DNA Ladder (NEB #N3232) as a marker. Results indicated DNA was of high-integrity and suitable for long range PCR.

B: Monarch-purified gDNA from human whole blood, HeLa cells and mouse tail was diluted to produce a five log range of input template concentrations. The results were generated using primers targeting gHEME (human whole blood) and gREL (HeLa, mouse tail) for qPCR assays with the Luna® Universal qPCR Master Mix (NEB #M3003) and cycled on a Bio-Rad® CFX Touch qPCR thermal cycler. Results indicated that DNA is highly pure and free from inhibitors, optimal for qPCR.



Monarch Total RNA Miniprep Kit

#T2010S 50 preps 248 €

Companion Products:

Monarch RNA Purification Columns #T2007L 100 columns + collection tubes 153 €
Monarch gDNA Removal Columns #T2017L 100 columns + collection tubes 128 €
Monarch Collection Tubes II #T2018L 100 collection tubes 27 €
Monarch DNA/RNA Protection Reagent #T2011L 56 ml 91 €
Monarch RNA Lysis Buffer #T2012L 100 ml 61 €
Monarch Total RNA Miniprep Enzyme Pack #T2019L 1 pack 61 €
Monarch RNA Priming Buffer #T2013L 56 ml61 €
Monarch RNA Wash Buffer #T2014L 50 ml 61 €
Nuclease-free Water #B1500S 25 ml 25 € #B1500L 100 ml 61 €

- Use with a wide variety of sample types
- Purify RNA of all sizes, including miRNA & small RNAs > 20 nucleotides
- Includes DNase I, gDNA removal columns, Proteinase K, and a stabilization reagent
- Efficiently remove contaminating genomic

 DNA
- Protocols available for RNA fractionation and RNA cleanup
- Save money with value pricing for an all-in-one kit

See our "Guidelines for Choosing Sample Input Amounts," including expected yields on page 357.


Description: The Monarch Total RNA Miniprep Kit is a comprehensive solution for sample preservation, cell lysis, gDNA removal, and purification of total RNA from a wide variety of biological samples, including cultured cells, blood, and mammalian tissues. Additionally, tough-to-lyse samples, such as bacteria, yeast, and plant, can be processed with additional steps that enhance lysis. Cleanup of enzymatic reactions or purification of RNA from TRIzol® -extracted samples is also possible using this kit. Purified RNA has high quality metrics, including $\rm A_{260/280}$ and $\rm A_{260/230}$ ratios > 1.8, high RIN scores, and minimal residual gDNA. Captured RNA ranges in size from full-length rRNAs down to intact miRNAs. Additionally, differential binding conditions allow selective capture or exclusion of the sub-200 nucleotide RNA pool that includes miRNA, 5S rRNA, and tRNA. Purified RNA is suitable for downstream applications, such as RT-qPCR, cDNA synthesis, RNA-seq, Northern blot analysis, etc.

The Monarch Total RNA Miniprep Kit Includes:

- Monarch gDNA Removal Columns
- Monarch RNA Purification Columns
- Monarch Collection Tubes II
- Monarch DNA/RNA Protection Reagent (2X)
- Monarch RNA Lysis Buffer
- Monarch Proteinase K
- Monarch Proteinase K Resuspension Buffer
- Monarch Proteinase K Reaction Buffer
- Monarch DNase I
- Monarch DNase I Reaction Buffer
- Monarch RNA Priming Buffer
- Monarch RNA Wash Buffer (5X)
- Monarch Nuclease-free Water

SPECIFICATIONS	
Binding Capacity	100 μg RNA
RNA Size	≥ 20 nt
Purity	A _{260/280} and A _{260/230} usually ≥ 1.8
Input Amount	up to 107 cells or 50 mg tissue*
Elution Volume	30–100 µl
Yield	varies depending on sample type
Compatible downsteam applications	RNA Library prep for NGS, RT-PCR, RT-qPCR, Northern blots

^{*}See "Guidelines for Choosing Sample Input Amounts" on page 357.



Monarch-purified RNA is high-quality and compatible with a wide variety of downstream applications. Total RNA from a broad array of sample types was purified using the Monarch Total RNA Miniprep Kit. Aliquots were run on an Agilent Bioanalyzer® 2100 using the Nano 6000 RNA chip (S. cerevisiae RNA was run using a plant Nano assay). RIN values and O.D. ratios confirm the overall integrity and purity of the RNA. To demonstrate compatibility with downstream applications, samples were subsequently used for RT-PCR (+/- RT) (A) for detection of 4 different RNA species using Protoscript® II Reverse Transcriptase/LongAmp Taq DNA Polymerase, NGS library prep (B) using NEBNext® Ultra* II RNA Library Prep Kit and RT-qPCR (C) using Luna® One-Step RT-qPCR Reagents.

NEW

Monarch RNA Cleanup Kits

Monarch RNA Cleanup Kit (10 μg) #T2030S 10 preps. . 52€ #T2030L 100 preps 284 € Monarch RNA Cleanup Kit (50 ug) #T2040S 10 preps 50 € #T2040L 100 preps 280 € Monarch RNA Cleanup Kit (500 μg) 10 preps 58 € #T2050S #T2050L 100 preps 440 €

Companion Products:

Monarch RNA Cleanup Columns (10 μg) #T2037L 100 columns + collection tubes 200 e
Monarch RNA Cleanup Columns (50 μg) #T2047L 100 columns + collection tubes 198 ·
Monarch RNA Cleanup Columns (500 μg) #T2057L 100 columns + collection tubes 355 ε
Monarch Collection Tubes II #T2018L 100 tubes27
Monarch RNA Cleanup Binding Buffer #T2041S 80 ml86

#T2041S 80 ml86 €

Monarch RNA Cleanup Wash Buffer

#T2042S 40 ml32 €

Nuclease-free Water

#B1500S 25 ml25 €

#B1500L

- Choose from 3 different binding capacities and flexible elution volumes
- Quickly and easily purify large quantities of highquality RNA from in vitro transcription (IVT) reactions

100 ml 61 €

 Efficiently remove unincorporated nucleotides from your RNA sample

Great for RNA cleanup following *in vitro* transcription with HiScribe™ Kits, see pages 185–187.

Description: The Monarch RNA Cleanup Kits provide a fast and simple silica spin column-based solution for RNA cleanup and concentration after any enzymatic reaction (including in vitro transcription, DNase I treatment, capping and labeling) and after other purification methods such as phenol/ chloroform extraction. The Monarch RNA Cleanup Kits are available in 3 different binding capacities: 10 μg , 50 μg and 500 μg . Each kit contains unique columns, all designed to prevent buffer retention and ensure no carryover of contaminants, enabling low-volume elution of highly-pure RNA. Following the standard protocol, RNA \geq 25 nt is purified with this kit; however, a modified protocol is available to enable the binding of RNA as small as 15 nt (including miRNAs).


Applications:

- RNA Cleanup and Concentration (including from the TRIzol aqueous phase)
- Enzymatic Reaction Cleanup
- In vitro Transcription Cleanup
- RNA Gel Extraction
- · RNA Fractionation

The Monarch RNA Cleanup Kits Include:

- Monarch RNA Cleanup Columns (10, 50 or 500 μg)
- Monarch RNA Cleanup Binding Buffer
- Monarch RNA Cleanup Wash Buffer
- Monarch Collection Tubes II
- Nuclease-free Water

Monarch RNA Cleanup Kit	NEB #T2030 (10 μg)	NEB #T2040 (50 μg)	NEB #T2050 (500 μg)
Binding Capacity	10 µg	50 μg	500 μg
RNA Size Range	≥ 25	nt (≥ 15 nt with modified prot	ocol)
Typical Recovery		70-100%	
Eluion Volume	6–20 µl	20–50 μΙ	50–100 μl
Purity		A _{260/280} > 1.8 and A _{260/230} > 1.8	
Protocol Time	5 minutes of spin a	and incubation time	10–15 minutes of spin and incubation time



Recovery of RNA from Monarch RNA Cleanup Kits with Varying Elution Volumes. rRNA (10, 50 or 500 μg, respectively of 16S and 23S Ribosomal Standard from E. coli, Sigma) was purified using a Monarch RNA Cleanup Kit (10 μg, NEB #72030) (50 μg, NEB #72040) (500 μg, NEB #72050). Nuclease-free water was used to elute the RNA. The percent recovery of the RNA was calculated from the resulting A₂₀₀ as measured using a Trinean® DropSense® 16. ~80% of RNA can be efficiently recovered in 6 μl from the Monarch RNA Cleanup Kit (10 μg, NEB #72030), 20 μl from the Monarch RNA Cleanup Kit (50 μg, NEB #72040), and 50 μl from the Monarch RNA Cleanup Kit (500 μg, NEB #72050).





Choking the Oceans with Plastic

By the year 2050, it is predicted that the Earth's oceans will be home, by weight, to more plastic than fish. Eight million metric tons of plastic trash enters the sea every year. This plastic collects in all the oceans on Earth in large, circular currents, or gyres. The largest gyre has been coined The Great Pacific Garbage Patch. Many estimates are made to convey the immensity of this patch — it is 1.6 million km² or 618,000 miles², or three times the size of France. In other words, it is big!

Ocean debris consists of large and small pieces of plastic; however, the predominant contaminant is tiny confetti-sized pieces of plastic that have been photodegraded by the sun. This fog of particulate easily enters the food chain, and approximately 12,000 to 14,000 tons are ingested by fish and invertebrates each year, threatening the biodiversity of ocean life.

Much of the media coverage focuses on sea turtles, birds, seals and other marine animals entangled in, or ingesting, the trash. However, the plastic problem also extends to chemical pollution and its worrying effect on the reproductivity of various marine organisms. Previous research on plasticizers, such as bisphenol A (BPA) and phthalates, has shown that they induce endocrine toxicity and reproductive alterations. Furthermore, plastic polymers attract polychlorinate biphenyls (PCBs) and other toxic chemicals.

Pacific Oysters (*Crassostrea gigas*) that were experimentally exposed to polystyrene microparticles at concentrations estimated for the water-sediment interface (a typical oyster habitat) produced fewer and smaller egg cells, less-mobile sperm and slow-growing offspring¹. Additionally, exposure to nanoplastics caused green algae (*Scenedesmus obliquus*) to exhibit reduced growth and chlorophyll concentrations, and a small planktonic crustacean (*Daphnia magna*) to present reduced neonate body size and increased malformations².

Many marine organisms, such as bryozoans and crustaceans, are known to colonize floating wood or seaweed, using them as a raft. Now, non-rafting, invasive species, such as coral pathogens, are using plastic debris to travel to, and damage, new ecosystems. Also "hitching a ride" are communities of microbes that are vastly different than those found in the surrounding seawater, for example, members of the genus *Vibrio*, which is a human pathogen. A marine insect, *Halobates sericeus*, that typically lays its eggs on natural rafts, such as shells and bird feathers, has now been observed laying eggs on plastic debris. The long-term effects of this have yet to be determined.

In the first decade of this century, more plastic was generated than ever before — and every piece produced is likely still here on the planet. Intervention at the source of this problem, specifically a reduction of single-use plastics, clearly needs to be promoted. However, we can find a glimmer of hope in ongoing efforts — massive trawling clean-ups, GPS tagging of trash to model its movement, materials scientists turning their attention to more environmentally-friendly packaging, large companies reducing their waste and consumers who are making more informed decisions every day.

- (1) Sussarellu, R., et al (2016) PNAS, 113, 2430-2435.
- (2) Besseling, E., et al (2014) Environ Sci Technol, 48, 12336–12343.

NEBNext® Reagents for Next Generation Sequencing



The leading reagents for sample preparation for next generation sequencing.

Library preparation is a critical part of the next generation sequencing workflow; successful sequencing requires the generation of high quality libraries of sufficient yield and quality.

As sequencing technologies improve and capacities expand, boundaries continue to be pushed on sample preparation: high performance is required from ever-decreasing input quantities and from samples of lower quality.

To meet these growing challenges, the NEBNext suite of products continues to evolve to support next generation sequencing with sample preparation tools that streamline workflows, minimize inputs, and improve library quality and yields.

NEBNext reagents are available for sample preparation for DNA, RNA, ChIP, FFPE, Small RNA, single cell and microbiome samples, for use with Illumina®, Ion Torrent™ and other sequencing platforms. Products are in user-friendly formats including kits and modules, with bulk or customized formats also available. Adaptors and primers are available separately, for maximized flexibility. Use of NEBNext products has been cited in thousands of peer-reviewed publications.

For additional convenience and cost-effectiveness in high-throughput workflows, NEBNext reagents are also available in bulk and customized formats. For more information, contact NEBSolutions@neb.com.

Featured Products

- NEBNext Ultra™ II DNA FS DNA Library Prep Kits
- 142 NEBNext Ultra II RNA Library Prep Kits
- 144 NEBNext Single Cell/Low Input RNA Kits
- 146 NEBNext Adaptors & Primers
- 151 NEBNext Direct® Custom Ready Panels
- 141 NEBNext Enzymatic Methyl-seq Kit

Featured Tools & Resources



Visit NEBNextSelector.neb.com for help with selecting products.



Publications Related to Sample Prep



Visit www.NEBNext.com to keep up to date on everything NEBNext.





REAGENTS FOR ILLUMINA SEQUENC	ING		Adaptors & Primers	
Product Selection Charts			NEBNext Multiplex Oligos	
DNA Product Selection Chart	1	36	(96 Unique Dual Index Primer Pairs)	146
RNA Product Selection Chart	1	37	NEBNext Multiplex Oligos (Dual Index Primers Set 1)	146
			NEBNext Multiplex Oligos	140
Ultra II DNA			(Dual Index Primers Set 2)	146
NEBNext Ultra II DNA Library Prep Kit			NEBNext Multiplex Oligos (Index Primers Set 1)	146
for Illumina	1	39	NEBNext Multiplex Oligos (Index Primers Set 2)	146
NEBNext Ultra II DNA Library Prep		00	NEBNext Multiplex Oligos (Index Primers Set 3)	146
with Sample Purification Beads	1	39	NEBNext Multiplex Oligos (Index Primers Set 4)	146
NEBNext Ultra II FS DNA Library Prep Kit for Illumina	1	39	NEBNext Multiplex Oligos (96 Index Primers)	146
NEBNext Ultra II FS DNA Library Prep	'	33	NEBNext Multiplex Oligos	
with Sample Purification Beads	1	39	(Methylated Adaptor, Index Primers Set 1)	146
			NEBNext Adaptor Dilution Buffer	146
METHYLATION ANALYSIS			Target Enrichment	
NEBNext Enzymatic Methyl-seg Kit	1	41	NEBNext Direct Cancer HotSpot Panel	150
NEBNext Enzymatic Methyl-seg Conversion Mode	ule 1	41	NEBNext Direct BRCA1/BRCA2	150
NEBNext Q5U Master Mix	1	41	NEBNext Direct Custom Ready Panels	151
NEBNext Multiplex Oligos for Enzymatic Methy	l-seq			
(Unique Dual Index Primer Pairs)	1	41	Library Quantification	
			NEBNext Library Quant Kit for Illumina	156
Ultra II RNA			NEBNext Library Dilution Buffer	156
NEBNext Ultra II Directional RNA Library Prep Kit for Illumina	4	42		
NEBNext Ultra II Directional RNA Library Prep	'	42	REAGENTS FOR ION TORRENT	
with Sample Purification Beads	1	42		
NEBNext Ultra II RNA Library Prep			NEBNext Fast DNA Library Prep Set for Ion Torrent	157
Kit for Illumina	1	42	NEBNext Fast DNA Fragmentation	
NEBNext Ultra II RNA Library Prep			& Library Prep Set for Ion Torrent	157
with Sample Purification Beads	1	42	2, · · ·, · · · · · · · · ·	
Cinale Cell DNA			REAGENTS FOR TRNA DEPLETION	
Single Cell RNA			NEBNext rRNA Depletion Kit (Human/Mouse/Rat)	147
NEBNext Single/Cell Low Input RNA Library Prep Kit for Illumina	1	44	NEBNext rRNA Depletion Kit (Human/Mouse/Rat)	
NEBNext Single/Cell Low Input cDNA Synthesis			with RNA Sample Purification Beads	147
& Amplification Module		44	NEBNext Poly(A) mRNA	
NEBNext Single Cell Lysis Module	1	44	Magnetic Isolation Module	147
			NEBNext Globin RNA & rRNA Depletion Kit	
Small RNA			(Human/Mouse/Rat)	148
NEBNext Multiplex Small RNA Library Prep Set	t		NEBNext Globin RNA & rRNA Depletion Kit	
for Illumina (Set 1)	1	45	(Human/Mouse/Rat) with RNA Sample	
NEBNext Multiplex Small RNA Library Prep Set			Purification Beads	148
for Illumina (Set 2)		45	NEBNext Magnetic Separation Rack	148
NEBNext Multiplex Small RNA Library Prep Kit		AE.		
for Illumina (Index Primers 1-48) NEBNext Small RNA Library Prep Set		45	REAGENTS SUITABLE FOR ANY PLATFO	RM
for Illumina (Multiplex Compatible)	1	45	DNA Separation	
ioi manina (marapiox compansis)			NEBNext Microbiome DNA Enrichment Kit	152
Modules & Enzymes				
NEBNext Ultra II FS DNA Module	- 1	40	DNA Repair	
Ultra II RNA First Strand Synthesis Module	1	43	NEBNext FFPE DNA Repair Mix	153
Ultra II Directional RNA Second Strand				
Synthesis Module	1	43	Modules & Enzymes	
Ultra II Non-directional RNA Second Strand		40	NEBNext Magnesium RNA Fragmentation Module	143
Synthesis Module	1	43	NEBNext dsDNA Fragmentase	154
NEBNext Ultra II End Repair/	40 4	12	NEBNext dsDNA Fragmentase Reaction Buffer v2	
	40, 1 40, 1		NEBNext Ultra II Q5 Master Mix	155
NEBNext Ultra II Ligation Module NEBNext Ultra II Q5 Master Mix		43 55		
WEDNOK! OILIA II QO WASICI IVIIA		00	Magnetic Separation	
			NEBNext Magnetic Separation Rack	148

One or more of these products are covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs. Inc. For more information, please email us at gbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications.


Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at https://www.neb.com/support/terms-of-sale. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB. ILLUMINA® is a registered trademark of Illumina, Inc.

ION TORRENT® is a trademark of Life Technologies, Inc.

Illumina® DNA Product Selection Chart

Use the following chart to determine the best NEBNext® products for your Illumina DNA library prep needs. For the most up-to-date product and pricing information, visit **NEBNext.com**.

For help selecting products, try our online product selection tool at NEBNextSelector.neb.com




^{*} Module and EM-seq Oligos also available. Reagents for original Ultra workflow and standard workflow are also available. See ordering information.

Illumina RNA Product Selection Chart

Use the following chart to determine the best NEBNext products for your Illumina RNA sequencing needs. For the most up-to-date product and pricing information, visit **NEBNext.com**.

For help selecting products, try our online product selection tool at NEBNextSelector.neb.com




^{*} For human/mouse/rat samples.


Even more $from less - NEBNext^{\otimes} Ultra^{TM} II$ DNA Library Prep Kits for Illumina

Are you challenged with trying to obtain higher library yields using ever-decreasing input amounts? Each component in the NEBNext Ultra II DNA Library Prep Kit from NEB has been reformulated, resulting in a several-fold increase in library yield with as little as 500 picograms of input DNA. These advances deliver unprecedented performance, while enabling lower inputs and fewer PCR cycles.

To learn more and to view performance data, visit NEBNextUltrall.com.



The NEBNext Ultra II DNA Library Prep Kit produces the highest yield libraries from a broad range of input amounts. Libraries were prepared from Human NA19240 genomic DNA using the input amounts and numbers of PCR cycles shown. Manufacturers' recommended protocols were followed, with the exception that size selection was omitted.



NEBNext Ultra II produces the highest rates of conversion to adaptor-ligated molecules from a broad range of input amounts. Libraries were prepared from Human NA19240 genomic DNA using the input amounts and library prep kits shown without an amplification step, and following manufacturers' recommendations. qPCR was used to quantitate adaptor-ligated molecules, and quantitation values were then normalized to the conversion rate for Ultra II. The Ultra II kit produces the highest rate of conversion to adaptor-ligated molecules, for a broad range of input amounts.



NEBNext Ultra II DNA and FS DNA Library Prep Kits for Illumina

NEBNext Ultra II DNA Library Prep Kit for Illumina

#E7645S 24 reactions535 € #E7645L 96 reactions2045 €

NEBNext Ultra II DNA Library Prep with Sample Purification Beads

#E7103S 24 reactions 615 € #E7103L 96 reactions 2352 €

NEBNext Ultra II FS DNA Library Kit

for Illumina #E7805S

#E7805L

24 reactions 665 € 96 reactions 2526 €

NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads

#E6177S 24 reactions 730 € #E6177L 96 reactions 2758 €

See ordering information for NEBNext Ultra II modules.

- Get more of what you need, with the highest library yields
- Use to generate high quality libraries even when you have only limited amounts of DNA, with inputs as low as 500 pg
- Prepare libraries from ALL of your samples, including GC-rich targets and FFPE DNA samples
- Improve yields for target enrichment applications
- Save time with streamlined workflows, reduced hands-on time, and automation compatibility, and enjoy the flexibility of kit or module format products
- Access reliable and easy-to-use, scalable enzymatic DNA fragmentation, integrated into the Ultra II DNA workflow with the FS version of the kit
- Enjoy the flexibility and reliability of the gold standard SPRIselect size selection and clean-up beads, supplied in just the amounts you need


Visit **NEBNextUltrall.com** for more information, including our technical note and protocol videos

Description: The NEBNext Ultra II DNA Library Prep Kits for Illumina meet the challenge of constructing high quality libraries from ever-decreasing input quantities. The reagents for each step in the library preparation workflow have been reformulated to enable high yield preparation of high quality libraries from 500 picograms to 1 microgram of input DNA. This new generation of NEBNext reagents uses a fast, streamlined, automatable workflow and enables use of fewer PCR cycles while


also improving GC coverage. The kit is also compatible with PCR-free workflows and is effective with challenging samples such as FFPE DNA.

The Ultra II FS DNA Library Prep Kit combines robust enzymatic DNA fragmentation with end repair and dA-tailing, integrated into a streamlined library prep workflow.

Both the Ultra II DNA and Ultra II FS DNA kits are available with or without SPRIselect® beads.



NEBNext Ultra II provides uniform GC coverage for microbial genomic DNA over a broad range of GC composition and input amounts. Libraries were made using 500 pg, 1 ng and 100 ng of the genomic DNAs shown and the Ultra II DNA Library Prep Kit and sequenced on an Illumina MiSeq®. Reads were mapped using Bowtie 2.2.4 and GC coverage information was calculated using Picard's CollectGCBiasMetrics (v1.117). Expected normalized coverage of 1.0 is indicated by the horizontal grey line, the number of 100 bp regions at each GC% is indicated by the vertical grey bars, and the colored lines represent the normalized coverage for each library.



${\it NEBNext\ Ultra\ II\ provides\ the\ highest\ and\ most\ uniform\ coverage\ of\ difficult\ sequence\ regions.}$

A: Indexed libraries were prepared from 100 ng of Human NA19240 genomic DNA using a PCR-free workflow or the library prep kits shown, following manufacturers' recommendations. The PCR-free library was prepared using NEBNext Ultra II. Libraries were sequenced on the Illumina NextSeq[®] 500. 420 million reads were randomly extracted from each dataset, to produce an average coverage of 10X. Reads were mapped to the GRCh37 reference genome using Bowtie 2.2.4. Reads on each region were counted using bedtools v2.19.1. The number of reads overlapping distinct difficult, low-coverage regions of the human genome (1) are shown for each library. Ultra II provides the highest and most uniform coverage of these difficult regions, and provides the coverage closest to that obtained with a PCR-free protocol.

B: From the 420 million 75 bp reads randomly extracted from each dataset, 10X coverage was expected. The percentage of difficult regions covered at $\geq 10X$ is shown for each library prep kit and for the PCR-free workflow. Ultra II provides the highest percentage of reads at $\geq 10X$ coverage and also provides the coverage closest to that obtained with a PCR-free protocol. (1) Aird, D. et al. (2011). Analyzing and minimizing PCR amplification bias in Illumina sequencing libraries. Genome Biology 12(2), R18.

NEBNext Ultra II DNA Reagents for Illumina Sequencing

NEBNext Ultra II Kits for DNA are available with or without integrated enzymatic DNA fragmentation. Note that adaptors and primers are supplied separately. In addition to stringent QCs on individual components, the NEBNext DNA kits are also functionally validated by preparation of a library, followed by Illumina sequencing.

 $\label{eq:linear_line$

Fragmentation	End Repair/dA-Tailing	Adaptor Ligation	Clean Up/ Size Selection	PCR Enrichment	Clean Up	Tota			
	NEBNext Ultra II DNA Lib	rary Prep Kit for Illumi	na (NEB #E7645)			6			
	Ultra II End Prep Enzyme Mix Ultra II End Prep Reaction Buffer (10X)	Ultra II Ligation Master Mix Ligation Enhancer		NEBNext Ultra II Q5 Master Mix		Ha (not fragr 12 -			
	NEBNext Ultra II Library Prep with Sample Purification Beads (NEB #E7103)								
	Ultra II End Prep Enzyme Mix Ultra II End Prep Reaction Buffer (10X)	Ultra II Ligation Master Mix Ligation Enhancer	NEBNext Sample Purification Beads (SPRIselect)	NEBNext Ultra II Q5 Master Mix	NEBNext Sample Purification Beads (SPRIselect)	1.7 -			
NEBNext Ultra II FS DNA Library Prep Kit for Illumina (NEB #E7805)									
Ultra II FS Enzyme Ultra II FS Reaction		Ultra II Ligation Master Mix Ligation Enhancer		NEBNext Ultra II Q5 Master Mix		Har (ind			
NEBNext Ultra II FS	NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads (NEB #E6177)								
Ultra II FS Enzyme Mix Ultra II FS Reaction Buffer		Ultra II Ligation Master Mix Ligation Enhancer	NEBNext Sample Purification Beads (SPRIselect)	NEBNext Ultra II Q5 Master Mix	NEBNext Sample Purification Beads (SPRIselect)	1.4 –			
NEBNext Ultra II FS DNA Module (NEB #E7810)									
	Ultra II FS Enzyme Mix Ultra II FS Reaction Buffer								
NEBNext dsDNA Fragmentase® (NEB #M0348)	NEBNext Ultra II End Repair/dA-Tailing Module (NEB #E7546)	NEBNext Ultra II Ligation Module (NEB #E7595)		NEBNext Ultra II Q5 Master Mix (NEB #M0544)					
NEBNext dsDNA Fragmentase Reaction Buffer v2 Magnesium Chloride	Ultra II End Prep Enzyme Mix Ultra II End Prep Reaction Buffer (10X)	Ultra II Ligation Master Mix Ligation Enhancer		NEBNext Ultra II Q5 Master Mix					

NEBNext Enzymatic Methyl-seq

NEW NEBNext Enzymatic Methyl-seq Kit #E7120S 24 reactions 895 € #E7120L 96 reactions 3360 €

NEBNext Enzymatic Methyl-seq Conversion

Module #E7125S 24 reactions 185 € #E7125L 96 reactions 675 €

NEW

NEBNext Q5U™ Master Mix

#M0597S 50 reactions 125 € #M0597L 250 reactions 495 €

NEW


NEBNext Multiplex Oligos for Enzymatic
Methyl-seq (Unique Dual Index Primer Pairs)
#E7140S 24 reactions 135 €
#E7140L 96 reactions 535 €

- Superior sensitivity of detection of 5-mC and 5-hmC
- Larger library insert sizes
- More uniform GC coverage
- Greater mapping efficiency
- High-efficiency library preparation


NEBNext Enzymatic Methyl-seq is an enzymatic alternative to bisulfite conversion with superior performance. For more information, including extensive performance data, visit NEBNext.com.

Description: While bisulfite sequencing has been the gold standard for the study of DNA methylation, this conversion treatment is damaging to DNA, resulting in DNA fragmentation, loss and GC bias. The NEBNext Enzymatic Methyl-seq Kit (EM-seq™) provides an enzymatic alternative to whole genome bisulfite sequencing (WGBS), combined with highefficiency streamlined library preparation suitable for Illumina sequencing.

The highly effective EM-seq enzymatic conversion minimizes damage to DNA and, in combination with the supplied NEBNext Ultra II library preparation workflow reagents, results in high quality libraries that enable superior detection of 5-mC and 5-hmC from fewer sequencing reads.



NEBNext Enzymatic Methyl-seq libraries have larger insert sizes 50 ng Human NA12878 genomic DNA was sheared to 300 bp using the Covaris® S2 instrument and used as input into EM-seq and WGBS protocols. For WGBS, NEBNext Ultra II DNA was used for library construction, followed by the Zymo Research EZ DNA Methylation-Gold™ kit for bisulfite conversion. Libraries were sequenced on an Illumina MiSeq (2 x 76 bases) and insert sizes were determined using Picard 2.18.14. The normalized frequency of each insert size was plotted, illustrating that library insert sizes are larger for EM-seq than for WGBS, and indicating that EM-seq does not damage DNA as bisulfite treatment does in WGBS.



EM-seq identifies more CpGs than WGBS, at lower sequencing coverage depth with superior uniformity of GC coverage. 10, 50 and 200 ng Human NA12878 genomic DNA was sheared to 300 bp using the Covaris S2 instrument and used as input into EM-seq and WGBS protocols. For WGBS, NEBNext Ultra II DNA was used for library construction, followed by the Zymo Research EZ DNA Methylation-Gold Kit for bisulfite conversion. Libraries were sequenced on an Illumina NovaSeq® 6000 (2 x 100 bases). Reads were aligned to hg38 using bwa-meth 0.2.2.

A: Coverage of CpGs with EM-seq and WGBS libraries was analyzed using 324 million paired end reads, and each top and bottom strand CpGs were counted independently, yielding a maximum of 56 million possible CpG sites. EM-seq identifies more CpGs at lower depth of sequencing.

B: GC coverage was analyzed using Picard 2.17.2 and the distribution of normalized coverage across different GC contents of the genome (0-100%) was plotted. EM-seq libraries have significantly more uniform GC coverage, and lack the AT over-representation and GC under-representation typical of WGBS libraries.

NEBNext Ultra II Library Prep Kits for RNA

NEBNext Ultra II Directional RNA Library Prep Kit for Illumina

#E7760S 24 reactions1028 € #E7760L 96 reactions3495 €

NEBNext Ultra II Directional RNA

Library Prep with Sample Purification Beads
#E7765S 24 reactions1140 €
#E7765L 96 reactions3880 €

NEBNext Ultra II RNA Library Prep Kit for Illumina

#E7770S 24 reactions 980 € #E7770L 96 reactions 3325 €

NEBNext Ultra II RNA Library Prep with Sample Purification Beads


#E7775S 24 reactions1080 € #E7775L 96 reactions3685 €

See ordering information for NEBNext Ultra II modules.


- Get more of what you need, with the highest library yields
- Generate high quality libraries with limited amounts of RNA:
 - 10 ng–1μg Total RNA (polyA mRNA workflow)
 - 5 ng-1μg (rRNA depletion workflow)
- Minimize bias, with fewer PCR cycles required
- Maximize the flexibility to order reagents for your specific workflow needs
 - Directional (strand-specific, using the "dUTP method") and non-directional workflow options available
 - rRNA Depletion and poly(A) mRNA Isolation reagents are available separately
 - Adaptors and primers for multiplexing, in 12- and 96-index formats, are available separately
- Save time with streamlined workflows, reduced hands-on time, and automation compatibility
- Enjoy the reliability of the gold standard SPRIselect size selection and clean-up beads, supplied in just the amounts you need
- Rely on robust performance, even with low quality RNA

Visit **UltraliRNA.com** to learn more and to view performance data

Do you need increased sensitivity and specificity from your RNA-seq experiments? Do you have ever-decreasing amounts of input RNA? To address these challenges, our next generation of RNA library prep kits have been reformulated at each step, resulting in several fold higher yields of high quality libraries and enabling use of lower input amounts and fewer PCR cycles. The kits have streamlined, automatable workflows and are available for directional (strand-specific, using the "dUTP method") and non-directional library prep, with the option of SPRISelect beads for size-selection and clean-up steps.




NEBNext Ultra II Directional RNA produces the highest yields, from a range of input amounts. Poly(A)-containing mRNA was isolated from 10 ng, 100 ng and 1 µg of Universal Human Reference RNA (Agilent #740000) and libraries were made using the NEBNext Ultra II Directional RNA kit, Kapa Stranded mRNA-Seq kit, Kapa mRNA HyperPrep kit and Illumina TruSeq Stranded mRNA Kit. The input RNA amount and number of PCR cycles are indicated. Library yields from an average of three replicates are shown. Error bars indicate standard deviation. Library yields were assessed using the Agilent® Bioanalyzer®.



Uniformity of Coverage across the DAM1 transcript. Poly(A)-containing mRNA was isolated from 10 ng, 100 ng and 1 µg of Universal Human Reference RNA (Agilent #740000) and libraries were made using the NEBNext Ultra II Directional RNA kit, Kapa Stranded mRNA-Seq kit, Kapa mRNA Hyperpep kit and Illumina TruSeq Stranded mRNA Kit. Coverage across transcript ENST00000369541.3 (DAM1) was assessed by mapping reads directly to the transcriptome (Hisat 2.0.3) and assessing coverage using bedtools cov in 100 bins along the transcript length. Libraries prepared using the NEBNext Ultra II Directional RNA Kit provided superior coverage across the transcript at 100 ng and 10 ng input amounts.

NEBNext Ultra II RNA Reagents for Illumina Sequencing

NEBNext Ultra II RNA Kits are available for directional (strand-specific) and non-directional library preparation, and for bulk RNA and single cell samples. These kits utilize streamlined workflows and have been designed for performance with input amounts as low as 5 ng. Note that reagents for rRNA depletion and poly(A) mRNA enrichment are supplied separately, as are adaptors and primers. In addition to stringent QCs on individual components, the NEBNext RNA kits are functionally validated by preparation of a library, followed by Illumina sequencing.



NEBNext Single Cell/Low Input RNA Library Prep

NEW

NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina

#E6420S 24 reactions1195 € #E6420L 96 reactions3965 €

NFW

NEBNext Single Cell/Low Input cDNA Synthesis & Amplification Module

#E6421S 24 reactions 650 € #E6421L 96 reactions 2210 €

NEW

NEBNext Single Cell Lysis Module


#E5530S 96 reactions90 €

- Generate the highest yields of highquality full-length transcript sequencing libraries from single cells, or as little as 2 pg-200 ng total RNA
- Experience unmatched detection of low abundance transcripts
- Rely on consistent transcript detection for a wide range of input amounts and sample types
- Obtain full-length, uniform transcript coverage, regardless of input amount or sample type
- Use with cultured or primary cells, or total RNA
- Save time with a fast, streamlined workflow, minimal handling steps and hands-on time
 - Single-tube protocol from cell lysis to cDNA
 - Enzymatic DNA fragmentation, end repair and dA-tailing reagents in a single enzyme mix, with a single protocol, regardless of GC content
- Available with or without library construction reagents


Description: The unique workflow of the NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina meets the demand for a highly sensitive, yet robust method that consistently generates high-quality, full-length transcript sequencing data from a single cell or ultra-low input RNA.

Optimized cDNA synthesis and amplification steps incorporate template switching, and utilize a unique protocol and suite of reagents.

cDNAs are generated directly from single cells or 2 pg–200 ng RNA, and even low-abundance transcripts are represented in the high yields of cDNA obtained. This is followed by library construction that incorporates the Ultra II FS enzymatic DNA fragmentation/end repair/dA-tailing mix in a simple and efficient workflow.



Increased transcript detection with the NEBNext Single Cell/Low Input RNA Library Prep Kit Sequencing libraries were generated from Jurkat single cells (6 replicates) using the NEBNext Single Cell/Low Input RNA Library Prep Kit, or the SMART-Seq v4 Ultra Low Input RNA Kit for Sequencing (Clontech # 634891) plus the Nextera XT DNA Library Prep Kit (Illumina #FC-131-1096). Libraries were sequenced on an Illumina NextSeq 500 using paired-end mode (2 x 76 bp). TPM = Transcripts per Kilobase Million. Each dot represents the number of transcripts identified at the given TPM range, and each box represents the median, first and third quartiles per replicate and method. Salmon 0.6 was used for read mapping and quantification of all GENCODE v25 transcripts. Panels show the number of transcripts detected within the following TPM ranges: 1-5, 5-10, 10-50 and > 50 TPM. Increased identification of low abundance transcripts is observed with the NEBNext libraries.



Higher library yields with the NEBNext Single Cell/Low Input RNA Library Prep Kit Sequencing libraries were generated from HeLa, Jurkat and M1 single cells or 10 pg of Universal Human Reference (UHR) RNA (Agilent #740000) with recommended amounts of ERCC RNA Spike-In Mix I (Thermo Fisher Scientific® #4456740). The NEBNext Single Cell/Low Input RNA Library Prep Kit, or the SMART-Seq v4 Ultra Low Input RNA Kit for Sequencing (Clontech #634891) plus the Nextera® XT DNA Library Prep Kit (Illumina #FC-131-1096) were used. Error bars indicate standard deviation for 6-11 replicates. For the NEBNext workflow ~80% of the cDNA was used as input into sequencing library preparation, and libraries were amplified with 8 PCR cycles. For the SMART-Seq v4/Nextera XT workflow, as recommended, 125 pg of cDNA was used as input in sequencing library preparation and 12 PCR cycles were used for amplification. Error bars indicate standard deviation for 6-11 replicates.

NEBNext Small RNA Library Prep Kits

NEBNext Multiplex Small RNA Library Prep Set for Illumina (Set 1)

#E7300S 24 reactions1470 € #E7300L 96 reactions4998 €

NEBNext Multiplex Small RNA Library Prep Set for Illumina (Set 2)

#E7580S 24 reactions1470 € #E7580L 96 reactions4998 €

NEBNext Multiplex Small RNA Library Prep Kit for Illumina (Index Primers 1-48)

#E7560S 96 reactions4998 €


NEBNext Small RNA Library Prep Set for Illumina (Multiplex Compatible)

#E7330S 24 reactions1200 € #E7330L 96 reactions4080 €

For adenylation of custom ssDNA adaptors, the 5´ DNA Adenylation Kit is available (NEB #E2610).

The novel NEBNext Small RNA workflow has been optimized to minimize adaptor-dimers while producing high-yield, high-diversity libraries. Adaptors and primers are included in the Small RNA kits, and multiplexing options are available. The Multiplex kit contains index primers, and the Multiplex-Compatible kit enables use with your own barcode primers.

In addition to stringent QCs on individual components, the NEBNext Small RNA kits are functionally validated by library preparation of a Small RNA library, followed by Illumina sequencing. Reagent lots are also reserved specifically for inclusion in NEBNext kits. Most of these components are available in master mix format, reducing the number of vials provided in the kits, and reducing pipetting steps.





Isabel has been a Technical Support Scientist at NEB for over 3 years. She is a member of the NEB Running Club and the Hostess of the NEB Internal Tech Support Blog.

NEBNext Adaptors & Primers for Illumina

NEW

NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs) #E6440S 96 rxns (96 indices) 535 € #E6440L 384 rxns (96 indices) 1925 € NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 1) #E7600S 96 rxns (8 x 12 indices) 461 € NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 2) #E7780S 96 rxns (8 x 12 indices) 461 € NEBNext Multiplex Oligos for Illumina (Index Primers Set 1) #E7335S 24 rxns (12 indices) 106 € #E7335L 96 rxns (12 indices) 384 €

NEBNext Multiplex Oligos for Illumina (Index Primers Set 2)

#E7500S 24 rxns (12 indices) 106 € #E7500L 96 rxns (12 indices) 384 €

- Increased ligation efficiency
- Minimized adaptor-dimer formation
- Increased library yields
- Increased sample identification specificity (dual barcodes)
- Unique dual index pairs enable detection of barcode hopping
- Large number of barcodes/indices available
- Index pooling guidelines and sample sheets are provided

NEBNext Multiplex Oligos for Illumina (Index Primers Set 3)

#E7710S 24 rxns (12 indices) 106 € #E7710L 96 rxns (12 indices) 384 €

NEBNext Multiplex Oligos for Illumina (Index Primers Set 4)

#E7730S 24 rxns (12 indices) 106 € #E7730L 96 rxns (12 indices) 384 €

NEBNext Multiplex Oligos for Illumina (96 Index Primers)

#E6609S 96 rxns (96 indices) 680 € #E6609L 384 rxns (96 indices)2720 €

NEBNext Multiplex Oligos for Illumina (Methylated Adaptor, Index Primers Set 1)

#E7535S 24 rxns (12 indices) 135 € #E7535L 96 rxns (12 indices) 486 €

NEW


NEBNext Adaptor Dilution Buffer

#B1430S 1 x 9.6 ml32 €

Description: Designed for use in library prep for DNA, ChIP DNA and RNA (but not Small RNA), the NEBNext Adaptors enable high-efficiency adaptor ligation and high library yields, with minimized adaptor-dimer formation. Incorporating a novel hairpin loop structure, the NEBNext Adaptor ligates with increased efficiency to end-repaired, dA-tailed DNA. The loop contains a U, which is removed by treatment with USER® Enzyme (a combination of UDG and Endo VIII), to open up the loop and make it available as a substrate for PCR. During PCR, barcodes can be incorporated by use of the NEBNext index primers, thereby enabling multiplexing. NEBNext Oligos can be used with NEBNext products, and with other standard Illumina-compatible library preparation protocols.

Single or dual barcode primer options are available. Unique dual index primer pairs are available to address the "index hopping" seen with certain Illumina sequencig instruments.

Functional Validation: Each set is functionally validated by construction of libraries, followed by Illumina sequencing.



NEBNext rRNA Depletion Kits (Human/Mouse/Rat)

NEBNext rRNA Depletion Kit (Human/Mouse/Rat)

#E6310S 6 reactions 330 € #E6310L 24 reactions 1200 € #E6310X 96 reactions 4320 €

NEBNext rRNA Depletion Kit (Human/Mouse/Rat) with RNA Sample Purification Beads

#E6350S 6 reactions 374 € #E6350L 24 reactions 1375 € #E6350X 96 reactions 4950 €


- A single kit that performs reliably well for all of your RNA samples:
 - FFPE (degraded) or high-quality (intact) RNA
 - 10 ng to 1 μg input amounts
- Remove > 95% of rRNA, and obtain more relevant sequence reads from your sample
- Obtain a more complete transcriptome picture through retention of noncoding & incomplete RNAs that are lost with oligo d(T) poly(A) mRNA enrichment methods
- Suitable for use with human, mouse or rat samples
- Easily integrated upstream of any downstream random-primed cDNA synthesis protocol
- Enjoy the reliability of the gold standard RNAClean beads, supplied in just the amounts you need

Description: The NEBNext rRNA Depletion Kit (Human/Mouse/Rat) employs an RNase H-based method (1,2) to deplete both cytoplasmic (5S rRNA, 5.8S rRNA, 18S rRNA and 28S rRNA) and mitochondrial ribosomal RNA (12S rRNA and 16S rRNA) from human total RNA preparations. This product is suitable for both intact and degraded RNA (e.g., FFPE RNA). The resulting rRNA-depleted RNA is suitable for RNA-Seq, random-primed cDNA synthesis, or other downstream RNA analysis applications. This kit is now available with or without RNAClean® beads.

- (1) Adiconis, X. et al (2013) Nature Methods, 10, 623-629.
- (2) Morlon, J.D. et al (2012) PLoS One, 77 e42882.

The rRNA Depletion Kit Includes:

- RNase H
- RNase H Reaction Buffer (10X)
- NEBNext rRNA Depletion Solution
- NEBNext Probe Hybridization Buffer
- DNase I (RNase-free)
- DNase I Reaction Buffer
- Nuclease-free Water
- NEBNext RNA Sample Purification Beads (#E6350)



NEBNext rRNA Depletion Kit (Human/Mouse/Rat) Workflow

NEBNext Poly(A) mRNA Magnetic Isolation Module

#E7490S 24 reactions67 € #E7490L 96 reactions242 € **Description:** The NEBNext Poly(A) mRNA Magnetic Isolation Module is designed to isolate intact poly(A)+ RNA from previously-isolated total RNA. The technology is based on the coupling of Oligo $d(T)_{25}$ to 1 μ m paramagnetic beads which are then used as the solid support for the direct binding of poly(A)+ RNA. Thus, the procedure permits the manual processing of multiple samples and can be adapted for automated high-throughput applications. Additionally, magnetic separation technology permits elution of intact mRNA in small volumes eliminating the need for precipitating the

poly(A)+ transcripts in the eluent. Intact poly(A)+ RNA which is fully representative of the mRNA population of the original sample can be obtained in less than one hour.

The mRNA Magnetic Isolation Module Includes:

- NEBNext Oligo d(T)₂₅ Beads
- NEBNext RNA Binding Buffer (2X)
- NEBNext Wash Buffer
- Nuclease-free Water
- NEBNext Tris Buffer

NEBNext Globin & rRNA Depletion Kits

NEW

NEBNext Globin & rRNA Depletion Kit (Human/Mouse/Rat)

#E7750S 6 reactions 340 € #E7750L 24 reactions 1235 € #E7750X 96 reactions 4450 €

NEV

NEBNext Globin & rRNA Depletion Kit (Human/Mouse/Rat) with RNA Sample Purification Beads

#E7755S 6 reactions 350 € #E7755L 24 reactions 1290 € #E7755X 96 reactions 4635 €

- Efficient, specific depletion of globin mRNA and rRNA
- Suitable for low-quality or high-quality RNA
- Compatible with a broad range of input amounts: 10 ng−1 µg
- Optional integration with poly(A) mRNA isolation workflows for removal of globin RNAs, rRNAs, and noncoding RNAs
- Fast workflow: 2 hours, with less than 10 minutes hands-on time
- Available with optional RNAClean beads

Description: The great majority of RNA in blood samples is comprised of globin mRNA as well as cytoplasmic and mitochondrial ribosomal RNAs (rRNA). These highly abundant RNA species can conceal the biological significance of less abundant transcripts, and so their efficient and specific removal is desirable.

The NEBNext Globin & rRNA Depletion Kit (Human/ Mouse/Rat) employs the NEBNext RNase H-based RNA depletion workflow to deplete the following:

- Globin mRNA (HBA1/2, HBB, HBD, HBM, HBG1/2, HBE1, HBQ1 and HBZ)
- Cytoplasmic rRNA (5S, 5.8S, 18S and 28S)
- . Mitochondrial rRNA (12S and 16S)


The kit is effective with human, mouse and rat total RNA preparations, both intact and degraded. The resulting depleted RNA is suitable for RNA-seq, random-primed cDNA synthesis, or other downstream RNA analysis.

This kit can also be used following poly(A) mRNA enrichment (e.g., using the NEBNext poly(A) mRNA Magnetic Isolation Module, NEB #E7490), so that the final depleted RNA contains only mRNA of interest and no non-coding RNA.

This kit is available with or without RNAClean beads.

The Globin & rRNA Depletion Kits Include:

- RNase H
- RNase H Reaction Buffer (10X)
- NEBNext Globin & rRNA Depletion Solution
- NEBNext Probe Hybridization Buffer
- DNase I (RNase-free)
- DNase I Reaction Buffer
- Nuclease-free Water
- NEBNext RNA Sample Purification Beads (#E7785)



NEBNext Globin & rRNA Depletion Kit efficiently removes Globin mRNA and rRNA. Ribosomal RNA (rRNA) and globin mRNA were depleted from Human, Mouse, and Rat Whole Blood Total RNA (100 ng) using the NEBNext Globin & rRNA Depletion Kit (Human/Mouse/Rat) (NEB #E7750) or the Globin-Zero® Gold rRNA Removal Kit (Illumina #GZG1224). Libraries were prepared from the depleted RNA using the NEBNext Ultra II Directional RNA Library Prep Kit for Illumina and sequenced on an Illumina NextSeq® instrument (2 x 75 bp). Reads were down sampled to 20 million reads per sample for analysis, and were identified as globin mRNA (A) or rRNA (B) using Mirabait (6 or more, 25-mers). The data represents an average of 3-4 replicates. Error bars indicate standard error.

NEBNext Magnetic Separation Rack

#S1515S 24 tubes (0.2 ml)475 €

- Fast separations in purification and size-selection steps in next generation sequencing workflows
- Small-scale separation of magnetic particles
- Anodized aluminum rack with Neodymium Iron Boron (NdFeB) rare earth magnets, the most powerful commercially available
- 24 tube capacity:
 8- and 12-strip 0.2 ml PCR tubes or individual 0.2 ml PCR tubes

Description: Next generation sequencing library preparation workflows include magnetic bead-based purification and size-selection steps and it is important for library yield and quality that bead separation be highly efficient and fast.


The NEBNext Magnetic Separation Rack was designed for this application and contains rare earth Neodymium Iron Boron (NdFeB) magnets, the most powerful commercially available magnets, in an anodized aluminium rack. The rack holds 24 0.2 ml tubes, and is compatible with single tubes or strip tubes.



NEBNext Direct* – Target Enrichment for NGS

NEBNext Direct enables highly specific target enrichment of genomic regions of interest. This innovative approach to target enrichment balances the speed and precision of multiplexed PCR-based approaches with the content scalability typical of hybridization-based methods.

This flexibility allows a singular workflow for assays ranging from single gene tests to comprehensive panels including hundreds of genes. Regardless of sample type or assay content, NEBNext Direct allows you to enrich your targets with precision.



NEBNext Direct employs a fast hybridization-based workflow that couples capture with library preparation



NEBNext Direct Cancer HotSpot Panel

#E7000S 8 reactions 1450 € #E7000L 24 reactions 4200 € #E7000X 96 reactions 11180 €


- Generate a higher percentage of your sequencing reads aligning to your targets
- Eliminate the need to over-sequence, reducing cost per sample
- Obtain uniform sequencing of all targets, regardless of GC content
- Save time with a 1-day workflow that combines enrichment with library preparation
- Generate high quality libraries with limited input amounts and degraded DNA samples, including FFPE and ctDNA
- Distinguish molecular duplicates, reducing false positive variants and improving sensitivity

Visit NEBNextDirect.com for more information, including additional performance data

Description: Using a novel approach to target enrichment, the NEBNext Direct Cancer Hotspot Panel enables highly specific hybridization-based capture of 190 common cancer targets from 50 genes. The NEBNext Direct technology offers significant advantages over both traditional in-solution hybridization and multiplex PCR protocols. Target enrichment is combined with library preparation, reducing processing time and minimizing sample loss. Ideal for automation, NEBNext Direct enables deep sequencing of genomic regions of interest for the discovery and identification of low frequency variants from challenging sample types.

TARGETS INCLUDE REGIONS FROM THE FOLLOWING CANCER-RELATED GENES				
ABL1	EGFR	GNAQ	KRAS	PTPN11
AKT1	ERBB2	GNAS	MET	RB1
ALK	ERBB4	HNF1A	MLH1	RET
APC	EZH2	HRAS	MPL	SMAD4
ATM	FBXW7	IDH1	NOTCH1	SMARCB1
BRAF	FGFR1	IDH2	NPM1	SM0
CDH1	FGFR2	JAK2	NRAS	SRC
CDKN2A	FGFR3	JAK3	PDGFRA	STK11
CSF1R	FLT3	KDR	PIK3CA	TP53
CTNNB1	GNA11	KIT	PTEN	VHL

For research use only, not intended for diagnostic use.



The NEBNext Direct Cancer HotSpot Panel demonstrates the ability to accurately detect a range of nucleic acid variants. This figure shows the expected versus observed variant allele frequencies (VAF) across the range of well-characterized variants present in a pool of 24 HapMap samples screened against the NEBNext Direct Cancer HotSpot Panel. 100 ng of input DNA was used, samples were sequenced on the Illumina MiSeq using 2 x 75 bp sequencing, and standard data analysis and variant calling algorithms were used. We were able to successfully detect 100% of the 168 truth variants present across a range of 2–100% VAF. The high degree of linearity across this broad dynamic range demonstrates the ability of the NEBNext Direct Cancer HotSpot Panel to accurately predict variant allele frequencies across a broad dynamic range.

NEBNext Direct BRCA1/BRCA2 Panel


#E6627S 8 reactions660 € #E6627L 24 reactions1900 € #E6627X 96 reactions6950 €

- Generate a higher percentage of your sequencing reads aligning to your targets
- Eliminate the need to over-sequence, reducing cost per sample
- Obtain uniform sequencing of all targets, regardless of GC content
- Save time with a 1-day workflow that combines enrichment with library preparation
- Generate high quality libraries with limited input amounts and degraded DNA samples, including FFPE and ctDNA
- Distinguish molecular duplicates, reducing false positive variants and improving sensitivity

Visit **NEBNextDirect.com** for more information, including additional performance data

Description: NEBNext Direct employs a unique hybridization-based enrichment workflow that hybridizes baits directly to genomic DNA, without the need for upfront library preparation. The BRCA1/BRCA2 panel

demonstrates extremely high specificity and unmatched coverage uniformity across a wide range of DNA inputs, allowing highly sensitive calling of germline and somatic variants while maximizing sequencer efficiency.



The NEBNext Direct BRCA1/BRCA2 Panel delivers highly efficient enrichment of BRCA1 and BRCA2 coding regions with a high percentage of reads mapping to targets. This histogram shows the percent of reads aligned to the human genome, and the percent of reads mapped to the targets included in the BRCA1/BRCA2 Panel across different input DNA amounts. 10 ng, 100 ng and 1 µg of purified genomic DNA was enriched using the NEBNext Direct BRCA1/BRCA2 Panel. Sequencing reads were generated on an Ilumina® MiSeq with 2 x 75 bp reads, 8 bp Sample ID and 12 bp unique molecular identifier. Sequencing read alignments were performed with BWA-MEM, and PCR duplicates were filtered using the UMIs.



NEBNext Direct Custom Ready Panels

#E6631S 8 reactions #E6631L 24 reactions #E6631X 96 reactions

Visit www.neb.com/E6631 to learn more and request a quote.

- Choose from a single gene to hundreds of genes.
- Experience unmatched specificity and coverage uniformity
- Eliminate synthesis and optimization steps, for a faster turnaround
- Improve sensitivity with our Unique Molecular Identifier (UMI)
- Generate results in one day with our automation-friendly workflow


Panels can be designed and ordered by visiting www.neb.com/ CustomReadyPanelForm

Description: Employing the unique NEBNext Direct hybridization-based enrichment method, NEBNext Direct Custom Ready Panels allow rapid customization of targeted gene panels for Illumina sequencing. Select from a list of human genes where baits have been carefully designed and optimized to provide complete coverage of the full coding (exon) regions. High quality panels can be designed by you and rapidly delivered, from any combination of genes.


The genes available through the NEBNext Direct Custom Ready offering are continually updated, and currently include those associated with a variety of translational research areas, including cancer, neurological disorders, cardiological disease, autism, severe combined immunodeficiency, cystic fibrosis and the recommended genes for incidental findings by the American College of Medical Genetics. The full list of genes currently available can be found at www.neb.com/CustomReadyPanelForm.

NEBNext Direct Custom Ready Panels can include anywhere from a single specific gene up to 1.5 megabases of total target territory. There are no limitations on genes that can be combined together in a panel. Each panel is tested prior to shipment, and sequencing results are returned through a custom Performance Report.


Bait sets for each gene included in the panel have undergone a rigorous development and optimization process to maximize specificity and target coverage uniformity.



NEBNext Direct Custom Ready Panels demonstrate optimum performance across a wide range of panel sizes. Key target enrichment metrics demonstrate consistent performance across a range of panel sizes. 100 ng of DNA was tested against panels of 1, 10, 25, 50 and 100 genes, and sequenced using Illumina paired-end 150 bp sequencing. Larger panels included all genes present in smaller panels.



Sensitivity in detection of variants across panel size and DNA input amount. 24 HapMap samples were blended to create a range of variant allele frequencies (VAF) down to 2%. 25, 50, 100, 200 and 500 ng of this blended DNA was enriched using NEBNext Direct Custom Ready Panels of 1, 10, 25, 50 and 100 genes. Larger panels were inclusive of the genes in smaller panels. Resulting libraries were sequenced using 2 x 150 bp Illumina sequencing and variants were called using Mutect and Vardict variant calling algorithms.



NEBNext Direct Custom Ready Panels demonstrate retention of target behavior across panel sizes. *IGV image of coverage profile for 4 BRAF exons included in panels of 1, 10, 25, 50 and 100 genes, demonstrates consistent target behavior with the addition of gene targets. 100 ng of DNA was used as input for NEBNext enrichment using the 5 panels, including the BRAF gene. Libraries were sequenced on an Illumina 2 x 150 basepair sequencing.*

NEBNext Microbiome DNA Enrichment Kit

#E2612S 6 reactions 210 € #E2612L 24 reactions 750 €


- Effective separation of microbial DNA from contaminating host DNA
- Fast, simple protocol
- Compatible with downstream applications including next generation sequencing on all platforms, qPCR and end point PCR
- Suitable for a wide range of sample types
- No requirement for live cells
- Captured host DNA can also be eluted and retained

Description: The NEBNext Microbiome DNA Enrichment Kit facilitates separation of microbial DNA from methylated host DNA (including human) by selective binding and removal of the CpG methylated host DNA (1).


Functional Validation: Each set of reagents are functionally validated by enriching *E. coli* DNA from a mixture of *E. coli* and human DNA. Enrichment is evaluated through library construction and sequencing of the enriched sample on an Illumina sequencer.


The Microbiome Enrichment Kit Includes:

- NEBNext MBD2-Fc Protein
- NEBNext Protein A Magnetic Beads
- NEBNext Bind/Wash Buffer (5X)
- 16S RNA Universal Bacteria Control Primers
- RPL30 Human DNA Control Primers
- (1) Feehery, G.R. et al. (2013) PLoS One, 8: e76096.
- (2) Chen, T., et al. (2010) *Database*, Vol. 2010, Article ID bag013, doi: 10.1093/database/bag013
- (3) Langmead, B., et al. (2009) *Genome Biol.* 10:R25 doi:10.1186/qb-2009-10-3-r25



Salivary Microbiome DNA Enrichment. DNA was purified from pooled human saliva DNA (Innovative Research) and enriched using the NEBNext Microbiome DNA Enrichment Kit. Libraries were prepared from unenriched and enriched samples and sequenced on the SOLiD 4 platform. The graph shows percentages of 500 M–537 M SOLiD™4 50 bp reads that mapped to either the Human reference sequence (hg19) or to a microbe listed in Human Oral Microbiome Database (HOMD)[2]. (Because the HOMD collection is not comprehensive, ~80% of reads in the enriched samples do not map to either database.) Reads were mapped using Bowtie 0.12.7[3] with typical settings (2 mismatches in a 28 bp seed region, etc.). SOLID™4 is a registered trademark of Life Technologies, Inc.





NEBNext Microbiome DNA Enrichment Kit workflow

Microbiome Diversity is Retained after Enrichment with the NEBNext Microbiome DNA Enrichment Kit. DNA was purified from pooled human saliva DNA (Innovative Research) and enriched using the NEBNext Microbiome DNA Enrichment Kit. Libraries were prepared from unenriched and enriched samples, followed by sequencing on the SOLiD 4 platform. The graph shows a comparison between relative abundance of each bacterial species listed in HOMD[2] before and after enrichment with the NEBNext Microbiome DNA Enrichment Kit. Abundance is inferred from the number of reads mapping to each species as a percentage of all reads mapping to HOMD. High concordance continues even to very low abundance species (inset). We compared 501M 50 bp SOliD 4 reads in the enriched dataset to 537M 50 bp SOLiD 4 reads in the unenriched dataset. Reads were mapped using Bowtie 0.12.7[3] with typical settings (2 mismatches in a 28 bp seed region, etc).

* Niesseria flavescens — This organism may have unusual methylation density, allowing it to bind the enriching beads at a low level. Other Niesseria species (N. mucosa, N. sicca and N. elognata) are represented, but do not exhibit this anomalous enrichment.

NEBNext FFPE DNA Repair Mix

#M6630S 24 reactions173 € #M6630L 96 reactions604 €

- Construct high-quality NGS libraries from FFPE DNA samples
- Use upstream of library prep for any NGS platform
- No alteration of DNA sequence
- Rely on NEB's NGS validation process for FFPE DNA library prep

Description: Archiving of clinical materials as Formalin-Fixed, Paraffin-Embedded (FFPE) samples is a common practice. However, the methods used for fixation and storage significantly damage and compromise the quality of nucleic acids from these samples. As a result, it can be challenging to obtain useful information, including high quality sequence data, especially when sample amounts are limited. The NEBNext FFPE DNA Repair Mix is a cocktail of enzymes formulated to repair DNA, and specifically optimized and validated for repair of FFPE DNA samples. The FFPE DNA Repair Mix increases library yield and overall library success rates, without introduction of bias.

Functional Validation: Each lot is functionally validated by repair of FFPE DNA followed by library construction and Illumina sequencing.

For the most up-to-date product and pricing information, visit **NEBNext.com**

FFPE DAMAGE TYPE	REPAIRED BY THE FFPE DNA ENZYME REPAIR MIX?
Deamination of cytosine to uracil	Yes
Nicks and gaps	Yes
Oxidized bases	Yes
Blocked 3' ends	Yes
DNA fragmentation	No
DNA-protein crosslinks	No

Table 1: Types of FFPE DNA damage and ability to be repaired by the NEBNext FFPE DNA Repair Mix.




Figure 1: Effect of FFPE DNA Repair Mix on library yields. An example of Agilent Bioanalyzer traces of libraries prepared from stomach tumor FFPE DNA that was treated with the FFPE DNA Repair Mix, or was untreated, before library construction. Yield improvements of 101% to 458% have been observed.

NEBNext dsDNA Fragmentase®

#M0348S 50 reactions 98 € #M0348L 250 reactions 392 €

Companion Product:

- Generation of dsDNA fragments for sequencing on next generation sequencing platforms
- Generation of dsDNA fragments for libraries

37° ₩

Description: NEBNext dsDNA Fragmentase generates dsDNA breaks in a time-dependent manner to yield 50-1,000 bp DNA fragments depending on reaction time. NEBNext dsDNA Fragmentase contains two enzymes, one randomly generates nicks on dsDNA and the other recognizes the nicked site and cuts the opposite DNA strand across from the nick, producing dsDNA breaks. The resulting DNA fragments contain short overhangs, 5´-phosphates, and 3´-hydroxyl groups. The random nicking activity of NEBNext dsDNA Fragmentase has been confirmed by preparing libraries for nextgeneration sequencing. A comparison of the sequencing results between gDNA prepared with NEBNext dsDNA fragmentase and with mechanical shearing demonstrates that the NEBNext dsDNA Fragmentase does not introduce any detectable bias during the sequencing library

Source: NEBNext dsDNA Fragmentase is composed of endonucleases isolated from two different *E. coli* sources: one construct expresses a fusion protein consisting of *E. coli* maltose binding protein and *Vibrio vulnificus* nuclease mutant protein; the other expresses a fusion protein consisting of maltose binding protein and T7 endonuclease mutant protein.

preparation and no difference in sequence coverage is

observed using the two methods.

Reaction Conditions: 1X NEBNext dsDNA Fragmentase Reaction Buffer v2, supplemented with 100 μM MgCl₂, when required. Incubate at 37°C.

1X NEBNext dsDNA Fragmentase Reaction Buffer v2:

20 mM Tris-HCl 10 mM MgCl_2 50 mM NaCl 0.15% Triton X-100 pH 7.5 @ $25^{\circ}\mathrm{C}$

Reagents Supplied with Enzyme:

10X NEBNext dsDNA Fragmentase Reaction Buffer v2 200 mM MgCl_a

Heat Inactivation: 65°C for 15 minutes in the presence of 50 mM DTT

MINELUTE® is a registered trademark of the Qiagen Group




Figure 1: Fragmentation of E. coli gDNA. E. coli gDNA was fragmented with NEBNext dsDNA Fragmentase for the indicated times and purified on MinElute® columns.



Chaithanya has been with NEB for 5 years, and is currently a Development Scientist. His work focuses on the development of new tools for NGS sample preparation. When he isn't busy at the bench, Chaithanya enjoys photography and a good cricket match.

NEBNext Ultra II Q5® Master Mix

Additional Products:

NEBNext Q5 I	HotStart HiFi PCR Master Mix
#M0543S	50 reactions 90 €
#M0543L	250 reactions 360 €
NEBNext High	n-Fidelity 2X PCR Master Mix
NEBNext High #M0541S	n-Fidelity 2X PCR Master Mix 50 reactions (50 µI) 90 €
0	,

- Next generation sequencing library preparation
- High-fidelity amplification
- Uniform GC coverage
- Improves sequencing library coverage of known difficult regions of the human genome


Description: The NEBNext Ultra II Q5 Master Mix is a new formulation of Q5 DNA Polymerase that has been optimized for robust, high-fidelity amplification of next-generation sequencing (NGS) libraries. This formulation further improves the uniformity of amplification of libraries, including superior performance with GC-rich regions.

The polymerase component of the master mix, Q5 High-Fidelity DNA Polymerase, is a novel thermostable DNA polymerase that possesses 3′→5′ exonuclease activity, and is fused to a processivity-enhancing Sso7d domain. Q5 also has the highest fidelity available (> 100-fold higher than that of *Taq* DNA Polymerase and ~12-fold higher than that of *Pyrococcus furiosus* (Pfu) DNA Polymerase), resulting in ultra-low error rates.

The NEBNext Ultra II Q5 Master Mix is an aptamerbased hot start formulation that allows convenient room temperature reaction set up. The convenient 2X master mix format contains dNTPs, Mg** and a proprietary buffer, and requires only the addition of primers and DNA template for robust amplification. NEBNext Ultra II Q5 Master Mix is also included in the NEBNext Ultra II DNA Library Prep Kit for Illumina.

Source: An *E. coli* strain that carries the Q5 High-Fidelity DNA Polymerase gene.

Reaction Conditions: NEBNext Ultra II Q5 Master Mix, DNA template and 1 μ M primers in a total reaction volume of 50 μ l.




NEBNext Ultra II Q5 Master Mix provides improved coverage of known low-coverage regions of the human genome. Libraries were prepared from Human NA19240 genomic DNA. One library was not amplified. The other two libraries were amplified using 5 cycles of PCR with NEBNext Q5 Hot Start HiFi PCR Master Mix (NEB #M0543) or with NEBNext Ultra II Q5 Master Mix (NEB #M0544). Libraries were sequenced on an Illumina NextSeq 500. 420 million 75 bp reads were randomly extracted from each dataset, representing an average coverage of 10X. Reads were mapped to the GRCh37 reference genome using Bowtie 2.2.4. Reads on each region were counted using bedtools v2.19.1.

A: The number of reads overlapping distinct low coverage regions of the human genome (1) are shown for each library.

B: From the 420 million 75 bp reads randomly extracted from each dataset, 10X coverage was expected. The % of difficult regions covered at > 10X are shown for each library. The NEBNext Ultra II Q5 Master Mix provides improved coverage of these known low coverage regions, without drop-outs, and shows similar coverage to the unamplified sample.

(1) Popatov, V. and Ong, J.L. (2017). Examining Sources of Error in PCR by Single-Molecule Sequencing. PLoS ONE. 12(1):e0169774.



NEBNext Ultra II Q5 Master Mix provides uniform GC coverage for microbial genomic DNA with a broad range of GC composition. Libraries were made using 100 ng of the genomic DNAs shown and the NEBNext Ultra II DNA Library Prep Kit. Libraries were amplified using the NEBNext Ultra II Q5 Master Mix, and sequenced on an Illumina MiSeq. GC coverage information was calculated using Picard's CollectGCBiasMetrics (v1.117). Expected normalized coverage of 1.0 is indicated by the horizontal grey line, the number of 100 bp regions at each GC% is indicated by the vertical grey bars, and the colored lines represent the normalized coverage for each library. NEBNext Ultra II Q5 Master Mix provides uniform GC coverage regardless of the GC content of the DNA.

NEBNext Library Quant Kit for Illumina

NEBNext Library Quant Kit for Illumina #E7630S 100 reactions105 € #E7630L 500 reactions450 €

Companion Product:

NEBNext Library Dilution Buffer #B6118S 7.5 ml $\dots 33 \in$

- Provides more accurate and reproducible quant values than alternative methods and kits
- Compatible with libraries with a broad range of insert sizes and GC content, made by a variety of methods
- Requires only 4 standards, allowing more libraries to be quantitated per kit
- Supplied with a convenient, Library Dilution Buffer
- The NEBNext Library Quant Master Mix requires only the addition of primers
- Utilizes a single extension time for all libraries, regardless of insert size
- Library quant values can be easily calculated using NEB's online tool, at NEBioCalculator.neb.com
- ROX is included in the kit, for use with qPCR instruments that require a reference dye for normalization

With NEBNext, optimal cluster density is achieved from quantitated libraries with a broad range of library size and GC content. Libraries of 310–963 bp from the indicated sources were quantitated using the NEBNext Library Quant Kit, then diluted to 8 pM and loaded onto a MiSeq (v2 chemistry; MCS v2.4.1.3). Library concentrations ranged from 7–120 nM, and resulting raw cluster density for all libraries was 965–1300 k/mm² (ave. =1199). Optimal cluster density was achieved using concentrations determined by the NEBNext Library Quant Kit for all library sizes.

Three 340–400 bp libraries were quantitated by 4 different users 2–4 times using either the NEBNext or Kapa Library Quantification Kit (Universal). A notable improvement in quantitation consistency was observed for concentrations determined by the NEBNext Kit (orange) versus those from the Kapa kit (gray).


Description: Accurate quantitation of nextgeneration sequencing libraries is essential for maximizing data output and quality from each sequencing run. For Illumina sequencing specifically, accurate quantitation of libraries is critical to achieve optimal cluster densities, a requirement for optimal sequence performance, gPCR is considered to be the most accurate and effective method of library quantitation, providing considerably higher consistency and reproducibility than electrophoresis or spectrophotometry, which measure total nucleic acid concentration. Amplification-based methods quantitate only those molecules that contain both adaptor sequences, thereby providing a more accurate estimate of the concentration of the library molecules that can be sequenced.


The NEBNext Library Quant Kit delivers significant improvements to qPCR-based library quantitation for next gen sequencing. The NEBNext Library Quant Kit


for Illumina contains components that are optimized for qPCR-based quantitation of libraries prepared for Illumina next-generation sequencing platforms. The NEBNext Library Quant Kit contains primers which target the P5 and P7 Illumina adaptor sequences and a set of high-quality, pre-diluted DNA standards to enable reliable quantitation of diluted DNA libraries between 150–1000 bp.

The Library Quant Kit Includes:

- NEBNext Library Quant Master Mix
- NEBNext Library Quant Primer Mix
- NEBNext Library Quant DNA Standards 1-4
- ROX (Low) and ROX (High)
- NEBNext Library Dilution Buffer









NEBNext Reagents for Ion Torrent™: DNA Library Preparation

NEBNext kits are available for DNA library preparation for lon Torrent, with or without enzymatic DNA fragmentation. In addition to stringent QCs on individual components, the NEBNext DNA kits are functionally validated by library preparation of a genomic DNA library, followed by lon Torrent sequencing. Reagent lots are reserved specifically for inclusion in NEBNext kits. Most of these reagents are provided in master mix format, reducing the number of vials provided in the kits, and reducing pipetting steps. Adaptors and primers for singleplex libraries are supplied in the kits. For multiplexed libraries, the lon XPress™ Barcode Adaptors from Thermo Fisher Scientific can be used.



^{*}Note that a minimum of 100 ng is recommended when used in conjunction with Ion Express Barcode Adaptors.



Varying GC Content Libraries. 0.5 µg of DNA from 3 different genomes with varying GC content were used to construct 200 bp and 400 bp libraries using the NEBNext Fast DNA Fragmentation and Library Prep Set for Ion Torrent, analyzed by the Agilent Bioanalyzer.

ION XPRESS" is a trademark of Life Technologies, Inc.

NEBNext Reagents for DNA Library Preparation – Ordering Information

Illumina Platform:

KITS FOR ILLUI	NINA DNA LIBRARY PREPARATION	NEB#	SIZE	PRICE
	NEBNext Ultra II DNA Library Prep Kit for Illumina	E7645S/L	24/96 rxns	535 €/2045 €
	NEBNext Ultra II DNA Library Prep with Sample Purification Beads	E7103S/L	24/96 rxns	615 €/2352 €
	NEBNext Ultra II FS DNA Library Prep Kit for Illumina	E7805S/L	24/96 rxns	665 €/2526 €
DNA & ChIP	NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads	E6177S/L	24/96 rxns	730 €/2758 €
	NEBNext Ultra DNA Library Prep Kit for Illumina	E7370S/L	24/96 rxns	488 €/1854 €
	NEBNext DNA Library Prep Master Mix Set for Illumina	E6040S/L	12/60 rxns	356 €/1424 €
	NEBNext ChIP-Seq Library Prep Master Mix Set for Illumina	E6240S/L	12/60 rxns	290 €/1160 €
	NEBNext Enzymatic Methyl-seq Kit	E7120S/L	24/96 rxns	895 €/3360 €
	NEBNext Enzymatic Methyl-seq Conversion Module	E7125S/L	24/96 rxns	185 €/675 €
MODULES & EN	IZYMES	NEB #	SIZE	PRICE
	NEBNext FFPE DNA Repair Mix	M6630S/L	24/96 rxns	173 €/604 €
	NEBNext Microbiome DNA Enrichment Kit	E2612S/L	6/24 rxns	210 €/750 €
	NEBNext Ultra II FS DNA Module	E7810S/L	24/96 rxns	290 €/1036 €
	NEBNext Ultra II End Repair/dA-Tailing Module	E7546S/L	24/96 rxns	262 €/835 €
	NEBNext Ultra II Ligation Module	E7595S/L	24/96 rxns	395 €/1270 €
	NEBNext Ultra II Q5 Master Mix	M0544S/L	50/250 rxns	99 €/395 €
	NEBNext Ultra End Repair/dA-Tailing Module	E7442S/L	24/96 rxns	228 €/731 €
DNA & ChIP	NEBNext Ultra Ligation Module	E7445S/L	24/96 rxns	365 €/11172 €
J	NEBNext dsDNA Fragmentase	M0348S/L	50/250 rxns	98 €/392 €
	NEBNext End Repair Module	E6050S/L	20/100 rxns	86 €/347 €
	NEBNext dA-Tailing Module	E6053S/L	20/100 rxns	102 €/408 €
	NEBNext Quick Ligation Module	E6056S/L	20/100 rxns	306 €/1224 €
	NEBNext Q5 Hot Start HiFi PCR Master Mix	M0543S/L	50/250 rxns	90 €/360 €
	NEBNext High-Fidelity 2X PCR Master Mix	M0541S/L	50/250 rxns	90 €/360 €
	NEBNext Q5U Master Mix	M0597S/L	50/250 rxns	125 €/495 €
	NEBNext dsDNA Fragmentase Reaction Buffer v2	B0349S	6 ml	29 €
ADAPTORS & P	RIMERS	NEB#	SIZE	PRICE
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 1)	E7335S/L	24/96 rxns	106 €/384 €
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 2)	E7500S/L	24/96 rxns	106 €/384 €
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 3)	E7710S/L	24/96 rxns	106 €/384 €
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 4)	E7730S/L	24/96 rxns	106 €/384 €
	NEBNext Multiplex Oligos for Illumina (96 Index Primers)	E6609S/L	96/384 rxns	680 €/2720 €
	NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs)	E6440S/L	96/384 rxns	535 €/1925 €
	NEBNext Multiplex Oligos for Illumina (Methylated Adaptor, Index Primers Set 1)	E7535S/L	24/96 rxns	135 €/486 €
	NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 1)	E7600S	96 rxns	461 €
	NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 2)	E7780S	96 rxns	461 €
	NEBNext Multiplex Oligos for Enzymatic Methyl-seq (Unique Dual Index Primer Pairs)	E7140S/L	24/96 rxns	135 €/535 €
	NEBNext Adaptor Dilution Buffer	B1430S	1 x 9.6 ml	32 €
ARGET ENRIC		NEB #	SIZE	PRICE
	NEBNext Direct Cancer HotSpot Panel	E7000S/L/X	8/24/96 rxns	1450 €/4200 €/11180 €
	NEBNext Direct BRCA1/BRCA2 Panel	E6627S/L/X	8/24/96 rxns	660 €/1900 €/6950 €
	NEBNext Direct Custom Ready Panels	E6631S/L/X	8/24/96 rxns	request quote
IBRARY QUAN		NEB #	SIZE	PRICE
	NEBNext Library Quant Kit for Illumina	E7630S/L	100/500 rxns	105 €/450 €
	NEBNext Library Dilution Buffer	B6118S	7.5 ml	33 €

Ion Torrent Platform:

PRODUCTS FOR D	ONA LIBRARY PREPARATION	NEB #	SIZE	PRICE
DNA	NEBNext Fast DNA Library Prep Set for Ion Torrent	E6270S/L	50 rxns	808 €
	NEBNext Fast DNA Fragmentation & Library Prep Set for Ion Torrent	E6285S/L	50 rxns	960 €

Suitable for Any Sequencing Platform:

DNA ENRIC	HMENT	NEB #	SIZE	PRICE
DNA	NEBNext Microbiome DNA Enrichment Kit	E2612S/L	6/24 rxns	210 €/750 €
DNA REPAII	R	NEB #	SIZE	PRICE
DNA	NEBNext FFPE DNA Repair Mix	M6630S/L	24/96 rxns	173 €/604 €
MODULES 8	& ENZYMES	NEB #	SIZE	PRICE
DNA	NEBNext dsDNA Fragmentase	M0348S/L	50/250 rxns	98 €/392 €
	NEBNext Ultra II Q5 Master Mix	M0544S/L	50/250 rxns	99 €/395 €
	NEBNext Q5 Hot Start HiFi PCR Master Mix	M0543S/L	50/250 rxns	90 €/360 €
	NEBNext High-Fidelity 2X PCR Master Mix	M0541S/L	50/250 rxns	90 €/360 €
	NEBNext dsDNA Fragmentase Reaction Buffer v2	B0349S	6 ml	29 €
MAGNETIC	SEPARATION	NEB #	SIZE	PRICE
	NEBNext Magnetic Separation Rack	S1515S	24 tubes	475 €



Bo is a Research Scientist and has been with NEB for more than 3 years. When she is not busy at the bench, Bo is an avid swimmer.

NEBNext Reagents for RNA Library Preparation – Ordering Information

Illumina Platform:

KITS FOR ILLUM	INA RNA LIBRARY PREPARATION	NEB #	SIZE	PRICE
Directional	NEBNext Ultra II Directional RNA Library Prep Kit for Illumina	E7760S/L	24/96 rxns	1028 €/3495
INA	NEBNext Ultra II Directional RNA Library Prep with Sample Purification Beads	E7765S/L	24/96 rxns	1140 €/3880 :
	NEBNext Ultra Directional RNA Library Prep Kit for Illumina	E7420S/L	24/96 rxns	1172 €/3801
Man dinastrant	NEBNext Ultra II RNA Library Prep Kit for Illumina	E7770S/L	24/96 rxns	980 €/3325
lon-directional NA	NEBNext Ultra II RNA Library Prep with Sample Purification Beads	E7775S/L	24/96 rxns	1080 €/3685
	NEBNext Ultra RNA Library Prep Kit for Illumina	E7530S/L	24/96 rxns	1071 €/3425
	NEBNext Multiplex Small RNA Library Prep Set for Illumina (Set 1)	E7300S/L	24/96 rxns	1470 €/4998
mall RNA	NEBNext Multiplex Small RNA Library Prep Set for Illumina (Set 2)	E7580S/L	24/96 rxns	1470 €/4998
IIIaii iiiva	NEBNext Multiplex Small RNA Library Prep Kit for Illumina (Index Primers 1-48)	E7560S	96 rxns	4998
	NEBNext Small RNA Library Prep Set for Illumina (Multiplex Compatible)	E7330S/L	24/96 rxns	1200 €/4080
ingle Cell	NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina	E6420S/L	24/96 rxns	1195 €/3965
IODULES & ENZ	ZYMES	NEB#	SIZE	PRICE
	NEBNext rRNA Depletion Kit (Human/Mouse/Rat)	E6310S/L/X	6/24/96 rxns	330 €/1200 €/4320
	NEBNext rRNA Depletion Kit (Human/Mouse/Rat) with RNA Sample Purification Beads	E6350S/L/X	6/24/96 rxns	374 €/1375 €/4950
	NEBNext Poly(A) mRNA Magnetic Isolation Module	E7490S/L	24/96 rxns	67 €/242
	NEBNext Magnesium RNA Fragmentation Module	E6150S	200 rxns	45
	NEBNext Ultra II RNA First Strand Synthesis Module	E7771S/L	24/96 rxns	165 €/528
	NEBNext Ultra II Directional RNA Second Strand Synthesis Module	E7550S/L	24/96 rxns	390 €/1245
NA	NEBNext Ultra II Non-directional RNA Second Strand Synthesis Module	E6111S/L	20/100 rxns	295 €/1180
	NEBNext RNA First Strand Synthesis Module	E7525S/L	24/96 rxns	152 €/487
	NEBNext Single Cell/Low Input cDNA Synthesis and Amplification Module	E6421S/L	24/96 rxns	650 €/2210
	NEBNext Single Cell Lysis Module	E5530S	96 rxns	90
	NEBNext Globin & rRNA Depletion Kit	E7750S/L/X	6/24/96 rxns	340 €/1235 €/4450
	NEBNext Globin & rRNA Depletion Kit (Human/Mouse/Rat) with RNA Sample Purification Beads	E7755S/L/X	6/24/96 rxns	350 €/1290 €/4635
	NEBNext Ultra End Repair/dA-Tailing Module	E7442S/L	24/96 rxns	228 €/731
	NEBNext Ultra Ligation Module	E7445S/L	24/96 rxns	365 €/1172
	NEBNext End Repair Module	E6050S/L	20/100 rxns	\$86 €/347
	NEBNext dA-Tailing Module	E6053S/L	20/100 rxns	102 €/408
NA	NEBNext Quick Ligation Module	E6056S/L	20/100 rxns	306 €/1224
	NEBNext Ultra II Q5 Master Mix	M0544S/L	50/250 rxns	99 €/395
	NEBNext Q5 Hot Start HiFi PCR Master Mix	M0543S/L	50/250 rxns	90 €/360
	NEBNext High-Fidelity 2X PCR Master Mix	M0541S/L	50/250 rxns	90 €/360
DAPTORS & PR		NEB#	SIZE	PRICE
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 1)	E7335S/L	24/96 rxns	106 €/384
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 2)	E7500S/L	24/96 rxns	106 €/384
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 3)	E7710S/L	24/96 rxns	106 €/384
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 4)	E7730S/L	24/96 rxns	106 €/384
	NEBNext Multiplex Oligos for Illumina (96 Index Primers)	E6609S/L	96/384 rxns	680 €/2720
	NEBNext Multiplex Oligos for Illumina (Oual Index Primers Set 1)	E7600S	96 rxns	461
	NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 1)	E7780S	96 rxns	461
	NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs)	E6440S/L	96/384 rxns	535 €/1925
	NEBNext Adaptor Dilution Buffer	B1430S	1 x 9.6 ml	32
IBRARY QUANT		NEB #	SIZE	PRICE
IBNAITT QUANT	NEBNext Library Quant Kit for Illumina	E7630S/L	100/500 rxns	105 €/450


Suitable for Any Sequencing Platform:


MODULES & ENZ	TYMES	NEB #	SIZE	PRICE
	NEBNext rRNA Depletion Kit (Human/Mouse/Rat)	E6310S/L/X	6/24/96 rxns	330 €/1200 €/4320 €
	NEBNext rRNA Depletion Kit (Human/Mouse/Rat) with RNA Sample Purification Beads	E6350S/L/X	6/24/96 rxns	374 €/1375 €/4950 €
	NEBNext Poly(A) mRNA Magnetic Isolation Module	E7490S/L	24/96 rxns	67 €/242 €
RNA	NEBNext Magnesium RNA Fragmentation Module	E6150S	200 rxns	45 €
IIIVA	NEBNext Ultra II RNA First Strand Synthesis Module	E7771S/L	24/96 rxns	165 €/528 €
	NEBNext Ultra II Non-directional RNA Second Strand Synthesis Module	E6111S/L	20/100 rxns	295 €/1180 €
	NEBNext Globin & rRNA Depletion Kit	E7750S/L/X	6/24/96 rxns	340 €/1235 €/4450 €
	NEBNext Globin & rRNA Depletion Kit (Human/Mouse/Rat) with RNA Sample Purification Beads	E7755S/L/X	6/24/96 rxns	350 €/1290 €/4635 €
MAGNETIC SEPA	RATION	NEB #	SIZE	PRICE
	NEBNext Magnetic Separation Rack	S1515S	24 tubes	475 €

Featured Online Tools

NEBNext Selector

Use this tool to guide you through selection of NEBNext reagents for next generation sequencing sample preparation. Try it out at **NEBNextSelector.neb.com**









In 2006, beekeepers were reporting huge losses in their honeybee colonies — 30-90% colony loss, up from 15-20% in previous years. This phenomenon was termed Colony Collapse Disorder (CCD), because worker bees would leave their colony, never to return, ultimately causing the entire colony to collapse.

It has since been generally agreed upon that many different environmental factors can cause CCD. Farming practices no longer include cover crops of clover and alfalfa — which are highly nutritious for bees — and the land is now used to grow large monocrops. Infection with Varroa mite leaves the bees immunocompromised and more susceptible to infection with other diseases. Overuse of neurotoxic pesticides (specifically, neonicotinoids) causes death, or in lower doses, disorientation that affects the bees' ability to fly back to their hive.

Furthermore, honeybees are often treated like livestock and are overworked. They are shipped from farm to farm to pollinate vast monocrops, without time to rejuvenate before they are moved on to the next crop that requires pollination. Additionally, like any managed livestock, there are high levels of disease among honeybee colonies, because of the density in which they are kept, and pathogen spillover to other species of wild pollinators is common.

The observation in the past couple of years is that these losses are beginning to plateau. However, this could be due to more controlled breeding of honeybees to meet agricultural demands rather than addressing the environmental triggers. Did we avoid a catastrophic extinction, or is there a broader threat that remains for the health of our pollinators?

There are 20,000 species of bees in the world responsible for pollinating one third of the world's crops, which equates to two to six billion dollars in global agriculture every year. Not only do our fruits and vegetables rely on pollination to eventually make it to our table, but bees also pollinate alfalfa hay, which feeds farm animals.

Pollination of 50% of crops worldwide requires the services of not just honeybees, but a whole host of wild pollinators, including butterflies, wasps, beetles and many other bee species. Truly efficient pollination requires all of these pollinators working together — at different times of the day and on crops best suited to their abilities. For example, tomato plants are either wind pollinated or require particular strong-winged bees, such as bumblebees, to perform "buzz pollination", or sonication, whereby the bee vibrates at a high frequency to release the pollen. Loss of these pollinators means that in some farming areas, tomato plants need to be hand pollinated.

While there has been much discussion about protecting honeybees, we should in fact be thinking about the well-being of ALL the crop pollinators. Being mindful about farming practices, pesticide use and habitat disruption can help to ensure the protection and proliferation of these important workers.

Markers & Ladders (DNA, RNA & Protein)



A wide range of ladders and markers to meet your macromolecule quantification needs.

New England Biolabs provides a wide range of ladders and markers, featuring exceptional quality, uniform band intensities, convenient band spacing and easy-to-identify reference bands.

Double-stranded DNA markers are available for conventional electrophoresis. Conventional electrophoresis markers (size range: ~10 to 2.3 x 10^4 bp) include: HindIII and BstEII digests of Lambda DNA, BstNI and MspI digests of pBR322 DNA, and a HaeIII digest of ϕ X174 DNA.

We also supply a series of DNA ladders ranging from 25 bp to 48.5 kb. Our 100 bp DNA ladder, 1 kb DNA Ladder and 1 kb Plus DNA Ladder are available in several formats: Conventional, Quick-Load® using either non-fluorescing purple dye or bromophenol blue as a tracking dye, and TriDye™ containing three dyes to track gel migration. Our Quick-Load Purple DNA ladders utilize our purple loading dye that improves band visibility by casting no UV shadow and are supplied with a vial of purple dye.

Our RNA ladders and markers have a size range of 17 to 9,000 bases. Both ssRNA ladders are supplied with 2X sample buffer and feature a higher intensity fragment to serve as a reference band. The dsRNA ladders are suitable for use as a size size standard in dsRNA and RNAi analysis on both denaturing polyacrylamide and agarose gels.

For your protein analysis needs, NEB offers a selection of highly pure protein standards. Sizes range from 10–250 kDa, which is ideal for accurate molecular weight determination for a wide range of expressed proteins.

Featured Products

168 Quick-Load Purple
1 kb Plus DNA Ladder

171 Color Prestained Protein Standard, Broad Range (10–250 kDa)

Featured Tools & Resources



Visit www.neb.com/DNAladders to find selection charts for NEB's DNA markers and ladders.




DNA MARKERS & LADDERS

DNA Ladders Size Ranges of DNA Ladders 166 Gel Loading Dye, Purple (6X) 166 Gel Loading Dye, Purple (6X), No SDS 166 1 kb DNA Ladder 167 100 bp DNA Ladder 167 1 kb Plus DNA Ladder 167 50 bp DNA Ladder 167 Fast DNA Ladder Low Molecular Weight DNA Ladder 167 PCR Marker Quick-Load Purple 1 kb Plus DNA Ladder 168 Quick-Load Purple 1 kb DNA Ladder 168 Quick-Load Purple 100 bp DNA Ladder 168 Quick-Load Purple 50 bp DNA Ladder 168 Quick-Load Purple Low Molecular Weight DNA Ladder 168 Quick-Load DNA Ladders 168 Quick-Load 1 kb Extend DNA Ladder 168 TriDye DNA Ladders 168 PFG Ladders 168 Supercoiled DNA Ladder 170

Conventional DNA Markers	
Lambda DNA-Mono Cut Mix	169
Lambda DNA-HindIII Digest	169
Lambda DNA-BstEII Digest	169
φX174 DNA-HaeIII Digest	169
pBR322 DNA-BstNI Digest	169
pBR322 DNA-Mspl Digest	169
RNA MARKERS & LADDERS	
dsRNA Ladder	170
microRNA Marker	170
ssRNA Ladder	170
Low Range ssRNA Ladder	170
PROTEIN STANDARDS	
Unstained Protein Standard,	
Broad Range (10-200 kDa)	171
Color Prestained Protein Standard,	
Broad Range (10-250 kDa)	171
Blue Prestained Protein Standard,	
Broad Range (11-250 kDa)	171

Size Ranges of DNA Ladders from NEB



Purple Loading Dye


Gel Loading Dye, Purple (6X)

#B7024S 4.0 ml 43 €

Gel Loading Dye, Purple (6X), no SDS

#B7025S 4.0 ml 43 €

Our Gel Loading Dye, Purple (6X) (with and without SDS) is now supplied with all unstained DNA Ladders, sharpens bands and eliminates the UV shadow seen with other dyes. These pre-mixed loading buffers contain a combination of two dyes, Dye 1 (pink/red) and Dye 2 (blue). The red dye serves as the tracking dye for both agarose and non-denaturing polyacrylamide gel electrophoresis. The two dyes separate upon electrophoresis; the red band is the major indicator and migrates similarly to Bromophenol Blue on agarose gels. Specifically chosen, this dye does not leave a shadow under UV light. EDTA is also included to chelate magnesium (up to 10 mM) in enzymatic reactions, thereby stopping the reaction. The dyes also contain Ficoll, which creates brighter and tighter bands when compared to glycerol loading dyes. Gel Loading Dye, Purple (6X) contains SDS, which often results in sharper bands, as some restriction enzymes are known to remain bound to DNA following cleavage.



DNA Ladders

1 kb DNA Ladder
#N3232S 200 gel lanes ... 51 €
#N3232L 1,000 gel lanes204 €

100 bp DNA Ladder
#N3231S 100 gel lanes56 €
#N3231L 500 gel lanes224 €

1 kb Plus DNA Ladder

#N3200S 100-200 gel lanes 52 € #N3200L 500-1,000 gel lanes 208 €

50 bp DNA Ladder

#N3236S 100-200 gel lanes 63 € #N3236L 500-1,000 gel lanes 252 € Low Molecular Weight DNA Ladder

#N3233S 100 gel lanes 63 € #N3233L 500 gel lanes 252 €

Fast DNA Ladder


#N3238S 50-200 gel lanes 55 €

PCR Marker

#N3234S 100 gel lanes 63 € #N3234L 500 gel lanes 252 € NEB offers a variety of DNA Ladders with sizes ranging from 25 bp to 48.5 kb for use in agarose gel electrophoresis.

- Stable at room temperature
- Sharp, uniform bands
- Easy-to-identify reference bands
- Supplied with 1 vial of Gel Loading Dye, Purple (6X), no SDS
- Can be used for sample quantification (see www.neb.com for mass values)

Concentration: 1 kb Plus and 50 bp DNA Ladders are supplied at 1,000 μ g/ml. 1 kb, 100 bp and Low Molecular Weight DNA Ladders are supplied at 500 μ g/ml. PCR Marker is supplied at 300 μ g/ml. Fast DNA Ladder is supplied at 25 μ g/ml in ready-to-load format.



DNA Ladders in Convenient Pre-mixed Formats

Quick-Load® Purple Formats:

	·
#N0550S 125	1 kb Plus DNA Ladder* -250 gel lanes 69 € -750 gel lanes
Quick-Load Purple #N0552S #N0552L	1 kb DNA Ladder* 125 gel lanes 60 € 375 gel lanes
Quick-Load Purple #N0551S #N0551L	100 bp DNA Ladder* 125 gel lanes 80 € 375 gel lanes
	50 bp DNA Ladder* i-250 gel lanes86 €
Quick-Load Purple Weight DNA Ladde #N0557S	

* Supplied with 1 vial of Gel Loading Dye,

TriDye™ Formats:

,	us DNA Ladder
#N3270S	125-250 gel lanes 73
TriDye 1 kb DI	NA Ladder
#N3272S	125 gel lanes 62
TriDye 100 bp	DNA Ladder

125 gel lanes 83 €

€.

€

Ouick-Load Formate

#N3271S


QUICK-LUAU	ruillais.
Quick-Load 1	kb Plus DNA Ladder
#N0469S	125-250 gel lanes71 €
Quick-Load 1	kb DNA Ladder
#N0468S	125 gel lanes 61 €
#N0468L	375 gel lanes 146 €
Quick-Load 1	kb Extend DNA Ladder
#N3239S	125 gel lanes 67 €

Quick-Load 100 bp DNA Ladder

#N0467S 125 gel lanes 80 € #N0467L 375 gel lanes 200 €

- Ready-to-load
- Stable at 25°C
- Uniform band intensities
- Easy-to-identify reference bands
- Defined mass profile for sample quantification

Our 1 kb Plus, 1 kb and 100 bp DNA Ladders are offered in four formats. Choose from the conventional ladder, the Quick-Load version using either non-fluorescing, purple dye or bromophenol blue as a tracking dye, or TriDye containing three dyes to facilitate monitoring of gel migration.





PFG Ladders

Lambda PFG Ladder

Purple (6X), no SDS.

#N0341S 50 gel lanes 160 €

NEW

Midrange PFG Marker

#N0342S 50 gel lanes 160 €

The Lambda PFG Ladder consists of one gelsyringe dispenser, sufficient for 50 gel lanes. Successively larger concatemers of lambda DNA (cl857 ind 1 Sam7) are embedded in 1% LMP agarose. This product is designed to be used as size markers for pulsed-field gel electrophoresis (PFG). Size range: 48.5–1,018 kb.

MidRange PFG Marker consists of concatemers of λ DNA isolated from the bacteriophage λ (cl857 ind1

Kilobases 242.5 -227.5 -209 0 -Kilobases 194.0 727.5 679.5 630.5 179.0 -160.5 -582.0 145 5 -485.0 — 130.5 -436.5 — 388.0 -112.0 -339.5 -97.0 -291.0 -82 0 -242.5 -194.0 63.5 -145.5 48.5 -33.5 -97.0 48.5 15.0 Lambda PFG Midrange PFG Ladder Marker 1% agarose gel, 4.5 V/cm, 1% agarose gel, 6 V/cm, 15°C for 48 hours. Switch 15°C for 24 hours. Switch times ramped from 5-120 seconds. times ramped from 1-25 seconds. Sam7) mixed with Xhol digested λ DNA embedded in 1% LMP agarose and supplied in a gelsyringe dispenser. Xhol produces fragments of 15.0 and 33.5 kb. These fragments anneal to and form concatemers with intact λ DNA. It is designed for use as a size marker for pulsed field gel electrophoresis (PFG). Size range: 15-291 kb.

Concentration: 50 µg/ml.

Usage Recommendations

Lambda PFG Ladder: The photograph represents the pulsed field gel separation of Lambda PFG Ladders using a CHEF apparatus, for 48 hours at 15°C in 0.5X TBE made with Milli-Q® water allowing resolution of 15 to 21 Lambda PFG Ladder bands.

Midrange PFG Marker: The photograph represents the pulsed field gel separation of Mid Range PFG Marker using a CHEF apparatus, for 24 hours at 15°C in 0.5X TBE made with Milli-Q water.

Usage Notes for PFG Ladders: Recommended plug sizes are from 5-10 µl. A 10 µl plug (one small graduation on the GelSyringe volume scale) contains approximately 0.5 µg of DNA. Each gelsyringe yields 50+ plugs.

Extrude agarose from gelsyringe carefully and slice plugs from the end with a sharp blade. One plug is sufficient for one lane of a gel. Place the plug at the front of the well and seal with molten agarose just above gelling temperature (-42–45°C). Allow no bubbles to form.

Melting plugs will cause denaturation of concatemers.

Never attach the agarose plugs to the gel comb before the gel is poured. Heat from the solidifying gel will cause the Lambda concatemers to denature.

MILLI-Q® is a registered trademark of Millipore Corporation.

Conventional DNA Markers

Lambda DNA-Mono Cut Mix
#N3019S 100 gel lanes70 €

Lambda DNA-HindIII Digest
#N3012S 150 gel lanes68 €
#N3012L 750 gel lanes272 €

Lambda DNA-BstEII Digest
#N3014S 150 gel lanes70 €

φX174 DNA-HaelII Digest

#N3026S 50 gel lanes76 € #N3026L 250 gel lanes304 €

pBR322 DNA-BstNI Digest

#N3031S 50 gel lanes73 € #N3031L 250 gel lanes292 €

pBR322 DNA-Mspl Digest

#N3032S 50 gel lanes73 € #N3032L 250 gel lanes292 € NEB offers a wide range of double-stranded DNA molecular weight markers for conventional agarose gel electrophoresis. These standards have a size range of approximately 10–23,000 base pairs.

The typical pattern generated by each of the conventional markers is shown below. The number of fragments generated by each marker, as well as the specific fragment sizes for each of the conventional markers can be found on the datacard, as well as our website, www.neb.com.


In addition, the conventional markers can be used for mass approximation. See the datacard or website for DNA mass information

The Lambda DNA-Mono Cut Mix is best separated by Pulsed Field Gel Electrophoresis. It is designed to be used as an RFLP marker on Southern blots because it provides an uncomplicated gel pattern.

Concentration: pBR322 DNA-BstNI Digest, pBR322 DNA-MspI Digest and φX174 DNA-HaeIII Digest are supplied at 1,000 μg/ml. Lambda DNA-BstEII Digest, Lambda DNA-HindIII Digest and Lambda DNA-Mono Cut Mix are supplied at 500 μg/ml.

Usage Recommendation: Dilution of these markers is recommended in TE or other buffer of minimal ionic strength. DNA may denature if diluted in dH_nO.


The cohesive ends of fragments 1 and 4 of the Lambda DNA-HindIII and Lambda DNA-BstEII Digests can be separated by heating to 60°C for 3 minutes.



Supercoiled DNA Ladder

#N0472S 100 gel lanes 105 €

Mass (ng) Kilobases



Ladder0.5 μg/lane.
0.8% TAE agarose gel.

The Supercoiled DNA ladder contains 9 proprietary supercoiled plasmids, ranging in size from 2 to 10 kb, that are suitable for use as supercoiled molecular weight standards for agarose electrophoresis. The 5 kb plasmid has an increased intensity to serve as a reference band.

Source: The 9 proprietary plasmids are purified, phenol extracted and equilibrated to 10 mM Tris-HCl (pH 8.0) and 1 mM EDTA.

Concentration: 500 µg/ml.

Notes: This ladder may contain some traces of nicked DNA and dimers above the 10 kb plasmid. To minimize nicking of the supercoiled DNA, always use sterile pipette tips and avoid multiple freeze-thaw cycles. The migration of supercoiled plasmids in agarose gels can change depending on agarose concentration, buffer and electrophoresis conditions. Dilute in TE or other buffer of minimal ionic strength. DNA may denature if diluted in dH_nO.

Usage Recommendation: Centrifuge briefly and mix gently before use. We recommend loading $0.5 \,\mu g \, (1 \,\mu l)$ of the Supercoiled DNA Ladder diluted in sample buffer. This ladder was not designed for precise quantification of DNA mass, but can be used for approximating the mass of DNA in comparably intense samples of similar size. The approximate mass of DNA in each of the bands in our Supercoiled DNA ladder is as follows (assuming a $0.5 \,\mu g$ loading):

Band	Base Pairs	DNA Mass
1	10,000	45 ng
2	8,000	45 ng
3	6,000	45 ng
4	5,000	136 ng
5	4,000	45 ng
6	3,500	45 ng
7	3,000	45 ng
8	2,500	45 ng
9	2,017	45 ng

RNA Markers & Ladders

dsRNA Ladder #N0363S

25 gel lanes 90 €

microRNA Marker

#N2102S 100 gel lanes 69 €

ssRNA Ladder

#N0362S 25 gel lanes 69 €

Low Range ssRNA Ladder


#N0364S 100 gel lanes 68 €

NEB offers several RNA Markers and Ladders with a size range from 17 to 9,000 bases. The Low Range ssRNA Ladder and the ssRNA Ladder are suitable for use as RNA size standards on denaturing or native gels. Both are supplied with RNA Loading Dye (2X) (NEB #B0363) and feature a higher intensity fragment to serve as a reference band. The microRNA Marker, provided in a ready-to-load denaturing solution, is ideally used as a size marker on denaturing polyacrylamide gels or northern blots and is best visualized stained with SYBR®-Gold. It is supplied

with a 3´-biotinylated 21-mer oligonucleotide probe that can be labeled with \(\gamma 32-P-ATP \) and T4 PNK (NEB #M0201). The ds RNA Ladder is suitable for use as a size standard in dsRNA and RNAi analysis on both polyacrylamide and agarose gels.

Concentration: Low Range ssRNA Ladder and dsRNA Ladder are supplied at $500 \mu g/ml$. ssRNA Ladder is supplied at $2,000 \mu g/ml$. MicroRNA Marker is supplied at $12 ng/\mu l$.

SYBR® is a registered trademark of Molecular Probes, Inc.



Protein Standards

NEW Unstained Protein Standard, Broad Range (10-200 kDa)

#P7717S 150 gel lanes 128 € #P7717L 750 gel lanes 498 €

NEW

Color Prestained Protein Standard, Broad Range (10-250 kDa)

150 gel lanes155 € #P7719S #P7719L 750 gel lanes 630 €

NEW

Blue Prestained Protein Standard, Broad Range (11-250 kDa)

#P7718S 150 gel lanes 140 € #P7718L 750 gel lanes 588 € Companion Products:


Blue Loading Buffer Pack 3X Blue Loading Buffer (8 ml) 30X Reducing Agent (1 ml)

#B7703S27 €

NEB offers a selection of highly pure protein standards available as unstained, blue prestained or color prestained (containing two colored reference bands for easy identification). Sizes range from 10 to 250 kDa which is ideal for calculating molecular weight determination for a wide range of expressed proteins. NEB protein standards provide uniform band intensities, convenient band spacing and easy-toidentify reference bands.

Recommended Load Volume: $3 \, \mu l$

Note: For calculating molecular weight determinations, use NEB's Unstained Protein Standard, Broad Range.







The Ripple Effect

The removal or extinction of one species in an ecosystem can create gaps in the food chain. This gap can cause a "ripple" throughout the ecosystem, sequentially threatening different species due to disturbances in their food source, habitat or specialized relationships with other species.

In the late 1970s, the Large Blue Butterfly (*Maculinear arion*) was declared extinct in Britain. This beautiful, rare species was highly sought after in the 19th century by butterfly enthusiasts, and consequently, its numbers declined. Conservation efforts to fend off collectors and preserve Large Blue colonies began in the 1920s, and part of this effort included the removal of grazing sheep.

The Large Blue has a complex, predatory relationship with the red ant, *Myrmica sabuleti*. The butterfly lays its eggs on the leaves of wild thyme plants, which is where the caterpillar feeds before eventually dropping to the ground. It deceptively hunches its body so that it appears to be the size of the ant larvae and produces a pheromone which mimics that of the red ant grub. The pheromone attracts the red ant, fooling it into adopting the caterpillar and transporting it to its underground nest where it feeds on ant larvae. The caterpillar will only survive in the nest of the *Myrmica sabuleti*, all other red ant species in this environment will recognize it as an impostor and kill it. The caterpillars go as far as to mimic the sound of the Queen ants so that the worker ants will feed and clean them. At least 230 ant larvae and 354 ant workers are required to guarantee the survival of just one Large Blue Butterfly.

Researcher Jeremy Thomas carefully observed the last surviving colony of Large Blue Butterflies from 1972 until it disappeared in 1977. He pieced together the butterfly's intricate, exploitative relationship with the red ant and also concluded that a ripple throughout the ecosystem began with the seemingly minor modification of removing grazing sheep. This habitat felt a ripple when the grass grew longer and the soil temperature decreased by just 1 or 2°C, resulting in an unsuitable habitat for the temperature- and humidity-sensitive red ant. The red ant population plummeted, and this left its predator, the very visible and appreciated Large Blue, without its food source. The Large Blue's lifecycle was no longer viable now that the incredibly specialized environment in which it had survived, was gone.

A habitat recovery project involved re-establishing grazing animals on over 100 sites. Large Blue larvae were imported from Sweden and placed in red ant nests in the 1980s and 1990s. A period of optimization of grazing conditions was required before the red ant, and then the Large Blue, thrived.

In recent years, the number of insects in Britain and Europe has declined more than plants and birds. Thomas' findings highlighted the need to focus on the ecological conditions that cause insect numbers to decrease and to integrate their needs with modern land use. The intensive research project to restore the habitat of the Large Blue was the first effort to reverse the decrease in numbers of an endangered insect species, which has helped other declining insect populations, and hopefully will help to prevent other such "ripple effects".

Genome Editing



Programmable nucleases for your applications.

Easily changing the sequence specificity of a DNA binding protein enables many new possibilities for the detection and manipulation of DNA, including for genome editing – creating targeted changes in the DNA of living cells. CRISPR (Clustered Regularly Interspaced Short Palindromic Repeats) -associated (Cas) nucleases have been adapted from bacterial immune systems into useful tools for biotechnology. Cas nucleases including Cas9 and Cas12a, are attractive for genome editing because they are easily programmable with a single guide RNA or crRNA/tracrRNA to introduce a double-stranded break at a specific target. The resultant breaks are repaired by cellular machinery, in some applications using an exogenous repair template.

The ease of programming Cas nucleases, and their conversion into nicking endonucleases, or DNA binding proteins without nuclease activity has expanded their use to include delivering a specific cargo to a locus for applications including visualization, activation, repression and base editing. Furthermore, Cas nuclease and their variants are useful tools for the detection and manipulation of DNA *in vitro*.

Featured Products

178 EnGen Mutation Detection Kit

179 EnGen sgRNA Synthesis Kit, S. pyogenes

177 EnGen Lba Cas12a (Cpf1)

177 EnGen Sau Cas9

Featured Tools & Resources



Visit www.neb.com/GenomeEditing for more information, including our feature article and latest brochure.





FEATURED PRODUCTS SUPPORTING CRISPR WORKFLOWS

EnGen Spy Cas9 NLS	177	™ Q5 Site-Directed Mutagenesis Kit	
EnGen Mutation Detection Kit	178	(with or without competent cells)	92
EnGen sgRNA Synthesis Kit, S. pyogenes	179	Q5 High-Fidelity DNA Polymerases	63
EnGen Spy Cas9 Nickase	177	■ NEBuilder HiFi DNA Assembly Master Mix	88
EnGen Spy dCas9 (SNAP-tag)	177	■ NEBuilder HiFi DNA Assembly Cloning Kit	88
EnGen Lba Cas12a (Cpf1)	177	HIScribe T7 ARCA mRNA Kit (with or without tailing	186
EnGen Sau Cas9	177	IIII HiScribe T7 High Yield RNA Synthesis Kit	185
Cas9 Nuclease, S. pyogenes	177	HiScribe T7 Quick High Yield RNA Synthesis Kit	185
Monarch Total RNA Miniprep Kit	130	I HiScribe SP6 RNA Synthesis Kit	185
		T7 Endonuclease I	113
	EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i> EnGen Spy Cas9 Nickase EnGen Spy dCas9 (SNAP-tag) EnGen Lba Cas12a (Cpf1) EnGen Sau Cas9 Cas9 Nuclease, <i>S. pyogenes</i>	EnGen Mutation Detection Kit 178 EnGen sgRNA Synthesis Kit, S. pyogenes 179 EnGen Spy Cas9 Nickase 177 EnGen Spy dCas9 (SNAP-tag) 177 EnGen Lba Cas12a (Cpf1) 177 EnGen Sau Cas9 177 Cas9 Nuclease, S. pyogenes 177	EnGen Mutation Detection Kit EnGen sgRNA Synthesis Kit, S. pyogenes EnGen Spy Cas9 Nickase EnGen Spy dCas9 (SNAP-tag) EnGen Spy dCas9 (SNAP-tag) EnGen Lba Cas12a (Cpf1) EnGen Sau Cas9 177 EnGen Sau Cas9 178 EnGen Sau Cas9 179 EnGen Sau Cas9 170 EnGen Sau Cas9 171 EnGen Sau Cas9 172 EnGen Sau Cas9 173 EnGen Sau Cas9 174 EnGen Sau Cas9 175 EnGen Sau Cas9 176 EnGen Sau Cas9 177 EnGen Sau Cas9 EnGen Sau Cas

Recombinant Enzyme

Featured NEB Products Supporting CRISPR Workflows

New England Biolabs provides reagents to support a broad variety of CRISPR/Cas genome editing approaches. From introduction of Cas and single guide RNA (sgRNA) on plasmids, to direct introduction of Cas ribonucleoprotein (RNP) and detection of edits using next generation sequencing or enzymatic mutation detection, NEB® provides reagents that simplify and shorten genome editing workflows.

PRODUCT	CRISPR/Cas9 APPLICATION	NEB #	SIZE	PRICE
EnGen® Spy Cas9 NLS	in vitro cleavage of dsDNA. Genome engineering by direct introduction of active ribonucleoproteins	M0646T/M	400/2,000 pmol	158 €/632 €
EnGen Mutation Detection Kit	Determination of the targeting efficiency of genome editing protocols	E3321S	25 rxns	210€
EnGen sgRNA Synthesis Kit, S. pyogenes	Generation of microgram quantities of custom sgRNA	E3322S	20 rxns	420€
<mark>NEW</mark> EnGen Spy Cas9 Nickase	in vitro nicking of dsDNA. Genome engineering by direct introduction of active nickase complexes.	M0650S/T	70/400 pmol	63 €/158 €
NEW EnGen Spy dCas9 (SNAP-tag°)	Programmable binding of DNA. Compatible with SNAP-tag substrates for visualization and enrichment.	M0652S/T	70/400 pmol	63 €/158 €
NEW EnGen Lba Cas12a (Cpf1)	in vitro cleavage of dsDNA. Genome engineering by direct introduction of active nuclease complexes. Recognizes 5´-TTTN PAM.	M0653S/T	70/2,000 pmol	74 €/264 €
<mark>NEW</mark> EnGen Sau Cas9	in vitro cleavage of dsDNA. Genome engineering by direct introduction of active nuclease complexes. Recognizes 5'-NNGRRT-3' PAM.	M0654S/T	70/400 pmol	60 €/150 €
Cas9 Nuclease, S. pyogenes	in vitro cleavage of dsDNA. Genome engineering by direct introduction of active ribonucleoproteins.	M0386S/T/M	70/400/ 2,000 pmol	53 €/140 €/ 550 €
Monarch® Total RNA Miniprep Kit	Purification of total RNA, with a binding capacity of up to 100 μg	T2010S	50 preps	248€
Q5° Site-Directed Mutagenesis Kit (with or without competent cells)	Insertion of target sequence into a Cas9-sgRNA construct and modification of HDR templates	E0554S/ E0552S	10 rxns	192 €/132 €
Q5 High-Fidelity DNA Polymerases	High-fidelity construct generation for use with CRISPR workflows and for sequencing	Multiple	Multiple	see page 63
NEBuilder® HiFi DNA Assembly Master Mix	Single-tube, isothermal generation of the Cas9-sgRNA construct and HDR templates	E2621S/L/X	10/50/250 rxns	162 €/648 €/ 2589 €
NEBuilder HiFi DNA Assembly Cloning Kit	Single-tube, isothermal generation of the Cas9-sgRNA construct and HDR templates	E5520S	10 rxns	190€
HiScribe T7 ARCA mRNA Kit (with or without tailing)	Generation of Cas9 mRNA with ARCA cap	E2060S/ E2065S	20 rxns	376 €/322 €
HiScribe T7 High Yield RNA Synthesis Kit	Generation of sgRNA and Cas9 mRNA	E2040S	50 rxns	223€
HiScribe T7 Quick High Yield RNA Synthesis Kit	Generation of sgRNA and Cas9 mRNA	E2050S	50 rxns	273€
T7 Endonuclease I	Determination of the editing efficiency of genome editing experiments	M0302S/L	250/1,250 units	69 €/276 €

Programmable Cas Nucleases

The highest efficiency strategy for genome editing with CRISPR/Cas nucleases is direct introduction of Cas/guide RNA complexes.

This method further simplifies CRISPR/Cas workflows and has been reported to increase on-target editing activity and reduce off-target events.

NEB provides purified Cas9 nucleases from *S. pyogenes* and *S. aureus*, and Cas12a (Cpf1) nuclease from *Lachnospiraceae* bacterium ND2006.

In addition, NEB provides variants of Cas9 from *S. pyogenes*, including nicking endonuclease and endonuclease deficient versions.

PRODUCT	NEB #	FEATURES	SIZE	PRICE
Cas9 Nuclease, S. pyogenes	M0386S/T/M	Ideal for in vitro digestion of dsDNA Compatible with EnGen sgRNA Synthesis Kit, S. pyogenes (NEB #E3322) and the EnGen Mutation Detection Kit (NEB #E3321) For help with oligo design, try EnGen sgbNA Template Oligo Designer – now included in the NEBioCalculator® Tool	70/300/600 pmol	53 €/140 €/550 €
EnGen Spy Cas9 NLS	M0646T/M	Ideal for direct introduction of Cas9/sgRNA complexes Dual NLS for improved transport to the nucleus Compatible with EnGen sgRNA Synthesis Kit, S. pyogenes (NEB #E3322) and the EnGen Mutation Detection Kit (NEB #E3321)	400/2,000 pmol	158 €/632 €
NEW EnGen Spy Cas9 Nickase	M0650S/T	 Variant of Cas9 nuclease differing by a point mutation (D10A) in the RuvC nuclease domain Capable of generating nicks, but not cleaving DNA DNA double strand breaks can be generated, with reduced off-target cleavage, by targeting two sites with EnGen Cas9 Nickase in close proximity Compatible with the EnGen sgRNA Synthesis Kit, S. pyogenes (NEB #E3322) 	70/400 pmol	63 €/158 €
NEW EnGen Spy dCas9 (SNAP-tag)	M0652S/T	An inactive mutant of Cas9 nuclease that retains programmable DNA binding activity The N-terminal SNAP-tag allows for covalent attachment of fluorophores, biotin, and a number of other conjugates useful for visualization and target enrichment Compatible with the EnGen sgRNA Synthesis Kit, S. pyogenes (NEB #E3322)	70/400 pmol	63 €/158 €
NEW EnGen Lba Cas12a (Cpf1)	M0653S/T	 T-rich TTTN PAM sequence opens up additional genomic regions for targeting Shorter, 40-44 base guide RNA Two nuclear localization signals for improved transport to the nucleus 5' overhanging termini on cleavage products Active from 16 to 48°C Maintains activity at lower temperatures than the Acidaminococcus orthologs, permitting editing in ectothermic organisms such as zebra fish and xenopus High concentration enzyme can be used for microinjection, electroporation and lipofection 	70/2,000 pmol	74€/264€
NEW EnGen Sau Cas9	M0654S/T	 5´-NNGRRT-3´ PAM Dual NLS for improved transport to nucleus High concentration enzyme can be used for microinjection electroporation and lipofection Cleaves 3 bases upstream of PAM, blunt-ended cleavage 	70/400 pmol	60 €/150 €

EnGen Mutation Detection Kit

25 reactions 210 €

500 reactions 720 €

#E3321S

#M0494L

#T1030S 50 reactions ... 88 €
#T1030L 250 reactions ... 398 €

EnGen Spy Cas9 NLS
#M0646T 400 pmol ... 158 €
#M0646M 2,000 pmol ... 632 €

Monarch PCR & DNA Cleanup Kit (5 µg)

EnGen Spy Cas9 Nickase
#M0650S 70 pmol63 €
for high (20X) concentrationn
#M0650T 400 pmol158 €

EnGen Spy dCas9 (SNAP-tag) #M0652S 70 pmol 63 € for high (20X) concentration #M0652T 400 pmol 158 € Cas9 Nuclease, S.pyogenes

#M0386S 70 pmol53 € for high (20X) concentration #M0386T 400 pmol140 € for high (20X) concentration #M0386M 2,000 pmol550 € EnGen Lba Cas12a (Cpf1)

#M0653S 70 pmol \dots 74 \in for high (100X) concentration #M0653T 2,000 pmol \dots 264 \in

EnGen Sau Cas9
#M0654S 70 pmol60 €
for high (20X) concentration
#M0654T 400 pmol150 €
T7 Endonuclease I

T7 Endonuclease-based detection of

genome editing events


Description: The EnGen Mutation Detection Kit provides reagents for detection of on-target genome editing events. In the first step, regions targeted for genome editing (i.e. CRISPR/Cas9, TALENs, Zinc-finger Nucleases) are amplified using Q5 Hot Start High-Fidelity 2X Master Mix. Upon denaturation and re-annealing, heteroduplexes are formed when mutations from insertions and deletions (indels) are present in the amplicon pool. In the second step, annealed PCR products are digested with EnGen T7 Endonuclease I, a structure-specific enzyme that will recognize mismatches larger than 1 base. Both strands of the DNA are cut when a mismatch is present, which results in the formation of smaller fragments. Analysis of the resulting fragments provides an estimate of the efficiency of the genome editing experiments.

The EnGen Mutation Detection Kit includes a Control Template and Primer Mix that can be used as a control for the PCR reaction and T7 Endonuclease I digestion. The Control Template and Primer Mix provided contains two plasmids and primers that when amplified, denatured and re-annealed will form heteroduplexes that contain a 10-base insertion. This structure is a substrate for T7 Endonuclease I. The digestion of the 600 bp heteroduplex containing amplicon yields products of 200 bp and 400 bp. 600 bp parental homoduplexes are uncleaved, and are easily distinguished from cleaved heteroduplexes when separated and visualized by agarose gel electrophoresis or fragment analysis instrument.

The protocol has been optimized so that PCR products generated by the Q5 Hot Start High-Fidelity 2X Master Mix can be introduced directly into the T7 Endonuclease I digestion without the need for purification. Digestion of the heteroduplex is complete in only 15 minutes, and Proteinase K is included to stop the reaction efficiently. Additional Q5 Hot Start High-Fidelity 2X Master Mix is also included to allow for optimization of target site amplification before digestion.

The EnGen Mutation Kit Includes:

- Q5 Hot Start High-Fidelity 2X Master Mix
- NEBuffer 2
- EnGen T7 Endonuclease I
- Control Template and Primer Mix
- Proteinase K., Molecular Biology Grade
- Quick-Load® Purple 1 kb Plus DNA Ladder
- Gel Loading Dye, Purple (6X), no SDS



Genomic DNA is amplified with primers bracketing the modified locus. PCR products are then denatured and re-annealed yielding three classes of possible structures. Duplexes containing a mismatch greater than one base are digested by T7 Endonuclease I. The DNA is then electrophoretically separated and fragment analysis is used to estimate targeting efficiency.

EnGen sgRNA Synthesis Kit, S. pyogenes

#E3322S 20 reactions420 €

Companion Products:

Johnpanion Products.						
EnGen Spy Cas9 NLS #M0646T #M0646M		158 €				
EnGen Spy Cas9 Nickase #M0650S for high (20X) concentration	70 pmol	63 €				
#M0650T	400 pmol	158 €				
EnGen Spy dCas9 (SNAF	0,					
#M0652S	70 pmol	63 €				
for high (20X) concentration #M0652T	400 pmol	158 €				
Cas9 Nuclease, <i>S.pyoge</i> #M0386S		53 €				
for high (20X) concentration #M0386T	400 pmol	140 €				
for high (20X) concentration #M0386M	2,000 pmol	550 €				
RNA Loading Dye, (2X) #B0363S	4 ml	49€				
Alkaline Phosphatase, Ca #M0290S #M0290L	1,000 units	(CIP) 72 € 288 €				
Monarch Total RNA Mini #T2010S		248 €				

 Rapid generation of microgram quantities of sgRNAs in less than one hour

NEBioCalculator®

Configure target-specific DNA oligos design for use with the EnGen sgRNA Synthesis Kit, *S. pyogenes* with our oligo design tool accessible through NEBioCalculator® at NEBioCalculator.neb.com


Description: The EnGen sgRNA Synthesis Kit, S. pyogenes provides a simple and quick method for transcribing high yields of sgRNA in a single 30 minute reaction, using the supplied reagents and target-specific DNA oligos designed by the user.

In nature, *S. pyogenes* Cas9 is programmed with two separate RNAs, the crRNA and tracrRNA. The crRNA, or CRISPR RNA sequence contains approximately 20 nucleotides of homology complementary to the strand of DNA opposite and upstream of a PAM (Protospacer Adjacent Motif) (NGG) sequence. The tracrRNA, or transactivating crRNA, contains partial complementary sequence to the crRNA as well as the sequence and secondary structure that is recognized by Cas9. These sequences have been adapted for use in the lab by combining the tracrRNA and crRNA into one long single guide RNA (sgRNA) species capable of complexing with Cas9 to recognize and cleave the target DNA.

The EnGen sgRNA Synthesis Kit, *S. pyogenes* combines an *S. pyogenes* Cas9-specific Scaffold Oligo (included in the EnGen 2X sgRNA Reaction Mix) that is partially complementary to the target-specific oligos designed by the user. The two oligos anneal at the overlapping region and are filled in by the DNA polymerase, creating a double-stranded DNA (dsDNA) template for transcription. Synthesis of the dsDNA template and transcription of RNA occur in a single reaction, resulting in the generation of a functional sgRNA.

The EnGen sgRNA Synthesis Kit, S. pyogenes Includes:

- EnGen sgRNA Enzyme Mix
- EnGen 2X sgRNA Reaction Mix, S. pyogenes
- DNase I (RNase-free)
- EnGen sgRNA Control Oligo, S. pyogenes
- DTT (0.1 M)



A. The target-specific oligo contains the T7 promoter sequence, ~20 nucleotides of target-specific sequence and a 14 nucleotide overlap sequence complementary to the S. pyogenes Cas9 Scaffold Oligo supplied in the reaction mix. Target-specific oligos (or EnGen sgRNA Control Oligo, S. pyogenes) are mixed with the EnGen 2X sgRNA Reaction Mix (NTPs, dNTPs, S. pyogenes Cas9 Scaffold Oligo), 0.1 M DTT and the EnGen sgRNA Enzyme Mix (DNA and RNA polymerases) at room temperature. B. At 37°C the two oligos anneal at the 14 nucleotide overlap region of complementarity. C. The DNA polymerase extends both oligos from their 3' ends creating a double-stranded DNA template. D. The RNA polymerase recognizes the double-stranded DNA of the T7 promoter and initiates transcription. The resulting sgRNA contains the target-specific/crRNA sequence as well as the tracrRNA. All steps occur in a single reaction during a 30 minute incubation at 37°C.





The Downstream Effects of Overfishing

Fishing has been a mainstay throughout the history of civilization. However, the methods, equipment and extent of fishing have drastically changed in order to feed an overpopulated world. As a result, we are rapidly depleting our oceans of what was once thought to be an endless supply of fish. What is not as readily understood is that gaps in the food web created by overfishing certain species, and the currently-used industrial fishing methods, such as bottom trawling, are destroying entire ecosystems.

The global fishing industry is worth nearly 250 billion dollars per year, which equates to 90 million tons of fish. Eighty percent of commercial fish stocks have been declared overexploited, yet this is not reflected in the variety of fish available in local markets. This is due to the fact that industrial-scale fishing efforts have moved further offshore and to deeper depths of the ocean to find new sources. "Super trawlers" the size of ocean liners use probes, radar, sonar, helicopters and spotter planes to hunt down marine life. Nets that can measure 40 km long catch targeted species of fish, as well as non-targeted species ("bycatch"), such as seabirds, turtles and dolphins, which are thrown back into the sea, often dead or dying.

Modern fishing techniques completely devastate ocean habitats. Bottom trawlers drag nets along the ocean floor, destroying rich, complex ecosystems and communities of invertebrates on the sea bed as they comb for scallops and shrimp. They leave behind a barren landscape of sand and gravel. The disregard for entire ecosystems also extends to coastal areas where mangroves are destroyed for shrimp fishing, leaving local coastal communities without storm protection, a natural form of water filtration, or nursery habitats for marine life.

Large deep-water fish, such as Bluefin tuna, are commercially very valuable and are aggressively targeted. They take longer to reach sexual maturity and do not reproduce as often. They are unable to replenish themselves fast enough to meet demand and consequently, younger fish that have not yet reached spawning age are targeted. Depleting the large and predatory fish significantly affects food web dynamics. The smaller, more resilient fish remain, and eventually, the range of fish observed are smaller and less diverse.

The Aleutian Islands, off the coast of northern Alaska, saw the effects of overfishing on an entire ecosystem when baleen and sperm whales were overfished, removing a food source for Orcas. This forced the Orcas to feed on less calorie-rich animals, including otters. The primary prey of otters is sea urchin, which in turn, feed on kelp. The decline in otters led to a drastic increase in kelp grazing by urchin, and the destruction of kelp forests. When nesting fish do not have the protection of kelp, their larvae are vulnerable to predation. In the case of the Aleutian Islands, this ultimately led to the collapse of the fishing industry.

Overfishing affects the balance of marine ecosystems and the livelihood of millions of people. As awareness grows, steps are being taken to protect our oceans. Marine reserves, or "no fishing" zones are being established so that fish stocks can recover. Organizations such as the Marine Stewardship Council certify sustainably caught fish, indicating to conscientious consumers the products that support restoration of marine life. Supermarket chains are taking on a vital role offering sustainably caught seafood to customers. These measures give hope that with even greater awareness we may just be able to restore this essential resource.

RNA Reagents



A broad portfolio of reagents to support RNA research.

RNA molecules play a multitude of cellular roles in all kingdoms of life, perhaps reflecting the hypothetical, prehistoric RNA world. RNA is an essential carrier of genetic information, can scaffold molecular interactions, catalyze chemical reactions, and influences gene expression.

In the last several years, our understanding of RNAs as regulatory molecules in the cell has dramatically changed. Many new classes of small and large non-coding RNAs, with largely unexplored functions, have been reported, ushering in a renaissance of RNA-focused research in biology.

Small RNAs play a major role in the post-transcriptional regulation of gene expression in eukaryotic and prokaryotic organisms. Small RNAs play central roles in CRISPR pathways of adaptive immunity. Researchers have capitalized on these pathways to enable analyses of gene function not previously possible.

Introducing RNA, and in particular mRNA, into cells is an exciting new platform for inducing phenotypic changes in cells with implications for therapy, vaccines, and research applications.


New England Biolabs continues its strong tradition of providing high quality reagents to support RNA research. Our expanding range of products includes tools for the synthesis, processing, cleanup, isolation, analysis, amplification, copying and cloning of RNA molecules. Further, all NEB products pass stringent quality control assays to ensure the highest level of functionality and purity.

Featured Products

- 186 HiScribe™ T7 ARCA mRNA Kit (with tailing)
- 185 HiScribe T7 High Yield RNA Synthesis Kits
- 193 Template Switching RT Enzyme Mix
- 188 Hi-T7™ RNA Polymerase

Featured Tools & Resources

- **184** Avoiding Ribonuclease Contamination
- 195 Reported Activities for RNA Ligases
- View our video tutorial describing high yield *in vitro* synthesis of both capped and uncapped mRNA.
- View our video for avoiding ribonuclease contamination.





	Avaiding Dihanualagas Contamination	104	-Di))	DNA E' Dyranhaanhabydralaaa (Danii)	100
	Avoiding Ribonuclease Contamination			RNA 5´ Pyrophosphohydrolase (RppH) 5´ Deadenylase	199 199
	RNA Synthesis			RNase I,	200
RR	HiScribe T7 Quick High Yield RNA Synthesis Kit	105		RNase H	200
RR	HiScribe T7 High Yield RNA Synthesis Kit	105		Thermostable RNase H	200
RR	HiScribe SP6 RNA Synthesis Kit	105		RNase HII	200
RR	HiScribe T7 ARCA mRNA Kit	100		Quick Dephosphorylation Kit	99
RR	HiScribe T7 ARCA mRNA Kit (with tailing)	100		Antarctic Phosphatase	100
RR	EnGen sgRNA Synthesis Kit, S. pyogenes	179		Alkaline Phosphatase Calf Intestinal (CIP)	100
	Recommended HiScribe Kits by Application	187	RR	Shrimp Alkaline Phosphatase (rSAP)	100
RR	T3 RNA Polymerase	188	RR	T4 Polynucleotide Kinase	98
RR	T7 RNA Polymerase	188	Rii	ShortCut RNase III	201
RR	SP6 RNA Polymerase	400		XRN-1	201
RR	Hi-T7 RNA Polymerase	188			
RR	E. coli Poly(A) Polymerase	188		RNase Control	
RR	Poly(U) Polymerase			RNase Inhibitor, Murine	202
			RR	RNase Inhibitor, Human Placenta	203
	Ribonucleotide Solution Set	189		Ribonucleoside Vanadyl Complex	203
	Ribonucleotide Solution Mix	189		RNA Library Preparation for	
	Pyrophosphatase, Inorganic (<i>E. coli</i> , yeast)	189		Next Generation Squencing	
	Thermostable Inorganic Pyrophosphatase	189		NEBNext Ultra II Directional RNA	
	Vaccinia Capping System	190		Library Prep Kit for Illumina	142
RR	mRNA Cap 2'-O-Methyltransferase	190		NEBNext Ultra II Directional RNA Library	
	RNA Cap Analogs	190		Prep with Sample Purification Beads	142
	3´-Desthiobiotin-GTP 3´-Biotin-GTP	191		NEBNext Ultra II RNA Library Prep Kit for Illumina	142
-Dill		191		NEBNext Ultra II RNA Library Prep	
hii	yDcpS	191		with Sample Purification Beads	142
	cDNA Synthesis, RT-PCR & RT-qPCR			NEBNext Single Cell/Low Input RNA	144
	cDNA Synthesis Selection Chart	191		Library Prep Kit for Illumina NEBNext Single Cell/Low Input cDNA	144
RR	ProtoScript II Reverse Transcriptase	192		Synthesis & Amplification Module	144
	LunaScript RT SuperMix Kit	70	Dill		
	M-MuLV Reverse Transcriptase	192	Kii	Exonuclease T	202
	AMV Reverse Transcriptase	192	D))	Nucleoside Digestion Mix DNase I (RNase-free)	202
RR	WarmStart RTx Reverse Transcriptase	193	m	Divase i (nivase-liee)	202
RR	Template Switching RT Enzyme Mix	193		RNA Cleanup & Isolation	
	Primers for cDNA Synthesis	193		EpiMark N6-Methyladenosine Enrichment Kit	204
RR	ProtoScript II First Strand cDNA Synthesis Kit	194	RR	p19 siRNA Binding Protein	204
RR	ProtoScript First Strand cDNA Synthesis Kit	194		Magnetic mRNA Isolation Kit	204
RR	Luna Universal One-Step RT-qPCR Kit	71		Oligo d(T) ₂₅ Magnetic Beads	205
RR	Luna Universal Probe One-Step RT-qPCR Kit	71		Oligo d(T) ₂₅ Cellulose Beads	205
RR	One Taq One-Step RT-PCR Kit	72		polyA Spin mRNA Isolation Kit	205
RR	One Tag RT-PCR Kit	72		Streptavidin Magnetic Beads	205
				Monarch Total RNA Miniprep Kit	130
	RNA Ligases & Modifying Enzymes			Monarch RNA Cleanup Kit (10 μg)	131
	Reported Activities for T4 RNA Ligases	195		Monarch RNA Cleanup Kit (50 μg)	131
	RNA Ligase Selection Chart	196		Monarch RNA Cleanup Kit (500 μg)	131
	T4 RNA Ligase 1	197		RNA Markers & Ladders	
	T4 RNA Ligase 2	197		dsRNA Ladder	207
	T4 RNA Ligase 2, truncated	197		microRNA Marker	207
	T4 RNA Ligase 2, truncated K227Q and KQ	198		ssRNA Ladder	207
	RtcB Ligase	198		Low Range ssRNA Ladder	207
	Thermostable 5' AppDNA/RNA Ligase	198		RNA Loading Dye (2X)	207
	5´ DNA Adenylation Kit	199		Universal miRNA Cloning Linker	207
RR	SplintR Ligase	199		Oniversal mintivA cioning Linker	207

Recombinant Enzyme

One or more of these products are covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at gbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications. Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at https://www.neb.com/support/terms-of-sale. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.



Avoiding Ribonuclease Contamination

RNA is more susceptible to degradation than DNA, due to the ability of the 2´ hydroxyl groups to act as nucleophiles. Many ribonucleases (RNases) bypass the need for metal ions by taking advantage of the 2´ hydroxyl group as a reactive species. These RNases are resistant to metal chelating agents and, some of them, like RNase A family enzymes, can survive prolonged boiling or autoclaving. RNase A-type enzymes rely on active site histidine residues for catalytic activity (1) and can be inactivated by the histidine-specific alkylating agent diethyl pyrocarbonate (DEPC).

Sources of RNase Contamination:

RNases are found in all cell types and organisms from prokaryotes to eukaryotes. These enzymes generally have very high specific activity, meaning tiny amounts of contamination in an RNA sample is sufficient to destroy the RNA. The major sources of RNase contamination in a typical laboratory include:

- · Aqueous solutions, reagents used in experiments
- · Environmental exposure RNases are in the air, on most surfaces and on dust
- · Contact with skin

Laboratory Precautions (2,3):

New England Biolabs' enzymes certified for RNA work have been purified to be free of ribonucleases. However, it is possible to reintroduce RNases during the course of experimentation from various sources. RNase contamination can be prevented by following a few common sense laboratory procedures:

- Always wear gloves during an experiment and change them often, especially after contact with skin, hair or
 other potentially RNase-contaminated surfaces, such as doorknobs, keyboards and animals.
- Use RNase-free solutions. Use RNase-free certified, disposable plasticware and filter tips whenever possible.
- · Maintain a separate area for RNA work. Carefully clean the surfaces.
- Decontaminate glassware by baking at 180°C or higher for several hours or by soaking in freshly prepared 0.1% (v/v) DEPC in water or ethanol for 1 hour, followed by draining and autoclaving. Autoclaving will destroy any unreacted DEPC which can otherwise react with other proteins and RNA.
- Decontaminate polycarbonate or polystyrene materials (e.g., electrophoresis tanks) by soaking in 3% hydrogen peroxide for 10 minutes. Remove peroxide by extensively rinsing with RNase-free water prior to use.

Preparation of Solutions (2,3):

Preparation of solutions using the following suggestions can help prevent RNase contamination:

- As an alternative to the historic use of DEPC, which can inhibit enzymatic reactions if not completely removed, we have found that Milli-Q[®] (Millipore) purified water is sufficiently free of RNases for most RNA work.
- DEPC treatment of solutions is accomplished by adding 1 ml DEPC (Sigma) per liter of solution, stirring for 1 hour, and autoclaving for one hour to remove any remaining DEPC. [Note: Compounds with primary amine groups (e.g., Tris) which will react with DEPC, cannot be DEPC-treated. Other compounds, which are not stable during autoclaving, cannot be DEPC-treated].
- Solutions and buffers (e.g. DTT, nucleotides, manganese salts) should be prepared by dissolving the solid (highest available purity) in autoclaved DEPC-treated or Milli-Q water and passing the solution through a 0.22 µm filter to sterilize.

Inhibitors of Ribonucleases:

RNA can also be protected from RNase activity by using one of the following RNase inhibitors:

- RNase Inhibitor, Murine, (NEB #M0314) is a recombinant protein RNase inhibitor of murine origin. Like
 human and porcine RNAse, it is a specific inhibitor of RNases A, B and C, but is more stable due to improved
 resistance to oxidation (4). The inhibitor requires low concentration of DTT (< 1 mM) to maintain activity,
 making it ideal for reactions where low DTT concentration is required (e.g., real-time RT-PCR).
- RNase Inhibitor, Human Placenta, (NEB #M0307), a recombinant protein of human placental origin, is a specific inhibitor for RNases A, B and C. Similar to the Murine RNase Inhibitor, it is compatible with many enzymatic reactions involving RNA (e.g., in vitro transcription, RT-PCR, ligation, etc.).
- Ribonucleoside Vanadyl Complex (NEB #S1402) is a transition-state analog inhibitor of RNase A-type enzymes with K_i = 1 X 10⁻⁵ M. This complex is compatible with many RNA isolation procedures, but it should not be used in the presence of EDTA. The complex also inhibits many other enzymes used in RNA work (5).

References

- (1) Fersht, A.R. (1977) Enzyme Structure and Mechanism Freeman, Reading, PA, 325-329.
- (2) Blumberg, D.D. (1987) Methods Enzymol., 152, 20-24.
- (3) Sambrook, J., Fritsch, E.F. and Maniatis, T. (1989). Molecular Cloning: A Laboratory Manual, (2nd ed.), (pp. 7.3–7.5). Cold Spring Harbor: Cold Spring Harbor Laboratory Press.
- (4) Kim, B.M. et al. (1999) Protein Science, 8, 430-434.
- (5) Berger, S.L. (1987) Methods Enzymol., 152, 227-234



€

€ €

€

₽ €

₽

₽ €

€ €

€

HiScribe[™] T7 High Yield RNA Synthesis Kits

HiScribe T7 High Yield RNA Synthesis Kit #E2040S 50 reactions 223 €

HiScribe T7 Quick High Yield RNA Synthesis Kit #E2050S 50 reactions 273 €

- Synthesis of long and short RNA transcripts
- Incorporation of modified nucleotides
- Incorporation of labeled nucleotides
- Generation of capped RNA using cap
- Synthesis of radioactively labeled probes with high or low specific activity

Description: NEB's HiScribe T7 High Yield RNA Synthesis Kits offer robust in vitro RNA transcription of many kinds of RNA, including internally labeled and co-transcriptionally capped transcripts. Utilizing T7 RNA Polymerase, these kits achieve efficient transcription with small amounts of template, and can generate up to 180 µg per reaction, or up to 30-40 µg of capped RNA using cap analog. RNA generated can be used in a variety of applications, including RNA structure/function studies, ribozyme biochemistry, probes for RNase protection assays and hybridization-based blots, anti-sense RNA and RNAi experiments, microarray analysis, microinjection, and in vitro translation and RNA vaccines.

The HiScribe T7 High Yield RNA Synthesis Kit is an extremely flexible system, with separate NTPs included for flexible reaction setup. The HiScribe T7 Quick High Yield RNA Synthesis Kit is provided in master mix format for fast reaction setup.

RX

The HiScribe T7 High Yield RNA Synthesis Kit Includes:

- T7 RNA Polymerase Mix
- T7 Reaction Buffer (10X)
- ATP, GTP, UTP, CTP (100 mM)
- FLuc Control Template

The HiScribe T7 Quick High Yield RNA Synthesis Kit Includes:

- T7 RNA Polymerase Mix
- NTP Buffer Mix
- FLuc Control Template
- DNase I (RNase-free)
- LiCl Solution

HiScribe SP6 RNA Synthesis Kit

#E2070S 50 reactions 273 €

- Synthesis of long and short RNA transcripts
- Incorporation of modified nucleotides
- Incorporation of labeled nucleotides
- Generation of capped RNA using cap
- Synthesis of radioactively labeled probes with high or low specific activity

Use Monarch RNA Cleanup Kits to purify your synthesized RNA, see page 130-131.

Description: The HiScribe SP6 RNA Synthesis Kit is designed for the in vitro transcription of RNA using SP6 RNA Polymerase. This kit is suitable for synthesis of high yield RNA transcripts and allows for incorporation of cap analogs (not included) or modified nucleotides (not included) to obtain capped, biotin-labeled or dye-labeled RNA. The kit is also capable of synthesizing high specific activity radiolabeled RNA for use as probes or targets.

RNA synthesized from this kit is suitable for many applications including RNA structure and function studies, ribozyme biochemistry, probes for RNase protection or gel shift assays, hybridization-based blots, anti-sense RNA or RNAi experiments, microarray analysis, microinjection and in vitro translation studies.

This kit contains sufficient reagents for 50 reactions of 25 μl each. Each standard reaction yields \ge 80 μg of RNA from 1 µg SP6 Control Template DNA. Each kit can yield ≥ 4 mg of RNA.

RX

The HiScribe SP6 RNA Synthesis Kit Includes:

- SP6 Reaction Buffer (10X)
- ATP (Tris), GTP (Tris), UTP (Tris), CTP (Tris), (50 mM)
- SP6 Control Template
- SP6 RNA Polymerase Mix
- DNase I (RNase-free)
- LiCl Solution

Companion Products:

	1,000 units71 €	Low Range ssRNA Ladder #N0364S 25 gel lanes68
#M0303L RNA Loading Dye (2) #B0363S		E. coli Poly(A) Polymerase #M0276S 100 units72 #M0276L 500 units288
Vaccinia Capping Sys #M2080S	stem 400 units140 €	mRNA Cap 2´-O-Methyltransferase #M0366S 2,000 units57
#M0493S	idelity DNA Polymerase 100 units112 € 500 units508 €	3´-O-Me-m ⁷ G(5´)ppp(5´)G RNA Cap Structure Analog #S1411S 1 µmol148
#T1030L	50 preps 88 € 250 preps 398 €	#\$1411L 5 \(\mu\)mol \(\ldots\)592 G(5')ppp(5')A RNA Cap Structure Analog #\$1406\$ 1 \(\mu\)mol \(\ldots\)138 #\$1406L 5 \(\mu\)mol \(\ldots\)552
#T2030L	10 preps 52 € 100 preps 284 €	G(5')ppp(5')G RNA Cap Structure Analog #S1407S 1 µmol119 #S1407L 5 µmol476
	10 preps 50 € 100 preps 280 €	m ⁷ G(5 ⁻)ppp(5 ⁻)G RNA Cap Structure Analog #S1404S 1 µmol155 #S1404L 5 µmol620
#T2050S	10 preps 58 € 100 preps 440 €	m'G(5')ppp(5')A RNA Cap Structure Analog #S1405S 1 µmol155 #S1405L 5 µmol620

RX

HiScribe T7 ARCA mRNA Kits

HiScribe T7 ARCA mRNA Kit

#E2065S 20 reactions 322 €

HiScribe T7 ARCA mRNA Kit (with tailing) #E2060S 20 reactions 376 €

Companion Products:

DNase I (RNase-Free)

#M0303S 1,000 units 71 € #M0303L 5,000 units 284 €

RNA Loading Dye (2X)

#B0363S 4 ml 49 €

- Synthesis of capped and tailed mRNA
- Incorporation of modified nucleotides
- Template removal and mRNA purification reagents included

Description: Most eukaryotic mRNAs require a 7-methyl guanosine (m7G) cap structure at the 5° end and a Poly(A) tail at the 3´ end to be efficiently translated. By using a DNA template encoding a poly(A) tail, the HiScribe T7 ARCA mRNA Kit can be used to synthesize capped and tailed mRNAs. The cap structure is added to the mRNA by co-transcriptional incorporation of Anti-Reverse Cap Analog (ARCA, NEB #S1411) using T7 RNA Polymerase. The transcription reaction can be set up easily by combining the ARCA/NTP Mix, T7 RNA Polymerase Mix and a suitable DNA template. The kit also allows for partial incorporation of 5mCTP, Pseudo-UTP and other modified nucleotides into mRNA. mRNAs synthesized with the kit can be used for cell transfection, microinjection, in vitro translation and RNA vaccines. Poly(A) tail is incorporated during the transcription reaction. The kit also includes DNase I and LiCI for DNA template removal and quick mRNA purification.

The HiScribe T7 ARCA mRNA Kit (with tailing) is designed for quick production of ARCA capped and polv(A) tailed mRNA in vitro from templates without encoded poly(A) tails.

The HiScribe T7 ARCA mRNA Kit Includes:

- T7 RNA Polymerase Mix
- ARCA/NTP Mix
- DNase I (RNase-free)
- LiCI Solution
- CLuc Control Template

The HiScribe T7 ARCA mRNA Kit (with tailing) Includes:

- T7 RNA Polymerase Mix
- ARCA/NTP Mix
- DNase I (RNase-free)
- E. coli Poly(A) Polymerase
- Poly(A) Polymerase Reaction Buffer
- LiCI Solution
- CLuc Control Template

Advantages:

- · Quick reaction setup and streamlined protocol
- Enables partial incorporation of 5mCTP, Pseudo-UTP and other modified CTP and UTP
- · Ultra-high quality components ensure mRNA integrity

EnGen® sgRNA Synthesis Kit, S. pyogenes

See page 179 for more information.

#E3322S 20 reactions 420 €

The EnGen sgRNA Synthesis Kit, S. pyogenes provides a simple and quick method for transcribing high yields of sgRNA in a single 30 minute reaction, using the supplied reagents and target-specific DNA oligos designed by the user.



NEB's RNA Research Division works to discover, understand and develop enzymes to streamline RNA workflows. Several of its team members are pictured here.

















Recommended HiScribe RNA Synthesis Kits by Application

A A A A A A A A A A A A A A A A A A A	2		T7 KIT	S		SP6 KITS
	APPLICATION	HiScribe T7 High Yield RNA Synthesis Kit (#E2040)	HiScribe T7 Quick High Yield RNA Synthesis Kit (#E2050)	HiScribe T7 ARCA mRNA Kit (#E2065)	HiScribe T7 ARCA mRNA (with Tailing) (#E2060)	HiScribe SP6 RNA Synthesis Kit (#E2070)
	Fluorescent labeling: FAM, Cyanine (Cy) dyes, etc. • Fluorescent <i>in situ</i> hybridization (FISH)		~			~
Probe labeling MRNA & RNA for transfection	Non-fluorescent labeling: Biotin, Digoxigenin In situ hybridization Blot hybridization with secondary detection Microarray		V			~
	High specific activity radiolabeling Blot hybridization RNase protection	~				~
	Streamlined mRNA synthesis with ARCA co-transcriptional capping and enzymatic poly(A) tailing Transfection Microinjection In vitro translation				~	
	Streamlined ARCA capped RNA synthesis Template encoded poly(A) tails Non polyadenylated transcripts Transfection Microinjection In vitro translation			V		
	Co-transcriptional capping with alternate cap analogs Transfection Microinjection In vitro translation		V			~
ti alisietiivii	Post-transcriptional capping with Vaccinia Capping System Transfection Microinjection In vitro translation	V	V			~
	Complete substitution of NTPs: 5-mC, pseudouridine, etc. Post-transcriptional capping with Vaccinia mRNA Capping System	V				~
	Partial substitution of NTPs: 5-mC, pseudouridine, etc.		V	V	~	~
	Unmodified RNA		V			~
	Hairpins, short RNA, dsRNA • Gene knockdown		~			~
	Complete substitution of NTPs • Aptamer selection • Isotopic labeling	~				~
Structure, function, & binding studies	Partial substitution of one or more NTPs • Aptamer selection • Structure determination		V			~
	Unmodified RNA • SELEX • Structure determination		V			~

RNA Polymerases

T3 RNA Polymerase

#M0378S 5,000 units72 €

T7 RNA Polymerase

#M0251S 5,000 units69 € 25,000 units 276 € #M0251L

SP6 RNA Polymerase

#M0207S 2,000 units67 € #M0207L 10,000 units 268 €

Hi-T7 RNA Polymerase

#M0658S 5,000 units83 €

- Radiolabeled RNA probe preparation
- RNA generation for in vitro translation
- RNA generation for studies of RNA structure, processing and catalysis

Description: Initiation of transcription with T3, T7 and SP6 RNA Polymerases is highly specific for the T3, T7 and SP6 phage promoters, respectively. Cloning vectors have been developed which direct transcription from the T3, T7 or SP6 promoter through polylinker cloning sites. These vectors allow in vitro synthesis of

defined RNA transcripts from a cloned DNA sequence.

Hi-T7 RNA Polymerase is an engineered, thermoactive T7 RNA Polymerase. Hi-T7 uses T7 RNA Polymerase Promoters. It increases capping efficiency and eliminates dsRNA by-product formation during synthesis.

R\\\ 37°

Reaction Conditions: 1X RNAPol Reaction Buffer. Supplement with 0.5 mM each ATP, UTP, GTP, CTP (not included) and DNA template containing the appropriate promoter. Incubate at 37°C (T3, T7 and SP6) or 50°C (Hi-T7).

Unit Definition: One unit is defined as the amount of enzyme required to incorporate 1 nmol ATP into an acidinsoluble material in 1 hour at 37°C or 50°C for Hi-T7. Unit assay conditions can be found at www.neb.com.

Concentration: T3 RNA Polymerase: 50,000 units/ ml. T7 RNA Polymerase: 50,000 units/ml. SP6 RNA Polymerase: 20,000 units/ml. Hi-T7 RNA Polymerase: 50,000 units/ml.

E. coli Poly(A) Polymerase

#M0276S 100 units72 € #M0276L 500 units 288 €

Companion Products:

Adenosine-5' Triphosphate (ATP) #P0756S 1 ml 35 € #P0756L 5 ml 138 €

RNase Inhibitor, Murine #M0314S 3,000 units 71 € 15.000 units 284 € #M0314L

- Labeling of RNA with ATP or cordycepin 5´-triphosphate
- Poly(A) tailing of RNA for cloning or affinity purification
- Enhances translation of RNA transferred into eukaryotic cells

Description: E. coli Poly(A) Polymerase catalyzes the template independent addition of AMP from ATP to the 3' end of RNA.

Reaction Conditions: 1X Poly(A) Polymerase Reaction Buffer. Supplement with 1 mM ATP. Incubate at 37°C. May heat inactivate at 65°C for 20 minutes.

Rii 37° 165

Reagents Supplied with Enzyme:

10X Poly(A) Polymerase Reaction Buffer 10 mM ATP

Unit Definition: One unit is defined as the amount of enzyme that will incorporate 1 nmol of AMP into RNA in a 20 µl volume in 10 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 5,000 units/ml

Poly(U) Polymerase

#M0337S 60 units88 €

Companion Products:

RNase Inhibitor, Murine

#M0314S 3,000 units 71 € #M0314L 15,000 units 284 €

Ribonucleotide Solution Set

#N0450S 10 µmol of each 69 € #N0450L 50 μmol of each 277 €

- Labeling of RNA with UTP
- Poly(U) tailing of RNA for cloning

Description: Poly(U) Polymerase catalyzes the template-independent addition of UMP from UTP or AMP from ATP to the 3´ end of RNA.

Reaction Conditions: 1X NEBuffer 2. Supplement with 0.5 mM UTP (not supplied). Incubate at 37°C. May heat inactivate at 65°C for 20 minutes.

R**%** 37° **₩**

Unit Definition: One unit is defined as the amount of enzyme that incorporates 1 nmol of UMP into RNA in a 50 μl volume in 10 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 2,000 units/ml

Notes: Poly(U) Polymerase in NEBuffer 2 will incorporate UMP or AMP from UTP or ATP into RNA. Tailing length of poly(U) varies with UTP and primer concentration. Poly(U) Polymerase is highly processive under low primer concentrations (< 100 pmol).



















RNA REAGENTS

E. coli RNA Polymerase, Core Enzyme & Holoenzyme

E. coli RNA Polymerase, Core Enzyme #M0550S 100 units 178 €

E. coli RNA Polymerase, Holoenzyme #M0551S 50 units 133 €

- RNA synthesis from E. coli promoter
- Transcription initiation studies
- In vitro translation with PURExpress

Description: *E. coli* RNA Polymerase Core Enzyme consists of 5 subunits designated α , α , β ′, β , and ω . The enzyme is free of sigma factor and does not initiate specific transcription from bacterial and phage DNA promoters. The enzyme retains the ability to transcribe RNA from nonspecific initiation sequences. Addition of sigma factors will allow the enzyme to initiate RNA synthesis from specific bacterial and phage promoters. The core enzyme has a molecular weight of approximately 400 kDa.

E. coli RNA Polymerase Holoenzyme is the core enzyme saturated with sigma factor 70. The Holoenzyme initiates RNA synthesis from sigma 70 specific bacterial and phage promoters.

37°

Reaction Conditions: 1X *E. coli* RNA Polymerase Reaction Buffer, 0.5 mM of each rNTP and DNA template. Incubate at 37°C.

Unit Definition: One unit is the amount of enzyme required to incorporate 1 nmol of NTP into RNA in 10 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 1.000 units/ml

Ribonucleotides

Ribonucleotide Solution Set

#N0450S 10 μmol of each 69 € #N0450L 50 μmol of each 277 €

Ribonucleotide Solution Mix

#N0466S 10 μmol of each72 € #N0466L 50 μmol of each288 € **Description:**

Ribonucleotide Solution Set:

The set consists of four separate solutions of ATP, GTP, CTP and UTP, pH 7.5, as sodium salts. Each nucleotide is supplied at 100 mM.

Ribonucleotide Solution Mix:

A buffered equimolar solution of ribonucleotide triphosphates: rATP, rCTP, rGTP and rUTP, pH 7.5, as sodium salts. Each nucleotide is supplied at a concentration of 25 mM (total rNTP concentration equals 100 mM).

Note: To ensure maximum activity upon long-term storage, aliquot and store at -80°C.

Pyrophosphatases

Pyrophosphatase, Inorganic (E. coli)

#M0361S 10 units65 €

#M0361L 50 units260 €

Pyrophosphatase, Inorganic (yeast)

#M2403S 10 units 67 € #M2403L 50 units 268 €

Thermostable Inorganic Pyrophosphatase
#M0296S 250 units 72 €
#M0296L 1,250 units 288 €

- Increasing RNA yield in transcription reactions
- Enhancing DNA replication

Description: Inorganic pyrophosphatase (PPase) catalyzes the hydrolysis of inorganic pyrophosphate to form orthophosphate.

$$P_2 O_7^{-4} + H_2 O \rightarrow 2HPO_4^{-2}$$

Source: Pyrophosphatase, Inorganic (*E. coli*) is prepared from a clone of the *E. coli* inorganic pyrophosphatase gene.

Pyrophosphatase, Inorganic (yeast) is an *E. coli* strain containing a genetic fusion of the *Saccharomyces cerevisiae ppa* gene and the gene coding for *Mycobacterium xenopi* GyrA intein. Developed by BioHelix Corporation, now a wholly owned subsidiary of Quidel Corporation, and produced at New England Biolabs.

RX W

Thermostable Inorganic Pyrophosphatase is an *E. coli* strain carrying a plasmid encoding a pyrophosphatase from the extreme thermophile *Thermococcus litoralis*.

Unit Definition: One unit is defined as the amount of enzyme that will generate 1 µmol of phosphate per minute from inorganic pyrophosphate under standard reaction conditions. Unit assay conditions can be found at www neb com

Concentration: Pyrophosphatase, Inorganic (*E. coli*) and Pyrophosphatase, Inorganic (yeast) 100 units/ml. Thermostable Inorganic Pyrophosphatase 2,000 units/ml.

Thermostable Inorganic Pyrophosphatase retains 100% activity after incubation at 100°C for 4 hours.

Vaccinia Capping System

#M2080S 400 units 140 €

Companion Product:

RNase Inhibitor, Murine

#M0314S 3,000 units 71 € #M0314L 15,000 units 284 €

- Capping mRNA prior to in vivo or in vitro translation
- Labeling 5´end of mRNA

Description: Based on the Vaccinia Virus Capping Enzyme, the Vaccinia Capping System provides the necessary components to add 7-methylguanylate cap structures (Cap 0) to the 5´ end of RNA. In eukaryotes, these terminal cap structures are involved in stabilization, transport and translation of mRNAs. Enzymatic production of capped RNA is an easy way to improve the stability and translational competence of RNA used for in vitro translation, transfection and microinjection. Alternatively, use of labeled GTP in a reaction provides a convenient way to label any RNA containing a 51

This single enzyme is composed of two subunits (D1 and D12) and has three enzymatic activities (RNA triphosphatase and guanylyltransferase by the D1 subunit and guanine methyltransferase by the D12 subunit), all

terminal triphosphate.

R**%** 37°

necessary for addition of a complete Cap 0 structure, m⁷Gppp(5´)N. All capped structures are added in the proper orientation, unlike co-transcriptional addition of some cap analogs.

Reaction Conditions: 1X Capping Buffer. Supplement with 0.5 mM GTP and 0.1 mM SAM. Incubate at 37°C.

Reagents Supplied:

Vaccinia Capping Enzyme Capping Buffer (10X) GTP Solution (10 mM) SAM Solution (32 mM)

Unit Definition: One unit is defined as the amount of enzyme required to incorporate 10 pmol of (α -32P) GTP into an 80 nucleotide transcript in 1 hour at 37°C.

Concentration: 10,000 units/ml

mRNA Cap 2'-O-Methyltransferase

2,000 units57 €

#M0366S

- Enhances translation of RNA
- Improving mRNA expression during microinjection and transfection

Description: mRNA Cap 2´-O-Methyltransferase adds a methyl group at the 2´-O position of the first nucleotide adjacent to the cap structure at the 5' end of the RNA. The enzyme utilizes S-adenosylmethionine (SAM) as a methyl donor to methylate capped RNA (Cap 0) resulting in a Cap 1 structure.

Reaction Conditions: 1X Capping Buffer. Supplement with 0.2 mM SAM (supplied). Incubate at 37°C.

RX 37°

Reagents Supplied:

Capping Buffer (10X) SAM (32 mM)

Unit Definition: One unit is defined as the amount of enzyme required to methylate 10 pmoles of 80 nt long capped RNA transcript in 1 hour at 37°C.

Concentration: 50,000 units/ml

RNA Cap Analog Selection Chart

The 5' terminal m⁷G cap present on most eukaryotic mRNAs promotes translation, *in vitro*, at the initiation level. For most RNAs, the cap structure increases stability, decreases susceptibility to exonuclease degradation, and promotes the formation of mRNA initiation complexes. Certain prokaryotic mRNAs with 5' terminal cap structures are translated as efficiently as eukaryotic mRNA in a eukaryotic cell-free protein synthesizing system. Splicing of certain eukaryotic substrate RNAs has also been observed to require a cap structure.

PRODUCT	APPLICATION	SIZE	PRICE
Anti-Reverse Cap Analog 3'-0-Me-m'G(5') ppp(5')G (#S1411S/L)	 Produces 100% translatable capped transcripts Co-transcriptional capping with T7 (NEB #M0251), SP6 (NEB #M0207) and T3 RNA polymerases Synthesis of m⁷G capped RNA for <i>in vitro</i> splicing assays Synthesis of m⁷G capped RNA for transfection or microinjection 	1/5 µmol	148 €/ 592 €
Standard Cap Analog m ⁷ G(5´)ppp(5´)G (#S1404S/L)	 Co-transcriptional capping with T7, SP6 and T3 RNA polymerases Synthesis of m⁷G capped RNA for <i>in vitro</i> splicing assays Synthesis of m⁷G capped RNA for transfection or microinjection 	1/5 µmol	155 €/ 620 €
Unmethylated Cap Analog G (5´)ppp(5´)G (#S1407S/L)	Co-transcriptional capping with T7, SP6 and T3 RNA polymerases Synthesis of unmethylated G capped RNA	1/5 µmol	119 €/ 476 €
Methylated Cap Analog for A +1 sites m ⁷ G(5´)ppp(5´)A (#S1405S/L)	 Co-transcriptional capping with T7 RNA polymerase from the phi2.5 promoter that contains an A at the transcription initiation site Synthesis of m⁷G capped RNA for <i>in vitro</i> splicing assays Synthesis of m⁷G capped RNA for transfection or microinjection 	1/5 µmol	155 €/ 620 €
Unmethylated Cap Analog for A +1 sites G(5')ppp(5')A (#S1406S/L)	Co-transcriptional capping with T7 RNA polymerase from the phi2.5 promoter that contains an A at the transcription initiation site Synthesis of unmethylated G capped RNA Synthesis of A capped RNA	1/5 µmol	138 €/ 552 €

















3'-Desthiobiotin-GTP & 3'-Biotin-GTP

3´-Desthiobiotin-GTP

#N0761S 0.5 μmol 318 €

NEW

3´-Biotin-GTP

#N0760S 0.5 μmol 318 €

Description: 3´-Desthiobiotin-GTP or 3´-Biotin GTP are guanosine triphosphate (GTP) analogs which are modified at their 3´ position with desthiobiotin or biotin, respectively. When used with the Vaccinia Capping System, (NEB #M2080) these reagents enable affinity tagging of RNA triphosphate ends. Tagged RNAs are enriched by binding to Hydrophilic Streptavidin Magnetic Beads

(NEB #\$1421). Desthiobiotin-tagged RNAs can be eluted with free biotin. This approach is used in Cappable-seq, a method developed at NEB for directly enriching the 5'-ends of primary transcripts (1).

Reference:

(1) Ettwiller, L. et al. (2016) BMC Genomics, 17,199.



#M0463S 4,000 units 250 €

- mRNA decapping, enabling recapping with tagged-GTP analogs
- Biotinylation of 5´ ends of primary transcripts
- Recappable-seq

Description: yDcpS decapping enzyme from *S. cerevisiae* hydrolyzes the phosphodiester bond between the gamma and beta phosphates of m⁷G capped mRNA, leaving behind a diphosphorylated 5' end and m⁷GMP. yDcpS is capable of decapping full length mRNAs and the diphosphorylated 5' end it leaves behind is suitable for recapping using Vaccinia Capping Enzyme.

RX NEBU 37° 📆

Reaction Conditions: 1X yDcpS Reaction Buffer. Incubate at 37°C.

Concentration: 200.000 units/ml

cDNA Synthesis Selection Chart

cDNA SYNTHESIS	FEATURES	SIZE	PRICE
KITS AND MIXES			
NEW LunaScript® RT SuperMix Kit (NEB #E3010)	Ideal for cDNA synthesis in a two-step RT-qPCR workflow Single tube supermix contains random hexamer and oligo-dT primers, dNTPs, Murine RNase Inhibitor, and Luna Reverse Transcriptase Visible blue tracking dye for easy reaction setup Fast 13-minute protocol	25/100 rxns	125 €/420 €
ProtoScript® II First Strand cDNA Synthesis Kit (NEB #E6560)	Generates cDNA at least 10 kb in length Contains ProtoScript II Reverse Transcriptase, an enzyme with increased thermostability and reduced RNase H activity Convenient 2-tube kit includes dNTPs, Oligo-dT primer and Random Primer Mix	30/150 rxns	162 €/648 €
ProtoScript First Strand cDNA Synthesis Kit (NEB #E6300)	Generates cDNA at least 5 kb in length Contains M-MuLV Reverse Transcriptase Convenient 2-tube kit includes dNTPs, Oligo-dT primer and Random Primer Mix	30/150 rxns	155 €/620 €
NEW Template Switching RT Enzyme Mix (NEB #M0466)	Incorporates a universal adaptor sequence at the 3' end of cDNA during the RT reaction Enzyme mix and buffer are optimzed for efficient template switching RT enzyme mix includes RNase Inhibitor High sensitivity for cDNA amplification — enables transcriptome analysis by RNA-seq from single cells or as low as 2 pg of human total RNA Robust and simple workflow for 5' Rapid Amplification of cDNA Ends (RACE) Retains the complete 5' end of transcripts for 2nd Strand cDNA Synthesis	20/100 rxns	88 €/352 €
STANDALONE REAGENTS			
ProtoScript II Reverse Transcriptase (NEB #M0368) An alternative to SuperScript® II	RNase H ⁻ mutant of M-MuLV Reverse Transcriptase with increased thermostability and reduced RNase H activity Increased reaction temperatures (37–50°C)	4,000/10,000/40,000 units	88 €/\$160 €/ 574 €
M-MuLV Reverse Transcriptase (NEB #M0253)	Robust reverse transcriptase for a variety of templates Standard reaction temperatures (37–45°C)	10,000/50,000 units	70 €/280 €
AMV Reverse Transcriptase (NEB #M0277)	Robust reverse transcriptase for a broad temperature range (37–52°C) Can be used for templates requiring higher reaction temperatures	200/1,000 units	72 €/288 €
WarmStart® RTx Reverse Transcriptase (NEB #M0380)	Permits room temperature reaction setup Increased reaction temperatures (50–65°C) Optimized for RT-LAMP isothermal detection	50/250 rxns	65 €/260 €

ProtoScript® II Reverse Transcriptase

#M0368S 4.000 units80 € #M0368L 10,000 units 160 € #M0368X 40,000 units 574 €

Companion Products:

RNase H

#M0297S 250 units 70 € #M0297L 1,250 units 280 €

Monarch Total RNA Miniprep Kit

#T2010S 50 preps 248 €

For a complete listing of Deoxynucleotide Solutions, see page 79.

- Efficient reverse transcription from different starting RNA amounts
- Increased thermostability
- Generates cDNA up to 10 kb or more

RX 42° 65

Description: ProtoScript II Reverse Transcriptase Incubate at 42°C for 50 minutes. If random primers are is a recombinant M-MuLV reverse transcriptase with used, a 10 minute incubation at room temperature is reduced RNase H activity and increased thermostability. recommended before transferring to 42°C. May heat It can be used to synthesize first strand cDNA at higher inactivate at 65°C for 20 minutes. temperatures than the wild type M-MuLV. The enzyme is

Unit Definition: One unit is defined as the amount of enzyme that will incorporate 1 nmol of dTTP into acid-insoluble material in a total reaction volume of 50 µl in 10 minutes at 37°C using poly(rA) • oligo(dT)18 as template. Unit assay conditions can be found at www. neb.com.

Concentration: 200,000 units/ml

RR PCR W WW

LunaScript® RT SuperMix Kit

#E3010S 25 reactions 125 € #E3010L 100 reactions 420 €

LunaScript RT SuperMix Kit is an optimized master mix containing all the necessary components for first strand cDNA synthesis in the context of a two-step RT-qPCR workflow. It features the thermostable Luna Reverse Transcriptase, which supports cDNA synthesis at elevated temperatures. Murine RNase Inhibitor is

active up to 50°C, providing higher specificity, higher

yield of cDNA and more full length cDNA product, up

Reaction Conditions: 1X ProtoScript II Reverse

M-MuLV (RNase H⁻), supplemented with 0.5 mM

Transcriptase Reaction Buffer, 10 mM DTT, 200 units

dNTPs (not included) and 5 µM dT23 VN (not included).

M-MuLV Reverse Transcriptase RNase H-.

to 12 kb in length. This product was formerly known as

also included to protect template RNA from degradation. The LunaScript RT SuperMix Kit contains random hexamer and poly-dT primers, allowing for even coverage across the length of the RNA targets.

See page 70 for more information

M-MuLV Reverse Transcriptase

#M0253S ..70 € 10,000 units #M0253L 50,000 units 280 €

Companion Product:

Monarch Total RNA Miniprep Kit #T2010S 50 preps 248 €

- cDNA synthesis
- RNA Sequencing
- RT-PCR

Description: Molonev Murine Leukemia Virus (M-MuLV) Reverse Transcriptase is an RNA-directed DNA polymerase. This enzyme can synthesize a complementary DNA strand initiating from a primer using either RNA (cDNA synthesis) or single-stranded DNA as a template. M-MuLV Reverse Transcriptase lacks 3 → 5 exonuclease activity.

RN 🙀

Reaction Conditions: 1X M-MuLV Reverse Transcriptase Reaction Buffer. Supplement with dNTPs (not included). Incubate at 37-42°C. May heat inactivate at 65°C for 20 minutes.

Unit Definition: One unit is the amount of enzyme required to incorporate 1 nmol of dTTP into an acidinsoluble material in 10 minutes at 37°C using poly(rA) • oligo(dT) as template primer. Unit assay conditions can be found at www.neb.com.

Concentration: 200.000 units/ml

AMV Reverse Transcriptase

#M0277S 200 units72 € #M0277L 1,000 units 288 €

Companion Product:

Monarch Total RNA Miniprep Kit #T2010S 50 preps 248 €

- cDNA synthesis
- RNA Sequencing
- RT-PCR

Description: Avian Myeloblastosis Virus (AMV) Reverse Transcriptase is an RNA-directed DNA polymerase. This enzyme can synthesize a complementary DNA strand initiating from a primer using RNA (cDNA synthesis) or single-stranded DNA as a template.

Source: Avian Myeloblastosis Virus (AMV)

Reaction Conditions: 1X AMV Reverse Transcriptase Reaction Buffer. Supplement with dNTPs (not included) Incubate at 37-42°C.

No

Unit Definition: One unit is defined as the amount of enzyme required to incorporate 1 nmol of dTTP into an acid-insoluble form in 10 minutes at 37°C using poly(rA) • oligo(dT) as template primer. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 and 25,000 units/ml

Storage Note: Once thawed, store at -20°C. Repeated freeze thaw cycles will inactivate the enzyme. Aliquots can be stored for longer periods at −70°C.



















WarmStart® RTx Reverse Transcriptase

#M0380S 50 reactions65 € #M0380L 250 reactions260 €

Companion Product:

Monarch Total RNA Miniprep Kit #T2010S 50 preps248 €

- RT-LAMP
- cDNA Synthesis
- RT reactions requiring room temperature setup

For more information on products for LAMP, see pages 74.

Description: WarmStart RTx Reverse Transcriptase is a unique *in silico*-designed, RNA-directed DNA polymerase coupled with a reversibly-bound aptamer that inhibits RTx activity below 40°C. This enzyme can synthesize a complementary DNA strand initiating from a primer using RNA (cDNA synthesis) or single-stranded DNA as a template. RTx is a robust enzyme for RNA detection in amplification reactions and is particularly well suited for use in Loop-mediated Isothermal Amplification (LAMP). The WarmStart property enables high-throughput applications, room-temperature setup, and increases the consistency and specificity of amplification reactions. RTx contains intact RNase H activity.

RR W W

Reaction Conditions: 1X Isothermal Amplification Buffer, template, primer, dNTPs and 0.25–0.5 μl of WarmStart RTx Reverse Transcriptase in a reaction volume of 25 μl. Incubate at 50–55°C for cDNA synthesis or directly at 65°C for One-step RT-LAMP. May heat inactivate at 80°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme that will incorporate 1 nmol of dTTP into acid-insoluble material in a total reaction volume of 50 µl in 20 minutes at 50°C using poly(rA) ●oligo(dT)18 as template.

Concentration: 15,000 units/ml

NEW

Template Switching RT Enzyme Mix

#M0466S 20 reactions88 € #M0466L 100 reactions352 €

Companion Products:

NEBNext High-Fidelity 2X PCR Master Mix #M0541S 50 reactions #M0541I 250 reactions 360 € Q5 Hot Start High-Fidelity 2X Master Mix #M0494S 100 rxns (2 x 1.25 ml) 180 € #M0494I 500 rxns (10 x 1.25 ml) 720 € 500 rxns (1 x 12.5 ml) 172 € #M0494X NEB PCR Cloning Kit #E1202S 20 reactions 295 € Monarch Total RNA Miniprep Kit #T2010S 50 preps 248 €

LongAmp Hot Start *Taq* 2X Master Mix
#M0533S 100 reactions 165 €
#M0533L 500 reactions 660 €

- Incorporates a universal adaptor sequence at the 3´ end of cDNA during the RT reaction
- High sensitivity for cDNA amplification enables transcriptome analysis by RNA-seq from single cells or as low as 2 pg of human total RNA
- Robust and simple workflow for 5´-Rapid Amplification of cDNA Ends (RACE)
- Retains the complete 5' end of transcripts for 2nd Strand cDNA Synthesis

Description: Template switching reverse transcription (RT) incorporates a universal adaptor sequence to the 3´-end of cDNA. This convenient feature can be utilized in several downstream applications:


- cDNA synthesis and amplification in a onetube reaction
- 5' Rapid Amplification of cDNA Ends (RACE)
- 2nd strand cDNA synthesis that keeps the 5´ end of the transcripts intact

RR W W

The Template Switching RT Enzyme Mix is optimized for efficient template switching during the RT reaction. The enzyme mix contains RNase Inhibitor in a specially formulated buffer, making reactions easy to setup with no additives needed. It is highly sensitive and specific and can generate RNA-seq libraries from as little as 2 pg of human total RNA or 5´-RACE from 10 ng of total RNA, both with minimal background.

Reagents Supplied:

Template Switching RT Enzyme Mix (10X) Template Switching RT Buffer (4X)



Primers for cDNA Synthesis

Oligo d(N)_n primers are used for the priming and sequencing of mRNA adjacent to the 3´-poly A tail or tailed cDNA. Note: #S1316 does not contain a 5´-phosphate.

Reference: References for properties and applications of these products can be found at www.neb.com.

PRODUCT	NEB #	SIZE	PRICE
Random Primer 6 (5´d(N ₆)3´) ~14.6 nmol	S1230S	1.0 A ₂₆₀ unit	97 €
Random Primer 9 (5´d(N ₉)3´) ~11.6 nmol	S1254S	1.0 A ₂₆₀ unit	97 €
Oligo d(T) ₂₃ VN	S1327S	1.0 A ₂₆₀ unit	97 €
Random Primer Mix	S1330S	100 μΙ (60 μΜ)	97 €
Oligo d(T) ₁₈ mRNA Primer	S1316S	5.0 A ₂₆₀ units	91 €

ProtoScript II First Strand cDNA Synthesis Kit

#E6560S 30 reactions 162 € #E6560L 150 reactions 648 €

Companion Products:

polyA Spin mRNA Is	solation Kit
#S1560S	8 isolations224 €
Magnetic mRNA Iso	lation Kit
#S1550S	25 isolations 321 €
RNase Inhibitor, Mu	rine
#M0314S	3,000 units 71 €
#M0314L	15,000 units 284 €

Monarch Total RNA Miniprep Kit #T2010S 50 preps 248 €

Random Primer Mix

100 μI (60 μM)97 € #S1330S

Oligo d(T)23VN*

#S1327S 1.0 A₂₆₀ unit 97 €

*Note: V = A, G or C, and N = A, G, C or T

- Enzyme and reaction mixes add flexibility to reaction setup
- Suitable for any PCR format
- Efficient reverse transcription from different starting RNA amounts
- Generates cDNA at least 10 kb

Description: ProtoScript II First Strand cDNA Synthesis Kit features two optimized mixes. ProtoScript II Enzyme Mix and ProtoScript II Reaction Mix. The enzyme mix combines ProtoScript II Reverse Transcriptase and Murine RNase Inhibitor, while the reaction mix contains dNTPs and an optimized buffer. ProtoScript II Reverse Transcriptase is a recombinant M-MuLV reverse transcriptase with reduced RNase H activity and increased thermostability. It can be used to synthesize first strand cDNA at higher temperatures than the wild type M-MuLV. The enzyme is active up to 50°C, providing higher specificity and higher yield of cDNA.

The kit also provides two optimized primers for reverse transcription and nuclease-free water. An anchored Oligo-dT primer $[d(T)_{23}VN]$ forces the primer to anneal to the beginning of the poly(A) tail. The optimized Random Primer Mix provides random and consistent priming sites covering the entire RNA templates including both mRNAs and non-polyadenylated RNAs. The first strand cDNA product generated is up to 10 kb.

RX

RX

The ProtoScript II First Strand cDNA Synthesis Kit Includes:

- 10X ProtoScript II Enzyme Mix
- 2X ProtoScript II Reaction Mix
- Random Primer Mix (60 μM), Oligo d(T)₂₃VN Primer (50 µM)** and Nuclease-free Water

For robust amplification of a wide range of DNA templates, we recommend One Taq® or Q5® High-Fidelity DNA Polymerases.

ProtoScript First Strand cDNA Synthesis Kit

#E6300S 30 reactions 155 € #E6300L 150 reactions 620 €

Companion Products:


polyA Spin mRNA Isola	ation Kit						
#S1560S	8 isolations224 €						
Magnetic mRNA Isolati	on Kit						
#S1550S	25 isolations321 €						
Monarch Total RNA Mi	Monarch Total RNA Miniprep Kit						
#T2010S	50 preps248 €						

- Enzyme and reaction mixes add flexibility to reaction setup
- Suitable for any PCR format
- Efficient reverse transcription from different starting RNA amounts
- Generates cDNA at least 5 kb

Description: ProtoScript First Strand cDNA Synthesis Kit features two optimized mixes, ProtoScript Enzyme Mix and Protoscript Reaction Mix. Protoscript Enzyme Mix combines M-MuLV Reverse Transcriptase and RNase Inhibitor, Murine, while Protoscript Reaction Mix contains dNTPs and an optimized buffer. The kit also contains two optimized primers for reverse transcription and nuclease-free water. An anchored oligo-dT primer [d(T)₂₂VN] forces the primer to anneal to the beginning of the polyA tail. The optimized Random Primer Mix provides random and consistent priming sites covering the entire RNA template including both mRNAs and nonpolyadenylated RNAs. The first strand cDNA product generated is over 13.0 kb. This product was formerly known as M-MuLV First Strand cDNA Synthesis Kit.

The ProtoScript First Strand cDNA Synthesis Kit Includes:

- 10X ProtoScript Enzyme Mix
- 2X ProtoScript Reaction Mix
- Random Primer Mix (60 μM), Oligo d(T)₂₂VN Primer (50 μM)** and Nuclease-free Water



First Strand cDNA Synthesis with the ProtoScript Kit.

Reactions were performed at 42°C using 2 µg of human spleen total RNA. Negative control reactions (-RT) were carried out with 1X ProtoScript Reaction Mix. A fraction of the first strand cDNA product was used to amplify sequences specific for three different messenger RNAs using 1X LongAmp® Taq 2X Master Mix (NEB #M0287). Lane 1: 1.1 kb of beta-actin gene. Lane 2: RT control of 1.1 kb of beta-actin gene. Lane 3: 4.7 kb of Xrn-1 gene. Lane 4: -RT control of 4.7 kb of Xrn-1 gene. Lane 5: 9.8 kb of guanine nucleotide exchange factor p532. Lane 6: -RT control of 9.8 kb of quanine nucleotide exchange factor p532. Marker M is 1 kb Plus DNA Ladder (NEB #N3200)











Enzymes for Innovation

37° Incubation Temperature



^{**}Oligo d(T)23 VN and Random Primer Mix contain 1 mM dNTP

^{**}Oligo d(T)22 VN and Random Primer Mix contain 1 mM dNTP

RT-PCR & RT-qPCR Kits

Luna® Universal (One-Step RT-qPCR Kit
#E3005S	200 reactions 218 €
#E3005L	500 reactions498 €
#E3005X	1,000 reactions869 €
#E3005E	2,500 reactions 1924 €
Luna Universal Pr	obe One-Step RT-qPCR Kit
#E3006S	200 reactions198 €
#E3006L	500 reactions448 €
#E3006X	1,000 reactions788 €

2,500 reactions 1738 €

#E3006E

See pages 71-72 for more information. One Tag® One-Step RT-PCR Kit

#E5315S 30 reactions 154 € One Taq RT-PCR Kit #E5310S 30 reactions 143 €

RX

RNA Ligase Activity Chart

NEB offers a variety of ligases for DNA and RNA research. Many of these enzymes are recombinant and all offer the quality and value you have come to expect from our products. The chart below highlights reported activities of our T4 ligases ranked by application. A substrate-based selection chart for DNA ligases can be found on page 94.

Reported Activities and Applications for T4 Ligases

T4 RNA Ligase 1 ①	T4 RNA Ligase 2 ① **********************************	T4 RNA Ligase 2 Truncated 14 RNA Ligase 2 Truncated K22TQ 14 RNA Ligase 2 Truncated KQ 1
T4 DNA Ligase ①H P ②	Thermostable 5' App DNA/RNA Ligase 1 OH App 2 OH App 3 OH App	T3 DNA Ligase ① OH P
SplintR Ligase ① OH P ② OH P	RtcB Ligase ① ————— P OH ② ~~~~~ P OH ② ———————————————————————————————————	5' Adenylation Kit 1 5' PSS DNA — AppSS DNA 2 5' PRNA — App RNA 3 SSDNA 3' p — SSDNA 3' ppA 4 RNA 3' p — RNA 2',3' cyclic p

The ligation activities depicted have been reported, but may require optimized reaction conditions.

RNA Ligase Selection Chart

RNA Ligase Selection Chart

	T4 RNA Ligase 1	T4 RNA Ligase 2	T4 RNA Ligase 2 Truncated	T4 RNA Ligase 2, Truncated K227Q	T4 RNA Ligase 2, Truncated KQ	Thermostable 5´ App DNA/RNA Ligase	5´ Adeny- lation Kit	SplintR® Ligase	RtcB Ligase
RNA APPLICATIONS									
Ligation of nicks in dsRNA		***							
Labeling of 3´ termini of RNA	***		*	*	*	*			
Ligation of ssRNA to ssRNA	***								
Ligation of preadenylated adaptors to RNA	**		**	**	***	**			
5´ Adenylation	•						***		
Ligation of 3´P and 5´OH of ssRNA									***
DNA APPLICATIONS									
Ligation of preadenylated adaptors to ssDNA						***			
DNA/RNA APPLICATIONS									
Joining of RNA & DNA in a ds-structure		**							
ssDNA Ligation with RNA Splint		**						***	
Ligation of RNA and DNA with 3´P and 5´OH									**
NGS APPLICATIONS									
NGS Library Prep ssRNA-ssDNA (ligation)	A		A	A	A				
NGS Library Prep ssRNA-ds-Adaptor splinted ligation		A							
FEATURES									
Thermostable						~	V		
Recombinant	V	V	~	V	V	V	V	~	~

- Optimal, recommended ligase for selected application
 - Works well for selected application
 - Will perform selected application, but is not recommended
 - ▲ Please consult the specific NGS protocol to determine the optimal enzyme for your needs



Harriet is the Chief IP Counsel for NEB and has been with the company for 17 years. When she is not working, Harriet enjoys traveling with a good history book.

















196

T4 RNA Ligase 1 (ssRNA Ligase)

T4 RNA Ligase 1 (ssRNA Ligase)

#M0204S 1,000 units63 € #M0204L 5,000 units252 €

T4 RNA Ligase 1 (ssRNA Ligase)

High Concentration

#M0437M 5,000 units 253 €

Companion Products:

Adenosine 5´ -Triphosphate (ATP) #P0756S 1 ml 35 € #P0756L 5 ml 138 €

Universal miRNA Cloning Linker

#S1315S 5 μg148 €

- Ligation of ss-RNA and DNA
- Labeling of 3´-termini of RNA with 5´-[³²P] pCp
- Inter- and intramolecular joining of RNA and DNA molecules
- Synthesis of ss-oligodeoxyribonucleotides
- Incorporation of unnatural amino acids into proteins

Description: Catalyzes ligation of a 5′ phosphoryl-terminated nucleic acid donor to a 3′ hydroxyl-terminated nucleic acid acceptor through the formation of a 3′→5′ phosphodiester bond with hydrolysis of ATP to AMP and PPi. Substrates include single-stranded RNA and DNA as well as dinucleoside pyrophos-

Reaction Conditions: 1X T4 RNA Ligase Reaction Buffer. Supplement with 1 mM ATP. Incubate at 25°C. May heat inactivate at 65°C for 15 minutes or boiling for 2 minutes.

RN 25° 1664

Notes on Use: Addition of DMSO to 10% (v/v) is required for pCp ligation.

Reagents Supplied with Enzyme:

10X T4 RNA Ligase Reaction Buffer

10 mM ATP (with NEB #M0204) or 100 mM ATP (with NEB #M0437) 50% PEG 8000

Unit Definition: One unit is defined as the amount of enzyme required to convert 1 nmol of 5´-[³²P] rA16 into a phosphate resistant form in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 or 30,000 units/ml

T4 RNA Ligase 2 (dsRNA Ligase)

#M0239S 150 units78 € #M0239L 750 units312 €

- Cohesive-end adaptor ligation
- Best choice for ligating nicks in dsRNA
- Suitable for ligating 3´ OH of RNA to 5´ phosphate of DNA in a DNA/ RNA hybrid

Description: T4 RNA Ligase 2, also known as T4 Rnl2 (gp24.1), has both intermolecular and intramolecular RNA strand-joining activity. Unlike T4 RNA Ligase 1 (NEB #M0204), T4 RNA Ligase 2 is much more active joining nicks on double stranded RNA than on joining the ends of single stranded RNA. The enzyme requires an adjacent 5′ phosphate and 3′ OH for ligation. The enzyme can also ligate the 3′ OH of RNA to the 5′ phosphate of DNA in a double stranded structure.

RX 37° 🐝

Reaction Conditions: 1X T4 RNA Ligase 2 Reaction Buffer. Incubate at 37°C. May heat inactivate at 80°C for 5 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to ligate 0.4 µg of an equimolar mix of a 23-mer and 17-mer RNAs in a total reaction volume of 20 µl in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

T4 RNA Ligase 2, truncated

#M0242S 2,000 units69 € #M0242L 10,000 units276 €

Companion Product:

Universal miRNA Cloning Linker #S1315S 5 µg148 €

- Ligate a pre-adenylated DNA or RNA sequence tag to any RNA 3 end
- Join a single-stranded adenylated primer to small RNAs for cDNA library creation
- Join a single-stranded adenylated primer to RNA for strand-specific cDNA library construction

Description: T4 RNA Ligase 2, truncated (T4 Rnl2tr) specifically ligates the pre-adenylated 5' end of DNA or RNA to the 3' end of RNA. The enzyme does not require ATP, but does need the pre-adenylated substrate. T4 Rnl2tr is expressed from a plasmid in *E. coli* which encodes the first 249 amino acids of the full-length T4 RNA Ligase 2. Unlike the full-length ligase, T4 Rnl2tr cannot ligate the phosphorylated 5' end of RNA or DNA to the 3' end of RNA. This enzyme, also known as Rnl2 (1–249), has been used for optimized linker ligation for the cloning of microRNAs. This enzyme reduces background ligation, because it can only use pre-adenylated linkers.

RN 25° (65)

Reaction Conditions: 1X T4 RNA Ligase Reaction Buffer. Incubate at 25°C. May heat inactivate at 65°C for 20 minutes.

Reagents Supplied with Enzyme:

10X T4 RNA Ligase Reaction Buffer 50% PEG 8000

Unit Definition: 200 units is defined as the amount of enzyme required to give 80% ligation of a 31-mer RNA to the pre-adenylated end of a 17-mer [DNA Universal miRNA Cloning Linker (NEB #S1315)] in a total reaction volume of 20 µl in 1 hour at 25°C. Unit assay conditions can be found at www.neb.com.

Concentration: 200,000 units/ml

T4 RNA Ligase 2, truncated K227Q and truncated KQ

T4 RNA Ligase 2, truncated K227Q #M0351S 2,000 units 69 € #M0351L 10,000 units 276 €

T4 RNA Ligase 2, truncated KQ

#M0373S 2,000 units 69 € #M0373L 10,000 units 276 €

Companion Product:

Universal miRNA Cloning Linker #S1315S 5 μg 148 \in

- Ligate a pre-adenylated DNA or RNA sequence tag to any RNA 3´ end
- Join a single stranded adenylated primer to small RNAs for cDNA library creation
- Join a single stranded adenylated primer to RNA for strand-specific cDNA library construction

Description: T4 RNA Ligase 2, truncated KQ (T4 Rnl2tr KQ) specifically ligates the pre-adenylated 5´ end of DNA or RNA to the 3´ OH end of RNA. The enzyme does not use ATP for ligation, but requires pre-adenylated linkers.

T4 Rnl2tr KQ is a double-point mutant of T4 RNA Ligase 2, truncated (NEB #M0242). Mutation of K227 in T4 RNA Ligase 2 reduces enzyme lysyl adenylation. K227Q reduces the formation of undesired ligation products (concatemers and circles) by T4 Rnl2tr, by reducing the trace activity of T4 Rnl2tr in transfer of adenylyl groups from linkers to the 5´-phosphates of input RNAs. Mutation of R55K in T4 Rnl2tr K227Q increases the ligation activity of the enzyme to levels similar to T4 Rnl2tr.

The exclusion of ATP, use of pre-adenylated linkers, and the reduced enzyme lysyl adenylation activity provide the lowest possible background in ligation reactions. This enzyme has been used for optimized linker ligation for high-throughput sequencing library construction of small RNAs.

R 25° **₩**

Source: Expressed as an MBP fusion from a plasmid in *E. coli* which encodes the first 249 amino acids of the full length T4 RNA Ligase 2. T4 Rnl2tr K227Q has a lysine to glutamine mutation at position 227. T4 Rnl2tr KQ has an arginine to lysine and lysine to glutamine mutation at positions 55 and 227, respectively.

Reaction Conditions: 1X T4 RNA Ligase Reaction Buffer. Incubate at 25°C. May heat inactivate at 65°C for 20 minutes.

Reagents Supplied with Enzyme:

10X T4 RNA Ligase Reaction Buffer 50% PEG 8000

Unit Definition: 200 units is defined as the amount of enzyme required to give 80% ligation of a 31-mer RNA to the pre-adenylated end of a 17-mer DNA [Universal miRNA Cloning Linker (NEB #S1315)] in a total reaction volume of 20 µl in 1 hour at 25°C. Unit assay conditions can be found at www.neb.com.

Concentration: 200,000 units/ml

RtcB Ligase

#M0458S 25 reactions73 €

- Ligate ssRNA or ssDNA with a 3´-phosphate or a 2´,3´-cyclic phosphate to the 5´-OH of ssRNA
- Circularization of ssRNA with compatible ends

This is an **Enzyme for Innovation** (EFI). EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit **www.neb.com/EnzymesforInnovation** to view the full list.

Description: RtcB Ligase from *E. coli* joins single stranded RNA with a 3´-phospate or 2´,3´-cyclic phosphate to another RNA with a 5´-hydroxyl. Ligation requires both GTP and MnCl₂ and proceeds through a 3´-guanylate intermediate. With substrates having a 2´,3´-cyclic phosphate end, hydrolysis to a 3´-phosphate precedes 3´ end activation with GMP and ligation.

Source: RtcB Ligase is expressed as His-tagged fusion in *E. coli*.

RN 😻 37°

Reaction Conditions: 1X RtcB Reaction Buffer. Supplement with 0.1 mM GTP and 1 mM MnCl₂ (supplied). Incubate at 37°C.

Reagents Supplied with Enzyme:

10X RtcB Reaction Buffer MnCl₂ (10 mM) GTP (10 mM)

Concentration: 15 µM

Thermostable 5' AppDNA/RNA Ligase

#M0319S 10 reactions70 € #M0319L 50 reactions280 €

Companion Product:

Universal miRNA Cloning Linker #\$1315\$ 5 μ g 148 \in

- Ligation of ssDNA to an adenylated DNA linker for NGS library construction
- Ligation of an adenylated linker to RNA at elevated temperatures for small RNA NGS library construction

Description: Thermostable 5´ App DNA/RNA Ligase is a point mutant of catalytic lysine of RNA ligase from *Methanobacterium thermoautotrophicum*. This enzyme is ATP independent. It requires a 5´ pre-adenylated linker for ligation to the 3´-OH end of either RNA or single stranded DNA (ssDNA). The enzyme is also active in ligation of RNA with 2´-O-methylated 3´ end to 5´-adenylated linkers. The optimal temperature for ligation reaction is 60–65°C. The mutant ligase is unable to adenylate the 5´-phosphate of RNA or ssDNA, which reduces the formation of undesired ligation products (concatemers and circles).

RX 65° W

The ability of the ligase to function at 65°C might reduce the constraints of RNA secondary structure in RNA ligation experiments.

Reaction Conditions: 1X NEBuffer 1. Incubate at 65°C.

Reagents Supplied with Enzyme:

10X NEBuffer 1 10X MnCl₂

Concentration: 20 µM

Usage Note: For optimal ligation of ssDNA to preadenylated linkers, we recommend using NEBuffer 1 supplemented with manganese (supplied).



















NA KEAGENI

5 DNA Adenylation Kit

#E2610S 10 reactions 116 € #E2610L 50 reactions 464 €

- Enzymatic 5' adenylation of ss-DNA linkers for next gen sequencing
- One-step reaction gives quantitative adenylation
- Simpler than existing chemical and enzymatic methods
- Reduces need for purification of reaction product

Description: The 5´ DNA Adenylation Kit is a simple and efficient enzymatic method for generating 5´-adenylated DNA oligonucleotides using *Mth* RNA ligase, ATP and single stranded 5´-phosphorylated DNA. The kit is optimized to produce the adenylated DNA intermediate with or without 3´ terminator. The 5´ DNA adenylation kit routinely generates greater than 95% conversion of pDNA to AppDNA. This highly efficient process eliminates the need for gel isolation of the product and increases overall yield.

R\\\ 65°

The 5' DNA Adenylation Kit Includes:

- Mth RNA Ligase (Recombinant)
- 5´ DNA Adenylation Reaction Buffer (10X)
- 1 mM ATP

Note: The low turnover of the enzyme requires an approximately equimolar concentration of the enzyme and the oligonucleotide substrate. Adenylated DNA linkers can be used for 3´-end ligation of RNA in cDNA library preparation for next generation sequencing protocols.

SplintR® Ligase

#M0375S 1,250 units95 € #M0375L 6,250 units428 €

- Ligation of adjacent, single-stranded DNA splinted by a complimentary RNA
- Characterization of miRNAs and mRNAs, including SNPs

Description: SplintR Ligase, also known as PBCV-1 DNA Ligase or *Chlorella* virus DNA Ligase, efficiently catalyzes the ligation of adjacent, single-stranded DNA oligonucleotides splinted by a complementary RNA strand. This previously unreported activity may enable novel approaches for characterization of miRNAs and mRNAs, including SNPs. SplintR is ideally suited for many target enrichment workflows with applications in next-generation sequencing and molecular diagnostics. The robust activity of the enzyme and its affinity for RNA-splinted DNA substrates (apparent Km = 1 nM) enable sub-nanomolar detection of unique RNA species within a complex mixture, making SplintR ligase a superior choice for demanding RNA detection technologies.

RR 25° KS

Reaction Conditions: 1X SplintR Ligase Reaction Buffer. Incubate at 25°C. May heat inactivate at 65°C for 20 minutes

Unit Definition: One unit is defined as the amount of enzyme needed to ligate (to 50% completion) 2 picomoles of a tripartite FAM-labeled DNA:RNA hybrid substrate in a 20 µl reaction at 25°C in 15 minutes in 1X SplintR Ligase Reaction Buffer. Unit assay conditions can be found at www.neb.com.

Concentration: 25,000 units/ml

RNA 5' Pyrophosphohydrolase (RppH)

#M0356S

200 units92 €

- Conversion of 5´-triphosphate RNA to monophosphate RNA
- Preparation of 5´-phosphate RNA for ligation
- Characterization of RNA 5´ ends

Description: The bacterial RNA 5´ Pyrophosphohydrolase (RppH) removes pyrophosphate from the 5´ end of triphosphorylated RNA to leave a 5´ monophosphate RNA. The RppH protein was also known as NudH/YgdP which can split diadenosine penta-phosphate to ADP and ATP.

Reaction Conditions: 1X NEBuffer 2. Incubate at 37°C.

R**%** 37°

10X NEBuffer 2

Unit Definition: One unit is defined as the amount of enzyme that converts 1 µg 300-mer RNA transcript into an XRN-1 digestible RNA in 30 minutes at 37°C.

Concentration: 5,000 units/ml

Reagents Supplied with Enzyme:

5' Deadenylase

#M0331S

1,000 units70 €

- Deadenylation of 5´ end of DNA and RNA
- Aprataxin-dependent DNA repair assay
- Analysis of dinucleoside tetraphosphate

Description: Yeast 5´ Deadenylase is a member of the HIT (histidine triad) family of proteins and specifically a member of the Hint branch. It is the yeast orthologue of aprataxin. Mutations in human aprataxin have been known to be involved in the neurological disorder known as ataxia oculomotor apraxia-1. The human protein has been shown to resolve abortive ligation intermediates by removing the AMP at the 5´ end of DNA (AMP-DNA hydrolase activity). It also repairs DNA damage at 3´ ends by removing 3´-phosphate and 3´-phosphoglycolate. Human aprataxin acts on small molecules, such as nucleotide polyphosphates diadenosine tetraphosphate (AppppA) and lysyl-AMP.

RX 30° ₩

The 5´ Deadenylase is encoded by the *HNT3* gene of *S. cerevisiae*. NEB has shown this protein is capable of deadenylation from 5´ end of DNA and RNA, leaving the phosphate at 5´ end. It also cleaves AppppA into ATP and AMP. Its activity on lysyl-AMP is not detectable.

Reaction Conditions: 1X NEBuffer 1 and 5–50 pmol adenylated DNA (AMP-DNA) in 20 μ l. Incubate at 30°C. May heat inactivate at 70°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to remove 10 pmoles of AMP from a 5'-adenylated DNA oligo in 10 minutes at 30°C.

Concentration: 50,000 units/ml

RNase I_f

#M0243S 5,000 units66 € #M0243L 25,000 units 264 €

- Eliminates RNA from DNA and protein preparations
- Degradation of single-stranded RNA to mono-, di- and trinucleotides
- Used in ribonuclease protection assays

Description: Ribonuclease I. (RNase I.) is an RNA endonuclease which will cleave at all RNA dinucleotide bonds leaving a 5' hydroxyl and 2', 3' cyclic monophosphate. It has a preference for single-stranded RNA over double-stranded RNA. RNase I, is a recombinant protein fusion of RNase I (from E. coli) and maltose-

Reaction Conditions: 1X NEBuffer 3. Incubate at 37°C. May heat inactivate at 70°C for 20 minutes.

binding protein. It has identical activity to RNase I.

R**%** 37° ₩

Notes on Use: RNase I, will not degrade DNA. It has a strong preference for single-stranded RNA over doublestranded RNA.

Unit Definition: One unit is defined as the amount of enzyme required to fully digest 1 picomole of synthetic ssRNA 33-mer in a total reaction volume of 10 µl in 15 minutes in 1X NEBuffer 3 as visualized on a 20% acrylamide gel (40:1 Bis) stained with SYBR® Gold. Unit assay conditions can be found at www.neb.com.

Concentration: 50.000 units/ml SYBR® is a registered trademark of Molecular Probes, Inc.

RNase H

#M0297S 250 units70 € 1,250 units 280 € #M0297L

- Removal of poly(A) tails of mRNA hybridized to poly(dT)
- Removal of mRNA during second strand cDNA synthesis

Description: Ribonuclease H (RNase H) is an endoribonuclease which specifically hydrolyzes the phosphodiester bonds of RNA which is hybridized to DNA. This enzyme does not digest single- or doublestranded DNA.

Reaction Conditions: 1X RNase H Reaction Buffer. Incubate at 37°C. May heat inactivate at 65°C for 20 minutes.

R**%** 37° **₩**

Unit Definition: One unit is defined as the amount of enzyme required to produce 1 nmol of ribonucleotides from 20 picomoles of a fluorescently labelled 50 base pair RNA-DNA hybrid in a total reaction volume of 50 µl in 20 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 5,000 units/ml

Thermostable RNase H

#M0523S 250 units 140 €

- Higher stringency RNA structure mapping and site-specific RNA cleavage
- Removal of poly(A) tails from mRNA hybridized to oligo(dT)
- Removal of mRNA during second strand cDNA synthesis
- Component of isothermal amplification methods

Description: Thermostable RNase H specifically recognizes and cleaves the phosphodiester bonds of an RNA strand in an RNA-DNA hybrid while leaving the DNA strand intact. This thermostable nuclease exhibits the same enzymatic properties as E. coli RNase H, but is active at much higher temperatures.

Reaction Conditions: 1X RNase H Reaction Buffer. Incubate at 50°C.

R% 50° 1464

Unit Definition: One unit is defined as the amount of enzyme required to produce 1 nmol of ribonucleotides from 40 picomoles of a fluorescently labelled 25 base pair RNA-DNA hybrid in a total reaction volume of 50 µl in 20 minutes at 50°C. Unit assay conditions can be found at www.neb.com.

Concentration: 5.000 units/ml

RNase HII

#M0288S 250 units72 € 1,250 units 288 € #M0288L

- Nicking of products generated with a polymerase that will incorporate ribonucleotides
- Generation of a double-stranded break at the site of an incorporated ribonucleotide when used with T7 Endo I
- Degradation of the RNA portion of Okazaki fragments or other RNA-DNA hybrids

Description: Ribonuclease HII (RNase HII) is an endoribonuclease that preferentially nicks 5° to a ribonucleotide within the context of a DNA duplex. The enzyme leaves 5´ phosphate and 3´ hydroxyl ends. RNase HII will also nick at multiple sites along the RNA portion of an Okazaki fragment.

Source: An *E. coli* strain containing a genetic fusion of the RNase HII gene (rnhB) from E. coli and the gene coding for maltose binding protein (MBP). Following affinity chromatography, RNase HII is cleaved from the fusion construct by Factor Xa and then purified away from both MBP and Factor Xa. RNase HII cleaved from MBP has four additional amino acids at its N-terminus (Ile-Ser-Glu-Phe).

RX 37° **166**

Reaction Conditions: 1X ThermoPol Reaction Buffer. Incubate at 37°C.

Unit Definition: One unit is defined as the amount of enzyme required to yield a fluorescence signal consistent with the nicking of 100 pmol of synthetic double-stranded DNA substrate containing a single ribonucleotide near the quencher of a fluorophore/quencher pair in 30 minutes at 37°C in 1X ThermoPol Reaction Buffer. Unit assay conditions can be found at www.neb.com.

Concentration: 5,000 units/ml

Note: Incubation with 0.1% SDS is sufficient to inactivate RNase HII.

















RNA REAGENT

Phosphorylation and Dephosphorylation

See pages 98-100 for more information.

Quick Dephosphorylation Kit

#M0508S 100 units 68 € #M0508L 500 units 272 €

Antarctic Phosphatase

#M0289S 1,000 units 68 € #M0289L 5,000 units 272 €

Alkaline Phosphatase Calf Intestinal (CIP)

#M0290S 1,000 units 72 €

#M0290L 5,000 units 288 €

Shrimp Alkaline Phosphatase (rSAP)

#M0371S 500 units 60 € #M0371L 2,500 units 240 €

T4 Polynucleotide Kinase

#M0201S 500 units 58 € #M0201L 2,500 units 232 €

ShortCut® RNase III

#M0245S 200 units 123 € #M0245L 1,000 units 492 €

- Generates siRNAs for RNA interference studies
- Gene silencing
- Target validation
- Removal of long dsRNAs

Description: ShortCut RNase III converts long double-stranded RNA into a heterogeneous mix of short (18–25 bp) interfering RNAs (siRNA) suitable for RNA interference in mammalian cells. 1.5 units (1 μ I) of ShortCut RNase III is sufficient to convert 1 μ g of dsRNA into siRNA.

Source: An *E. coli* strain containing a genetic fusion of the *E. coli* RNase III gene (rnc) and the gene coding for maltose binding protein (MBP).

Reaction Conditions: 1X ShortCut RNase III Reaction Buffer. Supplement with 20 mM MnCl₂ (supplied). Incubate at 37°C.


Unit Definition: One unit is the amount of enzyme required to digest 1 μ g of dsRNA to siRNA in a total reaction volume of 50 μ l in 20 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 2,000 units/ml

RN 37° W

Advantages of the ShortCut RNase III:

- Make effective siRNAs against any gene target
- Heterogeneous population of siRNA ensures silencing of target gene
- From DNA template to transfection in just 1 day
- Eliminates trial and error approach of synthetic siRNA



siRNA production by ShortCut RNase III: A) Varying amounts of ShortCut RNase III were incubated with 2 μg of a 500 bp dsRNA for 20 minutes. B) dsRNA fragments (1 kb and 175 bp) were digested with ShortCut RNase III. Digests were analyzed by 20% TBE polyacrylamide gel electrophoresis.

XRN-1

#M0338S 20 units88 € #M0338L 100 units352 €

 Removal of RNA containing 5' monophosphate from an RNA mixture

Description: XRN-1 is a highly processive 5′ to 3′ exoribonuclease, requiring a 5′ monophosphate. It also acts on 5′ monophosphate ssDNA with reduced efficiency.

Reaction Conditions: 1X NEBuffer 3. Incubate at 37°C. May heat inactivate at 70°C for 10 minutes.

R% 37° ₩

Unit Definition: One unit is defined as the amount of enzyme that digests 1 μg of phosphorylated yeast RNA in 60 minutes at 37°C.

Concentration: 1,000 units/ml

37°

Exonuclease T

#M0265S 250 units72 € 1,250 units 288 € #M0265L

Generate blunt ends in DNA or RNA with 3' overhangs

Description: Exonuclease T (Exo T), also known as RNase T, is a single-stranded RNA or DNA specific nuclease that requires a free 3´ terminus and removes nucleotides in the $3' \rightarrow 5'$ direction. Exo T can be used to generate blunt ends from RNA or DNA having 3' extensions.

Source: Exonuclease T is overexpressed and purified as a C-terminal fusion to maltose-binding protein (MBP). MBP is removed from Exo T by Factor Xa cleavage and then purified. When cleaved from MBP has an additional amino acid on the N-terminus and a Phe instead of a Met (i.e. Glu-Phe-Exo T instead of Met-Exo T).

RR 25° 65

Reaction Conditions: 1X NEBuffer 4. Incubate at 25°C. May heat inactivate at 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to produce 0.1 nmol of TCA soluble nucleotides from 1 nmol of [3H]-labeled polythymidine in a total reaction volume of 100 µl in 30 minutes at 25°C in 1X NEBuffer 4 with 1 nmol [3H]-labeled polythymidine DNA.

Concentration: 5,000 units/ml

Usage Note: Exo T is has different activity on RNA vs. DNA. For RNA, 1 unit of Exo T is required to completely digest 1.0 pmol of rA20 under standard reaction condition as measured by gel electrophoresis.

Nucleoside Digestion Mix

#M0649S 50 reactions 119 €

- Convenient one-step protocol
- Digests both DNA and RNA to single nucleosides
- Low-glycerol formulation significantly reduces glycerol-induced ion suppression during MS analysis

Description: The Nucleoside Digestion Mix is a mixture of enzymes that provides a convenient one-step method to generate single nucleosides from DNA or RNA. Optimized for quantitative analysis by liquid chromatographymass spectrometry (LC-MS), this reagent eliminates the need for sequential multi-step, time-consuming digestion protocols. The Nucleoside Digestion Mix digests ssDNA, dsDNA, DNA/RNA hybrids and RNA (except mRNA cap structures) containing epigenetically modified (m5C,

hm5C, f5C, ca5C, m4C, m6A, etc.), unnatural, or damaged bases. Moreover, the low-glycerol formulation (< 1%) significantly reduces glycerol-induced ion suppression during mass spectrometry analysis.

Reaction Conditions: 1X Nucleoside Digestion Mix Reaction Buffer. Incubate at 37°C. May heat inactivate at 65°C for 20 minutes.

DNase I (RNase-free)

1.000 units71 € #M0303L 5,000 units 284 €

- Degradation of DNA template in transcription reactions
- Removal of contaminating genomic DNA from RNA samples
- DNase I footprinting
- Nick Translation

Description: DNase I (RNase-free) is an endonuclease that nonspecifically cleaves DNA to release di-, tri- and oligonucleotide products with 5´-phosphorylated and 3´-hydroxylated ends. DNase I acts on single- and double-stranded DNA, chromatin and RNA:DNA hybrids.

Reaction Conditions: 1X DNase I Reaction Buffer. Incubate at 37°C. May heat inactivate at 75°C for 10 minutes.

RR 37° 🛣

RX

Unit Definition: One unit is defined as the amount of enzyme which will completely degrade 1 µg of pBR322 DNA in a total reaction volume of 50 µl in 10 minutes at 37°C.

Complete degradation is defined as the reduction of the majority of DNA fragments to tetranucleotides or smaller.

Concentration: 2.000 units/ml

Note: EDTA should be added to a final concentration of 5 mM to protect RNA from being degraded during enzyme inactivation.

RNase Inhibitor, Murine

3,000 units71 € #M0314S #M0314L 15,000 units 284 €

- Inhibits common eukaryotic RNases
- Compatible with Tag Polymerase, AMV or M-MuLV Reverse Transcriptases
- cDNA synthesis & RT-PCR
- In vitro transcription/translation
- Enzymatic RNA labeling reaction

RNase Inhibitor, Murine has significantly improved resistance to oxidation compared to human & porcine RNase inhibitors.

Description: RNase Inhibitor, Murine is a 50 kDa recombinant protein of murine origin. It specifically inhibits RNases A, B and C by binding noncovalently in a 1:1 ratio with high affinity. It is not effective against RNase 1, RNase T1, S1 Nuclease, RNase H or RNase from Aspergillus. No inhibition of polymerase activity is observed when used with Taq DNA Polymerase, AMV or M-MuLV Reverse Transcriptases, or Phage RNA Polymerases (SP6, T7, or T3).

Recombinant RNase Inhibitor, Murine does not contain the pair of cysteines identified in the human version that are very sensitive to oxidation and lead to inactivation of the inhibitor. As a result, RNase Inhibitor, Murine has significantly improved resistance to oxidation compared to the human/porcine RNase inhibitors, and is stable at low DTT concentrations (< 1 mM). This makes it ideal for reactions where high concentration DTT is adverse to the reaction (e.g., RT-qPCR).

Unit Definition: One unit is defined as the amount of RNase Inhibitor, Murine required to inhibit the activity of 5 ng of RNase A by 50%. Activity is measured by the inhibition of hydrolysis of cytidine 2´, 3´-cyclic monophosphate by RNase A.

Concentration: 40,000 units/ml







PCR PCR Enzyme













202

RNase Inhibitor, Human Placenta

#M0307S 2,000 units78 € #M0307L 10,000 units312 €

- Inhibits common eukaryotic RNases
- Compatible with Taq Polymerase, AMV or M-MuLV Reverse Transcriptases
- Active over a broad pH range (pH 5–8)
- cDNA synthesis reactions
- In vitro transcription/translation

Description: RNase Inhibitor, Human Placenta is a recombinant human placental protein which specifically inhibits ribonucleases (RNases) A, B and C. It is not effective against RNase 1, RNase T1, S1 Nuclease, RNase H or RNase from *Aspergillus*. In addition, no inhibition of polymerase activity is observed when RNase Inhibitor, Human Placenta is used with *Taq* DNA Polymerase, AMV or M-MuLV Reverse Transcriptases, or Phage RNA Polymerases (SP6, T7, or T3).

RX

The 50 kDa protein inhibits RNases by binding noncovalently in a 1:1 ratio with an association constant greater than 10¹⁴.

Unit Definition: One unit is defined as the amount of RNase Inhibitor, Human Placenta required to inhibit the activity of 5 ng of RNase A by 50%. Activity is measured by the inhibition of hydrolysis of cytidine 2´, 3´-cyclic monophosphate by RNase A.

Concentration: 40,000 units/ml

Ribonucleoside Vanadyl Complex

#S1402S 10 ml (200 mM)82 €

Ribonucleoside Vanadyl Complex is an equimolar mixture of all four ribonucleosides, complexed with oxovanadium IV by a modification of the procedures by Berger (1). Each lot of the complex is assayed for oxovanadium V content and inhibition of ribonuclease activity.

Vanadium complexes are used in mRNA purifications as exogenous ribonuclease inhibitors. They are compatible with cell lysis techniques and with sucrose gradient fractionation of cytoplasmic components.

Reference:

(1) Berger, S.L. and Birkenmeier, C.S. (1979) *Biochemistry* 18, 5143–5149.

NEBNext® Reagents for RNA Library Preparation

NEBNext Ultra II Directional RNA Library Prep Kit for Illumina

#E7760S 24 rxns1028 € #E7760L 96 rxns3495 €

NEBNext Ultra II Directional RNA

Library Prep with Sample Purification Beads
#E7765S 24 rxns1140 €
#E7765L 96 rxns3880 €

NEBNext Ultra II RNA Library Prep Kit

for Illumina

#E7770S 24 rxns 980 € #E7770L 96 rxns 3325 €

NEBNext Ultra II RNA Library Prep with Sample Purification Beads

#E7775S 24 rxns1080 € #E7775L 96 rxns3685 €

NEW

NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina

#E6420S 24 rxns1195 € #E6420L 96 rxns3965 €

NEW


NEBNext Single Cell/Low Input cDNA Synthesis & Amplification Module

#E6421S 24 rxns650 € #E6421L 96 rxns2210 €

rRNA depletion and poly(A) mRNA isolation reagents are available separately. More information on NEBNext reagents for RNA library preparation can be found on pages 142–148.

Do you need increased sensitivity and specificity from your RNA-seq experiments? Do you have ever-decreasing amounts of input RNA? To address these challenges, our next generation of RNA library prep kits have been reformulated at each step, resulting in several fold higher yields of high quality libraries and enabling use of lower

input amounts, including single cells, and fewer PCR cycles. The kits have streamlined, automatable workflows and are available for directional (strand-specific, using the "dUTP method") and non-directional library prep, with the option of SPRISelect beads for size-selection and clean-up steps.



NEBNext Ultra II Directional RNA produces the highest yields, from a range of input amounts. Poly(A)-containing mRNA was isolated from 10 ng, 100 ng and 1 µg of Universal Human Reference RNA (Agilent #740000) and libraries were made using the NEBNext Ultra II Directional RNA kit, Kapa™ Stranded mRNA-Seq kit, Kapa mRNA HyperPrep kit and Illumina® TruSeq® Stranded mRNA Kit. The input RNA amount and number of PCR cycles are indicated. Library yields from an average of three replicates are shown. Error bars indicate standard deviation. Library yields were assessed using the Agilent® Bioanalyzer®.

EpiMark® N6-Methyladenosine Enrichment Kit

#E1610S 20 reactions 420 €

- Enrichment for m6A modified RNA in immunoprecipitation protocols
- Enriched RNA can be used directly for next gen sequencing or RT-qPCR

Description: The EpiMark N6-Methyladenosine Enrichment Kit contains a rabbit monoclonal antibody specific for N6-Methyladenosine (m6A). The kit also contains two control RNAs, one with m6A modification (Gaussia luciferase) and one without (Cypridina luciferase) to monitor enrichment and depletion. The GLuc RNA control was transcribed in the presence of 20% m6ATP and 80% ATP.

This kit can be used to enrich m6A modified RNA in immunoprecipitation protocols for downstream analysis by next-generation RNA sequencing or RT-qPCR. Modified RNA is isolated from a fragmented RNA sample by binding to the N6-Methyladenosine

Epi

R_{{

antibody attached to Protein G Magnetic Beads. After multiple wash and clean-up steps, the enriched RNA is eluted in nuclease-free water and is ready for further analysis.

The EpiMark N6-Methyladenosine **Enrichment Kit Includes:**

- N6-Methyladenosine Antibody
- m6A Control RNA (100 nM)
- Unmodified Control RNA (100 nM)

N6-Methyladenosine Antibody is produced by Cell Signaling Technology, Inc. and sold by New England Biolabs, Inc

p19 siRNA Binding Protein

#M0310S 1.000 units82 €

Companion Products:

Chitin Resin #S6651S 20 ml 85 € #S6651L 100 ml340 € Chitin Magnetic Beads #E8036S 5 ml122 € #F8036I 25 ml 488 €

Description: The p19 siRNA Binding Protein (19 kDa) from the plant Carnation Italian Ringspot Virus (CIRV) binds siRNAs with nanomolar affinity. The dimeric protein preferentially binds 21 -nucleotide siRNAs with a 2 -nucleotide 3' extension and a 5' phosphate. The protein binds RNA in a size-dependent and sequenceindependent manner. If the siRNAs are 4 bases longer, the affinity for the protein is reduced about 100 fold. When p19 siRNA Binding Protein is expressed in plants it suppresses RNA interference.

Source: p19 siRNA Binding Protein is cloned and expressed in E. coli as a fusion protein with an amino terminal MBP (Maltose-binding protein) and a carboxy terminal CBD (Chitin binding domain).

Unit Definition: One unit is defined as amount of protein that binds to 10 ng of siRNA at 25°C in 1 hour. Unit assay conditions can be found at www.neb.com.

Molecular Weight: 67 kDa Concentration: 10.000 units/ml

- High affinity binding of siRNAs
- Affinity purification of siRNA with chitin magnetic beads

Magnetic mRNA Isolation Kit

#S1550S 25 isolations 321 €

Component Sold Separately:

Oligo d(T)₂₅ Magnetic Beads #S1419S 25 mg 262 €

Companion Products:


6-tube Magnetic Separation Rack #S1506S 6 tubes (1.5 ml) 206 € 12-tube Magnetic Separation Rack 12 tubes (1.5 ml) 326 € #S1509S

- Suitable for automated high-throughput applications
- Eliminates need for organic solvents
- No need to precipitate poly(A)+ transcripts in eluent
- Obtain intact poly(A)+ RNA in < 1 hour
- Negligible gDNA contamination

Description: The New England Biolabs Magnetic mRNA Isolation Kit is designed to isolate intact poly(A)+ RNA from cells and tissue without requiring phenol or other organic solvents. The technology is based on the coupling of Oligo $d(T)_{25}$ to 1 µm paramagnetic beads which are then used as the solid support for the direct binding of poly(A)+ RNA. Thus, the procedure permits the manual processing of multiple samples and can be adapted for automated high-throughput applications. Additionally, magnetic separation technology permits elution of intact mRNA in small volumes eliminating the need for precipitating the poly(A)+ transcripts in the eluent. Intact poly(A)+ RNA which is fully representative of the mRNA population of the original sample can be obtained in less than one hour. Oligo $d(T)_{25}$ Magnetic Beads can be reused up to three times and the researcher has the option of eluting the isolated mRNA or using the bound dT DNA as a primer in a first-strand cDNA reaction.

The Magnetic mRNA Isolation Kit Includes:

- Oligo d(T)₂₅ Magnetic Beads
- Lvsis/Binding Buffer
- Wash Buffer I, II, III, Elution Buffer



Consistency and wide isolation range are demonstrated by poly(A)+ RNA isolation from duplicate samples of decreasing numbers of HEPG2 cells (5 x 105 to 1 x 103) by direct lysis/ binding in microtiter plates followed by mRNA isolation with the magnetic method. 1/10th of isolated mRNA is converted to oligo (dT) primed cDNA using ProtoScript M-MuLV First Strand cDNA Synthesis Kit (NEB #E6300) and qPCR done with validated primers for the peptidylpropyl isomerase, a low-abundance housekeeping gene.















Oligo d(T), Magnetic Beads

#S1419S 25 mg 262 €

Companion Product:

96-well Microtiter Plate Magnetic Separation Rack

#S1511S 96 well 489 € Description: An affinity matrix for the small-scale isolation of mRNA from crude cell lysates and tissue. The isolation occurs through the hybridization of covalently coupled oligo $d(T)_{25}$ to the poly(A) region present in most eukaryotic mRNAs. Applications include direct mRNA isolation following lysis and second-round purification of previously isolated total RNA. The magnetic separation technology is scalable and permits elution of intact mRNA in small volumes eliminating the need for precipation of the isolated mRNA. Beads can be reused up to three times, and the researcher has the option of eluting the isolated mRNA or using the bound d(T)_{as} as a primer in a first-strand cDNA reaction.

Beads are supplied as a 5 mg/ml suspension in phosphate buffer (PBS) (pH 7.4), containing 0.05% Tween-20 and 0.05% NaN₃

Support Matrix: 1 µm nonporous superparamagnetic microparticles.

Binding Capacity: 1 mg of Oligo d(T)₂₅ Beads will bind 10 μg of poly(A)+ RNA.

Oligo (dT)₂₅ Cellulose Beads

#S1408S 250 mg 132 €

Description: An affinity matrix used for the isolation of mRNA containing polyadenylic (poly A) regions. This matrix consists of oligo $(dT)_{25}$ covalently coupled to a cross-linked cellulose bead.

Support Matrix: A cross-linked cellulose bead.

Binding Capacity: > 400 O.D. per gram of cellulose

polyA Spin™ mRNA Isolation Kit

#S1560S 8 isolations 224 €

Companion Product:

Oligo d(T) Cellulose Beads #S1408S 250 mg132€

Description: The polyA Spin Kit is a rapid and convenient alternative for the purification of full-length poly(A)+ eukaryotic messenger (mRNA) from samples of total RNA. Poly(A)+ RNA selection is made by affinity chromatography using spin columns prepackaged with NEB's Oligo (dT)₀₅-Cellulose Beads (NEB #S1408). Its high binding capacity and rapid hybridization kinetics make it an ideal support for chromatographic poly(A)+ RNA selection. Intact mRNA can be isolated in as little as forty minutes from multiple eukaryotic cell lysates or samples of total RNA from various sources, including fungi, plants and animal tissue.

Reagents sufficient for the isolation and subsequent precipitation of poly(A)+ RNA from eight samples of as much as 1 mg of total RNA are provided. The isolated RNA can be used for in vitro translation, preparation of cDNA libraries, northern analysis, subtractive hybridization or differential display.

The polyA Spin Kit Includes:

- Prepacked Oligo dT₂₅-Cellulose Beads
- Sterile microfuge recovery tubes
- Wash Buffer, Elution Buffer, Glycogen, NaCl and NaOAc

Streptavidin Magnetic Beads

#S1420S 5 ml (20 mg) 280 €

Companion Products:

6-Tube Magnetic Separation Rack #S1506S 6 tubes (1.5 ml) 206 €

96-well Microtiter Plate Magnetic Separation Rack #S1511S 96 well489 €

superparamagnetic particles covalently coupled to a highly pure form of streptavidin. The beads can be used to capture of biotin labeled substrates including antigens, antibodies and nucleic acids. The strength of the biotin-streptavidin interaction coupled with low non-specific binding permits captured substrates to be useful as ligands in subsequent experiments including mRNA isolation and the capture of primary or

Description: Streptavidin Magnetic Beads are 1 µm

Beads are supplied as a 4 mg/ml suspension in phosphate buffer (PBS) (pH 7.4) containing 0.1% BSA, 0.05% Tween-20 and 0.05% NaN_a.

secondary antibodies.

Support Matrix: 1 µm non-porous superparamagnetic microparticle.

Binding Capacity: The beads will bind greater than 1000 pmol of free biotin per mg and greater than 500 pmol of single-stranded 25 bp biotinylated oligonucleotide per mg.

NEW

Monarch Kits for Cleanup & Isolation

See pages 130-131 for more information.

Monarch Total RNA	A Miniprep Kit
#T2010S	50 preps248 €
Monarch RNA Clea	anup Kit (10 μg)
#T2030S	10 preps 52 €
#T2030L	100 preps 284 €
Monarch RNA Clea	anup Kit (50 μg)
#T2040S	10 preps 50 €
#T2040L	100 preps 280 €
Monarch RNA Clea	anup Kit (500 µg)
#T2050S	10 preps 58 €
#T2050L	100 preps 440 €

The Monarch Total RNA Miniprep Kit is a comprehensive solution for sample preservation, cell lysis, gDNA removal, and purification of total RNA from a wide variety of biological samples, including cultured cells, blood, and mammalian tissues. Additionally, tough-tolyse samples, such as bacteria, yeast, and plant, can be processed with additional steps that enhance lysis. Cleanup of enzymatic reactions or purification of RNA from TRIzol® -extracted samples is also possible using this kit. Purified RNA has high quality metrics, including $A_{260/280}$ and $A_{260/230}$ ratios \geq 1.8, high RIN scores, and minimal residual gDNA. Captured RNA ranges in size from full-length rRNAs down to intact miRNAs. Additionally, differential binding conditions allow selective capture or exclusion of the sub-200 nucleotide RNA pool that includes miRNA, 5S rRNA, and tRNA. Purified RNA is suitable for downstream applications, such as RT-qPCR, cDNA synthesis, RNA-seq, Northern blot analysis, etc.

The Monarch RNA Cleanup Kits provide a fast and simple silica spin column-based solution for RNA cleanup and concentration after any enzymatic reaction (including *in vitro* transcription, DNase I treatment, capping and labeling) and after other purification methods such as phenol/chloroform extraction. The Monarch RNA Cleanup Kits are available in 3 different binding capacities: 10 μ_{B} , 50 μ_{B} and 500 μ_{B} . Each kit contains unique columns, all designed to prevent buffer retention and ensure no carryover of contaminants, enabling low-volume elution of highly-pure RNA. Following the standard protocol, RNA \geq 25 nt is purified with this kit; however, a modified protocol is available to enable the binding of RNA as small as 15 nt (including miRNAs).

The Monarch Total RNA Miniprep Kit Includes:

- Monarch gDNA Removal Columns
- Monarch RNA Purification Columns
- Monarch Collection Tubes II
- Monarch DNA/RNA Protection Reagent (2X)
- Monarch RNA Lysis Buffer
- Monarch Proteinase K
- Monarch Proteinase K Resuspension Buffer
- Monarch Proteinase K Reaction Buffer
- Monarch DNase I
- Monarch DNase I Reaction Buffer
- Monarch RNA Priming Buffer
- Monarch RNA Wash Buffer (5X)
- Monarch Nuclease-free Water

The Monarch RNA Cleanup Kits Include:

- Monarch RNA Cleanup Columns (10, 50 or 500 μg)
- Monarch RNA Cleanup Binding Buffer
- Monarch RNA Cleanup Wash Buffer
- Monarch Collection Tubes II
- Nuclease-free Water



Claire has been with NEB for just over a year as a Quality Systems & Documentation Specialist. She is already an active member of the NEB community, having participated on the Raffle Committee, Education Committee, and in the Gloucester BioTech Academy/

RNA Markers & Ladders

dsRNA Ladder

#N0363S 25 gel lanes90 €

microRNA Marker

#N2102S 100 gel lanes 69 €

ssRNA Ladder

#N0362S 25 gel lanes 69 €


Low Range ssRNA Ladder

#N0364S 100 gel lanes 68 €

NEB offers several RNA Markers and Ladders with a size range from 17 to 9,000 bases. The ssRNA ladders are suitable for use as RNA size standards on denaturing or native gels. Both are supplied with 2X Loading Buffer and feature a higher intensity fragment to serve as a reference band. The microRNA Marker, supplied in a ready-to-load denaturing solution, is ideally used as a size marker on polyacrylamide gels or Northern blots and is best visualized stained with ssRNA fluorescent dyes. It is

supplied with a 3'-biotinylated 21-mer oligonucleotide probe that can also be labeled with $[\gamma^{-32}P]$ ATP and T4 PNK (NEB #M0201). The ds RNA Ladder is suitable for use as a size standard in dsRNA analysis on both polyacrylamide and agarose gels.

Concentration: Low Range ssRNA Ladder and dsRNA Ladder are supplied at $500 \, \mu g/ml$. ssRNA Ladder is supplied at $2,000 \, \mu g/ml$. MicroRNA Marker is supplied at $12 \, ng/\mu l$.



RNA Loading Dye (2X)

#B0363S 4 ml49 €

Description: The RNA Loading Dye, (2X) is a premixed loading dye for use with denaturing and non-denaturing PAGE/agarose gels.

RNA Loading Dye Composition:

1X RNA Loading Dye:

47.5% Formamide, 0.01% SDS, 0.01% Bromophenol Blue, 0.005% xylene cyanol and 0.5 mM EDTA.

Universal miRNA Cloning Linker

#S1315S 0.83 nmol 148 €

Companion Product:

T4 RNA Ligase 2, truncated KQ #M0373S 2,000 units 69 € #M0373L 10,000 units 276 € This 5-adenylated, 3´-blocked oligoribonucleotide can be used for cloning short RNAs according to the procedure of Bartel (1). RNA ligase recognizes the "activated" adenylated oligo and covalently links (ligates) its 5´ end to the 3´ OH of a second single stranded sequence in the absence of ATP. In a mixture of nucleic acids use of the 5´ adenylated, 3´ blocked oligo with T4 RNA Ligase 2.

truncated, T4 RNA Ligase 2, truncated K227Q or T4 RNA Ligase 1 (w/o ATP) results in ligation of the target oligo only.

The sequence of the adenylated DNA oligo is 5'-rAppCTGTAGGCACCATCAAT-NH_a 3'.

Reference:

(1) Lau et al. (2001) Science, 294, 858-856.





Microbes Cleaning Up Our Mess

In a world where we consume 90 million barrels of oil every day, spills at oil rigs are inevitable. In fact, according to the U.S. Environmental Protection Agency, there are 70 spills per day — most are insignificant. However, in 2010, 4.9 million barrels of crude oil spewed from the BP Macondo well, located in the Gulf of Mexico, over a period of 87 days. This unfortunate event was termed the Deepwater Horizon spill, and is the largest marine oil spill in history.

Extensive damage to marine ecosystems resulted from the spill, including heavy oiling of four of the five endangered, protected turtle species that live and breed in the Gulf of Mexico. Additionally, many marine mammals that live in the Gulf of Mexico experienced reproductive failure and organ damage following the spill.

Large fire booms were used to surround the slick and burn off approximately 5% of the oil. Chemical dispersants broke 40% of the oil into tiny droplets suspended under the ocean surface in a plume that measured 161 km (100 miles) long. Scientists were surprised to find that the clean-up occurred at a much faster rate than predicted, because it was aided by hydrocarbon-degrading microorganisms. Further, the use of chemical dispersants, while controversial due to their toxicity, actually aided this process by generating smaller size oil droplets (and a larger surface area) for what was estimated to be a bloom of 100 sextillion microbes.

Background levels of oil exist in many marine ecosystems due to slow seepage from underwater reservoirs. Luckily, a well-established community of bacteria capable of degrading the compounds contained in oil is also present. The short doubling time of these bacteria and their ability for horizontal gene transfer allow for the rapid proliferation of this diverse community following a spill. These hydrocarbon-degraders use the oil as an energy source. Each microbe population has a distinct set of hydrocarbon degradation genes with specificity for breaking down different constituents of the oil, and they proliferate at different time points following a spill, depending on their specificity. Following degradation of the oil, the microbes stop rapidly dividing and are consumed by other organisms further up the food chain.

Oil from the Deepwater Horizon spill washed into wetlands and marshes all along the northern Gulf of Mexico, and heavily oiled a coastal plant named *Spartina alterniflora*, commonly known as smooth cordgrass. Cordgrass is a foundational plant that stabilizes land and minimizes erosion. Scientists found that in addition to the presence of oil in the plant tissues, the composition of endophytes living symbiotically in cordgrass roots showed an increase in oil-degrading bacteria. It is well known that these bacteria can degrade oil in oceans and on beaches, but what researchers are trying to establish is whether the bacteria continue to process the oil inside the plant tissue, and if so, which bacteria are most efficient at accomplishing this. Scientists are also exploring whether the plant roots can deliver hydrocarbon-degrading bacteria to the buried oil in the marsh. This raises the possibility for bacterial inoculation of cordgrass, which can then be planted and used to help restore coastal areas exposed to an oil spill.

The Deepwater Horizon oil spill was catastrophic to marine habitats. What resulted from the close study of these habitats was the identification and understanding of a diverse community of microorganisms, which are highly evolved to deal with their specialized environment and when necessary, can clean up the mess we make. Yet another reason to appreciate microbes.

Protein Expression & Purification Technologies



NEB offers an array of solutions for robust expression of your target protein.

At first glance, recombinant protein expression looks quite simple. Essentially, DNA encoding a target protein is cloned downstream of a promoter in an expression vector. This vector is then introduced into a host cell, and the cell's protein synthesis machinery produces the desired protein. In practice, however, protein expression can be very challenging because so many factors may influence the process. For example, each protein folds in its own unique manner, a process that may be influenced by the choice of expression host. Similarly, some proteins require post-translational modifications or proper insertion into a biological membrane. Finally, some proteins may have an activity that is detrimental to the host. Thus, no single solution exists for successful production of all recombinant proteins. Instead, it is beneficial to have access to a wide range of expression tools, and a willingness to explore multiple approaches to better one's chances for success.

NEB offers an array of expression systems offering different advantages, enabling you to choose the strategy that best suits your protein expression and purification needs. Many share a compatible polylinker, enabling the gene of interest to be easily shuffled between systems. Additionally, a selection of competent cells is available for *in vitro* expression of difficult-to-expess proteins.

Featured Products

PURExpress® In Vitro Protein Synthesis Kit

217 Competent Cells for Protein Expression

Featured Tools & Resources

Protein Expression Selection Chart

361 Enhancing Efficiency

362 Protein Expression with T7 Express Strains



Did you know you can find application notes using our expression systems on our website? Choose the applications tab at www.neb.com to learn more.



Protein Expression Selection Chart	212	Cell-free Expression PURExpress <i>In Vitro</i> Protein Synthesis Kit	216	
EXPRESSION SYSTEMS		PURExpress Δ Ribosome Kit	216	
E. coli		PURExpress Δ (aa, tRNA) Kit	216	
pMAL Protein Fusion & Purification System	213	PURExpress Δ RF123 Kit	216	
pMAL Companion Products:		PURExpress Companion Products:		
Amylose Resin	213	PURExpress Disulfide Bond Enhancer	216	
Amylose Resin High Flow	213	E. coli Ribosome	216	
Amylose Magnetic Beads	213	Competent Cells		
Anti-MBP Magnetic Beads	213		9–233	
Anti-MBP Monoclonal Antibody	213	Competent Cell Expression Strain Selection Cha		
Factor Xa Protease	213		217	
pMAL Vectors	213	MAGNETIC MATRICES AND RACKS		
IMPACT Kit	214	Protein A Magnetic Beads	218	
IMPACT Vector Comparison	214	Protein G Magnetic Beads	218	
IMPACT Companion Products:		Streptavidin Magnetic Beads	218	
Anti-CBD Monoclonal Antibody	214	Hydrophilic Streptavidin Magnetic Beads	218	
Chitin Resin	214	Goat Anti-Rabbit IgG Magnetic Beads	219	
IMPACT Vectors	214	Goat Anti-Mouse IgG Magnetic Beads	219	
		Goat Anti-Rat IgG Magnetic Beads	219	
Yeast		Chitin Magnetic Beads Amylose Magnetic Beads	219 219	
K. lactis Protein Expression Kit	215	Ni-NTA Magnetic Beads	219	
K. lactis Companion Products:		Anti-MBP Magnetic Beads	213	
BstXI	32	6-Tube Magnetic Separation Rack	219	
SacII	47	12-Tube Magnetic Separation Rack	219	
Yeast Medium Pack	215	50 ml Magnetic Separation Rack	219	
		96-Well Microtiter Plate Magnetic Separation Rack	219	
		SNAP-Cantura Magnetic Reads	283	

Protein Expression & Purification Selection Chart

NEB offers an array of expression systems with different advantages, enabling you to choose the strategy that best suits your protein expression and purification needs. Many share a compatible polylinker, enabling the gene of interest to be easily shuffled between systems. Additionally, a selection of competent cells is available for *in vitro* expression of difficult-to-expess proteins.

APPLICATION	КІТ	ADVANTAGES		
	pMAL™ Protein Fusion and Purification System	Substantial yields (up to 100 mg/L) in more than 75% cases tested; uses the strong P _{tac} promoter		
High Yield Expression	K. lactis Protein Expression Kit	Uses the strong <i>LAC4</i> promoter; multiple integrations of plasmid results in higher yield		
	IMPACT™ Kit	Uses the T7 promoter for high level regulated expression		
Cell-free Expression	PURExpress® In Vitro Protein Synthesis Kit	Quickly generates analytical amounts of protein		
	K. lactis Protein Expression Kit	Easily co-express 2–4 proteins		
Co-expression of Multiple Proteins	PURExpress In Vitro Protein Synthesis Kits	Bicistronic constructs or multiple plasmids with appropriate transcription and translation regulatory elements can be used		
Enhanced Calubility	pMAL Protein Fusion and Purification System	Fusion to MBP enhances solubility of proteins in E.coli*		
Enhanced Solubility	K. lactis Protein Expression Kit	Utilizes <i>K. lactis</i> eukaryotic folding pathway		
	IMPACT Kit	Utilizes an intein-chitin binding domain (CBD) tag on either the N- or C- terminus, offering single-step purification		
Affinity Tag Chromatography	pMAL Protein Fusion and Purification System	Fusion to MBP allows for purification on amylose resin		
	K. lactis Protein Expression Kit	Vectors are sold separately that generate fusions to MBP allowing for purification on amylose resin		
Post-translational Modification	K. lactis Protein Expression Kit	Secretion of both N- and O- glycosylated proteins		
Periplasmic Expression	pMAL Protein Fusion and Purification System	Periplasmic expression enhances folding of proteins with disulfide bonds		
Secreted Expression	K. lactis Protein Expression Kit	Eliminates cell lysis step, simplifying purification		
	K. lactis Protein Expression Kit	Secretion of protein from the cell		
Tovio Dustoino	IMPACT Kit	Can express the toxic gene in two pieces and ligate proteins together using intein mediated protein ligation (IPL)		
Toxic Proteins	pMAL Protein Fusion and Purification System	Can export toxic proteins into periplasmic space		
	PURExpress In Vitro Protein Synthesis Kits	Cell-free environment not affected by "toxicity" of expressed protein		
Protein Modification	IMPACT Kit	Generates proteins with reactive ends (N-terminal cysteine and/or C-terminal thioester) allowing for labeling or ligation of proteins or peptides.		
	PURExpress In Vitro Protein Synthesis Kits	Allows introduction of modified, unnatural or labeled amino acids		
No Additional Amino Acid Decidues	IMPACT Kit	Start of native protein is fused adjacent to site of intein cleavage		
No Additional Amino Acid Residues	pMAL Protein Fusion and Purification System	Start of protein is fused adjacent to protease site		

^{*}Kapust and Waugh (1999) Protein Science, 8, 1668-1674.

pMAL™ Protein Fusion and Purification System

#E8200S 626 €

Components Sold Separately:

Amylose Resin #E8021S #E8021L	15 ml186 € 100 ml970 €
Factor Xa Protease #P8010S #P8010L	50 μg71 € 250 μg284 €
pMAL-p5X Vector #N8109S	10 μg117€
pMAL-c5X Vector #N8108S	10 µg120€

- Reliable expression: substantial yields (up to 100 mg/L) in more than 75% of the cases tested
- Expression in either the cytoplasm or periplasm: folding of proteins with disulfide bonds can be enhanced by expression in the periplasm or expression in the cytoplasm of SHuffle® cells
- Enhanced solubility: MBP fusion proteins demonstrate improved solubility when expressed in E. coli
- Gentle elution with maltose: no detergents or harsh denaturants

For a more detailed description and a restriction map of the pMAL-p5X vector, including the MCS, see page 374. See our website, www.neb.com, for pMAL sequence files.

Description: In the pMAL Protein Fusion and Purification System, the cloned gene is inserted into a pMAL vector downstream from the malE gene, which encodes maltose-binding protein (MBP). This results in the expression of a MBP-fusion protein. The technique uses the strong P_{tac} promoter and the translation initiation signals of MBP to express large amounts of the fusion protein. The fusion protein is then purified by a one-step affinity purification specific for MBP (Figure 1).

The system uses the pMAL vectors which include a sequence coding for the recognition site of a specific protease adjacent to the multiple cloning site. This allows the protein of interest to be cleaved from MBP after purification, without adding any vector-derived residues to the protein. For this purpose, the polylinker includes a restriction site superimposed on the sequence coding for the site of the specific protease. This is where the gene of interest is inserted.

Expression from the pMAL vectors yields up to 100 mg fusion protein from a liter of culture. In most cases, the expressed protein is soluble, as fusion to MBP has been proven to enhance the solubility of proteins expressed in *E. coli*. While no expression system works with every cloned gene, the pMAL Protein Fusion and Purification System gives substantial yields of protein in more than 75% of the cases tested.

In pMAL-c5X, the *malE* signal sequence is deleted, resulting in cytoplasmic expression of the fusion protein. The pMAL-p5X Vector contains the *malE* signal sequence directing the fusion protein through the cytoplasmic membrane into the periplasm.

References: References for properties and applications of this product can be found at www.neb.com.




Figure 1: Schematic of the pMAL System

The pMAL Protein Fusion and Purification System includes:

- pMAL-c5X, pMAL-p5X
- Amylose Resin (binding capacity > 4 mg/ml volume)
- Factor Xa Protease
- MBP5 (marker for SDS-polyacrylamide gels)
- MBP5-paramyosin-ΔSal (control protein for Factor Xa cleavage)
- E. coli ER2523 (NEB Express)

pMAL Companion Products

Amylose Resin #E8021S 15 ml 186 € #E8021L 100 ml 970 €

Amylose Resin High Flow

#E8022S 15 ml 246 € #E8022L 100 ml1474 € **Description:** Affinity matrix used for isolation of proteins fused to maltose-binding protein. It is a composite amylose/agarose bead. Amylose Resin High Flow is a more rigid bead, suitable for use in automated chromatography systems.

Binding Capacity: Amylose Resin and Amylose Resin High Flow: >4 mg MBP5-paramyosin Δ Sal fusion protein/ml of bed volume.

Amylose Magnetic Beads

#E8035S 25 mg 246 €

Anti-MBP Magnetic Beads

#E8037S 10 mg 264 €

Description: Affinity matrices for the small-scale isolation and purification of MBP-fusion proteins. Amylose or monoclonal Anti-MBP are covalently coupled to a paramagnetic particle through a linkage that is stable and leak resistant over a wide pH range.

Binding Capacity: 1 mg of Amylose Magnetic Beads will bind \geq 10 μ g of MBP-fusion protein.

1 mg of Anti-MBP Magnetic Beads will bind 5 μ g of MBP-paramyosin fusion protein.

Anti-MBP Monoclonal Antibody

#E8032S 0.05 ml 174 € #E8032L 0.25 ml 696 € **Description:** Anti-MBP Monoclonal Antibody is a murine anti-maltose-binding protein antibody, isotype

IgG2a. It is purified from tissue culture supernatant by protein A affinity chromatography.

IMPACT[™] Kit

#E6901S......378 €

Components Sold Se	parately:
pTXB1 Vector #N6707S	10 μg114 €
pTYB21 Vector #N6709S	10 μg114 €
Chitin Resin #S6651S #S6651L	20 ml 85 € 100 ml 340 €
Anti-CBD Monoclonal Anti #E8034S	body 0.05 ml70 €
pTWIN1 Vector #N6951S (not included with the kit)	10 μg110 €

- Single-column purification without the use of proteases
- Produce target protein without vector derived amino acids
- Fusion to either N- or C-terminus of target protein
- Ligation and labeling of recombinant proteins
- Isolation of proteins with or without N-terminal methionine

For a more detailed description and a restriction map of the pTXB1 and pTYB21 vectors, including the MCS, see page 378–379; visit www.neb.com for sequence files.

Description: The IMPACT (Intein Mediated Purification with an Affinity Chitin-binding Tag) Kit utilizes engineered protein splicing elements (inteins) to purify recombinant proteins by a single column (Figure 1). This kit distinguishes itself from other protein fusion systems by its ability to separate a recombinant protein from the affinity tag without the use of a protease.

The IMPACT Kit allows fusion of a tag consisting of an intein and the chitin binding domain (CBD) from *Bacillus circulans*, to either the C-terminus (pTXB1) or the N-terminus (pTYB21) of a target protein (Figure 2). In the presence of thiols, such as DTT, the intein undergoes specific self-cleavage which releases the target protein. The pTXB1 vector can also be used to express a protein with a C-terminal thioester for use in Intein-mediated Protein Ligation (IPL). The IPL reaction, also referred to as expressed protein ligation, allows for ligation of a peptide or a protein with a N-terminal cysteine to a bacterially expressed protein with a C-terminal thioester through a native peptide bond for use in protein labeling and semisynthesis. For more information on the IMPACT System and IPL, visit www.neb.com.

pTXB1 is a *E. coli* expression vector that utilizes a minintein from the *Mycobacterium xenopi gyrA* gene [*Mxe* GyrA intein; 22 kDa]. This intein has been modified and combined with the CBD to create an affinity tag which can be bound to chitin beads (NEB #S6651). Release of the target protein is induced by thiol-reagents such as DTT or 2-mercaptoethanesulfonic acid (for ligation).

The pTYB21 vector allows for the fusion of the intein tag containing the *Saccharomyces cerevisiae* (*Sce*) *VMA*1 intein and CBD to the N-terminus of the target protein.

pTWIN1 Vector is available separately and enables isolation of proteins with an N-terminal cysteine and/or a C-terminal thioester. The polylinker is designed for the in-frame fusion of a target gene between the modified *Ssp* DnaB and *Mxe* GyrA inteins. The presence of the CBD facilitates purification.

References: References for properties and applications of this product can be found at www.neb.com.

The IMPACT Kit Includes:

- E. coli expression vectors: pTXB1, pTYB21 and pMXB10 Control Vector
- Chitin Beads: Affinity matrix (binding capacity=2 mg/ml)
- Anti-CBD Monoclonal Antibody
- 1, 4-Dithiothreitol (DTT) & SDS Loading Buffer




Figure 1: Purification of Maltose Binding Protein (MBP) in a single affinity purification step: Lane 1: uninduced cell extract. Lane 2: induced cell extract showing expressed fusion protein. Lane 3: MBP fraction eluted after inducing cleavage overnight at 4°C. Lane 4: MBP ligated to a peptide containing an N-terminal cysteine. Marker M is the protein ladder.




Figure 2: Schematic of the IMPACT System.

Table 1: Guide to IMPACT vectors and applications.

VECTORS	SITE OF TARGET PROTEIN FUSION	INTEIN TAG (KDA)	RECOMMENDED CLONING SITES ^a	PREFERRED RESIDUES AT CLEAVAGE SITE ^b	METHOD OF CLEAVAGE ^{c,d}	APPLICATIONS
pTXB1	C-terminus	Mxe GyrA intein (28)	Ndel-Sapl/Spel	Y, F, O, N, T, K, A, H, M (Unfavorable residues- S, P, D, G)	DTT (or MESNA) pH 8.0-8.5, 4°C	Purification; C-terminal thioester for ligation and modification
pTYB21	N-terminus	Sce VMA1 intein (56)	Sapl/Bsml/Ndel- Pstl	A, Q, M, G, L, N, W, F, Y (Unfavorable residues- P, S, C, T, R)	DTT ^d pH 8.0-8.5, 25°C	Purification
pTWIN1	C-terminus (Intein 2)	Mxe GyrA intein (28)	Ndel-Sapi/Spel	M, Y, F, LEM (Unfavorable residues- S, P, E, D)	DTT (or MESNA) pH 8.0-8.5, 4°C	Purification; C-terminal thioester for ligation and modification

^a NEBuilder HiFi DNA Assembly Cloning Kit (NEB #E5520) can be used to generate construct without the use of restriction enzymes. ^bActual cleavage efficiency is dependent on the adjacent residues as well as the folding of the fusion protein. ^cDithiothreitol (DTT) is used only for protein purification. 2-mercaptoethanesulfonic acid (MESNA) is used for isolation of proteins possessing a C-terminal thioester for ligation, labeling and cyclization. ^dCysteine can be used in the place of DTT.

K. lactis Protein Expression Kit

#E1000S820 €

Components Sold Separately:

Sacii								
#R0157S	2,000 units	63€						
#R0157L	10,000 units	252€						
Yeast Med #B9017S		\$38€						
K. lactis GG799 Competent Cells								
#C1001S	5 transformation reactions	320€						

- Clone and express genes toxic to E. coli
- Simultaneous expression of multiple genes
- No expensive antibiotics or methanol required
- Easy-to-use protocols for those inexperienced with yeast systems
- Attractive commercial sublicensing

Restriction map for pKLAC2 can be found on page 373; sequence files are available at www.neb.com.

Description: The K. lactis Protein Expression Kit provides an easy method for expressing a gene of interest in the yeast Kluyveromyces lactis. The gene is cloned into the integrative pKLAC-series of vectors, and may be expressed intracellularly or secreted. The K. lactis system offers several advantages over other yeast and bacterial expression systems. Abundant overexpression of protein is achieved through high culture densities as well as the ability to integrate multiple copies of the vector. The pKLAC-series of vectors use a strong LAC4 promoter, which has been modified to lack expression in E. coli, making this system useful for expressing toxic genes. For the selection of transformants, no expensive antibiotics are required. In addition, no methanol is required in growth media. Finally, the K. lactis system can express post-translationally modified proteins, making it a useful alternative to bacterial expression systems.


pKLAC2 is a general purpose expression vector. Using this vector, proteins may be produced intracellularly or may be fused to the $\it K.$ lactis $\it cc$ -mating factor sequence for secreted expression. Vector pKLAC2 contains an MCS that is compatible with other expression systems available from NEB.

GG799 competent cells are provided in the *K. lactis* Protein Expression Kit. GG799 cells are characterized by very high cell density growth and efficient expression of foreign proteins. GG799 cells have no genetic modifications.


The K. lactis Protein Expression Kit Includes:

- pKLAC2 Vector and pKLAC1-malE
 Control Plasmid
- SacII
- Integration Primer Set
- CutSmart® Buffer
- K. lactis GG799 Competent Cells and Transformation Reagent
- Yeast Medium Powder & Acetamide Solution

References: References for properties and applications of these products can be found at www.neb.com.



Secreted protein processing. In the nucleus, an integrated expression vector encoding a fusion between the $\alpha\textsc{-MF}$ domain (blue) and a desired protein (black) is expressed. A signal peptide in the $\alpha\textsc{-MF}$ domain directs entry of the fusion protein into the endoplasmic reticulum (ER) and is removed by signal peptidase (SP). The fusion protein is transported to the Golgi apparatus where the Kex protease removes the $\alpha\textsc{-MF}$ domain. The protein of interest is then secreted from the cell.



Protein Expression in K. lactis. SDS-polyacrylamide gel electrophoresis separation of secreted recombinant maltosebinding protein (MBP) detected directly in peptone rich growth medium by Coomassie staining. Lane 1: Protein Molecular Weight Markers. Lane 2: spent culture medium (15 µl) from Wild-type K. lactis cells. Lane 3: spent culture medium (15 µl) from K. lactis cells harboring an integrated expression cassette containing the E. coli malE gene.

Companion Products:

Enterokinase, ligl	ht chain	One <i>Taq</i> DNA F	Polymerase
#P8070S	480 units112 €	#M0480S	200 units38 €
#P8070L	2,560 units448 €	#M0480L	1,000 units152 €
BstXI		#M0480X	5,000 units600 €
#R0113S	1,000 units68 €	One Taq Hot St	art DNA Polymerase
#R0113L	5,000 units272 €	#M0481S	200 units78 €
		#M0481L	1,000 units312 €
		#M0481X	5,000 units 1248 €
		Deoxynucleoti	de (dNTP) Solution Mix
		#N0447S #N0447L	8 μmol of each*61 € 40 μmol of each244 €
		#110441 L	το μποι οι σασίι244 C

^{*}Available in 4 x 0.2 ml aliquots

PURExpress® In Vitro Protein Synthesis Kits

PURExpress *In Vitro* Protein Synthesis Kit #E6800S 10 reactions261 € #E6800L 100 reactions 2292 €

PURExpress Δ Ribosome Kit

#E3313S 10 reactions 350 €

PURExpress ∆ (aa, tRNA) Kit

#E6840S 10 reactions350 €

PURExpress ∆ RF123 Kit

#E6850S 10 reactions 358 €

Companion Product:

PURExpress Disulfide B	Bond Enhancei	
#E6820S	50 reactions	242€
E. coli Ribosome		
#P0763S	1 mg	182 €

- Generation of analytical amounts of proteins for further characterization
- Confirmation of open reading frames
- Generation of truncated proteins to identify active domains and functional residues
- Introduction of modified, unnatural or labeled amino acids (NEB #E6840, #E6850)
- tRNA structure and function studies (NEB #E6840)
- Ribosome structure and function studies (NEB #E3313, #P0763)
- Release factor function studies/ ribosome display (NEB #E6850)
- Epitope mapping

PURExpress is based on the PURE system technology originally developed by Dr. Takuya Ueda at the University of Tokyo and commercialized as the **PURESYSTEM™** by Biocomber (Tokyo, Japan).

PURESYSTEM™ is a trademark of Post Genome Institute

Description: A rapid method for gene expression analysis, PURExpress is a novel cell-free transcription/ translation system reconstituted from the purified components necessary for *E. coli* translation. The nuclease-free and protease-free nature of the PURExpress platform preserves the integrity of DNA and RNA templates/complexes and results in proteins that are free of modification and degradation. Transcription and translation are carried out in a one-step reaction, and require the mixing of only two tubes. With results available in a few hours, PURExpress saves valuable laboratory time and is ideal for high throughput technologies.

Advantages:

- Suitable for circular or linear DNA template
- Visualize synthesized protein directly on a Coomassie stained gel
- Protein expression in approximately 2 hours
- Transcription/translation components can be removed by affinity chromatography


PURExpress Disulfide Bond Enhancer: This proprietary blend of proteins and buffer components is designed to correctly fold target proteins with multiple disulfide bonds produced in PURExpress reactions or *E. coli* S30 extracts. Added at the beginning of a reaction, the components assist with the oxidation of cysteine thiols and correcting mis-oxidized substrates, increasing yield of soluble and functionally active protein.




PURExpress Δ RF123 Kit: Release factors are involved in termination of protein translation by recognizing the stop codons in an mRNA sequence. In a ribosome display experiment using PURExpress, lack of release factors could stabilize the ternary complex of mRNA-ribosome-nascent protein. As a result, the cDNA recovery could be higher. In this kit, the three release factors are supplied separately, allowing users to perform a protein synthesis reaction/ribosome display experiment with/ without release factors of their choice.

PURExpress Δ **(aa, tRNA) Kit:** The tRNA and amino acids are supplied separately in this kit, allowing users to run a protein synthesis reaction by adding modified amino acids and tRNA mixtures to the reaction.

E. coli Ribosome: The 70S *E. coli* Ribosome consists of a small subunit (30S) and a large subunit (50S). This preparation of ribosomes is highly active in NEB's PURExpress Protein Synthesis Kit (NEB #E6800), and can be used in ribosome structure and function studies, as a target for drug screening and as starting material for isolation of native ribosomal RNAs (5S, 16S, 23S). *E. coli* Ribosome is supplied as a 33.3 mg/ml solution.



Protein expression using the PURExpress In Vitro Protein Synthesis Kit. 25 µI reactions containing 250 ng template DNA and 20 units RNase Inhibitor were incubated at 37°C for 2 hours. 2.5 µI of each reaction was analyzed by SDS-PAGE using a 10–20% Tris-glycine gel. The red dot indicates the protein of interest. Marker M is the Protein Ladder



PURExpress Disulfide Bond Enhancer. (PDBE) promotes proper folding of active vtPA. Reactions were set up according to PURExpress specifications with the vtPA template DNA. After a 2 hour incubation at 37°C, 5 µl of each reaction was used in an activity assay and cleavage of the chromogenic substrate was monitored for one hour

PURExpress Kit Components

<i>IN VITRO</i> PROTEIN SYNTHESIS (E6800S)	Δ RIBOSOME (E3313S)	Δ RF 123 (E6850S)	Δ (aa, tRNA) (E6840S)
Solution A	Solution A	Solution A	Solution A (minus aa and tRNA)
Solution B	Factor Mix	Solution B (minus RF1, RF2 and RF3)	Solution B
Control (DHFR) template	Control (DHFR) template	Control (DHFR) template	Control (DHFR) template
	Control Ribosomes	RF1, RF2 and RF3	Amino Acid Mixture
			E. coli tRNA

Competent Cell Expression Strain Selection Chart

- Free of animal products
- T1 phage resistance (fhuA2)
- B strains are deficient in proteases Lon and OmpT
- A variety of convenient formats
- Bulk sales capabilities with custom packaging

NEB offers a wide selection of competent cell strains designed for the expression of a variety of proteins. Proteins with multiple disulfide bonds are correctly oxidized to significantly higher yields with SHuffle® strains. Tunable T7 expression is achieved with Lemo21(DE3), an ideal strain for difficult targets including membrane proteins. NiCo21(DE3) can be used for expression and purification of His-tagged proteins. NEB Express and T7 Express are

offered with varying levels of expression control. Several strains are available with added control by lac^p . Only NEB offers exceptional control of T7 expression from the lysY gene, which is ideal for proteins that are difficult to express or toxic to the cell. Each strain is provided with a protocol for optimal expression.

For more information see pages 229-233.

Expression Strains

CHARACTERISTICS	STRAIN	NEB #	SIZE
Versatile non-T7 expression strain Protease deficient	NEB Express Competent <i>E. coli*</i>	C2523H/I	20 x 0.05 ml/6 x 0.2 ml
* Control of IPTG induced expression from P_{bc} , P_{bc} and P_{bc} * Protease deficient	NEB Express I ^q Competent E. coli	C3037I	6 x 0.2 ml
Most popular T7 expression strain Protease deficient	T7 Express Competent E. coli	C2566H/I	20 x 0.05 ml/6 x 0.2 ml
T7 expression Protease deficient Better reduction of basal expression	T7 Express <i>lysY</i> Competent <i>E. coli</i>	C3010I	6 x 0.2 ml
To expression Protease deficient Highest level of expression control To expression	T7 Express <i>lysY/l^q</i> Competent <i>E. coli</i>	C3013I	6 x 0.2 ml
Enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm Protease deficient/B strain	SHuffle® Express Competent E. coli	C3028J	12 x 0.05 ml
Enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm T7 expression Protease deficient/B strain	SHuffle T7 Express Competent <i>E. coli</i>	C3029J	12 x 0.05 ml
To expression Protease deficient/B strain Tightly controlled expression of toxic proteins Enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm	SHuffle T7 Express <i>lysY</i> Competent <i>E. coli</i>	C3030J	12 x 0.05 ml
Enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm T7 expression/K12 strain	SHuffle T7 Competent <i>E. coli</i>	C3026J	12 x 0.05 ml
Routine expression for non-T7 Vectors Protease deficient	BL21 Competent E. coli	C2530H	20 x 0.05 ml
Routine T7 Expression Protease deficient	BL21(DE3) Competent E. coli	C2527H/I	20 x 0.05 ml/6 x 0.2 ml
Tunable T7 Expression for difficult targets Protease deficient	Lemo21(DE3) Competent E. coli	C2528J	12 x 0.05 ml
Expression and purification of His-tagged proteins Protease deficient	NiCo21(DE3) Competent E. coli	C2529H	20 x 0.05 ml

Note: Store Competent Cells at -80° C. Once thawed, do not refreeze. Storage at -20° C will result in a significant decrease in transformation efficiency. Cells lose efficiency whenever they are warmed above -80° C, even if they do not thaw.

^{*} NEB Express is the recommended strain for the pMAL Protein Fusion and Purification System.

Magnetic Affinity Matrices



Magnetic particles are ideally suited for applications involving high-throughput proteomic screening, small-scale protein isolation, immunomagnetic isolations or cell separation experiments. Magnetic affinity purification of tagged proteins, antigens, antibodies and nucleic acids can be done conveniently and quickly, with minimal time necessary for separation of the solid-phase from solution. In addition, immobilized substrates remain biologically active and can be eluted in small volumes or serve as ligands in subsequent pull-down or target interaction experiments involving DNA or proteins.

- No centrifugation required; matrix can be regenerated without loss of binding capacity.
- Minimal sample loss during pipetting because magnetic beads concentrate at the side of the tube instead of the bottom.

Immunomagnetic Isolation

Protein A Magnetic Beads

#S1425S 1 ml 178 €

Protein G Magnetic Beads

#S1430S 1 ml 198 €

- Small-scale purification or immunoprecipitation of IgG species
- No centrifugation required
- Regenerate matrix without binding capacity loss

Description: Protein A and Protein G Magnetic Beads are affinity matrices for the small-scale isolation and purification of immunoglobulins. They exhibit high affinity for subclasses of IgG from many mammalian species including human, rabbit and mouse. Truncated recombinant forms of Protein A and Protein G are covalently coupled to a nonporous superparamagnetic particle. These truncated proteins exhibit higher unit binding for the IgG Fc region and lower non-specific binding. The affinity for IgG varies with species and subclasses of IgG within a species (see table).

The proteins are coupled through a linkage that is stable and leak resistant over a wide pH range. This permits the immunomagnetic purification of IgG from ascites, serum or cell culture supernatants; after which, the matrix can be re-generated without loss of binding capacity. Protein A and Protein G magnetic beads can also be used to immunoprecipitate target proteins from crude cell lysates using selected primary antibodies. In addition, specific antibodies can be chemically crosslinked to the Protein A or Protein G coated surface to create a reusable immunprecipitation bead, avoiding co-elution of antibody with target antigen.

Supplied as a 1 ml suspension in PBS Buffer.

Support Matrix: $2 \mu m$ nonporous superparamagnetic microparticle

Binding Capacity: 1 ml of Protein A or Protein G Magnetic Beads binds > 280 µg of Human IgG.

References: References for enzyme properties and applications for these products can be found at www.neb.com.

Affinity of Protein A/G for IgG Species:

Species	Immunoglobulin	Binding to Protein A	Binding to Protein G
Human	lgG (normal)	++++	++++
	lgG1	++++	++++
	lgG2	++++	++++
	lgG3	_	++++
	lgG4	++++	++++
	IgM	_	_
	IgA	_	-
	IgE	_	-
Mouse	lgG1	+	++++
	lgG2a	++++	++++
	lgG2b	+++	+++
	lgG3	++	+++
Rat	lgG1	_	+
	lgG2a	_	++++
	lgG2b	_	++
	lgG2c	+	++
Goat	lgG	+/-	++
Rabbit	lgG	++++	+++
Sheep	IgG	+/-	++

Magnetic Beads

Streptavidin Magnetic Beads

#S1420S 5 ml 280 €

Hydrophilic Streptavidin Magnetic Beads

#S1421S 5 ml 280 €

Description: Streptavidin Magnetic Beads are 1 µm superparamagnetic particles covalently coupled to a highly pure form of streptavidin. The beads can be used to capture biotin labeled substrates including antigens, antibodies and nucleic acids. The strength of the biotin-streptavidin interaction, an association constant (Ka) of 10¹⁵ M⁻¹, coupled with the low non-specific binding of streptavidin, permits captured substrates to be useful as ligands in subsequent experiments including mRNA isolation and the capture of primary or secondary antibodies.

Hydrophilic Streptavidin Magnetic Beads have a modified bead surface generated by using a unique combination of blocking reagent to give better consistency, much lower non-specific binding and improved handling properties resulting in superior signal to noise in applications involving DNA immunoprecipitations.

Beads are supplied as a 4 mg/ml suspension in phosphate buffer (PBS).

Support Matrix: 1 µm nonporous superparamagnetic microparticle (#S1420) or 2 µM nonporous superparamagnetic microparticle (#S1421).

Binding Capacity: \$1420 will bind > 1000 pmol of free biotin per mg or > 500 pmol of ss-25 bp biotinylated oligonucleotide per mg. \$1421 will bind > 800 pmol of free biotin per mg or > 400 pmol of ss-25 bp biotinylated oligonucleotide per mg.

References: References for enzyme properties and applications for these products can be found at www.neb.com.

Immobilized Secondary Antibodies

Goat Anti-Rabbit IgG Magnetic Beads #S1432S 20 mg 178 €

Goat Anti-Mouse IgG Magnetic Beads #S1431S 20 mg 178 €

Goat Anti-Rat IgG Magnetic Beads #S1433S 20 mg 176 € Goat Anti-Rabbit IgG Magnetic Beads: An affinity matrix for the small-scale immunomagnetic separation and purification of rabbit IgGs. Goat Anti-Rabbit IgG is covalently coupled to a nonporous superparamagnetic particle. This secondary antibody binds the heavy chain of all rabbit IgG subclasses and is suitable for immunoassays that employ a rabbit IgG primary polyclonal antibody. Cell separations and sorting can be accomplished using rabbit IgG antibody to defined cell surface antigens

Goat Anti-Mouse IgG Magnetic Beads: An affinity matrix for the small-scale immunomagnetic separation and purification of mouse IgGs. Anti-Mouse IgG is covalently coupled to a nonporous superparamagnetic particle. This secondary antibody binds the heavy chain of mouse IgG and is suitable for immunoassays that employ a mouse IgG primary monoclonal antibody. Cell separations and sorting can be accomplished using a mouse IgG antibody to defined cell surface antigens.

Goat Anti-Rat IgG Magnetic Beads: An affinity matrix for the small-scale immunomagnetic separation and purification of rat IgG's. Anti-Rat IgG is covalently coupled to a 1µm nonporous superparamagnetic particle. This secondary antibody binds the Fc portion of all monoclonal rat IgG subclasses and is suitable for immunoassays that employ a rat IgG primary monoclonal antibody. Cell separations and sorting can be accomplished using a rat IgG antibody to defined cell surface antigens.

Supplied as a 1 ml suspension in PBS Buffer.

Support Matrix: 1 µm nonporous superparamagnetic microparticle.

Binding Capacity: 1 mg will bind 5 µg of lgG.

Magnetic Protein Purification

Chitin Magnetic Beads

#E8036S 5 ml 122 € #E8036L 25 ml 488 €

Amylose Magnetic Beads

#E8035S 25 mg 246 €

NEW

Ni-NTA Magnetic Beads

#\$1423\$ 1 ml 190 € #\$1423L 5 ml 899 €

Companion Product:

Anti-MBP Magnetic Beads #E8037S 10 mg264 € Chitin Magnetic Beads: An affinity matrix for the small-scale isolation of target proteins fused to a chitin binding domain (CBD). Chitin beads have been prepared having a magnetite core. This permits the magnetic isolation of CBD-fusion proteins from cell culture supernatants; after which, the matrix can be regenerated without loss of binding capacity. Immobilized fusion proteins can be used in subsequent experiments to capture target proteins from crude cell lysates that interact with the immobilized fusion protein.

Chitin Magnetic Beads are supplied as a 50:50 (v/v) suspension in water containing 20% ethanol.

Amylose Magnetic Beads: An affinity matrix for the small-scale isolation and purification of MBP-fusion proteins. Amylose is covalently coupled to a superparamagnetic particle through a linkage that is stable and leak resistant over a wide pH range. This permits the isolation of MBP-fusion proteins from cell culture supernatants. Immobilized fusion proteins can be used in subsequent experiments to capture target proteins from crude cell lysates that interact with the immobilized MBP-fusion protein.

Amylose Magnetic Beads are supplied as a 25 mg/ml suspension in water containing 20% ethanol.

Ni-NTA Magnetic Beads: An affinity matrix for the small-scale isolation and purification of polyhistidine-tagged (His-tagged) fusion proteins in manual or

automated formats. Immobilized Metal Affinity Chromatography (IMAC) purifications employing Ni-NTA (nickel-nitrilotriacetic acid) magnetic beads can be performed under native or denaturing conditions, which permit efficient binding and purification of insoluble proteins, proteins that aggregate in inclusion bodies, or proteins that possess a tertiary structure that sequester the polyhistidine affinity tag. Immobilized His-tagged proteins can be used in subsequent experiments to pull-down proteins that may interact with the immobilized protein.

Support Matrix: Chitin Magnetic Beads are approximately 10–100 μm superparamagnetic microparticle.

Amylose Magnetic Beads are $1-10~\mu m$ superparamagnetic microparticle.

Ni-NTA Magnetic Beads are spherical, agarose-based super-paramagnetic microparticles ranging in size from $20{\text -}100~\mu m$.

Binding Capacity: 1 ml of Chitin Magnetic Beads will bind 2 mg of CBD fusion protein.

1 mg of Amylose Magnetic Beads will bind 10 μg of MBP-fusion protein.

Ni-NTA: Varies with target, typically ≥ 7.5 mg Histagged fusion protein/ml bed volume.

Magnetic Racks

6-Tube Magnetic Separation Rack

#S1506S 6 tubes (1.5 ml) 206 €

12-Tube Magnetic Separation Rack

#S1509S 12 tubes (1.5 ml) 326 €

50 ml Magnetic Separation Rack

#S1507S 4 tubes (50 ml) 292 €

96-Well Microtiter Plate Magnetic Separation Rack #S1511S 96-well 489 €

Description: The Magnetic Separation Racks are designed to be used for small-scale separations using magnetic particles. The magnets are located on the sides of the racks, resulting in minimal sample loss during pipetting.

Magnets: Neodymium rare earth permanent magnets.

Dimensions: 2-Tube Rack (1" x 2" x 1¾"), 6-Tube Rack (3" x 2" x 1¾"), 12-Tube Rack (5½" x 2" x 1¾"), 50 ml Rack (3½" x 4¼" x 3½") and 96-Well Microtiter Plate Rack (5½" x 1¼" x 3¾").

SNAP-Capture Magnetic Beads

See page 283 for more information.





The Effect of Dams on River Systems

Dams have substantially contributed to urban development and industrialization. Historically, they were considered a sign of a thriving civilization, providing hydroelectricity, irrigation, flood control and a supply of water to nearby towns. However, they come with an environmental cost, and aging dams have limited capability to meet the energy needs of a modern world.

Essentially, dams obstruct the flow of rivers — there are 84,000 dams in the USA that are three feet high or greater, blocking 600,000 miles of river. They also have a profound effect on the river ecosystem. Upstream of a dam, the water is stagnant. Sediment, rocks and wood that would normally flow downstream and shape the landscape, build up and affect the coastline ecosystem. Weeds and algae proliferate and reservoir depth results in colder water, which in turn reduces the amount of oxygen and changes the nutrient composition for marine life. The dam effectively turns the river into a lake.

The Elwha and Glines Canyon River Dams were built on the Elwha River, which mostly lies within Olympic National Park in Washington, USA, in 1913 and 1927, respectively. They were built in order to harness the power of the river to generate electricity. At the time, they energized economic growth in the region, but 100 years later they provided only a minimal amount of the electricity needs of the district. The cost of keeping the dams outweighed the benefits, and they were removed between 2011 and 2013. The removal of the dams was two decades in planning. This gave scientists the opportunity to document the surrounding ecosystem before and after the demolition, so that the feasibility of other dam removal projects could be assessed.

The two dams blocked wild salmon runs, and this had a profound effect, not only on salmon numbers but also on the surrounding ecosystem upstream of the dam. Salmon are a keystone species. They thrive and grow on marine nutrients in the sea, and when they return to the river and travel upstream to spawn, they transport valuable nutrients to the wildlife. If they die or are killed by other animals in the ecosystem, such as bears, otter and eagles, these nutrients are distributed into the surrounding vegetation. When the salmon disappeared, so did the animals that relied on them for survival.

The Elwha River was choked with 33 tons of sediment — dismantling the dams reshaped 13 miles of the Elwha River and expanded the river delta at the Pacific Ocean. Wood and sediment reshaped the shore — beaches, kelp beds and eelgrass beds flourished. Willows were soon thriving by the river. Moss grew and created a microclimate for other plants. The dams had prevented Pacific Salmon from reaching 90% of their habitat, and now they were returning to parts of the river that had not seen salmon in 100 years.

This massive project has provided a wealth of information for ecologists. Not only do we better understand the effects of putting dams in place, but we can now predict the outcome of future dam removal projects. Not every dam is obsolete, but as the outdated and unsafe dams are coming down, the flow of rivers and valuable ecosystems are being restored.

Competent Cells



NEB has a competent cell strain for your needs.

Choose the right cells for your cloning and protein expression applications from NEB's portfolio of high efficiency competent cell strains.

For cloning experiments, choose from a variety of formats, including chemical and electrocompetent. These *E. coli* strains are T1 phage resistant and are Endonucleasel-deficient for high-quality plasmid preparations. Additionally, all competent cells from NEB are free of animal products.

NEB also offers a wide variety of competent cell strains ideal for many protein expression applications. These strains address the needs of difficult protein expression control, toxic protein expression and cytoplasmic disulfide bond formation. NEB Express, T7 Express and SHuffle® strains are available with varying levels of control. I^q strains feature added control from increased supply of Lac repressor ($Iacl^9$). Only NEB offers the exceptional control of expression from the IysY gene that reduces basal expression from T7 strains without inhibiting IPTG-induced expression. Our Lemo21(DE3) strain features tunable T7 expression for difficult targets. Each strain is provided with a detailed protocol for optimal expression.





Featured Products

NEB Cloning Competent E. coli Sampler

NEB Stable Competent E. coli

229 BL21 Competent E. coli

232 SHuffle Express Competent E. coli

Featured Tools & Resources

361 Enhancing Transformation Efficiency

Protein Expression with T7 Express Strains

Troubleshooting Guide for Cloning



Visit www.neb.com to find additional online tools, including our Competitor Cross-reference Tool for comparing NEB strains to other commercially available strains.



Strain Properties Chart	224	Protein Expression Strains	
Convenient Formats	225	BL21 Competent <i>E. coli</i>	229
Cloning Strains NEB Cloning Competent <i>E. coli</i> Sampler NEB Turbo Competent <i>E. coli</i>	226	BL21(DE3) Competent <i>E. coli</i> Lemo21(DE3) Competent <i>E. coli</i> NiCo21(DE3) Competent <i>E. coli</i> NEB Express Competent <i>E. coli</i> (High Efficiency)	229 229 230 230
(High Efficiency) Electrocompetent	226 226	NEB Express I ^q Competent <i>E. coli</i> (High Efficiency)	230
NEB 10-beta Competent <i>E. coli</i> (High Efficiency) Electrocompetent	227 227	T7 Express Competent <i>E. coli</i> (High Efficiency)	231
NEB 5-alpha Competent <i>E. coli</i> (High Efficiency)	227	T7 Express <i>lysY</i> Competent <i>E. coli</i> (High Efficiency)	231
(Subcloning Efficiency) Electrocompetent	227 227	T7 Express <i>lysY/I^q</i> Competent <i>E. coli</i> (High Efficiency)	231
NEB 5-alpha F´I ^q Competent <i>E. coli</i> (High Efficiency)	227	SHuffle Express Competent <i>E. coli</i> SHuffle T7 Express Competent <i>E. coli</i>	232 232
NEB Stable Competent <i>E. coli</i> (High Efficiency)	228	SHuffle T7 Competent <i>E. coli</i> SHuffle T7 Express <i>lysY</i> Competent <i>E. coli</i> How do I express my protein in SHuffle?	233 233 233
dam ⁻ /dcm ⁻ Competent <i>E. coli</i>	228	Tion do roxproce my protein in citatile:	

Strain Properties

There are many properties to consider when choosing a strain for your experiments.

Requirements such as plasmid preparation, blue/white screening, proper disulfide bond formation and fast colony growth necessitate specific strain choices. The following selection chart highlights the characteristics

of NEB's strains to help select the optimal strain for a particular experiment.

CAUTION: Chemically Competent *E. coli* contain DMSO, a hazardous material. Review the MSDS before handling.

Cloning Strain Properties

STRAIN PROPERTIES	TRANSFORMATION EFFICIENCY (cfu/µg) ⁽¹⁾			BLUE/WHITE SCREENING	laciª	lysY	COLONIES VISIBLE AFTER 6.5 HRS.	endA ^{¬(2)}	PROTEASE DEFICIENT ⁽³⁾	F	recA ⁻	T7 RNA Polymerase	CYTOPLASMIC Disulfide Bond Formation(4)	DRUG Resistance ⁽⁵⁾
NEB Turbo	1–3 x 10 ⁹	K12	•	•	•	_	•	•	_	•	_	_	-	nit
NEB 5-alpha	1–3 x 10 ⁹	K12	•	•	-	-	_	•	_	-	•	_	-	none
NEB 5-alpha F´ I ^q	1–3 x 10 ⁹	K12	•	•	•	_	_	•	_	•	•	_	-	tet
NEB 10-beta	1–3 x 10 ⁹	K12	•	•	_	-	-	•	_	_	•	_	-	str
dam-/dcm-	1–3 x 10 ⁶	K12	•	_	_	_	_	•	_	_	_	_	-	cam, str, nit
NEB Stable	1–3 x 10 ⁹	K12	•	•	•	_	-	•	-	•	•	_	-	tet, str

Protein Expression Strain Properties

NEB Express	0.6-1 x 10 ⁹	В	•	_	_	_	-	•	•	_	_	_	_	nit
NEB Express I ^q	0.6-1 x 10 ⁹	В	•	_	•	_	_	•	•	-	_	_	_	cam, nit
T7 Express	0.6-1 x 10 ⁹	В	•	-	-	-	-	•	•	-	-	•	_	nit
T7 Express <i>lysY</i>	0.6-1 x 10 ⁹	В	•	_	_	•	_	•	•	-	_	•	_	cam, nit
T7 Express <i>lysY/I</i> ^q	0.6-1 x 10 ⁹	В	•	-	•	•	-	•	•	-	_	•	_	cam, nit
SHuffle® Express	1 x 10 ⁷	В	•	-	•	_	-	•	•	_	_	_	•	spec ⁽⁶⁾ , nit
SHuffle T7 Express	1 x 10 ⁷	В	•	-	•	-	-	•	•	_	_	•	•	spec ⁽⁶⁾ , nit
SHuffle T7 Express <i>lysY</i>	1 x 10 ⁷	В	•	_	•	•	-	•	•	-	_	•	•	cam, spec ⁽⁶⁾ , nit
SHuffle T7	1 x 10 ⁶	K12	•	-	•	_	-	_	-	•	_	•	•	str, spec, nit
BL21	1–5 x 10 ⁷	В	•	_	_	_	_	_	•	-	_	_	_	none
BL21(DE3)	1–5 x 10 ⁷	В	•	-	_	_	-	_	•	-	_	•	-	none
Lemo21(DE3)	1–3 x 10 ⁷	В	•	_	_	•	-	-	•	-	_	•	_	cam
NiCo21(DE3)	1–5 x 10 ⁷	В	•	-	-	_	-	-	•	_	-	•	_	none

- (1) TE are given high-quality for high efficiency chemically competent strains. TE for electrocompetent strains are 1-4 x 10 10 cfu/µg. TE for subcloning strains are >1 x 10 6 cfu/µg.
- (2) Important for high-quality plasmid preparation.
- (3) Lacks Lon and OmpT protease activity.
- (4) Constitutively expresses a chromosomal copy of the disulfide bond isomerase DsbC.
- (5) nit = nitrofurantoin, tet = tetracycline, cam = chloramphenicol, str = streptomycin, spec = spectinomycin
- (6) Resistance to low levels of streptomycin may be observed.



For help with choosing the right competent cell strain, try NEBcloner at NEBcloner.neb.com

Convenient Formats

NEB provides these superior competent cell strains in several formats for your convenience. Most are available as 50 µl single-use transformation tubes and many are available in larger, 200 µl tubes for multiple simultaneous reactions. Our most popular cloning strains are

available as electrocompetent cells. NEB 5-alpha is also available in a lower efficiency, subcloning format for substantial value, as well as in a 96-well plate, 384-well plate and striptube formats.

Cloning Formats

FORMATS	50 µI Single-USE (H,J-Formats)	200 µl TUBES (I-FORMAT)	ELECTROCOMPETENT (K-FORMAT)	SUBCLONING (F-FORMAT)	96-WELL PLATE (P-FORMAT)	384-WELL PLATE (R-FORMAT)	8-TUBE STRIPS (U-FORMAT)	OUTGROWTH MEDIUM & CONTROL PLASMID INCLUDED
NEB Turbo	•	•	•	_	_			•
NEB 5-alpha	•	•	•	•	•	•	•	•
NEB 5-alpha F´ Iª	•	•	_	_	-			•
NEB 10-beta	•	•	•	-	-			•
dam / dcm	•	•	_	_				•
NEB Stable	•	•	_	-	-			•

NEB Express	•	•	_	-	-		•
NEB Express I ^q	_	•	_	-	_		_
T7 Express	•	•	-	-	-		•
T7 Express <i>lysY</i>	_	•	_	_	_		_
T7 Express /ysY/I ^q	-	•	_	-	-		-
SHuffle® Express	•	_	_	_	_		_
SHuffle T7 Express	•	-	-	-	-		_
SHuffle T7 Express <i>lysY</i>	•	_	_	_	_		_
SHuffle T7	•	-	-	-	-		-
BL21	•	_	_	_	_		•
BL21(DE3)	•	•	_	-	-		•
Lemo21(DE3)	•	_	_	_	_		• (1)
NiCo21(DE3)	•	-	-	-	-		•

⁽¹⁾ Rhamnose solution is provided instead of SOC, control plasmid is included.



Using another competent cell strain? Try our competitor cross reference tool to find out which NEB strain is compatible.



NEB Cloning Competent E. coli Sampler

#C1010S 8 tubes 120 €

Companion Product:

SOC Outgrowth Medium #B9020S 4 x 25 ml 83 €

- Outgrowth medium and control plasmid included
- Value pricing
- Free of animal products

Description: A sample pack of four cloning strains of *E. coli* suitable for high efficiency transformation.

Please refer to the individual datacards for each reagent's recommended use and storage conditions.


Transformation Efficiency:

1-3 x 109 cfu/µg pUC19 DNA

Resistance: T1 phage (fhuA2)

The Sampler Includes:

- NEB Turbo Competent E. coli (High Efficiency)
- NEB 5-alpha Competent E. coli (High Efficiency)
- NEB 10-beta Competent E. coli (High Efficiency)
- NEB Stable Competent E. coli (High Efficiency)
- SOC Outgrowth Medium
- NEB 10-beta/Stable Outgrowth Medium
- pUC19 Vector



Benefit from high transformation efficiencies: Transformation efficiences were compared using manufacturers' recommended protocols. Values shown are the average of triplicate experiments.

NEB Turbo Competent E. coli

NEB Turbo Competent *E. coli* (High Efficiency)

#C2984H 20 x 0.05 ml 261 € #C2984I 6 x 0.2 ml 198 €

NEB Turbo Electrocompetent *E. coli* #C2986K 6 x 0.1 ml 244 €

Companion Product:

SOC Outgrowth Medium #B9020S 4 x 25 ml 83 €

- Tight expression control (lacl^q)
- Colonies visible after 6.5 hours
- Isolate DNA after 4 hrs growth
- 5 minute transformation protocol with Amp^R plasmids
- Clone toxic genes
- Free of animal products

Description: *E. coli* cells featuring fast colony growth (6.5 hours) and tight expression control.

Genotype: F´ $proA^*B^*$ $lacl^*\Delta lacZM15 / fhuA2 <math>\Delta (lacproAB)$ glnV galK16 galE15 R(zgb-210::Tn10)Tet S endA1 thi-1 $\Delta (hsdS-mcrB)5$

Features:

- · Suitable for blue/white screening
- Activity of nonspecific endonuclease I (endA1) eliminated for highest quality plasmid preparations


Transformation Efficiency:

High Efficiency: 1-3 x 109 cfu/µg pUC19 DNA

Electrocompetent: > 1 x 1010 cfu/µg pUC19 DNA

Resistance: T1 phage (fhuA2), Nit


Sensitivity: Amp, Cam, Kan, Spec, Str, Tet



Miniprep DNA can be prepared from a single overnight colony after inoculation and only 3 hours growth. DNA yield doubles after an additional hour of growth.

Reagents Supplied:

SOC Outgrowth Medium pUC19 Control DNA



NEB Turbo Transformation: With NEB Turbo, colonies are visible after only 8 hours. Ligation products were transformed into 50 µl of NEB Turbo Competent E. coli and plated on LB/Amp. Plates were incubated for 8 hours, 10 hours and 12 hours at 37°C. NEB Turbo features fast colony growth and blue/white selection to simplify cloning experiments.



NEB 10-beta Competent E. coli

NEB 10-beta Competent *E. coli* (High Efficiency)

#C3019H 20 x 0.05 ml 231 € #C3019I 6 x 0.2 ml 179 €

NEB 10-beta Electrocompetent *E. coli* #C3020K 6 x 0.1 ml 214 €

Companion Product:

NEB 10-beta/Stable Outgrowth Medium #B9035S 4 x 25 ml 82 €

- Clone large plasmids and BACs
- DH10B derivative
- Free of animal products

Description: A DH10B derivative suitable for a wide range of applications, including large plasmid and BAC cloning.

Genotype: Δ(ara-leu) 7697 araD139 fhuA ΔlacX74 galK16 galE15 e14- φ80dlacZΔM15 recA1 relA1 endA1 nupG rpsL (Str[®]) rph spoT1 Δ(mrr-hsdRMS-mcrBC)

Features:

- Efficient transformation of methylated DNA derived from eukaryotic sources or unmethylated DNA derived from PCR, cDNA and many other sources
- · Suitable for blue/white screening without IPTG
- Activity of nonspecific endonuclease I (endA1) eliminated for highest quality plasmid preparations
- Reduced recombination of cloned DNA (recA1)

Transformation Efficiency:

High Efficiency: 1-3 x 109 cfu/μg pUC19 DNA


Electrocompetent: > 2 x 1010 cfu/µg pUC19 DNA

Resistance: T1 phage (fhuA2), Str

Sensitivity: Amp, Cam, Kan, Nit, Spec, Tet

Reagents Supplied:

NEB 10-beta/Stable Outgrowth Medium pUC19 Control DNA



Effect of Plasmid Size on Transformation Efficiency:
NEB 10-beta chemically competent cells are more efficiently
transformed with large plasmids than NEB 5-alpha cells. The
difference in TE between the two cell lines increases with the size
of the plasmid being transformed.

NEB 5-alpha Competent E. coli

NEB 5-alpha Competent *E. coli* (High Efficiency)

#C2987H 20 x 0.05 ml197 €
#C2987I 6 x 0.2 ml153 €
#C2987P 1 x 96 well plate493 €
#C2987R 1 x 384 well plate ...1096 €
#C2987U 96 x 50 μl/tube742 €

NEB 5-alpha Competent *E. coli* (Subcloning Efficiency)*

#C2988J 6 x 0.4 ml71 €

NEB 5-alpha Electrocompetent *E. coli* #C2989K 6 x 0.1 ml 184 €

Companion Product:

- DH5α derivative
- Free of animal products

Description: A DH5 α derivative and versatile *E. coli* cloning strain.

Genotype: fhuA2Δ(argF-lacZ)U169 phoA glnV44 φ80Δ(lacZ)M15 gyrA96 recA1 relA1 endA1 thi-1 hsdR17

Features:

- Efficient transformation of unmethylated DNA derived from PCR, cDNA and many other sources (hsdR)
- Suitable for blue/white screening
- Activity of nonspecific endonuclease I (endA1) eliminated for highest quality plasmid preparations
- Reduced recombination of cloned DNA (recA1)

Transformation Efficiency:

High Efficiency: 1-3 x 109 cfu/μg pUC19 DNA

Subcloning Efficiency: > 1 x 106 cfu/µg pUC19 DNA

Electrocompetent: > 1 x 10¹⁰ cfu/μg pUC19 DNA

Resistance: T1 phage (fhuA2)

Sensitivity: Amp, Cam, Kan, Nit, Spec, Str, Tet

Reagents Supplied:

SOC Outgrowth Medium pUC19 Control DNA

NEB 5-alpha F'I^q Competent E. coli

NEB 5-alpha F' I q Competent E. coli (High Efficiency)

#C2992H 20 x 0.05 ml 197 € #C2992I 6 x 0.2 ml 154 €

Companion Product::

- Tight expression control (lacl^q)
- F´Strain with extremely high TE
- DH5α derivative
- Free of animal products

Description: An F´ *E. coli* strain with extremely high transformation efficiency suitable for toxic gene cloning.

Genotype: F´ $proA^+B^+$ $lacl^a$ $\Delta(lacZ)M15$ zzf::Tn10 (Tet^a) / $fhuA2\Delta(argF-lacZ)U169$ phoA glnV44 $\phi80$ $\Delta(lacZ)M15$ qyrA96 recA1 endA1 thi-1 hsdR17

Features:

- Efficient transformation of unmethylated DNA derived from PCR, cDNA and many other sources (hsdR)
- · Suitable for blue/white screening
- Activity of nonspecific endonuclease I (endA1) eliminated for highest quality plasmid preparations
- · Reduced recombination of cloned DNA (recA1)
- Suitable for propagation of M13 clones

Transformation Efficiency:

High Efficiency: 1–3 x 109 cfu/μg

Resistance: T1 phage (fhuA2), Tet

Sensitivity: Amp, Cam, Kan, Nit, Spec, Str

Reagents Supplied: SOC Outgrowth Medium

SOC Outgrowth Mediur pUC19 Control DNA

^{*} NEB 5-alpha Competent *E. coli* (Subcloning Efficiency) is not supplied with SOC Outgrowth Medium or pUC19 Control DNA.

NEB Stable Competent E. coli

NEB Stable Competent *E. coli* (High Efficiency)

#C3040H 20 x 0.05 ml 318 € #C3040I 6 x 0.2 ml 245 €

Companion Product

NEB 10-beta/Stable Outgrowth Medium #B9035S 4 x 25 ml82

- T1 phage resistance (fhuA)
- Free of animal products
- Carries endA mutation (isolated plasmids are free of Endol)

Description: Chemically competent *E. coli* cells suitable for high efficiency transformation and isolation of plasmid clones containing repeat elements.

Genotype: F´ $proA^*B^*$ $lacF^*\Delta(lacZ)M15$ zzf::Tn10 $(Tet^*)/\Delta(ara-leu)$ 7697 araD139 fhuA $\Delta lacX74$ galK16 galE15 e14- $\phi80dlacZ\Delta M15$ recA1 relA1 endA1 nupG rpsL (Str^*) rph spoT1 $\Delta(mrr-hsdRMS-mcrBC)$

Features:

- Activity of nonspecific endonuclease I (endA1) abolished for highest quality plasmid preparations
- Rapid growth recA strain

Applications:

- · Cloning unstable inserts
- · Isolating and propagating retroviral/lentiviral clones
- Compatible with Gibson Assembly® Reactions, as well as ligation reactions

Transformation Efficiency:

1-5 x 108 cfu/µg pUC19 DNA (NEB #C3040H)

> 1 x 108 cfu/µg pUC19 DNA (NEB #C3040I)

Resistance: T1 phage (fhuA), Str, Tet Sensitivity: Amp, Cam, Kan, Nit, Spec

Reagents Supplied:

NEB 10-beta/Stable Outgrowth Medium pUC19 Control DNA

dam / dcm Competent E. coli

#C2925H 20 x 0.05 ml 244 € #C2925I 6 x 0.2 ml 188 €

Companion Product:

SOC Outgrowth Medium #B9020S

4 x 25 ml 83 €

- Isolate plasmids free of Dam and Dcm methylation
- Free of animal products

Description: Methyltransferase deficient *E. coli* cells suitable for growth of plasmids free of Dam and Dcm methylation.

Genotype: ara-14 leuB6 fhuA31 lacY1 tsx78 glnV44 galK2 galT22 mcrA dcm-6 hisG4 rfbD1 R(zgb210::Tn10) Tet^s endA1 rspL136 (Str^R) dam13::Tn9 (Cam^R) xylA-5 mtl-1 thi-1 mcrB1 hsdR2

Features:

- Allows for propagation of plasmids free of Dam and Dcm methylation
- Activity of nonspecific endonuclease I (endA1) abolished for highest quality plasmid preparations

Transformation Efficiency:

1-3 x 106 cfu/µg pUC19 DNA

Resistance: T1 phage (fhuA31), Cam, Nit, Str

Sensitivity: Amp, Kan, Spec, Tet


Reagents Supplied: SOC Outgrowth Medium pUC19 Control DNA

NEB Tube Opener

#C1008S

2 units 20 €

Description: Use to open a variety of microcentrifuge tubes. Can be used for snap-on caps or screw-on caps.



BL21 Competent E. coli

#C2530H

20 x 0.05 ml 194 €

Companion Product:

SOC Outgrowth Medium #B9020S

4 x 25 ml83 €

- Ideal for P_{lac}, P_{tac}, P_{trc}, ParaBAD expression vectors
- Resistance to phage T1 (fhuA2)
- Protease deficient
- Free of animal products

Description: Widely used non-T7 expression *E. coli* strain. Suitable for transformation and protein expression. This strain does not express the T7 RNA Polymerase.

Genotype: fhuA2 [lon] ompT gal [dcm] ∆hsdS

Features:

. Deficient in proteases Lon and OmpT

Transformation Efficiency:

1-5 x 107 cfu/µg pUC19 DNA

Resistance: T1 phage (*fhuA2*)

Sensitivity: Amp, Cam, Kan, Nit, Spec, Str, Tet

Reagents Supplied: SOC Outgrowth Medium pUC19 Control DNA

BL21(DE3) Competent E. coli

#C2527H #C2527I 20 x 0.05 ml 194 € 6 x 0.2 ml 151 €

Companion Product:

SOC Outgrowth Medium #B9020S

4 x 25 ml 83 €

- Routine T7 expression
- Free of animal products
- Protease deficient B strain

Description: Widely used T7 expression *E. coli* strain.

Genotype: fhuA2 [lon] ompT gal (λ DE3) [dcm] $\Delta hsdS$ λ DE3 = λ sBamHlo $\Delta EcoRl-B$ int::(lacl::PlacUV5::T7 qene1) i21 $\Delta nin5$

Features:

- · Deficient in proteases Lon and OmpT
- Resistant to phage T1 (fhuA2)

Transformation Efficiency:

 $1-5 \times 10^7$ cfu/µg pUC19 DNA

Resistance: T1 phage (fhuA2)

Sensitivity: Amp, Cam, Kan, Nit, Spec, Str, Tet

Reagents Supplied:SOC Outgrowth Medium pUC19 Control DNA

Lemo21(DE3) Competent E. coli

#C2528J

12 x 0.05 ml 224 €

- Expression of difficult targets
- Membrane protein expression
- Ideal for periplasmic expression
- Expression of toxic proteins
- Proteins with solubility issues

Description: Lemo21(DE3) Competent *E. coli* is a tunable T7 expression strain designed for the expression of challenging proteins. A derivative of BL21(DE3), Lemo21(DE3) offers the host features of this popular expression strain, with the added benefit of being able to control expression levels by varying the level of T7 lysocyme (*lysY*), the natural inhibitor of T7 RNA Polymerase. The fine control of expression makes Lemo21(DE3) ideal for membrane proteins, toxic proteins, secreted proteins and proteins prone to insoluble expression.

Genotype: fhuA2 [lon] ompT gal (λ DE3) [dcm] $\Delta hsdS$ /pLemo(Cam[®]) λ $DE3 = \lambda$ sBamHlo $\Delta EcoRI-B$ int::(lacl::PlacUV5::T7 gene1) i21 $\Delta nin5$ pLemo = pACYC184-PrhaBAD-Iys Y

Features:

- Enhanced BL21(DE3) derivative
- · Fine control of expression
- Greatest range of expression of any T7 strain (0-2,000 µM rhamnose)
- Potential elimination of inclusion body formation

Transformation Efficiency:

High Efficiency: 1–3 x 10⁷ cfu/μg pUC19 DNA

Resistance: T1 phage (*fhuA2*), Cam Sensitivity: Amp, Kan, Nit, Spec, Str, Tet

Reagents Supplied: L-rhamnose solution pUC19 control DNA

NiCo21(DE3) Competent E. coli

#C2529H 20 x 0.05 ml 373 €

Companion Product:

SOC Outgrowth Medium #B9020S 4 x 25 ml83 €

- Superior alternative to BL21(DE3) for routine protein expression
- Improved purity of target proteins isolated by IMAC
- Free of animal products

Description: Poly-histidine tagged recombinant proteins that are isolated by immobilized metal affinity chromatography (IMAC) are often contaminated with significant amounts of endogenous *E. coli* metal binding proteins. The protein expression strain NiCo21(DE3) has been engineered to minimize *E. coli* protein contamination of IMAC fractions: GImS is mutated to eliminate binding to IMAC resins and three other proteins (SlyD, ArnA and Can) are tagged to enable rapid removal by chitin affinity chromatography.

Genotype: can::CBD fhuA2 [lon] ompT gal (λ DE3) [dcm] arnA::CBD slyD::CBD glmS6Ala Δ hsdS λ DE3 = λ sBamHlo Δ EcoRl-B int::(lacl::PlacUV5::T7 gene1) i21 Δ nin5

Features:

- · Identical growth characteristics as BL21(DE3)
- · Deficient in proteases Lon and OmpT

Transformation Efficiency:

 $1-5 \times 10^7$ cfu/µg pUC19 DNA

Resistance: T1 phage (fhuA2)

Sensitivity: Amp, Cam, Kan, Nit, Spec, Str, Tet

Reagents Supplied: SOC Outgrowth Medium pUC19 Control DNA

NEB Express Competent E. coli

NEB Express Competent *E. coli* (High Efficiency)

#C2523H 20 x 0.05 ml 194 € #C2523I 6 x 0.2 ml 151 €

Companion Product:

SOC Outgrowth Medium #B9020S 4 x 25 ml 83 €

- Enhanced BL21 derivative ideal for P_{lac}, P_{tac}, P_{tac}
 P_{tre} expression vectors
- Fast growth from colonies
- Free of animal products
- Protease deficient

Description: A versatile non-T7 expression *E. coli* strain. NEB Express is the recommended host strain for pMAL protein fusion and purification system.

Genotype: fhuA2 [lon] ompT gal sulA11 $R(mcr-73::miniTn10--Tet^s)2$ [dcm] $R(zgb-210::Tn10--Tet^s)$ $endA/\Delta(mcrC-mrr)114::IS10$

Features:

- · Deficient in proteases Lon and OmpT
- · Does not restrict methylated DNA

Transformation Efficiency:

High Efficiency: 0.6-1 x 10° cfu/µg pUC19 DNA

Resistance: T1 phage (fhuA2), Nit

Sensitivity: Amp, Cam, Kan, Tet, Spec, Str

Reagents Supplied: SOC Outgrowth Medium pUC19 Control DNA

NEB Express I^q Competent E. coli

NEB Express I^q Competent *E. coli* (High Efficiency)

#C3037I 6 x 0.2 ml 158 €

Companion Product:

- Enhanced BL21 derivative ideal for P_{lac} , P_{tac} , P_{trc} , P_{75} expression vectors
- Better control of IPTG induced expression with non-T7 plasmids
- Fast growth from colonies
- lacl^q reduces basal expression
- Protease deficient
- Free of animal products

Description: *E. coli* cells featuring control of IPTG induced expression with non-T7 plasmids.

Genotype: MiniF lacF (Cam^R) / fhuA2 [lon] ompT gal sulA11 R(mcr-73::miniTn10--Tet^S)2 [dcm] R(zgb-210::Tn10--Tet^S) endA/Δ(mcrC-mrr)114::IS10

Features:

- Deficient in proteases Lon and OmpT
- Does not restrict methylated DNA
- Ideal for controlled protein expression from pUC19 and pUC19 derivatives

Transformation Efficiency:

High Efficiency: 0.6-1 x 109 cfu/µg pUC19 DNA

Resistance: T1 phage (*thuA2*), Cam, Nit Sensitivity: Amp, Kan, Spec, Str, Tet

T7 Express Competent E. coli

T7 Express Competent E. coli (High Efficiency)

#C2566H 20 x 0.05 ml 194 € 6 x 0.2 ml 151 € #C2566I

Companion Product:

SOC Outgrowth Medium 4 x 25 ml83 € #B9020S

- Enhanced BL21 derivative
- Popular T7 expression strain
- Fast growth from colonies
- Free of animal products

Description: Enhanced BL21 E. coli derivative for T7 expression

Genotype: fhuA2 lacZ::T7 gene1 [lon] ompT gal sulA11 R(mcr-73::miniTn10--Tet^S)2 [dcm] R(zgb-210::Tn10--Tet^s) endA1 ∆(mcrC-mrr)114::IS10

Features:

- T7 RNA Polymerase in the lac operon no λ prophage
- Deficient in proteases Lon and OmpT
- · Does not restrict methylated DNA

Transformation Efficiency:

High Efficiency: 0.6-1 x 109 cfu/µg pUC19 DNA

Resistance: T1 phage (fhuA2), Nit

Sensitivity: Amp, Cam, Kan, Spec, Str, Tet

Reagents Supplied: SOC Outgrowth Medium pUC19 Control DNA

T7 Express *lysY* Competent *E. coli*

T7 Express lysY Competent E. coli (High Efficiency)

#C3010I 6 x 0.2 ml 158 €

Companion Product:

SOC Outgrowth Medium #B9020S

4 x 25 ml 83 €

- Enhanced BL21 derivative
- T7 Lysozyme for expression control
- Clone toxic genes
- Fast growth from colonies
- Free of animal products

Description: Enhanced BL21 E. coli derivative for T7 expression with enhanced reduction of basal expression.

Genotype: MiniF lysY (Cam^R) / fhuA2 lacZ::T7 gene1 [lon] ompT gal sulA11 R(mcr-73::miniTn10--Tet^s)2 [dcm] $R(zgb-210::Tn10--Tet^s)$ endA1 $\Delta(mcrC-mrr)$ 114::IS10

Features:

- T7 RNA Polymerase in the lac operon no λ prophage
- Control of T7 RNA Polymerase by T7 lysozyme allows potentially toxic genes to be expressed
- LysY is a variant of T7 lysozyme lacking amidase activity, thus cells are not susceptible to lysis during induction
- Deficient in proteases Lon and OmpT
- · Does not restrict methylated DNA
- · No Cam requirement

Transformation Efficiency:

High Efficiency: 0.6-1 x 109 cfu/µg pUC19 DNA

Resistance: T1 phage (fhuA2), Cam, Nit Sensitivity: Amp, Kan, Spec, Str, Tet

T7 Express $lysY/I^q$ Competent E. coli

T7 Express lysY/Iq Competent E. coli (High Efficiency)

#C3013I 6 x 0.2 ml 158 €

Companion Product:

SOC Outgrowth Medium #B9020S 4 x 25 ml 83 €

- Enhanced BL21 derivative
- Tight control of expression (lacl^q)
- Highest level of expression control
- Clone toxic genes
- Fast growth from colonies
- Free of animal products

Description: Enhanced BL21 E. coli derivative with highest level of T7 expression control.

Genotype: MiniF lysY laclq(CamR) / fhuA2 lacZ::T7 gene1 [lon] ompT gal sulA11 R(mcr-73::miniTn10--Tet s)2 [dcm] R(zgb-210::Tn10--Tet s) endA1 Δ (mcrCmrr)114::IS10

Features:

- T7 RNA Polymerase in the lac operon no λ prophage
- Tight control of expression by laclq allows potentially toxic genes to be cloned
- · Control of T7 RNA Polymerase by T7 lysozyme allows toxic genes to be expressed
- LysY is a variant of T7 lysozyme lacking amidase activity, thus cells are less susceptible to lysis during induction
- · Deficient in proteases Lon and OmpT
- · Does not restrict methylated DNA
- · No Cam requirement

Transformation Efficiency:

High Efficiency: 0.6-1 x 109 cfu/µg pUC19 DNA

Resistance: T1 phage (fhuA2), Cam, Nit

Sensitivity: Amp, Kan, Spec, Str, Tet

Disulfide Bonds

Features of Shuffle® Strains:

- Engineered E. coli K12 or B strains promote disulfide bond formation in the cytoplasm
- Constitutively expresses a chromosomal copy of the disufide bond isomerase DsbC
- DsbC promotes the correction of misoxidized proteins into their correct form
- The cytoplasmic DsbC is a chaperone that can also assist in the folding of proteins that do not require disulfide bonds
- Alternative expression strain for proteins that do not fold in wild-type E. coli, independent of redox state


SHuffle strains from NEB are engineered *E. coli* strains capable of expressing proteins with increasing disulfide bond complexity in the cytoplasm. SHuffle strains express the disulfide bond isomerase DsbC within the cytoplasm. DsbC isomerizes mis-oxidized substrates into their correctly folded state greatly enhancing the fidelity of disulfide bond formation. Cytoplasmic expression also results in significantly higher protein yields of

disulfide bonded proteins when compared to periplasmic expression. SHuffle strains are sensitive to kan, amp, tet and in most cases, cam, which makes them able to express proteins from a wide variety of expression vectors offering greater versatility in experimental design.

References:

References for properties and applications for these products can be found at www.neb.com.





Disulfide bond formation in the cytoplasm of wild type E. coli is not favorable, while SHuffle is capable of correctly folding proteins with multiple disulfide bonds in the cytoplasm.

SHuffle Express Competent E. coli

#C3028J

12 x 0.05 ml 224 €

- Folds disulfide bonded proteins in the cytoplasm
- Protease deficient
- Enhanced BL21 derivative
- Free of animal products

Description: *E. coli* cells with enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm.

Genotype: fhuA2 [lon] ompT ahpC $gal \lambda att::pNEB3-r1-cDsbC$ (Spec®, lac!®) $\Delta trxB$ sulA11 R(mcr-73-::miniTn10--Tet\$) 2 [dcm] R(zgb-210::Tn10--Tet\$) endA1 Δgor $\Delta (mcrC-mrr)114::IS10$

Transformation Efficiency:

1 x 107 cfu/µg pUC19 DNA

Resistance: T1 phage (fhuA2), Nit, Str*, Spec

Sensitivity: Amp, Cam, Kan, Tet

*Resistance to low levels of streptomycin may be observed.

SHuffle T7 Express Competent E. coli

#C3029J

12 x 0.05 ml 224 €

- Folds disulfide bonded proteins in the cytoplasm
- T7 expression
- Protease deficient B strain
- Enhanced BL21 derivative
- Free of animal products

Description: T7 Expression *E. coli* strain with enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm.


Genotype: fhuA2 lacZ::T7 gene1 [Ion] ompT ahpC gal λ att::pNEB3-r1-cDsbC (Spec^R, lact^B) Δ trxB sulA11 R(mcr-73-::miniTn10--Tet^S)2 [dcm] R(zgb-210::Tn10--Tet^S) endA1 Δ gor Δ (mcrC-mrr)114::IS10

Transformation Efficiency: 1 x 10⁷ cfu/µg pUC19 DNA

Resistance: T1 phage (fhuA2), Nit, Str*, Spec

Sensitivity: Amp, Cam, Kan, Tet

*Resistance to low levels of streptomycin may be observed.



COMPETENT CELLS

SHuffle T7 Competent E. coli

#C3026J

12 x 0.05 ml 224 €

- Folds disulfide bonded proteins in the cytoplasm
- T7 expression
- K12 strain
- Free of animal products

Description:T7 Expression *E. coli* strain with enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm.

Genotype: F´ lac, pro, lach' / Δ (ara-leu)7697 araD13 fhuA2 lacZ::T7 gene1 Δ (phoA)Pvull phoR ahpC* galE (or U) galK λ att::pNEB3-r1-cDsbC (SpecR, lach) Δ trxB rpsL150(StrR) Δ gor Δ (malF)3

Transformation Efficiency: 1 x 106 cfu/µg pUC19 DNA

Resistance: T1 phage (fhuA2), Nit, Spec, Str*

Sensitivity: Amp, Cam, Kan, Tet

*Resistance to low levels of streptomycin may be observed.

SHuffle T7 Express lysY Competent E. coli

#C3030J

12 x 0.05 ml 224 €

- Express toxic proteins (lysY)
- Folds disulfide bonded proteins in the cytoplasm
- T7 expression
- Protease deficient B strain
- Free of animal products
- Enhanced BL21 derivative

Description: *E. coli* strain with tight T7 Expression control and enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm.

Genotype: MiniF lysY(Cam^R) / fhuA2 lacZ::T7 gene1 [lon] ompT ahpC gal λatt::pNEB3-r1-cDsbC (Spec^R, lac^R) ΔtrxB sulA11 R(mcr-73-::miniTn10--Tet^S)2 [dcm] R(zgb-210::Tn10--Tet^S) endA1 Δgor Δ(mcrC-mrr)114::IS10

Transformation Efficiency:

1 x 107 cfu/µg pUC19 DNA

Resistance: T1 phage (fhuA2), Cam, Nit, Str*, Spec

Sensitivity: Amp, Kan, Tet

*Resistance to low levels of streptomycin may be observed.

How do I express my protein in SHuffle cells?

Currently there are two SHuffle cell lines available from NEB; SHuffle (NEB #C3026) based on *E. coli* K12, and SHuffle Express (NEB #C3028, #C3029, #C3030) based on *E. coli* B.

We recommend testing both B and K12 expression strains, as we do see variability in expression depending on the protein of interest (Table 1). If T7 expression is not necessary, then we recommend comparing NEB #C3026 and #C3028. If T7 expression is necessary, test NEB #C3026 and #C3029. If T7-driven expression of a protein is toxic, switch to a non-leaky lysY version (NEB #C3030). Once the strain is chosen expression conditions should be optimized. This can include temperature as well as auto expression (1).

View our online tutorial for tips on setting up reactions with SHuffle.

Table 1. Percentage of relative solubility of various proteins using SHuffle (K12 and B strains):

	RELATIVE '		
PROTEIN	K12	В	# CYSTEINES
Gluc	65	100	10
Urokinase	60	100	24
vtPA	5	100	12
BSA	100	0	35
Polymerase	100	0	0
Nuclease	100	10	4

Results are determined based on protein levels detected by SDS-PAGE (not shown)

Reference:

 Ke, N. and Berkmen, M. (2014) Current Protocols Molecular Biology 16.1B.21.





Taking Molecular Tools to the Jungle

Climate change, habitat destruction and pollution are causing accelerated species extinction, particularly in tropical ecosystems. Extinction rates that were one species/million/year have increased to 100–1,000 species/million/year. Unfortunately, most of these species have yet to be identified. With an incomplete understanding of species diversity, it is difficult to know where to direct conservation efforts and resources. This knowledge gap has prompted an urgency to catalog as much biodiversity as possible, as quickly as possible.

DNA barcoding is a standardized method of identification that utilizes a short region of the mitochondrial cytochrome c oxidase I (COI) gene in animals (and various other sequences in plants, fungi and protists) to document species quickly and inexpensively. The mutation rate of mitochondrial DNA over relatively short evolutionary periods reflects the diversity between species. It should be noted that barcoding is distinct from, and does not supersede, specialized taxonomic identification of subtle anatomical differences between species, but the combination of molecular and morphological data improves the characterization and delimitation of species.

Biologists have used barcoding in large projects, such as the Census of Marine Life, a 10-year study that assessed the biodiversity and distribution of the Earth's aquatic ecosystems. This study identified 190,000 species, including 6,000 potentially new species. Scientists discovered new habitats, symbiotic relationships and microbial biospheres. They found species that are in decline and new examples of ecosystem resilience.

The colossal effort to document all of the organisms on our planet before they become extinct is a race against time, and it presents researchers with logistical hurdles related to cost and sampling. First, the countries with the greatest biodiversity, such as tropical regions, are often the countries that do not have abundant research resources dedicated to conducting this work. Second, the laws that govern the international transportation of biological materials from biodiverse regions to resource-rich areas can cause delays that compromise sample integrity or prevent sample transportation altogether.

These limitations have necessitated the development of small, portable sequencing tools and technologies that can be used on-site. Conservation scientists can now transport a small DNA sequencing platform, PCR thermocycler, microcentrifuge, reagents and a laptop in a backpack to remote regions, where they can rapidly and cost-effectively extract DNA, amplify and sequence barcodes. This process can be carried out within 24 hours of sample collection, accelerating data acquisition significantly.

Collecting genetic information at the source allows easy documentation of information regarding species health, geographic distribution, hybrid zones, as well as the identification of new species. Rapid access to this information can help focus conservation efforts, and guide the allocation of appropriate resources when planning for species conservation. Geographical areas of highest biological value can be identified and protected, and laws can be implemented to preserve the most endangered species.

Documenting organisms in their habitats is a less invasive, expedient method that can assist conservationists in protecting our planet's immense, yet diminishing, biodiversity.

Glycobiology & Protein Tools



Trust NEB's expertise in enzymology when you need glycobiology reagents.

Glycobiology

Proteomics, the systematic study of proteins in biological systems, has expanded the knowledge of protein expression, modification, interaction and function. However, in eukaryotic cells, the majority of proteins are post-translationally modified (1). A common post-translational modification, essential for cell viability, is the attachment of glycans. Glycosylation defines the adhesive properties of glycoconjugates, and it is largely through glycan–protein interactions that cell–cell and cell–pathogen contacts occur, a fact that accentuates the importance of glycobiology.

Glycomics, the study of glycan expression in biological systems, relies on effective enzymatic and analytical techniques for correlation of glycan structure with function. Glycobiology is a small but rapidly growing field in biology, with relevance to biomedicine, biotechnology, biofuels and basic research. Glycan molecules modulate many other processes important for cell and tissue differentiation, metabolic and gene regulation, protein activity, protein clearance, transport and more (2-9).

Protein Tools

Not only are proteins a major structural component of living systems, they can also be effector molecules whose states determine downstream activities. Therefore, studying the protein complement within a cell can reveal the mechanisms behind many of the cell's responses to its environment. Given the vast number of applications for protein analysis, several tools and methods for its study exist; determining the correct method for your application is paramount to success.

Phage display technology is an *in vitro* screening technique for identifying ligands for proteins and other macromolecules. At the crux of phage display technology is the ability to express peptide or protein sequences as fusions to the coat proteins of a bacteriophage. Libraries of phage-displayed peptides or proteins are thereby physically linked to their encoding nucleic acid, allowing selection of binding partners for myriad target types by iterative rounds of *in vitro* panning and amplification, followed by DNA sequencing.

All NEB products pass stringent quality control assays to ensure the highest level of functionality and purity.

Featured Products

241 Rapid™ PNGase F

243 N-Glycan Sequencing Kit

254 Thermolabile Proteinase K

Featured Tools & Resources

Protein Tools & Glycomics Overview




Glycobiology Unit Conversion Chart



Visit www.NEBglycosidase.com to view our online tutorial on N- and O-linked glycosylation.

(1) Khoury G.A. et al. (2011) Scientific Reports 1: 90. PMID: 22034591.
(2) Varki A. (1993) Glycobiology, 3, 97–130. PMID: 8490246.
(3) Zhao Y.Y. et al. (2008) Cancer Sci. 99, 1304–1310. PMID: 18492092.
(4) Zhao Y. et al. (2008) FEBS J. 275. 1393–1948. PMID: 13834838.
(5) Skropeta D. (2009) Bioorg. Med. Chem. 17, 2645–2653. PMID: 19285412.
(6) Neu U. et al. (2011) J. Cilin. Immunol. 31, 10–21. PMID: 2118745.
(7) Cerliani J.P. et al. (2011) J. Cilin. Immunol. 31, 10–21. PMID: 21184154.
(8) Aarnoudse C.A. et al. (2006) Curr. Opin. Immunol. 18, 105–111. PMID: 16303292.
(9) Arnold J.N. (2006) Immunol. Lett. 106, 103–110. PMID: 16814399.



	Endoglycosidases			Heparin Lyases
RR	Endo F2	238	RR	Bacteroides Heparinase I, II, III
RR	Endo F3	238		
RR	Endo D	239		Phage Display Ph.D7 Phage Display Peptide Library Kit
RR	Endo H, Endo H,	239		
	Endo S	239		Ph.D12 Phage Display Peptide Library Kit
	PNGase A	240		Ph.DC7C Phage Display Peptide Library Kit
RR	PNGase F and PNGaseF, recombinant	240		Ph.D12 Library
	RNase B	240		Ph.D. Peptide Display Cloning System
RR	Rapid PNGase F and			Proteases & Proteome Analysis
	Rapid PNGase F (non-reducing format)	241	RX	IdeZ Protease
	Rapid PNGase F Antibody Standard	241		Trypsin-digested BSA MS Standard
	Remove-iT PNGase F	241		(CAM Modified)
RR	O-Glycosidase	242		Trypsin-ultra, Mass Spectrometry Grade
RR	Protein Deglycosylation Mix II	242	RX	Endoproteinase GluC
	Fetuin	242	RX	Endoproteinase AspN
	Exoglycosidases			Endoproteinase LysC
Dil	N-Glycan Sequencing Kit	243		Proteinase K, Molecular Biology Grade
	α -N-Acetylgalactosaminidase	243	Rit	Thermolabile Proteinase K
	β- <i>N</i> -Acetylglucosaminidase S	243		Factor Xa Protease
		243		Enterokinase, light chain
	β- <i>N</i> -Acetylhexosaminidase, α1-2 Fucosidase	244	Rit	Furin
	α1-3.4 Fucosidase	244	RX	TEV Protease
		244		Protein Phosphatases
	α1-2,3,4,6 Fucosidase	244 245		Properties of Protein Phosphatases from NE
	α 1-2,4,6 Fucosidase O α 1-3.6 Galactosidase	245 245	RR	Lambda Protein Phosphatase (Lambda PP)
				<i>p</i> -Nitrophenylphosphate (PNPP)
	α1-3,4,6 Galactosidase	245		Sodium Orthovanadate
	β1-3 Galactosidase	245		
	β1-3,4 Galactosidase	246		Protein Kinases
	β1-4 Galactosidase S	246		cAMP-dependent Protein Kinase (PKA)
	α1-2,3 Mannosidase	246		Casein Kinase II (CK2)
	α1-2,3,6 Mannosidase	246		Adenosine 5´ Triphosphate
	α1-6 Mannosidase	247		
	α2-3,6,8 Neuraminidase	247		Recombinar
RR	α2-3,6,8,9 Neuraminidase A	247		

Endoglycoceramidase I (EGCase I)

Glycosidases

- Enabling Novel Technologies
- Unique Specifications
- Exceptional Value
- High Purity

NEB offers a selection of endoglycosidases and exoglycosidases for glycobiology research. Many of these reagents are recombinant, and all undergo several quality control assays, enabling us to provide products with lower unit cost, high purity, and reduced lot-to-lot variation.

All of our glycosidases are tested for contaminants. Since p-nitrophenyl-glycosides are not hydrolyzed by some exoglycosidases, we use only fluorescentlylabeled oligosaccharides to screen for contaminating glycosidases.

NEB's glycosidases are provided with 10X buffer to ensure optimal activity. Using more than one glycosidase simultaneously is a common timesaving procedure. Selecting the best buffer to provide reaction conditions that optimize enzyme activity is an important consideration.

Reaction Buffer Compositions:

Visit www.neb.com for details.

Endo F2

#P0772S

480 units 208 €



- Removal of complex biantennary N-linked glycans from glycoproteins and glycopeptides
- Useful for determining N-alycosylation sites

Description: Endo F2 is a highly specific recombinant endoglycosidase which cleaves within the chitobiose core of asparagine-linked complex biantennary and high mannose oligosaccharides from glycoproteins and glycopeptides. Endo F2 cleaves biantennary glycans at a rate approximately 20 times greater than high mannose glycans. The activity of Endo F2 is identical on biantennary structures with and without core fucosylation. However, Endo F2 is not active on hybrid or tri- and tetra-antennary oligosaccharides. Endo F2 is tagged with a chitin binding domain (CBD) for easy removal from a reaction and is supplied glycerol free for optimal performance in HPLC and MS intensive methods.

RR 37° 65

Source: Cloned from Elizabethkingia miricola (formerly Flavobacterium meningosepticum) and expressed in E. coli.

Reagents Supplied:

10X GlycoBuffer 4

Molecular Weight: 39,800 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the carbohydrate from 10 µg Porcine Fibrinogen in 1 hour at 37°C in a total reaction volume of 10 µl.


Concentration: 8,000 units/ml

Heat Inactivation: 65°C for 10 minutes

Endo F3

#P0771S

240 units 187 €



- Removal of complex biantennary and triantennary N-linked glycans from glycoproteins and glycopeptides
- Useful for determining N-glycosylation sites

Description: Endo F3 is a highly specific recombinant endoglycosidase which cleaves within the chitobiose core of asparagine-linked fucosylated-biantennary and triantennary complex oligosaccharides from glycoproteins. Endo F3 is tagged with a chitin binding domain (CBD) for easy removal from a reaction and is supplied glycerol free for optimal performance in HPLC and MS intensive methods.

Source: Cloned from *Elizabethkingia miricola* (formerly Flavobacterium meningosepticum) and expressed in E. coli.

RR 37° 165

Reagents Supplied:

10X GlycoBuffer 4

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the carbohydrate

from 10 µg Porcine Fibrinogen in 1 hour at 37°C in a

total reaction volume of 10 µl.

Molecular Weight: 38,800 daltons

Concentration: 8,000 units/ml

Heat Inactivation: 65°C for 10 minutes












Endo D

#P0742S 1,500 units 187 € 7,500 units 748 € #P0742L



X= (H or oligosaccharide)

- Removal of paucimannose N-linked glycans from glycoproteins and glycopeptides
- Useful for determining N-glycosylation sites

RR 37° 165

Description: Endo D, also known as Endoglycosidase D, is a recombinant glycosidase, which cleaves within the chitobiose core of paucimannose N-linked glycans, with or without extensions in the antennae.

Endo D is tagged with a chitin binding domain (CBD) for easy removal from a reaction, and is supplied glycerol-free for optimal performance in HPLC and MS intensive methods.

Source: A truncated Endo D gene cloned from Streptococcus pneumoniae and expressed in E. coli as a fusion to chitin binding domain.

Reagents Supplied:

10X DTT

10X GlycoBuffer 2

Molecular Weight: 140,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to remove > 95% of the carbohydrate from 10 µg of glycosidase-trimmed (trimannosyl core) Fetuin in 1 hour at 37°C in a total reaction volume of 10 µl.


Concentration: 50,000 units/ml

Heat Inactivation: 65°C for 10 minutes

Endoglycosidase H

Endo H #P0702S 10,000 units72 € #P0702L 50,000 units 290 €

Endo H, #P0703S 100,000 units72 € #P0703L 500,000 units 290 €



Endo H and Endo H, cleave only high mannose structures (n = 2-150, x = $(Man)_{1-2}$, y = H) and hybrid structures (n = 2, x and/or y = AcNeu-Gal-GlcNAc).

 Removal of high mannose N-glycans from glycoproteins

Description: Endoglycosidase H is a recombinant glycosidase which cleaves within the chitobiose core of high mannose and some hybrid oligosaccharides from N-linked glycoproteins.

Endo H, is a recombinant protein fusion of Endoglycosidase H and maltose binding protein. It has identical activity to Endo H.

Source: Endo H and Endo H, have been cloned from Streptomyces plicatus and overexpressed in E. coli.

Reaction Conditions:

Denature glycoprotein in 1X Glycoprotein Denaturing Buffer at 100°C for 10 minutes. Incubate in 1X GlycoBuffer 3 at 37°C. Heat inactivation: 65°C for 10 minutes.

Reagents Supplied:

10X Glycoprotein Denaturing Buffer 10X GlycoBuffer 3

R\\\ 37° \\

Molecular Weight:

Endo H: 29.000 daltons Endo H.: 70,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to remove > 95% of the carbohydrate from 10 µg of denatured RNase B in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration:

Endo H concentration: 500,000 units/ml Endo H, concentration: 1,000,000 units/ml

Usage Notes: Enzymatic activity is not affected by SDS.

To dealy cosylate a native aly coprotein, longer incubation time as well as more enzyme may be required.

Endo S

#P0741S 6.000 units 187 € #P0741L 30,000 units 748 €



- Removal of N-glycans from native IgG
- Useful for determining N-glycosylation sites

Description: Endo S is an endoglycosidase with a uniquely high specificity for removing N-linked glycans from the chitobiose core of the heavy chain of native IgG Endo S is tagged with a chitin binding domain (CBD) for easy removal from a reaction and is supplied glycerol-free for optimal performance in HPLC- and MSintensive methods

Source: Endo S is cloned from *Streptococcus pyogenes* and overexpressed as a fusion to the chitin binding domain in E. coli.

RR 37° 155

Reagents Supplied: 10X GlycoBuffer 1

Molecular Weight: 136,000 daltons

Unit Definition: 5 µg of IgG in 1X GlycoBuffer 1 are incubated with two-fold dilutions of Endo S for 1 hour at 37°C. Separation of reaction products is visualized by

SDS-PAGE.


Concentration: 200.000 units/ml Heat Inactivation: 55°C for 10 minutes





PNGase A

#P0707S 150 units 262 € #P0707L 750 units 1048 €



PNGase A hydrolyzes N-glycan chains from glycoproteins/peptides regardless of the presence of xylose or fucose. [x = H or Man or GlcNAc].

 Removal of N-linked glycans from glycoproteins 37° 1664

Description: PNGase A is a recombinant amidase, which cleaves between the innermost GlcNAc and asparagine residues of high mannose, hybrid, and short complex oligosaccharides such as those found in plant and insect cells from N-linked glycoproteins and glycopeptides. PNGase A differs from PNGase F in that it cleaves N-linked glycans with or without $\alpha(1,3)$ -linked core fucose residues.

Source: Cloned from *Oryza sativa* (rice) and expressed in *Pichia pastoris*.

Reaction Conditions:

Denature 1 μ g of recombinant Avidin in 1X Glycoprotein Denaturing Buffer at 100°C for 10 minutes. After the addition of NP-40 and GlycoBuffer 3, two-fold dilutions of PNGase A are added and the reaction mix is incubated for 1 hour at 37°C. Heat inactivation: 65°C for 10 minutes.

Reagents Supplied:

10X Glycoprotein Denaturing Buffer 10X GlycoBuffer 3 10% NP-40

Molecular Weight: 63,800 daltons

Unit Definition: One unit is defined as the amount of enzyme required to remove > 95% of the carbohydrate from 1 μg of denatured recombinant Avidin produced in Maize in 1 hour at 37°C in a total reaction volume of 10 μl

Concentration: 5.000 units/ml

Usage Notes: PNGase A is active on both glycoproteins and glycopeptides.

PNGase A cannot cleave larger *N*-glycans such as those from Fetuin, Fibrinogen, IgG, Lactoferrin and Transferrin.

PNGase A is able to cleave high mannose *N*-glycan structures from Man 3 up to Man 9.

PNGase F & PNGase F, Recombinant

PNGase F

#P0704S 15,000 units 160 € #P0704L 75,000 units 640 €

PNGase F (Glycerol-free)

#P0705S 15,000 units 160 € #P0705L 75,000 units 640 €

PNGase F, Recombinant

#P0708S 15,000 units 160 € #P0708L 75,000 units 640 €

PNGase F (Glycerol-free), Recombinant #P0709S 15,000 units 160 €

#P0709L 75,000 units 640 €



PNGase F hydrolyzes nearly all types of N-glycan chains from glycopeptides/proteins [x = H or oligosaccharide].

 Removal of N-linked glycans from glycoproteins **Description:** Peptide-N-Glycosidase F, also known as PNGase F, is an amidase which cleaves between the innermost GlcNAc and asparagine residues of high mannose, hybrid and complex oligosaccharides from *N*-linked glycoproteins. A glycerol-free version of PNGase F is also offered for HPLC methods.

Source: NEB #P0704 and NEB #P0705 are purified from *Flavobacterium meningosepticum*.

NEB #P0708 and NEB #P0709 are cloned from *Elizabeth-kingia miricola* (formerly *Flavobacterium meningosepti-cum*) and expressed in *E. coli*.

Reaction Conditions:

Denature glycoprotein in 1X Glycoprotein Denaturing Buffer at 100°C for 10 minutes. Heat inactivation: 75°C for 10 minutes.

RX 37° 📆

Reagents Supplied:

10X Glycoprotein Denaturing Buffer 10X GlycoBuffer 2 10% NP-40

Molecular Weight: 36,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to remove > 95% of the carbohydrate from 10 μ g of denatured RNase B in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 500,000 units/ml

Usage Notes: Since PNGase F activity is inhibited by SDS, it is essential to have NP-40 present in the reaction mixture.

To deglycosylate a native glycoprotein, longer incubation time as well as more enzyme may be required.

After addition of NP-40 and Glycobuffer 2, two-fold dilutions of PNGase F are added and the reaction mix is incubated for 1 hour at 39°C.

Companion Products:

RNase B (control substrate) #P7817S 250 µg66 €

Endoglycosidase Reaction Buffer Pack 10X GlycoBuffer 2, 10X GlycoBuffer 3, 10X Glycoprotein Denaturing Buffer, 10% NP-40 (1 ml of each) #B0701S23 €











Rapid™ PNGase F & Rapid PNGase F (non-reducing format)

Rapid PNGase F

#P0710S 50 reactions 439 €

Rapid PNGase F (non-reducing format) #P0711S 50 reactions 439 €

Companion Product:

Rapid PNGase F Antibody Standard #P6043S 250 ng320 €



- Complete deglycosylation of antibodies and fusion proteins in minutes
- Release of all N-glycans rapidly and without hias
- Optimal activity is ensured for 12 months, if stored properly
- Purified to > 95% homogeneity, as determined by SDS-PAGE

Description: Rapid PNGase F is an improved reagent that allows the complete and rapid deglycosylation of antibodies and fusion proteins in minutes. All N-glycans are released rapidly and without bias, and are ready to be prepared for downstream chromatography or mass spectrometry analysis. Rapid PNGase F creates an optimized workflow which reduces processing time without compromising sensitivity or reproducibility.

Developed for proteomic applications, Rapid PNGase F (non-reducing format) enables complete and rapid dealycosylation while preserving disulfide bonds. This facilitates high throughput proteomics applications and methods for antibody characterization by mass spectrometry such as intact mass analysis. Rapid PNGase F (non-reducing format) combines the advantages of Rapid PNGase F (fast processing time), with non-reducing conditions, preserving quaternary structure.

RN 50° 16

Heat inactivation: 75°C for 10 minutes.

Reagents Supplied:


Rapid PNGase F


Rapid PNGase F Reaction Buffer (5X)

Rapid PNGase F (non-reducing format)


Rapid PNGase F (non-reducing format) Buffer (5X)

Specificity: Rapid PNGase F cleaves all complex. hybrid and high-mannose type glycans from antibodies and related proteins. Core α 1-3 fucosylation (found in immunoglobulins expressed in plant or insect cells) is resistant to both PNGase F and Rapid PNGase F.





ESI-TOF analysis of an antibody before and after treatment with Rapid PNGase F.



ESI-TOF analysis of an antibody before and after treatment with Rapid PNGase F (non-reducing format).

Remove-iT® PNGase F

#P0706S 6,750 units 190 € #P0706L 33,750 units 760 €

Companion Products:

Chitin Magnetic Beads

#E8036S 5 ml 122 € #F8036I 25 ml 488 €

6-Tube Magnetic Separation Rack

6 tubes 206 € #S1506S

12-Tube Magnetic Separation Rack

#S1509S 12 tubes 326 €



Remove-iT PNGase F hydrolyzes nearly all types of N-glycan chains from glycopeptides/proteins [x = H or]oligosaccharide].

 Removal of N-linked glycans from glycoproteins

Description: Remove-iT PNGase F is an amidase

which cleaves between the innermost GIcNAc and asparagine residues of high mannose, hybrid, and complex oligosaccharides from N-linked glycoproteins. Remove-iT PNGase F is tagged with a chitin binding domain (CBD) for easy removal from a reaction and is supplied glycerol free for optimal performance in HPLC and MS intensive methods.

Source: Purified from *Flavobacterium meningosepticum*.

Reaction Conditions:

Denature glycoprotein in 1X DTT (40 mM) at 55°C for 10 minutes. Incubate in 1X GlycoBuffer 2 for 1 hour at 37°C. Heat inactivation: 75°C for 10 minutes.

Reagents Supplied:

10X DTT

10X GlycoBuffer 2

Molecular Weight: 41,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to remove > 95% of the carbohydrate from 5 µg of DTT denatured RNase B in 1 hour at 37°C in a total reaction volume of 10 µl

Concentration: 225.000 units/ml

37° ₩

Usage Notes: Using typical RNase B denaturing conditions with NEB Glycoprotein Denaturing Buffer containing SDS and DTT, Remove-iT PNGase F yields a higher concentration of 500,000 U/ml.

If using Remove-iT PNGase F under typical PNGase F denaturing conditions, it is essential to have NP-40 in the reaction mixture, as Remove-iT PNGase F is inhibited by SDS. It is not known why this non-ionic detergent counteracts the SDS inhibition.

Removal of Remove-iT PNGase F from the deglycosylation reaction can be scaled up linearly with larger volumes of chitin magnetic beads.

O-Glycosidase


#P0733S 2,000,000 units 135 € #P0733L 10,000,000 units 540 €


Companion Products:

O-Glycosidase & Neuraminidase Bundle 1 bundle187 € #E0540S

α2-3,6,8 Neuraminidase

2,000 units 68 € #P0720S #P0720L 10,000 units 272 €





Removal of Core 1 and Core 3 O-linked disaccharide glycans from glycoproteins **Description:** O-Glycosidase, also known as Endo- α -N-Acetylgalactosaminidase, catalyzes the removal of Core 1 and Core 3 O-linked disaccharides from glycoproteins.

Source: Cloned from Enterococcus faecalis and expressed in E. coli.

Reagents Supplied:

10X Glycoprotein Denaturing Buffer 10X GlycoBuffer 2 10% NP-40

Molecular Weight: 147,000 daltons

RX 37° 65

Unit Definition: One unit is defined as the amount of enzyme required to remove 0.68 nmol of O-linked disaccharide from 5 mg of neuraminidase-digested, non-denatured fetuin in 1 hour at 37°C in a total reaction volume of 100 μ l (1 unit of both $\emph{O}\text{-Glycosidase}$ and PNGase F will remove equivalent molar amounts of O-linked disaccharides and N-linked oligosaccharides, respectively).

Concentration: 40,000,000 units/ml Heat Inactivation: 65°C for 10 minutes Usage Note: Remove sialic acids if present

Protein Deglycosylation Mix II

20 reactions 508 € #P6044S

- Fast reaction setup
- Enzyme mixture ensures effective deglycosylation of N-and O-linked glycans
- Can be used under native and reducing conditions
- Enzymatic deglycosylation leaves intact core structures suitable for mass spectrometry analysis

Description: The Protein Deglycosylation Mix II contains all of the enzymes, reagents and controls necessary to remove almost all N-linked and simple O-linked glycans, as well as some complex O-linked glycans. This mix contains enzymes sufficient for 20 reactions or the cleavage of as much as 2 mg of glycoprotein.


Reagents Supplied:

Deglycosylation Mix Buffer 1 (10X) Deglycosylation Mix Buffer 2 (10X)

Deglycosylation Enzyme Mix II:

PNGase F (Glycerol-free), Recombinant: 10,000 units/vial O-Glycosidase: 80,000 units/vial α2-3,6,8,9 Neuraminidase A: 400 units/vial β1-4 Galactosidase S: 960 units/vial β-N-Acetylhexosaminidase,: 300 units/vial

R\\\ 37°



Enzymatic Deglycosylation of Bovine Fetuin under both native (10X Deglycosylation Mix Buffer 1) and reducing (10X Deglycosylation Mix Buffer 2) conditions. 20 µg reactions were loaded onto a 10-20% Tris-glycine SDS-PAGE gel. Lane 1: Color Prestained Protein Standard, Broad Range (11-245 kDa), Lane 2: 20 µg untreated Fetuin control, Lane 3: 20 µg Fetuin deglycosylated under native conditions with Deglycosylation Mix Buffer 1, Lane 4: 20 µg Fetuin deglycosylated under reducing conditions with Deglycosylation Mix Buffer 2, Lane 5: 5 µl Protein Deglycosylation Mix II

Fetuin

#P6042S 500 μg73 € **Description:** Fetuin is a glycoprotein containing sialylated N-linked and O-linked glycans that can be used as a positive control for endoglycosidase enzymes.

Source: Fetal Calf Serum Molecular Weight: 64 kDa Concentration: 10 mg/ml

Note: 500 µg is enough for approximately 20 reactions. Due to heterogeneous glycosylation, Fetuin runs as a doublet on an SDS-PAGE gel.











N-Glycan Sequencing Kit

#E0577S

20 reactions 710 €

- Recombinant enzymes with no detectable endoglycosidase or other exoglycosidase contaminating activities
- Compatible with N-linked glycans released from a variety of CHO and murine derived antibodies, as well as N-linked glycans released from other glycoproteins
- Simultaneous digestion with other exoglycosidases
- Optimal activity and stability for up to 12 months

Description: The *N*-Glycan Sequencing Kit consists of seven well characterized, highly pure, recombinant exoglycosidase enzymes selected to simplify the process of characterizing typical N-linked glycan structures.

R**%** 37°

The N-Glycan Sequencing Kit Includes:

- α2-3,6,8,9 Neuraminidase A
- α2-3 Neuraminidase S
- β-N-Acetylglucosaminidase S
- β1-4 Galactosidase S
- α1-3,4,6 Galactosidase
- $-\alpha$ 1-2.4.6 Fucosidase O
- α1-2,3,6 Mannosidase
- Zinc
- GlycoBuffer 1



α-N-Acetylgalactosaminidase

#P0734S #P0734L

3,000 units 135 € 15,000 units 540 €



Description: α -N-Acetylgalactosaminidase is a highly specific exoglycosidase that catalyzes the hydrolysis of α -linked D-*N*-Acetylgalactosamine residues from oligosaccharides and N-glycans attached to proteins.

Source: Cloned from *Chryseobacterium* meningosepticum and expressed in E. coli.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 µg/ml BSA. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

RN BSA 37° KK

Reagents Supplied:

10X GlycoBuffer 1 100X BSA

Molecular Weight: 47,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal α -D-*N*-Acetylgalactosamine from 1 nmol of (GalNAc α 1-3) (Fucα1-2)Galβ1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 20,000 units/ml

β-N-Acetylglucosaminidase S

#P0744S #P0744L

100 units 129 €

500 units 516 €



Removal of bisecting β-GlcNAc residues

Description: β-*N*-Acetylglucosaminidase S is a highly specific exoglycosidase that catalyzes the hydrolysis of terminal, non-reducing β -N-Acetylglucosamine residues from oligosaccharides.

Source: Cloned from *Streptococcus pneumoniae* and expressed in E. coli.

Reaction Conditions: 1X GlycoBuffer 1. Incubate at 37°C.

Reagents Supplied:

10X GlycoBuffer 1

RR 37° Wb

Molecular Weight: 125,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal, non-reducing β-N-Acetylglucosamine from 1 nmol GlcNAcβ1-4GlcNAcβ1-4GlcNAc-AMC, in 1 hour at 37°C in a total reaction volume of 10 μl.

Concentration: 4,000 units/ml

β -N-Acetylhexosaminidase,

#P0721S 500 units66 € #P0721L 2,500 units264 €





Active only on linear substrates

Description: β-*N*-Acetylhexosaminidase, is a recombinant protein fusion of β -*N*-Acetylhexosaminidase and maltose binding protein with identical activity to β -*N*-Acetylhexosaminidase. It catalyzes the hydrolysis of terminal β -p-*N*-Acetylgalacto-samine and glucosamine residues from oligosaccharides.

Source: Cloned from *Streptomyces plicatus* and overexpressed in *E. coli.*

Reaction Conditions: 1X GlycoBuffer 1. Incubate at 37°C. Heat inactivation: 75°C for 10 minutes.

R**₹** 37° ₩

Reagents Supplied: 10X GlycoBuffer 1

Molecular Weight: 100,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal β-D-N-Acetylgalactosamine from 1 nmol of GalNAcβ1-4Galβ1-4Glc-AMC, in 1 hour at 37°C in 10 μ1 volume.

Concentration: 5.000 units/ml

α1-2 Fucosidase

#P0724S 1,000 units 130 € #P0724L 5,000 units 520 €



Active only on linear substrates

Description: α1-2 Fucosidase is a highly specific exoglycosidase that catalyzes the hydrolysis of linear α1-2 linked ι -fucopyranosyl residues from oligosaccharides. In this case, a linear substrate is defined as having no branching on the adjacent residue.

Source: Cloned from *Xanthomonas manihotis* and expressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 μg/ml BSA. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

RN BSA 37° V654

Reagents Supplied:

10X GlycoBuffer 1 100X BSA

Molecular Weight: 70,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the α -L-fucose from 1 nmol of Fuc α 1-2Gal β 1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 20,000 units/ml

α1-3,4 Fucosidase

#P0769S 200 units 129 € #P0769L 1,000 units 516 €



Description: α 1-3,4 Fucosidase, (also known as AMF) is a broad specificity exoglycosidase that catalyzes the hydrolysis of α 1-3 and α 1-4 linked fucose residues from oligosaccharides and glycoproteins.

Source: Cloned from the sweet almond tree (*Prunus dulcis*) and expressed in *Pichia pastoris*.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 μg/ml BSA. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

RN BSA 37° KK

Reagents Supplied: 10X GlycoBuffer 1 100X BSA

Molecular Weight: 56,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the α -fucose from 1 nmol of Gal β 1-4GlcNAc β 1-3(Fuc α 1-3)Gal β 1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 µl

Concentration: 4,000 units/ml

α 1-2,3,4,6 Fucosidase

#P0748S 400 units 133 € #P0748L 2,000 units 532 €



Description: α 1-2,3,4,6 Fucosidase is a broad specificity exoglycosidase that catalyzes the hydrolysis of α 1-2, α 1-3, α 1-4 and α 1-6 linked ι -fucopyranosyl residues from oligosaccharides.

Source: Cloned from bovine kidney and expressed in *E. coli*.

Reaction Conditions: 1X Glycobuffer 1. Supplement with 100 μg/ml BSA. Incubate at 37°C. Heat inactivation: 100°C for 10 minutes.

RN BSA 37° WW

Reagents Supplied: 10X Glycobuffer 1

100X BSA

Molecular Weight: 51,800 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the α -1-fucose from 1 nmol of Fuc α 1-2Gal β 1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 8,000 units/ml











α1-2,4,6 Fucosidase O

#P0749S #P0749L

400 units 131 € 2,000 units 524 €



Description: α 1-2,4,6 Fucosidase 0 is a broad specificity exoglycosidase that catalyzes the hydrolysis of terminal α 1-2, α 1-4 and α 1-6 linked fucose residues from oligosaccharides. α 1-2,4,6 Fucosidase O cleaves α 1-6 fucose residues more efficiently than other linkages.

Source: Cloned from *Omnitrophica* bacterium and expressed in E. coli.

Reaction Conditions: 1X Glycobuffer 1. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

Reagents Supplied: 10X Glycobuffer 1

RX 37° 1664

Molecular Weight: 49,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the fucose from 1 nmol of GOF from human IgG [GIcNAc_B1-2Manα1-6(GlcNAcβ1-2Manα1-3)Manβ1-4GlcNAcβ1-4GlcNAc(Fucα1-6)-AMAC], in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 2,000 units/ml

α1-3,6 Galactosidase

#P0731S #P0731L

100 units 128 € 500 units 512 €



Description: α 1-3, 6 Galactosidase is a highly specific exoglycosidase that catalyzes the hydrolysis of α 1-3, 6 linked p-galactopyranosyl residues from oligosaccharides.

Source: Cloned from Xanthomonas manihotis and expressed in E. coli.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 µg/ml BSA. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

Optimal incubation times and enzyme concentrations must be determined empirically for a particular substrate.

RN BSA 37° 1654

Reagents Supplied: 10X GlycoBuffer 1 100X BSA

Molecular Weight: 70,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal, α -p-galactose from 1 nmol Gal α 1-3Gal β 1-4Gal-AMC, in 1 hour at 37°C in a total reaction volume of 10 μl.

Concentration: 4,000 units/ml

α1-3,4,6 Galactosidase

#P0747S #P0747L

200 units 131 € 1,000 units 524 €



Description: α 1-3,4,6 Galactosidase is a highly specific exoglycosidase that catalyzes the hydrolysis of α 1-3, α 1-4 and α 1-6 linked D-galactopyranosyl residues from oligosaccharides.

Source: Cloned from green coffee bean and expressed in E. coli.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 μ g/ml BSA. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

Optimal incubation times and enzyme concentrations must be determined empirically for a particular substrate.

RN BSA 37° VIII

Reagents Supplied: 10X GlycoBuffer 1 100X BSA

Molecular Weight: 39,700 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal, α -D-galactose from 1 nmol Gal α 1-3Gal β 1-4Gal-AMC, in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 8,000 units/ml

β1-3 Galactosidase

#P0726S #P0726L

500 units 130 € 2,500 units 520 €



Description: β1-3 Galactosidase is a highly specific exoglycosidase that catalyzes the hydrolysis of β 1-3 and, at a much lower rate. 61-6 linked p-galactopyranosyl residues from oligosaccharides. The approximate kinetic data show > 100-fold preference for β 1-3 over β 1-6 linkages and > 500-fold preference from β 1-3 over β 1-4 linkages.

Source: Cloned from Xanthomonas manihotis and expressed in E. coli.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 μg/ml BSA. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

RN BSA 37° KK

Reagents Supplied:

10X GlycoBuffer 1 100X BSA

Molecular Weight: 66,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal β -Dgalactose from 1 nmol of Gal\u00e41-3GlcNAc\u00bb1-3Gal\u00bb1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 10,000 units/ml

β1-3,4 Galactosidase

#P0746S 400 units 129 € #P0746L 2,000 units 516 €



Description: β1-3,4 Galactosidase, cloned from bovine testis and also known as BTG, is a highly specific exoglycosidase that catalyzes the hydrolysis of terminal β1-3 and β1-4 linked galactose residues from oligosaccharides.

Source: Cloned from bovine testis and expressed in *Pichia pastoris*.

Reaction Conditions: 1X GlycoBuffer 4. Incubate at 37°C.Heat inactivation: 65°C for 10 minutes.

RR 37° 65

Reagents Supplied: 10X GlycoBuffer 4

Molecular Weight: 71.000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal, β -p-galactose from 1 nmol Gal β 1-4GlcNAc β 1-3Gal β 1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 III

Concentration: 8,000 units/ml

β1-4 Galactosidase S

#P0745S 400 units 129 € #P0745L 2,000 units 516 €



Description: β1-4 Galactosidase S is a highly specific exoglycosidase that catalyzes the hydrolysis of β1-4 linked galactose residues from oligosaccharides.

Source: Cloned from *Streptococcus pneumoniae* and expressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Incubate at 37°C. Heat inactivation: 75°C for 10 minutes.

RR 37° K

Reagents Supplied: 10X GlycoBuffer 1

Molecular Weight: 231,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal, β -D-galactose from 1 nmol Gal β 1-4GlcNAc β 1-3Gal β 1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 8,000 units/ml

α1-2,3 Mannosidase

#P0729S 640 units 130 € #P0729L 3,200 units 520 €



Description: α 1-2,3 Mannosidase is a highly specific exoglycosidase that catalyzes the hydrolysis of α 1-2 and α 1-3 linked p-mannopyranosyl residues from oligosaccharides.

Source: Cloned from *Xanthomonas manihotis* and expressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 μg/ml BSA. Incubate at 37°C. Heat inactivation: 55°C for 10 minutes.

RN BSA 37° 😽

Reagents Supplied: 10X GlycoBuffer 1 100X BSA

Molecular Weight: 90,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal α -D-Mannose from 1 nmol of Man α 1-3Man β 1-4GlcNAc-AMC, in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 32.000 units/ml

α1-2,3,6 Mannosidase

#P0768S 80 units 130 € #P0768L 400 units 520 €



Description: α 1-2,3,6 Mannosidase, cloned from Jack Bean, and also known as JBM, is a broad specificity exoglycosidase that catalyzes the hydrolysis of terminal α 1-2, α 1-3 and α 1-6 linked mannose residues from oligosaccharides. α 1-2,3,6 Mannosidase has a slight preference for α 1-2 mannose residues over α 1-3 and α 1-6 mannose residues.

Source: Cloned from *Canavalia ensiformis* (Jack Bean) and expressed in *Pichia pastoris*.

Reaction Conditions: 1X GlycoBuffer 4 and 1X Zinc. Incubate at 37°C. Heat inactivation: 95°C for 10 minutes.

RN 37° 1

Reagents Supplied:

10X GlycoBuffer 4 10X Zinc

Molecular Weight: 110,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal mannose from 1 nmol of Man(α 1,3)-Man(β 1,4)-GlcNAc-AMC, in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 2,000 units/ml

α1-6 Mannosidase

#P0727S #P0727L

800 units 130 € 4,000 units 520 €



Description: α 1-6 Mannosidase is a highly specific exoglycosidase that removes unbranched α 1-6 linked p-mannopyranosyl residues from oligosaccharides. When used in conjunction with α 1-2,3 Mannosidase, the α 1-6 Mannosidase will cleave α 1-6 Mannose residues from branched carbohydrate substrates.

Source: Cloned from Xanthomonas manihotis and expressed in E. coli.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 µg/ml BSA. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

RN BSA 37° WSF

Reagents Supplied:

10X GlycoBuffer 1 100X BSA

R**₹** 37° **★**

Molecular Weight: 58,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal α -Dmannose from 1 nmol of $Man\alpha 1$ -6 $Man\alpha 1$ -6Man-AMC, in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 40,000 units/ml

Note: p-nitrophenyl- α -D-mannopyranoside is NOT a substrate for this enzyme.

α 2-3,6,8 Neuraminidase

#P0720S #P0720L

2,000 units68 € 10,000 units 272 €



Active from pH 4.5 to 8.5

Description: Neuraminidase is the common name for Acetyl-neuraminyl hydrolase (Sialidase). This Neuraminidase catalyzes the hydrolysis of α 2-3, α 2-6 and α 2-8 linked N-acetylneuraminic acid residues from glycoproteins and oligosaccharides.

Source: Cloned from Clostridium perfringens and overexpressed in E. coli.

Reaction Conditions: 1X GlycoBuffer 1. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

Reagents Supplied: 10X GlycoBuffer 1

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal α-Neu5Ac from 1 nmol of Neu5Acα2-3Galβ1-3GlcNAcβ1-3Galβ1-4Glc-AMC, in 5 minutes at 37°C in a total reaction volume of 10 µl.

Concentration: 50.000 units/ml

Molecular Weight: 43,000 daltons

Note: This enzyme shows a preference for $\alpha 2.3$ and α 2,6 linkages over α 2,8 linkages.

α2-3,6,8,9 Neuraminidase A

#P0722S #P0722L

800 units69 € 4,000 units 276 €



 Removes branched sialic acid residues that are linked to an internal residue

Description: Neuraminidase is the common name for Acetyl-neuraminyl hydrolase (Sialidase). α 2-3,6,8,9 Neuraminidase A catalyzes the hydrolysis of all linear and branched non-reducing terminal sialic acid residues from glycoproteins and oligosaccharides. The enzyme releases α 2-3 and α 2-6 linkages at a slightly higher rate than α 2-8 and α 2-9 linkages.

Source: Cloned from Arthrobacter ureafaciens and expressed in E. coli.

Reaction Conditions: 1X GlycoBuffer 1. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

RX 37° 65

Reagents Supplied:

10X GlycoBuffer 1

Molecular Weight: 100,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal α -Neu5Ac from 1 nmol Neu5Ac α 2-3Gal β 1-3GlcNAc β 1-3Galβ1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 20,000 units/ml

α2-3 Neuraminidase S

#P0743S #P0743L

400 units68 € 2,000 units 272 €



Description: Neuraminidase is the common name for Acetyl-neuraminyl hydrolase (Sialidase). α2-3 Neuraminidase S is a highly specific exoglycosidase that catalyzes the hydrolysis of α 2-3 linked N-acetyl-neuraminic acid residues from glycoproteins and oligosaccharides.

Source: Cloned from Streptococcus pneumoniae and expressed in E. coli.

Reaction Conditions: 1X GlycoBuffer 1. Incubate at 37°C. Heat inactivation: 75°C for 10 minutes.

RX 37° ₩

Reagents Supplied:

10X GlycoBuffer 1

Molecular Weight: 74,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the α -Neu5Ac from 1 nmol of Neu5Acα2-3Galβ1-3GlcNAcβ1-3Galβ1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 8.000 units/ml

IdeZ Protease (lgG-specific)

#P0770S

4,000 units 420 €

human lgG1, lgG3, lgG4: CPAPELLG GPSVF human lgG2: CPAPPVA GPSVF murine IgG2a: CPAPNLLG GPSVF murine IgG3: CPPGNILG GPSVF

 Complete fragmentation of antibodies and immunoglobulin fusion proteins in 30 minutes under native conditions

Description: IdeZ Protease (IgG-specific) is a recombinant antibody specific protease cloned from Streptococcus equi subspecies zooepidemicus that recognizes all human, sheep, monkey, and rabbit IgG subclasses, specifically cleaving at a single recognition site below the hinge region, yielding a homogenous pool of F(ab´)2 and Fc fragments. IdeZ Protease more effectively cleaves murine IgG2a than IdeS.

Source: Cloned from Streptococcus equi subspecies zooepidemicus and expressed in E. coli.

RN 37° 👑

Reaction Conditions:

1X GlycoBuffer 2. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.


Reagents Supplied with Enzyme:

10X GlycoBuffer 2

Molecular Weight: 35,578 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of 1 µg of human IgG, in 15 minutes at 37°C in a total reaction volume of 10 µl.

Concentration: 80,000 units/ml



Digestion of IgG with IdeZ Protease (IgG-specific), followed by denaturation.



Angela, Alycia and Anne have worked at NEB for 9, 9 and 30 years, respectively. This mother/daughter trio has worked in various departments within the company and are well-loved by the NEB community.











Bacteroides Heparinase I

#P0735S 240 units 153 € #P0735L 600 units 306 €

HO NHSO₃H COOH

Denotes either glucuronic acid or iduronic acid.

All structural determinants for enzyme specificity are displayed in red.

 Degradation of heparin and heparan sulfate glycosaminoglycans **Description:** *Bacteroides* Heparinase I cloned from *Bacteroides eggerthii*, also called Heparin Lyase I, is active on heparin and the highly sulfated domains of heparan sulfate. The reaction yields oligosaccharide products containing unsaturated uronic acids which can be detected by UV spectroscopy at 232 nm.

Source: Cloned from *Bacteroides eggerthii* and expressed in *E. coli*.

Reaction Conditions:

10 μg heparin substrate, 10 μl *Bacteroides* Heparinase Reaction Buffer and H_2O in a total reaction volume of 100 μl. Incubate reaction at 30°C. Heat inactivation: 100°C for 1 minute.

RX 30° W/

Reagents Supplied:

10X Bacteroides Heparinase Reaction Buffer


Molecular Weight: 42,000 daltons

Unit Definition: One unit is defined as the amount of enzyme that will liberate 1.0 μ mol unsaturated oligosaccharides from porcine mucosal heparin per minute at 30°C and pH 7.0 in a total reaction volume of 100 μ l. For unit assay conditions, visit www.neb.com.

Concentration: 12,000 units/ml

Bacteroides Heparinase II

#P0736S 80 units 231 € #P0736L 200 units 462 €



Denotes either glucuronic acid or iduronic acid.

All structural determinants for enzyme specificity are displayed in red.

 Degradation of heparin and heparan sulfate glycosaminoglycans **Description:** *Bacteroides* Heparinase II cloned from *Bacteroides eggerthii*, also called Heparin Lyase II, is active on heparin and heparan sulfate. The reaction yields oligosaccharide products containing unsaturated uronic acids which can be detected by UV spectroscopy at 232 nm.

Source: Cloned from *Bacteroides eggerthii* and expressed in *E. coli*.

Reaction Conditions:

 $10 \mu g$ heparin substrate, $10 \mu l$ *Bacteroides* Heparinase Reaction Buffer and H₂O in a total reaction volume of $100 \mu l$. Incubate reaction at $30 \, ^{\circ}$ C. Heat inactivation: $100 \, ^{\circ}$ C for 1 minute.

RR 30° W

Reagents Supplied:

10X Bacteroides Heparinase Reaction Buffer


Molecular Weight: 86,000 daltons

Unit Definition: One unit is defined as the amount of enzyme that will liberate 1.0 μmol unsaturated oligosaccharides from porcine mucosal heparin per minute at 30°C and pH 7.0 in a total reaction volume of 100 μl. For unit assay conditions, visit www.neb.com.

Concentration: 4,000 units/ml

Bacteroides Heparinase III

#P0737S 14 units 310 € #P0737L 35 units 620 €



Denotes either glucuronic acid or iduronic acid.
All structural determinants for enzyme specificity are displayed in red.

 Degradation of heparan sulfate glycosaminoglycans **Description:** Bacteroides Heparinase III cloned from Bacteroides eggerthii, also called Heparin Lyase III, is active on heparan sulfate. The reaction yields oligosaccharide products containing unsaturated uronic acids which can be detected by UV spectroscopy at 232 nm.

Source: Cloned from *Bacteroides eggerthii* and expressed in *E. coli*.

Reaction Conditions:

10 μg heparan sulfate substrate, 10 μl *Bacteroides* Heparinase Reaction Buffer and H_2 0 in a total reaction volume of 100 μl. Incubate reaction at 30°C. Heat inactivation: 100°C for 1 minute.

RX 30° 👑

Reagents Supplied:

10X Bacteroides Heparinase Reaction Buffer

Molecular Weight: 75,000 daltons

Unit Definition: One unit is defined as the amount of enzyme that will liberate 1.0 µmol unsaturated oligosaccharides from heparan sulfate per minute at 30°C and pH 7.0 in a total reaction volume of 100 µl. For unit assay conditions, visit www.neb.com.

Concentration: 700 units/ml

Ph.D.™ Phage Display Peptide Library Kits

The Ph.D. Kits Include:

- Sufficient Phage Display Library for 10 separate panning experiments, complexity of 10⁹ clones
- 28 gIII Sequencing Primer (100 pmol)
- 96 gIII Sequencing Primer (100 pmol)
- Host E. coli strain ER2738
- Control Target (Streptavidin) and Elutant (Biotin)
- Detailed Protocols

Description: Phage display describes a selection technique in which a library of peptide or protein variants is expressed on the outside of a phage virion, while the genetic material encoding each variant resides on the inside. This creates a physical linkage between each variant protein sequence and the DNA encoding it, which allows rapid partitioning based on binding affinity to a given target molecule (antibodies, enzymes, cell-surface receptors, etc.) by an in vitro selection process called panning. In its simplest form (Figure 1), panning is carried out by incubating a library of phage-displayed peptides with a plate (or bead) coated with the target, washing away the unbound phage, and eluting the specifically-bound phage. The eluted phage is then amplified and taken through additional binding/amplification cycles to enrich the pool in favor of binding sequences. After 3-4 rounds, individual clones are characterized by DNA sequencing and ELISA.

New England Biolabs offers 3 pre-made random peptide libraries, as well as the cloning vector M13KE for construction of custom libraries. The pre-made libraries consist of linear heptapeptide (Ph.D.-7) and dodecapeptide (Ph.D.-12) libraries, as well as a disulfide-constrained heptapeptide (Ph.D.-C7C) library. All of the libraries have complexities in excess of 2 billion independent clones. The randomized peptide sequences in all three libraries are expressed at the N-terminus of the minor coat protein plll, resulting in a valency of 5 copies of the displayed peptide per virion.

The Ph.D. libraries have been used for myriad applications, including epitope mapping (Figure 2), identification of protein-protein contacts and enzyme inhibitors and discovery of peptide ligands for GroEL, HIV, semiconductor surfaces and small-molecule fluorophores and drugs.




Figure 1: Routine Phage Display Workflow. Round 1: Incubate 10¹¹ pfu Ph.D. library + target incubation, wash away non-binders, elute bound phage, enrich selected phage with amplification in E. coli. Carry out 3-4 rounds of selection and then proceed with sequencing and/or phage-ELISA.

β-endorphin: Y G G F M T S E K Q T P... 1st round sequences: V G W I S P P I H I P T

G W I Н Q S Υ W R L D S S Α S Ρ F S н н Р 2nd round sequences:

Υ	G	G	F	L	1	G	L	Q	D	Α	S
Υ	G	G	F	Н	Υ	K	Ε	Τ	G	Α	L
Υ	Q	Р	D	N	Р	S	R	Q	1	Α	N
٧	Υ	С	Υ	1	N	Q	S	M	1	G	N
Н	Н	D	T	Ε	Υ	R	T	Τ	Q	L	S
N	L	K	F	Р	Т	N	Р	K	Α	M	W
L	Р	N	L	Τ	W	Α	L	M	Р	R	Α
D	N	W	Р	Υ	R	Р	S	F	S	L	S
S	Н	N	T	Υ	S	Α	Р	R	Р	S	Α
S	L	L	Н	Υ	Α	S	S	L	S	L	M
V	Т	M	N	Т	K	Т	Р	G	Р	M	Р

3rd round sequences:




Figure 2: Epitope mapping of an anti- β -endorphin monoclonal antibody with the Ph.D.-12 library.

The Ph.D.-12 library was panned against anti- β -endorphin antibody 3-E7 in solution (10 nM antibody), followed by affinity capture of the antibody-phage complexes onto Protein A-agarose (rounds 1 and 3) or Protein G-agarose (round 2). Bound phage were eluted with 0.2 M glycine-HCl, pH 2.2. Selected sequences from each round are shown aligned with the first 12 residues of β -endorphin; consensus elements are boxed. The results clearly show that the epitope for this antibody spans the first 7 residues of β -endorphin, and that the bulk of the antibody-antigen binding energy is contributed by the first 4 residues (YGGF), with some flexibility allowed in the third position. Additionally, the conserved position of the selected sequences within the 12 residue window indicates that the free α -amino group of the N-terminal tyrosine is part of the epitope.


Ph.D. Peptide Display Cloning System

Ph.D. Peptide Display Cloning System 20 µg M13KE gIII Cloning Vector 16 µg Extension Primer

#E8101S174 €

Description: The Ph.D. Peptide Display Cloning System facilitates the display of custom peptide libraries on the surface of bacteriophage M13 as coat protein fusions, creating a physical linkage between each displayed peptide and its encoding DNA sequence. Peptide ligands for a variety of targets can then be selected by the straightforward method of panning. The supplied display vector M13KE is an M13 derivative with cloning sites engineered for N-terminal pllI fusion, resulting in a valency of 5 displayed peptides per virion. The use of a phage vector, rather than a phagemid, simplifies the intermediate amplification steps, since neither antibiotic selection nor helper phage superinfection are required. Since displayed proteins longer than 20–30 amino acids have a deleterious effect on the infectivity function of pllI

in phage vectors, **this vector is suitable only for the display of short peptides.** Included with the cloning vector is an insert extension primer as well as a detailed protocol for cloning a peptide library into M13KE.



Trypsin-digested BSA MS Standard (CAM Modified)


#P8108S

500 pmol78 €

Standardization range: 500 to 2400 Da

Description: A complex mixture of peptides produced by Trypsin digestion of Bovine Serum Albumin (BSA) that was reduced and alkylated with Iodoacetimide (CAM modified). This peptide mixture can be used to test a Matrix-Assisted Laser Desorption/Ionization Time-Of-Flight (MALDI-TOF) or Electrospray Ionization (ESI) mass spectrometer (TOF, Q-TOF or Ion Trap).

Source: BSA (GENBANK P02769) was digested using Trypsin-ultra, Mass Spectrometry Grade, (NEB #P8101).



One hundred fmol of resuspended peptide mix may be analyzed by reverse phase liquid chromatography with on-line MS/MS analysis, for example with a Proxeon EASY-nLC and by Orbitrap Mass Spectrometer. Both analytical methods reveal a range of peptides in the standard. At least sixty percent sequence coverage is seen after database search, with greater than 15 unique peptides being identified.

Trypsin-ultra[™], Mass Spectrometry Grade

#P8101S

100 μg81 €

Lys/Arg ▼XXX

- Digestion of proteins for proteomic analysis by mass spectrometry
- Protein and peptide identification
- TPCK treatment eliminates chymotryptic activity
- Free of contaminating proteases

Description: Trypsin-ultra, Mass Spectrometry Grade is a serine endopeptidase. It selectively cleaves peptide bonds C-terminal to lysine and arginine residues. Trypsin-ultra, Mass Spectrometry Grade is treated with L-(tosylamido-2-phenyl) ethyl chloromethyl ketone (TPCK) to inactivate any remaining chymotryptic activity. It is modified by acetylation of the ε -amino groups of lysine residues to prevent autolysis. Trypsin-ultra, Mass Spectrometry Grade cleaves at Lys-Pro and Arg-Pro bonds at a much slower rate than other amino acid residues.

Source: Isolated from bovine (*Bos taurus*) pancreas.

Reaction Conditions:

1X Trypsin Ultra Reaction Buffer. Incubate at 37°C.

Reagents Supplied with Enzyme:

2X Trypsin-ultra Reaction Buffer

Molecular Weight: 23,675 daltons

Reconstitution: Trypsin-ultra, Mass Spectrometry Grade should be reconstituted by the addition of 20-200 µl of high purity water. Rapid autolysis is a function of enzyme concentration.

Notes: Can be stored frozen in solution at -20°C for up to 1 week. A decrease in activity will occur if stored in solution. Use only freshly reconstituted protease for best results

Sue is a member of our NEBNext Production team, and has worked at NEB for over 25 years! Sue has many interests outside of NEB, including sewing, biking, fishing and spending time with her grandchildren.













GLYCOBIOLOGY & PROTEIN TOOLS

Endoproteinase GluC

#P8100S

50 μg77 €

XX-Glu[▼]XX

- Ideal for proteomic analysis by mass spectrometry
- Protein & peptide identification
- Free of contaminating proteases.
 Produced from a protease-deficient Bacillus subtilis strain

Description: Endoproteinase GluC (*Staphylococcus aureus* Protease V8) is a serine proteinase which selectively cleaves peptide bonds C-terminal to glutamic acid residues. Endoproteinase GluC also cleaves at aspartic acid residues at a rate 100–300 times slower than at glutamic acid residues.

Source: *Staphylococcus aureus* Protease V8 gene cloned and expressed in *Bacillus subtilis*.

Reaction Conditions:

1X GluC Reaction Buffer. Incubate at 37°C.



Reagents Supplied with Enzyme:

2X GluC Reaction Buffer

Molecular Weight: 29,849 daltons

Reconstitution: Endoproteinase GluC should be reconstituted by the addition of 50–500 µl of high purity water. Rapid autolysis is a function of enzyme concentration.

Note: Can be stored frozen in solution at –20°C for up to 2 weeks. A decrease in activity will occur if stored in solution. Use only freshly reconstituted protease for best results.

Endoproteinase AspN

#P8104S

50 μg76 €

XX▼Asp-XXX

- Ideal for proteomic analysis by mass spectrometry
- Free of contaminating proteases
- Best suited for peptide identification

Description: Endoproteinase AspN (flavastacin) is a zinc metalloendopeptidase which selectively cleaves peptide bonds N-terminal to aspartic acid residues.

Source: Purified from Flavobacterium menigosepticum.

Reaction Conditions:

1X AspN Reaction Buffer. Incubate at 37°C.

Reagents Supplied with Enzyme:

2X AspN Reaction Buffer

Molecular Weight: 40,089.9 daltons



Reconstitution: Endoproteinase AspN should be reconstituted by the addition of $50-500 \mu I$ of high purity water. Rapid autolysis is a function of enzyme concentration.

Notes: Can be stored frozen in solution at -20° C for up to 2 weeks. A decrease in activity will occur if stored in solution. Use only freshly reconstituted protease for best results.

Endoproteinase LysC

#P8109S

20 μg 149 €

XX-Lys[▼]XXX

- Ideal for proteomic analysis by mass spectrometry
- Free of contaminating proteases
- Best suited for peptide identification

Description: LysC is a serine endoproteinase, isolated from *Lysobacter enzymogenes*, that cleaves peptide bonds C-terminal to lysine residues. LysC is a sequencing grade enzyme and is suitable for proteomics and glycobiology applications.

Note: Resuspend in 200 μ l double-distilled water to make 100 ng/ μ l LysC solution in 10 mM Tris-HCl, pH 8.0. The solution can be stored at 4°C for several days or in single-use aliquots at -20°C for several months.

37°

Source: Isolated from *Lysobacter enzymogenes*.

Molecular Weight: 30,000 daltons

Reconstitution: Endoproteinase LysC should be reconstituted in 200 µl double-distilled water to make a 100 ng/µl solution in 10 mM Tris-HCl, pH 8.0. Rapid autolysis is a function of enzyme concentration.

Proteinase K, Molecular Biology Grade

#P8107S

2 ml78 €

- Isolation of plasmid and genomic DNA
- Isolation of RNA
- Inactivation of RNases, DNases and enzymes in reactions

Description: Proteinase K is a subtilisin-related serine protease that will hydrolyze a variety of peptide bonds.

Source: Engyodontium album (Tritirachium album)

Reaction Conditions: Proteinase K is active in a wide range of buffers including all NEB specific restriction endonuclease buffers. It is highly active between pH 7.5 and 12.0 and temperatures between 20 and 60°C. Proteinase K is also active in chelating agents such as EDTA and activity is stimulated in up to 2% SDS or 4 M urea.

Calcium is important for thermostability of Proteinase K but it is not required for catalysis, therefore Proteinase K is also active in buffers containing chelating agents such as EDTA.

Molecular Weight: 28,900 daltons

Unit Definition: One unit will digest urea-denatured hemoglobin at 37°C (pH 7.5) per minute to produce equal absorbance as 1.0 µmol of L-tyrosine using Folin & Ciocalteu's phenol reagent.

Concentration: 800 units/ml

Thermolabile Proteinase K

#P8111S

30 units 150 €

- Isolation of plasmid and genomic DNA
- Isolation of RNA
- Inactivation of RNases, DNases and enzymes in reactions
- Removal of enzymes from DNA to improve cloning efficiency
- PCR purification

Description: Thermolabile Proteinase K is an engineered, subtilisin-related serine protease that will hydrolyze a variety of peptide bonds.

Source: Cloned from *Engyodontium album (Tritirachium* album), mutagenized to increase thermolability of the enzyme and expressed in K. lactis.

Molecular Weight: 29,000 daltons



Reaction Conditions: Thermolabile Proteinase K is active in a wide range of buffers. It is highly active between pH 7.0-9.5 and temperatures 20-40°C. It is active in chelating agents such as EDTA and activity is stimulated in up to 1% SDS.

Unit Definition: One unit is defined as the amount of enzyme required to release 1.0 µmol of 4-nitroaniline per minute from N-Succinyl-Ala-Ala-Pro-Phe-p-nitroanilide at 25°C, in a total reaction volume of 105 µl.

Concentration: 120 units/ml

Factor Xa Protease

#P8010S #P8010L

50 μg71 € 250 μg 284 €

Ile-Glu/Asp-Gly-Arg▼

Description: Factor Xa cleaves after the arginine residue in its preferred cleavage site Ile-(Glu or Asp)-Gly-Arg. It will sometimes cleave at other basic residues, depending on the substrate conformation. The most common secondary site, among those that have been sequenced, is Gly-Arg. There seems to be a correlation between proteins that are unstable in E. coli and those that are cleaved by Factor Xa at secondary sites; this may indicate that these proteins are in a partially unfolded state. Factor Xa will not cleave a site followed by proline or arginine.

Source: Purified from bovine plasma and activated by treatment with the activating enzyme from Russell's viper venom.

Molecular Weight: The predominant form of Factor Xa has a molecular weight of approximately 43,000 daltons, consisting of two disulfide-linked chains of approximately 27,000 daltons and 16,000 daltons. On SDS-PAGE, the reduced chains have apparent molecular weights of 30,000 daltons and 20,000 daltons.

Unit Definition: 1 µg of Factor Xa will cleave 50 µg of test substrate to 95% completion in 6 hours or less. Unit assay conditions can be found at www.neb.com.

Concentration: 1 mg/ml

Removal: Factor Xa will bind specifically to benzamidine-agarose (e.g., GE Life Sciences #17-5143-02).











Enterokinase, light chain

#P8070S 480 units 112 € #P8070L 2,560 units 448 €

Asp-Asp-Asp-Asp-Lys[▼]

Description: Enterokinase is a specific protease that cleaves after the lysine at its cleavage site, Asp-Asp-Asp-Asp-Lys. It will sometimes cleave at other basic residues, depending on the conformation of the protein substrate.

Source: Purified from *K. lactis* containing a clone of the light chain of the bovine enterokinase gene.

Enterokinase will not cleave a site followed by proline.

Molecular Weight: The molecular weight of the light chain of enterokinase is 26,300 daltons. Its apparent molecular weight on SDS-PAGE is 31 kDa.

R{{

Unit Definition: One unit is defined as the amount of enzyme required to cleave 25 μg of a MBP-EK-paramyosin-ΔSal substrate to 95% completion in 16 hours at 25°C in a total reaction volume of 25 μl. Unit assay conditions can be found at www.neb.com.

Concentration: 16,000 U/ml

Removal: Enterokinase will bind specifically to trypsin inhibitor agarose (e.g., Sigma T-0637).

Furin

#P8077S 50 units 138 € #P8077L 250 units 552 €

Arg-X-X-Arg▼

Description: Furin is an ubiquitous subtilisin-like proprotein convertase. It is the major processing enzyme of the secretory pathway and is localized in the trans-golgi network. Substrates of Furin include blood clotting factors, serum proteins and growth factor receptors such as the insulin-like growth factor receptor. The minimal cleavage site is Arg-X-X-Arg[▼]. However, the enzyme prefers the site Arg-X-(Lys/Arg)-Arg[▼]. An additional arginine at the P6 position appears to enhance cleavage. Furin is inhibited by EGTA, α1-Antitrypsin Portland and polyarginine compounds.

Note: The ability to cleave a particular substrate appears to depend on its tertiary structure as well as on the amino acids immediately surrounding the cleavage site.

RX

Source: Isolated from *Spodoptera frugiperda* (Sf9) cells infected with recombinant baculovirus carrying truncated human furin (kindly provided by R. Fuller).

Molecular Weight: The calculated molecular weight of truncated human furin is 52,700 daltons. Its apparent molecular weight in SDS-PAGE gels is 57,000 daltons.

Unit Definition: One unit is defined as the amount of furin required to cleave 25 μ g of a MBP-FN-paramyosin- Δ Sal substrate to 95% completion in 6 hours at 25°C in a total reaction volume of 25 μ l. Unit assay conditions can be found at www.neb.com.

Concentration: 2,000 units/ml

NEW

TEV Protease

#P8112S 1,000 units 105 €

ENLYFQ*(G/S)

Description: TEV Protease, also known as Tobacco Etch Virus (TEV) Protease, is a highly specific cysteine protease that recognizes the amino-acid sequence Glu-Asn-Leu-Tyr-Phe-Gln-(Gly/Ser) and cleaves between the Gln and Gly/Ser residues. It is often used for the removal of affinity purification tags such as maltose-binding protein (MBP) or poly-histidine from fusion proteins. TEV Protease has a 7xHis-tag for easy removal from a reaction using nickel affinity resins and has been engineered for greater performance.

RX 30° \

Source: Cloned from Tobacco Etch Virus and expressed in *E. coli*.

Molecular Weight: 28,000 daltons

Unit Definition: 1 unit of TEV Protease will cleave 2 μg of MBP-fusion protein, MBP5-TEV-paramyosin ΔSal , to 95% completion in a total reaction volume of 10 μl in 1 hour at 30°C in 50 mM Tris-HCl (pH 7.5 @ 25°C) with 0.5 mM EDTA and 1 mM DTT.

Concentration: 10,000 units/ml

Properties of Protein Phosphatases from NEB

The significance of protein phosphorylation in regulating the function and activity of protein factors and enzymes is now well established. Analysis of the presence of such phosphorylation, and its attendant effects, is often aided by removal of the protein phosphate groups by phosphatases.

Lambda Protein Phosphatase (Lambda PP)

#P0753S 20,000 units 144 € #P0753L 100,000 units 576 €

Companion Products:

p-Nitrophenylphosphate (PNPP)

Non-specific substrate for protein, alkaline and acid phosphatases

#P0757S 1 ml30 € #P0757L 5 ml120 €

Sodium Orthovanadate (Vanadate) General inhibitor for protein phosphotyrosyl specific phosphatases

#P0758S 1 ml29 € #P0758L 5 ml116 € **Description:** Lambda Protein Phosphatase (Lambda PP) is a Mn²-dependent protein phosphatase with activity towards phosphorylated serine, threonine and tyrosine residues. Lambda-PPase is active on phosphorylated histidine residues.

Source: Isolated from a strain of *E. coli* that carries the bacteriophage lambda ORF221 open reading frame (kindly provided by Dr. D. Barford).

Reaction Conditions:

1X NEBuffer for Protein MetalloPhosphatases (PMP). Supplement with 1 mM MnCl₂. Incubate at 30°C. Heat inactivation: 65°C for 1 hour in the presence of 50 mM EDTA.

RN 30° 65

Reagents Supplied with Enzyme:

10X NEBuffer for Protein MetalloPhosphatases (PMP) 10X MnCl₂ (10 mM)

Molecular Weight: 25,000 daltons

Unit Definition: One unit is defined as the amount of enzyme that hydrolyzes 1 nmol of ρ -Nitrophenyl Phosphate (50 mM, NEB #P0757) in 1 minute at 30°C in a total reaction volume of 50 μ l.

Concentration: 400,000 units/ml



Lydia has been a member of the Marketing Communications Team for over 3 years. She is also the host of the NEB podcast, Lessons from Lab & Life.











256

GLYCOBIOLOGY & PROTEIN TOOLS

Protein Kinases

The reversible addition of phosphate groups to proteins is important for the transmission of signals within eukaryotic cells and, as a result, protein phosphorylation and dephosphorylation regulate many diverse cellular processes. As the number of known protein kinases has increased at an ever-accelerating pace, it has become more challenging to determine which protein kinases interact with which substrates in the cell. The determination of consensus phosphorylation site motifs by amino acid sequence alignment of known substrates has proven useful in this pursuit. These motifs can be helpful for predicting phosphorylation sites for specific protein kinases within a potential protein substrate.

Since the determinants of protein kinase specificity involve complex 3-dimensional interactions, these motifs, short amino-acid sequences describing the primary structure around the phosphoacceptor residue, are a significant oversimplification of the issue. They do not take into account possible secondary and tertiary structural elements, or determinants from other polypeptide chains or from distant locations within the same chain. Furthermore, not all of the residues described in a particular specificity motif may carry the same weight in determining recognition and phosphorylation by the kinase. As a consequence, they should be used with some caution.

On the other hand, many of the residues within these consensus sequences have in fact proven to be crucial recognition elements, and the very simplicity of these motifs has made them useful in the study of protein kinases and their substrates. In addition to the prediction of phosphorylation sites, short synthetic oligopeptides based on consensus motifs are often excellent substrates for protein kinase activity assays.

The table below summarizes the specificity motifs for protein kinases that are available from NEB. Phosphoacceptor residue is indicated in red, amino acids which can function interchangeably at a particular residue are separated by a slash (/), and residues which do not appear to contribute strongly to recognition are indicated by an "X".

Visit www.neb.com for more information on protein kinases from NEB, including detailed information on recognition determinants.

PROTEIN KINASE	NEB#	RECOGNITION DETERMINANT	SIZE	PRICE
cAMP-dependent Protein Kinase (PKA), Catalytic Subunit	P6000S/L	R-R-X-S/T Y Y= hydrophobic residue	100,000/500,000 units	145 €/580 €
Casein Kinase II (CK2)	P6010S/L	S-X-X-E/D	10,000/50,000 units	145 €/580 €

D = aspartic acid, E = glutamic acid, R = arginine, S = serine, T = threonine, Y = tyrosine, X = any amino acid

Note: More specific information on recognition determinants for each kinase can be found on the corresponding product page at www.neb.com.


Companion Product:

Adenosine-5' Triphosphate (ATP)

#P0756S 1.0 ml35 € #P0756L 5.0 ml138 €

Endoglycoceramidase I (EGCase I)

#P0773S 150 mU73 €



This is an Enzyme for Innovation (EFI). EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit www.neb.com/ EnzymesforInnovation to view the full list.

Description: Endoglycoceramidase I (ECGase I) catalyzes the hydrolysis of the β-glycosidic linkage between oligosaccharides and ceramides in various glycosphingolipids. One unit of R. triatomea EGCase I is defined as the amount of enzyme required to hydrolyze 1 µmol of ganglioside GM1a per minute at 37°C.

Source: EGCase I is isolated from a strain of *E. coli*, which contains the cloned EGCase I gene from Rhodococcus triatomea.

RN 😻 37° 🛞

Reaction Conditions: 1X EGCase I Reaction Buffer. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Reagents Supplied: 10X EGCase I Reaction Buffer

Molecular Weight: 50,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to hydrolyze 1 µmol of ganglioside GM1a per minute at 37°C.

Concentration: 6 units/ml



The Importance of the Top Predator

The collapse of an ecosystem, if we are not observant, can take decades to fully realize and understand. In the mid-1990s, the banks of the Lamar River in Yellowstone National Park, USA, were eroded, streams had widened, willow shoots that had previously flourished by the river in the 1920s were nonexistent, aspen trees were not regenerating and beavers were scarce. How is it that the collapse of this ecosystem could be the result of an event that occurred in the 1920s?

The case in Yellowstone National Park began as an investigation into the disappearance of aspen trees. Scientists examined the age of the remaining trees by drilling into the cores and establishing age versus diameter relationships. They discovered that the aspen tree had not regenerated in the past 70 years, and by back-dating, it was realized that this 70-year hiatus correlated with the disappearance of the wild wolf from Yellowstone. The downstream effects that occurred as a result of the sudden absence of the wild wolves highlights the complexity of a balanced ecosystem.

A top predator, also known as an apex predator, is at the head of the food chain. In the early 1900s, Yellowstone's top predator, the wild wolf, was viewed as a threat and was hunted until it was completely eliminated from the park by the 1920s.

Pieces of the puzzle started to fit together when scientists examined the effect of unchecked population growth of elk, the primary prey of the wild wolf. Vegetation, including aspen and willow tree saplings, were overgrazed by the elk. The riverside willow had provided material for beavers to build their dams, and in turn, the dams provided more water and nutrients for growth of the riverside vegetation. Additionally, the beaver's protective dams had tempered the seasonal changes in the river flow. With no material to build dams, the beavers also disappeared.

In 1995, in a highly controversial move, 31 wild wolves were relocated from Canada to Yellowstone. Their movement and behavior were observed as they hunted elk and deer. The elk carcasses not only fed the wolves, but also coyotes, ravens, magpies, eagles and finally, insects.

The elk were also tracked, and it was observed that they avoided the gorges and valleys where they were easy prey. With less elk and deer grazing on willow and aspen, the long-gone vegetation started to regenerate; more berries and insects followed, and then various bird species. Beaver families moved back into the area and used the willow to build their dams, which created a habitat for otter, muskrat and reptiles. The wolves also killed coyotes, and so the rabbit and mice populations grew, which fed hawks, foxes, badgers and weasels. The regeneration of the riverside vegetation stabilized the river banks against erosion, and subsequently, the rivers narrowed and became more fixed in their course. Pools began to form, creating habitats for other organisms.

The Yellowstone example gives insight into the trophic cascade that can result from removing a key organism from a balanced landscape and ecosystem, and while a top predator kills certain species, the downstream effects give life to many other species. The regeneration of this amazing ecosystem shows the incredible ability of a vast biome to restore itself, and the lessons learned here could lead to better predator management decisions in other locations.

Epigenetics



Simplify your epigenetics research with EpiMark® validated products.

Epigenetics is the study of heritable changes in the phenotype of a cell or organism that are not encoded in the DNA of the genome. The molecular basis of an epigenetic profile arises from covalent modifications of the protein and DNA components of chromatin. The epigenetic profile of a cell often dictates cell memory and cell fate and, thus influences mammalian development.

The epigenetic code is hypothesized to be the combined effects of histone modifications and DNA methylation on gene expression. While the genetic code for an individual is the same in every cell, the epigenetic code is tissue- and cell-specific, and may change over time as a result of aging, disease or environmental stimuli (e.g., nutrition, life style, toxin exposure). Cross-talk between histone modifications, DNA methylation or RNAi pathways are being studied in such areas as cancer, X chromosome inactivation, and imprinting.

For over 40 years, New England Biolabs has been committed to understanding the mechanisms of restriction and methylation of DNA. This expertise in enzymology has led to the development of a suite of validated products for epigenetics research. These unique solutions to study DNA and histone modifications are designed to address some of the challenges of the current methods. All NEB products pass stringent quality control assays to ensure the highest level of functionality and purity.

Featured Products

- 262 EpiMark 5-hmC and 5-mC Analysis Kit
- 263 EpiMark Methylated DNA Enrichment Kit
- 267 EpiMark Nucleosome Assembly Kit
- 273 NEBNext Methyl-seq

Featured Tools & Resources

- Videos of NEB Scientists
 Discussing Epigenetics
- Epigenetics-related FAQs
- Feature Articles
 - Visit www.EpiMark.com to view an interactive tutorial explaining the phenomenon of epigenetics at the molecular level.





269

Histone H4 Human, Recombinant

Recombinant Enzyme

273

Methyl-seq (Unique Dual Index Primer Pairs)

EpiMark® 5-hmC and 5-mC Analysis Kit

#E3317S 20 reactions 248 €

- Reproducible quantitation of 5-hmC and 5-mC within a specific loci
- Easy-to-follow protocols
- Compatible with existing techniques
- Amenable to high throughput

Visit EpiMark.com to view a video tutorial for this kit.


Description: The EpiMark 5-hmC and 5-mC Analysis Kit is a simple and robust method for the identification and quantitation of 5-hydroxymethylcytosine (5-hmC) and 5-methylcytosine (5-mC) within a specific DNA locus. This enzymatic approach utilizes the differential methylation sensitivity of the isoschizomers Mspl and Hpall in a simple 3-step protocol.

Briefly, genomic DNA is treated with T4 phage β-glucosyltransferase and UDP-glucose, which glucosylates all 5-hmC present. DNA is digested with Mspl and Hpall, two isoschizomers with different methylation sensitivity. Endpoint or real time PCR can then be used to identify and quantitate the different methylation states. Designed to simplify methylation analysis, the EpiMark Kit expands the potential for new biomarker discovery.

RX Epi NEB4

The EpiMark 5-hmC and 5-mC **Analysis Kit Includes:**

- T4 Phage β-glucosyltransferase
- UDP-Glucose
- Mspl
- Hpall
- Proteinase K
- Control DNA (unmodified, 5-mC and 5-hmC)
- Forward and reverse control primer mix
- NEBuffer 4



Analysis of the different methylation states in Balb/C mouse tissue samples (locus 12) using the EpiMark 5-hmC and 5-mC Analysis Kit. A) Endpoint PCR of the 6 different reactions needed for methylation analysis. The boxed lanes indicate the presence of 5-hmC. B) Real time PCR data was used to determine amounts of 5-hmC and 5-mC present. The results demonstrate a variation in 5-hmC levels in the tissue sources indicated.

T4 Phage β -glucosyltransferase

#M0357S 500 units #M0357L 2,500 units 352 €

- Glucosylation of 5-hydroxymethylcytosine in DNA
- Immunodetection of 5-hydroxymethylcytosine in DNA
- Labeling of 5-hydroxymethylcytosine residues by incorporation of [3H]- or [14C]-glucose into 5-hmC-containing DNA acceptor after incubation with [3H]- or [14C]-UDP-GIC
- Detection of 5-hydroxymethylcytosine in DNA by protection from endonuclease cleavage

Description: T4 Phage β-glucosyltransferase specifically transfers the glucose moiety of uridine diphosphoglucose (UDP-Glc) to the 5-hydroxymethylcytosine (5-hmC) residues in double-stranded DNA, making beta-glucosyl-5-hydroxymethylcytosine.

Reaction Conditions: 1X NEBuffer 4 and 40 μM UDP-Glucose. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

RX Epi NEB4 dil B 37° 165

Reagents Supplied with Enzyme:

10X NEBuffer 4

50X UDP-Glucose (2 mM)

Unit Definition: One unit is defined as the amount of enzyme required to protect 0.5 µg T4gt-DNA against cleavage by Mfel restriction endonuclease.

Concentration: 10,000 units/ml















EpiMark Methylated DNA Enrichment Kit


#E2600S

25 reactions 451 €

- High-affinity binding provides greater sensitivity
- Elution in a small volume simplifies downstream applications
- Easy-to-follow protocol yields enriched fractions in less than 2 hours
- Enriched methylated DNA fragments can be easily ligated to doublestranded adaptors for next generation sequencing
- Highly pure product from a wide range of input DNA concentrations

Description: The EpiMark Methylated DNA Enrichment Kit enables the enrichment of double-stranded CpG methylated DNA based on CpG methylation density. It utilizes the methyl-CpG binding domain of human MBD2a protein as a capture agent. The protein is fused to the Fc tail of human IgG1 (MBD2a-Fc), which is coupled to Protein A Magnetic Beads (MBD2a-Fc/Protein A Bead). This stable complex will selectively bind double-stranded methylated CpG containing DNA. The high binding affinity of the MBD2a-Fc coupled beads and optimized reagents increases sensitivity and accuracy. This kit contains all the individual components necessary to achieve enrichment in less than two hours using a four step process.

Enrichment Workflow



Epi

The EpiMark Methylated DNA **Enrichment Kit Includes:**

- MBD2-Fc protein
- Protein A Magnetic Beads
- Bind/wash buffer
- High-salt elution buffer
- Fragmented HeLa DNA
- Line element primers for methylated-locus controls
- RPL30 primers for unmethylated locus contols
- MirA locus control primers

EpiMark N6-Methyladenosine Enrichment Kit

#E1610S

20 reactions 420 €

Description: The EpiMark N6-Methyladenosine Enrichment Kit can be used to enrich m6A modified RNA in immunoprecipitation protocols for downstream analysis by next-generation RNA sequencing or RT-qPCR. Modified RNA is isolated from a fragmented RNA sample by

binding to the N6-Methyladenosine antibody attached to Protein G Magnetic Beads. After multiple wash and clean-up steps, the enriched RNA is eluted in nuclease-

free water and is ready for further analysis.

See page 204 for more information.

See page 68 for more information.

EpiMark Bisulfite Conversion Kit

#E3318S

48 reactions 160 €

- Complete conversion of unmodified cytosine to uracil
- Easy-to-follow protocol
- Reliable and consistent results
- Purification columns included

Description: This technique can reveal the methylation status of every cytosine residue, and it is amenable to massively parallel sequencing methods. Bisulfite conversion involves the conversion of unmodified cytosines to uracil, leaving the modified bases 5-mC and 5-hmC. The EpiMark Bisulfite Conversion Kit is designed for the detection of methylated cytosine, using a series of alternating cycles of thermal denaturation, followed by incubation with sodium bisulfite. This kit includes all the reagents necessary for complete bisulfite conversion, including spin columns. Amplification of bisulfite-treated samples can then be performed using EpiMark Hot Start Taq DNA Polymerase.

The EpiMark Bisulfite Conversion Kit Includes:

Sodium metabisulfite

Epi

Epi

- Solubilization buffer
- Desulphonation reaction buffer
- EpiMark spin columns with 2 ml collection tubes
- Binding buffer
- Wash buffer
- Elution buffer

R\\\ Epi

EpiMark Hot Start Taq DNA Polymerase

#M0490S #M0490L

100 reactions 61 € 500 reactions 244 €

dilA Diluent Buffer



Description: EpiMark Hot Start *Tag* DNA Polymerase is a mixture of *Taq* DNA Polymerase and a temperature sensitive, aptamer-based inhibitor. The inhibitor binds reversibly to the enzyme, inhibiting polymerase activity at temperatures below 45°C, but releases the enzyme during normal PCR cycling conditions. This permits assembly of reactions at room temperature. An advantage of the

aptamer-based hot start mechanism is that it does not converted DNA

require a separate high temperature incubation step to activate the enzyme. The advanced aptamer-based hot-start activity coupled with the supplied optimized reaction buffer makes the EpiMark Hot Start Tag DNA Polymerase an excellent choice for use on bisulfite-

Methylation-Dependent Restriction Enzymes

- Specificity to epigenetically-relevant DNA modifications (5-mC and 5-hmC)
- Easy-to-follow protocols
- Less harsh than bisulfite conversion
- Simplified data analysis

The EpiMark Suite of products has been validated for use in epigenetics applications. Visit **EpiMark.com** for more information

Many restriction enzymes are sensitive to DNA methylation states. Cleavage can be blocked or impaired when a particular base in the recognition site is modified. The MspJI family of restriction enzymes are dependent on methylation and hydroxymethylation for cleavage to occur (1). These enzymes excise 32-base pair DNA fragments containing a centrally located 5-hmC or 5-mC modified residue that can be extracted and sequenced. Due to the known position of this epigenetic modification, bisulfite conversion is not required prior to downstream analysis.

These EpiMark validated, methylation-dependent restriction enzymes expand the potential for mapping epigenetic modifications and simplify the study of DNA methylation. Additionally, they provide an opportunity to better understand the role of 5-hydroxymethylcytosine in the genome.

Reference:

(1) Cohen-Karni, D. et al. (2011) *PNAS*, 108, 11040–11045.

AbaSI

#R0665S 1.000 units113 €

 ${}^{x}C = {}^{ghm}C, {}^{hm}C, {}^{m}C \text{ or } C$

Description: AbaSI is a DNA modification-dependent endonuclease that recognizes 5-glucosylhydroxymethylcytosine (ghmC) in double-stranded DNA and cleaves 11–13 bases 3′ from the modified C leaving a 2–3 base 3′ overhang. The enzyme only cleaves if there is a G residue 20–23 nucleotides 3′ from the modified C. AbaSI also recognizes 5-hydroxymethylcytosine (hmC) at a much lower efficiency. It does not recognize DNA with 5-methylcytosine (mC) or unmodified cytosine.

CutSmart Rii Epi dii C 25° 165

Reaction Conditions: 1X CutSmart Buffer + 1 mM DTT (supplied), 25°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

FspEI

#R0662S 200 units 106 €

 $5' \dots G^m G(N)_{12}^{\vee} \dots 3'$ $3' \dots G(N)_{16}^{\vee} \dots 5'$ **Description:** FspEI is a modification-dependent endonuclease which recognizes C^mC sites and generates a double-stranded DNA break on the 3' side of the modified cytosine at N_{12}/N_{16} . Recognized cytosine modifications include C5-methylation (5-mC) and C5-hydroxymethylation (5-hmC).

CutSmart Ri Epi dil B 37° 👑

Reaction Conditions: 1X CutSmart Buffer + 1X Enzyme Activator Solution (supplied), 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 5.000 units/ml

Note: Star activity may result from extended digestion.

LpnPI

#R0663S 200 units108 €

5′... C^mC D G (N)₁₀... 3′ 3′... G G H C (N)₁₄... 5′ **Description:** LpnPI is a modification-dependent endonuclease which recognizes C^mCDG sites and generates a double-stranded DNA break on the 3' side of the modified cytosine at N₁₀/N₁₄. Recognized cytosine modifications include C5-methylation (5-mC) and C5-hydroxymethylation (5-hmC).

CutSmart Rill Epi dil B 37° 📆

Reaction Conditions: 1X CutSmart Buffer + 1X Enzyme Activator Solution (supplied), 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

Note: Star activity may result from extended digestion.

MspJI

#R0661S 200 units108 € #R0661L 1,000 units432 €

5′...[™]C N N R (N)₉ ♥...3′ 3′... G N N Y (N)₁₃...5′ **Description:** MspJI is a modification-dependent endonuclease that recognizes m CNNR sites and generates a double-stranded DNA break on the 3′ side of the modified cytosine at N_g/N₁₃. The recognized cytosine modifications include C5-methylation (5-mC) and C5-hydroxymethylation (5-hmC).

CutSmart RR Epi dil B 37° 65

Reaction Conditions: 1X CutSmart Buffer + 1X Enzyme Activator Solution (supplied), 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

Note: Star activity may result from extended digestion.

McrBC

#M0272S 500 units 70 € #M0272L 2,500 units 280 €

5´...Pu^mC (N₄₀₋₃₀₀₀) Pu^mC... 3´

- Determination of the methylation state of CpG dinucleotides
- Detection of cytosine methylated DNA

Description: McrBC is an endonuclease that cleaves DNA containing methylcytosine* on one or both strands. McrBC will not act upon unmethylated DNA. Sites on the DNA recognized by McrBC consist of two half-sites of the form (G/A)mC. These half-sites can be separated by up to 3 kb, but the optimal separation is 55–103 base pairs. McrBC requires GTP for cleavage, but in the presence of a non-hydrolyzable analog of GTP, the enzyme will bind to methylated DNA specifically, without cleavage. McrBC will act upon a pair of PumCG sequence elements, thereby detecting a high proportion of methylated CpGs, but will not recognize Hpall/Mspl sites (CCGG) in which the internal cytosine is methylated.

Reaction Conditions: 1X NEBuffer $2 + 200 \mu g/ml$ BSA + 1 mM GTP, 37°C. Heat inactivation: 65°C for 20 minutes.

RN Epi NEB 2 BSA dil B 37° 165

Reagents Supplied with Enzyme:

10X NEBuffer 2 100X BSA 100X GTP (100 mM)

Control Plasmid DNA

Unit Definition: One unit is defined as the amount of enzyme required to cleave $0.5~\mu g$ of a plasmid containing multiple McrBC sites in 1 hour at $37^{\circ}C$ in a total reaction volume of $50~\mu l$.

Concentration: 10,000 units/ml

Note: McrBC makes one cut between each pair of halfsites, cutting close to one half-site or the other, but cleavage positions are distributed over several base pairs approximately 30 base pairs from the methylated base. Therefore, when multiple McrBC half-sites are present in DNA (as is the case with cytosine-methylated genomic DNA), the flexible nature of the recognition sequence results in an overlap of sites, producing a smeared, rather than a sharp, banding pattern.

Additional Restriction Enzymes for Epigenetic Analysis

Dnnl

See pages 34, 38, 41 for more information.

Visit pages 334–336 for a complete list of methylation-sensitive restriction enzymes.

Methylation sensitive restriction enzymes can be used to generate fragments for further analysis. When used in conjunction with an isoschizomer that has the same recognition site, but is methylation insensitive, information about methylation status can be obtained.

#R0176S #R0176L	1,000 units 65 € 5,000 units 260 €
DpnII #R0543S #R0543L	1,000 units 72 € 5,000 units 288 €
for high (5X) concentration #R0543T #R0543M	1,000 units 72 € 5,000 units 288 €

Hpall #R0171S #R0171L	2,000 units 65 € 10,000 units 260 €
for high (5X) concentration #R0171M	on 10,000 units 260 €
MspI #R0106S #R0106L	5,000 units 64 € 25,000 units 256 €
for high (5X) concentration #R0106T #R0106M	5,000 units 64 € 25,000 units 256 €

Single Letter Code: R = A or G Y = C or T M = A or C K = G or T S = C or G W = A or T H = A or C or T (not G) B = C or G or T (not A)

V = A or C or G (not T) D = A or G or T (not C) N = A or C or G or T

^{* 5-}methylcytosine, 5-hydroxymethylcytosine or N4-methylcytosine


Histones

- Purification and characterization of enzymes that modify histone proteins
- Octamer and nucleosome modification
- Carrier chromatin immunoprecipitation (CChIP)
- High throughput studies

In eukaryotes, nuclear DNA is assembled into chromatin by nucleoprotein complexes. The primary unit of chromatin, a nucleosome core particle (NCP), is an octamer complex made up of two molecules each of Histone H2A, H2B, H3 and H4 and approximately 147 base pairs of nuclear DNA. Histone H1 further condenses the DNA by binding the linker segments between NCP complexes (1,2). Histones undergo diverse post-translational modification including acetylation, phosphorylation, mono-, di- or tri-methylation, ubiquitination, isomerization and ADP-ribosylation. Through their potential combinatorial sequences on a given histone and their reversibility, these modifications dynamically restrict or recruit numerous other proteins or protein complexes onto chromatin (3). The study of their roles in gene regulation (4), cellular stress events (4), aging and DNA repair (5) is revealing the multiple functions of histone modifications in the fate of a cell. Additional variability is incorporated into the system by histone variants. Acting individually or combinatorially in conjunction with DNA modification, histone modifications and histone variants are thought to establish an epigenetic code or epigenetic mechanism of gene regulation (3).

In total, seven human histones, including three histone H3 variants (see alignment below), have been individually cloned and expressed in E. coli and then highly purified from cell extracts at NEB. Mass spectrometry analysis demonstrates that these recombinant histones are free of post-translational modifications. These histones are ideal substrates for the purification and characterization of histone modifying enzymes.

To aid in studying intact nucleosomes, we now offer the EpiMark Nucleosome Assembly Kit. The precise mixing of a 2:1 ratio of Histone H2A/H2B Dimer to Histone H3.1/H4 Tetramer generates a recombinant human histone octamer, and in the presence of DNA forms nucleosomes (7,8). Enzymes that are unable to modify individual histones or DNA may be active on these nucleosome core particles, the histone dimer, or the histone tetramer (9,10). The NCPs also may be used as carrier chromatin in CChIP (carrier chromatin immunoprecipitation) assays (11). The recombinant human histone dimer and recombinant human histone tetramer are also available as separate products.



Experience the purity of Histones from NEB, SDS-PAGE analysis of the histones available from NEB.

- Histone H1° (NEB #M2501) 1 μg
- Histone H2A (NEB #M2502) 1 μg
- Histone H2B (NEB #M2505) 1 μg
- Histone H3.1 (NEB #M2503) 1 μg
- Histone H3.2 (NEB #M2506) 1 μg
- 6: Histone H3.3 (NEB #M2507) 1 μg
- Histone H4 (NEB #M2504) 1 μg
- Histone H2A/H2B Dimer (NEB #M2508) 2 μg
- Histone H3.1/H4 Tetramer (NEB #M2509) 2 µg
- 10: NFB Protein Ladder

References:

- Kornberg, R.D. (1977) Annu. Rev. Biochem., 46, 931-954
- van Holde, K.E. (1989) Chromatin, 1-497.
- Kim, J.K., Samaranayake, M. and Pradhan, S. (2009) Cell. Mol. Life Sci., 66, 596-612.
- Huang, J. et al. (2006) Nature, 444, 629-632
- Pahlich, S., Zakaryan, R.P. and Gehring, H. (2006) Biochim. Biophys. Acta., 1764, 1890-1903.
- Hake, S.B. et al. (2006) J. Biol. Chem., 281, 559-568
- Luger et al. (1999) Methods in Mol. Biol., 119,
- Comb, D. and Mersha, F. unpublished observa-
- Li, Yan et al. (2009) J. Biol. Chem., 284, 34283-34295.
- Fung, J. et al. (2002) Current Biology, 12, 1086-1099.
- (11) O'Neill et al. (2006) Nature Genetics, 38, 835-841.

10	20	30	40	50	60	70		
ARTKQTARKS	TGGKAPRKQL	ATKAARKSAP	ATGGVKKPHR	YRPGTVALRE	IRRYQKSTEL	LIRKLPFQRL	histone	H3.1
ARTKQTARKS	TGGKAPRKQL	ATKAARKSAP	ATGGVKKPHR	YRPGTVALRE	IRRYQKSTEL	LIRKLPFQRL	histone	Н3.2
ARTKQTARKS	TGGKAPRKQL	ATKAARKSAP	${\color{red}{\bf S}}{\color{blue}{\bf T}}{\color{blue}{\bf G}}{\color{blue}{\bf V}}{\color{blue}{\bf K}}{\color{blue}{\bf F}}{\color{blue}{\bf H}}{\color{blue}{\bf R}}$	YRPGTVALRE	IRRYQKSTEL	LIRKLPFQRL	histone	H3.3
80	90	100	110	120	130			
VREIAQDFKT	DLRFQSSAVM	ALQEACEAYL	${\tt VGLFEDTNLC}$	AIHAKRVTIM	PKDIQLARRI	RGERA	histone	H3.1
VREIAQDFKT	DLRFQSSAVM	${\tt ALQEASEAYL}$	${\tt VGLFEDTNLC}$	${\tt AIHAKRVTIM}$	PKDIQLARRI	RGERA	histone	H3.2
				3 7 11 3 11 D 1 1 1 1 1 1	DEDICT ADDI	DOED !	histone	*** 2
VREIAQDFKT	DLRFQSAAIG	ALQEASEAYL	VGLFEDTNLC	ATHAKKVTIM	PKDIQLARRI	RGERA	nistone	H3.3

Sequence alignment of Human Histone variants H3.1, H3.2 and H3.3. Human Histone H3.1, 3.2 and 3.3 vary by only a few amino acids (changes are highlighted in red), but are associated with different biological functions (6)





EpiMark Nucleosome Assembly Kit

0.2 nmol.....115 €

#E5350S 20 reactions 596 €

Histone H3.1/H4 Tetramer Human, Recombinant

Components Sold Separately:

#M2509S 1 nmol.....234 €

Histone H2A/H2B Dimer Human, Recombinant
#M2508S 2 nmol.....234 €

Nucleosome Control DNA

Highly pure, recombinant system

#N1202S

- Pre-formed histone dimer and tetramer complexes simplify octamer formation
- Components stable for one year
- Ideal for ChIP Assay, HAT Assay and enzyme modification assays (e.g., methylation studies)

Description: This kit contains the components necessary to form an unmodified recombinant human nucleosome using your own target DNA or the supplied control DNA. The protocol requires the mixing of already formed and purified recombinant human Histone H2A/H2B Dimer and Histone H3.1/H4 Tetramer in the presence of DNA in high salt, followed by dialysis down to low salt, to make nucleosomes. One tetramer associates with two dimers to form the histone octamer on the DNA, generating a nucleosome. A method for assaying nucleosome formation by gel shift assay is also provided. These nucleosomes may serve as a better substrate for enzymes that are inactive on the DNA or one of the core histones alone. Each described reaction creates nucleosomes from ~50 pmol of a 208 bp DNA and may be scaled depending on the experiment.


Histone H2A/H2B Dimer Human, Recombinant is generated by refolding the denatured, purified subunits H2A and H2B, followed by gel filtration. Histone H3.1/H4 Tetramer Human, Recombinant is generated by refolding the denatured, purified subunits H3.1 and H4, followed by gel filtration. Both the dimer and tetramer are highly pure and are available separately for histone modification studies.

The Nucleosome Control DNA is a 208 base pair fragment from *Lytechinus variegates* 5SrDNA, and can be used for mononucleosome formation.

The EpiMark Nucleosome Assembly Kit Includes:

- Histone H2A/H2B Dimer
- Histone H3.1/H4 Tetramer
- Control DNA

R\ Epi



Gel shift assay to visualize nucleosome formation.

Samples from nucleosome assembly reactions were run on 6% polyacrylamide gel in 0.5X TBE. Lane 1: Low Molecular Weight DNA Ladder (NEB #N3233). Lane 2: Nucleosome Control DNA. Lane 3: 0.5:1 ratio of Octamer* to DNA. Lane 4: 1:1 ratio of Octamer* to DNA.

*Octamer = 2:1 mix of Histone H2A/H2B Dimer and Histone H3.1/H4 Tetramer.

Histone H1º Human, Recombinant

#M2501S


100 μg80 €

Description: Histone H1 acts on the linker region of polynucleosome DNA to condense the chromatin into structures of ~30 nm. It is not necessary for octamer or nucleosome core particle formation. Eight different Histone H1 proteins have been identified in the human genome. Histone H10 is a non replication-dependent histone that is highly expressed in terminally differentiated cells.

Synonyms: Histone H1.0, Histone H1(0), Histone H1'

Concentration: 1 mg/ml





Mass Spectrometry Analysis of Histone H1^o Human, Recombinant. The average mass calculated from primary sequence is 20731.53 Da.

Histone H2A Human, Recombinant

#M2502S


100 μg80 €

Description: Histone H2A combines with Histone H2B to form the H2A-H2B heterodimer. Two H2A/H2B heterodimers interact with an H3/H4 tetramer to form the histone octamer. Histone H2A is also modified by various enzymes and can act as a substrate for them. These modifications have been shown to be important in gene regulation.

Concentration: 1 mg/ml



Ri Epi



Mass Spectrometry Analysis of Histone H2A Human, Recombinant. The average mass calculated from primary sequence is 13990.28 Da.

Histone H2B Human, Recombinant

#M2505S


100 μg80 €

Description: Histone H2B combines with Histone H2A to form the H2A-H2B heterodimer. Two H2A/H2B heterodimers interact with an H3/H4 tetramer to form the histone octamer. Histone H2B is also modified by various enzymes and can act as a substrate for them. These modifications have been shown to be important in gene regulation.

Synonyms: Histone H2B/8, Histone H2B.1, Histone

H2B-GL105

Concentration: 1 mg/ml



Mass Spectrometry Analysis of Histone H2B Human, Recombinant. The average mass calculated from primary sequence is 13788.97 Da.

16.000

18.000

12.000

Histone H3.1 Human, Recombinant


#M2503S

100 μg80 €

Description: Histone H3 combines with Histone H4 to form the H3/H4 tetramer. Two H2A/H2B heterodimers interact with an H3/H4 tetramer to form the histone octamer. It is also modified by various enzymes and can act as a substrate for them. These modifications have been shown to be important in gene regulation.

Histone H3.1, an H3 variant that has thus far only been found in mammals, is replication dependent and is associated with gene activation and gene silencing.

Synonyms: Histone H3/a Concentration: 1 mg/ml



Mass Spectrometry Analysis of Histone H3.1 Human, Recombinant. The average mass calculated from primary sequence is 15272.89 Da.

Epi EpiMark Validated

NEB4 Optimal Buffer

Heat Inactivation

Cloned at NEB

Recombinant Enzyme

BSA Requires BSA

37° Incubation Temperature

dilA Diluent Buffer

Histone H3.2 Human, Recombinant

#M2506S

100 μg80 €


Description: Histone H3 combines with Histone H4 to form the H3/H4 tetramer. Two H2A/H2B heterodimers interact with an H3/H4 tetramer to form the histone octamer. It is also modified by various enzymes and can act as a substrate for them. These modifications have been shown to be important in gene regulation.

Histone H3.2, an H3 variant that is found in all eukaryotes except budding yeast, is replication dependent and is associated with gene silencing.

Synonyms: Histone H3/m, H3/o

Concentration: 1 mg/ml

Ri Epi



Mass Spectrometry Analysis of Histone H3.2 Human, Recombinant. The average mass calculated from primary sequence is 15258.06 Da.

Histone H3.3 Human, Recombinant

#M2507S

100 μg80 €


Description: Histone H3 combines with Histone H4 to form the H3/H4 tetramer. Two H2A/H2B heterodimers interact with an H3/H4 tetramer to form the histone octamer. It is also modified by various enzymes and can act as a substrate for them. These modifications have been shown to be important in gene regulation.

Histone H3.3, an H3 variant that is found in all eukaryotes from yeast to human, is replication and cell cycle phase independent and is the most common H3 in non-dividing cells. It has been shown to be enriched in covalent modifications associated with gene activation.

Synonyms: Histone H3.3A, H3F3, H3.3B

Concentration: 1 mg/ml

Ri Epi



Mass Spectrometry Analysis of Histone H3.3 Human, Recombinant. The average mass calculated from primary sequence is 15176.96 Da.

Histone H4 Human, Recombinant

#M2504S


100 μg80 €

Description: Histone H3 combines with Histone H4 to form the H3/H4 tetramer. Two H2A/H2B heterodimers interact with an H3/H4 tetramer to form the histone octamer. Histone H4 is also modified by various enzymes and can act as a substrate for them. These modifications have been shown to be important in gene regulation.

Synonyms: For HIST2H4 gene: H4/N, H4F2, H4FN

Concentration: 1 mg/ml

Ri Epi




Mass Spectrometry Analysis of Histone H4 Human, Recombinant. The average mass calculated from primary sequence is 11236.15 Da.

RX

See pages 114-116 for more information.

DNA Methyltransferases



CpG Methyltransferas	se (M.Sssl)	EcoGII Methyltra	insferase
#M0226S #M0226L		#M0603S	200 units 73 €
for high (5X) concentration		EcoRI Methyltra	nsferase
#M0226M	500 units 284 €	#M0211S	10,000 units 65 €
GpC Methyltransferas	se (M.CviPI)	Haelll Methyltra	nsferase
#M0227S	200 units 74 €	#M0224S	500 units 73 €
#M0227L	1,000 units 296 €		
		Hhal Methyltran:	sferase
Alul Methyltransferas	е	#M0217S	1,000 units 73 €
#M0220S	100 units 73 €		
		Hpall Methyltran	sferase
BamHI Methyltransfe	rase	#M0214S	100 units 69 €
#M0223S	100 units 74 €		
#M0223L	500 units 296 €	Mspl Methyltran	sferase
		#M0215S	100 units 73 €
dam Methyltransferas	se		
#M0222S	500 units 75 €	Taq I Methyltrans	sferase

Human DNA (cytosine-5) Methyltransferase (Dnmt1)

#M0222L

#M0230S 50 units 123 € #M0230L 250 units 492 €

5'... CG...3' 3′... G C ... 5′ Human Dnmt1 | CH₃ 5'... ĊG...3' 3'... G C ... 5' CH₃

CH₃

Description: Dnmt1 methylates cytosine residues in hemimethylated DNA at 5´...CG...3´. Mammalian Dnmt1 is believed to be involved in carcinogenesis, embryonic development and several other biological functions. The bulk of the methylation takes place during DNA replication in the S-phase of the cell cycle.

2,500 units 300 €

Reaction Conditions: 1X Dnmt1 Reaction Buffer. Supplement with 100 μg/ml BSA and 160 μM S-adenosylmethionine. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

RX Epi SAM 37° K

#M0219S

Reagents Supplied with Enzyme:

10X Dnmt1 Reaction Buffer 100X BSA

32 mM S-adenosylmethionine (SAM)

Unit Definition: One unit is the amount of enzyme required to catalyze the transfer of 1 pmol of methyl group to poly dl.dC substrate in a total reaction volume of 25 µl in 30 minutes at 37°C.

1,000 units 73 €

Concentration: 2,000 units/ml

Note: For DNA modification and protection applications, M.SssI (NEB #M0226) is preferred because it efficiently methylates both unmethylated and hemimethylated DNA substrates.

















EPIGENETICS

HeLa Genomic DNA

#N4006S

15 μg 110 €

- PCR, SNP analysis and southern blotting
- Genomic DNA library construction
- Control DNA for Methylation-specific PCR (MSP), Bisulfite sequencing, Methylationsensitive Single-Nucleotide Primer Extension (MS-SNuPE,) Combined Bisulfite Restriction Analysis (COBRA), Bisulfite treatment and PCR-single strand Confirmation Polymorphism Analysis (Bisulfite-PCR-SSCP/BiPS)

Description: HeLa (cervix adenocarcinoma) cells were grown to confluency in DMEM plus 10% fetal bovine serum. Genomic DNA was isolated by a standard genomic purification protocol, phenol extracted and equilibrated to 10 mM Tris-HCI (pH 7.5) and 1 mM EDTA.

Note: NEB is one of the only suppliers of this product.

A₂₆₀/₂₈₀ Ratio: 1.87

5-methyl-dCTP

#N0356S

1 μmol in 0.1 ml70 €

Generation of fully methylated cytosinesubstituted DNA

Description: Cytosine modification at carbon 5 (C5) represents an important epigenetic modification. It is also believed to be the starting substrate for the Ten-Eleven Translocation (TET) family of enzymes and their associated oxidation pathways. 5-methyl-dCTP offers the ability to enzymatically make defined fully methylated cytosine-substituted DNA, which can be used for a variety of biochemical and cellular applications.

5-methyl-dCTP (2´-deoxy-5-methylcytidine 5'-triphosphate) is supplied as a 10 mM solution at pH 7. Nucleotide concentration is determined by measurements of absorbance at 260 nm.

Formula: $C_{10}H_{15}N_3O_{13}P_3$ (free acid) Concentration: 10 mM solution

Molecular weight: 481.1 (acid form)

Diluent Compatibility: Can be diluted using sterile distilled water, preferably Milli-Q® water, or can be diluted using sterile TE [10 mM Tris-HCI, 1 mM EDTA

(pH 7.5)]



NEB's Golf Committee has been organizing our charity golf tournament for over 10 years, with the proceeds benefiting Ipswich High School students. Pictured are committee members Kari, Deana, Ted, Tanya and Karen.

NEBNext® Reagents for ChIP-Seq Library Preparation

Epi

NEBNext reagents are a series of highly pure reagents that facilitate library preparation of DNA or RNA for downstream applications, such as next generation sequencing and expression library construction. These reagents undergo stringent quality controls and functional validation, ensuring maximum yield, convenience and value.

For sample preparation of a ChIP-Seq DNA library, NEB offers kits, oligos and modules that support standard or fast workflows. To decide which products to choose, use the selection chart below. For more information on NEBNext reagents for library preparation, see pages 134-161.





Katrina is a Production Planner and has been at NEB for over 4 years. In this role, she works with many departments within NEB to ensure timely product release. She is quite active in the NEB community, and is a member of the NEB Soccer Club, Running Club and Holiday Raffle Committee.





NEBNext Enzymatic Methyl-seq

NEW

NEBNext Enzymatic Methyl-seq Kit
#E7120S 24 reactions 895 €
#E7120L 96 reactions 3360 €

NEW

NEBNext Enzymatic Methyl-seq Conversion Module

#E7125S 24 reactions 185 € #E7125L 96 reactions 675 €

NEW

NEBNext Q5U Master Mix

#M0597S 50 reactions 125 € #M0597L 250 reactions 495 €

NEW


NEBNext Multiplex Oligos for Enzymatic
Methyl-seq (Unique Dual Index Primer Pairs)
#E7140S 24 reactions 135 €
#E7140L 96 reactions 535 €

- Superior sensitivity of detection of 5-mC and 5-hmC
- Larger library insert sizes
- More uniform GC coverage
- Greater mapping efficiency
- High-efficiency library preparation


NEBNext Enzymatic Methyl-seq is an enzymatic alternative to bisulfite conversion with superior performance. For more information, including extensive performance data, visit NEBNext.com.

Description: While bisulfite sequencing has been the gold standard for the study of DNA methylation, this conversion treatment is damaging to DNA, resulting in DNA fragmentation, loss and GC bias. The NEBNext Enzymatic Methyl-seq Kit (EM-seq™) provides an enzymatic alternative to whole genome bisulfite sequencing (WGBS), combined with highefficiency streamlined library preparation suitable for Illumina sequencing.

The highly effective EM-seq enzymatic conversion minimizes damage to DNA and, in combination with the supplied NEBNext Ultra II library preparation workflow reagents, results in high quality libraries that enable superior detection of 5-mC and 5-hmC from fewer sequencing reads.



NEBNext Enzymatic Methyl-seq libraries have larger insert sizes 50 ng Human NA12878 genomic DNA was sheared to 300 bp using the Covaris® S2 instrument and used as input into EM-seq and WGBS protocols. For WGBS, NEBNext Ultra II DNA was used for library construction, followed by the Zymo Research EZ DNA Methylation-Gold™ kit for bisulfite conversion. Libraries were sequenced on an Illumina MiSeq (2 x 76 bases) and insert sizes were determined using Picard 2.18.14. The normalized frequency of each insert size was plotted, illustrating that library insert sizes are larger for EM-seq than for WGBS, and indicating that EM-seq does not damage DNA as bisulfite treatment does in WGBS.




EM-seq identifies more CpGs than WGBS, at lower sequencing coverage depth with superior uniformity of GC coverage. 10, 50 and 200 ng Human NA12878 genomic DNA was sheared to 300 bp using the Covaris S2 instrument and used as input into EM-seq and WGBS protocols. For WGBS, NEBNext Ultra II DNA was used for library construction, followed by the Zymo Research EZ DNA Methylation-Gold Kit for bisulfite conversion. Libraries were sequenced on an Illumina NovaSeq® 6000 (2 x 100 bases). Reads were aligned to hg38 using bwa-meth 0.2.2.

A: Coverage of CpGs with EM-seq and WGBS libraries was analyzed using 324 million paired end reads, and each top and bottom strand CpGs were counted independently, yielding a maximum of 56 million possible CpG sites. EM-seq identifies more CpGs at lower depth of sequencing.

B: GC coverage was analyzed using Picard 2.17.2 and the distribution of normalized coverage across different GC contents of the genome (0-100%) was plotted. EM-seq libraries have significantly more uniform GC coverage, and lack the AT over-representation and GC under-representation typical of WGBS libraries.





The Importance and Fragility of the Amazon Biome

The expanse and biodiversity of the world's largest rainforest is hard to fathom. The Amazon covers almost 7 million km² (2.7 million miles²) in the northern region of South America. The Amazon River flows for 6,000 km (3,700 miles) from Peru, across Brazil, and into the Atlantic. The uniqueness of the biodiversity in the Amazon is matched nowhere else on Earth. It is home to one in 10 species on our planet, and 25% of all terrestrial species. There are 400 billion trees that belong to 16,000 different species. In Brazil alone, there are 100 different species of monkeys; nine new species have been discovered in the last decade. Dozens of plant and animal species are still being discovered every year.

Plant and animal species form intricate, cooperative relationships with each other that have evolved over millions of years. In the Amazon, adaption is the key to survival. For example, the Brazil nut tree cannot be cultivated outside of the Amazon because of its incredibly specialized, mutualistic relationship with large-bodied bees strong enough to open its petal and pollinate other flowers, the agouti who disperses the Brazil nut and sustains a healthy population of Brazil nut trees, and the scent of a specific forest-dwelling orchard that male bees use to attract females.

The Amazon is often referred to as the lungs of our planet — the trees help to regulate Earth's atmosphere by "breathing" in carbon dioxide and releasing oxygen, thereby keeping the Earth cooler. It is estimated that 40% of human-made carbon dioxide is cleared by rainforests.

The Amazon regulates local and global weather patterns by absorbing heat and releasing water vapor via photosynthesis into the atmosphere, seeding the clouds with rain. The Amazon River carries 20% of the world's water to the sea.

Unfortunately, the Amazon is under constant threat by humans who are logging and farming the land, mining its resources, damming its rivers, and destroying indigenous lands and cultures. Deforestation has the Amazon at the "tipping point", with just over 80% of trees remaining. This drastic reduction in the number of trees reduces the amount of moisture released into the atmosphere, thereby jeopardizing the rainfall patterns and the replenishment of the rivers. The act of cutting down the trees releases stored carbon into the atmosphere, contributing to global warming. The immense plant diversity of the rainforest is being replaced with monocrops, such as soy or palm, and the animal diversity replaced by a single species, such as grazing cattle.

Substantial efforts are being made to protect large swaths of forest. In Brazil, armed guards protect against illegal logging, and while this activity has decreased from 25,000 km²/year (15,500 miles²) to just over 6,000 km²/year (3,700 miles²), it is almost impossible to guard an area that is larger than India. Governments are expanding protected areas — for example, the Colombian government expanded protection around Chiribiquete National Park to incorporate an area that is home to three uncontacted and isolated tribes and more than 200,000 paintings of pre-Columbian art. Further, technology is being carefully introduced to contacted tribes to give them the control to map and manage over 70 million acres of ancestral rainforest.

Protecting the entire forest from human destruction is a seemingly impossible feat. However, global awareness of the significance of this rainforest, its biological, cultural and economic riches, can turn the tide and ensure its survival for generations to come.

Cellular Analysis



Novel tools to study expression & function of proteins.

Cell imaging analysis can use fluorescent dyes, fluorophore-labeled molecules or recombinant protein plasmid systems. Recombinant protein labeling systems and bioluminescent reporter systems are among the most sensitive fluorescence methods for imaging expression, transport, co-localization and degradation in either fixed or living cells. Protein labeling systems offer many advantages. For example, color changes can be easily implemented by using different substrates. Protein labeling systems can involve the use of tag-specific antibodies or antibodies to separate epitopes engineered into a plasmid tag system for detection. Protein labeling systems can be used with non-cell permeable substrates to enable the specific imaging of cell surface targets. This strategy is not possible with bioluminescent recombinant systems. In living cells, protein labeling substrates can be introduced and followed in cells over time. Two separate cellular targets can also be imaged simultaneously, using protein labeling systems with mutually exclusive, tag-specific fluorescent substrates.

Studies of protein expression, interactions and structure, often use reporter systems to introduce and select for gene targets in cells. Reporter genes confer drug resistance, bioluminescence or fluorescence properties in the cells into which they are introduced. Typical reporter studies link reporter genes directly to a promoter region of interest, the function of which can be monitored by the reporter activity. Protein fusion tagging is used to detect subcellular localization, degradation, protein-protein interactions, etc. Typical fusion tags are fluorescent proteins (e.g., eGFP) or small protein epitopes (e.g., FLAG, Myc HA) which can be detected by fluorescence FACS or western blots. New generations of reporter gene systems expand the range of applications and enhance experimental possibilities.

Featured Products

279 SNAP-Cell® Starter Kit

280 SNAP-Cell 647-SiR

280 SNAP-Surface® 649

Featured Tools & Resources

Troubleshooting Guide for SNAP-tag Technology

Application Notes

Videos and Tutorials: SNAP-tag Technology




Cellular Imaging & Analysis

Overview of SNAP- and CLIP-tag Technologies	278
Comparison of SNAP-tag/CLIP-tag	
Technology to GFP	279
SNAP-tag Starter Kit	279
Fluorescent Substrates	280
Blocking Agents	281
Anti-SNAP-tag Antibody, Polyclonal	281
SNAP-tag Purified Protein	281
Cloning Vectors	282
Biotin Labels	282
Vista Label	282
SNAP-Capture	283
Building Blocks	283



Cellular Imaging & Analysis



Live HeLa cell transfected with pSNAP_i-tubulin and pCLIP_i-Cox8A (mitochondrial cytochrome oxidase 8A). Cells were labeled with 3 µM SNAP-Cell TMR-Star (red) and 5 µM CLIP-Cell 505 (green) for 25 minutes and counterstained with Hoechst 33342 (blue) for nuclei.

Features of SNAP-tag and CLIP-tag:


- Clone and express once, then use with a variety of substrates
- Non-toxic to living cells
- Wide selection of fluorescent substrates
- · Highly specific covalent labeling
- · Simultaneous dual labeling

Applications of SNAP-tag and CLIP-tag:

- Simultaneous dual protein labeling inside or on the surface of live cells
- Protein localization and translocation
- Pulse-chase experiments
- Receptor internalization studies
- Selective cell surface labeling
- Protein pull down assays
- Protein detection in SDS-PAGE
- Flow cytometry
- High throughput binding assays in microtiter plates
- Biosensor interaction experiments
- FRET-based binding assays
- Single molecule labeling
- Super-resolution microscopy

New England Biolabs offers an innovative technology for studying the function and localization of proteins in live and fixed cells. Covalent protein labeling brings simplicity and versatility to the imaging of mammalian proteins in live cells, as well as the ability to capture proteins *in vitro*. The creation of a single genetic

construct generates a fusion protein which, when covalently attached to a variety of fluorophores, biotin or beads, provides a powerful tool for studying proteins. For added flexibility, NEB offers two systems in which the protein is labeled by a self-labeling fusion protein (SNAP-tag® and CLIP-tag™).



Protein labeling with SNAP-tag (gold) and CLIP-tag (purple). The SNAP- or CLIP-tag is fused to the protein of interest (blue). Labeling occurs through covalent attachment to the tag, releasing either a guanine or a cytosine moeity.

SNAP-tag and CLIP-tag – Self-Labeling Tag Technology

The SNAP- and CLIP-tag protein labeling systems enable the specific, covalent attachment of virtually any molecule to a protein of interest. There are two steps to using this system: cloning and expression of the protein of interest as a SNAP-tag fusion, and labeling of the fusion with the SNAP-tag substrate of choice. The SNAP-tag is based on the human O⁶-alkylguanine-DNA-alkyltransferase (hAGT), a DNA repair protein. SNAP-tag substrates are fluorophores, biotin or beads conjugated to guanine or chloropyrimidine leaving groups via a benzyl linker. In the labeling reaction, the substituted benzyl group of the substrate is covalently attached to the SNAP-tag. CLIP-tag is a modified version of SNAP-tag, engineered to react with benzylcytosine rather than benzylguanine derivatives. When used in conjunction with SNAP-tag, CLIP-tag enables the orthogonal and complementary labeling of two proteins simultaneously in the same cells.

SNAP-Cell®: SNAP-Cell labels are cell-permeant and uniquely suited for the labeling of SNAP-tag fusion proteins inside living or fixed cells, on cell surfaces or *in vitro*. These labels are spread across the visible spectrum, ranging from blue to red. Non-fluorescent cell-permeable blocking agent is also available.

SNAP-Surface®: SNAP-Surface labels are non-cell-permeant and routinely used to label SNAP-tag fusion proteins on the surface of living cells, in fixed cells or *in vitro*. These labels are spread across the visible spectrum and include the photostable AlexaFluor® dyes and a variety of other commonly used fluorophores. Non-fluorescent, non-cell-permeable blocking agent is also available.

CLIP-CeII™: CLIP-Cell labels are cell-permeant and uniquely suited for the labeling of CLIP-tag fusion proteins inside living or fixed cells, on cell surfaces or *in vitro*. The CLIP-tag is a derivative of the SNAP-tag that reacts with orthogonal substrates, allowing simultaneous labeling of two expressed proteins with different fluorophores. Non-fluorescent cell-permeable blocking agent is also available.

CLIP-Surface™: CLIP-Surface labels are non-cell permeant and routinely used to label CLIP-tag fusion proteins on the surface of living cells, in fixed cells or in vitro. The CLIP-tag is a derivative of the SNAP-tag that reacts with orthogonal substrates, allowing simultaneous labeling of two expressed proteins with different fluorophores. The labels include fluorophores at commonly used areas of the visible spectrum, such as 488, 547 and 647 nm.



Find an overview of SNAP-tag labeling.

Comparison of SNAP-tag/CLIP-tag Technologies to GFP

While SNAP/CLIP-tag technologies are complementary to GFP, there are several applications for which SNAP- and CLIP-tag self-labeling technologies are advantageous.

APPLICATION	SNAP-tag/CLIP-tag	GFP AND OTHER FLUORESCENT PROTEINS
Time-resolved fluorescence	Fluorescence can be initiated upon addition of label	Color is genetically encoded and always expressed. Photoactivatable fluorescent proteins require high intensity laser light, which may activate undesired cellular pathways (e.g., apoptosis)
Pulse-chase analysis	Labeling of newly synthesized proteins can be turned off using available blocking reagents (e.g., SNAP-Cell® Block)	Fluorescence of newly synthesized proteins cannot be specifically quenched to investigate dynamic processes
Ability to change colors	A single construct can be used with different fluorophore substrates to label with multiple colors	Requires separate cloning and expression for each color
Surface specific labeling	Can specifically label subpopulation of target protein expressed on cell surface using non-cell permeant substrates	Surface subpopulation cannot be specifically visualized
Single molecule detection	Conjugation with high quantum yield and photostable fluorophores	Fluorescent proteins are generally less bright and photobleach quicker than most organic fluorophores
Visualizing fixed cells	Resistant to fixation; strong labeling	Labile to fixation; weak labeling
Pull-down studies	"Bait" proteins can be covalently captured on BG beads	Requires anti-GFP antibody to non-covalently capture "bait" protein, complicating downstream analysis
Live animal imaging	Cell permeable far-red dye available, permitting deep tissue visualization	Signal is easily quenched by fixation (whole-mount specimens or thick sections); limited spectral flexibility and weaker fluorescence

SNAP-tag Starter Kit

SNAP-Cell® Starter Kit #E9100S

1 set 278 €

- Detailed yet easy to follow labeling protocol
- Robust, well characterized fluorophores included
- Control plasmid included; expressed proteins have well defined subcellular localization
- Examine protein localization either inside cells or on the surface of cells
- Intracellular labeling
- Cell surface labeling
- In vitro analysis

To enable researchers to quickly and easily begin using our protein labeling system, our Starter Kit includes all the components necessary to covalently attach either a red or a green fluorophore to SNAP-tag fusion proteins in living cells, fixed cells or in vitro. The SNAP-Cell Starter Kit contains a mammalian expression plasmid (pSNAP,) encoding the SNAP-tag flanked by restriction sites for cloning a gene of interest, and two cell-permeable fluorescent SNAP-tag substrates. A positive control plasmid (pSNAP,-Cox8A), encoding a SNAP-tagged protein (cytochrome c oxidase 8A) with a well-characterized mitochondrial localization, is also included. Lastly, a negative control "blocking agent" (SNAP-Cell Block) is included that interacts with the SNAP-tag, but is not fluorescent. There are two steps to using this system: subcloning and expression of the protein of interest as a SNAP, fusion, and labeling of the fusion with the SNAP-tag substrate of choice.

Fluorophores


Each of the fluorophores have been extensively validated and selected for their brightness and stability. Furthermore, they have been assessed for cell permeability.

SNAP-Cell TMR-Star is a photostable red fluorescent substrate that can be used to label SNAP-tag fusion proteins inside living cells or fixed cells, on cell surfaces, or *in vitro*. This cell-permeant substrate is based on tetramethylrhodamine and suitable for imaging with standard rhodamine filter sets. When covalently bound to SNAP-tag proteins, it has an excitation maximum at 554 nm and an emission maximum at 580 nm.

SNAP-Cell 505-Star is a photostable green fluorescent substrate that can be used to label SNAP-tag fusion proteins inside living cells or fixed cells, on cell surfaces, or *in vitro*. This cell-permeant substrate is suitable for imaging with standard fluorescein filter sets. When covalently bound to SNAP-tag proteins, it has an excitation maximum at 504 nm and an emission maximum at 532 nm.


The SNAP-Cell Starter Kit Includes:

- pSNAP, Vector
- pSNAP,-Cox8A Control Plasmid
- SNAP-Cell 505-Star
- SNAP-Cell TMR-Star
- SNAP-Cell Block



SNAP-Cell: Live CHO-K1 cells transiently transfected with pSNAP_T-H2B. Cells were labeled with SNAP-Cell 505-Star (green) for 15 minutes at 37°C, 5% CO₂,

Fluorescent Substrates for Protein Labeling



Live HeLa cell transfected with pSNAP_r-ER (endoplasmic reticulum) and pCLIP_r-tubulin. Cells were labeled with 3 µM SNAP-Cell TMR-Star (red) and 5 µM CLIP-Cell 505 (green) for 25 minutes and counterstained with Hoechst 33342 (blue) for nuclei.

- Fluorescently label SNAP-tag or CLIPtag fusions for cellular imaging
- Labels span fluorescent imaging spectrum from aqua (430 nm) to far-red (647+ nm) wavelengths
- Cell-permeable and non-cell-permeable labels available

NEB offers a large selection of fluorescent labels (substrates) for SNAP-tag and CLIP-tag fusion proteins. SNAP-tag substrates consist of a fluorophore conjugated to guanine or chloropyrimidine leaving groups via a benzyl linker, while CLIP-tag substrates consist of a fluorophore conjugated to a cytosine leaving group via a benzyl linker. These substrates

will label their respective tags without the need for additional enzymes. Cell-permeant substrates (SNAP- and CLIP-Cell) are suitable for both intracellular and cell-surface labeling, whereas non-cell-permeant substrates (SNAP- and CLIP-Surface) are specific for fusion proteins expressed on the cell surface only.

SELF-LABELING TAG							
	APPLICATIONS				SIZE	PRICE	
	Cell-Permeable						
	SNAP-Cell 430	S9109S	421	444,484	50 nmol	326 €	
	SNAP-Cell 505-Star	S9103S	504	532	50 nmol	326 €	
	SNAP-Cell Oregon Green®	S9104S	490	514	50 nmol	360 €	
	SNAP-Cell TMR-Star	S9105S	554	580	30 nmol	326 €	
	SNAP-Cell 647-SiR	S9102S	645	661	30 nmol	326 €	
SNAP-tag	Non-cell-permeable						
	SNAP-Surface Alexa Fluor® 488	S9129S	496	520	50 nmol	360 €	
	SNAP-Surface 488	S9124S	506	526	50 nmol	315 €	
	SNAP-Surface Alexa Fluor 546	S9132S	558	574	50 nmol	348 €	
	SNAP-Surface 549	S9112S	560	575	50 nmol	326 €	
	SNAP-Surface 594	S9134S	606	626	50 nmol	315 €	
	SNAP-Surface Alexa Fluor 647	S9136S	652	670	50 nmol	360 €	
	SNAP-Surface 649	S9159S	655	676	50 nmol	326 €	
	APPLICATIONS	NEB #	EXCITATION*	EMISSION*(1)	SIZE	PRICE	
	Cell-Permeable						
	CLIP-Cell 505	S9217S	504	532	50 nmol	326 €	
01.10.1	CLIP-Cell TMR-Star	S9219S	554	580	30 nmol	326 €	
CLIP-tag	Non-cell-permeable						
	CLIP-Surface 488	S9232S	506	526	50 nmol	326 €	
	CLIP-Surface 547	S9233S	554	568	50 nmol	326 €	
	CLIP-Surface 647	S9234S	660	673	50 nmol	326 €	

^{*} Excitation and emission values determined experimentally for labeled protein tag.

This table lists all currently available fluorescent substrates for SNAP-tag and CLIP-tag, along with excitation and emission wavelengths (determined from a labeled fusion tag, rather than the unreacted substrate).

⁽¹⁾ Colors are based on the electromagnetic spectrum. Actual color visualization may vary.

CELLULAR ANALYS

Blocking Agents

- Irreversible blocking
- Ideal for pulse-chase applications

Blocking agents are non-fluorescent substrates that block the reactivity of the SNAP- or CLIP-tag intracellularly (SNAP-Cell Block and CLIP-Cell Block) or on the surface of live cells (SNAP-Surface Block and CLIP-Cell Block). They can be used to generate inactive controls in live cell and *in vitro* labeling experiments performed with SNAP- or CLIP-tag fusion proteins.

SNAP- and CLIP-Cell Block are highly membrane permeant and once inside the cell react with the SNAP- or CLIP-tag, irreversibly inactivating them for subsequent labeling steps.

SNAP-Surface Block also reacts with the SNAP-tag irreversibly, inactivating it for subsequent labeling steps. This blocker is largely membrane impermeant essentially limiting blocking to cell surface-exposed SNAP-tags.

PRODUCT	NEB #	APPLICATION	SIZE	PRICE
SNAP-Cell Block	S9106S	Block SNAP-tag inside live cells and <i>in vitro</i>	100 nmol	128 €
CLIP-Cell Block	S9220S	Block CLIP-tag inside or on the surface of live cells and <i>in vitro</i>	100 nmol	128 €
SNAP-Surface Block	S9143S	Block SNAP-tag on the surface of live cells and in vitro	200 nmol	128€

Anti-SNAP-tag Antibody, Polyclonal

#P9310S

100 μl 256 €

Description: The Anti-SNAP-tag Antibody (Polyclonal) can be used in Western blots with SNAP-tag and CLIP-tag proteins. Polyclonal antibodies are produced from the immunization of rabbit with purified recombinant SNAP-tag protein and affinity purified using SNAP-BG resin.

Sensitivity: 5 ng of SNAP-tag per load in Western blotting.

Recommended Dilution: 1:1000

SNAP-tag Purified Protein

#P9312S

50 μg88 €

Description: SNAP-tag Purified Protein can be used as a positive control for *in vitro* labeling with various SNAP-tag fluorescent substrates. The coding sequence of SNAP-tag was cloned into a pTXB1 derived *E. coli* T7 expression vector. SNAP-tag protein was expressed and purified according to the instructions in the IMPACT™ kit manual (NEB #E6901). The purified SNAP-tag protein was dialyzed into


1X phosphate buffered saline (PBS) solution containing 1 mM DTT at 1 mg/ml (50 μ M) and stored at -80°C.

Molecular Weight: 19,694



Chris has been at NEB for over 12 years and is currently a member of the Production Team. Chris is a fitness enthusiast and enjoys running, obstacle course racing, endurance events—NEB is lucky to have him as a member of the Gym Committee.

Cloning Vectors



Live HeLa cell transfected with pSNAP_i-tubulin and pCLIPf-H2B constructs generated using pSNAP_i and pCLIP, vectors. Cells were labeled with 3 µM SNAP-Cell TMR-Star (red) and 5 µM CLIP-Cell 505 (green) for 25 minutes and counterstained with Hoechst 33342 (blue) for nuclei.

 Vectors for mammalian and bacterial expression available Vectors are available for SNAP-tag and CLIP-tag fusion protein expression and labeling in mammalian and bacterial systems as part of the starter kits. The mammalian SNAP, and CLIP, vectors express faster-reacting variants of the SNAP- and CLIP-tags than previously available vectors. Improved polylinker sequences both upstream and downstream from the tag allow expression of the tag on either end of the protein of interest, under control of the CMV promoter, SNAP,-tag and CLIP,-tag expression vectors contain a neomycin resistance (NeoR) gene for selection of stable transfectants, together with an IRES element for efficient expression of both the fusion protein and NeoR. Codon usage has been optimized for mammalian expression. Control plasmids encoding fusion proteins that are localized to the nucleus (H2B), mitochondria (Cox8A) and cell surface (ADRβ2, NK1R) are also available through Addgene.

The bacterial expression vector pSNAP-tag(T7)-2 includes cloning sites both upstream and downstream from the SNAP-tag, which is under control of the T7 promoter. Codon usage in the SNAP-tag gene has been optimized for *E. coli* expression.

Source: Isolated from an *E. coli* strain by a standard plasmid purification procedure. Plasmids have been purified free of endotoxins for efficient transfection

Concentration: 500 µg/ml

Restriction Map: For a more detailed description and restriction map of pSNAP, Vector, see page 389. See www.neb.com for sequence and map files for all expression and control plasmids.

PRODUCT	NEB #	FEATURES	SIZE	PRICE
pSNAP _f Vector	N9183S	stable and transient mammalian expression	20 μg	161 €
pSNAP-tag(T7)-2 Vector	N9181S	bacterial expression under T7 control	20 μg	161 €
pCLIP _f Vector	N9215S	stable and transient mammalian expression	20 µg	161 €

Biotin Labels


SNAP-Biotin® #S9110S 50 nmol 258 €

CLIP-Biotin

#S9221S 50 nmol 258 €

- Label SNAP-tag and CLIP-tag fusions with biotin
- Compatible with a variety of streptavidin conjugates
- Attach to streptavidin surfaces on microtiter plates

For optimal flexibility with existing technologies, biotinylated labels are available for studies using streptavidin platforms. Cell-permeant (SNAP-Biotin and CLIP-Biotin) substrates are based on biotin with an amidocaproyl linker. Biotin labels are suitable for applications such as biotinylation of fusion proteins for detection with streptavidin fluorophore conjugates or labeling in solution for analysis by SDS-PAGE/ Western blot. Biotin labels are also used for capture with streptavidin surfaces for binding and interaction studies.



Western blot analysis of biotin labeling reactions using anti-Biotin Antibody (CST #7075). SNAP-tag, CLIP-tag, and ACP-tag-MBP (5 μΜ) labeled with a biotin-containing substrate (10 μΜ). Marker M is Biotinylated protein ladder (CST #7727).


Vista Label

SNAP-Vista® Green

#S9147S 50 nmol 326 €

- Label protein fusions in cell lysates, or as purified proteins
- Detect proteins on SDS-PAGE
- Use with standard gel documentation equipment

SNAP-Vista Green fluorescent substrate can be used to label SNAP-tag fusions in cell lysates or as purified proteins for detection by SDS-PAGE. The substrate is optimal for visualization using a laser based gel scanner.



Typical SDS-PAGE of SNAP-Vista Green labeled proteins visualized using a gel scanner (Tyhoon 9400).


ELLULAR ANALYSI

SNAP-Capture

- Selectively capture SNAP-tag fusion proteins from solution
- Ideal for protein pull-down experiments or proteomic analysis

SNAP-Capture products are agarose or magnetic agarose beads coupled to a benzylguanine substrate, used to selectively capture and immobilize SNAP-tag fusion proteins from solution. These beads have a

high loading capacity for SNAP-tag fusion proteins and show very low non-specific adsorption of proteins from a complex lysate, making them especially suitable for pull-down applications.



Substrate structure on SNAP-Capture Pull-Down Resin

PRODUCT	NEB #	SIZE	PRICE
SNAP-Capture Pull-Down Resin	S9144S	2 ml	218 €
SNAP-Capture Magnetic Beads	S9145S	2 ml	210 €

Building Blocks

- Synthesize new SNAP-tag and CLIP-tag substrates
- Make surfaces for protein immobilization
- Attach novel molecules or ligands to proteins
- Create custom substrates for protein labeling

For advanced users with novel probes interested in working with SNAP-tag and CLIP-tag labeling technologies, building blocks are available for linkage of the core benzylguanine (BG) and benzylcytosine (BC) moieties to activated esters, primary amines and thiol groups. A variety of functional groups allows the choice of chemical coupling approaches to suit the molecule or surface to be coupled. Couple onto surfaces such as Biacore® chips or microarrays for specific protein immobilization. Couple onto peptides, proteins and DNA oligomers.

Couple onto new fluorophores or affinity reagents for specific protein labeling. Labeling is gentle, precise, and versatile: one label is covalently bound under biological conditions in a defined position.

References:

References for enzyme properties and applications for this product can be found at www.neb.com.

PRODUCT	NEB #	STRUCTURE	APPLICATION	SIZE	PRICE
BG-NH2	S9148S	N NH ₂	SNAP-tag substrate. Suitable for linkage to NHS esters and other activated carboxylic esters.	2 mg	384 €
BG-PEG-NH2	S9150S	N N NH2 H N O O O NH2	SNAP-tag substrate. PEG-linker gives superior flexibility. Particularly suited for immobilization on solid surfaces.	2 mg	384 €
BG-GLA-NHS	\$9151\$	N N NH ₂	SNAP-tag substrate. Activated as NHS ester. Reacts with primary amines.	2 mg	384 €
BG-Maleimide	S9153S	N N NH ₂	SNAP-tag substrate. Activated as maleimide. Reacts with thiols.	2 mg	384 €
BC-NH2	S9236S	NH ₂ NH ₂	CLIP-tag substrate. Suitable for linkage to NHS esters and other activated carboxylic esters.	2 mg	384 €





Backyard Biodiversity

Habitat loss due to human development is the primary cause of diminishing biodiversity. Living in urban areas furnished with exotic plant species and manicured lawns fragments ecosystems, and it leads to a disconnect between our everyday activities and the unintended consequences that they have on nature. As a result, biodiversity is often treated as a commodity. But, there are simple things that can be done to help preserve biodiversity — starting with learning how to share our living space with the organisms that are so essential to our existance.

Whether you live in a house with a yard or an apartment with a balcony, you can contribute to creating continuity in nature by bringing back the plants that were once naturally found there — this supports the safe travel of animals between core habitats as they search for cover and forage for food and water. A "wildlife corridor" that incorporates many vegetative layers and provides food, water and shelter will be more resilient to perturbations.

Native plants are adapted to local conditions and are easier to maintain, particularly in arid regions, so leave the native plant species undisturbed. Landscape using native trees and vegetation, and remove invasive plant species. Also, plants protect themselves by producing distasteful or toxic chemicals, and native insects that have evolved with specific native plant lineages develop a tolerance for and only eat, these plants. Non-native plants produce different chemicals, which can be detrimental to the native insects. Systemic pesticides such as neonicotinoids should always be avoided, as these pesticides persist in all parts of the plant and can poison the pollinators.

Fragmented habitats and pesticides leave pollinators malnourished. However, nectar from a variety of native flowering plants that bloom throughout the season attracts many pollinators — birds and beneficial insects — which keep the pests at bay without the use of pesticides. In addition to nectar-producing flowers, plants that feed butterfly larvae are also important.

Leaving wooded areas to age and decompose on their own offers significant benefit to many species. A dead tree may provide shelter or a perch for woodpeckers and other birds, frogs and lizards. A pile of rocks or logs can serve as a home for chipmunks or toads. Decomposing logs provide a habitat for insects and worms to thrive, and nutrients for the soil, encouraging plant growth. Insects and worms then pass the energy from plants to non-plant eating animals further up the food chain, such as spiders, birds and amphibians.

A tree hollow takes up to 150 years to naturally develop and is essential for nesting and breeding. Nesting boxes for birds, bees and bats can help alleviate the shortage of hollows. A birdbath, pond, or a carefully planned rain garden will attract birds and aquatic wildlife, such as frogs and dragonflies.

In addition to our own backyards, one can consider getting involved with regional ecological restoration efforts. Good places to look for opportunities are land trusts, wildlife foundations, native plant societies, government agencies (e.g., Forest Service, Fish & Wildlife), and environmental organizations. Protecting habitats before they have been damaged is the best form of biodiversity conservation and is most successfully implemented by government regulations.

As human development and urban expansion continues to increase, we are outcompeting other species for space on Earth. Educating ourselves on how we can co-exist and provide wildlife with water, food, cover and a safe place to raise their young in our own surroundings can go a long way in conserving biodiversity.

Reference Appendix



Technical Support – for scientists, by scientists

As a partner to the scientific community, New England Biolabs is committed to providing top quality tools and scientific expertise. This philosophy still stands, and has led to long-standing relationships with many of our fellow scientists. NEB's commitment to scientists is the same regardless of whether or not they purchase product from NEB: their ongoing research is supported by our catalog, website and technical staff.

NEB's technical support model is unique as it utilizes most of the scientists at NEB. Several of our product lines have designated technical support scientists assigned to servicing customers in those application areas. Any questions regarding a product could be dealt with by one of the technical support scientists, the product manager who manufactures it, the product development scientist who optimizes it, or a researcher who uses the product in their daily research. As such, customers are supported by scientists and often experts in the product or its application.


To access technical support:

- Call 1-800-632-7799 (Monday Friday: 9:00 am 6:00 pm EST)
- Submit an online form at www.neb.com/techsupport
- Email info@neb.com
- International customers can contact a local NEB subsidiary or distributor.
 For more information see inside back cover.

Featured Tools & Resources

- 290 Optimizing Restriction Enzyme Reactions
- Performance Chart for Restriction Enzymes
- 349 Troubleshooting Guide for Cloning
- **334** Methylation Sensitivity
- 337 Guidelines for PCR Optimization
- Cleavage Close to the End of DNA Fragments
- 288 Online Interactive Tools
- View NEB TV Episode #22 to learn more about our Technical Support program.
- Visit the Tools & Resources tab at www.neb.com to find additional online tools, video tech tips and tutorials to help you in your research.





Online Interactive Tools Databases, Mobile Apps & Freezer Program Locator		288 289	PCR, qPCR & RT-qPCR General Guidelines for PCR Optimization PCR Troubleshooting Guide		337 338
Restriction Endonucleases Optimizing Restriction Enzyme Reactions Double Digestion Types of Restriction Enzymes Troubleshooting Guide	2	290 291 291 292	Optimization Tips for Luna qPCR Optimization Tips for Luna One-Step RT-qPCR Luna qPCR Troubleshooting Guide Luna One-Step RT-qPCR Troubleshooting Guide Cleavage Close to the End of DNA Fragments Activity of Restriction Enzymes in PCR Buffers	;	339 340 341 342 343 345
Performance Performance Chart for Restriction Enzymes Activity of Enzymes at 37°C Activity of DNA Modifying Enzymes in CutSmart Buffer Tips for Avoiding Star Activity	2	298 299 299 300	Troubleshooting Guide for Cloning	347–3 349–3 352–3	351
High-Fidelity (HF) Restriction Enzymes Reduced Star Activities of HF Enzymes Time-Saver Qualified Restriction Enzymes		301 301 303	Nucleic Acid Purification Troubleshooting Guide for DNA Cleanup & Plasmid Purification Using Monarch Kits		354
Isoschizomers/Recognition Sequences Cross Index of Recognition Sequences Palindromic Penta-nucleotide Recognition Sequences Specificities Greater Than 6 Bases Interrupted Palindromes Multiple & Nonpalindromic Recognition Sequences	304–3	305 305 306	Choosing Sample Input Amounts When Using the Monarch Genomic DNA Purification Kit Troubleshooting Guide for Genomic DNA Purification Using Monarch Kits Choosing Sample Input Amounts When Using the Monarch Total RNA Miniprep Kit Troubleshooting Guide for Total RNA Extraction & Purification using Monarch Kits	e ;	355 356 357 358
Alphabetized List of NEB Recognition Sequences Isoschizomers	309–3 311–3		Transformation Genetic Markers Enhancing Transformation Efficiency	359–3	360 361
Properties of Restriction Enzymes Survival in a Reaction Cleavage of Supercoiled DNA	328–3 3	329 329	Protein Expression with T7 Express Strains NGS Library Preparation DNA/RNA Input Guidelines for NGS Library Prep	(362 363
Generating New Cleavage Sites Recleavable Filled-in 5' Overhangs Recleavable Blunt Ends Compatible Cohesive Ends Methylation Sensitivity		330 331 333	Protein Expression Labeling with SNAP-tag Technology Troubleshooting Guide Frequencies of Restriction Sites in	;	364
Dam, Dcm and CpG Methylation	334–3	336	DNA Maps	365–: 369–:	
			Odds and Ends General Molecular Biology Information	382–	385

One or more of these products are covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at gbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications. Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at https://www.neb.com/support/terms-of-sale. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

Use the Tools & Resources tab at www.neb.com to access our growing selection of interactive technical tools. These tools can also be accessed directly in the footer of every web page.

NEB scientists are often involved in the development of online tools that will aid in their research. We are now making these tools and in some circumstances, the source code, available for you to evaluate. To learn more, visit www.neb.com/NEBetaTools.

Online Tools

Competitor Cross-Reference Tool



Use this tool to select another company's product and find out which NEB product is compatible. Choose either the product name or catalog number from the available selections, and this tool will identify the recommended NEB product.

DNA Sequences and Maps Tool



With the DNA Sequences and Maps Tool, find the nucleotide sequence files for commonly used molecular biology tools, including plasmid, viral and bacteriophage vectors.

Double Digest Finder



Use this tool to guide your reaction buffer selection when setting up double digests, a common timesaving procedure. Choosing the right buffers will help you to avoid star activity and loss of product.

Enzyme Finder



Use this tool to select restriction enzymes by name, sequence, overhang or type. Enter your sequence using single letter code nomenclature, and Enzyme Finder will identify the right enzyme for the job.

Glycan Analyzer



Use this tool to interpret ultra or high pressure liquid chromatography (UPLC/HPLC) N-glycan profiles following exoglycosidase digestions.

NEB Golden Gate Assembly Tool



Use this tool to assist with in silico DNA construct design for Golden Gate DNA assembly. It enables the accurate design of primers with appropriate Type IIS restriction sites and overlaps, quick import of sequences in many formats and export of the final assembly, primers and settings.

NEBaseChanger®



NEBaseChanger can be used to design primers specific to the mutagenesis experiment you are performing using the Q5® Site-Directed Mutagenesis Kit. This tool will also calculate a recommended custom annealing temperature based on the sequence of the primers by taking into account any mismatches.

NEBNext® Selector



Use this tool to guide you through the selection of NEBNext reagents for next generation sequencing sample preparation.

NEBcutter® V2.1



Identify the restriction sites within your DNA sequence using NEBcutter. Choose between Type II and commercially available Type III restriction enzymes to digest your DNA. NEBcutter will indicates cut frequency and methylation-state sensitivity.

NEBioCalculator®



Use this tool for your scientific calculations and conversions for DNA and RNA. Options include conversion of mass to moles, ligation amounts, conversion of OD to concentration, dilution and molarity. Additional features include sgRNA template oligo design and qPCR library quantification.

NEBcloner®



Use this tool to find the right products and protocols for each step (digestion, end modification, ligation, transformation and mutagenesis) of your next traditional cloning experiment. Also, find other relevant tools and resources to enable protocol optimization.

NEBuilder® Assembly Tool



Use this tool to design primers for your DNA assembly reaction, based on the entered fragment sequences and the polymerase being used for amplification.

PCR Fidelity Estimator



Estimate the percentage of correct DNA copies (those without base substitution errors) per cycle of PCR for selected DNA polymerases.

PCR Selector



Use this tool to help select the right DNA polymerase for your PCR setup. Whether your amplicon is long, complex, GC-rich or present in a single copy, the PCR selection tool will identify the perfect DNA polymerase for your reaction.

Tm Calculator



Determine the optimal annealing temperature for your amplicon with our Tm Calculator. Simply input your DNA polymerase, primer concentration and your primer sequence, and the Tm Calculator will guide you to successful reaction conditions.

Thermostable Ligase Reaction Temperature Calculator



This tool will help you estimate an optimal reaction temperature to minimize mismatch for thermostable ligation of two adjacent ssDNA probes annealed to a template.

Online Tools (continued)

Read Coverage Calculator



This tool allows for easy calculation of values associated with read coverage in NGS protocols.

Additional Databases

Polbase®



Polbase is a repository of biochemical, genetic and structural information about DNA Polymerases.

REBASE®



Use this tool as a guide to the ever-changing landscape of restriction enzymes. REBASE, the Restriction Enzyme DataBASE, is a dynamic, curated database of restriction enzymes and related proteins.

Mobile Apps

NEB Tools for iPhone®, iPad® or Android®



NEB Tools brings New England Biolabs' most popular web tools to your iPhone, iPad or Android devices.

- Use Enzyme Finder to select a restriction enzyme by category or recognition sequence, or search by name to find information on any NEB
 enzyme. Sort your results so they make sense to you, then email them to your inbox or connect directly to www.neb.com.
- Use Double Digest Finder to determine buffer and reaction conditions for experiments requiring two restriction enzymes.
- Use Tm Calculator to calculate annealing temperatures for your PCR reaction.
- · Also included are several popular calculators from the NEBioCalculator web app.

When using either of these tools, look for CutSmart*, HF* and Time-Saver* enzymes for the ultimate in convenience. NEB Tools enables quick and easy access to the most requested restriction enzyme information, and allows you to plan your experiments from anywhere.


IPHONE® and IPAD® are registered trademarks of Apple Computers, Inc. ANDROID® is a registered trademark of Google. Inc.


Looking for a Freezer Program?

NEBnow® Locator



NEBnow Freezer Programs are ideally suited for researchers in academics and industry looking for on-site access to the world's finest restriction enzymes and related products. NEB freezers offer you convenience, flexibility and value.





Optimizing Restriction Enzyme Reactions

While standard recommended reaction conditions are a good place to start, in some cases, optimization may be necessary to achieve the best results. Depending on the enzyme(s) being used, variables such as incubation time, number of enzyme units used, and reaction temperature should be tested to find the optimal reaction conditions for your substrate DNA and enzyme(s) of choice.

Protocol: Restriction Enzyme Reactions

	STANDARD Protocol	TIME-SAVER Protocol
DNA	up to 1 μg	up to 1 μg
10X NEBuffer	5 μl (1X)	5 μl (1X)
Restriction Enzymes	10 units*	1 μΙ
Total Volume	50 μΙ	50 μΙ
Incubation Temperature	Enzyme-dependent	Enzyme-dependent
Incubation Time	60 minutes	5–15 minutes**

TOOLS & RESOURCES

Visit NEBRestrictionEnzymes.com to find:

- · Online tutorials for setting up restriction enzyme digests
- · Tips to avoid star activity
- Restriction Enzyme Performance Chart
- · Troubleshooting guide
- · Access to NEB's online tools, including: **Enzyme Finder, Double Digest Finder** and **NEBcloner**



TIPS FOR OPTIMIZATION

Enzyme

- · Keep on ice when not in the freezer
- Should be the last component added to reaction
- Mix components by pipetting the reaction mixture up and down, or by "flicking" the reaction tube. Follow with a quick ("touch") spin-down in a microcentrifuge. Do not vortex the reaction.
- In general, we recommend 5-10 units of enzyme per µg DNA, and 10-20 units per µg of genomic DNA in a 1 hour digest

Star Activity

- · Unwanted cleavage that can occur when enzyme is used under sub-optimal conditions, such as:
- Too much enzyme present
- Too long of an incubation time
- Using a non-recommended buffer
- Glycerol concentrations above 5%
- · Star activity can be reduced by using a High-Fidelity (HF®) enzyme, reducing incubation time, using a Time-Saver™ enzyme or increasing reaction volume

DNA

- · Should be free of contaminants such as phenol, chloroform, alcohol, EDTA, detergents and salts. Spin column purification readily accomplishes this; extra washes during purification can also
- Methylation of DNA can affect digestion with certain enzymes. For more information about methylation visit www.neb.com/methylation

Buffer

- Use at a 1X concentration
- BSA is included in NEBuffer 1.1, 2.1, 3.1 and CutSmart® Buffer. No additional BSA is needed.
- Restriction enzymes that do not require BSA for optimal activity are not adversely affected if BSA is present in the reaction

Reaction Volume

- A 50 µl reaction volume is recommended for digestion of up to 1 µg of substrate. This helps maintain salt levels introduced by miniprepped DNA low enough that they don't affect enzyme activity.
- Enzyme volume should not exceed 10% of the total reaction volume to prevent star activity due to excess glycerol
- · Additives in the restriction enzyme storage buffer (e.g., glycerol, salt), as well as contaminants found in the substrate solution (e.g., salt, EDTA or alcohol), can be problematic in smaller reaction volumes

	RESTRICTION ENZYME*	DNA	10X Nebuffer
10 μl rxn**	1 unit	0.1 μg	1 µl
25 µl rxn	5 units	0.5 μg	2.5 µl
50 µl rxn	10 units	1 µg	5 μΙ

- * Restriction enzymes can be diluted using the recommended
- ** 10 µl rxns should not be incubated for longer than 1 hour to avoid evaporation

Incubation Time

- Incubation time for the Standard Protocol is 1 hour. Incubation for the Time-Saver Protocol is 5-15 minutes.
- · Visit www.neb.com/timesaver for list of Time-Saver qualified enzymes
- · It is possible, with many enzymes, to use fewer units and digest for up to 16 hours. For more information, visit www.neb.com

Storage

- Storage at -20°C is recommended for most restriction enzymes. For a few enzymes, storage at -80°C is recommended. Visit www.neb.com for storage information.
- 10X NEBuffers should be stored at -20°C

Stability

- · The expiration date is found on the label
- Long term exposure to temperatures above –20°C should be minimized whenever possible

^{*}Sufficient to digest all types of DNAs.

**Time-Saver qualified enzymes can also be incubated overnight with no star activity.

Double Digestion

Digesting a DNA substrate with two restriction enzymes simultaneously (double digestion) is a common timesaving procedure. Over 210 restriction enzymes are 100% active in CutSmart Buffer, making double digestion simple. If you are using an enzyme that is not supplied with CutSmart Buffer, the Performance Chart for Restriction Enzymes (pages 293–298) rates the percentage activity of each restriction endonuclease in the four standard NEBuffers.

Setting up a Double Digestion

- Double digests with CutSmart restriction enzymes can be set up in CutSmart Buffer. Otherwise, choose an NEBuffer that results in the most activity for both enzymes. If star activity is a concern, consider using one of our High-Fidelity (HF) enzymes.
- Set up reaction according to recommended protocol (see page 290). The final concentration of glycerol in any reaction should be less than 5% to minimize the possibility of star activity (see page 300). For example, in a 50 µl reaction, the total amount of enzyme added should not exceed 5 µl.
- If two different incubation temperatures are necessary, choose the optimal reaction buffer and set up reaction accordingly. Add the first enzyme and incubate at the desired temperature. Then, heat inactivate the first enzyme, if it can be heat inactivated, add the second enzyme and incubate at the recommended temperature.
- Depending on an enzyme's activity rating in a non-optimal NEBuffer, the number of units or incubation time may be adjusted to compensate for the slower rate of cleavage.

Setting up a Double Digestion with a Unique Buffer (designated "U")

 NEB currently supplies three enzymes with unique buffers: EcoRI, Sspl and DpnII. In most cases, DpnII requires a sequential digest. Note that EcoRI and Sspl have HF versions (NEB #R3101 and NEB #R3132, respectively) which is supplied with CutSmart Buffer.

Setting up a Sequential Digestion

- If there is no buffer in which the two enzymes exhibit > 50% activity, a sequential digest can be performed.
- Set up a reaction using the restriction endonuclease that has the lowest salt concentration in its recommended buffer and incubate to completion.
- Adjust the salt concentration of the reaction (using a small volume of a concentrated salt solution) to approximate the reaction conditions of the second restriction endonuclease.
- Add the second enzyme and incubate to complete the second reaction.
- Alternatively, a spin column can be used to isolate the DNA prior to the second reaction.

TOOLS & RESOURCES

Visit www.neb.com/nebtools for:

 Help choosing double digest conditions using NEB's, Double Digest Finder and NEBcloner®



Types of Restriction Enzymes

Restriction enzymes are traditionally classified into four types on the basis of subunit composition, cleavage position, sequence specificity and cofactor requirements. However, amino acid sequencing has uncovered extraordinary variety among restriction enzymes and revealed that at the molecular level there are many more than four different types.

Type I Enzymes are complex, multisubunit, combination restriction-and-modification enzymes that cut DNA at random far from their recognition sequences. Type I enzymes are of considerable biochemical interest, but they have little practical value since they do not produce discrete restriction fragments or distinct gel-banding patterns.

Type II Enzymes cut DNA at defined positions close to or within their recognition sequences. They produce discrete restriction fragments and distinct gel banding patterns, and they are the only class used in the laboratory for DNA analysis and gene cloning. Rather than forming a single family of related proteins, Type II enzymes are a collection of unrelated proteins of many different sorts. Type II enzymes frequently differ so utterly in amino acid sequence from one another, and indeed from every other known protein, that they exemplify the class of rapidly evolving proteins that are often indicative of involvement in host-parasite interactions.

Type III Enzymes are also large combination restriction-and-modification enzymes. They cleave outside of their recognition sequences and require two such sequences in opposite orientations within the same DNA molecule to accomplish cleavage; they rarely yield complete digests.

Type IV Enzymes recognize modified, typically methylated DNA and are exemplified by the McrBC and Mrr systems of *E. coli*.

TOOLS & RESOURCES

Visit the video library at www.neb.com to find:

Tutorials on Type I, II and III restriction enzymes



TYPE I, II AND III RESTRICTION ENZYMES



Restriction Enzyme Troubleshooting Guide

PROBLEM	CAUSE	SOLUTION
		Check the methylation sensitivity of the enzyme(s) to determine if the enzyme is blocked by methylation of the recognition sequence
	Restriction enzyme(s)	* Use the recommended buffer supplied with the restriction enzyme
Few or no transformants	didn't cleave completely	* Clean up the DNA to remove any contaminants that may inhibit the enzyme
		When digesting a PCR fragment, make sure to have at least 6 nucleotides between the recognition site and the end of the DNA molecule
	The restriction enzyme(s) is bound	Lower the number of units
The digested DNA ran as a	to the substrate DNA	* Add SDS (0.1–0.5%) to the loading buffer to dissociate the enzyme from the DNA
smear on an agarose gel	Nucleon contemination	Use fresh, clean running buffer and a fresh agarose gel
	Nuclease contamination	Clean up the DNA
		DNA isolated from a bacterial source may be blocked by Dam and Dcm methylation
	Classage in blooked	* DNA isolated from eukaryotic source may be blocked by CpG methylation
	Cleavage is blocked by methylation	Check the methylation sensitivity of the enzyme(s) to determine if the enzyme is blocked by methylation of the recognition sequence
		* If the enzyme is inhibited by Dam or Dcm methylation, grow the plasmid in a dam-/dcm- strain (NEB #C2925)
	Coll in his little	Enzymes that have low activity in salt-containing buffers (NEBuffer 3.1) may be salt sensitive, so clean up the DNA prior to digestion
	Salt inhibition	• DNA purification procedures that use spin columns can result in high salt levels, which inhibit enzyme activity. To prevent this, DNA solution should be no more than 25% of total reaction volume.
Incomplete restriction	Inhibition by PCR components	Clean up the PCR fragment prior to restriction digest
enzyme digestion	Using the wrong buffer	* Use the recommended buffer supplied with the restriction enzyme
	Too few units of enzyme used	* Use at least 3–5 units of enzyme per μg of DNA
	Incubation time was too short	Increase the incubation time
	Digesting supercoiled DNA	* Some enzymes have a lower activity on supercolled DNA. Increase the number of enzyme units in the reaction.
	Presence of slow sites	* Some enzymes can exhibit slower cleavage towards specific sites. Increase the incubation time, 1–2 hours is typically sufficient.
	Two sites required	* Some enzymes require the presence of two recognition sites to cut efficiently
	DNA is contaminated with an inhibitor	Assay substrate DNA in the presence of a control DNA. Control DNA will not cleave if there is an inhibitor present. Miniprep DNA is particularly susceptible to contaminants.
	an minului	Clean DNA with a spin column, resin or drop dialysis, or increase volume to dilute contaminant
	If larger bands than expected are seen	Lower the number of units in the reaction
	in the gel, this may indicate binding of the enzyme(s) to the substrate	* Add SDS (0.1–0.5%) to the loading buffer to dissociate the enzyme from the substrate
		Use the recommended buffer supplied with the restriction enzyme
		Decrease the number of enzyme units in the reaction
	Star activity	• Make sure the amount of enzyme added does not exceed 10% of the total reaction volume. This ensures that the total glycerol concentration does not exceed 5% v/v.
	•	• Decrease the incubation time. Using the minimum reaction time required for complete digestion will help prevent star activity.
Extra bands in the gel		• Try using a High-Fidelity (HF) restriction enzyme. HF enzymes have been engineered for reduced star activity.
		• Enzymes that have low activity in salt-containing buffers (e.g., NEBuffer 3.1) may be salt sensitive. Make sure to clean up the DNA prior to digestion.
	Partial restriction, enzyme dinest	• DNA purification procedures that use spin columns can result in high salt levels, which inhibit enzyme activity. To prevent this, DNA solution should be no more than 25% of total reaction volume.
	Partial restriction enzyme digest	* Clean-up the PCR fragment prior to restriction digest
		* Use the recommended buffer supplied with the restriction enzyme
		* Use at least 3–5 units of enzyme per µg of DNA and digest the DNA for 1–2 hours

New England Biolabs supplies > 210 restriction enzymes that are 100% active in a single buffer, CutSmart. This results in increased efficiency, flexibility and ease-of-use, especially when performing double digests.

This performance chart summarizes the activity information of NEB restriction enzymes. To help select the best conditions for double digests, this chart shows the optimal (supplied) NEBuffer and approximate activity in the four standard NEBuffers for each enzyme. Note that BSA is included in all NEBuffers, and is not provided as a separate tube. In addition, this performance chart shows recommended reaction temperature, heat-inactivation temperature, recommended diluent buffer, methylation sensitivity and whether the enzyme is Time-Saver qualified (e.g., cleaves substrate in 5–15 minutes under recommended conditions, and can be used overnight without degradation of DNA).

Chart Legend

- U Supplied with a unique reaction buffer that is different from the four standard NEBuffers. The compatibility with the four standard NEBuffers is indicated in the chart.
- RR Recombinant
- *e* Engineered enzyme for maximum performance
- Time-Saver qualified
- Indicates that the restriction enzyme requires two or more sites for cleavage
- SAM Supplied with a separate vial of S-adenosylmethionine (SAM). To obtain 100% activity, SAM should be added to the 1X reaction mix as specified on the product data card.
- dcm dcm methylation sensitivity
- dam methylation sensitivity
- CpG methylation sensitivity

Activity Notes (see last column)

FOR STAR ACTIVITY

- Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.
- 2. Star activity may result from extended digestion.
- 3. Star activity may result from a glycerol concentration of > 5%.
- * May exhibit star activity in this buffer.

FOR LIGATION AND RECUTTING

- a. Ligation is less than 10%
- b. Ligation is 25% 75%
- c. Recutting after ligation is < 5%
- d. Recutting after ligation is 50% 75%
- e. Ligation and recutting after ligation is not applicable since the enzyme is either a nicking enzyme, is affected by methylation, or the recognition sequence contains variable sequences.

TIE				o/ 10-				INCUB.	INACTIV.					
12		ENZYME	SUPPLIED Nebuffer	% ACT 1.1	2.1	EBUFFERS 3.1	CUTSMART	TEMP. (°C)	TEMP. (°C)	DIL.	SUBSTRATE	METHYLAT SENSITIV		NOTE(S)
RR	0	AatII	CutSmart	< 10	50*	50	100	37°	80°	В	Lambda		CpG	
RX		AbaSI	CutSmart	25	50	50	100	25°	65°	С	T4 wt Phage			е
RR	•	Accl	CutSmart	50	50	10	100	37°	80°	Α	Lambda		CpG	
RX	0	Acc65I	3.1	10	75*	100	25	37°	65°	Α	pBC4	dcm	CpG	
RR	•	Acil	CutSmart	< 10	25	100	100	37°	65°	Α	Lambda		CpG	d
RX	0	AcII	CutSmart	< 10	< 10	< 10	100	37°	No	В	Lambda		CpG	
RR	•	Acul	CutSmart + SAM	50	100	50	100	37°	65°	В	Lambda			1, b, d
RX		Afel	CutSmart	25	100	25	100	37°	65°	В	pXba		CpG	
RR	•	AfIII	CutSmart	50	100	10	100	37°	65°	Α	phiX174			
RX		AfIIII	3.1	10	50	100	50	37°	80°	В	Lambda			
RX		Agel	1.1	100	75	25	75	37°	65°	С	Lambda		CpG	
RX	9 <i>e</i>	Agel-HF	CutSmart	100	50	10	100	37°	65°	Α	Lambda		CpG	
RR	•	Ahdl	CutSmart	25	25	10	100	37°	65°	Α	Lambda		CpG	a
RX		Alel-v2	CutSmart	< 10	< 10	< 10	100	37°	80°	В	Lambda		CpG	
RX	0	Alul	CutSmart	25	100	50	100	37°	80°	В	Lambda			b
RX		Alwl	CutSmart	50	50	10	100	37°	No	Α	Lambda dam-	dam		1, b, d
RR	•	AlwNI	CutSmart	10	100	50	100	37°	80°	Α	Lambda	dcm		
RX	0	Apal	CutSmart	25	25	< 10	100	25°	65°	Α	pXba	dcm	CpG	
RR	•	ApaLI	CutSmart	100	100	10	100	37°	No	Α	Lambda HindIII		CpG	
RX	0	ApeKI	3.1	25	50	100	10	75°	No	В	Lambda		CpG	
RR	•	Apol	3.1	10	75	100	75	50°	80°	Α	Lambda			
RX	9 <i>e</i>	Apol-HF	CutSmart	10	100	10	100	37°	80°	В	Lambda			
RR	•	Ascl	CutSmart	< 10	10	10	100	37°	80°	Α	Lambda		CpG	
RX	0	Asel	3.1	< 10	50*	100	10	37°	65°	В	Lambda			3
RX		AsiSI	CutSmart	100	100	25	100	37°	80°	В	pXba (Xho digested)		CpG	2, b
RX	0	Aval	CutSmart	< 10	100	25	100	37°	80°	Α	Lambda		CpG	
RR	9	Avall	CutSmart	50	75	10	100	37°	80°	Α	Lambda	dcm	CpG	
RX	0	AvrII	CutSmart	100	50	50	100	37°	No	В	Lambda HindIII			
RX	•	Bael	CutSmart + SAM	50	100	50	100	25°	65°	Α	Lambda		CpG	е
RR	•	BaeGI	3.1	75	75	100	25	37°	80°	Α	Lambda			
RX	0	BamHI	3.1	75*	100*	100	100*	37°	No	Α	Lambda			3
RR	9 <i>e</i>	BamHI-HF	CutSmart	100	50	10	100	37°	No	Α	Lambda			
RX		Banl	CutSmart	10	25	< 10	100	37°	65°	Α	Lambda	dcm	CpG	1

Performance Chart for Restriction Enzymes (continued)

1		V			SUPPLIED	% AC	TIVITY IN N	FRUFFFRS		INCUB. TEMP.	INACTIV. TEMP.			METHYL	TION	
				ENZYME	NEBUFFER	1.1	2.1	3.1	CUTSMART	(°C)	(°C)	DIL.	SUBSTRATE	SENSITI		NOTE(S)
RX				Banll	CutSmart	100	100	50	100	37°	80°	Α	Lambda			2
R₩	0			Bbsl	2.1	100	100	25	75	37°	65°	В	Lambda			
R₩		e		BbsI-HF	CutSmart	10	10	10	100	37°	65°	В	Lambda			
R\lambda	•	2	?*site	Bbvl	CutSmart	100	100	25	100	37°	65°	В	pBR322			3
RX				BbvCl	CutSmart	10	100	50	100	37°	No	В	Lambda		CpG	1, a
RX				Bccl	CutSmart	100	50	10	100	37°	65°	Α	pXba			3, b
RX				BceAl	3.1	100*	100*	100	100*	37°	65°	Α	pBR322		CpG	1
RX		2	?*site	Bcgl	3.1 + SAM	10	75*	100	50*	37°	65°	А	Lambda	dam	CpG	е
RX	•			BciVI	CutSmart	100	25	< 10	100	37°	80°	С	Lambda			b
RX	•			Bcll	3.1	50	100	100	75	50°	No	А	Lambda dam-	dam		
RX		e		BcII-HF	CutSmart	100	100	10	100	37°	65°	В	Lambda dam-	dam	_	
RX	•			BcoDI	CutSmart	50	75	75	100	37°	No	В	Lambda		CpG	
R₩				Bfal	CutSmart	< 10	10	< 10	100	37°	80°	В	Lambda		_	2, b
RX	9	2	?+site	BfuAl	3.1	< 10	25	100	10	50°	65°	В	Lambda		CpG	3
RX	0			BgII	3.1	10	25	100	10	37°	65°	В	Lambda		CpG	
RX	•			BgIII	3.1	10	10	100	< 10	37°	No	A	Lambda			
RX	0			Blpl	CutSmart	50	100	10	100	37°	No oso	A	Lambda		00	d
RX	0			BmgBl	3.1	< 10	10	100	10	37°	65°	В	Lambda		CpG	3, b, d
RX RX				Bmrl	2.1	75	100	75	100*	37°	65°	В	Lambda HindIII			b
RX	•	e		Bmtl UE	3.1	100	100	100	100	37°	65° 65°	B B	pXba			2
RX			?+site	Bmtl-HF Bpml	CutSmart 3.1	50 75	100 100	100	100 100	37° 37°	65°	В	pXba			2
RX		<u> </u>	2 3110	Bpu10l	3.1	10	25	100	25	37°	80°	В	Lambda Lambda			3, b, d
RX	•			BpuEl	CutSmart + SAM	50*	100	50*	100	37°	65°	В	Lambda			d d
RX				Bsal	CutSmart	75*	75	100	100	37°	65°	В	pXba	dcm	CpG	3
RX	•	e		Bsal-HFv2	CutSmart	100	100	100	100	37°	80°	В	pXba	dcm	CpG	3
RR	9			BsaAl	CutSmart	100	100	100	100	37°	No	С	Lambda		CpG	
RR				BsaBl	CutSmart	50	100	75	100	60°	80°	В	Lambda dam-	dam	CpG	2
RR	0			BsaHl	CutSmart	50	100	100	100	37°	80°	С	Lambda	dcm	CpG	-
RX				BsaJI	CutSmart	50	100	100	100	60°	80°	A	Lambda			
RX	0			BsaWI	CutSmart	10	100	50	100	60°	80°	A	Lambda			
	0			BsaXI	CutSmart	50*	100*	10	100	37°	No	С	Lambda			е
RX	•			BseRI	CutSmart	100	100	75	100	37°	80°	A	Lambda			d
R\lambda				BseYI	3.1	10	50	100	50	37°	80°	В	Lambda		CpG	d
R₩	0	2	?+site	Bsgl	CutSmart + SAM	25	50	25	100	37°	65°	В	Lambda			d
R\\	0			BsiEl	CutSmart	25	50	< 10	100	60°	No	Α	Lambda		CpG	
RX				BsiHKAI	CutSmart	25	100	100	100	65°	No	Α	Lambda			
RX	0			BsiWl	3.1	25	50*	100	25	55°	65°	В	phiX174		CpG	
RX	0	e		BsiWI-HF	CutSmart	50	100	10	100	37°	No	В	phiX174		CpG	
RX	0			BsII	CutSmart	50	75	100	100	55°	No	Α	Lambda	dcm	CpG	b
RX	0			Bsml	CutSmart	25	100	< 10	100	65°	80°	Α	Lambda			
R₩	•			BsmAl	CutSmart	50	100	100	100	55°	No	В	Lambda		CpG	
R\lambda	•			BsmBI	3.1	10	50*	100	25	55°	80°	В	Lambda		CpG	
RX				BsmFl	CutSmart	25	50	50	100	65°	80°	Α	pBR322	dcm	CpG	1
RX	•			BsoBl	CutSmart	25	100	100	100	37°	80°	Α	Lambda			
RX	•			Bsp1286I	CutSmart	25	25	25	100	37°	65°	А	Lambda			3
	•			BspCNI	CutSmart + SAM	100	75	10	100	25°	80°	А	Lambda	_	_	b
RX				BspDI	CutSmart	25	75	50	100	37°	80°	А	Lambda	dam	CpG	
RR	9			BspEl	3.1	< 10	10	100	< 10	37°	80°	В	Lambda dam-	dam	CpG	
RX	0			BspHI	CutSmart	< 10	50	25	100	37°	80°	A	Lambda	dam		
RX			?*site	BspMI	3.1	10	50*	100	10	37°	65°	В	Lambda			

a. Ligation is less than 10% b. Ligation is 25% – 75%

c. Recutting after ligation is <5% d. Recutting after ligation is 50%-75%

e. Ligation and recutting after ligation is not applicable since the enzyme is either a nicking enzyme, is affected by methylation, or the recognition sequence contains variable sequences.

				ENZYME	SUPPLIED Nebuffer	% ACTI 1.1	VITY IN N 2.1	EBUFFERS 3.1	CUTSMART	INCUB. TEMP. (°C)	INACTIV. TEMP. (°C)	DIL.	SUBSTRATE	METHYLATION SENSITIVITY	
RX	•			BspQl	3.1	100*	100*	100	100*	50°	80°	В	Lambda		3
	0			Bsrl	3.1	< 10	50	100	10	65°	80°	В	phiX174		b
RX	0			BsrBl	CutSmart	50	100	100	100	37°	80°	A	Lambda	CpG	d
RX	0			BsrDI	2.1	10	100	75	25	65°	80°	Α	Lambda		3, d
RX	0	e		BsrFI-v2	CutSmart	25	25	0	100	37°	No	С	pBR322	CpG	
RX	•			BsrGl	2.1	25	100	100	25	37°	80°	A	Lambda		
RX	0	e		BsrGI-HF	CutSmart	10	100	100	100	37°	80°	Α	Lambda		
RX	•			BssHII	CutSmart	100	100	100	100	50°	65°	В	Lambda	CpG	
RX	0	e		BssSI-v2	CutSmart	10	25	< 10	100	37°	No	В	Lambda		
R₩				BstAPI	CutSmart	50	100	25	100	60°	80°	А	Lambda	CpG	b
RX	•			BstBl	CutSmart	75	100	10	100	65°	No	Α	Lambda	CpG	
RX	•			BstEII	3.1	10	75*	100	75*	60°	No	Α	Lambda		3
RX	0	e		BstEII-HF	CutSmart	< 10	10	< 10	100	37°	No	Α	Lambda		
RX	•			BstNI	3.1	10	100	100	75	60°	No	Α	Lambda		a
	•			BstUI	CutSmart	50	100	25	100	60°	No	Α	Lambda	CpG	b
RX	0			BstXI	3.1	< 10	50	100	25	37°	80°	В	Lambda	dem	3
R₩	•			BstYI	2.1	25	100	75	100	60°	No	Α	Lambda		
RX	0	e		BstZ17I-HF	CutSmart	100	100	10	100	37°	No	Α	Lambda	CpG	
RX	0			Bsu36l	CutSmart	25	100	100	100	37°	80°	С	Lambda HindIII		b
RX	0			Btgl	CutSmart	50	100	100	100	37°	80°	В	pBR322		
RX		_		BtgZl	CutSmart	10	25	< 10	100	60°	80°	Α	Lambda	CpG	3, b, d
RX	0	e		BtsI-v2	CutSmart	100	100	25	100	55°	No	Α	Lambda		
RX		e		BtsIMutI	CutSmart	100	50	10	100	55°	80°	A	pUC19		b
RX	0			BtsCI	CutSmart	10	100	25	100	50°	80°	В	Lambda	_	
- DW	0			Cac8I	CutSmart	50	75	100	100	37°	65°	В	Lambda	CpG	b
RX	0		Otoles	Clal	CutSmart	10	50	50	100	37°	65°	A	Lambda dam-	dam CpG	
RX	0		2*site	CspCl	CutSmart + SAM	10	100	10	100	37°	65°	A	Lambda		е
RX RX	V			CviAll	CutSmart	50	50	10	100	25° 37°	65°	C	Lambda		1, b
RX	•			CviKI-1 CviQI	CutSmart 3.1	25 75	100 100*	100	100 75*	25°	No No	A C	pBR322 Lambda		1, υ b
RX	0			Ddel	CutSmart	75 75	100	100	100	37°	65°	В	Lambda		D
RX	•			Dpnl	CutSmart	100	100	75	100	37°	80°	В	pBR322	CpG	b
	9			DpnII	U	25	25	100*	25	37°	65°	В	Lambda dam-	dam	b
RX	0			Dral	CutSmart	75	75	50	100	37°	65°	A	Lambda		
RX	0	e		DrallI-HF	CutSmart	< 10	50	10	100	37°	No	В	Lambda	CpG	b
	•			Drdl	CutSmart	25	50	10	100	37°	65°	A	pUC19	СрС	
RX				Eael	CutSmart	10	50	< 10	100	37°	65°	A	Lambda	dcm CpG	
RX	•			Eagl	3.1	10	25	100	10	37°	65°	В	pXba	CpG	
RX	•	e		Eagl-HF	CutSmart	25	100	100	100	37°	65°	В	pXba	CpG	
RX	•			Earl	CutSmart	50	10	< 10	100	37°	65°	В	Lambda	CpG	b, d
RX				Ecil	CutSmart	100	50	50	100	37°	65°	Α	Lambda	CpG	2
RX	•			Eco53kl	CutSmart	100	100	< 10	100	37°	65°	Α	pXba	CpG	3, b
RX	0			EcoNI	CutSmart	50	100	75	100	37°	65°	Α	Lambda		b
RX	•			Eco0109I	CutSmart	50	100	50	100	37°	65°	Α	Lambda HindIII	dcm	3
RX	•		2*site	EcoP15I	3.1 + ATP	75	100	100	100	37°	65°	Α	pUC19		е
R₩				EcoRI	U	25	100*	50	50*	37°	65°	С	Lambda	CpG	
RX	0	e		EcoRI-HF	CutSmart	10	100	< 10	100	37°	65°	С	Lambda	CpG	
RX	0			EcoRV	3.1	10	50	100	10	37°	80°	Α	Lambda	CpG	
R₩	0	e		EcoRV-HF	CutSmart	25	100	100	100	37°	65°	В	Lambda	CpG	
RX	0			Esp3I	CutSmart	100	100	< 10	100	37°	65°	В	Lambda	СрС	
RX				Fatl	2.1	10	100	50	50	55°	80°	Α	pUC19		

^{1.} Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

^{2.} Star activity may result from extended digestion.3. Star activity may result from a glycerol concentration of > 5%.

^{*} May exhibit star activity in this buffer.

Performance Chart for Restriction Enzymes (continued)

								INCUB.	INACTIV.					
1		ENZYME	SUPPLIED Nebuffer	% ACT 1.1	IVITY IN N 2.1	EBUFFERS 3.1	CUTSMART	TEMP. (°C)	TEMP. (°C)	DIL.	SUBSTRATE	METHYLA SENSITI		NOTE(S)
RX		Faul	CutSmart	100	50	10	100	55°	65°	Α	Lambda		CpG	3, b, d
RX	0	Fnu4HI	CutSmart	< 10	< 10	< 10	100	37°	No	Α	Lambda		CpG	а
RX	2+site	Fokl	CutSmart	100	100	75	100	37°	65°	Α	Lambda	dcm	CpG	3, b, d
RX	0	Fsel	CutSmart	100	75	< 10	100	37°	65°	В	pBC4	dcm	CpG	
RX	0	Fspl	CutSmart	10	100	10	100	37°	No	С	Lambda		CpG	b
RX		FspEl	CutSmart	< 10	< 10	< 10	100	37°	80°	В	pBR322	dcm		1, e
RX	•	Haell	CutSmart	25	100	10	100	37°	80°	Α	Lambda		CpG	
RX	0	Haelll	CutSmart	50	100	25	100	37°	80°	Α	Lambda			
RX		Hgal	1.1	100	100	25	100	37°	65°	Α	phiX174		CpG	1
RX	0	Hhal	CutSmart	25	100	100	100	37°	65°	Α	Lambda		CpG	
RX	•	HincII	3.1	25	100	100	100	37°	65°	В	Lambda		CpG	
RX		HindIII	2.1	25	100	50	50	37°	80°	В	Lambda			2
RX	• e	HindIII-HF	CutSmart	10	100	10	100	37°	80°	В	Lambda			
RX	0	Hinfl	CutSmart	50	100	100	100	37°	80°	Α	Lambda		CpG	
RX	•	HinP1I	CutSmart	100	100	100	100	37°	65°	Α	Lambda		CpG	
RX		Hpal	CutSmart	< 10	75*	25	100	37°	No	Α	Lambda		CpG	1
RX	•	Hpall	CutSmart	100	50	< 10	100	37°	80°	Α	Lambda		CpG	
RX	•	Hphl	CutSmart	50	50	< 10	100	37°	65°	В	Lambda	dam	CpG	b, d
RX		Hpy99I	CutSmart	50	10	< 10	100	37°	65°	Α	Lambda		CpG	
RX	0	Hpy166II	CutSmart	100	100	50	100	37°	65°	С	pBR322		CpG	
RX		Hpy188I	CutSmart	25	100	50	100	37°	65°	Α	pBR322	dam		1, b
RX		Hpy188III	CutSmart	100	100	10	100	37°	65°	В	pUC19	dam	CpG	3, b
RX	0	HpyAV	CutSmart	100	100	25	100	37°	65°		Lambda		CpG	3, b, d
RX		HpyCH4III	CutSmart	100	25	< 10	100	37°	65°	Α	Lambda			b
RX	0	HpyCH4IV	CutSmart	100	50	25	100	37°	65°	Α	pUC19		CpG	
RX	0	HpyCH4V	CutSmart	50	50	25	100	37°	65°	Α	Lambda			
RX		I-Ceul	CutSmart	10	10	10	100	37°	65°	В	pBHS Scal-linearized			
RX		I-Scel	CutSmart	10	50	25	100	37°	65°	В	pGPS2 NotI-linearized]	0.0	2
RX RX		Kasl	CutSmart	50	100	50	100	37°	65°	В	pBR322		CpG	3
RX	9 <i>e</i>	Kpnl	1.1 CutSmart	100	75	< 10	50	37°	No No	A	pXba			
RX		Kpnl-HF LpnPl	CutSmart	100 < 10	25 < 10	< 10 < 10	100 100	37° 37°	No 65°	A B	pXba			1, e
RX	0	Mbol	CutSmart	< 10 75	100	100	100	37°	65°	А	pBR322 Lambda dam-	dam	CpG	1, 5
RX	2*site	Mboll	CutSmart	100*	100	50	100	37°	65°	C	Lambda dam-	dam	opu	b
RX	2 01.0	Mfel	CutSmart	75	50	10	100	37°	No	A	Lambda	dum		2
RX	9 <i>e</i>	Mfel-HF	CutSmart	75 75	25	< 10	100	37°	No	A	Lambda			
RX	0	Mlul	3.1	10	50	100	25	37°	80°	A	Lambda		CpG	
RX	9 e	Mlul-HF	CutSmart	25	100	100	100	37°	No	A	Lambda		CpG	
RX	0	MluCl	CutSmart	100	100	100	100	37°	No	A	Lambda			
RX	0	Mlyl	CutSmart	50	50	10	100	37°	65°	A	Lambda			b, d
RX	2+site	Mmel	CutSmart + SAM	50	100	50	100	37°	65°	В	phiX174		CpG	b, c
RX	9	MnII	CutSmart	75	100	50	100	37°	65°	В	Lambda			b
RX		Mscl	CutSmart	25	100	100	100	37°	80°	С	Lambda	dcm		
RX	0	Msel	CutSmart	75	100	75	100	37°	65°	A	Lambda			
RX	0	MsII	CutSmart	50	50	< 10	100	37°	80°	A	Lambda			
RX	•	Mspl	CutSmart	75	100	50	100	37°	No	Α	Lambda			
RX	0	MspA1I	CutSmart	10	50	10	100	37°	65°	В	Lambda		CpG	
RX		MspJI	CutSmart	< 10	< 10	< 10	100	37°	65°	В	pBR322			1, e
RX	0	Mwol	CutSmart	< 10	100	100	100	60°	No	В	Lambda		CpG	
RX	2+site	Nael	CutSmart	25	25	< 10	100	37°	No	Α	pXba		CpG	b
RX	2+site	Narl	CutSmart	100	100	10	100	37°	65°	Α	pXba		CpG	
			Julion 101								p			

a. Ligation is less than 10% b. Ligation is 25% – 75%

c. Recutting after ligation is $<\!5\%$ d. Recutting after ligation is 50%-75%

e. Ligation and recutting after ligation is not applicable since the enzyme is either a nicking enzyme, is affected by methylation, or the recognition sequence contains variable sequences.

1		V	ENZYME	SUPPLIED Nebuffer	% ACT	TIVITY IN N 2.1	IEBUFFERS 3.1	CUTSMART	INCUB. TEMP. (°C)	INACTIV. TEMP. (°C)	DIL.	SUBSTRATE	METHYLATION SENSITIVITY	
RX		Sales and	Nb.BbvCl	CutSmart	25	100	100	100	37°	80°	Α	pUB		е
R₩			Nb.Bsml	3.1	< 10	50	100	10	65°	80°	Α	pBR322		е
RX			Nb.BsrDI	CutSmart	25	100	100	100	65°	80°	Α	pUC19		е
RX			Nb.BssSI	3.1	10	100	100	25	37°	No	В	pUC19		е
RX			Nb.Btsl	CutSmart	75	100	75	100	37°	80°	Α	phiX174		е
RX	0		Ncil	CutSmart	100	25	10	100	37°	No	Α	Lambda	CpG	l b
RX	0		Ncol	3.1	100	100	100	100	37°	80°	Α	Lambda		
RX	0	e	Ncol-HF	CutSmart	50	100	10	100	37°	80°	В	Lambda		
RX	•		Ndel	CutSmart	75	100	100	100	37°	65°	Α	Lambda		
RX	•	2+site	NgoMIV	CutSmart	100	50	10	100	37°	No	Α	pXba	CpG	1 1
RX	•		Nhel	2.1	100	100	10	100	37°	65°	С	Lambda HindIII	CpG	I
RX	•	e	Nhel-HF	CutSmart	100	25	< 10	100	37°	80°	С	Lambda HindIII	CpG	
RX	•		NIaIII	CutSmart	< 10	< 10	< 10	100	37°	65°	В	phiX174		
RX			NIaIV	CutSmart	10	10	10	100	37°	65°	В	pBR322	dcm	
RX		2+site	NmeAIII	CutSmart + SAM	10	10	< 10	100	37°	65°	В	phiX174		С
RX	•		Notl	3.1	< 10	50	100	25	37°	65°	С	pBC4	CpG	
RX	•	e	NotI-HF	CutSmart	25	100	25	100	37°	65°	Α	pBC4	CpG	
RX	•		Nrul	3.1	< 10	10	100	10	37°	No	Α	Lambda	dam	l b
RX	•	e	Nrul-HF	CutSmart	0	25	50	100	37°	No	Α	Lambda	dam	
RX	•		Nsil	3.1	10	75	100	25	37°	65°	В	Lambda		
RX	•	e	Nsil-HF	CutSmart	< 10	20	< 10	100	37°	80°	В	Lambda		
RX	•		Nspl	CutSmart	100	100	< 10	100	37°	65°	Α	Lambda		
RX			Nt.Alwl	CutSmart	10	100	100	100	37°	80°	Α	pUC101 dam-dcm-	dam	е
RX			Nt.BbvCI	CutSmart	50	100	10	100	37°	80°	Α	pUB	CpG	l e
RX			Nt.BsmAl	CutSmart	100	50	10	100	37	65°	Α	pBR322	CpG	l e
RX			Nt.BspQI	3.1	< 10	25	100	10	50°	80°	В	pUC19		е
RX			Nt.BstNBI	3.1	0	10	100	10	55°	80°	Α	T7		е
RX			Nt.CviPII	CutSmart	10	100	25	100	37°	65°	Α	pUC19	CpG	l e
RX	•		Pacl	CutSmart	100	75	10	100	37°	65°	Α	pNEB193		
RX	•		PaeR7I	CutSmart	25	100	10	100	37°	No	Α	Lambda HindIII	CpG	
RX			Pcil	3.1	50	75	100	50*	37°	80°	В	pXba		
RX	•		PfIFI	CutSmart	25	100	25	100	37°	65°	Α	pBC4		b
RX	•		PfIMI	3.1	0	100	100	50	37°	65°	Α	Lambda	dcm	3, b, d
RX			PI-PspI	U	10	10	10	10	65°	No	В	pAKR XmnI		
RX			PI-Scel	U	10	10	10	10	37°	65°	В	pBSvdeX XmnI		
RX		2+site	Plel	CutSmart	25	50	25	100	37°	65°	Α	Lambda	CpG	b, d
RX		2+site	PluTl	CutSmart	100	25	< 10	100	37°	65°	Α	pXba	CpG	l b
RX	•		Pmel	CutSmart	< 10	50	10	100	37°	65°	Α	Lambda	CpG	
RX	•		PmII	CutSmart	100	50	< 10	100	37°	65°	Α	Lambda HindIII	CpG	
RX	0		PpuMI	CutSmart	< 10	< 10	< 10	100	37°	No	В	Lambda HindIII	dcm	
RX	•		PshAl	CutSmart	25	50	10	100	37°	65°	Α	Lambda	CpG	
RX			Psil	CutSmart	10	100	10	100	37°	65°	В	Lambda		3
RX			PspGI	CutSmart	25	100	50	100	75°	No	Α	T7	dcm	3
RX			PspOMI	CutSmart	10	10	< 10	100	37°	65°	В	pXba	dcm CpG	
RX			PspXI	CutSmart	< 10	100	25	100	37°	No	В	Lambda HindIII	CpG	
RX	0		Pstl	3.1	75	75	100	50*	37°	80°	С	Lambda		
RX	0	e	PstI-HF	CutSmart	10	75	50	100	37°	No	С	Lambda		
RX	0		Pvul	3.1	< 10	25	100	< 10	37°	No	В	pXba	CpG	
RX	•	e	Pvul-HF	CutSmart	25	100	100	100	37°	No	В	pXba	CpG	
RX	0		Pvull	3.1	50	100	100	100*	37°	No	В	Lambda		
RX	0	e	PvuII-HF	CutSmart	< 10	< 10	< 10	100	37°	No	В	Lambda		

 $^{1. \ \, \}text{Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of $>5\%.}$

Star activity may result from extended digestion.
 Star activity may result from a glycerol concentration of > 5%.

^{*} May exhibit star activity in this buffer.

Performance Chart for Restriction Enzymes (continued)

Th	1	M							INCUB.	INACTIV.					
	-	F	ENZYME	SUPPLIED Nebuffer	% ACT 1.1	IVITY IN N 2.1	NEBUFFERS 3.1	CUTSMART	TEMP. (°C)	TEMP. (°C)	DIL.	SUBSTRATE	METHYLAT SENSITIV		NOTE(S)
RX	0		Rsal	CutSmart	25	50	< 10	100	37°	No	А	Lambda		CpG	
R\{		2*site	Rsrll	CutSmart	25	75	10	100	37°	65°	С	Lambda		CpG	
RX	0		Sacl	1.1	100	50	10	100	37°	65°	Α	Lambda HindIII			
RX	0	e	SacI-HF	CutSmart	10	50	< 10	100	37°	65°	Α	Lambda HindIII		CpG	
RX	0	2*site	SacII	CutSmart	10	100	10	100	37°	65°	Α	pXba		CpG	
RX	0		Sall	3.1	< 10	< 10	100	< 10	37°	65°	Α	Lambda HindIII		CpG	
RX	0	e	Sall-HF	CutSmart	10	100	100	100	37°	65°	Α	Lambda HindIII		CpG	
RX	0		Sapl	CutSmart	75	50	< 10	100	37°	65°	В	Lambda			
RX			Sau3AI	1.1	100	50	10	100	37°	65°	Α	Lambda		CpG	b
RX			Sau96I	CutSmart	50	100	100	100	37°	65°	Α	Lambda	dcm	CpG	
R₩	0		Sbfl	CutSmart	50	25	< 10	100	37°	80°	Α	Lambda			3
RX	0	e	Sbfl-HF	CutSmart	50	25	< 10	100	37°	80°	В	Lambda			
RX	0	e	Scal-HF	CutSmart	100	100	10	100	37°	80°	В	Lambda			
RX			ScrFI	CutSmart	100	100	100	100	37°	65°	С	Lambda	dcm	CpG	2, a
RX			SexAl	CutSmart	100	75	50	100	37°	65°	Α	pBC4 dcm-	dcm		3, b, d
RX			SfaNI	3.1	< 10	75	100	25	37°	65°	В	phiX174		CpG	3, b
RX			SfcI	CutSmart	75	50	25	100	37°	65°	В	Lambda			3
RX	0	2*site	Sfil	CutSmart	25	100	50	100	50°	No	С	pXba	dcm	CpG	
RX	0		Sfol	CutSmart	50	100	100	100	37°	No	В	Lambda HindIII	dcm	CpG	
RX		2*site	09.7	CutSmart	100	100	10	100	37°	65°	Α	Lambda		CpG	1
RX	•		Smal	CutSmart	< 10	< 10	< 10	100	25°	65°	В	Lambda HindIII		CpG	b
RX			Smll	CutSmart	25	75	25	100	55°	No	Α	Lambda			b
RX			SnaBl	CutSmart	50	50	10	100	37°	80°	Α	T7		CpG	1
RX	0	_	Spel	CutSmart	75	100	25	100	37°	80°	С	Adenovirus-2			
RX	0	e	Spel-HF	CutSmart	25	50	10	100	37°	80°	C	pXba			
RX			Sphl	2.1	100	100	50	100	37°	65°	В	Lambda			2
RX	0		SphI-HF	CutSmart	50	25	10	100	37°	65°	В	Lambda			
RX DW	0	e	Srfl	CutSmart	10	50	0	100	37°	65°	В	pNEB193-SrFI		CpG	
RX	0		Sspl	U	50	100	50	50	37°	65°	С	Lambda			
RX	0	e	SspI-HF	CutSmart	25	100	< 10	100	37°	65°	В	Lambda	r men		
RX DD	•		Stul	CutSmart	50	100	50	100	37°	No	A	Lambda	dcm dcm	CpG	
R* R*	9		StyD4I	CutSmart	10	100	100	100	37° 37°	65°	В	Lambda	uciii	Сра	h
RX	•	e	Styl-HF	3.1	10	25	100	10 100	37°	65° 65°	A	Lambda Lambda			b
RX	0		Swal	CutSmart 3.1	25 10	100 10	25 100	100	25°	65°	A B	pXba			b, d
RX	0		Taq ^a l	CutSmart	50	75	100	100	65°	80°	В	Lambda	dam		b, u
RX	0		Tfil	CutSmart	50	100	100	100	65°	No	С	Lambda	uum	CpG	
-4115	0		Tsel	CutSmart	75	100	100	100	65°	No No	В	Lambda		СрС	3
			Tsp45I	CutSmart	100	50	< 10	100	65°	No	A	Lambda			J
	0		TspMI	CutSmart	50*	75*	50*	100	75°	No	В	pUCAdeno		CpG	d
RX	0		TspRI	CutSmart	25	50	25	100	65°	No	В	Lambda			u
RX	0		Tth1111	CutSmart	25	100	25	100	65°	No	В	pBC4			b
R:	0		Xbal	CutSmart	< 10	100	75	100	37°	65°	A	Lambda Hindlll dam	dam		U
RX			Xcml	2.1	10	100	25	100	37°	65°	C	Lambda			2
R::	0		Xhol	CutSmart	75	100	100	100	37°	65°	A	Lambda Hindlll			b
RX	0		Xmal	CutSmart	25	50	< 10	100	37°	65°	A	pXba		CpG	3
RX	0		XmnI	CutSmart	50	75	< 10	100	37°	65°	A	Lambda			b
RX			Zral	CutSmart	100	25	10	100	37°	80°	В	Lambda		CpG	5
=11115			Liai	Outoman	100	20	10	100	01	00	D	Lambua			

a. Ligation is less than 10% b. Ligation is 25% – 75%

c. Recutting after ligation is <5% d. Recutting after ligation is 50%-75%

e. Ligation and recutting after ligation is not applicable since the enzyme is either a nicking enzyme, is affected by methylation, or the recognition sequence contains variable sequences.

Activity of Enzymes at 37°C

Listed below is the percentage of activity exhibited at 37°C for enzymes that have an optimal incubation temperature higher (thermophiles) or lower (25°C) than 37°C.

ENZYME	OPTIMAL TEMP. (°C)	% ACTIVITY AT 37°C
AbaSI	25°	0
Apal	25°	100*
ApeKI	75°	10
Apol	50°	50
Bael	25°	20
BcII	50°	50
BfuAI	50°	50
BsaBI	60°	20
BsaJI	60°	20
BsaWl	60°	20
BsiEl	60°	30
BsiHKAI	65°	5
BsiWI	55°	50
BsII	55°	30
BsmAl	55°	50
BsmBI	55°	20
BsmFI	65°	50
Bsml	65°	20

ENZYME	OPTIMAL TEMP. (°C)	% ACTIVITY AT 37°C
BspCNI	25°	75
BspQI	50°	10
Bsrl	65°	20
BsrDI	65°	30
BssHII	50°	75
BstAPI	60°	10
BstBI	65°	10
BstEII	60°	50
BstNI	60°	30
BstUI	60°	20
BstYI	60°	30
BtgZl	60°	75
BtsI-v2	55°	75
BtsCI	50°	50
BtsIMutI	55°	N/A
CviAII	25°	20
CviQI	25°	10
Fatl	55°	20

ENZYME	OPTIMAL TEMP. (°C)	% ACTIVITY At 37°C
Faul	55°	20
Mwol	60°	10
Nb.Bsml	65°	25
Nb.BsrDI	65°	75
Nt.BspQI	50°	80
Nt.BstNBI	55°	10
PI-PspI	65°	5
PspGI	75°	10
Sfil	50°	10
Smal	25°	50
SmII	55°	10
Swal	25°	50
TaqαI	65°	10
Tfil	65°	10
Tsel	65°	20
Tsp45I	65°	10
TspMI	75°	20
TspRI	65°	10
Tth111I	65°	10

Activity of DNA Modifying Enzymes in CutSmart Buffer

A selection of DNA modifying enzymes were assayed in CutSmart Buffer, in lieu of their supplied buffers. Functional activity was compared to the activity in its supplied buffer, plus required supplements. Reactions were set up according to the recommended reaction conditions, with CutSmart Buffer replacing the supplied buffer.

ENZYME	ACTIVITY IN CUTSMART	REQUIRED SUPPLEMENTS
Alkaline Phosphatase (CIP)	+++	
Antarctic Phosphatase	+++	Requires Zn ²⁺
Bst DNA Polymerase	+++	
CpG Methyltransferase (M. SssI)	+++	
DNA Polymerase I	+++	
DNA Polymerase I, Large (Klenow) Fragment	+++	
DNA Polymerase Klenow Exo-	+++	
DNase I (RNase-free)	+++	Requires Ca2+
E. coli DNA Ligase	+++	Requires NAD
Endonuclease III (Nth), recombinant	+++	
Endonuclease VIII	+++	
Exonuclease I	+++	
Exonuclease III	+++	
Exonuclease VII	+++	
Exonuclease V (Rec BCD)	+++	Requires ATP
GpC Methyltransferase (M. CviPI)	+	Requires DTT
Lambda Exonuclease	++	
McrBC	+++	

- +++ full functional activity
- + + 50-100% functional activity
- + 0-50% functional activity

ENZYME	ACTIVITY IN CUTSMART	REQUIRED Supplements
Micrococcal Nuclease	+++	Requires Ca ²⁺
Nuclease Bal-31	+++	
phi29 DNA Polymerase	+++	
Quick Dephosphorylation Kit	+++	
RecJ _f	+++	
Shrimp Alkaline Phosphatase (rSAP)	+++	
T3 DNA Ligase	+++	Requires ATP + PEG
T4 DNA Ligase	+++	Requires ATP
T4 DNA Polymerase	+++	
T4 Phage β-glucosyltransferase (T4-BGT)	+++	
T4 Polynucleotide Kinase	+++	Requires ATP + DTT
T4 PNK (3´ phosphatase minus)	+++	Requires ATP + DTT
T5 Exonuclease	+++	
T7 DNA Ligase	+++	Requires ATP + PEG
T7 DNA Polymerase (unmodified)	+++	
T7 Exonuclease	+++	
Thermolabile Exol	+++	
USER Enzyme, recombinant	+++	

^{*}Apal has 100% activity at 37°C, however the half-life of this enzyme at 37°C is only 30 minutes.

Tips for Avoiding Star Activity

Under non-standard reaction conditions, some restriction enzymes are capable of cleaving sequences which are similar but not identical to their defined recognition sequence. This altered specificity has been termed "star activity". Although the propensity for star activity varies, the vast majority of enzymes from New England Biolabs will not exhibit star activity when used under recommended conditions in their supplied NEBuffers. If an enzyme has been reported to exhibit star activity, it will be indicated in the product entry found in the catalog, on the supplied card and on our website.

CONDITIONS THAT Contribute to Star activity	STEPS THAT CAN BE TAKEN TO INHIBIT STAR ACTIVITY
High gluggral appeartment (s. E9/ y/h)	Restriction enzymes are stored in 50% glycerol, therefore the amount of enzyme added should not exceed 10% of the total reaction volume.
High glycerol concentration (> 5% v/v)	Use the standard 50 μl reaction volume to reduce evaporation during incubation.
High concentration of enzyme/µg of DNA ratio (varies with each enzyme, usually 100 units/µg)	Use the fewest units possible to achieve digestion. This avoids overdigestion and reduces the final glycerol concentration in the reaction.
Non-optimal buffer	Whenever possible, set up reactions in the recommended buffer. Buffers with differing ionic strengths and pHs may contribute to star activity.
Prolonged reaction time	Use the minimum reaction time required for complete digestion. Prolonged incubation may result in increased star activity, as well as evaporation.
Presence of organic solvents [DMSO, ethanol (4), ethylene glycol, dimethylacetamide, dimethylformamide, sulphalane (5)]	Make sure the reaction is free of any organic solvents, such as alcohols, that might be present in the DNA preparation.
Substitution of Mg ²⁺ with other divalent cations (Mn ²⁺ , Cu ²⁺ , Co ²⁺ , Zn ²⁺)	Use Mg ²⁺ as the divalent cation. Other divalent cations may not fit correctly into the active site of the restriction enzyme, possibly interfering with proper recognition.

Note: The relative significance of each of these altered conditions will vary from enzyme to enzyme.

New England Biolabs recommends setting up restriction enzyme digests in a 50 µl reaction volume. However, different methods may require smaller reaction volumes. When performing restriction enzyme digests in smaller reaction volumes, extra care must be taken to follow the steps listed above to avoid star activity. Alternatively, using our line of **High Fidelity (HF) restriction enzymes** will allow greater flexibility in reaction setup. Please visit **www.neb.com/HF** frequently to learn about new additions to the HF restriction enzyme product line.

Reference:

(1) Nasri, M. and Thomas, D. (1986) Nucleic Acids Res. 14, 811.

TOOLS & RESOURCES

Visit NEBRestrictionEnzymes.com to find:

- Online tutorials on how to avoid star activity, and for setting up restriction enzyme digests
- The full list of HF enzymes available
- · Troubleshooting guides



Giron has been with NEB for over 2 years as a Development Scientist. Giron's family lives in Germany and he frequently travels to spend time with them. Giron is also an experienced spiritual teacher and feng shui consultant, as well as a skilled violinist.





APPENDIX

High-Fidelity (HF) Restriction Enzymes

As part of our ongoing commitment to the study and improvement of restriction enzymes, NEB offers a line of High-Fidelity (HF) restriction enzymes. These engineered enzymes have the same specificity as the native enzymes, are all active in CutSmart Buffer and have reduced star activity. Star activity, or off-target cleavage, is an intrinsic property of restriction enzymes. Most restriction enzymes will not exhibit star activity when used under recommended reaction conditions. However, for enzymes that have reported star activity, extra caution must be taken to set up reactions according to the recommended conditions to avoid unwanted cleavage.

Many techniques such as cloning, genotyping, mutational analysis, mapping, probe preparation, sequencing and methylation detection employ a wide range of reaction conditions and require the use of enzymes under suboptimal conditions. HF enzymes with reduced star activity offer increased flexibility to reaction setup and help maximize results under a wide range of conditions.

In addition to reduced star activity, all of these engineered enzymes work optimally in CutSmart Buffer, which has the highest level of enzyme compatibility and will simplify double digest reactions. They are all Time-Saver qualified and digest substrate DNA in 5–15 minutes, and can also be incubated overnight without degredation of DNA. HF enzymes are available at the same price as the native enzymes and are supplied with purple loading dye.

Visit www.neb.com/HF to learn more about HF enzymes.

TOOLS & RESOURCES

Visit NEBRestrictionEnzymes.com to find:

- The full list of HF enzymes available
- Online tutorials on how to avoid star activity and setting up digests using the Time-Saver protocol



Reduced Star Activities of HF Enzymes

The following table indicates the number of units of HF enzyme that can be used compared to the native enzyme before any significant star activity is detected. The HF Factor refers to the X-fold increase in fidelity that is achieved by choosing an HF enzyme. This data clearly illustrates the flexibility that is offered by using an HF restriction enzyme.

PRODUCT Name	PRODUCT Number	BUFFER†	MAXIMUM UNITS WITH NO STAR ACTIVITY*	HF FACTOR
Agel-HF	#R3552	CutSmart	≥ 250	≥ 8
Agel	#R0552	1.1	32	
Apol-HF	#R3566	CutSmart	500	25
Apol	#R0566	3.1	20	
BamHI-HF	#R3136	CutSmart	≥ 4,000	≥ 125
BamHI	#R0136	3.1	32	
BbsI-HF	#R3539	CutSmart	≥ 500	≥ 4
Bbsl	#R0539	2.1	120	
BcII-HF	#R3160	CutSmart	500	16
BcII	#R0160	3.1	32	
Bmtl-HF	#R3658	CutSmart	1,000,000	62,500
Bmtl	#R0658	3.1	32	
Bsal-HFv2	#R3733	CutSmart	500	16
Bsal	#R0535	CutSmart	32	
BsiWI-HF	#R3553	CutSmart	100	1
BsiWl	#R0553	3.1	100	
BsrGI-HF	#R3575	CutSmart	≥ 2,000	≥ 62
BsrGl	#R0575	2.1	16	
BstEII-HF	#R3162	CutSmart	> 2,000	> 125
BstEII	#R0162	3.1	16	
BstZ17I-HF	#R3594	CutSmart	500	25
BstZ17I**	N/A	CutSmart	20	
DrallI-HF	#R3510	CutSmart	≥ 2,000	≥ 1,000
DrallI**	N/A	3.1	2	
Eagl-HF	#R3505	CutSmart	500	2
Eagl	#R0505	3.1	250	
EcoRI-HF	#R3101	CutSmart	16,000	64
EcoRI	#R0101	U	250	
EcoRV-HF	#R3195	CutSmart	≥ 64,000	≥ 64
EcoRV	#R0195	3.1	1,000	
HindIII-HF	#R3104	CutSmart	≥ 500,000	≥ 2,000
HindIII	#R0104	2.1	250	
Kpnl-HF	#R3142	CutSmart	≥ 1,000,000	≥ 62,500
Kpnl	#R0142	1.1	16	
Mfel-HF	#R3589	CutSmart	≥ 500	≥ 16
Mfel	#R0589	CutSmart	32	

MIUI #F NCOI-HF #F NCOI #F Nhel-HF #F Nhel #F NotI-HF #F NotI #F NruI-HF #F NsiI-HF #F	80198 3 80193 C 80193 3 831311 C 80131 2 80131 2 83189 C 80189 3 83192 C 80192 3	CutSmart 3.1 CutSmart CutSmart CutSmart CutSmart CutSmart CutSmart CutSmart	≥ 4,000 ≥ 2,000 ≥ 64,000 120 ≥ 32,000 120 ≥ 64,000 4,000	≥ 530
Ncol-HF	33193 C 30193 3 33131 C 30131 2 33189 C 30189 3 33192 C 30192 3	CutSmart 3.1 CutSmart 2.1 CutSmart 3.1 CutSmart 3.1	≥ 64,000 120 ≥ 32,000 120 ≥ 64,000	≥ 266
Ncol	R0193 3 R3131 C R0131 2 R3189 C R0189 3 R3192 C R0192 3	B.1 CutSmart P.1 CutSmart	120 ≥ 32,000 120 ≥ 64,000	≥ 266
Nhel-HF #F Nhel #F Notl-HF #F Notl #F Nrul-HF #F Nrul #F Nsil-HF #F Nsil #F	R3131 C R0131 2 R3189 C R0189 3 R3192 C R0192 3	CutSmart 2.1 CutSmart 3.1	≥ 32,000 120 ≥ 64,000	
Nhel #F NotI-HF #F NotI #F NruI-HF #F NruI #F NsiI-HF #F NsiI #F	R0131 2 R3189 0 R0189 3 R3192 0 R0192 3	2.1 CutSmart 3.1	120 ≥ 64,000	
NotI-HF #F NotI #F NruI-HF #F NruI #F NsiI-HF #F NsiI #F	R3189 C R0189 3 R3192 C R0192 3	CutSmart 3.1	≥ 64,000	\ 1C
NotI #F Nrul-HF #F Nrul #F Nsil-HF #F Nsil #F	R0189 3 R3192 C	3.1		× 10
Nrul-HF #F Nrul #F Nsil-HF #F Nsil #F	R3192 C		4.000	≥ 16
Nrul #F Nsil-HF #F Nsil #F	R0192 3	CutSmart	4,000	
Nsil-HF #F			≥ 32,000	64
Nsil #F	20407	8.1	≥ 500	
	R3127 C	CutSmart	≥ 8,000	2
PstI-HF #F	R0127 3	8.1	≥ 4,000	
	R3140 C	CutSmart	4,000	33
PstI #F	R0140 3	8.1	120	
Pvul-HF #F	R3150 C	CutSmart	≥ 16,000	≥ 32
Pvul #F	R0150 3	8.1	500	
PvuII-HF #F	R3151 C	CutSmart	500	32
Pvull #F	R0151 3	3.1	16	
SacI-HF #F	R3156 C	CutSmart	≥ 32,000	≥ 266
Sacl #F	R0156 1	.1	120	
Sall-HF #F	R3138 C	CutSmart	≥ 32,000	≥ 8,000
Sall #F	R0138 3	3.1	4	
Sbfl-HF #F	R3642 C	CutSmart	250	32
Sbfl #F	R0642 C	CutSmart	8	
Scal-HF #F	R3122 (CutSmart	250	62
Scal** N/	/A 3	3.1	4	
Spel-HF #F	R3133 (CutSmart	≥ 8,000	≥ 16
Spel #F	R0133 (CutSmart	500	
SphI-HF #F	R3182 (CutSmart	8,000	250
SphI #F	R0182 2	2.1	32	
SspI-HF #F	R3132 (CutSmart	500	16
Sspl #F	R0132 L	J	32	
Styl-HF #F	R3500 C	CutSmart	4,000	125
Styl #F		3.1		

- [†] Wild type enzymes were tested in supplied buffer for comparisons.
- * Wei, H. et al (2008) Nucleic Acids Reseach 36, e50.
- ** No longer available

Time-Saver Qualified Restriction Enzymes

Whether you are quickly screening large numbers of clones or setting up overnight digests, you will benefit from the high quality of our enzymes. Typically, a restriction digest involves the incubation of 1 μ l of enzyme with 1 μ g of purified DNA in a final volume of 50 μ l for 1 hour. However, to speed up the screening process, choose one of NEB's enzymes that are Time-Saver qualified. Over 185 of our enzymes will digest 1 μ g of substrate DNA in 5-15 minutes using 1 μ l of enzyme under recommended reaction conditions, and can also be used safely in overnight digestions. Unlike other suppliers, there is no special formulation, change in concentration or need to buy more expensive, new lines of enzymes to achieve digestion in 5-15 minutes. Nor do you have to worry if you incubate too long.

In an effort to provide you with as much information as possible, NEB has tested all of its enzymes on unit assay substrate as well as plasmid substrate and PCR Fragments. We recommend that this data be used as a guide, as it is not definitive for all plasmids. Restriction enzymes can often show site preference, presumably determined by the sequence flanking the recognition site. In addition, supercoiled DNA may have varying rates of cleavage. For more information, visit www.neb.com/TimeSaver. Note that there are some enzymes indicated below that can cut in 5-15 minutes, but cannot be incubated overnight. These are not Time-Saver qualified.

Since all of our enzymes are rigorously tested for nuclease contamination, you can also safely set up digests for long periods of time without sample degradation. Only NEB Time-Saver qualified enzymes offer power and flexibility — the power to digest in 5-15 minutes and the flexibility to withstand overnight digestions with no loss of substrate.

TOOLS & RESOURCES

Visit www.neb.com/TimeSaver to find:

- The full list of Time-Saver qualified restriction enzymes available
- Online tutorials on using Time-Saver qualified enzymes to speed up restriction enzyme digests

Chart Legend

- digests in 5 minutes
- digests in 15 minutes
- not completely digested in 15 minutes

ENZYME	UNIT ASSAY	SUBSTI PLASMID	RATE PCR
AatII	•	A	•
Accl	•	A	A
Acc65I	•	A	•
Acil	•	•	•
AcII	•	-	A
Acul	•	A	A
AfIII	•	•	•
Agel-HF	•	•	•
Ahdl	•	•	-
Alul	•	A	•
AlwNI	•	•	A
Apal	•	•	•
ApaLl	•	•	A
ApeKI	•	•	A
Apol	•	•	•
Apol-HF	•	•	A
Ascl	•	•	NT
Asel	•	•	NT
Aval	•	A	A
Avall	•	•	•
AvrII	•	NT	NT
Bael	•	•	A
BaeGI	•	A	A
BamHI	•	•	A
BamHI-HF	•	•	•
Bbsl	•	A	A
BbsI-HF	•	A	A
Bbvl	•	A	A
Bccl	•	A	A
BceAl	-	•	A
BciVI	•	•	A
BcII	•	A	A
BcII-HF	•	A	A
BcoDI	•	•	A

ENZYME	UNIT ASSAY	SUBST Plasmid	TRATE PCR
BfuAl	•	•	A
BfuCl	•	A	•
BgII	•	•	A
BgIII	•	•	A
Blpl	•	•	•
BmgBl	•	•	A
Bmrl	-	A	•
BmtI-HF	•	•	A
BpuEl	•	•	A
Bsal	•	•	•
Bsal-HFv2	-	•	A
BsaAl	•	•	
BsaHI	-	-	•
BsaWI	•	A	A
BsaXI	•	A	A
BseRI	•	•	
Bsgl	•	•	A
BsiEl	•	A	A
BsiHKAI	•	•	A
BsiWI	•	•	A
BsiWI-HF	-	•	A
BsII	•	•	•
Bsml	•	•	A
BsmAl	•	A	•
BsmBI	-	A	A
BsmFI	•	•	A
BsoBI	•	-	•
Bsp1286I	•	•	A
BspCNI	-	A	A
BspEl	•	A	A
BspHI	-	•	•
BspQI	•	•	A
Bsrl	•	•	A .
BsrBI	•	•	A
BsrDI	•	•	A

ENZYME	UNIT ASSAY	SUBS [*] PLASMID	TRATE PCR
BsrFI-v2		A	A
BsrGI		A	A
BsrGI-HF		•	A
BssHII	•	A	A
BssSI-v2		A	A
BstBI	•	•	A
BstEII	•	•	A
BstEII-HF	•	•	•
BstNI	•	•	A
BstUI	•	•	A
BstXI	•	•	A
BstYI	•	•	A
BstZ17I-HF	•	•	A
Bsu36I	•	A	•
BtgI	•	•	•
BtsI-v2	•	•	•
BtsCI	•	•	A
Cac8I	•	A	A
Clal	•	•	A
CspCI	•	•	A
CviAII	•	•	•
CviQI	•	•	•
Ddel	•	•	•
DpnI	•	•	A
DpnII	•	A	•
Dral	•	•	•
DrallI-HF	•	•	A
Drdl	•	•	•
Eagl	•	A	A
Eagl-HF	•	•	A
Earl	-	•	A
Eco53KI	•	•	•
EcoNI	•	•	•
Eco0109I	•	A	A
EcoP15I	•	A	A

ENZYME	UNIT ASSAY	SUBS [*] Plasmid	TRATE PCR
EcoRI	•	•	A
EcoRI-HF	•	•	•
EcoRV	•	•	A
EcoRV-HF	•	•	A
Esp3I		•	A
Fnu4HI	•	-	-
Fokl	•	•	•
Fsel	•	•	A
FspI	-	A	-
Haell		A	A
HaellI	•	•	•
Hgal	•	A	A
Hhal	•	-	A
HincII	•	A	•
HindIII-HF	•	•	•
Hinfl	•	•	•
HinP1I	•	A	•
Hpall	•	•	A
Hphl	•	A	A
Hpy166II	•	•	•
HpyAV	•	•	NT
HpyCH4IV	•	•	•
HpyCH4V	•	•	•
KpnI	•	•	•
KpnI-HF	•	•	•
Mbol	•	A	•
Mboll	•	•	•
Mfel	•	•	•
Mfel-HF	•	•	•
Mlul	•	•	•
Mlul-HF	•	•	A
MluCl	•	•	A
Mlyl	•	A	•
Mmel	•	•	A
MnII	•	•	•
Msel			•

ENZYME	UNIT Assay	SUBS [*] Plasmid	TRATE PCR
MsII	•	•	•
Mspl	•	•	•
MspA1I	•	•	•
Mwol		A	A
Ncil	•	•	•
Ncol	•	-	A
Ncol-HF	•	•	•
Ndel	•	•	A
NgoMIV	•	•	A
Nhel	•	•	A
Nhel-HF	•	•	-
NIaIII		A	
NmeAIII	•	A	A
Notl	•	•	A
Notl-HF	•	•	•
Nrul	•	•	A
Nrul-HF	•	-	A
Nsil	•	•	•
Nsil-HF	•	•	-
Nspl	•	•	A
Pacl	•	•	•
PaeR7I	•	A	A
Pflfl	•	-	A
PfIMI	•	A	A
Pmel	•	-	NT
PmII	•	A	A
PpuMI	•	A	A
PshAl		-	-
Pstl	•	•	•
PstI-HF	•	•	•
Pvul	•	A	•
Pvul-HF	•	•	•
Pvull	•	•	A
PvuII-HF	•	•	A
Rsal	•	•	•
Sacl	•	•	A

	UNIT SUBSTRATE		TRATE
ENZYME	ASSAY	PLASMID	PCR
SacI-HF	•	•	•
SacII	•	A	A
Sall	•	•	A
Sall-HF	•	•	A
Sapl	•	A	A
Sbfl	•	•	A
Sbfl-HF	•	•	A
Scal-HF	•	•	A
SfaNI	A	A	•
Sfil	•	A	A
Sfol	•	•	•
Smal	•	-	-
Spel	•	•	•
Spel-HF	-	-	A
Sphl	•	•	A
SphI-HF	•	•	A
Srfl	•	•	A
Sspl	•	•	A
SspI-HF	•	•	A
Stul	•	A	A
Styl	•	A	A
Styl-HF	•	•	A
StyD4I	•	A	A
Swal	•	A	A
Taq∝l	•	•	A
Tfil	•	•	A
Tsel	•	A .	A
TspMI	•	•	A
TspRI	•	•	A
Tth111I	•	•	A
Xbal	•	•	A
Xhol	•	•	A
Xmal	•	A .	-
Xmnl	•	•	A

Cross Index of Recognition Sequences

Sequences at the top of each column are written 5´ to 3´ according to convention. Open squares at the left of each row are place holders for nucleotides within a restriction enzyme recognition sequence; arrowheads indicate the point of cleavage.

Sequences of complementary strands and their cleavage sites are implied.

blue type = enzymes that recognize only one sequence

black type = enzymes that recognize multiple sequences (degenerate)

Palindromic Tetra- and Hexa-Nucleotide Recognition Sequences

	AATT	ACGT	AGCT	ATAT	CATG	CCGG	CGCG	CTAG	GATC	GCGC	GGCC	GTAC	TATA	TCGA	TGCA	TTAA
Y 0000	MluCl				Fatl				BfuCl Dpnll Mbol Sau3Al							
-		HpyCH2IV			CviAll	Mspl Hpall		Bfal		HinP1I		Csp6l CviQl		Taql		Msel
			Alul CviKI-1				BstUI		Dpnl		HaellI Phol CviKI-1	Rsal			HpyCH4V	
									BstKTI	Hhal						
		Tail			NIalli											
A ▼□□□□ T	Apol •		HindIII •		Pcil AfIIII	Agel • BsrFl BsaWl	Mlul ● AfIIII	Spel •	BgIII BstYII			Tatl				
A 🕇 T		AcII												Clal BspDl		Asel
A				Sspl •						Afel	Stul	Scal •				
A																
A Y T					Nspl					Haell					Nsile	
C [▼] □□□□ G	Mfel •				Ncol • Styl • Btgl	TspMI Xmal Acol Aval BsaJI BsoBI	BsaJI BtgI	Avrll BsaJI Styl •			Eagl ● Eael	BsiWI •	SfcI	PaeR7I Tlil Xhol Aval BsoBI Smll	SfcI	AfIII Smll
C O G				Ndel	BsăJI	BsaJI BsoBI BmeT110I								BsoBI Smll BmeT110I		
C - To G		PmII BsaAI	Pvull • MspA1I			Smal	MspA1I									
C G							Sacil		Pvul • BsiEl		BsiEI					
C G															Pstl •	
G C	EcoRI • Apol •					NgoMIV BsrFI	BssHII	Nhel •	BamHI • BstYI	Kasl Banl	PspOMI	Acc65I Banl		Sall •	ApaLI	
G OCOC		BsaHI								Narl BsaHl			AccI	Accl		
G and C		Zral	Ecl136II Eco53KI	EcoRV •	Cac8I	Nael Cac8l	Cac8I	Cac8I	NIaIV	Sfol NlaIV	NIaIV	NIaIV	BstZ17I • Hpy8I Hpy166II	HincII Hpy8I Hpy166II	Hpy81 Hpy16611	Hpal HincII Hpy8I
G C																Hpy166II
G □□□□▼C		AatII	Sacl • Banll BsiHKAI Bsp1286I		Sphl • Nspl			Bmtl •		PluTI Bbel Haell	Apal Banll BaeGl Bsp1286l	Kpnl ●			BaeGI Bsp1286I BsiHKAI	
T* A					BspHI	BspEI BsaWI Acol		Xbal	BcII •		Eael	BsrGI • Tatl				
T - A					Hpy188III	Hpy188III	Hpy188III	Hpy188III						BstBI		
T		SnaBI BsaAI					Nrul •			Fspl	MscI		Psil			Dral
T A																
T																

[•] HF (High-Fidelity) versions of these enzymes are available for simplified reactions and reduced star activity, at no additional cost. See page 301.

Palindromic Penta-Nucleotide Recognition Sequences

	AA 🗆 TT	AC □ GT	AG 🗆 CT	AT 🗆 AT	CA 🗆 TG		CG 🗆 CG	CT 🗆 AG	GA □ TC	GC □ GC	GG □ CC	GT 🗆 AC	TA 🗆 TA	TC □ GA	TG 🗆 CA	TT 🗆 AA
▼ □□ N □□						BssKI StyD4I						MaeIII				
						-		Ddel	Hinfl		Sau96I					
						ScrFI				Fnu4HI	BmgT120I					
N~		HpyCH4III												Hpy188I		
N 💑											Fmul					
N																
Y A						PspGI										
A									Tfil	ApeKI Tsel	Avall					
^ A						BstNI										
A ▼																
A																
A							Hpy99I									
₹ G												Tsp45I				
G C																
G G G G G G G G G G G G G G G G G G G						Ncil										
G T																
C										Taul BspUl						
G ~																

Single Letter Code:

In addition, see homing endonucleases on pages 55–56.

Note:

Enzymes marked with a "_A" are available from NEB.

 HF (High-Fidelity) versions of these enzymes are available for simplified reactions and reduced star activity, at no additional cost. See page 301.

Enzymes marked with a " \otimes " are not currently commercially available.

	DECIFICITIES OD	EATER THAN 6 BASES
ა		
	Aarl	CACCTGC(4/8)
	Aba6411II⊗	CRRTAAG
	AbaCIII⊗	CTATCAV
	Absl	CC/TCGAGG
	AcoY31II⊗	TAGCRAB
	AhyRBAHI⊗	GCYYGAC
	AhyYL17l⊗	YAAMGAG
	Ajul	(7/12)GAANNNNNNNTTGG(11/6)
	Alol	(7/12)GAACNNNNNNTCC(12/7)
	AlwFl⊗	GAAAYNNNNRTG
	AquIV⊗	GRGGAAG(19/17)
	Arsl	(8/13)GACNNNNNNTTYG(11/6)
	Ascl	GG/CGCGCC
_	710101	GCGAT/CGC
	Asp103l⊗	CGRAGGC
	AspJHL3II⊗	CGCCCAG
	AspNIH4III⊗	AAGAACB
	Asp114pII⊗	AGCABCC
_	Bael	(10/15)ACNNNNGTAYC(12/7)
	Barl	(7/12)GAAGNNNNNNTAC(12/7)
	Bbr57III⊗	GTRAAYG
_	BbvCI	CCTCAGC(-5/-2)
	BkrAM31DI⊗	RTTAAATM
	Ble402II⊗	GRAGCAG
	Bsp460III⊗	CGCGCAG
_	BspQI	GCTCTTC(1/4)
	CalB3II⊗	GRTTRAG
	Cbo67071IV⊗	GCRGAAG
	CcrNAIII⊗	CGACCAG
	Cdi81III⊗	GCMGAAG
	Cgl13032II⊗	ACGABGG
	Cly7489II⊗	AAAAGRG

S	PECIFICITIES GR	EATER THAN 6 BASES (CONT.)
A		(11/13)CAANNNNNGTGG(12/10)
	Ecl35734I⊗	GAAAYTC
	Eco4465II⊗	GAAABCC
	Eco43896II⊗	CRARCAG
\blacktriangle	Fsel	GGCCGG/CC
	FspAl	RTGC/GCAY
	FspPK15I⊗	GARGAAG
	GauT27I⊗	CGCGCAGG
	Jma19592II⊗	GRGCRAC
	KfII	GG/GWCCC
	Kpn156V⊗	CRTGATT
	Lmo370I⊗	AGCGCCG
	Lsp6406VI⊗	CRAGCAC
	Magl⊗	CRTTGAC(21/19)
	MauBl	CG/CGCGCG
	Mcr10l⊗	GAAGNNNNNCTC
	MkaDII⊗	GAGAYGT
	Mrel	CG/CCGGCG
	MspSC27II⊗	CCGCGAC
	Mtel	GCGC/NGCGC
	MtuHN878II⊗	CACGCAG
	NhaXl⊗	CAAGRAG
lacksquare		GC/GGCCGC
	NpeUS61II⊗	GATCGAC
$\overline{}$	Pacl	TTAAT/TAA
ī	Pal408l⊗	CCRTGAG
	Pasl	CC/CWGGG
	PfIPt14I⊗	RGCCCAC
	PfrJS12V⊗	GGCGGAG
	PinP23II⊗	CTRKCAG
	PliMl⊗	CGCCGAC
<u> </u>	Pmel	GTTT/AAAC
_	Ppil⊗	(7/12)GAACNNNNNCTC(13/8)
	PpiP13II⊗	CGCRGAC
	PpuMI	RG/GWCCY
-	Pse18267I⊗	RCCGAAG
	Psp0MII⊗	CGCCCAR(20/18)
_	PspXI	VC/TCGAGB
	Psrl	(7/12)GAACNNNNNNTAC(12/7)
	Pst145l⊗	CTAMRAG
	Pst273I⊗	GATCGAG
	Rba2021I⊗	CACGAGH
		CATCGAC(20/18)
	Rcel⊗ Pool⊙	
	Rpal⊗	GTYGGAG(11/9)
	RpaBl⊗	CCCGCAG(20/18)
	RpaB5l⊗	CGRGGAC(20/18)
	RpaTl⊗	GRTGGAG
	RspPBTS2III⊗	CTTCGAG

S	PECIFICITIES GR	EATER THAN 6 BASES (CONT.)
lack	RsrII	CG/GWCCG
A	Sapl	GCTCTTC(1/4)
\blacktriangle	Sbfl •	CCTGCA/GG
	SdeOSI⊗	(11/13)GACNNNNRTGA(12/10)
	SexAl	A/CCWGGT
▲	Sfil	GGCCNNNN/NGGCC
▲	SgrAl	CR/CCGGYG
	SgrDI	CG/TCGACG
	SmaUMH8I⊗	GCGAACB
	Sno506l⊗	GGCCGAG
	SpoDl⊗	GCGGRAG
	Srfl	GCCC/GGGC
	Sse8647I⊗	AG/GWCCT
	Ssp714II⊗	CGCAGCG
	SstE37I⊗	CGAAGAC(20/18)
	Sth20745III⊗	GGACGAC
lack	Swal	ATTT/AAAT
	TspARh3l⊗	GRACGAC
	UbaF9I⊗	TACNNNNNRTGT
	UbaF12I⊗	CTACNNNGTC
	UbaF13I⊗	GAGNNNNNNCTGG
	Vtu19109l⊗	CACRAYC

IN	TERRUPTED F	PALINDROMES
	Agsl	TTS/AA
_	Ahdl	GACNNN/NNGTC
A	Alel	CACNN/NNGTG
	Alfl⊗	(10/12)GCANNNNNNTGC(12/10)
▲	AlwNI	CAGNNN/CTG
	ApaBl⊗	GCANNNNN/TGC
A	ApeKI	G/CWGC
_	Avall	G/GWCC
	Bdal⊗	(10/12)TGANNNNNNTCA(12/10)
_	Bgll	GCCNNNN/NGGC
	Bisl	GC/NGC
_	Blpl	GC/TNAGC
	BIsI	GCN/GC
	BpII	(8/13)GAGNNNNNCTC(13/8)
	BsaBI	GATNN/NNATC
_	BsaJI	C/CNNGG
	BsaWI	W/CCGGW
_	BsiHKAI	GWGCW/C
A	BsII	CCNNNNN/NNGG
_	BstAPI	GCANNNN/NTGC
	BstEII •	G/GTNACC
_	BstNI	CC/WGG
_	BstXI	CCANNNNN/NTGG

IN	TERRUPTED PALIF	NDROMES (CONT.)
lack	Bsu36l	CC/TNAGG
	BthCl⊗	GCNG/C
_	Cac8I	GCN/NGC
	Cjul⊗	CAYNNNNRTG
_	Ddel	C/TNAG
	Dde51507I⊗	CCWGG
_	Dralll •	CACNNN/GTG
_	Drdl	GACNNNN/NNGTC
	EcoHI⊗	/CCSGG
_	EcoNI	CCTNN/NNNAGG
_	Eco0109I	RG/GNCCY
	Fall	(8/13)AAGNNNNNCTT(13/8)
	Fmul⊗	GGNC/C
_	Fnu4HI	GC/NGC
	Hael⊗	WGG/CCW
	HgiEll⊗	ACCNNNNNNGGT
_	Hinfl	G/ANTC
_	Hpy99I	CGWCG/
_	Hpy166II	GTN/NAC
_	Hpy188I	TCN/GA
_	Hpy188III	TC/NNGA
_	HpyCH4III	ACN/GT

IN	TERRUPTED PALIN	IDROMES (CONT.)
	HpyUM032XIII⊗	CYANNNNNNTRG
	Hsoll⊗	(8/14)CAYNNNNNRTG(14/8)
	KfII	GG/GWCCC
	MaeIII	/GTNAC
	MjalV⊗	GTNNAC
\blacktriangle	MsII	CAYNN/NNRTG
	Mtel	GCGC/NGCGC
	Mwol	GCNNNNN/NNGC
\blacktriangle	Ncil	CC/SGG
	Nhol⊗	GCWGC
\blacktriangle	NIaIV	GGN/NCC
	Pasl	CC/CWGGG
	Pcsl	WCGNNNN/NNNCGW
	Pfl8569l⊗	GCN/NGC
\blacktriangle	PfIFI	GACN/NNGTC
A	PfIMI	CCANNNN/NTGG
	PfoI	T/CCNGGA
	PpuMI	RG/GWCCY
A	PshAl	GACNN/NNGTC
	Psp03l⊗	GGWC/C
A	PspGI	/CCWGG
	Pssl⊗	RGGNC/CY

II	NTERRUPTED	PALINDROMES (CONT.)
\blacksquare	Rsrll	CG/GWCCG
A	Sau96I	G/GNCC
A	ScrFI	CC/NGG
	SetI	ASST/
A	SexAl	A/CCWGGT
_	Sfil	GGCCNNNN/NGGCC
	Sse8647I⊗	AG/GWCCT
A	Styl	C/CWWGG
A	StyD4I	/CCNGG
	Tatl	W/GTACW
	Taul	GCSG/C
A	Tfil	G/AWTC
A	Tsel	G/CWGC
A	Tsp45I	/GTSAC
A	TspRI	CASTGNN/
	Tssl⊗	GAGNNNCTC
A	Tth111I	GACN/NNGTC
	Unbl⊗	/GGNCC
	VpaK11Al⊗	/GGWCC
A	Xcml	CCANNNNN/NNNNTGG
_	Xmnl	GAANN/NNTTC

Multiple Recognition Sequences

TOOLS & RESOURCES

Visit the Tools & Resources tab at NEB.com to find:

 Access to our online tool NEBcutter, for help with restriction enzyme mapping

Single Letter Code:

Note:

Enzymes marked with a "^ are available from NEB. Enzymes marked with a "\otimes" are not currently commercially available.

M	ULTIPLE RECOGNI	TION SEQUENCES
	Aba6411II⊗	CRRTAAG
	AbaCIII⊗	CTATCAV
	AbaUMB2I⊗	YCCGSS
▲	Accl	GT/MKAC
	Aco12261II⊗	CCRGAG
	AcoY31II⊗	TAGCRAB
A	AfIIII	A/CRYGT
	Agsl	TTS/AA
	AhyRBAHI⊗	GCYYGAC
	AhyYL17l⊗	YAAMGAG
	AlwFl⊗	GAAAYNNNNNRTG
_	ApeKI	G/CWGC
▲	Apol	R/AATTY
	AquIV⊗	GRGGAAG(19/17)
	Arsl	(8/13)GACNNNNNNTTYG(11/6)
	Asp103l⊗	CGRAGGC
	AspBHI⊗	YSCNS(8/12)
	AspNIH4III⊗	AAGAACB
	Asp114pII⊗	AGCABCC
	Asu14238IV⊗	CGTRAC
	AteTl⊗	GGGRAG
	Aval	C/YCGRG

IMI	JLTIPLE RECOGNI	TION SEQUENCES (CONT.)
lack	Avall	G/GWCC
	Awo1030IV⊗	GCCRAG
_	Bael	(10/15)ACNNNNGTAYC(12/7)
A	BaeGI	GKGCM/C
_	Banl	G/GYRCC
A	Banll	GRGCY/C
	BanLl⊗	RTCAGG
	Bbr11l⊗	GGRCAG
	Bbr57III⊗	GTRAAYG
	BkrAM31DI⊗	RTTAAATM
	Ble402II⊗	GRAGCAG
	Bmgl⊗	GKGCCC
_	BsaAl	YAC/GTR
_	BsaHI	GR/CGYC
A	BsaWI	W/CCGGW
▲	BsiEl	CGRY/CG
A	BsiHKAI	GWGCW/C
_	BsoBI	C/YCGRG
_	Bsp1286I	GDGCH/C
_	BsrFI	R/CCGGY
_	BstNI	CC/WGG
A	BstYI	R/GATCY
_	Btgl	C/CRYGG
	CalB3II⊗	GRTTRAG
	Cba16038I⊗	CCTNAYNC
	Cbo67071IV⊗	GCRGAAG
	CchII⊗	GGARGA(11/9)
	Cch467III⊗	GNGAAAY
	Cco14983V⊗	GGGTDA
	Cco14983VI⊗	GCYGA
	Cdi81III⊗	GCMGAAG
	Cfupf3II⊗	GARCAG
	Cgl13032II⊗	ACGABGG
	Cje265V⊗	GKAAGC
	Cje54107III⊗	GKAAYC
	CjeFV⊗	GGRCA
	CjeNIII⊗	GKAAYG(19/17)
	CjeNV⊗	CCYGA
	Cjul⊗	CAYNNNNRTG
	Cjull⊗	CAYNNNNCTC
	Cly7489II⊗	AAAAGRG

M	ULTIPLE RECOGNIT	TION SEQUENCES (CONT.)
A	CviKI-1	RG/CY
	Dde51507I⊗	CCWGG
		Y/GGCCR
	Ecl35734I⊗	GAAAYTC
	Eco4465II⊗	GAAABCC
	Eco43896II⊗	CRARCAG
	EcoBLMcrX⊗	RCSRC(-3/-2)
	EcoE1140I⊗	ACCYAC
	EcoHI⊗	/CCSGG
	Eco57MI⊗	CTGRAG(16/14)
\blacktriangle	Eco0109I	RG/GNCCY
	Fail	YA/TR
	Fco1691IV⊗	GCVGAG
	FspAl	RTGC/GCAY
	FspPK15I⊗	GARGAAG
	Gdill⊗	CGGCCR(-5/-1)
	Hael⊗	WGG/CCW
\blacktriangle	Haell	RGCGC/Y
	HaelV⊗	(7/13)GAYNNNNNRTC(14/9)
	Hin4l⊗	(8/13)GAYNNNNNVTC(13/8)
\blacktriangle	HincII	GTY/RAC
A	Hpy99I	CGWCG/
	Hpy99XIV⊗	GGWTAA
	Hpy99XIV-mut1⊗	GGWCNA
	Hpy99XXII⊗	TCANNNNNTRG
	Hpy300XI⊗	CCTYNA
	HpyAXVI-mut1⊗	CRTTAA
	HpyAXVI-mut2⊗	CRTCNA
	HpyUM032XIII⊗	CYANNNNNNTRG
	HpyUM032XIII-mut1⊗	CYANNNNNNTTC
	Hsoll⊗	(8/14)CAYNNNNNRTG(14/8)
	Jma19592II⊗	GRGCRAC
	Jsp2502II⊗	GRNGAAT
	KfII	GG/GWCCC
	Kor51II⊗	RTCGAG
	Kpn156V⊗	CRTGATT
	KpnNH25III⊗	CTRGAG
	KpnNIH50I⊗	GCYAAG
	Lba2029III⊗	CYAAANG
	Lmo911II⊗	TAGRAG
	LlaG50l⊗ Lmo911ll⊗	CCGTKA TAGRAG

Multiple Recognition Sequences (continued)

MULTIPLE RECO	OGNITION SEQUENCES (CONT.)
Lpl1004ll⊗	AGGRAG
Lpnl⊗	RGC/GCY
▲ LpnPI	CCDG(10/14)
Lsp6406VI⊗	CRAGCAC
Maql⊗	CRTTGAC(21/19)
MkaDII⊗	GAGAYGT
Mmel	TCCRAC(20/18)
▲ MsII	CAYNN/NNRTG
MspA1I	CMG/CKG
Mspl7ll⊗	ACGRAG
MspJI	CNNR(9/13)
Ncil	CC/SGG
NhaXl⊗	CAAGRAG
Nhol+	GCWGC
NIi3877I⊗	CYCGR/G
NmeDl⊗	(12/7)RCCGGY(7/12)
Nspl	RCATG/Y
OspHL35III⊗	YAGGAG
Pal408l⊗	CCRTGAG
Pasl	CC/CWGGG
PcsI	WCGNNNN/NNNCGW
PflPt14l⊗	RGCCCAC
Pin17FIII⊗	GGYGAB
PinP23II⊗	CTRKCAG
PpiP13II⊗	CGCRGAC

M	ULTIPLE RECOGN	ITION SEQUENCES (CONT.)
A	PpuMI	RG/GWCCY
	Pse18267I⊗	RCCGAAG
	Psp03l⊗	GGWC/C
\blacktriangle	PspGI	/CCWGG
	Psp0MII⊗	CGCCCAR(20/18)
	PspPRI⊗	CCYCAG(15/13)
\blacktriangle	PspXI	VC/TCGAGB
	Pssl⊗	RGGNC/CY
	Pst145I⊗	CTAMRAG
	Pst14472I⊗	CNYACAC
	PsuGl⊗	BBCGD
	Rba2021I⊗	CACGAGH
	RdeGBIII⊗	(9/11)TGRYCA(11/9)
	Rlal⊗	VCW
	Rmu369III⊗	GGCYAC
	Rpal⊗	GTYGGAG(11/9)
	RpaB5l⊗	CGRGGAC(20/18)
	RpaTl⊗	GRTGGAG
\blacktriangle	RsrII	CG/GWCCG
	Sba460II⊗	GGNGAYG
	SdeAl⊗	CAGRAG(21/19)
	SdeOSI⊗	(11/13)GACNNNNRTGA(12/10)
	SenSARA26III⊗	ACRCAG
	SetI	ASST/
\blacktriangle	SexAl	A/CCWGGT

D/I	III TIDI E DECOCNI	TION SEQUENCES (CONT.)
		TION SEQUENCES (CONT.)
	SfcI	C/TRYAG
A	SgrAl	CR/CCGGYG
	SgrTl⊗	CCDS(10/14)
	SmaUMH8I⊗	GCGAACB
A	SmII	C/TYRAG
	SpoDl⊗	GCGGRAG
	Sse8647I⊗	AG/GWCCT
_	Styl	C/CWWGG
	SurP32all⊗	ACRGAG
	Tatl	W/GTACW
	Taul	GCSG/C
_	Tfil	G/AWTC
A	Tsel	G/CWGC
	Tsol⊗	TARCCA(11/9)
A	Tsp45I	/GTSAC
	TspARh3l⊗	GRACGAC
A	TspRI	CASTGNN/
	Tth111II⊗	CAARCA(11/9)
	UbaF9l⊗	TACNNNNNRTGT
	Van9116l⊗	CCKAAG
	Vdi96Ⅱ⊗	GNCYTAG
	VpaK11Al⊗	/GGWCC
	Vtu19109l⊗	CACRAYC
	Wvil⊗	CACRAG(21/19)

Nonpalindromic Recognition Sequences

Single Letter Code:

R = A or G K = G or T

Y = C or T S = C or GM = A or C W = A or TH = A or C or T V = A or C or G N = A or C or G or T B = C or G or T D = A or G or T

Enzymes marked with a "^" are available from NEB.

 HF (High-Fidelity) versions of these enzymes are available for simplified reactions and reduced star activity, at no additional cost. See page 301.

Enzymes marked with a " \otimes " are not currently commercially available.

N	ONPALINDROMIC	SEQUENCES
	Aarl	CACCTGC(4/8)
	Aba6411II⊗	CRRTAAG
	AbaB8342IV⊗	CATTAG
	AbaCIII⊗	CTATCAV
\blacktriangle	AbaSI	C(11/9)
	AbaUMB2I⊗	YCCGSS
	Acc65V⊗	GACGCA
	AceIII⊗	CAGCTC(7/11)
	AchA6III⊗	AGCCAG
\blacktriangle	Acil	CCGC(-3/-1)
	Aco12261II⊗	CCRGAG
	AcoY31II⊗	TAGCRAB
\blacktriangle	Acul	CTGAAG(16/14)
	Adh6U21I⊗	GAANCAG
	AhyRBAHI⊗	GCYYGAC
	AhyYL17I⊗	YAAMGAG
	Ajul	(7/12)GAANNNNNNNTTGG(11/6)
	Alol	(7/12)GAACNNNNNNTCC(12/7)
\blacktriangle	Alwl	GGATC(4/5)
	AlwFl⊗	GAAAYNNNNNRTG
	AmaCSI⊗	GCTCCA(11/9)
	ApyPl⊗	ATCGAC(20/18)
	AquII⊗	GCCGNAC(20/18)
	AquIII⊗	GAGGAG(20/18)
	AquIV⊗	GRGGAAG(19/17)
	Arsl	(8/13)GACNNNNNNTTYG(11/6)
	Asp103l⊗	CGRAGGC
	AspBHI⊗	YSCNS(8/12)
	AspDUT2V⊗	GNGCAAC
	AspJHL3II⊗	CGCCCAG
	AspNIH4III⊗	AAGAACB

N	ONPALINDROMI	C SEQUENCES (CON'T)
	AspSLV7III⊗	GTCTCA
	Asp114pII⊗	AGCABCC
	Asu14238IV⊗	CGTRAC
	AteTI⊗	GGGRAG
	Awo1030IV⊗	GCCRAG
\blacktriangle	Bael	(10/15)ACNNNNGTAYC(12/7)
	Bag18758I⊗	CCCGAG
	BanLl⊗	RTCAGG
	Barl	(7/12)GAAGNNNNNNTAC(12/7)
	Bbr11l⊗	GGRCAG
	Bbr52II⊗	GGCGAG
	Bbr57III⊗	GTRAAYG
\blacktriangle	Bbsl	GAAGAC(2/6)
\blacktriangle	Bbvl	GCAGC(8/12)
\blacktriangle	BbvCl	CCTCAGC(-5/-2)
\blacktriangle	Bccl	CCATC(4/5)
	Bce3081I⊗	TAGGAG
▲	BceAl	ACGGC(12/14)
	BceSIV⊗	(7/5)GCAGC(9/11)
	Bcefl⊗	ACGGC(12/13)
\blacktriangle	Bcgl	(10/12)CGANNNNNNTGC(12/10)
A	BciVI	GTATCC(6/5)
\blacktriangle	BcoDI	GTCTC(1/5)
	BfaSII⊗	GANGGAG
\blacktriangle	BfuAl	ACCTGC(4/8)
	BkrAM31DI⊗	RTTAAATM
	Ble402II⊗	GRAGCAG
	BloAll⊗	GAGGAC
	BmeDl⊗	C(2/0)
	Bmgl⊗	GKGCCC
\blacktriangle	BmgBl	CACGTC(-3/-3)

1	IONPALINDROMI	C SEQUENCES (CON'T)
A	Bmrl	ACTGGG(5/4)
▲	Bpml	CTGGAG(16/14)
	Bpu10I	CCTNAGC(-5/-2)
▲	BpuEl	CTTGAG(16/14)
▲	Bsal •	GGTCTC(1/5)
▲	BsaXI	(9/12)ACNNNNNCTCC(10/7)
	Bsbl⊗	CAACAC(21/19)
	BscAl⊗	GCATC(4/6)
	BscGl⊗	CCCGT
	BseMII	CTCAG(10/8)
	BseRI	GAGGAG(10/8)
▲	BseYl	CCCAGC(-5/-1)
	Bsgl	GTGCAG(16/14)
▲	Bsml	GAATGC(1/-1)
	BsmAl	GTCTC(1/5)
▲	BsmBl	CGTCTC(1/5)
▲	BsmFI	GGGAC(10/14)
	Bsp24l⊗	(8/13)GACNNNNNNTGG(12/7)
	Bsp460III⊗	CGCGCAG
	Bsp3004IV⊗	CCGCAT
▲	BspCNI	CTCAG(9/7)
	BspD6l⊗	GAGTC(4/6)
	BspGl⊗	CTGGAC
A	BspMI	ACCTGC(4/8)
	BspNCl⊗	CCAGA
A	BspQI	GCTCTTC(1/4)
A	Bsrl	ACTGG(1/-1)
A	BsrBI	CCGCTC(-3/-3)
A	BsrDI	GCAATG(2/0)
A	BssSI	CACGAG(-5/-1)
A	BtgZl	GCGATG(10/14)

M	ONDALINDROMI	C SEQUENCES (CON'T)
IN A		GCAGTG(2/0)
_	BtsIMutl	CAGTG(2/0)
	BtsCl	GGATG(2/0)
	Cal14237I⊗	GGTTAG
	CalB3II⊗	GRTTRAG
	Cau10061II⊗	GTTAAT
	Cba13II⊗ Cba16038I⊗	AGGAAT CCTNAYNC
	Cbo67071IV⊗	GCRGAAG
	CchII⊗	GGARGA(11/9)
	CchIII⊗	CCCAAG(20/18)
	Cch467III⊗	GNGAAAY
	Cco14983V⊗	GGGTDA
	Cco14983VI⊗ CcrNAIII⊗	GCYGA CGACCAG
	Cdil⊗	CATCG(-1/-1)
	Cdi81III⊗	GCMGAAG
	Cdi11397l⊗	GCGCAG
	Cdpl⊗	GCGGAG(20/18)
	Cdu23823II⊗	GTGAAG
	Cfupf3II⊗ Cgl13032I⊗	GARCAG GGCGCA
	Cgl13032ll⊗	ACGABGG
	Cjel⊗	(8/14)CCANNNNNNGT(15/9)
	Cje265V⊗	GKAAGC
	Cje54107III⊗	GKAAYC
	CjeFIII⊗	GCAAGG
	CjeFV⊗ CjeNII⊗	GGRCA GAGNNNNNGT
	CjeNIII⊗	GKAAYG(19/17)
	CjeNV⊗	CCYGA
	CjePl⊗	(7/13)CCANNNNNNNTC(14/8)
	CjeP659IV⊗	CACNNNNNNGAA
	Cjull⊗	CAYNNNNCTC
	Cla11845III⊗ Cly7489II⊗	GCGAA AAAAGRG
	Cma23826l⊗	CGGAAG
	Csp2014I⊗	GGAGGC
_	CspCl	(11/13)CAANNNNNGTGG(12/10)
	CstMI⊗	AAGGAG(20/18)
	DraRl⊗ Drdll⊗	CAAGNAC(20/18) GAACCA
	Earl	CTCTTC(1/4)
	Ecil	GGCGGA(11/9)
	Ecl234l⊗	CGGNAAG
	Ecl35734I⊗	GAAAYTC
	Eco4465II⊗ Eco43896II⊗	GAAABCC CRARCAG
	EcoBLMcrX⊗	RCSRC(-3/-2)
	EcoE1140I⊗	ACCYAC
	Eco57MI⊗	CTGRAG(16/14)
	EcoMVII⊗	CANCATC
	EcoNIH6II⊗	ATGAAG
	Eli8509II⊗	CCGGAG GACCAC
<u> </u>	EsaSSI⊗ Esp3I	CGTCTC(1/5)
	Esp3007I⊗	CAGAAG
	Exi27195I⊗	GCCGAC
_	Faul	CCCGC(4/6)
	Fco1691IV⊗	GCVGAG
	Finl⊗	GGGAC
	Fokl FspEl	GGATG(9/13) CC(12/16)
	FspPK15l⊗	GARGAAG
	FtnUV⊗	GAAACA
	GauT27I⊗	CGCGCAGG
	Gba708II⊗	ATGCAC
	Gdill⊗	CGGCCR(-5/-1)
	Gsal	CCCAGC(-1/-5)
	HaelV⊗	(7/13)GAYNNNNNRTC(14/9)

IV	ONPALINDROMIC	C SEQUENCES (CON'T)
	Haull⊗	TGGCCANNNNNNNNNNN/
	Hball⊗	GCCCAG
	HdeNY26I⊗	CGANNNNNTCC
	HdeZA17l⊗	GCANNNNNTCC
A	Hgal	GACGC(5/10)
	Hin4l⊗	(8/13)GAYNNNNNVTC(13/8)
A	Hphl	GGTGA(8/7)
	Hpy99XIII⊗	GCCTA
	Hpy99XIV⊗	GGWTAA
	Hpy99XIV-mut1⊗	GGWCNA
	Hpy99XXII⊗	TCANNNNNTRG
	Hpy300XI⊗	CCTYNA CCTTC(6/5)
•	HpyAVIV	GCGTA
	HpyAXIV⊗ HpyAXVI-mut1⊗	CRTTAA
	HpyAXVI-mut2⊗	CRTCNA
	HpyUM032XIII-mut1⊗	CYANNNNNNTTC
	HpyUM032XIV⊗	GAAAG
	HpyUM037X⊗	GTGGNAG, TNGGNAG
	Jma19592I⊗	GTATNAC
	Jma19592I⊗	GRGCRAC
	Jsp2502II⊗	GRNGAAT
	Kor51II⊗	RTCGAG
	Kpn156V⊗	CRTGATT
	KpnNH25III⊗	CTRGAG
	KpnNlH30lll⊗	GTTCNAC
	KpnNlH50l⊗	GCYAAG
	Lba2029III⊗	CYAAANG
	Lde4408II⊗	ACAAAG
	LlaG50l⊗	CCGTKA
	LmnI	GCTCC(1/-1)
	Lmo370l⊗	AGCGCCG
	Lmo911II⊗	TAGRAG
	Lpl1004ll⊗	AGGRAG
A	LpnPI	CCDG(10/14)
	Lra68l⊗	GTTCNAG
	LsaDS4I⊗	TGGAAT
	Lsp48III⊗	AGCACC
	Lsp6406VI⊗	CRAGCAC
	Magl⊗	CRTTGAC(21/19)
	Mba11I⊗	AGGCGA
\blacktriangle	Mboll	GAAGA(8/7)
	Mcr10l⊗	GAAGNNNNNCTC
	MkaDII⊗	GAGAYGT
\blacktriangle	Mlyl	GAGTC(5/5)
\blacktriangle	Mmel	TCCRAC(20/18)
\blacktriangle	MnII	CCTC(7/6)
	Mspl7II⊗	ACGRAG
\blacktriangle	MspJI	CNNR(9/13)
	MspSC27II⊗	CCGCGAC
	MtuHN878II⊗	CACGCAG
	NaI45188II⊗	ACCAGC
	Nbr128II⊗	ACCGAC
	NgoAVII⊗	GCCGC(7/7)
	NgoAVIII⊗	(12/14)GACNNNNNTGA(13/11)
	NhaXl⊗	CAAGRAG
	NIaCl⊗	CATCAC(19/17)
\blacktriangle	NmeAIII	GCCGAG(21/19)
	NpeUS61II⊗	GATCGAC
	OspHL35III⊗	YAGGAG
	PacIII⊗	GTAATC
	Pac19842II⊗	CCTTGA
	Pal408l⊗	CCRTGAG
	Pba2294I⊗	GTAAG
	Pcall⊗	GACGAG
	Pcr308II⊗	CCAAAG
	Pdi8503III⊗	CCGGNAG
	Pdu1735l⊗	CACCAC
		CACCAC GCAGT TCGTAG

N		C SEQUENCES (CON'T)
	PfIPt14I⊗	RGCCCAC
	PfrJS12IV⊗	TANAAG
	PfrJS12V⊗	GGCGGAG CTTCNAC
	PfrJS15III⊗ Pin17FIII⊗	GGYGAB
	PinP23II⊗	CTRKCAG
	PinP59III⊗	GAAGNAG
	PlaDl⊗	CATCAG(21/19)
A	Plel	GAGTC(4/5)
	PliMl⊗	CGCCGAC
	Ppil⊗	(7/12)GAACNNNNNCTC(13/8)
	PpiP13II⊗	CGCRGAC
	Pse18267I⊗ Psp0357II⊗	RCCGAAG GCGAAG
	Psp0MII⊗	CGCCCAR(20/18)
	PspPRI⊗	CCYCAG(15/13)
	Psrl	(7/12)GAACNNNNNNTAC(12/7)
	Pst145l⊗	CTAMRAG
	Pst273I⊗	GATCGAG
	Pst14472I⊗	CNYACAC
	PsuGl⊗	BBCGD
	Rba2021I⊗	CACGAGH
	Rcel⊗ RdeGBl⊗	CATCGAC(20/18) CCGCAG
	RdeGBII⊗	ACCCAG(20/18)
	RfIFIII®	CGCCAG
	Rlal⊗	VCW
	RIaII⊗	ACACAG(20/18)
	RleAl⊗	CCCACA(12/9)
	Rmu369III⊗	GGCYAC
	Rpal⊗	GTYGGAG(11/9)
	RpaBl⊗	CCCGCAG(20/18)
	RpaB5I⊗ RpaTI⊗	CGRGGAC(20/18) GRTGGAG
	Rsp008IV⊗	ACGCAG
	Rsp008V⊗	GCCCAT
	RspPBTS2III⊗	CTTCGAG
	Rtr1953I⊗	TGANNNNNTGA
	Saf8902III⊗	CAATNAG
A	Sapl	GCTCTTC(1/4)
	Sba460II⊗	GGNGAYG
	Sbo46l⊗ ScoDS2ll⊗	TGAAC GCTAAT
	SdeAl⊗	CAGRAG(21/19)
	SdeOSI®	(11/13)GACNNNNRTGA(12/10)
	Sen17963III⊗	CCAAAC
	SenA1673III⊗	GNGGCAG
	SenSARA26III⊗	ACRCAG
	SenTFIV⊗	GATCAG
A	SfaNI	GCATC(5/9)
	Sgel	CNNGNNNNNNNNN/
	SgrTl⊗ Siml⊗	CCDS(10/14) GGGTC(-3/0)
	SmaUMH5I⊗	CTTGAC
	SmaUMH8I⊗	GCGAACB
	Sno506l⊗	GGCCGAG
	SpnRII⊗	TCGAG
	SpoDl⊗	GCGGRAG
	Ssp714II⊗	CGCAGCG
	Ssp6803IV⊗	GAAGGC
	SspD5l⊗ SstE37l⊗	GGTGA(8/8) CGAAGAC(20/18)
	Sth132l⊗	CCCG(4/8)
	Sth20745III⊗	GGACGAC
	SthSt3II⊗	GAAGT
	Stsl⊗	GGATG(10/14)
	SurP32all⊗	ACRGAG
	Taqll	GACCGA(11/9)
	TaqIII⊗	CACCCA(11/9)
	Tsol⊗	TARCCA(11/9)

APPENDIX

Nonpalindromic Recognition Sequences (continued)

NONPALINDROMIC SEQUENCES (CON'T)		
TspARh3l⊗	GRACGAC	
TspDTI	ATGAA(11/9)	
TspGWI	ACGGA(11/9)	
Tstl⊗	(8/13)CACNNNNNNTCC(12/7)	
Tsul⊗	GCGAC	
Tth111II⊗	CAARCA(11/9)	
UbaF9l⊗	TACNNNNNRTGT	
UbaF11I⊗	TCGTA	
UbaF12l⊗	CTACNNNGTC	
UbaF13I⊗	GAGNNNNNNCTGG	

NONPALINDRO	MIC SEQUENCES (CON'T)
UbaF14l⊗	CCANNNNTCG
UbaPl⊗	CGAACG
Van9116l⊗	CCKAAG
Vdi96II⊗	GNCYTAG
Vtu19109I⊗	CACRAYC
Wvil⊗	CACRAG(21/19)
Xca85IV⊗	TACGAG
Ykrl⊗	C(10/9)
Yps3606l⊗	CGGAAG

Alphabetized List of NEB Recognition Sequences

All restriction enzyme recognition specificities and recommended enzymes available from New England Biolabs are listed below. For enzymes that recognize non-palindromic sequences, the complementary sequence of each strand is listed. For example, CCTC(7/6) and (6/7)GAGG both represent an MnII site. New entries are listed in **bold** type.

AA/CGTT AcII A/AGCTT HindIII-HF AAT/ATT SspI-HF /AATT MluCl A/CATGT Pcil A/CCGGT Agel-HF ACCTGC(4/8) BfuAl ACCTGC(4/8) BspMI A/CCWGGT SexAl Mlul-HF A/CGCGT ACGGC(12/14) BceAl A/CGT HpyCH4IV ACN/GT HpyCH4III (10/15)ACNNNNGTAYC(12/7) Bael (9/12)ACNNNNNCTCC(10/7) BsaXI A/CRYGT AfIIII A/CTAGT Spel-HF ACTGG(1/-1) Bsrl ACTGGG(5/4) Bmrl A/GATCT BgIII AGC/GCT Afel AG/CT Alul AGG/CCT Stul Scal-HF AGT/ACT AT/CGAT **BspDI** AT/CGAT Clal ATGCA/T Nsil-HF AT/TAAT Asel ATTT/AAAT Swal C(11/9) AbaSI (11/13)CAANNNNNGTGG(12/10) CspCI C/AATTG Mfel-HF CACGAG(-5/-1) BssSI-v2 CACGTC(-3/-3) BmgBI CAC/GTG PmII CACNNN/GTG DrallI-HF CACNN/NNGTG Alel-v2 (0/2)CACTG BtsIMutl (0/2)CACTGC BtsI-v2 CAG/CTG PvuII-HF CAGNNN/CTG AlwNI CAGTG(2/0) BtslMutl CASTGNN/ TspRI CA/TATG Ndel (0/2)CATCC BtsCI (13/9)CATCC Fokl (14/10)CATCGC BtgZI C/ATG CviAII

/CATG	Fatl
CATG/	NIaIII
(0/2)CATTGC	BsrDI
CAYNN/NNRTG	MsII
CC(12/16)	FspEl
(10/12)CCACNNNNNTTG	CspCl
(-1/1)CCAGT	Bsrl
CCANNNNN/NNNNTGG	Xcml
CCANNNNN/NTGG	BstXI
CCANNNN/NTGG	PfIMI
CCATC(4/5)	Bccl
C/CATGG	Ncol-HF
CCCAGC(-5/-1)	BseYI
(4/5)CCCAGT	Bmrl
CCCGC(4/6)	Faul
CCC/GGG	Smal
C/CCGGG	TspMI
C/CCGGG	Xmal
CCDG(10/14)	LpnPI
CCGC(-3/-1)	Acil
CCGC/GG	SacII
CCGCTC(-3/-3)	BsrBI
C/CGG	Hpall
C/CGG	Mspl
CC/NGG	ScrFI
/CCNGG	StyD4I
C/CNNGG	BsaJI
CCNNNNN/NNGG	BsII
C/CRYGG	Btgl
CC/SGG	Ncil
C/CTAGG	AvrII
CCTC(7/6)	MnII
CCTCAGC(-5/-2)	BbvCI
CCTGCA/GG	SbfI-HF
CCTNAGC(-5/-2)	Bpu10I
CC/TNAGG	Bsu36I
CCTNN/NNNAGG	EcoNI
CCTTC(6/5)	HpyAV
CC/WGG	BstNI
/CCWGG	PspGI
C/CWWGG	Styl-HF
(10/12)CGANNNNNNTGC(12/10)	Bcgl
CGAT/CG	Pvul-HF
CG/CG	BstUI
C/GGCCG	Eagl-HF
CG/GWCCG	RsrII
CGRY/CG	BsiEl
C/GTACG	BsiWl

All recognition sequences are written 5° to 3° using the single letter code nomenclature with the point of cleavage indicated by a "/".

Numbers in parentheses indicate point of cleaveage for non-palindromic enzymes. For example, GGTCTC(1/5) indicates cleavage at: 5´...GGTCTCN/...3´ 3´...CCAGAGNNNN/...5´

CGTCTC(1/5)	BsmBI
CGTCTC(1/5)	Esp3I
CGWCG/	Hpy99I
(14/10)CHGG	LpnPl
CMG/CKG	MspA1I
CNNR(9/13)	MspJI
CR/CCGGYG	SgrAl
C/TAG	Bfal
(14/16)CTCAAG	BpuEl
CTCAG(9/7)	BspCNI
(14/16)CTCCAG	Bpml
(8/10)CTCCTC	BseRI
C/TCGAG	PaeR7I
C/TCGAG	Xhol
(19/21)CTCGGC	NmeAIII
(-1/-5)CTCGTG	BssSI-v2
CTCTTC(1/4)	Earl
CTGAAG(16/14)	Acul
(7/9)CTGAG	BspCNI
(14/16)CTGCAC	Bsgl
CTGCA/G	PstI-HF
CTGGAG(16/14)	Bpml
C/TNAG	Ddel
C/TRYAG	SfcI
C/TTAAG	AfIII
(14/16)CTTCAG	Acul
CTTGAG(16/14)	BpuEl
C/TYRAG	Smll
C/YCGRG	Aval
C/YCGRG	BsoBl
(9/11)G	AbaSI
GAAGA(8/7)	Mboll
GAAGAC(2/6)	BbsI-HF
(4/1)GAAGAG	Earl
(4/1)GAAGAGC	BspQI
(4/1)GAAGAGC	Sapl
(5/6)GAAGG	НруАV
GAANN/NNTTC	Xmnl
GAATGC(1/-1)	Bsml
G/AATTC	EcoRI-HF
GACGC(5/10) GACGT/C	Hgal AatII
GAC/GTC	Zral
(-3/-3)GACGTG	BmgBl
GACN/NNGTC	PfIFI THA 1 1 1
GACN/NNGTC	Tth1111
GACNN/NNGTC	PshAl
GACNNN/NNGTC	Ahdl

GACNNNN/NNGTC	Drdl
(5/5)GACTC	Mlyl
(5/4)GACTC	Plel
(5/1)GAGAC	BcoDI
(5/1)GAGAC	BsmAl
(5/1)GAGACC	Bsal
(5/1)GAGACG	BsmBl
(5/1)GAGACG	Esp3I
(-3/-3)GAGCGG	BsrBl
GAG/CTC	Eco53kl
GAGCT/C	SacI-HF
(6/7)GAGG	MnII
GAGGAG(10/8)	BseRI
GAGTC(5/5)	Mlyl
GAGTC(4/5)	Plel
G/ANTC	Hinfl
GAT/ATC	EcoRV-HF
GA/TC	Dpnl
/GATC	DpnII
/GATC	Mbol
/GATC	Sau3Al
(5/4)GATCC	Alwl
(9/5)GATGC	SfaNI
(5/4)GATGG	Bccl
GATNN/NNATC	BsaBl
G/AWTC	Tfil
GCAATG(2/0)	BsrDI
GCAGC(8/12)	Bbvl
(8/4)GCAGGT	BfuAl
(8/4)GCAGGT	BspMI
GCAGTG(2/0)	BtsI-v2
(10/12)GCANNNNNNTCG	Bcgl
GCANNNN/NTGC	BstAPI
GCATC(5/9)	SfaNI
GCATG/C	SphI-HF
(-1/1)GCATTC	Bsml
GCCC/GGGC	Srfl
GCCGAG(21/19)	NmeAIII
GCC/GGC	Nael
G/CCGGC	NgoMIV
(14/12)GCCGT	BceAl
GCCNNNN/NGGC	Bgll
GCGAT/CGC	AsiSI
GCGATG(10/14)	BtgZI
GCG/C	Hhal
G/CGC	HinP1I
G/CGCGC	BssHII
(-1/-3)GCGG	Acil
(1/ 0)0000	71011

Alphabetized List of NEB Recognition Sequences (continued)

GC/GGCCGC	NotI-HF
(6/4)GCGGG	Faul
(10/5)GCGTC	Hgal
GC/NGC	Fnu4HI
GCN/NGC	Cac8I
GCNNNNN/NNGC	Mwol
GCTAG/C	BmtI-HF
G/CTAGC	Nhel-HF
GCTCTTC(1/4)	BspQI
GCTCTTC(1/4)	Sapl
(-2/-5)GCTGAGG	BbvCl
(12/8)GCTGC	Bbvl
(-1/-5)GCTGGG	BseYI
GC/TNAGC	Blpl
(-2/-5)GCTNAGG	Bpu10I
G/CWGC	ApeKI
G/CWGC	Tsel
GDGCH/C	Bsp1286I
(16/12)GG	FspEl
(7/10)GGAGNNNNNGT	BsaXI
(5/6)GGATAC	BciVI
GGATC(4/5)	Alwl
G/GATCC	BamHI-HF
GGATG(2/0)	BtsCl
GGATG(9/13)	Fokl
GG/CC	Haelll

GGCCGG/CC	Fsel
GGCCNNNN/NGGCC	Sfil
G/GCGCC	Kasl
GG/CGCC	Narl
GGCGC/C	PluTl
GGC/GCC	Sfol
GG/CGCGCC	Ascl
GGCGGA(11/9)	Ecil
GGGAC(10/14)	BsmFI
GGGCC/C	Apal
G/GGCCC	Psp0MI
G/GNCC	Sau96I
GGN/NCC	NIaIV
G/GTACC	Acc65I
GGTAC/C	KpnI-HF
GGTCTC(1/5)	Bsal-HFv2
GGTGA(8/7)	Hphl
G/GTNACC	BstEII-HF
G/GWCC	Avall
G/GYRCC	Banl
GKGCM/C	BaeGI
GR/CGYC	BsaHI
GRGCY/C	Banll
(7/12)GRTACNNNNGT	Bael
G/TAC	CviQI
GT/AC	Rsal

GTA/TAC	BstZ17I
GTATCC(6/5)	BciVI
(14/10)GTCCC	BsmFI
G/TCGAC	Sall-HF
GTCTC(1/5)	BcoDI
GTCTC(1/5)	BsmAl
(6/2)GTCTTC	BbsI-HF
G/TGCAC	ApaLI
GTGCAG(16/14)	Bsgl
GT/MKAC	Accl
GTN/NAC	Hpy166II
/GTSAC	Tsp45I
GTT/AAC	Hpal
GTTT/AAAC	Pmel
(18/20)GTYGGA	Mmel
GTY/RAC	HincII
GWGCW/C	BsiHKAI
R/AATTY	Apol
RCATG/Y	Nspl
R/CCGGY	BsrFI-v2
R/GATCY	BstYl
RGCGC/Y	Haell
RG/CY	CviKI-1
RG/GNCCY	Eco0109I
RG/GWCCY	PpuMI
TAC/GTA	SnaBl

(7/8)TCACC	Hphl
T/CATGA	BspHI
(9/11)TCCGCC	Ecil
T/CCGGA	BspEl
TCCRAC(20/18)	Mmel
T/CGA	Taqαl
TCG/CGA	Nrul-HF
TCN/GA	Hpy188I
TC/NNGA	Hpy188III
T/CTAGA	Xbal
(7/8)TCTTC	Mboll
T/GATCA	BcII
TG/CA	HpyCH4V
TGC/GCA	Fspl
TGG/CCA	Mscl
T/GTACA	BsrGI-HF
T/TAA	Msel
TTAAT/TAA	Pacl
TTA/TAA	Psil
TT/CGAA	BstBI
TTT/AAA	Dral
VC/TCGAGB	PspXI
W/CCGGW	BsaWl
YAC/GTR	BsaAl
Y/GGCCR	Eael
(13/9)YNNG	MspJI



NEB is the primary sponsor of the LabCentral
Learning Lab located in Kendall Square in
Cambridge, Massachusetts. We use this space
to provide hands-on training for our sales team
members, as well as our customers in the Boston
area. We also share the space with BioBuilder®, a
STEM education program focused on experiential
learning. Pictured here are several NEB members
participating in hands-on training.

ATTENDIX

Isoschizomers

Restriction enzymes that recognize the same sequence are isoschizomers. The first example discovered is called a prototype, and all subsequent enzymes that recognize the same sequence are isoschizomers of the prototype. The list below contains isoschizomers for commercially-available restriction endonucleases. It also specifies which isoschizomer is available from New England Biolabs.

All recognition sequences are written 5´ to 3´ using the single letter code nomenclature with the point of cleavage indicated by a "/".

Numbers in parentheses indicate point of cleavage for non-palindromic enzymes. For example, GGTCTC(1/5) indicates cleavage at: 5′...GGTCTCN/...3′

3´...CCAGAGNNNNN/...5´

Isoschizomers with alternative cleavage sites (neoschizomers) are indicated with a "^". Enzymes that are not currently commercially available are indicated with a " \otimes ". For more information on isoschizomers, visit **REBASE.neb.com**

Neoschizomers are a subset of isoschizomers that recognize the same sequence, but cleave at different positions from the prototype. Thus, AatII (recognition sequence: GACGT↓C) and ZraI (recognition sequence: GACGT↓C) are neoschizomers of one another, while HpaII (recognition sequence: C↓CGG) and MspI (recognition sequence: C↓CGG) are isoschizomers. Analogous designations are not appropriate for methyltransferases, where the differences between enzymes are not so easily defined and usually have not been well characterized.

FAITWAR	OFOURNOE	NEB	NED #	OFOLIENOF	OTHER ISOSOHITOMERS
ENZYME	SEQUENCE	ENZYME	NEB#	SEQUENCE	OTHER ISOSCHIZOMERS
Aanl	TTA/TAA	Psil	R0657	TTA/TAA	Psil
Aarl	CACCTGC(4/8)	D !!	Doroo	0.4.04.14.14.14.10.7.0	D D
Aasl	GACNNNN/NNGTC	Drdl	R0530	GACNNNN/NNGTC	Drdl, DseDl
AatII	GACGT/C	AatII	R0117	GACGT/C	Zral^
_		Zral^	R0659	GAC/GTC	
Aba6411II ⊗	CRRTAAG				
AbaB8342IV ⊗	CATTAG				
AbaCIII ⊗	CTATCAV				
AbaSI	C(11/9)	AbaSI	R0665	C(11/9)	
AbaUMB2I ⊗	YCCGSS				
Absl	CC/TCGAGG				
AccI	GT/MKAC	Accl	R0161	GT/MKAC	FbII, Xmil
AccII	CG/CG	BstUI	R0518	CG/CG	Bsh1236I, BspFNI, BstFNI, BstUI, MvnI
AccIII	T/CCGGA	BspEl	R0540	T/CCGGA	Aor13HI, BseAl, Bsp13I, BspEI, Kpn2I, MroI
Acc16I	TGC/GCA	Fspl	R0135	TGC/GCA	Fspl, Nsbl
Acc36I	ACCTGC(4/8)	BfuAl	R0701	ACCTGC(4/8)	BfuAl, BspMl, Bvel
		BspMI	R0502	ACCTGC(4/8)	
Acc65I	G/GTACC	Acc65I	R0599	G/GTACC	Asp718I, KpnI^, KpnI-HF^
		KpnI-HF^	R3142	GGTAC/C	
Acc65V ⊗	GACGCA				
AccB1I	G/GYRCC	Banl	R0118	G/GYRCC	Banl, BshNl, BspT107l
AccB7I	CCANNNN/NTGG	PfIMI	R0509	CCANNNN/NTGG	PfIMI, Van91I
AccBSI	CCGCTC(-3/-3)	BsrBl	R0102	CCGCTC(-3/-3)	BsrBl, Mbil
AceIII ⊗	CAGCTC(7/11)			(. ,	,
AchA6III ⊗	AGCCAG				
Acil	CCGC(-3/-1)	Acil	R0551	CCGC(-3/-1)	BspACI, Ssil
AcII	AA/CGTT	AcII	R0598	AA/CGTT	Psp1406l
AcIWI	GGATC(4/5)	Alwl	R0513	GGATC(4/5)	Alwl, BspPI
Acol	Y/GGCCR	Eael	R0508	Y/GGCCR	Eael
Aco12261Ⅱ ⊗	CCRGAG	Laui	110000	1/440011	Laui
AcoY31II ⊗	TAGCRAB				
Acsl	R/AATTY	Apol-HF	R3566	R/AATTY	Apol, Apol-HF, Xapl,
		Apol-Fir	R0641		
Acul	CTGAAG(16/14)		R0532	CTGAAG(16/14)	Eco57I
Acvi	CAC/GTG	PmII		CAC/GTG	BbrPI, Eco72I, PmaCI, PmII, PspCI
Acyl	GR/CGYC	BsaHl	R0556	GR/CGYC	BsaHl, BssNl, BstACl, Hin1l, Hsp92l
Adel	CACNNN/GTG	DrallI-HF	R3510	CACNNN/GTG	Dralli, Dralli-HF
Adh6U21I⊗	GAANCAG	0 :014	D0000	0.510	0. 014 0 :014 D. I.D. 1114
Afal	GT/AC	CviQI^	R0639	G/TAC	Csp61 [^] , CviQI [^] , Rsal, RsaNI [^]
		Rsal	R0167	GT/AC	
Afel	AGC/GCT	Afel	R0652	AGC/GCT	Aor51HI, Eco47III
Afil	CCNNNNN/NNGG	BsII	R0555	CCNNNNN/NNGG	Bsc4l, BseLl, BsII
AfIII	C/TTAAG	AfIII	R0520	C/TTAAG	Bfrl, BspTl, BstAFl, MspCl, Vha464l
AfIIII	A/CRYGT	AfIIII	R0541	A/CRYGT	
Agel	A/CCGGT	Agel-HF	R3552	A/CCGGT	Agel, Agel-HF, AsiGl, BshTl, CspAl, PinAl
Agsl	TTS/AA				
AhaIII ⊗	TTT/AAA	Dral	R0129	TTT/AAA	Dral
Ahdl	GACNNN/NNGTC	Ahdl	R0584	GACNNN/NNGTC	BmeRI, Dril, Eam1105I
Ahll	A/CTAGT	Spel-HF	R3133	A/CTAGT	Bcul, Spel, Spel-HF
AhyRBAHI⊗	GCYYGAC				
AhyYL17I⊗	YAAMGAG				
Ajil	CACGTC(-3/-3)	BmgBl	R0628	CACGTC(-3/-3)	BmgBl, Btrl
Ajnl	/CCWGG	BstNI^	R0168	CC/WGG	BciT130I^, BseBI^, BstNI^, Bst2UI^, EcoRII, MvaI^, Psp6I, PspGI
		PspGI	R0611	/CCWGG	
Ajul	(7/12)GAANNNNNNNTTGG(11/6)				
**	CACNN/NNGTG	Alel	R0634	CACNN/NNGTG	Olil

ENZYME	SEQUENCE	NEB Enzyme	NEB #	SEQUENCE	OTHER ISOSCHIZOMERS
Alfl ⊗	(10/12)GCANNNNNTGC(12/10)	ENETHIE		OEGOENOE'	- THE THE STORY OF
Alol	(7/12)GAACNNNNNNTCC(12/7)				
Alul	AG/CT	Alul	R0137	AG/CT	AluBl
					Alul
AluBl	AG/CT	Alul	R0137	AG/CT	
Alwi	GGATC(4/5)	Alwl	R0513	GGATC(4/5)	AcIWI, BspPI
Alw21I	GWGCW/C	BsiHKAI	R0570	GWGCW/C	Bbv12I, BsiHKAI
Alw26I	GTCTC(1/5)	BcoDI	R0542	GTCTC(1/5)	BcoDI, BsmAI, BstMAI
		BsmAl	R0529	GTCTC(1/5)	
Alw44I	G/TGCAC	ApaLI	R0507	G/TGCAC	ApaLI, Vnel
AlwFl⊗	GAAAYNNNNNRTG				
AlwNI	CAGNNN/CTG	AlwNI	R0514	CAGNNN/CTG	Cail, PstNI
Ama87I	C/YCGRG	Aval	R0152	C/YCGRG	Aval, BmeT110I, BsiHKCI, BsoBI, Eco88I
		BsoBI	R0586	C/YCGRG	
AmaCSI⊗	GCTCCA(11/9)				
Aor13HI	T/CCGGA	BspEl	R0540	T/CCGGA	Accill, BseAl, Bsp13l, BspEl, Kpn2l, Mrol
Aor51HI	AGC/GCT	Afel	R0652	AGC/GCT	Afel, Eco47III
		Alti	NU032	Add/dd1	Alei, Eco47III
Aoxl	/GGCC		D0444	0000010	D 400H D ONH
Apal	GGGCC/C	Apal	R0114	GGGCC/C	Bsp120I [^] , Psp0MI [^]
		Psp0MI [^]	R0653	G/GGCCC	
ApaBI⊗	GCANNNN/TGC	BstAPI^	R0654	GCANNNN/NTGC	BstAPI^
ApaLI	G/TGCAC	ApaLI	R0507	G/TGCAC	Alw44I, Vnel
\peKI	G/CWGC	ApeKI	R0643	G/CWGC	Tsel
		Tsel	R0591	G/CWGC	
Apol	R/AATTY	Apol-HF^	R3566	R/AATTY	Acsl, Xapl, Apol, Apol-HF^
ApyPI⊗	ATCGAC(20/18)	7.001 111	7.0000	1,70.111	7,001,7,401,7,401.11
Aqull⊗	GCCGNAC(20/18)				
	, ,				
\quIV ⊗	GRGGAAG(19/17)				
Arsl	(8/13)GACNNNNNNTTYG(11/6)				
Ascl	GG/CGCGCC	Ascl	R0558	GG/CGCGCC	PalAI, SgsI
Asel	AT/TAAT	Asel	R0526	AT/TAAT	PshBI, VspI
AsiGI	A/CCGGT	Agel-HF	R3552	A/CCGGT	Agel, Agel-HF, BshTl, CspAl, PinAl
AsiSI	GCGAT/CGC	AsiSI	R0630	GCGAT/CGC	Rgal, SfaAl, Sgfl
Asp103I⊗	CGRAGGC				
Asp700I	GAANN/NNTTC	XmnI	R0194	GAANN/NNTTC	MroXI, Pdml, Xmnl
Asp718I	G/GTACC	Acc65I	R0599	G/GTACC	Acc65I, KpnI^, KpnI-HF^
1307 101	d/d1A00	KpnI-HF^	R3142	GGTAC/C	Accool, replii , replii i li
A = = A O I	CICTACC				Audi Dial Vessil
AspA2I	C/CTAGG	AvrII	R0174	C/CTAGG	AvrII, BInI, XmaJI
AspBHI⊗	YSCNS(8/12)				
AspDUT2V⊗	GNGCAAC				
AspJHL3II⊗	CGCCCAG				
AspLEI	GCG/C	Hhal	R0139	GCG/C	BstHHI, Cfol, Hhal, Hin6I [^] , HinP1I [^] , HspAI [^]
		HinP1I^	R0124	G/CGC	
AspNIH4III ⊗	AAGAACB				
AspS9I	G/GNCC	Sau96I	R0165	G/GNCC	BmgT120I, Cfr13I, PspPI, Sau96I
AspSLV7III ⊗	GTCTCA	Ouusoi	110100	d/divoo	5111g 1 1201, 011 101, 1 3pt 1, 044301
Asp114pII⊗	AGCABCC	0- 001	DO10=	0.00100	AssOCI Desertions Of the D. D. O. CO.
Asul ⊗	G/GNCC	Sau96l	R0165	G/GNCC	AspS9I, BmgT120I, Cfr13I, PspPI, Sau96I
Asull	TT/CGAA	BstBl	R0519	TT/CGAA	Bpu14I, Bsp119I, BspT104I, BstBI, NspV, Sful
Asu14238IV ⊗	CGTRAC				
AsuC2I	CC/SGG	Ncil	R0196	CC/SGG	Bcnl, BpuMl, Ncil
AsuHPI	GGTGA(8/7)	Hphl	R0158	GGTGA(8/7)	Hphl
AsuNHI	G/CTAGC	Bmtl-HF^	R3658	GCTAG/C	Bmtl^, Bmtl-HF^, BspOI^, Nhel, Nhel-HF
		Nhel-HF	R3131	G/CTAGC	
AteTI⊗	GGGRAG			5/011.00	
Aval	C/YCGRG	Aval	R0152	C/YCGRG	Ama87I, BmeT110I, BsiHKCI, BsoBI, Eco88I
wai	o/ round				AIIIao 11, DIIIETTTUI, DSITINOI, DSUDI, EUU001
	0.101110.0	BsoBl	R0586	C/YCGRG	2 40 5 40 0 1 1 1 1 1 1 1
Avall	G/GWCC	Avall	R0153	G/GWCC	Bme18I, Eco47I, SinI, VpaK11BI
AvaIII ⊗	ATGCAT	Nsil-HF	R3127	ATGCA/T	EcoT22I, Mph1103I, Nsil, Nsil-HF, Zsp2I
AvrII	C/CTAGG	AvrII	R0174	C/CTAGG	AspA2I, BInI, XmaJI
Awo1030IV⊗	GCCRAG				
Axyl	CC/TNAGG	Bsu36I	R0524	CC/TNAGG	Bse21I, Bsu36I, Eco81I
3				.,	
	(40/4E)ACNININIOTAVO(40/7)	Bael	R0613	(10/15)ACNNNNGTAYC(12/7)	
Rael					
Bael BaeGl	(10/15)ACNNNNGTAYC(12/7) GKGCM/C	BaeGI	R0708	GKGCM/C	BseSI, BstSLI

		NEB			
ENZYME	SEQUENCE	ENZYME	NEB #	SEQUENCE	OTHER ISOSCHIZOMERS
Ball	TGG/CCA	Mscl	R0534	TGG/CCA	MISI, MIuNI, Mox20I, MscI, Msp20I
BamHI	G/GATCC	BamHI-HF	R3136	G/GATCC	BamHI, BamHI-HF
Banl	G/GYRCC	Banl	R0118	G/GYRCC	AccB1I, BshNI, BspT107I
Banll	GRGCY/C	Banll	R0119	GRGCY/C	Eco24I, EcoT38I, FriOI
BanLI ⊗	RTCAGG	Dami	110110	and on the	2002 11, 2001001, 11101
Barl	(7/12)GAAGNNNNNNTAC(12/7)				
Baul	CACGAG(-5/-1)	BssSI-v2	R0680	CACGAG(-5/-1)	BssSI, BssSI-v2, Bst2BI
Bbr11I⊗	GGRCAG	D3301 VZ	110000	0/10d/1d(5/ 1)	55501, 55501 V2, 551251
Bbr52II⊗	GGCGAG				
Bbr57III ⊗	GTRAAYG				
BbrPI	CAC/GTG	PmII	R0532	CAC/GTG	And Eco701 PmoCl Pmil PonCl
Bbsl		BbsI-HF	R3539		AcvI, Eco72I, PmaCI, PmII, PspCI
	GAAGAC(2/6)			GAAGAC(2/6)	Bbsl, Bbsl-HF, Bpil, BstV2l
Bbvl	GCAGC(8/12)	Bbvl	R0173	GCAGC(8/12)	BseXI, BstV1I, Lsp1109I
BbvII ⊗	GAAGAC(2/6)	BbsI-HF	R3539	GAAGAC(2/6)	Bbsl, Bbsl-HF, Bpil, BstV2I
Bbv12I	GWGCW/C	BsiHKAI	R0570	GWGCW/C	Alw21I, BsiHKAI
BbvCl	CCTCAGC(-5/-2)	BbvCl	R0601	CCTCAGC(-5/-2)	
Bccl	CCATC(4/5)	Bccl	R0704	CCATC(4/5)	D 51
Bce83I⊗	CTTGAG(16/14)	BpuEl	R0633	CTTGAG(16/14)	BpuEl
Bce3081I⊗	TAGGAG				
BceAl	ACGGC(12/14)	BceAl	R0623	ACGGC(12/14)	
Bcefl ⊗	ACGGC(12/13)	BceAI^	R0623	ACGGC(12/14)	BceAI^
Bcgl	(10/12)CGANNNNNNTGC(12/10)	Bcgl	R0545	(10/12)CGANNNNNNTGC(12/10)	
BciT130I	CC/WGG	BstNI	R0168	CC/WGG	AjnI^, BseBI, BstNI, Bst2UI, EcoRII^, MvaI, Psp6I^, PspGI^
		PspGI [^]	R0611	/CCWGG	
BciVI	GTATCC(6/5)	BciVI	R0596	GTATCC(6/5)	Bful, Bsul
BcII	T/GATCA	BcII-HF	R3160	T/GATCA	Fbal, Ksp22I, BcII, BcII-HF
Bcnl	CC/SGG	Ncil	R0196	CC/SGG	AsuC2I, BpuMI, Ncil
BcoDI	GTCTC(1/5)	BcoDI	R0542	GTCTC(1/5)	Alw26I, BsmAI, BstMAI
		BsmAl	R0529	GTCTC(1/5)	
Bcul	A/CTAGT	Spel-HF	R3133	A/CTAGT	Ahll, Spel, Spel-HF
Bdal ⊗	(10/12)TGANNNNNNTCA(12/10)	·			
Betl ⊗	W/CCGGW	BsaWl	R0567	W/CCGGW	BsaWI
Bfal	C/TAG	Bfal	R0568	C/TAG	FspBI, Mael, SspMI, XspI
BfaSII ⊗	GANGGAG			., .	and the second of the
Bfil ⊗	ACTGGG(5/4)	Bmrl	R0600	ACTGGG(5/4)	Bmrl, Bmul
Bfml	C/TRYAG	Sfcl	R0561	C/TRYAG	BstSFI, SfcI
Bfol	RGCGC/Y	Haell	R0107	RGCGC/Y	BstH2I, Haell
BfrI	C/TTAAG	AfIII	R0520	C/TTAAG	AfIII, BspTI, BstAFI, MspCI, Vha464I
Bful	GTATCC(6/5)	BciVI	R0596	GTATCC(6/5)	BciVI, Bsul
BfuAl	ACCTGC(4/8)	BfuAl	R0701	ACCTGC(4/8)	Acc36I, BspMI, Bvel
DiuAi	A00100(4/0)	BspMI	R0502	ACCTGC(4/8)	Accour, Dapivii, Dvci
Bgll	GCCNNNN/NGGC	Bgll	R0143	GCCNNNN/NGGC	
Bglll	A/GATCT	BgIII	R0144	A/GATCT	AcIMI Alud DooDI
Binl ⊗	GGATC(4/5) GC/NGC	Alwl	R0513	GGATC(4/5)	Aciwi, Alwi, BspPi Bisi^, Giui. Pkri^
Bisl	1				BISI", GIUI, PKII"
BkrAM31DI⊗	RTTAAATM				
Ble402II ⊗	GRAGCAG	Acodi	D0474	0.074.00	Ann AOL Avell Venn II
Bini	C/CTAGG	AvrII	R0174	C/CTAGG	AspA2I, AvrII, XmaJI
BloAll ⊗	GAGGAC				
Blpl	GC/TNAGC	Blpl	R0585	GC/TNAGC	Bpu1102I, Bsp1720I
BISI	GCN/GC				Bisl^, Glul^, Pkrl
BmcAl	AGT/ACT	Scal-HF	R3122	AGT/ACT	Scal, Scal-HF, Zrml
Bme18I	G/GWCC	Avall	R0153	G/GWCC	Avall, Eco47I, Sinl, VpaK11BI
Bme1390I	CC/NGG	ScrFl	R0110	CC/NGG	BmrFI, BstSCI [^] , MspR9I, ScrFI, StyD4I [^]
		StyD4I^	R0638	/CCNGG	
BmeRI	GACNNN/NNGTC	Ahdl	R0584	GACNNN/NNGTC	Ahdl, Dril, Eam11051
BmeT110I	C/YCGRG	Aval	R0152	C/YCGRG	Ama87I, Aval, BsiHKCI, BsoBI, Eco88I
		BsoBI	R0586	C/YCGRG	
Bmgl⊗	GKGCCC				
BmgBl	CACGTC(-3/-3)	BmgBl	R0628	CACGTC(-3/-3)	Ajil, Btrl
BmgT120I	G/GNCC	Sau96I	R0165	G/GNCC	AspS9I, Cfr13I, PspPI, Sau96I
Bmil	GGN/NCC	NIaIV	R0126	GGN/NCC	BspLI, NIaIV, PspN4I
Bmrl	ACTGGG(5/4)	Bmrl	R0600	ACTGGG(5/4)	Bmul
BmrFI	CC/NGG	ScrFl	R0110	CC/NGG	Bme1390I, BstSCI [^] , MspR9I, ScrFI, StyD4I [^]
2.1111	oo,naa	StyD4I^	R0638	/CCNGG	Sind toos, botton, morning out i, otypen
BmsI	GCATC(5/9)	SfaNI	R0172	GCATC(5/9)	Lwel, SfaNI
ופווטו	donio(J/3)	JIAINI	110172	UUATU(3/3)	LYVOI, OIGINI

		NEB			
ENZYME	SEQUENCE	ENZYME	NEB #	SEQUENCE	OTHER ISOSCHIZOMERS
Bmtl	GCTAG/C	BmtI-HF^	R3658	GCTAG/C	Bmtl, Bmtl-HF [^] , AsuNHI [^] , BspOI, NheI [^] , NheI-HF [^]
		Nhel-HF [^]	R3131	G/CTAGC	
Bmul	ACTGGG(5/4)	Bmrl	R0600	ACTGGG(5/4)	Bmrl
Boxl	GACNN/NNGTC	PshAl	R0593	GACNN/NNGTC	BstPAI, PshAI
Bpil	GAAGAC(2/6)	BbsI-HF	R3539	GAAGAC(2/6)	Bbsl, Bbsl-HF, BstV2I
3pH	(8/13)GAGNNNNNCTC(13/8)				
Bpml	CTGGAG(16/14)	Bpml	R0565	CTGGAG(16/14)	Gsul
Bpu10I	CCTNAGC(-5/-2)	Bpu10I	R0649	CCTNAGC(-5/-2)	
Bpu14I	TT/CGAA	BstBI	R0519	TT/CGAA	Asull, Bsp1191, BspT1041, BstBI, NspV, Sful
Bpu1102I	GC/TNAGC	Blpl	R0585	GC/TNAGC	Blpl, Bsp1720l
3puEl	CTTGAG(16/14)	BpuEl	R0633	CTTGAG(16/14)	51p1, 53p17201
ВриМI	CC/SGG	Ncil	R0196	CC/SGG	AsuC2I, BcnI, Ncil
Bsal	GGTCTC(1/5)	Bsal-HFv2	R3733	GGTCTC(1/5)	Bsal, Bsal-HFv2, Bso31l, BspTNI, Eco31l
Bsa29I	AT/CGAT	BspDI	R0557	AT/CGAT	BseCI, BshVI, BspDI, Bsu15I, BsuTUI, Clal
		Clal	R0197	AT/CGAT	
BsaAl	YAC/GTR	BsaAl	R0531	YAC/GTR	BstBAI, Ppu21I
BsaBI	GATNN/NNATC	BsaBI	R0537	GATNN/NNATC	Bse8l, BseJl
BsaHI	GR/CGYC	BsaHl	R0556	GR/CGYC	Acyl, BssNI, BstACI, Hin1I, Hsp92I
BsaJI	C/CNNGG	BsaJI	R0536	C/CNNGG	BseDI, BssECI
BsaWl	W/CCGGW	BsaWI	R0567	W/CCGGW	
BsaXI	(9/12)ACNNNNNCTCC(10/7)	BsaXI	R0609	(9/12)ACNNNNNCTCC(10/7)	
Bsbl ⊗	CAACAC(21/19)			,,, , :	
Bsc4l	CCNNNNN/NNGG	BsII	R0555	CCNNNNN/NNGG	Afil, BseLl, Bsll
BscGl⊗	CCCGT	Doll	110000	oommining mad	7till, Dooll, Doll
Bse1I		Dorl	D0E07	ACTCC(1/ 1)	DooMI Dorl
	ACTGG(1/-1)	Bsrl	R0527	ACTGG(1/-1)	BseNI, BsrI
Bse8I	GATNN/NNATC	BsaBI	R0537	GATNN/NNATC	BsaBI, BseJI
Bse21I	CC/TNAGG	Bsu36I	R0524	CC/TNAGG	Axyl, Bsu36l, Eco81l
Bse118I	R/CCGGY	BsrFI-v2	R0682	R/CCGGY	BsrFI-v2, BssAI, Cfr10I
BseAl	T/CCGGA	BspEl	R0540	T/CCGGA	AccIII, Aor13HI, Bsp13I, BspEI, Kpn2I, Mrol
BseBl	CC/WGG	BstNI	R0168	CC/WGG	AjnI^, BciT130I, BstNI, Bst2UI, EcoRII^, Mval, Psp6I^, PspGI^
		PspGI^	R0611	/CCWGG	
BseCl	AT/CGAT	BspDI	R0557	AT/CGAT	Bsa29I, BshVI, BspDI, Bsu15I, BsuTUI, Clal
	,	Clal	R0197	AT/CGAT	
BseDI	C/CNNGG	BsaJI	R0536	C/CNNGG	BsaJI, BssECI
Bse3DI	GCAATG(2/0)	BsrDI	R0574	GCAATG(2/0)	BseMI, BsrDI
BseGl	. ,		R0647	, ,	
DSEGI	GGATG(2/0)	BtsCI		GGATG(2/0)	BstF5I, BtsCI, FokI^
D 11	O ATNINI MINIATO	Fokl^	R0109	GGATG(9/13)	D DI D 01
BseJI	GATNN/NNATC	BsaBI	R0537	GATNN/NNATC	BsaBI, Bse8I
BseLI	CCNNNNN/NNGG	BsII	R0555	CCNNNNN/NNGG	Afil, Bsc4l, Bsll
BseMI	GCAATG(2/0)	BsrDI	R0574	GCAATG(2/0)	Bse3DI, BsrDI
BseMII	CTCAG(10/8)	BspCNI^	R0624	CTCAG(9/7)	BspCNI [^]
BseNI	ACTGG(1/-1)	Bsrl	R0527	ACTGG(1/-1)	Bse1I, BsrI
BsePI	G/CGCGC	BssHII	R0199	G/CGCGC	BssHII, Paul, Ptel
BseRI	GAGGAG(10/8)	BseRI	R0581	GAGGAG(10/8)	,,
BseSI	GKGCM/C	BaeGI	R0708	GKGCM/C	BaeGI, BstSLI
BseXI	GCAGC(8/12)	Bbvl	R0173	GCAGC(8/12)	Bbvl, BstV11, Lsp11091
BseX3I	C/GGCCG		R3505	C/GGCCG	
		Eagl-HF			BstZI, Eagl, Eagl-HF, EclXI, Eco52I
BseYl	CCCAGC(-5/-1)	BseYI	R0635	CCCAGC(-5/-1)	Gsal^, PspFl
Bsgl	GTGCAG(16/14)	Bsgl	R0559	GTGCAG(16/14)	
Bsh1236l	CG/CG	BstUI	R0518	CG/CG	Accil, BspFNI, BstFNI, BstUI, MvnI
Bsh1285I	CGRY/CG	BsiEl	R0554	CGRY/CG	BsiEl, BstMCl
BshFl	GG/CC	HaellI	R0108	GG/CC	Bsnl, BspANI, BsuRI, HaellI
BshNI	G/GYRCC	Banl	R0118	G/GYRCC	AccB1I, BanI, BspT107I
BshTl	A/CCGGT	Agel-HF	R3552	A/CCGGT	Agel, Agel-HF, AsiGl, CspAl, PinAl
BshVI	AT/CGAT	BspDI	R0557	AT/CGAT	Bsa29I, BseCl, BspDl, Bsu15I, BsuTUI, Clal
		Clal	R0197	AT/CGAT	,, , ,, ,,
Bsil ⊗	CACGAG(-5/-1)	BssSI-v2	R0680	CACGAG(-5/-1)	Baul, BssSI-v2, Bst2BI
BsiEl	, ,			. ,	
	CGRY/CG	BsiEl	R0554	CGRY/CG	Bsh1285I, BstMCI
BsiHKAI	GWGCW/C	BsiHKAI	R0570	GWGCW/C	Alw21I, Bbv12I
BsiHKCI	C/YCGRG	Aval	R0152	C/YCGRG	Ama87I, Aval, BmeT110I, BsoBI, Eco88I
		BsoBl	R0586	C/YCGRG	
BsiSI	C/CGG	Hpall	R0171	C/CGG	Hapll, Hpall, Mspl
		Mspl	R0106	C/CGG	
		IMShi	110100		
BsiWl	C/GTACG	BsiWI-HF	R3553	C/GTACG	BsiWI, BsiWI-HF, Pfl23II, PspLI

ENZYME	SEQUENCE	NEB Enzyme	NEB#	SEQUENCE	OTHER ISOSCHIZOMERS
	CCNNNNN/NNGG	Bsll	R0555	CCNNNNN/NNGG	
BsII					Afil, Bsc4l, BseLl
BsIFI	GGGAC(10/14)	BsmFl	R0572	GGGAC(10/14)	BsmFl, Faql
Bsml	GAATGC(1/-1)	Bsml	R0134	GAATGC(1/-1)	Mva1269I, PctI
BsmAl	GTCTC(1/5)	BcoDI BsmAI	R0542 R0529	GTCTC(1/5) GTCTC(1/5)	Alw26I, BcoDI, BstMAI
BsmBI	CGTCTC(1/5)	BsmBI Esp3I	R0580 R0734	CGTCTC(1/5) CGTCTC(1/5)	Esp3I
BsmFI	GGGAC(10/14)	BsmFl	R0572	GGGAC(10/14)	BsIFI, Fagl
BsnI	GG/CC	Haelll		GG/CC	7 1
Bso31I	GGTCTC(1/5)	Bsal-HFv2	R0108 R3733	GGTCTC(1/5)	BshFI, BspANI, BsuRI, HaellI
	, ,			, ,	Bsal, Bsal-HFv2, BspTNI, Eco311
BsoBl	C/YCGRG	Aval BsoBl	R0152 R0586	C/YCGRG C/YCGRG	Ama87I, Aval, BmeT110I, BsiHKCI, Eco88I
Don 101	T/CCGGA		R0540	T/CCGGA	Applit Apri 2111 Booki Booki Kooki Mroj
Bsp13I		BspEl			AccIII, Aor13HI, BseAI, BspEI, Kpn2I, Mrol
Bsp19I	C/CATGG (8/13)GACNNNNNNTGG(12/7)	Ncol-HF	R3193	C/CATGG	Ncol, Ncol-HF
Bsp24l ⊗	, ,	New LUE	D0100	T00/00A	DivMI Mad Mad HE Dad
Bsp68I	TCG/CGA	Nrul-HF	R3192	TCG/CGA	BtuMI, NruI, NruI-HF, RruI
Bsp119I	TT/CGAA	BstBI	R0519	TT/CGAA	Asull, Bpu14I, BspT104I, BstBI, NspV, Sful
Bsp120I	G/GGCCC	Apal^	R0114	GGGCC/C	Apal^, PspOMI
		Psp0MI	R0653	G/GGCCC	
Bsp143I	/GATC	DpnII	R0543	/GATC	BssMl, BstKTI^, BstMBl, Dpnll, Kzo9l, Mbol, Ndell, Sau3Al
		Mbol	R0147	/GATC	
		Sau3Al	R0169	/GATC	
Bsp460III ⊗	CGCGCAG				
Bsp1286I	GDGCH/C	Bsp1286I	R0120	GDGCH/C	MhII, Sdul
Bsp1407I	T/GTACA	BsrGI-HF	R3575	T/GTACA	BsrGI, BsrGI-HF, BstAUI
Bsp1720I	GC/TNAGC	Blpl	R0585	GC/TNAGC	Blpl, Bpu1102l
Bsp3004IV ⊗	CCGCAT				
BspACI	CCGC(-3/-1)	Acil	R0551	CCGC(-3/-1)	Acil, Ssil
BspANI	GG/CC	HaeIII	R0108	GG/CC	BshFi, BsnI, BsuRi, Haelli
BspCNI	CTCAG(9/7)	BspCNI	R0624	CTCAG(9/7)	BseMII^
BspDI	AT/CGAT	BspDI	R0557	AT/CGAT	Bsa29I, BseCI, BshVI, Bsu15I, BsuTUI, Clal
		Clal	R0197	AT/CGAT	
BspEl	T/CCGGA	BspEl	R0540	T/CCGGA	AccIII, Aor13HI, BseAI, Bsp13I, Kpn2I, MroI
BspFNI	CG/CG	BstUI	R0518	CG/CG	Accll, Bsh1236I, BstFNI, BstUI, MvnI
BspGl ⊗	CTGGAC				
BspHI	T/CATGA	BspHI	R0517	T/CATGA	Ccil, Pagl
BspLI	GGN/NCC	NIaIV	R0126	GGN/NCC	Bmil, NIaIV, PspN4I
BspLU11I ⊗	A/CATGT	Pcil	R0655	A/CATGT	Pcil, Pscl
BspMI	ACCTGC(4/8)	BfuAl	R0701	ACCTGC(4/8)	Acc36l, BfuAl, Bvel
Борічіі	7,00140(1,0)	BspMI	R0502	ACCTGC(4/8)	7,00001, 510711, 5701
BspMII ⊗	T/CCGGA	BspEl	R0540	T/CCGGA	Accili, Aor13Hi, BseAl, Bsp13i, BspEi, Kpn2i, Mroi
BspMAI	CTGCA/G	PstI-HF	R3140	CTGCA/G	Pstl, Pstl-HF
BspNCI⊗	CCAGA	1 30 111	110110	OTODAYO	1 30, 1 30 111
		Bmtl-HF	R3658	CCTAC/C	Acualitia Desti Desti LIE Nibala Nibal LIEA
BspOI	GCTAG/C	Nhel-HF^	R3131	GCTAG/C G/CTAGC	AsuNHI^, Bmtl, Bmtl-HF, NheI^, NheI-HF^
BspPI	GGATC(4/5)	Alwl	R0513	GGATC(4/5)	AcIWI, AlwI
BspQI	GCTCTTC(1/4)	BspQl	R0712	GCTCTTC(1/4)	Lgul, PciSl, Sapl
Бэрці	4010110(1/4)	Sapl	R0569	GCTCTTC(1/4)	Lydi, i didi, dapi
BspTI	C/TTAAG		R0520	C/TTAAG	AfIII, BfrI, BstAFI, MspCI, Vha464I
		AfIII			
BspT104I	TT/CGAA	BstBl	R0519	TT/CGAA	Asull, Bpu14I, Bsp119I, BstBI, NspV, Sful
BspT107I	G/GYRCC	Banl	R0118	G/GYRCC	AccB1I, BanI, BshNI
BspTNI	GGTCTC(1/5)	Bsal-HFv2	R3733	GGTCTC(1/5)	Bsal, Bsal-HFv2, Bso31I, Eco31I
Bsrl	ACTGG(1/-1)	Bsrl	R0527	ACTGG(1/-1)	Bse1I, BseNI
BsrBI	CCGCTC(-3/-3)	BsrBI	R0102	CCGCTC(-3/-3)	AccBSI, Mbil
BsrDI	GCAATG(2/0)	BsrDI	R0574	GCAATG(2/0)	Bse3DI, BseMI
BsrFl	R/CCGGY	BsrFI-v2	R0682	R/CCGGY	Bse118I, BsrFI-v2, BssAI, Cfr10I
BsrGI	T/GTACA	BsrGI-HF	R3575	T/GTACA	BsrGI, BsrGI-HF, Bsp1407I, BstAUI
BssAl	R/CCGGY	BsrFI-v2	R0682	R/CCGGY	Bse118I, BsrFI-v2, Cfr10I
BssECI	C/CNNGG	BsaJI	R0536	C/CNNGG	BsaJI, BseDI
BssHII	G/CGCGC	BssHII	R0199	G/CGCGC	BsePl, Paul, Ptel
BssMI	/GATC	DpnII	R0543	/GATC	Bsp143I, BstKTI^, BstMBI, DpnII, Kzo9I, MboI, NdeII, Sau3AI
		Mbol	R0147	/GATC	
		Sau3Al	R0169	/GATC	
		JauJAI	110100		
BssNI	GR/CGYC	BsaHl	R0556	GR/CGYC	Acyl, BsaHl, BstACl, Hin1l, Hsp92l
BssNI BssNAI	GR/CGYC GTA/TAC				Acyl, BsaHl, BstACl, Hin1l, Hsp92l Bst1107l, BstZ17l, BstZ17l-HF
		BsaHl	R0556	GR/CGYC	

ENZYME	SEQUENCE	NEB Enzyme	NEB#	SEQUENCE	OTHER ISOSCHIZOMERS
Bst6l	CTCTTC(1/4)	Earl	R0528	CTCTTC(1/4)	Eam1104I, Earl
Bst1107I	GTA/TAC	BstZ17I-HF	R3594	GTA/TAC	BssNAI, BstZ17I, BstZ17I-HF
BstACI	GR/CGYC	BsaHl	R0556	GR/CGYC	Acyl, BsaHl, BssNl, Hin1l, Hsp92l
BstAFI	C/TTAAG	AfIII	R0520	C/TTAAG	AfIII, BfrI, BspTI, MspCI, Vha464I
BstAPI	GCANNNN/NTGC	BstAPI	R0654	GCANNNN/NTGC	Tim, Bill, Bopti, Mopol, Vilato II
BstAUI	T/GTACA	BsrGI-HF	R3575	T/GTACA	Bsp1407I, BsrGI, BsrGI-HF
BstBI	TT/CGAA	BstBl	R0519	TT/CGAA	Asull, Bpu14I, Bsp119I, BspT104I, NspV, Sful
Bst2BI	CACGAG(-5/-1)	BssSI-v2	R0680	CACGAG(-5/-1)	Baul, BssSl, BssSl-v2
BstBAI	YAC/GTR	BsaAl	R0531	YAC/GTR	BsaAl, Ppu21I
Bst4CI	ACN/GT	HpyCH4III	R0618	ACN/GT	HpyCH4III, Taal
BstC8I	GCN/NGC	Cac8l	R0579	GCN/NGC	Cac8l
BstDEI	C/TNAG	Ddel	R0175	C/TNAG	Ddel, HpyF3I
BstDSI	C/CRYGG	Btgl	R0608	C/CRYGG	Btgl
BstEll	G/GTNACC	BstEII-HF	R3162	G/GTNACC	BstEII, BstEII-HF, BstPI, Eco91I, Eco065I, PspEI
BstENI	CCTNN/NNNAGG	EcoNI	R0521	CCTNN/NNNAGG	EcoNI, XagI
BstF5I	GGATG(2/0)	BtsCl	R0647	GGATG(2/0)	BseGI, BtsCI, FokI^
Dott Of	da/11d(2/0)	Fokl^	R0109	GGATG(9/13)	556di, 5660i, i okt
BstFNI	CG/CG	BstUI	R0518	CG/CG	Accil, Bsh1236i, BspFNI, BstUl, Mvni
BstH2I	RGCGC/Y	Haell	R0107	RGCGC/Y	Bfol, Haell
BstHHI	GCG/C	Hhal	R0139	GCG/C	AspLEI, CfoI, HhaI, Hin6I [^] , HinP1I [^] , HspAI [^]
Dourin	404/0	HinP1I^	R0124	G/CGC	7.6pcc1, 6101, 111101 , 11111 11 , 115p711
BstKTI	GAT/C	DpnII^	R0543	/GATC	Bsp143I^, BssMI^, BstMBI^, DpnII^, Kzo9I^, MboI^, NdeII^, Sau3AI/
Douver	uni/o	Mbol^	R0147	/GATC	DSP 1401, DSSIVII, DSLIVIDI, DPIIII, NZOSI, IVIDOI, NUCII, SUUSAI
		Sau3AI^	R0169	/GATC	
BstMAI	GTCTC(1/5)	BcoDI	R0542	GTCTC(1/5)	Alw26I, BcoDI, BsmAI
DSUVIAI	d1010(1/3)	BsmAl	R0529	GTCTC(1/5)	AIW20I, DCUDI, DSIIIAI
Do+MDI	ICATO		R0543	` ,	Pont 421 PooMI PotVTIA Ponti Kroßi Mhol Ndoll Cou 2Al
BstMBI	/GATC	DpnII Mbol	R0147	/GATC /GATC	Bsp1431, BssMI, BstKTI^, DpnII, Kzo9I, MboI, NdeII, Sau3AI
			R0169		
Do+MCI	CCDV/CC	Sau3Al		/GATC	Dob100EL DoiEL
BstMCI	CGRY/CG	BsiEl	R0554	CGRY/CG	Bsh1285I, BsiEl
BstMWI	GCNNNNN/NNGC	Mwol	R0573	GCNNNNN/NNGC	HpyF10VI, Mwol
BstNI	CC/WGG	BstNI	R0168	CC/WGG	AjnI^, BciT130I, BseBI, Bst2UI, EcoRII^, MvaI, Psp6I^, PspGI^
D-4NO	DOATO N/	PspGI^	R0611	/CCWGG	New L Vest
BstNSI	RCATG/Y	Nspl	R0602	RCATG/Y	Nspl, Xcel
BstPI	G/GTNACC	BstEII-HF	R3162	G/GTNACC	BstEII, BstEII-HF, Eco91I, Eco065I, PspEI
BstPAI	GACNN/NNGTC	PshAl	R0593	GACNN/NNGTC	Boxl, PshAl
BstSCI	/CCNGG	ScrFI [^]	R0110	CC/NGG	Bme1390I^, BmrFI^, MspR9I^, ScrFI^, StyD4I
D 1051	0.FD\40	StyD4I	R0638	/CCNGG	D(1.0()
BstSFI	C/TRYAG	Sfcl	R0561	C/TRYAG	Bfml, Sfcl
BstSLI	GKGCM/C	BaeGI	R0708	GKGCM/C	BaeGI, BseSI
BstSNI	TAC/GTA	SnaBl	R0130	TAC/GTA	Eco105I, SnaBI
BstUI	CG/CG	BstUI	R0518	CG/CG	Accil, Bsh1236i, BspFNi, BstFNi, Mvni
Bst2UI	CC/WGG	BstNI	R0168	CC/WGG	AjnI^, BciT130I, BseBI, BstNI, EcoRII^, Mval, Psp6I^, PspGI^
		PspGI^	R0611	/CCWGG	
BstV1I	GCAGC(8/12)	Bbvl	R0173	GCAGC(8/12)	Bbvl, BseXI, Lsp1109I
BstV2I	GAAGAC(2/6)	BbsI-HF	R3539	GAAGAC(2/6)	Bbsl, Bbsl-HF, Bpil
BstXI	CCANNNN/NTGG	BstXI	R0113	CCANNNNN/NTGG	
BstX2I	R/GATCY	BstYI	R0523	R/GATCY	BstYI, MfII, Psul
BstYI	R/GATCY	BstYI	R0523	R/GATCY	BstX2I, MfII, Psul
BstZI	C/GGCCG	Eagl-HF	R3505	C/GGCCG	BseX3I, EagI, EagI-HF, EcIXI, Eco52I
BstZ17I	GTA/TAC	BstZ17I-HF	R3594	GTA/TAC	BssNAI, Bst1107I, BstZ17I, BstZ17I-HF
Bsul	GTATCC(6/5)	BciVI	R0596	GTATCC(6/5)	BciVI, Bful
Bsu15I	AT/CGAT	BspDI	R0557	AT/CGAT	Bsa29I, BseCI, BshVI, BspDI, BsuTUI, Clal
		Clal	R0197	AT/CGAT	
Bsu36I	CC/TNAGG	Bsu36I	R0524	CC/TNAGG	Axyl, Bse21I, Eco81I
BsuRI	GG/CC	HaellI	R0108	GG/CC	BshFI, BsnI, BspANI, HaeIII
BsuTUI	AT/CGAT	BspDI	R0557	AT/CGAT	Bsa29I, BseCl, BshVI, BspDI, Bsu15I, Clal
		Clal	R0197	AT/CGAT	
Btgl	C/CRYGG	Btgl	R0608	C/CRYGG	BstDSI
BtgZl	GCGATG(10/14)	BtgZl	R0703	GCGATG(10/14)	
Btrl	CACGTC(-3/-3)	BmgBl	R0628	CACGTC(-3/-3)	Ajil, BmgBl
Btsl	GCAGTG(2/0)	Btsl-v2	R0667	GCAGTG(2/0)	. ₇ SgS.
BtsIMutl	CAGTG(2/0)	BtsIMutl	R0664	CAGTG(2/0)	
BtsCl	GGATG(2/0)	BtsCl	R0647	GGATG(2/0)	BseGI, BstF5I, FokI^
D.001	uu/11u(2/0)	Fokl [^]	R0109	GGATG(9/13)	Doddi, Doti Oi, i Oiti
		1 UNI	110103	durita(3/10)	

ENTWAL	OFOURNOR	NEB	NED "	OFOLIFNOE	OTHER ISOSOHITANIERS
ENZYME	SEQUENCE	ENZYME	NEB#	SEQUENCE	OTHER ISOSCHIZOMERS
BtuMI	TCG/CGA	Nrul-HF	R3192	TCG/CGA	Bsp68I, Nrul, Nrul-HF, Rrul
Bvel	ACCTGC(4/8)	BfuAl	R0701	ACCTGC(4/8)	Acc36I, BfuAI, BspMI
C		BspMI	R0502	ACCTGC(4/8)	
Cac8I	GCN/NGC	Cac8l	R0579	GCN/NGC	BstC8I
Cail	CAGNNN/CTG	AlwNI	R0514	CAGNNN/CTG	AlwNI, PstNI
Cal14237I ⊗	GGTTAG	AIWINI	NUJ 14	GAGININI/GTG	Alwivi, FSuvi
CalB3II⊗	GRTTRAG				
Caull ⊗	CC/SGG	Ncil	R0196	CC/SGG	AsuC2I, BcnI, BpuMI, Ncil
Cau10061II ⊗	GTTAAT	NGII	110130	00/300	ASUOZI, DOIII, DPUIVII, NOII
Cba13II ⊗	AGGAAT				
Cba16038I⊗	CCTNAYNC				
Cba100301⊗ Cbo67071IV⊗	GCRGAAG				
Cchll ⊗	GGARGA(11/9)				
CchIII ⊗	CCCAAG(20/18)				
Cch467III ⊗	GNGAAAY				
Ccil	T/CATGA	BspHI	R0517	T/CATGA	BspHI, PagI
CciNI	GC/GGCCGC	NotI-HF	R3189	GC/GGCCGC	Notl, Notl-HF
Cco14983V ⊗	GGGTDA	1100 111	110100	dojudoodo	100, 100 111
Cco14983VI ⊗	GCYGA				
CcrNAIII ⊗	CGACCAG				
Cdil⊗	CATCG(-1/-1)				
Cdi81III ⊗	GCMGAAG				
Cdi11397I⊗	GCGCAG				
Cdpl⊗	GCGGAG(20/18)				
Cdu23823II ⊗	GTGAAG				
Cfol	GCG/C	Hhal	R0139	GCG/C	AspLEI, BstHHI, Hhal, Hin6I^, HinP1I^, HspAI^
0.01	464,6	HinP1I^	R0124	G/CGC	7.69221, 2001111, 111101 , 11111 11 , 1169711
CfrI ⊗	Y/GGCCR	Eael	R0508	Y/GGCCR	Acol, Eael
Cfr9I	C/CCGGG	Smal^	R0141	CCC/GGG	Smal^, TspMI, Xmal
	0,00000	TspMI	R0709	C/CCGGG	,,
		Xmal	R0180	C/CCGGG	
Cfr10I	R/CCGGY	BsrFI-v2	R0682	R/CCGGY	Bse118I, BsrFI-v2, BssAI
Cfr13I	G/GNCC	Sau96I	R0165	G/GNCC	AspS9I, BmgT120I, PspPI, Sau96I
Cfr42I	CCGC/GG	SacII	R0157	CCGC/GG	Kspl, Sacil, Sfr303I, SgrBl
Cfupf3II ⊗	GARCAG				
Cgl13032l⊗	GGCGCA				
Cgl13032II ⊗	ACGABGG				
Cjel ⊗	(8/14)CCANNNNNNGT(15/9)				
Cie265V ⊗	GKAAGC				
Cje54107III ⊗	GKAAYC				
CjeFIII ⊗	GCAAGG				
CjeFV ⊗	GGRCA				
CjeNII ⊗	GAGNNNNNGT				
CjeNIII ⊗	GKAAYG(19/17)				
CjeNV ⊗	CCYGA				
CjePI ⊗	(7/13)CCANNNNNNNTC(14/8)				
CjeP659IV ⊗	CACNNNNNNGAA				
Cjul ⊗	CAYNNNNRTG				
CjuII ⊗	CAYNNNNCTC				
Clal	AT/CGAT	BspDI Clal	R0557 R0197	AT/CGAT AT/CGAT	Bsa29I, BseCI, BshVI, BspDI, Bsu15I, BsuTUI
Cla11845III ⊗	GCGAA	Oiui	110101	mpount	
Cly7489II ⊗	AAAAGRG				
Cma23826l ⊗	CGGAAG				
Cpol	CG/GWCCG	RsrII	R0501	CG/GWCCG	Cspl, Rsrll, Rsr2l
Csel	GACGC(5/10)	Hgal	R0154	GACGC(5/10)	Hgal
Csil	A/CCWGGT	SexAl	R0605	A/CCWGGT	Mabl, SexAl
Cspl	CG/GWCCG	RsrII	R0501	CG/GWCCG	Cpol, Rsrll, Rsr2l
Csp6l	G/TAC	CviQI	R0639	G/TAC	Afal^, CviQl, Rsal^, RsaNl
	-,	Rsal^	R0167	GT/AC	, , , , , , , , , , , , , , , , , , , ,
	GGAGGC				
Csp2014I ⊗					
Csp2014I⊗ CspAI	A/CCGGT	Agel-HF	R3552	A/CCGGT	Agel, Agel-HF, AsiGl, BshTl, PinAl
· ·		Agel-HF CspCl	R3552 R0645	A/CCGGT (11/13)CAANNNNNGTGG(12/10)	Agel, Agel-HF, AsiGl, BshTl, PinAl

ENZYME	SEQUENCE	NEB Enzyme	NEB#	SEQUENCE	OTHER ISOSCHIZOMERS
CviAll	C/ATG	CviAll	R0640	C/ATG	Fael^, Fatl^^, Hin1II^, Hsp92II^, NIaIII^
07// 111	o//iid	FatI^^	R0650	/CATG	radi , rati , riirrii , riopozii , maiii
		NIaIII^	R0125	CATG/	
Cvi II	DC/CV				CviVI 1
CviJI CviKI-1	RG/CY RG/CY	CviKI-1 CviKI-1	R0710 R0710	RG/CY	CviKl-1 CviJl
CviQl	G/TAC	CviQI	R0639	RG/CY G/TAC	
CVIQI	G/ IAC	Rsal^	R0167	GT/AC	Afal^, Csp6I, Rsal^, RsaNI
CviRI⊗	TG/CA	HpyCH4V	R0620	TG/CA	HpyCH4V
D					
Ddel	C/TNAG	Ddel	R0175	C/TNAG	BstDEI, HpyF3I
Dde51507I ⊗	CCWGG				
Dinl	GGC/GCC	Kasl^	R0544	G/GCGCC	Egel, Ehel, Kasl^, Mly113l^^, Narl^^, PluTl^^^, Sfol, SspDl^
		Narl^^	R0191	GG/CGCC	
		PluTI^^^	R0713	GGCGC/C	
D .	0.1.770	Sfol	R0606	GGC/GCC	
Dpnl	GA/TC	Dpnl	R0176	GA/TC	Mall
DpnII	/GATC	DpnII	R0543	/GATC	Bsp143I, BssMI, BstKTI^, BstMBI, Kzo9I, MboI, NdeII, Sau3AI
		Mbol	R0147	/GATC	
D 1	TTT /A A A	Sau3AI	R0169	/GATC	
Dral	TTT/AAA	Dral	R0129	TTT/AAA	Fac01001
Drall ⊗	RG/GNCCY	EcoO109I	R0503	RG/GNCCY	EcoO1091
Dralll Drall	CACNNN/GTG	DrallI-HF	R3510	CACNNN/GTG	Adel, Dralli, Dralli-HF
DraRI ⊗	CAAGNAC(20/18)	Dedl	DOEGO	C A CAINIAIAI (AIAICTC	And DooDI
Drdl	GACNNNN/NNGTC	Drdl	R0530	GACNNNN/NNGTC	Aasl, DseDl
DrdII ⊗	GAACCA	المما ٨	DOE04	CACNINIALAINICTO	Abd Decal Food1051
Dril	GACNNN/NNGTC	Ahdl	R0584	GACNNN/NNGTC	Ahdl, BmeRl, Eam1105l
Dsal ⊗	C/CRYGG GACNNNN/NNGTC	Btgl Drdl	R0608 R0530	C/CRYGG	BstDSI, Btgl
DseDI E	GACINININ/ININGTO	Diai	กบออบ	GACNNNN/NNGTC	Aasl, Drdl
Eael	Y/GGCCR	Eael	R0508	Y/GGCCR	Acol
Eagl	C/GGCCG	Eagl-HF	R3505	C/GGCCG	BseX3I, BstZI, Eagl, Eagl-HF, EciXI, Eco52I
Eam1104I	CTCTTC(1/4)	Earl	R0528	CTCTTC(1/4)	Bst6l, Earl
Eam1105I	GACNNN/NNGTC	Ahdl	R0584	GACNNN/NNGTC	Ahdl, BmeRl, Dril
Earl	CTCTTC(1/4)	Earl	R0528	CTCTTC(1/4)	Bst6I, Eam1104I
Ecil	GGCGGA(11/9)	Ecil	R0590	GGCGGA(11/9)	DS(O), Laitt 1041
Ecl136II	GAG/CTC	Eco53kl	R0116	GAG/CTC	EcolCRI, Eco53kl, Psp124Bl^, Sacl^, Sacl-HF^, Sstl^
20.100.1	arta, o ro	SacI-HF [^]	R3156	GAGCT/C	200.011, 200.0011, 1 0p. 2.12. 1 000. 1 1 1 000.
Ecl234I⊗	CGGNAAG				
Ecl35734I ⊗	GAAAYTC				
EclXI	C/GGCCG	Eagl-HF	R3505	C/GGCCG	BseX3I, BstZI, Eagl, Eagl-HF, Eco52I
Eco24I	GRGCY/C	Banll	R0119	GRGCY/C	Banll, EcoT38I, FriOI
Eco31I	GGTCTC(1/5)	Bsal-HFv2	R3733	GGTCTC(1/5)	Bsal, Bsal-HFv2, Bso31I, BspTNI
Eco32I	GAT/ATC	EcoRV-HF	R3195	GAT/ATC	EcoRV, EcoRV-HF
Eco47I	G/GWCC	Avall	R0153	G/GWCC	Avall, Bme18l, Sinl, VpaK11Bl
Eco47III	AGC/GCT	Afel	R0652	AGC/GCT	Afel, Aor51HI
Eco52I	C/GGCCG	Eagl-HF	R3505	C/GGCCG	BseX3I, BstZI, EagI, EagI-HF, EcIXI
Eco53KI	GAG/CTC	Eco53kl	R0116	GAG/CTC	Ecl136II, EcolCRI, Psp124BI^, SacI^, SacI-HF^, SstI^
		SacI-HF [^]	R3156	GAGCT/C	
Eco57I	CTGAAG(16/14)	Acul	R0641	CTGAAG(16/14)	Acul
Eco72I	CAC/GTG	PmII	R0532	CAC/GTG	Acvl, BbrPI, PmaCl, PmII, PspCl
Eco81I	CC/TNAGG	Bsu36I	R0524	CC/TNAGG	Axyl, Bse21l, Bsu36l
Eco88I	C/YCGRG	Aval	R0152	C/YCGRG	Ama87I, Aval, BmeT110I, BsiHKCI, BsoBI
		BsoBl	R0586	C/YCGRG	
Eco91I	G/GTNACC	BstEII-HF	R3162	G/GTNACC	BstEII, BstEII-HF, BstPI, EcoO65I, PspEI
Eco105I	TAC/GTA	SnaBl	R0130	TAC/GTA	BstSNI, SnaBI
Eco130I	C/CWWGG	Styl-HF	R3500	C/CWWGG	BssT1I, EcoT14I, Erhl, Styl, Styl-HF
Eco147I	AGG/CCT	Stul	R0187	AGG/CCT	Pcel, SseBl, Stul
Eco4465II ⊗	GAAABCC				
Eco43896Ⅱ ⊗	CRARCAG				
EcoBLMcrX ⊗	RCSRC(-3/-2)				
EcoE1140I⊗	ACCYAC				
EcolCRI	GAG/CTC	Eco53kl	R0116	GAG/CTC	Ecl136II, Eco53kI, Psp124BI^, SacI^, SacI-HF^, SstI^
		SacI-HF [^]	R3156	GAGCT/C	
Eco57MI⊗	CTGRAG(16/14)				
EcoMVII ⊗	CANCATC				

		NEB			
ENZYME	SEQUENCE	ENZYME	NEB#	SEQUENCE	OTHER ISOSCHIZOMERS
EcoNI	CCTNN/NNNAGG	EcoNI	R0521	CCTNN/NNNAGG	BstENI, Xagl
EcoNIH6II⊗	ATGAAG				
EcoO65I	G/GTNACC	BstEII-HF	R3162	G/GTNACC	BstEII, BstEII-HF, BstPI, Eco911, PspEI
Eco0109I	RG/GNCCY	Eco0109I	R0503	RG/GNCCY	E DI E DI IIE
EcoRI	G/AATTC	EcoRI-HF	R3101	G/AATTC	EcoRI, EcoRI-HF
EcoRII	/CCWGG	BstNI^	R0168	CC/WGG	Ajnl, BciT130I^, BseBI^, BstNI^, Bst2UI^, MvaI^, Psp6I, PspGI
FDV	O AT /ATO	PspGI	R0611	/CCWGG	F001 FDV FDV HE
EcoRV	GAT/ATC	EcoRV-HF	R3195	GAT/ATC	Eco32I, EcoRV, EcoRV-HF
EcoT14I	C/CWWGG	Styl-HF	R3500	C/CWWGG	BssT1I, Eco130I, Erhl, Styl, Styl-HF
EcoT22I	ATGCA/T	Nsil-HF	R3127	ATGCA/T GRGCY/C	Mph1103l, Nsil, Nsil-HF, Zsp2l
EcoT38I Eco53kI	GRGCY/C GAG/CTC	BanII Eco53kI	R0119 R0116	GAG/CTC	Banll, Eco24l, FriOI Ecl136ll, EcolCRI, Psp124BI^, SacI^, SacI-HF^, SstI^
ECOSSKI	dAd/CTC	SacI-HF [^]	R3156	GAGCT/C	EUT3011, EUT0111, FSp124bin, SdUn, SdUnin, SdUnin
Egel	GGC/GCC	Kasl^	R0544	G/GCGCC	Dinl, Ehel, Kasl^, Mly113l^^, Narl^^, PluTl^^^, Sfol, SspDI^
Lyci	440/400	Narl^^	R0191	GG/CGCC	טווו, בוופו, מפט , ועווץ דוטו , מפור , דועדו , טוטו, טשטו
		PluTI^^^	R0713	GGCGC/C	
		Sfol	R0606	GGC/GCC	
Ehel	GGC/GCC	Kasl^	R0544	G/GCGCC	Dinl, Egel, Kası^, Mly113I^^, Narl^^, PluTI^^^, Sfol, SspDI^
Liloi	440/400	Narl^^	R0191	GG/CGCC	Dilli, Egol, Nasi , Mily 1101 , Nati , 1 latti , 0101, 03pD1
		PluTI^^^	R0713	GGCGC/C	
		Sfol	R0606	GGC/GCC	
Eli8509II ⊗	CCGGAG	0.01		000,000	
Erhl	C/CWWGG	Styl-HF	R3500	C/CWWGG	BssT1I, Eco130I, EcoT14I, Styl, Styl-HF
EsaSSI ⊗	GACCAC	3tj. 1ti	110000	0,0111100	200111, 2001001, 2001111, 00,1, 00,1
Espl ⊗	GC/TNAGC	Blpl	R0585	GC/TNAGC	Blpl, Bpu1102l, Bsp1720l
Esp3l	CGTCTC(1/5)	BsmBI	R0580	CGTCTC(1/5)	BsmBl
		Esp3I	R0734	CGTCTC(1/5)	
Esp3007I ⊗	CAGAAG	·		()	
Exi27195I⊗	GCCGAC				
F					
Fael	CATG/	CviAII^	R0640	C/ATG	CviAll^, Fatl^^, Hin1ll, Hsp92ll, Nlalll
		FatI^^	R0650	/CATG	
		NIaIII	R0125	CATG/	
Fail	YA/TR				
Fall	(8/13)AAGNNNNNCTT(13/8)				
Faql	GGGAC(10/14)	BsmFI	R0572	GGGAC(10/14)	BsIFI, BsmFI
FatI	/CATG	CviAII^	R0640	C/ATG	CviAII^, FaeI^^, Hin1II^^, Hsp92II^^, NIaIII^^
		Fatl	R0650	/CATG	
		NlallI^^	R0125	CATG/	
Faul	CCCGC(4/6)	Faul	R0651	CCCGC(4/6)	
FauNDI	CA/TATG	Ndel	R0111	CA/TATG	Ndel
Fbal	T/GATCA	BcII-HF	R3160	T/GATCA	BcII, BcII-HF, Ksp22I
FbII	GT/MKAC	Accl	R0161	GT/MKAC	Accl, Xmil
Fco1691IV ⊗	GCVGAG	D 51	D. 6	00010/10/10	
Finl⊗	GGGAC	BsmFI	R0572	GGGAC(10/14)	BsIFI, BsmFI, Faql
FnuDII ⊗	CG/CG	BstUI	R0518	CG/CG	Accil, Bsh1236i, BspFNI, BstFNI, BstUI, MvnI
Fnu4HI	GC/NGC	Fnu4HI	R0178	GC/NGC	Fsp4HI, SatI
Fokl	GGATG(9/13)	BtsCI^	R0647	GGATG(2/0)	BseGI^, BstF5I^, BtsCI^
EriOI	CDCCV/C	Fokl	R0109	GGATG(9/13)	Panii Econdi EcoToni
FriOI	GRGCY/C	Banll	R0119	GRGCY/C	Banll, Eco24l, EcoT38l
Fsel	GGCCGG/CC	Fsel	R0588 R0135	GGCCGG/CC	Rigl Acc16I, Nsbl
Fspl	TGC/GCA RTGC/GCAY	Fspl	MUIJO	TGC/GCA	AUGTUI, NOUI
FspAl		Pfol	DUECO	C/TAG	Rfal Mael SchMl York
FspBI FspEI	C/TAG CC(12/16)	Bfal FspEl	R0568 R0662	CC(12/16)	Bfal, Mael, SspMl, Xspl
Fsp4HI	GC/NGC	Fnu4HI	R0178	GC/NGC	Fnu4HI, Satl
FspPK15I⊗	GARGAAG	i iiu4fii	110170	au/Nau	i nuti II, Jali
FtnUV ⊗	GAAACA				
G	UNNON				
GauT27I⊗	CGCGCAGG				
Gba708II⊗	ATGCAC				
Gdill⊗	CGGCCR(-5/-1)				
Glal	GC/GC				
Glul	GC/NGC				Bisl, Blsl^, Pkrl^
Gsal	CCCAGC(-1/-5)	BseYI^	R0635	CCCAGC(-5/-1)	BseYI^, PspFI^
Gsul	CTGGAG(16/14)	Bpml	R0565	CTGGAG(16/14)	Bpml
	0.00.00(10/11/	2piiii	.10000	5 . GG/10(10/17)	=

ENZYME	SEQUENCE	NEB Enzyme	NEB#	SEQUENCE	OTHER ISOSCHIZOMERS
H Hael ⊗	WGG/CCW				
Haell	RGCGC/Y	Haell	R0107	RGCGC/Y	Bfol, BstH2I
Haelll	GG/CC	Haelll	R0108	GG/CC	BshFI, BsnI, BspANI, BsuRI
HaeIV⊗	(7/13)GAYNNNNNRTC(14/9)	Haom	110100	du/00	Both I, Both, Bop IIII, Boutt
	C/CGG	Hpall	R0171	C/CGG	BsiSI, Hpall, Mspl
Hapll	6/666	пран Mspl	R0106	C/CGG	osioi, ripali, ivispi
Hball ⊗	GCCCAG				
HdeNY26I⊗	CGANNNNNTCC				
HdeZA17I⊗	GCANNNNNTCC				
Hgal	GACGC(5/10)	Hgal	R0154	GACGC(5/10)	Csel
HgiAl ⊗	GWGCW/C	BsiHKAI	R0570	GWGCW/C	Alw21I, Bbv12I, BsiHKAI
lgiCl⊗	G/GYRCC	Banl	R0118	G/GYRCC	AccB1I, Banl, BshNI, BspT107I
HgiEll ⊗	ACCNNNNNNGGT	Dum	110110	0/011100	7,000 H, Daill, Dollin, Dop 11071
-	GRGCY/C	Donll	R0119	GRGCY/C	Banll, Eco24l, EcoT38l, FriOl
HgiJII ⊗		Banll			
Hhal	GCG/C	Hhal HinP1I^	R0139 R0124	GCG/C G/CGC	AspLEI, BstHHI, CfoI, Hin6I^, HinP1I^, HspAI^
Hin1I	GR/CGYC	BsaHl	R0556	GR/CGYC	Acyl, BsaHl, BssNl, BstACl, Hsp92l
Hin1II	CATG/	CviAII^	R0640	C/ATG	CviAll^, Fael, Fatl^^, Hsp92II, NIaIII
	- ,	FatI^^	R0650	/CATG	, , , , ,
		NIaIII	R0125	CATG/	
lin 41 🛇	(0/12)\CV\NININININI\/T\C/12\0\	INIDINI	110120	UATU/	
Hin4l⊗	(8/13)GAYNNNNVTC(13/8)	Llov AV	DOCO4	CCTTC/C/E	UnvAV
Hin4II ⊗	CCTTC(6/5)	HpyAV	R0621	CCTTC(6/5)	HpyAV
Hin6I	G/CGC	Hhal^	R0139	GCG/C	AspLEI^, BstHHI^, CfoI^, HhaI^, HinP1I, HspAI
		HinP1I	R0124	G/CGC	
HinP1I	G/CGC	Hhal^	R0139	GCG/C	AspLEI^, BstHHI^, CfoI^, HhaI^, Hin6I, HspAI
		HinP1I	R0124	G/CGC	
HincII	GTY/RAC	HincII	R0103	GTY/RAC	Hindll
Hindll	GTY/RAC	HincII	R0103	GTY/RAC	HincII
HindIII	A/AGCTT	HindIII-HF	R3104	A/AGCTT	HindIII, HindIII-HF
Hinfl	G/ANTC	Hinfl	R0155	G/ANTC	
Hpal	GTT/AAC	Hpal	R0105	GTT/AAC	KspAl
- -Ipall	C/CGG	Hpall	R0171	C/CGG	BsiSI, HapII, MspI
	5,555	Mspl	R0106	C/CGG	,
Hphl	GGTGA(8/7)	HphI	R0158	GGTGA(8/7)	AsuHPI
Hpy8I	GTN/NAC	Hpy166II	R0616	GTN/NAC	Hpy166II
Hpy99I	CGWCG/	Hpy99I	R0615	CGWCG/	прутооп
	GCCTA	Пруээг	nuuio	Cawca/	
Hpy99XIII ⊗					
Hpy99XIV⊗	GGWTAA				
Hpy99XIV-mut1 ⊗	GGWCNA				
Hpy99XXII⊗	TCANNNNNTRG				
Hpy166II	GTN/NAC	Hpy166II	R0616	GTN/NAC	Hpy8l
Hpy178III ⊗	TC/NNGA	Hpy188III	R0622	TC/NNGA	Hpy188III
Hpy188I	TCN/GA	Hpy188I	R0617	TCN/GA	
Hpy188III	TC/NNGA	Hpy188III	R0622	TC/NNGA	
Hpy300XI ⊗	CCTYNA				
HpyAV	CCTTC(6/5)	HpyAV	R0621	CCTTC(6/5)	
lpyAXIV ⊗	GCGTA		110021	00110(0,0)	
HpyAXVI-mut1 ⊗	CRTTAA				
-	CRTCNA				
HpyAXVI-mut2 ⊗		HpyCH4III	R0618	ACN/GT	Bst4CI, Taal
HpyCH4III	ACN/GT	- 17			
HpyCH4IV	A/CGT	HpyCH4IV	R0619	A/CGT	HpySE526I, MaeII, TaiI^
HpyCH4V	TG/CA	HpyCH4V	R0620	TG/CA	D IDEL D I I
HpyF3I	C/TNAG	Ddel	R0175	C/TNAG	BstDEI, Ddel
lpyF10VI	GCNNNNN/NNGC	Mwol	R0573	GCNNNNN/NNGC	BstMWI, Mwol
HpySE526I	A/CGT	HpyCH4IV	R0619	A/CGT	HpyCH4IV, MaeII, TaiI^
HpyUM032XIII⊗	CYANNNNNNTRG				
HpyUM032XIII-mut1⊗	CYANNNNNNTTC				
HpyUM032XIV ⊗	GAAAG				
HpyUM037X⊗	GTGGNAG, TNGGNAG				
Hsp92I	GR/CGYC	BsaHI	R0556	GR/CGYC	Acyl, BsaHl, BssNl, BstACl, Hin1I
Hsp92II	CATG/	CviAII^	R0640	C/ATG	CviAII^, Fael, FatI^^, Hin1II, NIaIII
		FatI^^	R0650	/CATG	
				,	

ENZYME	SEQUENCE	NEB Enzyme	NEB#	SEQUENCE	OTHER ISOSCHIZOMERS
HspAl	G/CGC	Hhal^ HinP1I	R0139 R0124	GCG/C G/CGC	AspLEI [^] , BstHHI [^] , CfoI [^] , HhaI [^] , Hin6I, HinP1I
J					
Jma19592I ⊗	GTATNAC				
Jma19592II ⊗	GRGCRAC				
Jsp2502II ⊗ K	GRNGAAT				
Kasl	G/GCGCC	Kasl	R0544	G/GCGCC	Dinl^, Egel^, Ehel^, Mly113l^^, Narl^^, PluTl^^^, Sfol^, SspDl
		Narl^^	R0191	GG/CGCC	, , , , , , , , , , , , , , , , , , , ,
		PluTI^^^	R0713	GGCGC/C	
		Sfol^	R0606	GGC/GCC	
KfII	GG/GWCCC				
Kor51II ⊗	RTCGAG				
KpnI	GGTAC/C	Acc65I^	R0599	G/GTACC	Acc65I^, Asp718I^, KpnI, KpnI-HF
		Kpnl-HF	R3142	GGTAC/C	
Kpn2l	T/CCGGA	BspEl	R0540	T/CCGGA	AccIII, Aor13HI, BseAI, Bsp13I, BspEI, MroI
Kpn156V ⊗	CRTGATT				
KpnNH25III ⊗	CTRGAG				
KpnNlH30III ⊗	GTTCNAC				
KpnNlH50l ⊗	GCYAAG G/CCGGC				
Krol Kspl	CCGC/GG	SacII	R0157	CCGC/GG	Cfr421, SacII, Sfr3031, SgrBI
Ksp22I	T/GATCA	Sacii BcII-HF	R3160	T/GATCA	Bell, Bell-HF, Fbal
Ksp632I ⊗	CTCTTC(1/4)	Earl	R0528	CTCTTC(1/4)	Bst6I, Eam1104I, Earl
KspAl	GTT/AAC	Hpal	R0105	GTT/AAC	Hpal
Kzo9l	/GATC	DpnII	R0543	/GATC	Bsp143I, BssMI, BstKTI^, BstMBI, DpnII, Mbol, Ndell, Sau3Al
		Mbol	R0147	/GATC	· , · , · , · , · , · , · , · , · , · ,
		Sau3Al	R0169	/GATC	
L					
Lba2029III ⊗	CYAAANG				
Lde4408II⊗	ACAAAG				
Lgul	GCTCTTC(1/4)	BspQl Sapl	R0712 R0569	GCTCTTC(1/4) GCTCTTC(1/4)	BspQI, PciSI, SapI
LlaG50I⊗	CCGTKA				
Lmnl	GCTCC(1/-1)				
Lmo370I⊗	AGCGCCG				
Lmo911II⊗	TAGRAG				
Lpl1004II ⊗	AGGRAG	LanDI	DOCCO	0000(40/44)	
LpnPl Lra68l ⊗	CCDG(10/14) GTTCNAG	LpnPl	R0663	CCDG(10/14)	
Liaooi ⊗ LsaDS4l ⊗	TGGAAT				
Lsp48III ⊗	AGCACC				
Lsp1109I	GCAGC(8/12)	Bbvl	R0173	GCAGC(8/12)	Bbvl, BseXl, BstV1I
Lsp6406VI ⊗	CRAGCAC	5511		a 07 (a 0 (a) 12)	5511, 53511, 551111
Lwel	GCATC(5/9)	SfaNI	R0172	GCATC(5/9)	Bmsl, SfaNl
M					
Mabl	A/CCWGGT	SexAl	R0605	A/CCWGGT	Csil, SexAl
Mael	C/TAG	Bfal	R0568	C/TAG	Bfal, FspBl, SspMl, Xspl
Maell	A/CGT	HpyCH4IV	R0619	A/CGT	HpyCH4IV, HpySE526I, TaiI^
MaeIII	/GTNAC	5 /	Pa. 1 = -	0.4.550	
Mall	GA/TC	Dpnl	R0176	GA/TC	Dpnl
Maql ⊗	CRTTGAC(21/19)				
MauBl	CG/CGCGCG				
Mba11I⊗ MbiI	AGGCGA CCGCTC(-3/-3)	BsrBI	R0102	CCGCTC(-3/-3)	AccBSI, BsrBI
Mbol	/GATC	DpnII	R0102 R0543	/GATC	Bsp143I, BssMI, BstKTI^, BstMBI, DpnII, Kzo9I, NdeII, Sau3AI
MIDUI	Janu	Mbol	R0147	/GATC	סטף ויוסו, סטטויוו, סטנויווין עירווון, וענטטו, ווענטוו, טמטטאו
		Sau3Al	R0169	/GATC	
Mboll	GAAGA(8/7)	Mboll	R0148	GAAGA(8/7)	
Mcrl ⊗	CGRY/CG	BsiEl	R0554	CGRY/CG	Bsh1285I, BsiEI, BstMCI
Mcr10I⊗	GAAGNNNNNCTC				
Mfel	C/AATTG	Mfel-HF	R3589	C/AATTG	Munl, Mfel, Mfel-HF
MfII	R/GATCY	BstYI	R0523	R/GATCY	BstX2I, BstYI, Psul
Mhll	GDGCH/C	Bsp1286I	R0120	GDGCH/C	Bsp1286I, Sdul
MjaIV ⊗	GTNNAC	Hpy166II	R0616	GTN/NAC	Hpy8I, Hpy166II
MkaDII ⊗	GAGAYGT				
MIsI	TGG/CCA	Mscl	R0534	TGG/CCA	Ball, MluNl, Mox20l, Mscl, Msp20l

ENZVME	SECHENCE	NEB Enzyme	NED#	SECTIENCE	OTHER ISOSCHIZOMERS
ENZYME	SEQUENCE	Mlul-HF	NEB #	SEQUENCE	OTHER ISOSCHIZOMERS
/lul	A/CGCGT		R3198	A/CGCGT	Mlul, Mlul-HF
MuCl	/AATT	MluCl	R0538	/AATT	Sse9I, Tasl
/IluNI	TGG/CCA	Mscl	R0534	TGG/CCA	Ball, Mlsl, Mox20l, Mscl, Msp20l
Vilyl	GAGTC(5/5)	Mlyl	R0610	GAGTC(5/5)	Plel^, Ppsl^, Schl
		Plel^	R0515	GAGTC(4/5)	
Mly113I	GG/CGCC	KasI^	R0544	G/GCGCC	DinI^^, EgeI^^, EheI^^, KasI^, NarI, PluTI^^^, SfoI^^, SspDI^
		Narl	R0191	GG/CGCC	
		PluTl^^^	R0713	GGCGC/C	
		Sfol^^	R0606	GGC/GCC	
Mmel	TCCRAC(20/18)	Mmel	R0637	TCCRAC(20/18)	
Viniei	CCTC(7/6)	MnII	R0163	CCTC(7/6)	
Mox20I	TGG/CCA	Mscl	R0534	TGG/CCA	Ball, Misi, MiuNi, Msci, Msp201
	ATGCA/T		R3127	ATGCA/T	
Mph1103I		Nsil-HF	N3121	AIGUA/I	EcoT22I, Nsil, Nsil-HF, Zsp2I
Virel	CG/CCGGCG				
Mrol	T/CCGGA	BspEl	R0540	T/CCGGA	Accill, Aor13Hi, BseAl, Bsp13i, BspEl, Kpn2i
MroNI	G/CCGGC	Nael^	R0190	GCC/GGC	Nael^, NgoMIV, Pdil^
		NgoMIV	R0564	G/CCGGC	
MroXI	GAANN/NNTTC	XmnI	R0194	GAANN/NNTTC	Asp700I, PdmI, XmnI
Viscl	TGG/CCA	Mscl	R0534	TGG/CCA	Ball, Mlsl, MluNl, Mox20l, Msp20l
/Isel	T/TAA	Msel	R0525	T/TAA	SagAl, Tru11, Tru91
Visil	CAYNN/NNRTG	MsII	R0571	CAYNN/NNRTG	Rsel, SmiMI
Visil	C/CGG	Hpall	R0171	C/CGG	BsiSI, Hapll, Hpall
siohi	0/000				υσιοι, παριι, πραιι
Man 201	TOC/OCA	Mspl	R0106	C/CGG	Dall Mial Michil MacON M1
VIsp20I	TGG/CCA	Mscl	R0534	TGG/CCA	Ball, Misl, MiuNi, Mox20l, Mscl
MspA1I	CMG/CKG	MspA1I	R0577	CMG/CKG	
MspCI	C/TTAAG	AfIII	R0520	C/TTAAG	Afili, Bfri, BspTi, BstAFi, Vha464i
VIspI7II ⊗	ACGRAG				
MspJI	CNNR(9/13)	MspJI	R0661	CNNR(9/13)	
MspR9I	CC/NGG	ScrFI	R0110	CC/NGG	Bme1390I, BmrFI, BstSCI [^] , ScrFI, StyD4I [^]
		StyD4I^	R0638	/CCNGG	
MspSC27II ⊗	CCGCGAC				
Mssl	GTTT/AAAC	Pmel	R0560	GTTT/AAAC	Pmel
	TGC/GCA	Fspl	R0135	TGC/GCA	
VIstI ⊗		rspi	HU135	TGU/GUA	Acc16I, FspI, NsbI
Mtel	GCGC/NGCGC				
MtuHN878II ⊗	CACGCAG				
Munl	C/AATTG	Mfel-HF	R3589	C/AATTG	Mfel, Mfel-HF
Mval	CC/WGG	BstNI	R0168	CC/WGG	AjnI^, BciT130I, BseBI, BstNI, Bst2UI, EcoRII^, Psp6I^, PspGI^
		PspGI [^]	R0611	/CCWGG	
Mva1269I	GAATGC(1/-1)	Bsml	R0134	GAATGC(1/-1)	Bsml, Pctl
MvnI	CG/CG	BstUI	R0518	CG/CG	Accil, Bsh1236i, BspFNi, BstFNi, BstUi
Mwol	GCNNNN/NNGC	Mwol	R0573	GCNNNNN/NNGC	BstMWI, HpyF10VI
V	CONTINUE	IVIWOI	110070	domininini	Battivii, ripyr 10 vi
	CCC/CCC	MacI	D0100	GCC/GCC	MroNIA NgoMIVA Pdil
Nael	GCC/GGC	Nael	R0190	GCC/GGC	MroNI^, NgoMIV^, Pdil
	100100	NgoMIV^	R0564	G/CCGGC	
Val45188II ⊗	ACCAGC				
VarI	GG/CGCC	KasI^	R0544	G/GCGCC	Dinl^^, Egel^^, Ehel^^, Kasl^, Mly113I, PluTl^^^, Sfol^^, SspDl/
		Narl	R0191	GG/CGCC	
		PluTl^^^	R0713	GGCGC/C	
		Sfol^^	R0606	GGC/GCC	
Nbr128II⊗	ACCGAC				
Vcil	CC/SGG	Ncil	R0196	CC/SGG	AsuC2I, BcnI, BpuMI
Vcol	C/CATGG	Ncol-HF	R3193	C/CATGG	Bsp19l, Ncol, Ncol-HF
Vdel	CA/TATG	Ndel	R0111	CA/TATG	FauNDI
Ndell	/GATC	DpnII	R0543	/GATC	Bsp143I, BssMI, BstKTI^, BstMBI, DpnII, Kzo9I, MboI, Sau3AI
		Mbol	R0147	/GATC	
		Sau3Al	R0169	/GATC	
NgoAVII ⊗	GCCGC(7/7)				
•					
NgoAVIII ⊗	(12/14)GACNNNNNTGA(13/11)	Nael^	R0190	GCC/GGC	MroNl. Nael^. Pdil^
NgoAVIII ⊗		Nael^ NaoMIV	R0190 R0564	GCC/GGC	MroNI, Nael^, Pdil^
NgoAVIII ⊗ NgoMIV	(12/14)GACNNNNNTGA(13/11) G/CCGGC	Nael^ NgoMIV	R0190 R0564	GCC/GGC G/CCGGC	MroNI, Nael^, Pdil^
NgoAVII ⊗ NgoAVIII ⊗ NgoMIV	(12/14)GACNNNNNTGA(13/11) G/CCGGC CAAGRAG	NgoMIV	R0564	G/CCGGC	
NgoAVIII ⊗ NgoMIV	(12/14)GACNNNNNTGA(13/11) G/CCGGC				MroNI, Nael^, Pdil^ AsuNHI, Bmtl^, Bmtl-HF^, BspOl^, Nhel, Nhel-HF

		NEB			
ENZYME	SEQUENCE	ENZYME	NEB#	SEQUENCE	OTHER ISOSCHIZOMERS
NIalli	CATG/	CviAII^	R0640	C/ATG	CviAll^, Fael, Fatl^^, Hin1ll, Hsp92ll
		FatI^^	R0650	/CATG	,, , , , , , , , , , , , , , ,
		NIaIII	R0125	CATG/	
NIaIV	GGN/NCC	NIaIV	R0126	GGN/NCC	Bmil, BspLl, PspN4l
NIaCI ⊗	CATCAC(19/17)	111411	110120	danifino	511111, 500-21, 1 001111
NmeAIII	GCCGAG(21/19)	NmeAIII	R0711	GCCGAG(21/19)	
NmuCl	/GTSAC	Tsp45I	R0583	/GTSAC	TseFI, Tsp45I
Notl	GC/GGCCGC	NotI-HF	R3189	GC/GGCCGC	CciNI, NotI, NotI-HF
	GATCGAC	11011-111	110103	do/ddoodd	GGINI, NOU, NOU-I II
NpeUS61II ⊗	TCG/CGA	Nrul-HF	R3192	TCG/CGA	Bsp681, BtuM1, Nrul, Nrul-HF, Rrul
Nrul					
Nsbl	TGC/GCA	Fspl	R0135	TGC/GCA	Acc16l, Fspl
Nsil	ATGCA/T	Nsil-HF	R3127	ATGCA/T	EcoT22I, Nsil, Nsil-HF, Mph1103I, Zsp2I
Nspl	RCATG/Y	Nspl	R0602	RCATG/Y	BstNSI, Xcel
NspV	TT/CGAA	BstBI	R0519	TT/CGAA	Asull, Bpu14l, Bsp119l, BspT104l, BstBl, Sful
NspBII ⊗	CMG/CKG	MspA1I	R0577	CMG/CKG	MspA1I
0					
Olil	CACNN/NNGTG	Alel-v2	R0685	CACNN/NNGTG	Alel, Alel-v2
OspHL35III ⊗	YAGGAG				
P					
Pacl	TTAAT/TAA	Pacl	R0547	TTAAT/TAA	
PacIII ⊗	GTAATC				
Pac19842II ⊗	CCTTGA				
Pael	GCATG/C	SphI-HF	R3182	GCATG/C	Sphl, Sphl-HF
PaeR7I	C/TCGAG	PaeR7I	R0177	C/TCGAG	Sfr274I, Slal, Xhol
		Xhol	R0146	C/TCGAG	
Pagl	T/CATGA	BspHI	R0517	T/CATGA	BspHI, Ccil
Pal408I ⊗	CCRTGAG			.,	
PalAl	GG/CGCGCC	Ascl	R0558	GG/CGCGCC	Ascl, Sgsl
Pasl	CC/CWGGG	71001	110000	ad, ododoo	7.001, 0g01
Paul	G/CGCGC	BssHII	R0199	G/CGCGC	BsePI, BssHII, Ptel
Pba2294I ⊗	GTAAG	DSSLIII	110133	a/ododo	D361 1, D331 III, 1 (61
_					
Pcall ⊗	GACGAG	Chul	D0107	ACC/CCT	Food 471, Coopl, Chul
Pcel	AGG/CCT	Stul	R0187	AGG/CCT	Eco147I, SseBI, Stul
Pcil	A/CATGT	Pcil	R0655	A/CATGT	Pscl
PciSI	GCTCTTC(1/4)	BspQI	R0712	GCTCTTC(1/4)	BspQl, Lgul, Sapl
-		Sapl	R0569	GCTCTTC(1/4)	
Pcr308II ⊗	CCAAAG				
Pcsl	WCGNNNN/NNNCGW				
PctI	GAATGC(1/-1)	Bsml	R0134	GAATGC(1/-1)	Bsml, Mva1269I
Pdil	GCC/GGC	Nael	R0190	GCC/GGC	MroNI^, Nael, NgoMIV^
		NgoMIV^	R0564	G/CCGGC	
Pdi8503III ⊗	CCGGNAG				
Pdml	GAANN/NNTTC	XmnI	R0194	GAANN/NNTTC	Asp700I, MroXI, XmnI
Pdu1735l ⊗	CACCAC				
Penl ⊗	GCAGT				
Pfel	G/AWTC	Tfil	R0546	G/AWTC	Tfil
Pfl23II	C/GTACG	BsiWI-HF	R3553	C/GTACG	BsiWI, BsiWI-HF, PspLI
Pfl1108l ⊗	TCGTAG				
Pfl8569I ⊗	GCN/NGC				
PfIFI	GACN/NNGTC	PfIFI	R0595	GACN/NNGTC	Psyl, Tth111I
1 1111 1	artory reverse	Tth1111	R0185	GACN/NNGTC	1 3/1, 14111111
PfIMI	CCANNNN/NTGG	PfIMI	R0509	CCANNNN/NTGG	AccB7I, Van91I
PfIPt14I ⊗	RGCCCAC	1 HIVII	110000	OUAINININ/NTUU	1, vano 11
Pfol Pfol	T/CCNGGA				
PfrJS12IV ⊗	TANAAG				
PfrJS12V ⊗	GGCGGAG				
PfrJS15III ⊗	CTTCNAC		DOES-	A /0.000=	A 1 A 115 A 101 B 171 C 11
PinAl	A/CCGGT	Agel-HF	R3552	A/CCGGT	Agel, Agel-HF, AsiGl, BshTl, CspAl
Pin17FIII ⊗	GGYGAB				
	CTRKCAG				
PinP23II ⊗					
PinP59III ⊗	GAAGNAG				
	GAAGNAG GCN/GC				Bisl^, BIsl, Glul^
PinP59III ⊗		Mlyl^	R0610	GAGTC(5/5)	Bisl^, BIsl, Glul^ Mlyl^, Ppsl, Schl^
PinP59III ⊗ PkrI	GCN/GC	Miyi^ Plel	R0610 R0515	GAGTC(5/5) GAGTC(4/5)	
PinP59III ⊗ PkrI	GCN/GC				

Isoschizomers (continued)

FN7VME	SEQUENCE	NEB ENZYME	NER#	SECTIONS	OTHER ISOSCHIZOMERS
ENZYME PluTl	GGCGC/C	ENZYME Kasl^	NEB # R0544	SEQUENCE G/GCGCC	OTHER ISOSCHIZOMERS Dinl^^, Egel^^, Ehel^^, Kasl^, Mly113l^^^, Narl^^^, Sfol^^, SspDl/
PIUII	GGCGC/C				Diffirm, Egerm, Efferm, Nasin, Mily 1131mm, Nationn, Stolm, Sspur
		Narl^^^	R0191	GG/CGCC	
		PluTl	R0713	GGCGC/C	
		Sfol^^	R0606	GGC/GCC	
PmaCl	CAC/GTG	PmII	R0532	CAC/GTG	Acvl, BbrPl, Eco72I, PmII, PspCl
Pmel	GTTT/AAAC	Pmel	R0560	GTTT/AAAC	MssI
PmII	CAC/GTG	PmII	R0532	CAC/GTG	Acvl, BbrPl, Eco72l, PmaCl, PspCl
Ppil ⊗	(7/12)GAACNNNNNCTC(13/8)				
PpiP13II ⊗	CGCRGAC				
Ppsl	GAGTC(4/5)	Mlyl^	R0610	GAGTC(5/5)	Mlyl^, Plel, Schl^
		Plel	R0515	GAGTC(4/5)	
Ppu21I	YAC/GTR	BsaAl	R0531	YAC/GTR	BsaAl, BstBAl
PpuMI	RG/GWCCY	PpuMI	R0506	RG/GWCCY	Psp5II, PspPPI
Pscl	A/CATGT	Pcil	R0655	A/CATGT	Pcil
Pse18267I ⊗	RCCGAAG	1 611	110033	Noniai	1 011
PshAl		DobAl	DOEO2	C A CNINI /NINICTO	Dovl. DotDA1
	GACNN/NNGTC	PshAl	R0593	GACNN/NNGTC	Boxl, BstPAI
PshBl	AT/TAAT	Asel	R0526	AT/TAAT	Asel, Vspl
Psil	TTA/TAA	Psil	R0657	TTA/TAA	Aanl
Psp5II	RG/GWCCY	PpuMI	R0506	RG/GWCCY	PpuMI, PspPPI
Psp6I	/CCWGG	BstNI^	R0168	CC/WGG	Ajnl, BciT130I^, BseBI^, BstNI^, Bst2UI^, EcoRII, MvaI^, PspGI
		PspGI	R0611	/CCWGG	
Psp0357II ⊗	GCGAAG				
Psp1406I	AA/CGTT	AcII	R0598	AA/CGTT	AcII
Psp124BI	GAGCT/C	Eco53kI^	R0116	GAG/CTC	Ecl136II^, EcolCRI^, Eco53kI^, Sacl, Sacl-HF, Sstl
. 00 12 131	a. 100 1,70	Sacl-HF	R3156	GAGCT/C	20170011
PspCI	CAC/GTG	PmII	R0532	CAC/GTG	Acvl, BbrPl, Eco72I, PmaCl, PmII
PspEl	G/GTNACC	BstEII-HF	R3162	G/GTNACC	BstEII, BstEII-HF, BstPI, Eco91I, Eco065I
PspFI	CCCAGC(-5/-1)	BseYI	R0635	CCCAGC(-5/-1)	BseYI, Gsal^
PspGI	/CCWGG	BstNI^	R0168	CC/WGG	Ajnl, BciT130I^, BseBI^, BstNI^, Bst2UI^, EcoRII, MvaI^, Psp6I
		PspGI	R0611	/CCWGG	
PspLI	C/GTACG	BsiWI-HF	R3553	C/GTACG	BsiWI, BsiWI-HF, Pfl23II
PspN4I	GGN/NCC	NIaIV	R0126	GGN/NCC	Bmil, BspLl, NlaIV
Psp0MI	G/GGCCC	Apal^	R0114	GGGCC/C	Apal^, Bsp120I
		PspOMI	R0653	G/GGCCC	
Psp0MII⊗	CGCCCAR(20/18)				
PspPI	G/GNCC	Sau96I	R0165	G/GNCC	AspS9I, BmgT120I, Cfr13I, Sau96I
PspPPI	RG/GWCCY	PpuMI	R0506	RG/GWCCY	PpuMI, Psp5II
PspPRI ⊗	CCYCAG(15/13)	i paivii	110000	na, amoor	r paint, r opon
PspXI	VC/TCGAGB	PspXI	R0656	VC/TCGAGB	
•		гэрлі	nuusu	VC/TCUAUD	
PsrI	(7/12)GAACNNNNNNTAC(12/7)	5	50110	07001/0	D 144 D 4 D 4 15
PstI	CTGCA/G	PstI-HF	R3140	CTGCA/G	BspMAI, PstI, PstI-HF
Pst145I⊗	CTAMRAG				
Pst273I⊗	GATCGAG				
Pst14472I ⊗	CNYACAC				
PstNI	CAGNNN/CTG	AlwNI	R0514	CAGNNN/CTG	AlwNI, Cail
Psul	R/GATCY	BstYI	R0523	R/GATCY	BstX2I, BstYI, MfII
PsuGI ⊗	BBCGD				
Psyl	GACN/NNGTC	PfIFI	R0595	GACN/NNGTC	PfIFI, Tth111I
i oyi	G/101y/111UTO	Tth111I	R0185	GACN/NNGTC	. 111 1, 10111111
Dtol	CICCCC				ReaDI RecHII Paul
Ptel	G/CGCGC	BssHII	R0199	G/CGCGC	BsePI, BssHII, Paul
Pvul	CGAT/CG	Pvul-HF	R3150	CGAT/CG	Ple19I, Pvul, Pvul-HF
Pvull	CAG/CTG	PvuII-HF	R3151	CAG/CTG	Pvull, Pvull-HF
R					
Rba2021I ⊗	CACGAGH				
Rcel ⊗	CATCGAC(20/18)				
RdeGBI⊗	CCGCAG				
RdeGBII ⊗	ACCCAG(20/18)				
RdeGBIII ⊗	(9/11)TGRYCA(11/9)				
RfIFIII ⊗	CGCCAG				
Rgal	GCGAT/CGC	AsiSI	R0630	GCGAT/CGC	AsiSI, SfaAI, Sgfl
Rigl	GGCCGG/CC	Fsel	R0588	GGCCGG/CC	Fsel
Rlal ⊗	VCW				
RIaII ⊗	ACACAG(20/18)				
	CCCACA(12/9)				
RleAl ⊗	GGGAGA(12/9)				

ENZYME	SEQUENCE	NEB Enzyme	NEB#	SEQUENCE	OTHER ISOSCHIZOMERS
	GTYGGAG(11/9)	ENZTIVIE	NED#	SEQUENCE	OTHER ISOSCHIZOWERS
Rpal ⊗	. ,				
RpaBI⊗	CCCGCAG(20/18)				
RpaB5l⊗	CGRGGAC(20/18)				
RpaTI⊗	GRTGGAG		50.400	T00/004	0.001.01.01.01.01.01.01
Rrul	TCG/CGA	Nrul-HF	R3192	TCG/CGA	Bsp68I, BtuMI, Nrul, Nrul-HF
Rsal	GT/AC	CviQI^	R0639	G/TAC	Afal, Csp6I^, CviQI^, RsaNI^
		Rsal	R0167	GT/AC	
RsaNI	G/TAC	CviQI	R0639	G/TAC	Afal^, Csp6l, CviQl, Rsal^
		Rsal^	R0167	GT/AC	
Rsel	CAYNN/NNRTG	MsII	R0571	CAYNN/NNRTG	MsII, SmiMI
Rsp008IV ⊗	ACGCAG				
Rsp008V ⊗	GCCCAT				
RspPBTS2III ⊗	CTTCGAG				
RsrII	CG/GWCCG	RsrII	R0501	CG/GWCCG	Cpol, Cspl, Rsr2l
Rsr2I	CG/GWCCG	RsrII	R0501	CG/GWCCG	Cpol, Cspl, Rsrll
Rtr1953I ⊗	TGANNNNNTGA				- F - 1 - 2 F 1 - 2
S					
Sacl	GAGCT/C	Eco53kI^	R0116	GAG/CTC	Eci136II^, EcolCRI^, Eco53ki^, Psp124Bi, Saci, Saci-HF, Ssti
Jaci	dAd01/0	SacI-HF	R3156	GAGCT/C	E013011 , E0010111 , E0033N1 , 1 3p124b1, 3d01, 3d01 111, 33t1
Cooli	0000/00				Cfr401 Vanl Cfr0001 CarDI
Sacil	CCGC/GG	SacII	R0157	CCGC/GG	Cfr42I, KspI, Sfr303I, SgrBI
Saf8902III ⊗	CAATNAG	0.11.115	D0400	0./T0040	Call Call LIE
Sall	G/TCGAC	Sall-HF	R3138	G/TCGAC	Sall, Sall-HF
SanDI⊗	GG/GWCCC		_		KfII
Sapl	GCTCTTC(1/4)	BspQI	R0712	GCTCTTC(1/4)	BspQI, LguI, PciSI
		Sapl	R0569	GCTCTTC(1/4)	
SaqAl	T/TAA	Msel	R0525	T/TAA	Msel, Tru1l, Tru9l
SatI	GC/NGC	Fnu4HI	R0178	GC/NGC	Fnu4HI, Fsp4HI
Saul ⊗	CC/TNAGG	Bsu36l	R0524	CC/TNAGG	Axyl, Bse21I, Bsu36I, Eco81I
Sau96I	G/GNCC	Sau96I	R0165	G/GNCC	AspS9I, BmgT120I, Cfr13I, PspPI
Sau3AI	/GATC	DpnII	R0543	/GATC	Bsp143I, BssMI, BstKTI^, BstMBI, DpnII, Kzo9I, MboI, Ndell
		Mbol	R0147	/GATC	
		Sau3AI	R0169	/GATC	
Sba460II⊗	GGNGAYG			7	
Sbfl	CCTGCA/GG	SbfI-HF	R3642	CCTGCA/GG	Sbfl, Sbfl-HF, Sdal, Sse83871
Sb046I⊗	TGAAC	ODII III	113042	oo raan,aa	obii, obii fii, oddi, oscootfi
-		Scal-HF	R3122	AGT/ACT	DmoAl Cool Cool LIE 7rml
Scal	AGT/ACT				BmcAl, Scal, Scal-HF, Zrml
Schl	GAGTC(5/5)	Mlyl	R0610	GAGTC(5/5)	Mlyl, Plel^, Ppsl^
0 00011 0	007447	PleI^	R0515	GAGTC(4/5)	
ScoDS2II ⊗	GCTAAT	0 5	50110	00.000	D. 10001 D. El D.10011 M. DOL O. D.111
ScrFl	CC/NGG	ScrFI	R0110	CC/NGG	Bme1390I, BmrFI, BstSCI [^] , MspR9I, StyD4I [^]
		StyD4I^	R0638	/CCNGG	
Sdal	CCTGCA/GG	Sbfl-HF	R3642	CCTGCA/GG	Sbfl, Sbfl-HF, Sse8387I
SdeAl ⊗	CAGRAG(21/19)				
SdeOSI⊗	(11/13)GACNNNNRTGA(12/10)				
Sdul	GDGCH/C	Bsp1286I	R0120	GDGCH/C	Bsp1286I, MhII
Secl ⊗	C/CNNGG	BsaJI	R0536	C/CNNGG	BsaJI, BseDI, BssECI
Sen17963III ⊗	CCAAAC				
SenA1673III ⊗	GNGGCAG				
SenSARA26III ⊗	ACRCAG				
SenTFIV ⊗	GATCAG				
Setl	ASST/				
SexAl	A/CCWGGT	SexAl	R0605	A/CCWGGT	Csil, Mabl
SfaAl	GCGAT/CGC	AsiSI	R0630	GCGAT/CGC	AsiSI, Rgal, Sgfl
SfaNI	GCATC(5/9)	SfaNI	R0172	GCATC(5/9)	Bmsl, Lwel
Sfcl			R0561		Bfml, BstSFI
	C/TRYAG	Sfcl		C/TRYAG	,
Sfel ⊗	C/TRYAG	SfcI	R0561	C/TRYAG	Bfml, BstSFl, Sfcl
Sfil	GGCCNNNN/NGGCC	Sfil	R0123	GGCCNNNN/NGGCC	
Sfol	GGC/GCC	Kasl^	R0544	G/GCGCC	Dinl, Egel, Ehel, Kası [^] , Mly113l [^] , Narı [^] , PluTl ^{^^} , SspDl [^]
		Narl^^	R0191	GG/CGCC	
		PluTI^^^	R0713	GGCGC/C	
		Sfol	R0606	GGC/GCC	
Sfr274I	C/TCGAG	PaeR7I	R0177	C/TCGAG	PaeR7I, Slal, Xhol
		Xhol	R0146	C/TCGAG	
Sfr303I	CCGC/GG	SacII	R0157	CCGC/GG	Cfr42I, Kspl, SacII, SgrBI
Sful	TT/CGAA	BstBI	R0519	TT/CGAA	Asull, Bpu14I, Bsp119I, BspT104I, BstBI, NspV
Sgel	CNNGNNNNNNNN/				
Sgfl	GCGAT/CGC	AsiSI	R0630	GCGAT/CGC	AsiSI, Rgal, SfaAl
Syll	UUUAI/UUU	ASISI	นกดุวก	GUGAI/UGU	MOIOI, NYAI, OIAMI

Isoschizomers (continued)

		NEB			
ENZYME	SEQUENCE	ENZYME	NEB #	SEQUENCE	OTHER ISOSCHIZOMERS
SgrAl	CR/CCGGYG	SgrAl	R0603	CR/CCGGYG	
SgrBI	CCGC/GG	SacII	R0157	CCGC/GG	Cfr42I, KspI, SacII, Sfr303I
SgrDI	CG/TCGACG				
SgrTl⊗	CCDS(10/14)				
Sgsl	GG/CGCGCC	Ascl	R0558	GG/CGCGCC	Ascl, PalAl
		ASU	ทบววช	dd/CdCdCC	ASUI, FaiMi
Siml ⊗	GGGTC(-3/0)				
Sinl	G/GWCC	Avall	R0153	G/GWCC	Avall, Bme18l, Eco47l, VpaK11Bl
Slal	C/TCGAG	PaeR7I	R0177	C/TCGAG	PaeR7I, Sfr274I, XhoI
		Xhol	R0146	C/TCGAG	
Smal	CCC/GGG	Smal	R0141	CCC/GGG	Cfr9I^, TspMI^, XmaI^
	000,000	TspMI^	R0709	C/CCGGG	,
		· ·			
		Xmal^	R0180	C/CCGGG	
SmaUMH5I⊗	CTTGAC				
SmaUMH8I ⊗	GCGAACB				
Smil	ATTT/AAAT	Swal	R0604	ATTT/AAAT	Swal
SmiMI	CAYNN/NNRTG	MsII	R0571	CAYNN/NNRTG	MsII, Rsel
SmII	C/TYRAG	SmII	R0597	C/TYRAG	Smol
Smol	C/TYRAG	SmII	R0597	C/TYRAG	SmII
Snal ⊗	GTATAC	BstZ17I-HF	R3594	GTA/TAC	BssNAI, Bst1107I, BstZ17I, BstZ17I-HF
SnaBl	TAC/GTA	SnaBl	R0130	TAC/GTA	BstSNI, Eco105I
Sno506l⊗	GGCCGAG				
Spel	A/CTAGT	Spel-HF	R3133	A/CTAGT	Ahll, Bcul, Spel, Spel-HF
Sphl	GCATG/C	SphI-HF	R3182	GCATG/C	Pael, Sphl, Sphl-HF
SpII ⊗	C/GTACG	BsiWI-HF	R3553	C/GTACG	BsiWl, BsiWl-HF, Pfl23ll, PspLl
		III-IAAIen	110000	ojuinou	Dolivi, Dolivi III, I IIZOII, I apli
SpnRII ⊗	TCGAG				
SpoDI⊗	GCGGRAG				
Srfl	GCCC/GGGC	Srfl	R0629	GCCC/GGGC	
Sse9I	/AATT	MluCl	R0538	/AATT	MluCl, Tasl
Sse232I ⊗	CG/CCGGCG				Mrel
Sse8387I	CCTGCA/GG	Sbfl-HF	R3642	CCTGCA/GG	Sbfl, Sbfl-HF, Sdal
Sse8647I ⊗	AG/GWCCT	ODII III	1100-12	oordorydd	obii, obii Tii, oddi
		0.1	50107	100/007	5 447 5 4 6 4
SseBI	AGG/CCT	Stul	R0187	AGG/CCT	Eco147I, Pcel, Stul
Ssil	CCGC(-3/-1)	Acil	R0551	CCGC(-3/-1)	Acil, BspACI
Sspl	AAT/ATT	SspI-HF	R3132	AAT/ATT	Sspl, Sspl-HF
Ssp714II⊗	CGCAGCG				
Ssp6803IV ⊗	GAAGGC				
SspDI	G/GCGCC	Kasl	R0544	G/GCGCC	Dinl^, Egel^, Ehel^, Kasl, Mly113l^^, Narl^^, PluTl^^^, Sfol^
ουμοι	d/dCdCC				Dilli , Lyci , Lifei , Masi, Mily 1131 , Maii , 11011 , 3101
		Narl^^	R0191	GG/CGCC	
		PluTI^^^	R0713	GGCGC/C	
		Sfol^	R0606	GGC/GCC	
SspMI	C/TAG	Bfal	R0568	C/TAG	Bfal, FspBl, Mael, Xspl
Sstl	GAGCT/C	Eco53kI^	R0116	GAG/CTC	Ecl136II^, EcolCRI^, Eco53kI^, Psp124BI, SacI, SacI-HF
	u	SacI-HF	R3156	GAGCT/C	
CalF07I 🔿	0044040(00(40)	Saul-TIF	113130	UAUU1/U	
SstE37I ⊗	CGAAGAC(20/18)				
Sth132I⊗	CCCG(4/8)				
Sth20745III ⊗	GGACGAC				
SthSt3II ⊗	GAAGT				
Stul	AGG/CCT	Stul	R0187	AGG/CCT	Eco147I, Pcel, SseBI
Styl	C/CWWGG	Styl-HF	R3500	C/CWWGG	BssT1I, Eco130I, EcoT14I, Erhl, Styl, Styl-HF
StyD4I	/CCNGG	ScrFI [^]	R0110	CC/NGG	Bme1390I^, BmrFI^, BstSCI, MspR9I^, ScrFI^
		StyD4I	R0638	/CCNGG	
SurP32all ⊗	ACRGAG				
Swal	ATTT/AAAT	Swal	R0604	ATTT/AAAT	Smil
T.					
Taal	ACN/GT	HpyCH4III	R0618	ACN/GT	Bst4Cl, HpyCH4III
Tail	ACGT/	HpyCH4IV [^]	R0619	A/CGT	HpyCH4IV [^] , HpySE526I [^] , MaeII [^]
Taql	T/CGA	Taql	R0149	T/CGA	
TaqII	GACCGA(11/9)				
TagIII ⊗	CACCCA(11/9)				
Tasl	/AATT	MluCl	R0538	/AATT	MluCl, Sse9l
iuoi		IVIIUOI	110000	/AATT	Wildelf, Octob
Total	W/GTACW				
	GCSG/C				
Taul	GCSG/C G/AWTC	Tfil	R0546	G/AWTC	Pfel
Tatl Taul Tfil Tru11		Tfil Msel	R0546 R0525	G/AWTC T/TAA	Pfel Msel, SaqAl, Tru9l

		NEB			
ENZYME	SEQUENCE	ENZYME	NEB #	SEQUENCE	OTHER ISOSCHIZOMERS
TscAl	CASTGNN/	TspRI	R0582	CASTGNN/	TspRI
Tsel	G/CWGC	ApeKI	R0643	G/CWGC	ApeKI
1301	d/owdo	Tsel	R0591	G/CWGC	προιτί
TseFI	/GTSAC	Tsp45I	R0583	/GTSAC	NmuCl, Tsp45l
Tsol ⊗	TARCCA(11/9)	тэртэт	110000	JUTONO	Minuol, 13p4ol
-	/GTSAC	Top 4E1	R0583	/GTSAC	New CL Tool
Tsp45I		Tsp45I	KU383	/G15AC	NmuCI, TseFI
TspARh3I ⊗	GRACGAC		50010	AONIOT	D 1401 II OHAW T I
Tsp4Cl⊗	ACN/GT	HpyCH4III	R0618	ACN/GT	Bst4CI, HpyCH4III, Taal
TspDTI	ATGAA(11/9)				
TspEl ⊗	/AATT	MluCl	R0538	/AATT	MluCl, Sse9l, Tasl
TspGWI	ACGGA(11/9)				
TspMI	C/CCGGG	Smal^	R0141	CCC/GGG	Cfr9I, SmaI [^] , XmaI
		TspMI	R0709	C/CCGGG	
		Xmal	R0180	C/CCGGG	
TspRI	CASTGNN/	TspRI	R0582	CASTGNN/	TscAl
Tssl⊗	GAGNNNCTC				
Tstl ⊗	(8/13)CACNNNNNNTCC(12/7)				
Tsul ⊗	GCGAC				
Tth1111	GACN/NNGTC	PfIFI	R0595	GACN/NNGTC	PfIFI, Psyl
		Tth111I	R0185	GACN/NNGTC	,,
Tth111Ⅱ ⊗	CAARCA(11/9)		.10100	anonymior o	
U	0,0,0,0,0,0,0				
UbaF9I⊗	TACNNNNNRTGT				
UbaF11I⊗	TCGTA				
-	CTACNNNGTC				
UbaF12I⊗					
UbaF13l ⊗	GAGNNNNNCTGG				
UbaF14I ⊗	CCANNNNTCG				
UbaPl⊗	CGAACG				
V					
Van91I	CCANNNN/NTGG	PfIMI	R0509	CCANNNN/NTGG	AccB7I, PfIMI
Van9116I⊗	CCKAAG				
Vdi96II⊗	GNCYTAG				
Vha464I	C/TTAAG	AfIII	R0520	C/TTAAG	AfIII, BfrI, BspTI, BstAFI, MspCI
Vnel	G/TGCAC	ApaLI	R0507	G/TGCAC	Alw44I, ApaLI
VpaK11BI	G/GWCC	Avall	R0153	G/GWCC	Avall, Bme18I, Eco47I, SinI
Vspl	AT/TAAT	Asel	R0526	AT/TAAT	Asel, PshBI
Vtu19109I ⊗	CACRAYC				
W					
Wvil⊗	CACRAG(21/19)				
χ					
Xagl	CCTNN/NNNAGG	EcoNI	R0521	CCTNN/NNNAGG	BstENI, EcoNI
Xapl	R/AATTY	Apol-HF	R3566	R/AATTY	Acsl, Apol, Apol-HF
Xbal	T/CTAGA	Xbal	R0145	T/CTAGA	7,001, 7,001, 7,001 711
Xca85IV ⊗		Abui	110110	1/01/14/1	
ACCOUNTY (S)					
	TACGAG	Menl	R0602	RCATC/V	RetNCI Neni
Xcel	RCATG/Y	Nspl	R0602	RCATG/Y	BstNSI, NspI
Xcel Xcml	RCATG/Y CCANNNNN/NNNNTGG	Xcml	R0533	CCANNNNN/NNNNTGG	
Xcel	RCATG/Y	Xcml PaeR7I	R0533 R0177	CCANNNNN/NNNNTGG C/TCGAG	BstNSI, Nspl PaeR7I, Sfr274I, Slal
Xcel Xcml Xhol	RCATG/Y CCANNNNN/NNNNTGG C/TCGAG	Xcml PaeR7I Xhol	R0533 R0177 R0146	CCANNNNN/NNNNTGG C/TCGAG C/TCGAG	PaeR7I, Sfr274I, Slal
Xcel Xcml Xhol Xholl ⊗	RCATG/Y CCANNNNN/NNNNTGG C/TCGAG R/GATCY	Xcml PaeR7I Xhol BstYI	R0533 R0177 R0146 R0523	CCANNNNN/NNNNTGG C/TCGAG C/TCGAG R/GATCY	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfII, Psul
Xcel Xcml Xhol	RCATG/Y CCANNNNN/NNNNTGG C/TCGAG	Xcml PaeR7I XhoI BstYI Smal^	R0533 R0177 R0146 R0523 R0141	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG	PaeR7I, Sfr274I, Slal
Xcel Xcml Xhol Xholl ⊗	RCATG/Y CCANNNNN/NNNNTGG C/TCGAG R/GATCY	Xcml PaeR7I Xhol BstYI Smal^ TspMI	R0533 R0177 R0146 R0523 R0141 R0709	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG C/CCGGG	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfII, Psul
Xcel Xcml Xhol Xholl ⊗	RCATG/Y CCANNNNN/NNNNTGG C/TCGAG R/GATCY C/CCGGG	Xcml PaeR7I XhoI BstYI Smal^ TspMI Xmal	R0533 R0177 R0146 R0523 R0141 R0709 R0180	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG C/CCGGG C/CCGGG	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfII, Psul
Xcel Xcml Xhol Xholl ⊗	RCATG/Y CCANNNNN/NNNNTGG C/TCGAG R/GATCY	Xcml PaeR7I Xhol BstYI Smal^ TspMI	R0533 R0177 R0146 R0523 R0141 R0709	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG C/CCGGG	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfII, Psul
Xcel Xcml Xhol Xholl ⊗ Xmal	RCATG/Y CCANNNNN/NNNNTGG C/TCGAG R/GATCY C/CCGGG	Xcml PaeR7I XhoI BstYI Smal^ TspMI Xmal	R0533 R0177 R0146 R0523 R0141 R0709 R0180	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG C/CCGGG C/CCGGG	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfII, Psul Cfr9I, Smal^, TspMI
Xcel Xcml Xhol Xholl ⊗ Xmal Xmall ⊗	RCATG/Y CCANNNNN/NNNNTGG C/TCGAG R/GATCY C/CCGGG C/GGCCG	Xcml PaeR7I Xhol BstYI Smal^ TspMI Xmal Eagl-HF	R0533 R0177 R0146 R0523 R0141 R0709 R0180 R3505	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG C/CCGGG C/CCGGG C/CCGGG	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfII, Psul Cfr9I, Smal^, TspMI BseX3I, BstZI, Eagl, Eagl-HF, EclXI, Eco52I
Xcel Xcml Xhol Xholl ⊗ Xmal Xmall ⊗ XmaJl	RCATG/Y CCANNNNN/NNNNTGG C/TCGAG R/GATCY C/CCGGG C/GGCCG C/CTAGG	Xcml PaeR7I Xhol BstYI Smal^ TspMI Xmal Eagl-HF AvrII	R0533 R0177 R0146 R0523 R0141 R0709 R0180 R3505 R0174	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG C/CCGGG C/CCGGG C/CCGGG C/CCGGG	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfII, Psul Cfr9I, Smal^, TspMI BseX3I, BstZI, Eagl, Eagl-HF, EclXI, Eco52I AspA2I, AvrII, BinI
Xcel Xcml Xhol Xholl Xmall Xmall Xmall Xmall Xmil	RCATG/Y CCANNNNN/NNNNTGG C/TCGAG R/GATCY C/CCGGG C/GGCCG C/CTAGG GT/MKAC	Xcml PaeR7I Xhol BstYI Smal^ TspMI Xmal Eagl-HF AvrII Accl	R0533 R0177 R0146 R0523 R0141 R0709 R0180 R3505 R0174 R0161	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG C/CCGGG C/CCGGG C/CCGGG C/CTAGG GT/MKAC	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfII, Psul Cfr9I, Smal^, TspMI BseX3I, BstZI, Eagl, Eagl-HF, EclXI, Eco52I AspA2I, AvrII, BlnI Accl, FbII
Xcel Xcml Xhol Xhol ⊗ Xmal Xmall ⊗ XmaJl Xmil Xmnl	RCATG/Y CCANNNNN/NNNTGG C/TCGAG R/GATCY C/CCGGG C/GGCCG C/CTAGG GT/MKAC GAANN/NNTTC	Xcml PaeR7I Xhol BstYI Smal^ TspMI Xmal Eagl-HF AvrII Accl XmnI	R0533 R0177 R0146 R0523 R0141 R0709 R0180 R3505 R0174 R0161 R0194	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG C/CCGGG C/CCGGG C/CCGGG C/CTAGG GT/MKAC GAANN/NNTTC	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfII, Psul Cfr9I, Smal^, TspMI BseX3I, BstZI, Eagl, Eagl-HF, EclXI, Eco52I AspA2I, AvrII, BInI Accl, FbII Asp700I, MroXI, PdmI
Xcel Xcml Xhol Xholl ⊗ Xmal Xmall ⊗ XmaJl Xmil Xmnl Xspl Y	RCATG/Y CCANNNNN/NNNTGG C/TCGAG R/GATCY C/CCGGG C/GGCCG C/CTAGG GT/MKAC GAANN/NNTTC C/TAG	Xcml PaeR7I Xhol BstYI Smal^ TspMI Xmal Eagl-HF AvrII Accl XmnI	R0533 R0177 R0146 R0523 R0141 R0709 R0180 R3505 R0174 R0161 R0194	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG C/CCGGG C/CCGGG C/CCGGG C/CTAGG GT/MKAC GAANN/NNTTC	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfII, Psul Cfr9I, Smal^, TspMI BseX3I, BstZI, Eagl, Eagl-HF, EclXI, Eco52I AspA2I, AvrII, BInI Accl, FbII Asp700I, MroXI, PdmI
Xcel Xcml Xhol Xholl ⊗ Xmal Xmall ⊗ XmaJl Xmil Xmnl Xspl Y Yps36061 ⊗	RCATG/Y CCANNNNN/NNNTGG C/TCGAG R/GATCY C/CCGGG C/GGCCG C/CTAGG GT/MKAC GAANN/NNTTC	Xcml PaeR7I Xhol BstYI Smal^ TspMI Xmal Eagl-HF AvrII Accl XmnI	R0533 R0177 R0146 R0523 R0141 R0709 R0180 R3505 R0174 R0161 R0194	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG C/CCGGG C/CCGGG C/CCGGG C/CTAGG GT/MKAC GAANN/NNTTC	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfII, Psul Cfr9I, Smal^, TspMI BseX3I, BstZI, Eagl, Eagl-HF, EclXI, Eco52I AspA2I, AvrII, BInI Accl, FbII Asp700I, MroXI, PdmI
Xcel Xcml Xhol Xhol ⊗ Xmal Xmall ⊗ XmaJl Xmil Xmnl Xspl Y Yps36061 ⊗ Z	RCATG/Y CCANNNNN/NNNTGG C/TCGAG R/GATCY C/CCGGG C/GGCCG C/CTAGG GT/MKAC GAANN/NNTTC C/TAG CGGAAG	Xcml PaeR7I Xhol BstYI Smal^ TspMI Xmal Eagl-HF AvrII Accl XmnI Bfal	R0533 R0177 R0146 R0523 R0141 R0709 R0180 R3505 R0174 R0161 R0194 R0568	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG C/CCGGG C/CCGGG C/CCGGG C/CTAGG GT/MKAC GAANN/NNTTC C/TAG	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfiI, Psul Cfr9I, Smal^, TspMI BseX3I, BstZI, EagI, EagI-HF, EcIXI, Eco52I AspA2I, AvrII, BInI AccI, FbII Asp700I, MroXI, PdmI BfaI, FspBI, Mael, SspMI
Xcel Xcml Xhol Xholl ⊗ Xmal Xmall ⊗ XmaJl Xmil Xmnl Xspl Y Yps36061 ⊗	RCATG/Y CCANNNNN/NNNTGG C/TCGAG R/GATCY C/CCGGG C/GGCCG C/CTAGG GT/MKAC GAANN/NNTTC C/TAG	Xcml PaeR7I Xhol BstYI Smal^ TspMI Xmal Eagl-HF AvrII Accl Xmnl Bfal	R0533 R0177 R0146 R0523 R0141 R0709 R0180 R3505 R0174 R0161 R0194 R0568	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG C/CCGGG C/CCGGG C/CCGGG C/CTAGG GT/MKAC GAANN/NNTTC C/TAG	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfII, Psul Cfr9I, Smal^, TspMI BseX3I, BstZI, Eagl, Eagl-HF, EclXI, Eco52I AspA2I, AvrII, BInI Accl, FbII Asp700I, MroXI, PdmI
Xcel Xcml Xhol Xhol Xholl ⊗ Xmal Xmall ⊗ XmaJl Xmil Xmnl Xspl Y Yps36061 ⊗ Z Zral	RCATG/Y CCANNNNN/NNNTGG C/TCGAG R/GATCY C/CCGGG C/GGCCG C/CTAGG GT/MKAC GAANN/NNTTC C/TAG CGGAAG GAC/GTC	Xcml PaeR7I Xhol BstYI Smal^ TspMI Xmal Eagl-HF AvrII Accl Xmnl Bfal Aatll^ Zral	R0533 R0177 R0146 R0523 R0141 R0709 R0180 R3505 R0174 R0161 R0194 R0568	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG C/CCGGG C/CCGGG C/CCGGG C/CTAGG GT/MKAC GAANN/NNTTC C/TAG GACGT/C GAC/GTC	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfiI, Psul Cfr9I, Smal^, TspMI BseX3I, BstZI, EagI, EagI-HF, EclXI, Eco52I AspA2I, AvrII, BlnI AccI, FbII Asp700I, MroXI, PdmI BfaI, FspBI, MaeI, SspMI
Xcel Xcml Xhol Xhol ⊗ Xmal Xmall ⊗ XmaJl Xmil Xmnl Xspl Y Yps36061 ⊗ Z	RCATG/Y CCANNNNN/NNNTGG C/TCGAG R/GATCY C/CCGGG C/GGCCG C/CTAGG GT/MKAC GAANN/NNTTC C/TAG CGGAAG	Xcml PaeR7I Xhol BstYI Smal^ TspMI Xmal Eagl-HF AvrII Accl Xmnl Bfal	R0533 R0177 R0146 R0523 R0141 R0709 R0180 R3505 R0174 R0161 R0194 R0568	CCANNNNN/NNNTGG C/TCGAG C/TCGAG R/GATCY CCC/GGG C/CCGGG C/CCGGG C/CCGGG C/CTAGG GT/MKAC GAANN/NNTTC C/TAG	PaeR7I, Sfr274I, Slal BstX2I, BstYI, MfiI, Psul Cfr9I, Smal^, TspMI BseX3I, BstZI, EagI, EagI-HF, EcIXI, Eco52I AspA2I, AvrII, BInI AccI, FbII Asp700I, MroXI, PdmI BfaI, FspBI, Mael, SspMI

Survival in a Reaction

Restriction enzymes vary with respect to their ability to maintain activity in a reaction over an extended period of time.

+++ Enzyme is active > 8 hours

N/A Not Available

Enzyme is active 4-8 hours

ENZYME

Enzyme is active 2-4 hours No benefit from digesting over 1 hour

SURVIVAL

While most routine restriction digests are incubated for one hour or less at 37°C , there are certain applications that require the addition of less than 1 unit/µg of DNA and increasing the reaction time beyond one hour. The table below can be used as a guide when low levels of enzyme and extended reaction times are needed.

ENZTIVIE	SUNVIVAL
AatII	+++
AbaSI @25°C	N/A
Accl	+++
Acc65I	+
Acil	_
AcII	
	+
Acul	-
Afel	++
AfIII	+++
AfIIII	+++
Agel	+
Agel-HF	++
Ahdl	+++
Alel-v2	N/A
Alul	+++
Alwi	+++
AlwNI	
	+++
Apal @25°C	+++
ApaLI	+++
ApeKI @75°C	_
Apol @50°C	+++
Apol-HF	+++
Ascl	_
Asel	+++
AsiSI	+++
Aval	+
Avall	+++
AvrII	+++
Bael @25°C	+
BaeGI	+++
BamHI	+
BamHI-HF	+
Banl	+++
Banll	+++
Bbsl	++
BbsI-HF	_
Bbvl	++
BbvCl	
	+++
Bccl	+
BceAl	+++
Bcgl	++
BciVI	++
BcII @50°C	+
BcII-HF	N/A
BcoDI	+++
Bfal	+
BfuAl @50°C	++
Bgll	+++
-	
Bglll	++
Blpl	+
BmgBl	-
Bmrl	++
Bmtl	+++
Bmtl-HF	+++
Bpml	-

ENZYME	SURVIVAL
Bpu10I	+
BpuEl @25°C	+
Bsal	-
Bsal-HFv2	++
BsaAl	++
BsaBl @60°C	+++
BsaHI	+++
BsaJI @60°C	+++
BsaWI @60°C	+++
BsaXI	++
BseRI	+
BseYI	+
Bsql	+
BsiEl @60°C	++
BsiHKAI @65°C	_
BsiWI @55°C	_
BsiWI-HF	+++
BsII @55°C	+
Bsml @65°C	+++
BsmAl @55°C	+++
BsmBl @55°C	+
BsmFl @65°C	+++
BsoBl	+++
Bsp1286I	+
BspCNI @25°C	т _
BspDI	_
BspEl	+++
BspHI	++
BspMI	+++
BspQI @50°C	+++
Bsrl @65°C	++
BsrBI	+
BsrDI @65°C	
BsrFI-v2	+
BsrGl	+++
	+++
BsrGI-HF BssHII @50°C	+++
BssSI-v2	++
BstAPI @60°C	+++
BstBI @65°C	++
BstEll @60°C	+++
BstEII-HF	_
	_
BstNI @60°C BstUI @60°C	_
	+++
BstXI @55°C BstYI @60°C	+
Bst71@60°C BstZ17I-HF	++
	++
Bsu36I	+++
Btgl	+
BtgZI @60°C	_
BtsI-v2 @55°C	+++
BtsMutl @55°C	N/A
BtsCI @50°C	+
Cac8I	+++

Clal

For example, 1 unit of Aatll can be used to digest 8 μ g of DNA in a 16 hour digest (+ + +).

Extended activity was determined by performing the restriction endonuclease unit assay, using a 16 hour incubation in place of the standard 1 hour digestion. After the 16 hour digestion, extended activity enzymes (+ + +) required only 0.13 units to completely digest 1 µg of DNA. Intermediate activity enzymes required either 0.25 (++) or 0.50 (+) units for complete digestion over this extended incubation time. Finally, enzymes marked (–) required 1.0 unit for complete digestion, the same amount of enzyme required for a 1 hour digestion.

ENZYME

SURVIVA

Note: Reaction temperature is 37°C, unless otherwise noted.

ENZYME	SURVIVAL
CspCI	+++
CviAII @25°C	-
CviKI-1	++
CviQI @25°C	++
Ddel	+++
Dpnl	+++
DpnII	+++
Dral	+
DrallI-HF	+++
Drdl	+++
Eael	+++
Eagl	+++
Eagl-HF	+++
Earl	++
Ecil	_
Eco53kl	++
EcoNI	+++
Eco0109I	+++
EcoP15I	_
EcoRI	+++
EcoRI-HF	++
EcoRV	+
EcoRV-HF	+++
Esp3I	N/A
Fatl @55°C	_
Faul @55°C	_
Fnu4HI	+++
Fokl	++
Fsel	_
Fspl	+++
FspEl	+++
Haell	
HaellI	+++
Hgal	-
Hhal	++
Hincll	+++
HindIII	+++
HindIII-HF	+++
Hinfl	+++
HinP1I	+++
Hpal	++
Hpall	+++
Hphl	+++
Hpy99I	-
Hpy166II	+
Hpy188I	+++
Hpy188III	++
HpyAV	-
HpyCH4III	+++
HpyCH4IV	+++
HpyCH4V	+++
I-Ceul	++
I-Scel	++
Kasl	-
Kpnl	++
[h	1 T

ENZYME	SURVIVAL
KpnI-HF	+
LpnPI	-
Mbol	+++
Mboll	+
Mfel	+ +
Mfel-HF	++
Mlul	+++
Mlul-HF	+++
MluCl	_
Mlyl	-
Mmel	_
MnII	++
Mscl	+++
Msel	+++
MsII	+ +
Mspl	+
MspA1I	+ +
MspJI	+++
Mwol @60°C	+++
Nael	_
Narl	+++
Nb.BbvCI	+++
Nb.Bsml @65°C	++
Nb.BsrDI @65°C	++
Nb.BssSI	+++
Nb.BtsI	++
Ncil	+
Ncol	+++
Ncol-HF	+++
Ndel	++
NgoMIV	++
Nhel	++
Nhel-HF	+++
NIalli	+
NIaIV	+
NmeAIII	_
Notl	++
Notl-HF	+++
Nrul	+++
Nrul-HF	+++
Nsil	+++
Nsil-HF	+++
Nspl	++
Nt.Alwl	+++
Nt.BbvCl	+++
Nt.BsmAl	+++
Nt.BspQI @50°C	++
Nt.BstNBI @55°C	+
Nt.CviPII	_
Pacl	+++
PaeR7I	+++
Pcil	+++
PfIFI	+++
PfIMI	
PI-PspI @65°C	+
1 1-L 2h1 @00 P	+++

Survival in a Reaction (continued)

ENZYME	SURVIVAL
PI-Scel	+++
Phol @75°C	-
Plel	+++
PluTl	N/A
Pmel	-
PmII	+
PpuMI	+++
PshAl	-
Psil	++
PspGI @75°C	++
Psp0MI	+++
PspXI	+++
Pstl	+
PstI-HF	+++
Pvul	+++
Pvul-HF	+++
Pvull	+++

ENZYME	SURVIVAL
Pvull-HF	+++
Rsal	++
RsrII	+++
Sacl	+++
SacI-HF	+++
SacII	+++
Sall	+++
Sall-HF	+++
Sapl	+++
Sau3AI	+
Sau96I	+++
Sbfl	+
Sbfl-HF	+++
Scal-HF	+ +
ScrFI	+++
SexAl	++
SfaNI	++

ENZYME	SURVIVAL
SfcI	-
Sfil @50°C	++
Sfol	-
SgrAl	+++
Smal @25°C	-
SmII @55°C	++
SnaBl	+
Spel	++
Spel-HF	N/A
Sphl	+
SphI-HF	+
Srfl	+++
Sspl	++
SspI-HF	+++
Stul	+++
Styl	++
Styl-HF	+++

ENZYME	SURVIVAL
StyD4I	+
Swal @25°C	+ +
Taqal @65°C	+
Tfil @65°C	++
Tlil @16°C	++
Tsel @65°C	+
Tsp45I @65°C	+
TspMI @75°C	+
TspRI @65°C	+++
Tth1111 @65°C	++
Xbal	+++
XcmI	+++
Xhol	+++
Xmal	+++
XmnI	++
Zral	-

Cleavage of Supercoiled DNA

Restriction enzymes cleave different DNA substrates with varying efficiency. Restriction enzymes were tested for their ability to cleave various plasmids (pBR322, pUC19 and pLITMUS) under standard reaction conditions. Single sites were tested on each of these plasmids, depending on availability, and average values were taken when there was more than one data point available. Lambda DNA was used as the standard (1 unit to cleave in all cases).

ENZYME	UNITS TO CLEAVE PLASMID
AatII	3
AfIIII	1
Ahdl	1
Accl	4
Acc65I	1
AfIII	2
Agel	1
AlwNI	2
Apal	1
Apol	1
Asel	0.3
Aval	10
AvrII	1
Bael	3
BamHI	3
Banll	1
BgIII	8
Bpml	1
Bsal	2
BsaAl	20
BsaXI	2
BsiWl	3
Bsgl	1
Bsml	1
BspDI	1
BspEl	1
BspMI	**
BspQI	3
BsrFI	2

**	Requires two copies of its recognition
	sequence for cleavage to occur.

ENZYME	UNITS TO CLEAVE PLASMID
BsrGl	1
BssHI	4
Btgl	5
Clal	5
Eagl	10
Eco0109I	8
EcoNI	3
EcoRI	3
EcoRV	1
HincII	4
HindIII	5
Kasl	4
Kpnl	2
Mlul	2
Narl	10
Ncol	1
Ndel	3
NgoMIV	2
Nhel	5
Nrul	1
Nsil	1
Pcil	3
Psil	3
PstI	1
Pvul	2
Pvull	2
Sacl	5
Sall	10
Sapl	1

ENZYME	UNITS TO CLEAVE PLASMID
Scal	20
Smal	1
SnaBl	1
Spel	1
Sphl	3
Sspl	4
Stul	3
Styl	4
Tlil	2
TspMI	1
Tth111I	2
Xbal	2
Xhol	10
XmnI	5

New restriction sites can be generated by ligation of DNA fragments with compatible ends. These ends may be generated by:

- 1. Cleavage followed by fill-in of 5' overhangs to generate blunt ends.
- 2. Cleavage with two restriction enzymes that produce blunt ends.
- 3. Cleavage with two restriction enzymes that produce compatible overhangs.

Compatible ends, generated by each of the above methods, can be ligated to produce DNA sequences that often contain useful restriction endonuclease sites. Generation of these sites is detailed in the following tables.

Recleavable Filled-in 5' Overhangs

The table below lists a selection of restriction enzymes that generate 5´ overhangs which, if filled-in and ligated, result in new restriction sites. The combinations listed were identified by computer analysis, and have not necessarily been confirmed by experimentation. For a more complete listing visit our website, www.neb.com.

Restriction enzymes that have degenerate recognition specificities (e.g., recognize more than one sequence) have been excluded from this list. Where isoschizomers exist, only one member of each set is listed. Only commercially available enzymes have been listed.

30	TIAA G5		3UTTAATTAAG5
ENZYME	CLEAVAGE SITE	AFTER FILL-IN/ LIGATION	RECLEAVED BY
Acc65I	G/GTACC	GGTACGTACC	BsaAI, HpyCH4IV, RsaI, SnaBI ⁶
Acil	C/CGC	CCGCGC	(Acil), BstUl, Hhal
AcII	AA/CGTT	AACGCGTT	AfIIII, BstUI, MIuI ⁶
AfIII	C/TTAAG	CTTAATTAAG	Msel ² , Pacl ⁸ , MluCl
Agel	A/CCGGT	ACCGGCCGGT	BsiEI, (BsrFI)2, EaeI, EagI 6, HaeIII, HpaII
ApaLI	G/TGCAC	GTGCATGCAC	Cac8l, Nlalli, Nspl, Sphl ⁶
Ascl	GG/CGCGCC	GGCGCGCGCCC	(BssHII) ² , BstUI, Cac8I, Hhal
AvrII	C/CTAGG	CCTAGCTAGG	Alul, (Bfal) ²
BamHI	G/GATCC	GGATCGATCC	Alwl, Clal ⁶ , (DpnII ²), Taql
Bcll	T/GATCA	TGATCGATCA	Clal ⁶ , (DpnII ²), TaqI
Bfal	C/TAG	CTATAG	SfcI
BgIII	A/GATCT	AGATCGATCT	Clal ⁶ , DpnII ² , Taql
BsiWl	C/GTACG	CGTACGTACG	BsaAI, (BsiWI), HpyCH4IV, RsaI, SnaBI ⁶
BspDI/Clal	AT/CGAT	ATCGCGAT	BstUI, NruI ⁶
BspEl	T/CCGGA	TCCGGCCGGA	BsiEl, Eael, Eagl ⁶ , Haelll, (Hpall) ²
BspHI	T/CATGA	TCATGCATGA	(NIaIII) ² , NsiI ⁶
BsrGI	T/GTACA	TGTACGTACA	BsaAl, (Rsal) ² , SnaBl ⁶
BssHII	G/CGCGC	GCGCGCGCGC	(BssHII), BstUI, Cac8I, Hhal
BstBI	TT/CGAA	TTCGCGAA	BstUI, NruI ⁶
DpnII/MboI/Sau3AI	/GATC	GATCGATC	Clal ⁶ , (DpnII), TaqI
Eagl	C/GGCCG	CGGCCGGCCG	BsiEl, BsrFl, Cac8l, Eael ² , (Eagl ²), Fsel ⁸ , Haelll, Hpall, Nael
EcoRI	G/AATTC	GAATTAATTC	Asel ⁶ , Msel, MluCl, Xmnl ⁶
Fatl	/CATG	CATGCATG	BrfBI ⁶ , HpyCH4V, (Fatl ²)
HinP1I	G/CGC	GCGCGC	BssHII ⁶ , BstUI, Cac8I, (HhaI)
HindIII	A/AGCTT	AAGCTAGCTT	Alul, Bfal, Cac8l, Nhel ⁶
Hpall/Mspl	C/CGG	CCGCGG	Acil, BsaJl, BstUl, Btgl, MspA1l, SacII ⁶
HpyCH4IV	A/CGT	ACGCGT	AfIIII, BstUI, Mlul ⁶
Kasl	G/GCGCC	GGCGCGCGCC	(BssHII6)2, (BstUI)2, Cac8I, (HhaI)2
Mfel	C/AATTG	CAATTAATTG	Asel ⁶ , MluCl ²
Mlul	A/CGCGT	ACGCGCGCGT	BssHII ⁶ , BstUI, Cac8I, (HhaI) ²
MluCl	/AATT	AATTAATT	Asel ⁶ , Msel, (MluCl ²)
Narl	GG/CGCC	GGCGCGCC	Ascl ⁸ , BssHII, BstUI, Cac8I, Hhal
Ncol	C/CATGG	CCATGCATGG	NIaIII, NsiI ⁶
NgoMIV	G/CCGGC	GCCGGCCGGC	BsiEl, BsrFl, Cac8l, Eael, Eagl ⁶ , Haelll, Hpall, (NgoMIV ²)
Nhel	G/CTAGC	GCTAGCTAGC	Alul, Bfal, Cac8l, (Nhel)
Notl	GC/GGCCGC	GCGGCCGGCCGC	Acil, BsiEl, BsrFl, Cac8l, Eael, (Eagl ⁶) ² , Fnu4Hl, Fsel ⁸ , Haelll, Hpall, Nael ⁶
PaeR7I/XhoI	C/TCGAG	CTCGATCGAG	BsiEl, Dpnll, Pvul ⁶ , (Tagl ²)
Pcil	A/CATGT	ACATGCATGT	HpyCH4V , (NIaIII ²), NsiI ⁶
Psp0MI	G/GGCCC	GGGCCGGCCC	BsrFI, Cac8I, FseI ⁸ , HaeIII, HpaII, NaeI, Sau96I
PspXI	VC/TCGAGB	VCTCGATCGAGB	Pvul ⁶ , (Taql ²)
Sall	G/TCGAC	GTCGATCGAC	BsiEI, DpnII, PvuI ⁶ , TaqI
Spel	A/CTAGT	ACTAGCTAGT	Alul, (Bfal) ²
Tagl	T/CGA	TCGCGA	BstUI, NruI ⁶
Xbal	T/CTAGA	TCTAGCTAGA	Alul, Bfal
Xmal	C/CCGGG	CCCGGCCGGG	BsiEl, Eael, Eagl ⁶ , Haelll, Hpall, Ncil, ScrFl
			,, , ,

Table Notes

Enzymes in **bold** have recognition sequences of 6 or 8 bases. Sequence length is indicated by superscript (e.g., **Ascl*** = 8-base cutter).

Enzymes in parentheses indicate that the new sequence is still a substrate for the original enzyme.

A superscript 2 indicates that two identical sites have been generated within the filled-in/ligated sequence. For example, fill-in/ligation of AfIII generates the sequence CTTAATTAAG which contains two Msel sites (TTAA).

APPENDI

Recleavable Blunt Ends

The table below lists a selection of blunt-end cutters that produce recleavable ligation products. The combinations listed were identified by computer analysis, and although we have tried to ensure their accuracy, they have not necessarily been confirmed by experimentation. For a more complete listing visit our website, www.neb.com

RECLEAVED BY			
AGC/GCT Bst2171	ENZYME	LIGATED TO	RECLEAVED BY
AG/CT Bst2171		BstZ17I EcoRV SfoI FspI Nael	Alul SfaNI Haell, Hhal Hhal Acil, Fnu4HI
TYAC/GTR		BstZ17I MspA1I (CMG/CTG), PvuII EcoRV	Rsal Alul Mbol
CG/CG BstZ171	(YAC/GTR) (CAC/GTR) (TAC/GTR)	PmII PmII, SnaBI	BsaAl, Pmll BsaAl
GTA/TAC MspA11 (CMG/CKG), Nrul, Pvull, Stul Atel Alul Atel Hincll (GTY/GAC) Accl Sspl MluCl		BstZ17I EcoRV SfoI	Rsal Mbol Haelli
Taql		MspA1I (CMG/CKG), NruI, PvuII, StuI AfeI HincII (GTY/GAC)	Alul Acci
GAT/ATC		Nrul	Taql
TGC/GCA)		MspA1I (CMG/CKG), Pvull HaellI, MscI, Stul Afel, SfoI, Fspl	Alwl, Mbol SfaNl
GG/CC BstZ17I Mscl, Stul HaellI EcoRV AlwI, Mbol Stol Hincil (GTY/GAC) BsmFl Hincil (GTY/GAC) Hincil (GTC/RAC) Hincil (GTC/RAC) BstZ17I Accl (GTC/RAC) Haelli, Mscl, Stul BsmFl (GTC/RAC) Haelli, Mscl, Stul BsmFl (GTC/RAC) Hincil (GTC/RAC)		Afel, Sfol EcoRV Nael	Hhal SfaNI Acil, Fnu4HI
GTT/AAC		BstZ171 Mscl, Stul EcoRV Sfol	Rsal Haelll Alwl, Mbol Haelll, Sau96l
(GTC/RAC) BsrBI, MspA1I (CMG/CGG) Hpall (GTC/RAC) BstZ171 Accl (GTC/RAC) Haelli, Mscl, Stul BsmFI (GTT/RAC) Dral, Pmel Msel (GTC/RAC) Hpal HincII (GTT/RAC) Hpal HincII, Hpal, Msel (GTT/RAC) Nrul Taql (GTC/RAC) Rsal, Scal Tsp45I Mscl BsrBI, MspA1I (CMG/CGG) Acil KGG/CCA) BstZ171 Rsal Haelli, Stul Haelli Haelli EcoRV Alwi, Mbol Sfol Haelli, Sau96I		HincII (GTY/AAC) HincII (GTY/GAC)	HincII, Hpal, Msel HincII
(TGG/CCA) BstZ17I Rsal HaellI, Stul HaelII EcoRV AlwI, Mbol Sfol HaelII, Sau96I	(GTC/RAC) (GTC/RAC) (GTC/RAC) (GTT/RAC) (GTC/RAC) (GTT/RAC) (GTT/RAC)	BstZ17I Haelli, Msci, Stul Dral, Pmel Hpal Hpal Nrul	Accl BsmFl Msel Hincll Hincll, Hpal, Msel Taql
		BstZ17I HaellI, Stul EcoRV Sfol	Rsal Haelll Alwl, Mbol Haelll, Sau96l

Enzymes that have degenerate recognition sequences (e.g., recognize more than one sequence) are followed by a specific sequence in parentheses and are only listed if a non-degenerate equivalent does not exist. Be aware that these degenerate enzymes will cleave sequences in addition to the one listed. Where isoschizomers exist, only one member of each set is listed. Only commercially available enzymes are shown.

ENZYME	LIGATED TO	RECLEAVED BY
MspA11 (CAG/CKG) (CCG/CKG) (CCG/CKG) (CAG/CKG) (CMG/CKG) (CMG/CKG) (CMG/CKG) (CMG/CKG) (CMG/CKG) (CMG/CKG) (CMG/CKG) (CMG/CKG) (CMG/CKG) (CCG/CKG) (CCG/CKG) (CCG/CKG) (CCG/CKG) (CCG/CKG)	Alul Alul, HaellI, MscI, StuI BsaAI, FspI, HincII (GTY/GAC), PmII, SnaBI BsrBI, PvuII BsrBI BstZ17I BstUI, NruI Afel EcoRV SfoI SfoI Nael PvuII SmaI	Alul Acil Hpall Acil, MspA11 Acil, BsaJI, BstUI, MspA1I, Sac II Rsal Acil, BstUI Hpall Mbol HaellI, Hpall Hpall, Ncil, ScrFI Alul, MspA11, Pvull BsaJI, Hpall, Ncil, ScrFI
Nael (GCC/GGC)	BsrBI, MspA1I (CMG/CGG), Smal Afel, SfoI, FspI	Hpall, Ncil, ScrFl Acil, Fnu4Hl
Nrul (TCG/CGA)	BsrBI, MspA1I (CMG/CGG) BstZ17I BstUI Dral, HincII (GTY/AAC), HpaI, PmeI, RsaI, ScaI, SspI EcoRV SfoI	Acil, BstUI Rsal BstUI Taql Mbol, Taql Haelli
Pmel (GTTT/AAAC)	Dral, Swal HincII (GTY/AAC), Hpal Nrul	Dral, Msel Msel Taql
PmII (CAC/GTG)	BsaAI (YAC/GTA), SnaBI BsaAI (YAC/GTG) BsrBI, MspA1I (CMG/CGG) BmgBI	BsaAl BsaAl, Pmll Hpall Pmll
Pvull (CAG/CTG)	Alul BsrBI, MspA1I (CMG/CGG) BstZ17I EcoRV SfoI MspA1I (CMG/CTG)	Alul Acil, MspA1I Rsal Mbol HaeIII Alul, MspA1I, PvuII
Rsal (GT/AC)	HincII (GTY/GAC) Nrul Scal	Tsp45I TaqI RsaI
Scal (AGT/ACT)	HincII (GTY/GAC) NruI RsaI	Tsp45I Taql Rsal
Sfol (GGC/GCC)	Alul, BstUI, MspA1I (CMG/CKG), NruI, PvuII BsrBI, MspA1I (CMG/CGG) HaeIII, MscI, StuI AfeI EcoRV FspI NaeI SmaI	Haelli Haelli, Hpali Haelli, Sau96i Haeli, Hhai SfaNi Hhai Acii, Fnu4Hi Acii
Smal (CCC/GGG)	BsrBI, MspA1I (CMG/CGG) Afel, Sfol, FspI Nael	BsaJl, Hpall, Ncil, ScrFl Acil Hpall, Ncil, ScrFl
SnaBI (TAC/GTA)	BsaAI (YAC/GTA) BsaAI (YAC/GTG), PmII BsrBI, MspA1I (CMG/CGG)	BsaAl, SnaBl BsaAl Hpall
Sspl (AAT/ATT)	BstZ17I Nrul	MluCl Taql
Stul (AGG/CCT)	BsrBI, MspA1I (CMG/CGG) BstZ17I HaeIII, MscI EcoRV StoI HincII (GTY/GAC)	Acil Rsal HaellI Alwl, Mbol HaellI, Sau96I BsmFI
Swal (ATTT/AAAT)	Dral, Pmel	Dral, Msel

Compatible Cohesive Ends

Restriction enzymes that produce compatible cohesive ends often produce recleavable ligation products. The combinations listed were identified by computer analysis, and have not necessarily been confirmed by experimentation.

Where isoschizomers exist, only one member of each set is listed. A selection of enzymes available from New England Biolabs has been listed. For a more complete listing visit our website, **www.neb.com**

ENZYME LIGATED TO RECLEAVED BY Acc65I Banl (G/GTACC) Acc65I, Banl, Kpnl, NlaIV, Rsal (G/GTACC) BsiWI, BsrGI Rsal Accl (GT/CGAC) Acil, Acil, BsaHl (GR/CGYC), HinP1I, Hpall, Narl (GT/CGAC) Clal, BstBl, Tagl Taql Accl (GT/CGAC), AcII, Clal, BstBI, Taql Acil BsaHl (GR/CGCC), HinP1I, Narl (C/CGC) Acil Hpall AcII Accl (GT/CGAC), Acil, Clal, BstBl, (AA/CGTT) HinP1I, Hpall, Narl, Tagl Aval (C/CCGGG), Xmal Hpall, Ncil, ScrFl Agel (A/CCGGT) BsaWI, BspEI BsaWl, Hpall BsrFI (A/CCGGT), SgrAI (CA/CCGGTG) Agel, BsaWl, BsrFl, Hpall BsrFI, Hpall Apal, Banll, Bsp120l, Bsp1286l, Apal Banll (GGGCC/C), Bsp1286l (GGGCC/C) (GGGCC/C) Haelli, NialV, Sau961 ApaLl SfcI (C/TGCAG) (G/TGCAC) Apol (A/AATTY) EcoRI Apol, MluCl (G/AATTY) EcoRI Apol, EcoRI, MluCI (R/AATTY) Mfel, MluCl BstUI, Hhal Ascl Mlul (GG/CGCGCC) BssHII, BstUI, Cac8I, Hhal BssHII Bfal, Csp6l, Ndel Asel (AT/TAAT) Msel Msel AsiSI BsiEI (CGAT/CG) DpnII, PvuI (GCGAT/CGC) Pacl Msel Pvul DpnII, Pvul Aval Agel, BsaWl, BspEl, BsrFl (R/CCGGY), (C/CCGGG) NgoMIV, SgrAI (CR/CCGGYG) Hpall, Ncil, ScrFl (C/TCGAG) Aval, Taql, Xhol (C/TCGAG) Sall Tagl Aval, BsaJl, Hpall, Ncil, (C/CCGGG) Xmal ScrFI, Smal Avall, NlalV, Sau961 Avall PpuMI (RG/GACCY) (G/GWCC) Avall, Sau961 PpuMI (RG/GTCCY) Avall, BsmFl, NlaIV, Sau961 AvrII Nhel, Spel, Xbal (C/CTAGG) AvrII, Bfal, BsaJI, Styl Styl (C/CTAGG) Bcll, Dpnll Alwl, DpnII BamHI BgIII, BstYI (R/GATCY) Alwl, BstYI, DpnII (G/GATCC) BstYI (G/GATCC) Alwl, BamHI, BstYI, DpnII, NIalV Banl Acc65I, Banl, Kpnl, NlaIV, Rsal (G/GTACC) Acc65I (G/GCGCC) Banl, BsaHl, Haell, Hhal, Kasl Kasl, Narl, NlalV (G/GTACC) BsiWI, BsrGI Rsal (GGGCC/C) Apal, Bsp1286I (GGGCC/C) Apal, Banll, Bsp1286l, Haelll, NIaIV Sau96I Alul, Banll, BsiHKAI, (GAGCT/C) Bsp1286I (GAGCT/C), SacI Bsp1286I, SacI BamHI, BstYI (R/GATCY) Alwl, DpnII (T/GATCA) BgIII, Mbol DpnII

Enzymes that have degenerate recognition sequences (e.g., recognize more than one sequence) are followed by a specific sequence in parentheses and are only listed if a non-degenerate equivalent does not exist. Be aware that these degenerate enzymes will cleave sequences in addition to the one listed. A "-" denotes a ligation product that cannot be recleaved.

ENZVIAE	LICATED TO	DECLEAVED BY
ENZYME	LIGATED TO	RECLEAVED BY
Bfal (C/TAG)	Asel, Csp6l, Msel, Ndel	_
BgIII (A/GATCT)	BamHI, BstYI (R/GATCY) BcII, DpnII	Alwl, BstYI, DpnII DpnII
BsaHI (GR/CGYC) (GA/CGYC) (GG/CGYC) (GG/CGYC) (GA/CGYC) (GG/CGYC)	Accl (GT/CGAC), Clal, BstBl, Taql Acil, HinP11 Acil, HinP11 Hpall Narl Narl	Hgal Hhal Acil BsaHI, Hgal Banl, BsaHI, Haell, Hhal, Narl, NlaIV
BsaWI (W/CCGGW)	Agel, BsrFI (R/CCGGY), SgrAI (CR/ CCGGYG) AvaI (C/CCGGG), XmaI BspEI BsrFI (R/CCGGY), NgoMIV NgoMIV	Agel, BsaWl, BsrFl, Hpall Hpall, Ncil, ScrFl BsaWl, BspEl, Hpall BsrFl, Hpall Hpall
BsiEI (CGAT/CG) (CGAT/CG) (CGGC/CG)	Pacl Pvul Sacll	Msel BsiEl, Dpnll, Pvul Acil
BsiHKAI (GTGCA/C)	Bsp1286I (GTGCA/C) Bsp1286I (GAGCA/C) Bsp1286I (GAGCT/C), Sacl	BsiHKAI, Bsp1286I BsiHKAI, Bsp1286I Alul, BanII, BsiHKAI, Bsp1286I, SacI
	Pstl, Sbfl	Bsgl
BsiWI (C/GTACG)	Acc65I, Banl (G/GTACC), BsrGI	Rsal
Bsp1286l (GGGCC/C)	Apal, Banll (GGGCC/C)	Apal, Banll, Bsp1286l, Haelll, NlaIV, Sau96l
(GTGCA/C) (GGGCC/C) (GAGCT/C)	BsiHKAI BanII (GGGCC/C) BanII (GAGCT/C), BsiHKAI, SacI	ApaLi, BsiHKAI, Bsp1286I Banli, Bsp1286I Alui, Banli, BsiHKAI, Bsp1286I, Saci
(GWGCW/C) (GTGCA/C) (GTGCA/C)	BsiHKAI Nsil Pstl, Sbfl	BsiHKAI, Bsp1286I — BsgI
BspEI (T/CCGGA)	Agel, BsaWI, BsrFI (R/CCGGY), SgrAI (CR/CCGGYG) AvaI (C/CCGGG), XmaI BsaWI BsrFI (R/CCGGY), NgoMIV	BsaWl, Hpall Hpall, Ncil, ScrFl BsaWl, BspEl, Hpall Hpall
BspHI (T/CATGA)	Fatl, Ncol, Pcil	Fatl, NIallI
BsrFI (A/CCGGY) (G/CCGGY) (R/CCGGY) (A/CCGGY) (R/CCGGY) (G/CCGGY) (CR/CCGGYG)	Agel, BsaWI Agel, BsaWI, NgoMIV Aval (C/CCGGG), Xmal BsaWI, BspEI BsaWI, BspEI NgoMIV SgrAI	Agel, BsaWl, BsrFl, Hpall BsrFl, Hpall Hpall, Ncil, ScrFl BsaWl, Hpall Hpall BsrFl, Cac8l, Hpall, Nael BsrFl, Hpall
BsrGI (T/GTACA)	Acc651, Banl (G/GTACC), BsiWl	Rsal
BssHII (G/CGCGC)	Mlul Ascl	BstUI, Hhal BssHII, BstUI, Cac8I, Hhal
BstBI (TT/CGAA)	Accl (GT/CGAC), Clal, Taql Acil, Acll, BsaHI (GR/CGYC), HinP1I, HpaII, Narl	Taql

ENZYME	LIGATED TO	RECLEAVED BY
BstYI (A/GATCY) (G/GATCY) (R/GATCY) (G/GATCY) (A/GATCY)	BamHI, BgIII BamHI BcII, DpnII BcII, DpnII BgIII	Alwi, BstYl, Dpnll Alwi, BamHl, BstYl, Dpnll, NlaIV Dpnll Alwi, Dpnll Bglll, BstYl, Dpnll
Clal (AT/CGAT)	Accl (GT/CGAC), BstBI, Taql Acil, Acil, BsaHI (GR/CGYC), HinP1I, HpaII, Narl	Taql
DpnII/Mbol/ Sau3AI (/GATC)	BamHI, BstYI (R/GATCC) BcII, BgIII, BstYI (R/GATCY)	Alwl, DpnII DpnII
Eael (Y/GGCCR) (C/GGCCR) (T/GGCCR) (C/GGCCR) (T/GGCCR)	PspOMI Eagl Eagl Notl Notl	Haelli, Sau96l BsiEl, Eael, Eagl, Haelli Eael, Haelli Acil, BsiEl, Eael, Eagl, Fnu4Hl, Haelli Acil, Eael, Fnu4Hl, Haelli
Eagl (C/GGCCG)	PspOMI Eael (Y/GGCCR) Eael (C/GGCCG) Notl	Haelli, Sau96l Eael, Haelll BsiEl, Eael, Eagl, Haelll Acil, BsiEl, Eael, Eagl, Fnu4Hl, Haelll
EcoRI (G/AATTC)	Apol (G/AATTC) Apol (R/AATTY) Mfel, MluCl	Apol, EcoRI, MIuCI Apol, MIuCI MIuCI
FatI (/CATG)	BspHI, Ncol, Pcil	Fatl, NIalli
HinP1I (G/CGC)	Accl (GT/CGAC), Acll, Clal, BstBl, Taql Acil, BsaHl (GR/CGCC), Narl BsaHl (GR/CGTC) Hpall	— Hhal Hgal Acil
Hpall/Mspl (C/CGG)	Accl (GT/CGAC), Acll, Clal, BstBl, Taql Acil, BsaHl (GR/CGCC), HinP11, Narl	— Acil
Kasl (G/GCGCC)	Banl (G/GCGCC)	Banl, BsaHl, Haell, Hhal, Kasl, Narl, NlalV
Mfel (C/AATTG)	Apol (R/ATTTY), EcoRI, MluCI	MluCl
Mlul (A/CGCGT)	Ascl, BssHII	BstUI, Hhal
MluCl (/AATT)	Apol (R/AATTY), EcoRI, Mfel	MluCl
Msel (T/TAA)	Asel Bfal, Csp6l, Ndel	Msel —
Narl (GG/CGCC)	Accl (GT/CGAC), Acll, Clal, BstBl, Taql Acil, HinP11 BsaHl (GR/CGCC) BsaHl (GR/CGTC) Hpall	— Hhal Banl, BsaHl, Haell, Hhal, Narl, NlaIV BsaHl, Hgal Acil
Ncol (C/CATGG)	BspHl, Fatl, Pcil	Fatl, NIaIII
Ndel (CA/TATG)	Asel, Bfal, Csp6l, Msel	_
NgoMIV (G/CCGGC)	Agel, BsaWl, BsrFl (R/CCGGY), SgrAl Aval (C/CCGGG), Xmal BsaWl, BspEl BsrFl (R/CCGGC), SgrAl	BsrFl, Hpall Hpall, Ncil, ScrFl Hpall BsrFl, Cac8l, Hpall, Nael
Nhel (G/CTAGC)	AvrII, Spel, Styl (C/CTAGG), Xbal	Bfal
NIaIII (CATG/)	Sphl, Nspl	NIalli
Notl (GC/GGCCGC)	PspOMI Eagl	Acil, Eael, Fnu4HI, HaelII Acil, BsiEl, Eael, Eagl, Fnu4HI, HaelII
	Eael (Y/GGCCR)	Acil, BsiEl, Eael, Fnu4Hl, Haelll

ENZYME	LIGATED TO	RECLEAVED BY
Nsil (ATGCA/T)	BsiHKAI (GTGCA/C), Bsp1286I (GTGCA/C), Pstl, Sbfl	_
NspI (RCATG/Y)	Nialli, Sphi	NIaIII
Pacl (TTAAT/TAA)	AsiSI BsiEI (CGAT/CG), Pvul	Msel
Pcil (A/CATGT)	BspHI, Fatl, Ncol	Fatl, NIaIII
PluTI (GGCGC/C)	Haell	Haell
PpuMI (RG/GWCCY) (GG/GTCCY) (GG/GACCY)	Avall, Rsrll Avall, Rsrll Avall, Rsrll	Avall, Sau96l Avall, BsmFl, NialV, Sau96l Avall, NialV, Sau96l
PspOMI (G/GGCCC)	Eael (Y/GGCCR), Eagl Notl	Haelli, Sau96l Acil, Fnu4Hl, Haelli, Sau96l
PspXI (VC/TCGAGB)	Xhol, Tlil Sall	Xhol, Tlil Taql
PstI (CTGCA/G)	BsiHKAI, Bsp1286I (GTGCA/C) NsiI	Bsgl —
Pvul (CGAT/CG)	Sbfl AsiSI Pacl BsiEl (CGAT/CG)	Pstl DpnII, Pvul Msel BsiEI, DpnII, Pvul
RsrII (CG/GWCCG)	Avall, PpuMI (RG/GACCY) PpuMI (RG/GACCY) PpuMI (RG/GTCCY)	Avall, Sau96l Avall, NlaIV, Sau96l Avall, BsmFI, NlaIV, Sau96l
Sacl (GAGCT/C)	Banll (GAGCT/C), BsiHKAI, Bsp1286l (GAGCT/C)Alul, Banll, BsiHKAI, Bsp1286I, Sacl
SacII (CCGC/GG)	BsiEI (CGGC/CG)	Acil
Sall (G/TCGAC)	PspXI, XhoI	Taql
Sbfl (CCTGCA/GG)	BsiHKAI, Bsp1286I (GTGCA/C) NsiI PstI	Bsgl — Pstl
SfcI (C/TGCAG)	ApaLl	Bsgl
SgrAI (CR/CCGGYG)	See BsrFI	
Spel (A/CTAGT)	AvrII, Nhel, Styl (C/CTAGG), Xbal	Bfal
SphI (GCATG/C)	NIaIII, Nspl	Nialli
Styl (C/CTAGG) (C/CATGG)	AvrII Nhel, Spel, Xbal BspHI Ncol	Avrll, Bfal, BsaJl, Styl Bfal Nlalll BsaJl, Ncol, Nlalll, Styl
Taql (T/CGA)	Accl (GT/CGAC), Clal, BstBl Acil, Acll, BsaHl (GR/CGYC), HinP1I, Hpall, Narl	Taql
Xbal (T/CTAGA)	AvrII, Nhel, Spel, Styl (C/CTAGG)	Bfal
Xhol (Tlil) (C/TCGAG)	PspXI Sall	Xhol, Tlil Taql
Xmal (C/CCGGG)	Agel, BsaWl, BspEl, BsrFl, NgoMlV, SgrAl Aval (C/CCGGG)	Hpall, Ncil, ScrFl Aval, BsaJl, Hpall, Ncil, ScrFl, Smal, Xmal

Dam (G^mATC), Dcm (C^mCWGG) and CpG (^mCG) Methylation

DNA methyltransferases (MTases) that transfer a methyl group from S-adenosylmethionine to either adenine or cytosine residues are found in a wide variety of prokaryotes and eukaryotes. Methylation should be considered when digesting DNA with restriction endonucleases because cleavage can be blocked or impaired when a particular base in the recognition site is methylated.

Prokaryotic Methylation

In prokaryotes, MTases have most often been identified as elements of restriction/modification systems that act to protect host DNA from cleavage by the corresponding restriction endonuclease. Most laboratory strains of *E. coli* contain three site-specific DNA methyltransferases.

- Dam methyltransferases—methylation at the N⁶ position of the adenine in the sequence GATC (1,2).
- Dcm methyltransferases— methylation at the C⁵ position of cytosine in the sequences CCAGG and CCTGG (1,3).
- EcoKl methylase—methylation of adenine in the sequences AAC(N⁶A)GTGC and GCAC(N⁶A)GTT.

Some or all of the sites for a restriction endonuclease may be resistant to cleavage when isolated from strains expressing the Dam or Dcm MTase if the methylase recognition site overlaps the endonuclease recognition site. For example, plasmid DNA isolated from dam* E. coli is completely resistant to cleavage by Mbol, which cleaves at GATC sites.

Not all DNA isolated from $\it E. coli$ is methylated to the same extent. While pBR322 DNA is fully modified (and is therefore completely resistant to Mbol digestion), only about 50% of $\it \lambda$ DNA Dam sites are methylated, presumably because the methylase does not have the opportunity to methylate the DNA fully before it is packaged into the phage head. As a result, enzymes blocked by Dam or Dcm modification will yield partial digestion patterns with $\it \lambda$ DNA.

Restriction sites that are blocked by Dam or Dcm methylation can be un-methylated by cloning your DNA into a *dam*-, *dcm*- strain of *E. coli*, such as *dam*-/*dcm*- Competent *E. coli* (NEB #C2925).

Restriction sites can also be blocked if an overlapping site is present. In this case, part of the Dam or Dcm sequence is generated by the restriction enzyme sequence, followed by the flanking sequence. This situation should also be considered when designing restriction enzyme digests.

Eukaryotic Methylation

CpG MTases, found in higher eukaryotes (e.g., Dnmt1), transfer a methyl group to the C⁵ position of cytosine residues. Patterns of CpG methylation are heritable, tissue specific and correlate with gene expression. Consequently, CpG methylation has been postulated to play a role in differentiation and gene expression (4).

Note: The effects of CpG methylation are mainly a concern when digesting eukaryotic genomic DNA. CpG methylation patterns are not retained once the DNA is cloned into a bacterial host.

Methylation Sensitivity

The table below summarizes methylation sensitivity for NEB restriction enzymes, indicating whether or not cleavage is blocked or impaired by Dam, Dcm or CpG methylation if or when it overlaps each recognition site. This table should be viewed as a guide to the behavior of the enzymes listed rather than an absolute indicator. Consult REBASE (http://rebase.neb.com/rebase/), the restriction enzyme database, for more detailed information and specific examples upon which these guidelines are based.

References

- (1) Marinus, M.G. and Morris, N.R. (1973) *J. Bacteriol.*, 114, 1143–1150.
- (2) Geier, G.E. and Modrich, P. (1979) *J. Biol. Chem.*, 254, 1408–1413.
- (3) May, M.S. and Hattman, S. (1975) *J. Bacteriol.*, 123, 768–770.
- (4) Siegfried, Z. and Cedar, H. (1997) *Curr. Biol.*, 7, r305–307.

Legend

- Not Sensitive
- Blocked
- □ ol Blocked by Overlapping
- □ scol Blocked by Some Combinations of Overlapping
- Impaired
- ol Impaired by Overlapping
- scol Impaired by Some Combinations of Overlapping

Single Letter Code

ENZYME	SEQUENCE	Dam	Dcm	CpG
Aatll	GACGT/C	•	•	
AbaSI	^m C(11/9)	•	•	•
Accl	GT/MKAC	•	•	□ ol
Acc65I	G/GTACC	•	□ scol	□ scol
Acil	CCGC(-3/-1)	•	•	
AcII	AA/CGTT	•	•	
Acul	CTGAAG(16/14)	•	•	•
Afel	AGC/GCT	•	•	
AfIII	C/TTAAG	•	•	•
AfIIII	A/CRYGT	•	•	•
Agel	A/CCGGT	•	•	
Agel-HF	A/CCGGT	•	•	
Ahdl	GACNNN/NNGTC	•	•	
Alel-v2	CACNN/NNGTG	•	•	
Alul	AG/CT	•	•	•
Alwl	GGATC(4/5)		•	•
AlwNI	CAGNNN/CTG	•	□ ol	•
Apal	GGGCC/C	•	□ ol	□ ol
ApaLI	G/TGCAC	•	•	□ ol
ApeKI	G/CWGC	•	•	□ ol
Apol	R/AATTY	•	•	•
Apol-HF	R/AATTY	•	•	•
Ascl	GG/CGCGCC	•	•	•
Asel	AT/TAAT	•	•	•
AsiSI	GCGAT/CGC	•	•	•
Aval	C/YCGRG	•	•	
Avall	G/GWCC	•	□ ol	□ ol
AvrII	C/CTAGG	•	•	•
Bael	(10/15)ACNNNNGTAYC(12/7)	•	•	□ scol

ENZYME	SEQUENCE	Dam	Dcm	CpG
BaeGI	GKGCM/C	•	•	•
BamHI	G/GATCC	•	•	•
BamHI-HF	G/GATCC	•	•	•
Banl	G/GYRCC	•	□ scol	□ scol
Banll	GRGCY/C	•	•	•
Bbsl	GAAGAC(2/6)	•	•	•
BbsI-HF	GAAGAC(2/6)	•	•	•
Bbvl	GCAGC(8/12)	•	•	•
BbvCl	CCTCAGC(-5/-2)	•	•	♦ ol
Bccl	CCATC(4/5)	•	•	•
BceAl	ACGGC(12/14)	•	•	
Bcgl	(10/12)CGANNNNNNTGC(12/10)	♦ 0 l	•	□ scol
BciVI	GTATCC(6/5)	•	•	•
BcII	T/GATCA	•	•	•
BcII-HF	T/GATCA	•	•	•
BcoDI	GTCTC(1/5)	•	•	
Bfal	C/TAG	•	•	•
BfuAl	ACCTGC(4/8)	•	•	♦ ol
BgII	GCCNNNN/NGGC	•	•	□ scol
BgIII	A/GATCT	•	•	•
Blpl	GC/TNAGC	•	•	•
BmgBl	CACGTC(-3/-3)	•	•	
Bmrl	ACTGGG(5/4)	•	•	•
Bmtl	GCTAG/C	•	•	•
Bmtl-HF	GCTAG/C	•	•	•
Bpml	CTGGAG(16/14)	•	•	•
Bpu10I	CCTNAGC(-5/-2)	•	•	•
BpuEl	CTTGAG(16/14)	•	•	•
Bsal	GGTCTC(1/5)	•		□ scol

ENZYME	SEQUENCE	Dam	Dcm	CpG
Bsal-HFv2	GGTCTC(1/5)	•		□ scol
BsaAl	YAC/GTR	•	•	•
BsaBl	GATNN/NNATC	□ ol	•	□ scol
BsaHI	GR/CGYC	•	□ scol	
BsaJI	C/CNNGG	•	•	•
BsaWl	W/CCGGW	•	•	•
BsaXI	(9/12)ACNNNNNCTCC(10/7)	•	•	•
BseRI	GAGGAG(10/8)	•	•	•
BseYI	CCCAGC(-5/-1)	•	•	□ ol
Bsal	GTGCAG(16/14)	•	•	•
BsiEl	CGRY/CG	•	•	
BsiHKAI	GWGCW/C	•	•	•
BsiWl	C/GTACG	•	•	
BsiWI-HF	C/GTACG	•	•	
Bsll	CCNNNNN/NNGG	•	□ scol	_ scol
Bsml	GAATGC(1/-1)	•	0 3001	0 3001
BsmAl	GTCTC (1/5)			□ scol
	, ,	•	•	□ SCOI
BsmBl	CGTCTC(1/5)			_
BsmFl	GGGAC(10/14)	•	□ ol	□ 0l
BsoBl	C/YCGRG	•	•	•
Bsp1286I	GDGCH/C	•	•	•
BspCNI	CTCAG(9/7)	•	•	•
BspDI	AT/CGAT	□ ol	•	•
BspEl	T/CCGGA	□ ol	•	•
BspHI	T/CATGA	◇ 0l	•	•
BspMI	ACCTGC(4/8)	•	•	•
BspQI	GCTCTTC(1/4)			
Bsrl	ACTGG(1/-1)	•	•	•
BsrBI	CCGCTC(-3/-3)			□ scol
BsrDI	GCAATG(2/0)	•	•	•
BsrFI-v2	R/CCGGY	•	•	
BsrGI	T/GTACA	•	•	•
BsrGI-HF	T/GTACA	•	•	•
BssHII	G/CGCGC	•	•	
BssSI-v2	CACGAG(-5/-1)	•	•	•
BstAPI	GCANNN/NTGC	•	•	□ scol
BstBI	TT/CGAA	•	•	
BstEII	G/GTNACC	•	•	•
BstEII-HF	G/GTNACC	•	•	•
BstNI	CC/WGG	•	•	•
BstUI	CG/CG	•		
BstXI	CCANNNN/NTGG	•	□ scol	•
BstYI	R/GATCY	•	0 3001	
BstZ17I-HF	GTA/TAC	•	•	□ scol
Bsu36I	CC/TNAGG	•	•	■ SCOI
	C/CRYGG	•	•	•
Btgl				
BtgZl	GCGATG(10/14)	•	•	•
Btsl-v2	GCAGTG(2/0)	•	•	•
BtsCl	GGATG(2/0)	•	•	•
BtslMutl	CAGTG(2/0)	•	•	•
Cac8I	GCN/NGC	•	•	□ scol
Clal	AT/CGAT	□ ol	•	
CspCl	(11/13)CAANNNNNGTGG(12/10)	•	•	•
CviAII	C/ATG	•	•	•
CviKI-1	RG/CY	•	•	•
CviQI	G/TAC	•	•	•
Ddel	C/TNAG	•	•	•
DpnI	GA/TC	•	•	□ ol
DpnII	/GATC	•	•	•
Dral	TTT/AAA	•	•	•
DraIII-HF	CACNNN/GTG	•	•	♦ 0 l
Drdl	GACNNNN/NNGTC	•	•	□ scol
	Y/GGCCR	•	□ ol	□ ol

ENZYME	SEQUENCE	Dam	Dcm	CpG
Eagl	C/GGCCG	•	•	•
Eagl-HF	C/GGCCG	•	•	•
Earl	CTCTTC(1/4)	•		♦ 01
Ecil	GGCGGA(11/9)	•	•	□ scol
Eco53kl	GAG/CTC	•	•	□ scol
EcoNI	CCTNN/NNNAGG	•	•	•
Eco0109I	RG/GNCCY	•	□ ol	•
EcoP15I	CAGCAG(25/27)	•	•	•
EcoRI	G/AATTC	•	•	□ scol
EcoRI-HF	G/AATTC	•	•	□ scol
EcoRV	GAT/ATC	•	•	
EcoRV-HF	GAT/ATC	•	•	
Esp3I	CGCCTC(1/5)	•	•	• 555.
Fatl	/CATG	•	•	•
Faul	CCCGC(4/6)	•	•	•
Fnu4HI	GC/NGC	•	•	- □ 0l
Fokl		•	♦ ol	□ 01
-	GGATG(9/13)			
Fsel	GGCCGG/CC	•	♦ scol	•
Fspl	TGC/GCA	•	•	•
FspEl	C5mCNNNNNNNNNNNN	•	•	•
Haell	RGCGC/Y	•	•	
Haelll	GG/CC	•	•	•
Hgal	GACGC(5/10)	•		
Hhal	GCG/C	•	•	•
HincII	GTY/RAC	•	•	□ scol
HindIII	A/AGCTT	•	•	•
HindIII-HF	A/AGCTT	•	•	•
Hinfl	G/ANTC	•	•	□ scol
HinP1I	G/CGC	•	•	
Hpal	GTT/AAC	•	•	□ scol
Hpall	C/CGG	•	•	=
Hphl	GGTGA(8/7)	•		•
Hpy99I	CGWCG/	-	-	
Hpy166II	GTN/NAC	•	•	- □ 0l
Hpy188I	TCN/GA	□ ol	•	•
Hpy188III	TC/NNGA		•	□ ol
1 7	** *	□ ol		
HpyAV	CCTTC(6/5)	•	•	♦ 0 l
HpyCH4III	ACN/GT	•	•	•
HpyCH4IV	A/CGT	•	•	
HpyCH4V	TG/CA	•	•	•
Kasl	G/GCGCC	•	•	
KpnI	GGTAC/C	•	•	•
KpnI-HF	GGTAC/C	•		
LpnPI	C5mCDGNNNNNNNNNN	•	•	•
Mbol	/GATC	•		♦ 01
Mboll	GAAGA(8/7)	□ ol	•	•
Mfel	C/AATTG	•	•	•
Mfel-HF	C/AATTG	•	•	•
Mlul	A/CGCGT	•	•	
Mlul-HF	A/CGCGT	•	•	
MluCl	/AATT	•	•	•
Mlyl	GAGTC(5/5)	•	•	•
Mmel	TCCRAC(20/18)	•	•	□ ol
MnII	CCTC(7/6)	•	•	•
Mscl	TGG/CCA	•	□ ol	•
		•	□ 01	•
Msel	T/TAA			
MsII	CAYNN/NNRTG	•	•	•
Mspl	C/CGG	•	•	•
MspA1I	CMG/CKG	•	•	□ ol
MspJI	5mCNNRNNNNNNNNN	•	•	•
Mwol	GCNNNNN/NNGC	•	•	□ scol
Nael	GCC/GGC	•	•	•
Narl	GG/CGCC	•		

Dam (G^mATC), Dcm (C^mCWGG) and CpG (^mCG) Methylation (continued)

ENZYME	SEQUENCE	Dam	Dcm	CpG
Nb.BbvCl	CCTCAGC(none/-2)	•	•	•
Nb.Bsml	GAATGC(none/-1)	•	•	•
Nb.BsrDI	GCAATG(none/0)	•	•	•
Nb.BssSI	CACGAG(none/-1)	•	•	•
Nb.BtsI	GCAGTG(none/0)	•	•	•
Ncil	CC/SGG	•	•	♦ 0 l
Ncol	C/CATGG	•	•	•
Ncol-HF	C/CATGG	•	•	•
Ndel	CA/TATG	•	•	•
NgoMIV	G/CCGGC	•	•	•
Nhel	G/CTAGC	•	•	□ scol
Nhel-HF	G/CTAGC	•	•	□ scol
NIaIII	CATG/	•	•	•
NIaIV	GGN/NCC	•	□ ol	□ ol
NmeAIII	GCCGAG(21/19)	•	•	•
Notl	GC/GGCCGC	•	•	•
NotI-HF	GC/GGCCGC	•	•	
Nrul	TCG/CGA	□ ol	•	•
Nrul-HF	TCG/CGA		•	
Nsil	ATGCA/T	•	•	•
Nsil-HF	ATGCA/T	•	•	•
Nspl	RCATG/Y	•	•	•
Nt.Alwl	GGATC(4/none)		•	•
Nt.BbvCI	CCTCAGC(-5/none)	•	•	□ scol
Nt.BsmAl	GTCTCN(1/none)	•	•	
Nt.BspQI	GCTCTTC(1/none)	•	•	•
Nt.BstNBI	GAGTC(4/none)	•	•	•
Nt.CviPII	CCD(-3/none)	•	•	
Pacl	TTAAT/TAA	•	•	•
PaeR7I	C/TCGAG	•	•	
Pcil	A/CATGT	•	•	•
PfIFI	GACN/NNGTC	•	•	•
PfIMI	CCANNNN/NTGG	•	□ ol	•
Plel	GAGTC(4/5)	•	•	□ scol
PluTl	GGCGC/C	•	•	
Pmel	GTTT/AAAC	•	•	□ scol
PmII	CAC/GTG	•	•	
PpuMI	RG/GWCCY	•	□ ol	•
PshAl	GACNN/NNGTC	•	•	□ scol
Psil	TTA/TAA	•	•	•
PspGI	/CCWGG	•	•	•
Psp0MI	G/GGCCC	•	♦ scol	□ ol
PspXI	VC/TCGAGB	•	•	•
PstI	CTGCA/G	•	•	•
PstI-HF	CTGCA/G	•	•	
Pvul	CGAT/CG	•	•	
Pvul-HF	CGAT/CG	•	•	
Pvull	CAG/CTG	•	•	•
Pvull-HF	CAG/CTG	•	•	•

ENZYME	SEQUENCE	Dam	Dem	CpG
Rsal	GT/AC	•	•	□ scol
RsrII	CG/GWCCG	•	•	
Sacl	GAGCT/C	•	•	•
SacI-HF	GAGCT/C	•	•	□ scol
SacII	CCGC/GG	•	•	•
Sall	G/TCGAC	•	•	
Sall-HF	G/TCGAC	•	•	
Sapl	GCTCTTC(1/4)	•	•	•
Sau3AI	/GATC	•	•	□ ol
Sau96l	G/GNCC	•	□ ol	□ ol
Sbfl	CCTGCA/GG	•	•	•
Sbfl-HF	CCTGCA/GG	•	•	•
Scal-HF	AGT/ACT			
ScrFI	CC/NGG	•	□ ol	□ ol
SexAl	A/CCWGGT			
SfaNI	GCATC(5/9)	•	•	
SfcI	C/TRYAG			
Sfil	GGCCNNNN/NGGCC	•	♦ 01	□ scol
Sfol	GGC/GCC		□ scol	
SgrAl	CR/CCGGYG	•	•	
Smal	CCC/GGG	•	•	
SmII	C/TYRAG	•	•	•
SnaBl	TAC/GTA	•	•	•
Spel	A/CTAGT	•	•	•
Spel-HF	A/CTAGT	•	•	•
Sphl	GCATG/C	•	•	•
SphI-HF	GCATG/C	•	•	•
Srfl	GCCC/GGGC	•	•	•
Sspl	AAT/ATT	•	•	•
SspI-HF	AAT/ATT	•	•	•
Stul	AGG/CCT	•	□ ol	•
Styl	C/CWWGG	•	•	•
Styl-HF	C/CWWGG	•	•	•
StyD4I	/CCNGG	•	□ ol	♦ 0l
Swal	ATTT/AAAT	•	•	•
Taq¤l	T/CGA	□ 0l	•	•
Tfil	G/AWTC	•	•	□ scol
Tsel	G/CWGC	•	•	□ scol
Tsp45I	/GTSAC	•	•	•
TspMI	C/CCGGG	•	•	•
TspRI	NNCASTGNN/	•	•	•
Tth1111	GACN/NNGTC		•	•
Xbal	T/CTAGA CCANNNN/NNNNTGG	□ ol	•	•
Xcml Xhol				•
Xmal	C/TCGAG C/CCGGG	•	•	•
Xmai	GAANN/NNTTC	•	•	•
		•	_	•
Zral	GAC/GTC		•	•

PPENDIX

General Guidelines for PCR Optimization

New England Biolabs offers a diverse group of DNA Polymerases for PCR-based applications. Specific recommendations for PCR optimization can be found in the product literature or on the individual product webpages. However, these general guidelines will help to ensure success using New England Biolabs' PCR enzymes.

SETUP GUIDELINES

DNA Template

- Use high quality, purified DNA templates whenever possible. Please refer to specific product information for amplification from unpurified DNA (e.g., colony PCR or direct PCR).
- For low complexity templates (e.g., plasmid, lambda, BAC DNA), use 1 pg-10 ng of DNA per 50 µl reaction
- For higher complexity templates (e.g., genomic DNA), use 1 ng–1 μg of DNA per 50 μl reaction
- Higher DNA concentrations tend to decrease amplicon specificity, particularly for high numbers of cycles

Primers

- Primers should typically be 20–30 nucleotides in length, with 40–60% GC Content
- Primer Tm values should be determined with NEB's Tm Calculator (TmCalculator.neb.com)
- Primer pairs should have Tm values that are within 5°C
- Avoid secondary structure (e.g., hairpins) within each primer and potential dimerization between the primers
- Higher than recommended primer concentrations may decrease specificity
- When engineering restriction sites onto the end of primers, 6 nucleotides should be added 5' to the site
- Annealing temperatures should be determined according to specific enzyme recommendations.
 Please note that Q5® and Phusion® * annealing temperature recommendations are unique.
- Final concentration of each primer should be 0.05–1 µM in the reaction. Please refer to the more detailed recommendations for each specific enzyme.
- When amplifying products > 20 kb in size, primers should be ≥ 24 nucleotides in length with a GC content above 50% and matched Tm values above 60°C
- To help eliminate primer degradation and subsequent non-specific product formation, use a hot-start enzyme (e.g., One Taq® Hot Start DNA Polymerase or Q5 Hot Start High-Fidelity DNA Polymerase)

Magnesium Concentration

- Optimal Mg²⁺ concentration is usually 1.5–2.0 mM for most PCR polymerases
- Most PCR buffers provided by NEB already contain sufficient levels of Mg²⁺ at 1X concentrations.
- NEB offers a variety of Mg-free reaction buffers to which supplemental Mg²⁺ can be added for applications that require complete control over Mg²⁺ concentration
- Further optimization of Mg²⁺ concentration can be done in 0.2-1 mM increments, if necessary.
 For some specific applications, the enzyme may require as much as 6 mM Mg²⁺ in the reaction.
- Excess Mg²⁺ may lead to spurious amplification; Insufficient Mg²⁺ concentrations may cause reaction failure

Deoxynucleotides

- Ideal dNTP concentration is typically 200 µM of each, however, some enzymes may require as much as 400 µM each. Please refer to specific product literature for more detailed recommendations.
- Excess dNTPs can chelate Mg²⁺ and inhibit the polymerase
- Lower dNTP concentration can increase fidelity, however, yield is often reduced
- The presence of uracil in the primer, template, or deoxynucleotide mix will cause reaction failure when using archaeal PCR polymerases. Use One Taq or Taq DNA Polymerases for these applications.

Enzyme Concentration

- Optimal enzyme concentration in the reaction is specific to each polymerase. Please see the product literature for specific recommendations.
- In general, excess enzyme can lead to amplification failure, particularly when amplifying longer fragments

Starting Reactions

- Unless using a hot start enzyme (e.g., One Taq Hot Start DNA Polymerase or Q5 Hot Start High-Fidelity DNA Polymerase), assemble all reaction components on ice
- Add the polymerase last, whenever possible
- Transfer reactions to a thermocycler that has been pre-heated to the denaturation temperature.
 Please note that pre-heating the thermocycler is not necessary when using a hot start enzyme (e.g., One Taq Hot Start DNA Polymerase or Q5 Hot Start High-Fidelity DNA Polymerase).

CYCLING GUIDELINES

Denaturation

- Optimal denaturation temperature ranges from 94°–98°C and is specific to the polymerase in the reaction. Please refer to product information for recommended conditions.
- Avoid longer or higher temperature incubations unless required due to high GC content of the template
- For most PCR polymerases, denaturation of 5–30 seconds is recommended during cycling
- NEB's aptamer-based hot start enzymes do not require additional denaturation steps to activate the enzymes

Annealing

- Primer Tm values should be determined using the NEB Tm Calculator (TmCalculator.neb.com)
- For PCR polymerases other than Q5 High-Fidelity DNA Polymerase or Phusion High-Fidelity DNA Polymerase*, annealing temperatures are usually set at 2°-5°C below the lowest Tm of the primer pair
- When using Q5 High-Fidelity DNA Polymerase or Phusion High-Fidelity DNA Polymerase*, annealing temperatures should be set at 0°-3°C above the lowest Tm of the primer pair. Please refer to the product literature for detailed recommendations.
- Non-specific product formation can often be avoided by optimizing the annealing temperature or by switching to a hot start enzyme (e.g., One Taq Hot Start DNA Polymerase or Q5 Hot Start High-Fidelity DNA Polymerase)
- Annealing temperatures can be optimized by doing a temperature gradient PCR, starting at 5°C below the lowest Tm of the primer pair
- Ideally, primer Tm values should be less than the extension temperature. However, if Tm values are calculated to be greater than the extension temperature, a two-step PCR program (combining annealing and extension into one step) can be employed.

Extension

- Extension temperature recommendations range from 65°-72°C and are specific to each PCR polymerase. Please refer to the product literature for specific recommendations.
- Extension rates are specific to each PCR polymerase. In general, extension rates range from 15–60 seconds per kb. Please refer to the recommendations for each specific product.
- Longer than recommended extension times can result in higher error rates, spurious banding patterns and/or reduction of amplicon yields

^{*} Phusion DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific. Phusion® is a registered trademark and property of Thermo Fisher Scientific.

PCR Troubleshooting Guide

The following guide can be used to troubleshoot PCR reactions. Additional tips for optimizing reactions can be found in the technical reference section of our website, **www.neb.com**.

PROBLEM	POSSIBLE CAUSE	SOLUTION
	Low fidelity polymerase	Choose a higher fidelity polymerase such as Q5 High-Fidelity (NEB #M0491) or Phusion (NEB #M0530)* DNA Polymerases
	Suboptimal reaction conditions	Reduce number of cycles Decrease extension time
Sequence errors	Unbalanced nucleotide concentrations	Prepare fresh deoxynucleotide mixes
Sequence errors	Template DNA has been damaged	Start with a fresh template Try repairing DNA template with the PreCR® Repair Mix (NEB #M0309) Limit UV exposure time when analyzing or excising PCR product from the gel
	Desired sequence may be toxic to host	Clone into a non-expression vector Use a low-copy number cloning vector
	Incorrect annealing temperature	Recalculate primer Tm values using the NEB Tm calculator (TmCalculator.neb.com)
Incorrect product size	Mispriming	Verify that primers have no additional complementary regions within the template DNA
incorrect product size	Improper Mg ²⁺ concentration	Adjust Mg ²⁺ concentration in 0.2–1 mM increments
	Nuclease contamination	Repeat reactions using fresh solutions
	Incorrect annealing temperature	Recalculate primer Tm values using the NEB Tm calculator (TmCalculator.neb.com) Test an annealing temperature gradient, starting at 5°C below the lower Tm of the primer pair
	Poor primer design	Check specific product literature for recommended primer design Verify that primers are non-complementary, both internally and to each other Increase length of primer
	Poor primer specificity	Verify that oligos are complementary to proper target sequence
No product	Insufficient primer concentration	 Primer concentration can range from 0.05–1 μM in the reaction. Please see specific product literature for ideal conditions
	Missing reaction component	Repeat reaction setup
	Suboptimal reaction conditions	Optimize Mg ²⁺ concentration by testing 0.2–1 mM increments Thoroughly mix Mg ²⁺ solution and buffer prior to adding to the reaction Optimize annealing temperature by testing an annealing temperature gradient, starting at 5°C below the lower Tm of the primer pair
	Poor template quality	 Analyze DNA via gel electrophoresis before and after incubation with Mg²⁺ Check 260/280 ratio of DNA template
	Presence of inhibitor in reaction	Further purify starting template by alcohol precipitation, drop dialysis or commercial clean up kit Decrease sample volume
	Insufficient number of cycles	Rerun the reaction with more cycles
	Incorrect thermocycler programming	Check program, verify times and temperatures
	Inconsistent thermocycler block temperature	Test calibration of heating block
	Contamination of reaction tubes or solutions	* Autoclave empty reaction tubes prior to use to eliminate biological inhibitors * Prepare fresh solutions or use new reagents
	Complex template	Use Q5 High-Fidelity (NEB #M0491) or One Taq DNA Polymerase (NEB #M0482) For GC-rich templates, use One Taq DNA Polymerase (NEB #M0480) with One Taq GC Reaction Buffer (plus One Taq High GC Enhancer, if necessary) or Q5 High-Fidelity DNA Polymerase (NEB #M0491) with the High GC Enhancer For longer templates, we recommend LongAmp® Taq DNA Polymerase (NEB #M0323), Q5 or Q5 Hot Start High Fidelity DNA Polymerase (NEB #M0493)
	Premature replication	Use a hot start polymerase, such as Q5 Hot Start High-Fidelity (NEB #M0493) or One Taq Hot Start (NEB #M0481) DNA Polymerases Set up reactions on ice using chilled components and add samples to thermocycler preheated to the denaturation temperature
	Primer annealing temperature too low	Recalculate primer Tm values using the NEB Tm Calculator (TmCalculator.neb.com) Increase annealing temperature
	Incorrect Mg ²⁺ concentration	• Adjust Mg ²⁺ in 0.2–1 mM increments
Multiple or non-specific products	Poor primer design	Check specific product literature for recommended primer design Verify that primers are non-complementary, both internally and to each other Increase length of primer Avoid GC-rich 3´ ends
	Excess primer	• Primer concentration can range from 0.05–1 μM in the reaction. Please see specific product literature for ideal conditions.
	Contamination with exogenous DNA	Use positive displacement pipettes or non-aerosol tips Set-up dedicated work area and pipettor for reaction setup Wear gloves during reaction setup
	Incorrect template concentration	 For low complexity templates (e.g., plasmid, lambda, BAC DNA), use 1 pg—10 ng of DNA per 50 μl reaction For higher complexity templates (e.g., genomic DNA), use 1 ng—1 μg of DNA per 50 μl reaction

^{*} Phusion DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific. Phusion® is a registered trademark and property of Thermo Fisher Scientific.

Optimization Tips for Luna® qPCR

TIPS FOR OPTIMIZATION

New England Biolabs provides Luna products for your qPCR and RT-qPCR experiments. For more information on these products, visit **LUNAqPCR.com**. The following tips can be used to help optimize qPCR. For RT-qPCR guidelines, please see page 340.

Target Selection

- Short PCR amplicons, ranging from 70 to 200 bp, are recommended for maximum PCR efficiency
- Target sequences should ideally have a GC content of 40–60%
- Avoid highly repetitive sequences when possible

DNA Template

- Use high quality, purified DNA templates whenever possible. Luna qPCR is compatible with DNA samples prepared through typical nucleic acid purification methods.
- Template dilutions should be freshly prepared in either TE or water for each gPCR experiment
- Generally, useful concentrations of standard and unknown material will be in the range of 10^e copies to 1 copy. For gDNA samples from large genomes, (e.g., human, mouse) a range of 50–1 pg of gDNA is typical. For small genomes, adjust as necessary using 10^e–1 copy input as an approximate range. Note that for dilutions in the single-copy range, some samples will contain multiple copies and some will have none, as defined by the Poisson distribution.
- To generate cDNA, use of the LunaScript® RT SuperMix Kit (NEB #E3010) is recommended.
 Up to 1 µg total RNA, 1 µg mRNA or 100 ng specific RNA can be used in a 20 µl reaction.
- cDNA does not need to be purified before addition to the Luna reaction but should be diluted at least 1:20 before addition to qPCR

Primers

- Primers should typically be 15–30 nucleotides in length
- · Ideal primer content is 40-60% GC
- · Primer Tm should be approximately 60°C
- Primer Tm calculation should be determined with NEB's TmCalculator (TmCalculator.neb.com) using the Hot Start Tag setting
- For best results in qPCR, primer pairs should have Tm values that are within 3°C
- Avoid secondary structure (e.g., hairpins) within each primer and potential dimerization between primers
- G homopolymer repeats ≥ 4 should be avoided

- Optimal primer concentration for dye-based experiments (250 nM) is lower than for probebased experiments (400 nM). If necessary, the primer concentration can be optimized between 100–500 nM for dye-based qPCR or 200–900 nM for probe-based experiments.
- Higher primer concentrations may increase secondary priming and create spurious amplification products
- When using primer design software, enter sufficient sequence around the area of interest to permit robust primer design and use search criteria that permit cross-reference against relevant sequence databases to avoid potential off-target amplification.
- For cDNA targets, it is advisable to design primers across known exon-exon junctions in order to prevent amplification from genomic DNA
- Primers designed to target intronic regions can ensure amplification exclusively from genomic DNA

Hydrolysis Probes

- Probes should typically be 15–30 nucleotides in length to ensure sufficient quenching of the fluorophore
- The optimal probe concentration is 200 nM but may be optimized between 100 to 500 nM
- Both single or double-quenched probes may be used
- In general, non-fluorescence quenchers result in better signal-to-noise ratio than fluorescence quenchers
- Ideal probe content is 40-60% GC
- The probe Tm should be 5–10°C higher than the Tm of the primers to ensure all targeted sequences are saturated with probe prior to amplification by the primers
- Probes may be designed to anneal to either the sense or antisense strand
- Generally, probes should be designed to anneal in close proximity to either the forward or reverse primer without overlapping
- Avoid a 5´-G base which is known to quench 5´-fluorophores

Multiplexing

- Avoid primer/probe combinations that contain complementary sequences, and ensure target sequences do not overlap
- Probes should be designed such that each amplicon has a unique fluorophore for detection
- Select fluorophores based on the detection capabilities of the available real-time PCR instrument

- The emission spectra of the reporter fluorophores should not overlap
- Test each primer/probe combination in a singleplex reaction to establish a performance baseline. Ensure C_q values are similar when conducting the multiplex qPCR.
- Pair dim fluorescence dyes with high abundance targets and bright dyes with low abundance targets
- Optimization may require lower primer/probe concentrations to be used for high copy targets along with higher concentrations for low copy targets

Cycling Conditions

- Generally, best performance is achieved using the cycling conditions provided in the manual
- Longer amplicons (> 400 bp) can be used but may require optimization of extension times
- Due to the hot start nature of the polymerase, it is not necessary to preheat the thermocycler prior to use
- Select the "Fast" ramp speed where applicable (e.g., Applied Biosystems QuantStudio®)
- Amplification for 40 cycles is sufficient for most applications, but for very low input samples 45 cycles may be used

Reaction Setup

- For best results, keep reactions on ice prior to thermocycling
- A reaction volume of 20 μl is recommended for 96-well plates while a reaction volume of 10 μl is recommended for 384-well plates
- Reactions should be carried out in triplicate for each sample
- For each amplicon, ensure to include no template controls (NTC)
- To prevent carry-over contamination, treat reactions with 0.2 units/µl Antarctic Thermolabile UDG (NEB #M0372) for 10 minutes at room temperature prior to thermocycling
- The Luna reference dye supports broad instrument compatibility (High-ROX, Low-ROX, ROX-independent) so no additional ROX is required for normalization

Assay Performance

- Ensure 90–110% PCR efficiency for the assay over at least three log₁₀ dilutions of template
- Linearity over the dynamic range (R²) should ideally be ≥ 0.99
- Target specificity should be confirmed by product size, sequencing or melt-curve analysis

Optimization Tips for Luna One-Step RT-qPCR

TIPS FOR OPTIMIZATION

New England Biolabs provides Luna products for your qPCR and RT-qPCR experiments. For more information on these products, visit **LUNAqPCR.com**. The following tips can be used to help optimize your one-step RT-qPCR. For qPCR guidelines (DNA/cDNA starting material), please see page 339.

Target Selection

- Short PCR amplicons, ranging from 70 to 200 bp, are recommended for maximum PCR efficiency
- Target sequences should ideally have a GC content of 40–60%
- Avoid highly repetitive sequences when possible
- Target sequences containing significant secondary structure should be avoided

RNA Template

- Use high quality, purified RNA templates whenever possible. Luna qPCR is compatible with RNA samples prepared through typical nucleic acid purification methods.
- Prepared RNA should be stored in an EDTAcontaining buffer (e.g., 1X TE) for longterm stability
- Template dilutions should be freshly prepared in either TE or water for each qPCR experiment
- Treatment of RNA samples with DNase I (NEB #M0303) may minimize amplification from genomic DNA contamination
- Generally, useful concentrations of standard and unknown material will be in the range of 10⁸ copies to 10 copies. Note that for dilutions in the single-copy range, some samples will contain multiple copies and some will have none, as defined by the Poisson distribution. For total RNA, Luna One-Step Kits can provide linear quantitation over an 8-order input range of 1 µg−0.1 pg. For most targets, a standard input range of 100 ng−10 pg total RNA is recommended. For purified mRNA, input of ≤ 100 ng is recommended. For in vitrotranscribed RNA, input of ≤ 10⁹ copies is recommended.

Primers

- Primers should typically be 15–30 nucleotides in length
- Ideal primer content is 40–60% GC
- Primer Tm should be approximately 60°C
- Primer Tm calculation should be determined with NEB's TmCalculator. (TmCalculator.neb.com) using the Hot Start Tag setting.
- For best results in qPCR, primer pairs should have Tm values that are within 3°C
- Avoid secondary structure (e.g., hairpins) within each primer and potential dimerization between primers
- G homopolymer repeats ≥ 4 should be avoided

- The optimal primer concentration for dye-based experiments and probe-based experiments is 400 nM. If necessary, the primer concentration can be optimized between 100–900 nM.
- Higher primer concentrations may increase secondary priming and create spurious amplification products
- When using primer design software, enter sufficient sequence around the area of interest to permit robust primer design and use search criteria that permit cross-reference against relevant sequence databases to avoid potential off-target amplification
- It is advisable to design primers across known exon-exon junctions in order to prevent amplification from genomic DNA

Hydrolysis Probes

- Probes should typically be 15–30 nucleotides in length to ensure sufficient quenching of the fluorophore
- The optimal probe concentration is 200 nM but may be optimized between 100 to 500 nM
- Both single or double-quenched probes may
 be used.
- In general, non-fluorescence quenchers result in better signal-to-noise ratio than fluorescence quenchers
- Ideal probe content is 40-60% GC
- The probe Tm should be 5–10°C higher than the Tm of the primers to ensure all targeted sequences are saturated with probe prior to amplification by the primers
- Probes may be designed to anneal to either the sense or antisense strand
- Generally, probes should be designed to anneal in close proximity to either the forward or reverse primer without overlapping
- Avoid a 5´-G base which is known to quench 5´-fluorophores

Multiplexing

- Avoid primer/probe combinations that contain complementary sequences, and ensure target sequences do not overlap
- Probes should be designed such that each amplicon has a unique fluorophore for detection
- Select fluorophores based on the detection capabilities of the available real-time PCR instrument
- The emission spectra of the reporter fluorophores should not overlap
- Test each primer/probe combination in a singleplex reaction to establish a performance baseline. Ensure C_q values are similar when conducting the multiplex qPCR.
- Pair dim fluorescence dyes with high abundance targets and bright dyes with low abundance targets

 Optimization may require lower primer/probe concentrations to be used for high copy targets along with higher concentrations for low copy targets

Reverse Transcription

- The default reverse transcription temperature is 55°C
- For difficult targets, the temperature of reverse transcription may be increased to 60°C for 10 minutes
- Due to the WarmStart feature of the Luna RT, reverse transcription temperatures lower than 50°C are not recommended

Cycling Conditions

- Generally, best performance is achieved using the cycling conditions provided in the manual
- Longer amplicons (> 400 bp) can be used but may require optimization of extension times
- Due to the dual WarmStart/Hot Start feature of the Luna kits, it is not necessary to preheat the thermocycler prior to use
- Select the "Fast" ramp speed where applicable (e.g., Applied Biosystems QuantStudio).
- Amplification for 40 cycles is sufficient for most applications, but for very low input samples 45 cycles may be used

Reaction Setup

- For best results, keep reactions on ice prior to thermocycling
- A reaction volume of 20 µl is recommended for 96-well plates while a reaction volume of 10 µl is recommended for 384-well plates
- Reactions should be carried out in triplicate for each sample
- For each amplicon, ensure to include no template controls (NTC)
- A no Luna RT control should be conducted to guarantee amplification is specific for RNA input and not due to genomic DNA contamination
- To prevent carry-over contamination, treat reactions with 0.2 units/µl Antarctic Thermolabile UDG (NEB #M0372) for 10 minutes at room temperature prior to thermocycling
- The Luna reference dye supports broad instrument compatibility (High-ROX, Low-ROX, ROX-independent) so no additional ROX is required for normalization

Assay Performance

- Ensure 90–110% PCR efficiency for the assay over at least three log₁₀ dilutions of template.
- Linearity over the dynamic range (R²) should ideally be ≥ 0.99
- Target specificity should be confirmed by product size, sequencing or melt-curve analysis

Luna qPCR Troubleshooting Guide

PROBLEM	PROBABLE CAUSE(S)	SOLUTION(S)	
qPCR traces show low	Reagent omitted from qPCR assay	Medical short the control of the control of	
or no amplification	Reagent added improperly to qPCR assay	Verify all steps of the protocol were followed correctly	
	Incorrect cycling protocol	Refer to the proper qPCR cycling protocol in product manual	
	Incorrect channel selected for the qPCR thermal cycler	Verify correct optical settings on the qPCR instrument	
		Confirm the expiration dates of the kit reagents	
	DNA template or reagents are contaminated or degraded	Verify proper storage conditions provided in this user manual	
	211/1 complate of roagonic are contaminated of dograded	Rerun the qPCR assay with fresh reagents	
		Confirm template input amount	
Inconsistent qPCR traces for triplicate data	Improper pipetting during qPCR assay set-up	Ensure proper pipetting techniques	
ioi iripiicate data	qPCR plate film has lost its seal, causing evaporation in the well. The resulting qPCR trace may show significantly different	Ensure the qPCR plate is properly sealed before inserting into the qPCR thermal cycler.	
	fluorescence values relative to its replicates	Exclude problematic trace(s) from data analysis.	
	Poor mixing of reagents during qPCR set-up	Make sure all reagents are properly mixed after thawing them	
		Avoid bubbles in the qPCR plate	
	Bubbles cause an abnormal qPCR trace	Centrifuge the qPCR plate prior to running it in the thermal cycler	
		Exclude problematic trace(s) from data analysis	
DNA standard curve has a poor correlation coefficient/efficiency of the DNA standard curve falls outside the 90–110% range	Presence of outlying qPCR traces	Omit data produced by qPCR traces that are clearly outliers caused by bubbles, plate sealing issues, or other experimental problems	
	Improper pipetting during qPCR assay set-up	Ensure that proper pipetting techniques are used	
	Reaction conditions are incorrect	Verify that all steps of the protocol were followed correctly	
	Bubbles cause an abnormal qPCR trace	Avoid bubbles in the qPCR plate	
	Dubbles cause all abilibilital qr on trace	Centrifuge the qPCR plate prior to running it in the thermal cycler	
	Poor mixing of reagents	After thawing, make sure all reagents are properly mixed	
		• Ensure the threshold is set in the exponential region of qPCR traces	
	Threshold is improperly set for the qPCR traces	Refer to the real-time instrument user manual to manually set an appropriate threshold	
Melt curve shows different peaks for	Non-template amplification is occurring	Compare melt curve of NTC to samples	
low input samples	Infrequently, denaturation of a single species can occur in a	Redesign primers with a Tm of 60°C or use our Tm calculator to determine the optimal annealing temperature of the primers	
	biphasic manner, resulting in two peaks	Perform a primer matrix analysis to determine optimal primer concentrations	
No template control qPCR trace		Replace all stocks and reagents	
shows amplification, NTC C _q is close to or overlapping lower	Reagents are contaminated with carried-over products of previous qPCR (melt curve of NTC matches melt curve of higher	Clean equipment and setup area with a 10% chlorine bleach	
copy standards	input standards)	• Consider use of 0.2 U/µl Antarctic Thermolabile UDG to eliminate carryover products	
	Primers produce non-specific amplification (melt curve of NTC does not match melt curve of higher input standards)	• Redesign primers with a Tm of 60°C or use qPCR primer design software	



Mark has been with NEB for over 35 years and currently serves as our Senior Network Engineer, keeping our communications running smoothly.

Luna One-Step RT-qPCR Troubleshooting Guide

PROBLEM	PROBABLE CAUSE(S)	SOLUTION(S)
	Incorrect RT step temperature or RT step omitted	• For typical use, a 55°C RT step temperature is optimal for the Luna WarmStart Reverse Transcriptase.
	Incorrect cycling protocol	Refer to the proper RT-qPCR cycling protocol in product manual
	Reagent omitted from RT-qPCR assay	
qPCR traces show low or no amplification	Reagent added improperly to RT-qPCR assay	Verify all steps of the protocol were followed correctly
	Incorrect channel selected for the qPCR thermal cycler	Verify correct optical settings on the qPCR instrument
		Prepare high quality RNA without RNase/DNase contamination
		Confirm template input amount
	RNA template or reagents are contaminated or degraded	Confirm the expiration dates of the kit reagents
		Verify proper storage conditions provided in product manual
		Rerun the RT-qPCR assay with fresh reagents
	Improper pipetting during RT-qPCR assay set-up	Ensure proper pipetting techniques
	qPCR plate film has lost its seal, causing evaporation in the well. The resulting qPCR trace may show significantly different	Ensure the qPCR plate is properly sealed before inserting into the qPCR thermal cycler
Inconsistent qPCR traces	fluorescence values relative to its replicates.	Exclude problematic trace(s) from data analysis
for triplicate data	Poor mixing of reagents during RT-qPCR set-up	Make sure all reagents are properly mixed after thawing them
		Avoid bubbles in the qPCR plate
	Bubbles cause an abnormal qPCR trace	Centrifuge the qPCR plate prior to running it in the thermal cycler
		Exclude problematic trace(s) from data analysis
		Refer to the proper RT-qPCR cycling protocol in product manual
	Cycling protocol is incorrect	Use a 55°C RT step temperature
		• For ABI instruments, use a 1 minute 60°C annealing/extension step
	Presence of outlying qPCR traces	Omit data produced by qPCR traces that are clearly outliers caused by bubbles, plate sealing issues, or other experimental problems
Standard curve has a poor correlation	Improper pipetting during RT-qPCR assay set-up	Ensure that proper pipetting techniques are used
coefficient/efficiency of the standard	Reaction conditions are incorrect	Verify that all steps of the protocol were followed correctly
curve falls outside the 90–110% range	Bubbles cause an abnormal gPCR trace	Avoid bubbles in the qPCR plate
	bubbles cause an abhornial qi on trace	Centrifuge the qPCR plate prior to running it in the thermal cycler
	Poor mixing of reagents	After thawing, make sure all reagents are properly mixed
		Ensure the threshold is set in the exponential region of qPCR traces
	Threshold is improperly set for the qPCR traces	Refer to the real-time instrument user manual to manually set an appropriate threshold
	Non-template amplification is occurring	Compare melt curve of NTC to samples
Melt curve shows different peaks for low input samples	Non-template amplification is occurring Infrequently, denaturation of a single species can occur in a	Redesign primers with a Tm of 60°C or use our Tm calculator to determine the optimal annealing temperature of the primers
	biphasic manner, resulting in two peaks	Perform a primer matrix analysis to determine optimal primer concentrations
	Reagents are contaminated with carried-over products of previous	Replace all stocks and reagents
No template control qPCR trace shows	qPCR (Melt curve of NTC matches melt curve of higher input	Clean equipment and setup area with a 10% chlorine bleach
amplification/NTC C _q is close to or overlapping lower copy standards	standards)	Consider use of 0.2 U/µl Antarctic Thermolabile UDG to eliminate carryover products
	Primers produce non-specific amplification (Melt curve of NTC does not match melt curve of higher input standards)	Redesign primers with a Tm of 60°C or use qPCR primer design software
Amplification in No-RT control	RNA is contaminated with genomic DNA	Treat sample with DNase I
Ampinication in No-N1 Control	The Contaminated With goldening Divi	Redesign amplicon to span exon-exon junction

Cleavage Close to the End of DNA Fragments

Annealed 5´ FAM-labeled oligos were incubated with the indicated enzyme (10 units/1pmol oligo) for 60 minutes at the recommended incubation temperature and NEBuffer. The digest was run on a TBE acrylamide gel and analyzed by fluorescent imaging. The double stranded oligos were designed to have the indicated number of base pairs from the end followed by the recognition sequence and an additional 12 bases. In some cases asymmetric cleavage was observed and interpreted as a negative result. Asymmetric cleavage decreased with increasing base pairs from the end.

	BASE PAIRS FROM END				
ENZYME	1 bp	2 bp	3 bp	4 bp	5 bp
Acil		+	+	++	+++
Agel	+++	+++	+++	+++	+++
Agel-HF	+++	+++	+++	+++	+++
Alul	++	+++	+++	+++	+++
Apal Ascl	+++	+++	+++	+++	+++
AvrII	+++	+++	+++	+++	+++
	++	++	+++	+++	+++
BamHI	+	++	+++	+++	+++
BamHI-HF	+	+	+++	+++	+++
BbsI-HF	+++	+++	+++	+++	+++
BcII-HF	-	_	+++	+++	+++
BgIII	++	+++	+++	+++	+++
Bmtl	+++	+++	+++	+++	+++
Bmtl-HF	+++	+++	+++	+++	+++
Bsal	+++	+++	+++	+++	+++
Bsal-HFv2	+++	+++	+++	+++	+++
BsiWI	++	+++	+++	+++	+++
BsiWI-HF	+++	+++	+++	+++	+++
BsmBl	+++	+++	+++	+++	+++
BsrGI	+++	+++	+++	+++	+++
BssHII	+	+++	+++	+++	+++
BstZ17I-HF	+	+++	+++	+++	+++
Clal	_	-	+	+++	+++
Ddel	+++	+++	+++	+++	+++
Dpnl	-	++	++	NT	NT
DrallI-HF	+++	+++	+++	+++	+++
Eagl	++	+++	+++	+++	+++
Eagl-HF	+	+++	+++	+++	+++
EcoRI	+	+	++	++	+++
EcoRI-HF	+	+	++	+++	+++
EcoRV	++	++	++	++	+++
EcoRV-HF	+	++	++	++	+++
Esp3I	+++	+++	+++	+++	+++
Fsel	+	++	+++	+++	+++
HindIII	_	+	+++	+++	+++
HindIII-HF	-	+	+++	+++	+++
Hpal	+++	+++	+++	+++	+++
Kpnl	+	+++	+++	+++	+++
KpnI-HF	+	+++	+++	+++	+++
Mfel	+	++	+++	+++	+++
Mfel-HF	+	++	+++	+++	+++
Mlul	+	++	+++	+++	+++
Msel	+++	+++	+++	+++	+++
IVIOUI	T T T	T T T	+++	T T T	+++

Note: As a general rule and for enzymes not listed below, 6 base pairs should be added on on either side of the recognition site to cleave efficiently.

The extra bases should be chosen so that palindromes and primer dimers are not formed. In most cases there is no requirement for specific bases.

Chart Legend

- 0% + 0-20% ++ 20-50% +++ 50-100%

NT not tested

	BASE PAIRS FROM END					
ENZYME	1 bp	2 bo	3 bo	4 bp	5 bp	
Ncol		++	+++	+++	+++	
Ncol-HF	+	++	+++	+++	+++	
Ndel	+	+	+++	+++	+++	
Nhel	+	++	+++	+++	+++	
Nhel-HF	++	++	+++	+++	+++	
NIaIII	++	+++	+++	+++	+++	
Notl	++	++	++	++	++	
NotI-HF	++	++	++	++	++	
Nsil	+	+	+++	+++	+++	
Nspl	-	-	+	+	+++	
Pacl	+++	+++	+++	+++	+++	
Pcil	+++	+++	+++	+++	+++	
Pmel	+++	+++	+++	+++	+++	
Pstl	+	+++	+++	+++	+++	
PstI-HF	++	+++	+++	+++	+++	
Pvul	+++	+++	+++	+++	+++	
Pvul-HF	+++	+++	+++	+++	+++	
Pvull	++	++	++	+++	+++	
Pvull-HF	-	++	++	+++	+++	
Rsal	+	+++	+++	+++	+++	
Sacl	-	++	+++	+++	+++	
SacI-HF	-	+	+++	+++	+++	
SacII	+++	+++	+++	+++	+++	
Sall	-	++	+++	+++	+++	
Sall-HF	-	++	+++	+++	+++	
Sapl	+++	+++	+++	+++	+++	
Sau3AI	+++	+++	+++	+++	+++	
Sbfl	++	+++	+++	+++	+++	
Sbfl-HF	++	+++	+++	+++	+++	
Scal-HF	+	+++	+++	+++	+++	
Sfil	+++	+++	+++	+++	+++	
Smal	+++	+++	+++	+++	+++	
Spel	+	++	++	++	++	
Spel-HF	+	++	++	++	++	
Sphl	+++	+++	+++	+++	+++	
SphI-HF	++	++	+++	+++	+++	
Sspl	+	+++	+++	+++	+++	
SspI-HF	+	+++	+++	+++	+++	
Stul	+++	+++	+++	+++	+++	
Styl	+	++	+++	+++	+++	
Styl-HF	+	+++	+++	+++	+++	
Xbal	++	++	++	++	++	
Xhol	++	++	++	+++	+++	
Xmal	+++	+++	+++	+++	+++	

Activity of Restriction Enzymes in PCR Buffers

Frequently, a PCR product must be digested with restriction enzymes. For convenience, digestion can be performed directly in the PCR mix without any purification of the DNA. This table summarizes the activity of restriction enzymes on the DNA in *Taq*, Phusion*, One *Taq* and LongAmp *Taq* PCR mixes. 50 µl reactions containing 5 units of restriction enzyme were incubated at the appropriate temperature for 1 hour in a PCR mix containing the following: 1 µg DNA, 1 unit of DNA Polymerase and 1X ThermoPol Reaction Buffer, Standard *Taq* Reaction Buffer, Phusion HF Buffer, One *Taq* Standard Reaction Buffer or LongAmp *Taq* Reaction Buffer. Reactions were supplemented with 200 µM dNTPs. Enzyme activity was analyzed by gel electrophoresis.

Notes: The polymerase is still active and can alter the ends of DNA fragments after cleavage, affecting subsequent ligation. Primers containing the restriction

Taq IN Thermopol Q5 IN IN LONGAME FN7YMF **RXN BUFFEF Q5 BUFFER HF BUFFER RXN BUFFER** Taa RXN BUFFEI AatII <++ < + ++ Accl <++ <+ <+ + + ++ + +Acc65I <+ +++ Acil ++ ++ +++ +++ AcII +++ <+ <+ +++ +++ Acul +++ <+ ++ +++ +++ Afel +++ < + ++ +++ +++ AfIII <+ <+ <+ AfIIII <+ +++ <+ <+ Agel +++ +++ +++ <+ Agel-HF +++ <+ +++ +++ AhdI <+ _ <+ <+ Alel-v2 Alul +++ ++++++ +++ Alwl <+ <+ <+ <+ AlwNl <++ +++ <+ + Anal +++ < + < + < + ApaLI +++ <+ <+ +++ ApeKI <++ ++ +++ <+ Apol +++ ++ +++ ++ +++ Apol-HF + + +++ + + ++++ AscI +++ <+ <+ <+ Asel +++ <+ ++ AsiSI + + +<+ ++ +++ + + +Aval Avall +++ <+ ++ +++ +++ AvrII BaeGI +++ <+ +++ +++ +++ Bael < + ++ <+ <+ BamHI +++ <+ +++ +++ +++ BamHI-HF +++ < + <+ ++ Banl +++ <+ +++ +++ +++ Banll +++ +++ +++ +++ <+ Bbsl +++ <+ <+ +++ +++ BbsI-HF BbvCI +++ <+ <+ Bbvl Bccl <+ <+ <+ <+ <+ BceAl <+ <+ Bcgl <+ <+ + ++ ++ **BciVI** <+ BcII +++ ++ +++ +++ BcII-HF +++ + BcoDI <+ + Bfal <+ BfuAl <++ <+ BgII ++ < + <+ <+ BgIII <+ ++ <+ <+ Blpl <++ <+ <+ <+ BmgBI ++ + <+ <+ Bmrl <++ <+ +++ +++ +++ Bmtl +++ <+ ++ +++ +++ Bmtl-HF ++ <+ +++ ++ Bpml <+ <+ +++ <++ <++ BnuEl +++ ++ <++ <++ Bpu10I <+ <+ +++ ++ +++ BsaAl +++ ++ +++ + + ++++ BsaBl ++

enzyme recognition site can act as competitive inhibitors in the cleavage reaction. The use of restriction enzymes under non-optimal conditions may increase the likelihood of star activity. If any problems are encountered, the DNA should be purified by spin column or phenol/chloroform extraction followed by alcohol precipitation.

* Phusion DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Blolabs, linc. under agreement with, and under the performance specifications of Thermo Fisher Scientific. Pulsion® is a registered trademark and property of Thermo Fisher Scientific.

Chart Legend

Cleavage in extension mix with 5 units of enzyme:

+ + + complete cleavage + + ~ 50% cleavage + ~ 25% cleavage - no cleavage

** It has been shown that the addition of 1X Restriction Enzyme Buffer may help to improve the ability of some enzymes to cleave.

ENZYME	<i>Taq</i> IN Thermopol RXN Buffer	Q5 IN Q5 Buffer**	PHUSION In Phusion HF Buffer	ONE <i>Taq</i> In one <i>taq</i> RXN Buffer	LONGAMP <i>Taq</i> In Longamp <i>Taq</i> RXN Buffer
BsaHI	+++	+	+++	+++	+++
Bsal-HFv2	+	<+	+	+	++
BsaJI	+++	<+	++	+++	+++
BsaWl	<++	<+	++	+	+
BsaXI	<++	<+	<+	<+	<+
BseRI	+++	<+	++	++	+
BseYI	+++	++	++	+++	+++
Bsgl	<+	<+	+	<+	<+
BsiEl	+++	<+	++	++	++
BsiHKAI	-	++	+	-	-
BsiWl	+++	<+	+++	+++	+++
BsiWI-HF	-	-	-	-	-
BsII	+++	++	+++	+++	+++
BsmAl	+++	++	+++	<+	<+
BsmBI	<++	+	++	<+	<+
BsmFI	<+	+++	++	+	+
Bsml	+++	+	<+	+++	+
BsoBl	+++	+++	+++	++	+++
BspCNI	<+	<+	+	-	-
BspDI	<++	<+	++	+++	+++
BspEl	_	<+	<+	-	-
BspHI	+++	<+	+++	+++	+++
Bsp1286I	<+	<+	<+	<+	<+
BspMI	+++	<+	++	<+	<+
BspQI	+	++	+++	+++	+++
BsrBI	+++	<+	+	+++	+++
BsrDI	<+	<+	+	<+	<+
BsrFI-v2	<+	_	_	_	-
BsrGI	<+	+	+++	<+	+++
Bsrl	+++	<+	+++	++	+++
BssHII	+++	<+	+	+++	+++
BssSI-v2	+++	_	+	+++	+++
BstAPI	+++	<+	++	+++	+++
BstBI	+++	++	+++	+++	+++
BstEII	+++	<+	<+	+++	+++
BstEII-HF	+++	<+	<+	++	++
BstNI	+++	<+	<+	<+	<+
BstUI	+++	<+	<+	+++	+
BstXI	<++	+	+	+	<+
BstYI	+++	<+	<+	++	+
BstZ17I-HF	+++	_	+	+++	+++
Bsu36I	<+	<+	<+	<+	+
Btgl	+++	<+	+	<+	<+
BtgZI	+++	+	++	++	++
BtsI-v2	+++	_	+	+++	+++
BtsCl	+++	<+	<+	+++	+++
Cac8I	+++	<+	<+	+++	
Clal		<+			++
CspCl	++	_	<+ +	<+ <+	<+
CviAII	+++				+++
CviKI-1		<+	+	+++	
CviQI	+++	<+	++	+++	+++
Ddel	+++	+	+++	++	+++
	+++	++	+	+++	+++
Dpnl	+++	++	+++	++	++
DpnII	+++	++	+++	+++	++
Drall HE	+++	<+	+++	+++	+++
DrallI-HF	++	++	+++	++	++
Drdl	+++	<+	+++	+++	+++

ENZYME	<i>Taq</i> IN Thermopol RXN Buffer	Q5 IN Q5 BUFFER**	PHUSION In Phusion HF Buffer		LONGAMP <i>Taq</i> In Longamp <i>Taq</i> RXN BUFFER
Eael	+++	<+	_	<+	<+
Eagl Eagl-HF	<+	+++	+++	+++	+++
Earl	+	<+	+	++	++
Ecil	+++	< + + +	+++	+	<+
Eco53kl	+++	+ + <+	+++	<++ +++	<++ +++
EcoNI	+++	<+	+	+++	+++
EcoO109I	+++	<+	-	<+	+ + +
EcoP15I	<+	<+	+	<+	+
EcoRI	+	<+	+++	_	_
EcoRI-HF	+++	<+	+	+++	+++
EcoRV	<+	<+	+	-	<+
EcoRV-HF	+	<+	<+	+	++
Esp3I	+++	-	+++	+	+++
Fatl	++	<+	+++	<+	+++
Faul	+	<+	++	+++	++
Fnu4HI	+++	<+	<+	++	+
Fokl	+++	+	+	+++	+++
Fsel	+	<+	++	+++	-
Fspl	<++	<+	+	+	+
Haell	+++	<+	+++	+++	+++
HaellI	+++	<+	+++	+++	+++
Hgal	<+	<+	+	<++	<++
Hhal	+++	<+	+++	+++	+++
Hincll	+++	<+	<+	+++	+++
HindIII	+++	<+	+	++	+++
HindIII-HF	+++	<+	<+	+++	+++
Hinfl	+++	+++	+++	+	+++
HinP1I	+++	+	+++	+++	+++
Hpal	+++	<+	+++	+++	+++
Hpall	+++	<+	<+	<+	<+
Hphl	<++	<+	<+	<+	<+
HpyAV	+++	_	++	+	++
HpyCH4III	<++	<+	+	<++	<++
HpyCH4IV	+++	<+	<+	+++	+++
HpyCH4V	+++	<+	<+	+++	+++
Hpy99I	+++	_	+	<+	<+
Hpy188I	+++	<+	+	++	++
Hpy166II	+++	+	++	+++	+++
Hpv188III	+	<+	<+	+	<+
Kasl	+++	<+	+++	+++	_
Kpnl	+++	++	+	++	<+
Kpnl-HF	++	-	++	<+	<+
Mbol	+++	<+	+++	+++	+++
Mboll	+++	+	++	+	+
Mfel	+++	<+	<+	+++	+
Mfel-HF	+	_	_	+++	<+
MluCl	+	<+	<+	++	+
Mlul	+++	++	++	++	++
Mlul-HF	++	_	++	++	++
Mlyl	+++	+	++	<+	+
Mmel	<+	_	++	<+	<+
MnII	+++	+	+	+	+
Mscl	<+	<+	+	<+	<+
Msel	<+	<+	<+	<+	<+
MsII	+++	<+	+	+++	++
MspA1I	+++	<+	+++	++	+++
Mspl	+++	<+	+++	++	+++
Mwol	+++	+++	+++	++	+++
Nael	<+	<+	+	<+	<+
Narl	_	<+	++	+++	+++
Ncil	+++	<+	<+	+	<+
Ncol	+++	<+	+	++	++
Ncol-HF	+++	<+	_	++	+
Ndel	<++	++	+++	++	<+
NgoMIV	-	<+	+	<+	<+
Nhel	+++	<+	<+	+++	+++
Nhel-HF	+++	<+	-	++	++
NIalli	<+	<+	+	++	<+
NIaIV	+++	<+	+++	+++	+++
NmeAIII	<+	_	+++	<+	<+
Notl	++	<+	+	<+	<+
Notl-HF	+++	<+	<+	<+	+
1100 111		× r	<u> </u>	<u> </u>	т

	<i>Taq</i> IN		PHUSION	ONE <i>Taq</i>	LONGAMP Tag
ENZYME	THERMOPOL RXN BUFFER	Q5 IN Q5 Buffer**	IN PHUSION He buffer	IN ONE <i>Taq</i>	IN LONGAMP Tag RXN BUFFER
Nrul	++	+	+	++	++
Nrul-HF	++	-	-	+	-
Nsil	+++	+	+++	++	+
Nsil-HF	+++	++	+++	+++	+++
Nspl Pacl	+++	<+	<+	+++	++
Paci PaeR7I	+++	<+ <+	<+ <+	++	+++
Pcil	<+	<+	-	-	-
PfIFI	+++	<+	<+	<+	+
PfIMI	+	<+	+++	++	+++
Plel	+++	<+	<+	<+	<+
PluTI	+++	<+	+	+++	+++
Pmel Pmll	+++	<+	<+	+++	+++
PpuMI	+++	<+	+++	+++	+++
PshAl	+++	<+	<+	<+	<+
Psil	+++	<+	<+	<+	+++
PspGI	+++	+++	+++	+++	+++
Psp0MI	+++	<+	+	+++	+++
PspXI	+++	<+	++	+++	+++
Pstl-HF	++	+ <+	++	<+ ++	<+ +
Pvul	+++ <+	<+ <+	++	++	+ <+
Pvul-HF	+++	<+	+++	++	+++
Pvull	+++	<+	+	+++	+++
Pvull-HF	+	-	-	<+	<+
Rsal	+++	<+	++	+++	+++
RsrII	<++	-	-	<+	<+
Sacl-HF	+++	<+ <+	+ <+	++	++
SacII	+++	<+	+++	++	++
Sall	<+	+	++	-	-
Sall-HF	+	<+	+++	+	+++
Sapl	<++	<+	++	++	++
Sau3AI	+++	<+	<+	<+	<+
Sau96l	<++	+	+	+++	+++
Sbfl-HF	<++	<+	+	<+ <+	+++
Scal-HF	+	<+	<+	_	_
ScrFI	+++	+++	+++	+++	+++
SexAl	+++	<+	+++	+++	+++
SfaNI	-	<+	++	<++	<++
Sfcl	+++	<+	<+	+	+
Sfil Sfol	+++	-	-	+++	+++
SgrAl	+++	<+ <+	+++	+	+++
Smal	+++	<+	++	+++	+++
SmII	<+	<+	+	+	+
SnaBI	<+	<+	<+	+++	+++
Spel	+++	+	<+	+++	+++
Spel-HF	+++	-	<+	+++	+++
Sphl UE	+++	+	++	<+	<+
SphI-HF SrfI	+++	<+ <+	+++	+++	+++
Sspl-HF	++	<+	+	+++	+++
Stul	+++	<+	<+	+++	+++
StyD4I	<++	<+	+	<+	<+
Styl	<+	+	<+	<+	<+
Styl-HF	+	<+	<+	++	+++
Swal Taq ^a l	<+	<+	< + +	<+	+++
Tfil	+++	<+ <+	+ <+	+++	+++
Tsel	+++	+++	+++	+++	+++
Tsp45I	+++	-	-	+	<+
TspMI	+++	<+	+	+++	+++
TspRI	+	<+	<+	+++	+++
Tth1111	+++	<+	++	<+	+
Xbal	+++	-	<+	++	++
Xcml Xhol	+++	<+	+	+++	+++
Xmal	<+ +++	<+ <+	+++	++	+++
XmnI	+++	<+	<+	++	+++
Zral	+++	<+	<+	++	+

Getting Started with Molecular Cloning

Molecular cloning has traditionally used restriction enzymes to excise a fragment from source DNA, and to linearize a plasmid vector, while creating compatible ends. After purification, insert and vector are ligated to form a recombinant vector, which is transformed into an *E. coli* host. Alternatively, PCR can been used to generate both the vector and insert, which can be joined using a variety of techniques, such as standard DNA ligation, enzymatic joining using a recombinase or topoisomerase, or homologous recombination.

Regardless of the method chosen, the process can be made more efficient and successful by following good practices in the lab. The following tips will help improve the success of your cloning experiments.

1. Take the time to plan your experiments

Pay attention to the junction sequences and the effect on reading frames of any translated sequences. Check both the vector and insert for internal restriction sites (we recommend NEBcutter at **NEBcutter.neb.com**) prior to designing PCR primers that contain similar sites to those used for cloning. Verify that the antibiotic selective marker in the vector is compatible with the chosen host strain.

2. Start with clean DNA at the right concentration

Ensure that your source DNA is free of contaminants, including nucleases and unwanted enzymatic activities. Use commercially-available spin columns to purify starting DNA, (e.g., Monarch Plasmid Miniprep Kit, NEB #T1010 for DNA plasmids, Monarch PCR & DNA Cleanup Kit, NEB #T1030 for DNA Fragments). Completely remove solvents, such as phenol, chloroform and ethanol, prior to manipulation of the DNA. Elute DNA from the spin columns with salt-free buffer to prevent inhibition of the downstream steps, either restriction digestion or PCR amplification. Use a sufficient amount of DNA for the technique being used. Preparative restriction digests often require between 0.2 –2.0 µg, while single nanogram amounts are usually sufficient for DNA being used as a PCR template.

3. Perform your restriction digests carefully

The reaction volume should be compatible with the downstream step (e.g., smaller than the volume of the well of an agarose gel used to resolve the fragments). For a typical cloning reaction, this is often between $20-50\,\mu$ l. The volume of restriction enzyme(s) added should be no more than 10% of the total reaction volume, to ensure that the glycerol concentration stays below 5%; this is an important consideration to minimize star activity (unwanted cleavage).

4. Mind your ends

DNA ends prepared for cloning by restriction digest are ready for ligation without further modification, assuming the ends to be joined are compatible. If the ends are non-compatible, they can be modified using blunting reagents, phosphatases, etc.

DNA ends prepared by PCR for cloning may have a 3´ addition of a single adenine (A) residue following amplification using *Taq* DNA Polymerase (NEB #M0273). High-fidelity DNA polymerases, such as Q5 (NEB #M0491), leave blunt ends. PCR using standard commercial primers produces non-phosphorylated fragments, unless the primers were 5´ phosphorylated. The PCR product may need to be kinase treated to add a 5´ phosphate prior to ligation with a dephosphorylated vector.

5. Clean up your DNA prior to vector:insert joining

This can be done with gel electrophoresis or column purification (e.g., Monarch PCR & DNA Cleanup Kit, NEB #T1030). Isolating the desired DNA from unwanted parent vectors and/or other DNA fragments can dramatically improve your cloning results.

Confirm digested DNA on an agarose gel prior to ligation. For a single product, run a small amount of the digest, and then column purify to capture the remainder (e.g., Monarch PCR & DNA Cleanup Kit, NEB #T1030). When multiple fragments are produced and only one is to be used, resolve the fragments on a gel and excise the desired fragment under UV light. Using longwave (365 nm) UV light will minimize any radiation-induced DNA damage to the fragment. Recover the DNA fragment from the agarose slice using a gel extraction kit (e.g., Monarch DNA Gel Extraction Kit, NEB #T1020) or β -Agarase I (NEB #M0392).

6. Quantitate your isolated material

Simple quantitation methods, such as gel electrophoresis with mass standards or spectroscopic quantitation on low-input spectrophotometers (such as a NanoSpec®), ensure that the proper amount of material is used for the downstream joining reaction.

7. Follow the manufacturer's guidelines for the joining reaction

For traditional cloning, follow the guidelines specified by the ligase supplier. If a 3:1 molar ratio of insert to vector is recommended, try this first for best results. Using a 3:1 mass ratio is not the same thing (unless the insert and vector have the same mass). Ligation usually proceeds quickly and, unless your cloning project requires the generation of a high-complexity library that benefits from the absolute capture of every possible ligation product, long incubation times are not necessary.

Follow the manufacturer's guidelines for the joining reactions in PCR cloning and seamless cloning. If you are performing a cloning protocol for the first time, adhere to the recommended protocol for optimal results.

8. Use competent cells that are suited to your needs

While some labs prepare their own competent cells "from scratch" for transformations, the levels of competence achieved rarely matches the high levels attained with commercially-available competent cells. Commercially-available competent cells save time and resources, and make cloning more reproducible.

NANOSPEC® is a registered trademark of Nanometrics, Inc.



APPENDIX

Traditional Cloning Quick Guide

PREPARATION OF INSERT AND VECTORS

Insert From a Plasmid Source

 Digest plasmid with the appropriate restriction enzymes to produce a DNA fragment that can be cloned directly into a vector. Unidirectional cloning is achieved with restriction enzymes that produce non-compatible ends.

Insert From a PCR Product

- Design primers with appropriate restriction sites to clone unidirectionally into a vector
- Addition of 6 bases upstream of the restriction site is sufficient for digestion with most enzymes
- If fidelity is a concern, choose a proofreading polymerase such as Q5 High-Fidelity DNA Polymerase (NEB #M0491)
- Visit www.NEBPCRPolymerases.com for additional guidelines for PCR optimization
- Purify PCR product by running the DNA on an agarose gel and excising the band or by using a spin column (e.g., Monarch® DNA Gel Extraction Kit, NEB #T1020, Monarch PCR & DNA Cleanup Kit, NEB #T1030)
- · Digest with the appropriate restriction enzyme

Standard Restriction Enzyme Protocol

DNA	1 μg
10X NEBuffer	5 μl (1X)
Restriction Enzyme	10 units is sufficient, generally 1 µl is used
Nuclease-free Water	To 50 µI
Incubation Time	1 hour*
Incubation Temperature	Enzyme dependent

^{*} Can be decreased by using a Time-Saver qualified enzyme

Time-Saver Restriction Enzyme Protocol

DNA	1 µg
10X NEBuffer	5 μl (1X)
Restriction Enzyme	1 µl
Nuclease-free Water	Το 50 μΙ
Incubation Time	5–15 minutes*
Incubation Temperature	Enzyme dependent

^{*} Time-Saver qualified enzymes can also be incubated overnight with no star activity

Insert from Annealed Oligos

- Annealed oligos can be used to introduce a fragment (e.g., promoter, polylinker, etc.)
- Anneal two complementary oligos that leave protruding 5' or 3' overhangs for ligation into a vector cut with appropriate enzymes
- Non-phosphorylated oligos can be phosphorylated using T4 Polynucleotide Kinase (NEB #M0201)

Typical Annealing Reaction

Primer	1 μg
10X T4 Ligase Buffer	5 μl
Nuclease-free Water	To 50 μl
Incubation	85°C for 10 minutes, cool slowly (30-60 min.)

Vector

 Digest vector with appropriate restriction enzymes. Enzymes that leave non-compatible ends are ideal as they prevent vector self-ligation

DEPHOSPHORYLATION

- Dephosphorylation is sometimes necessary to prevent self-ligation. NEB offers four products for dephosphorylation of DNA:
- The Quick Dephosphorylation Kit (NEB #M0508), Shrimp Alkaline Phosphatase (rSAP) (NEB #M0371) and Antarctic Phosphatase (AP) (NEB #M0289) are heat-inactivatable phosphatases. They work in all NEBuffers, but AP requires supplementation with Zn²⁺
- Calf Intestinal Phosphatase (CIP) (NEB #M0290) will function under many different conditions and in most NEBuffers. However, CIP cannot be heat inactivated and requires a purification step (e.g., Monarch PCR & DNA Cleanup Kit, NEB #T1030) before ligation.

Dephosphorylation of 5' ends of DNA Using the Quick Dephosphorylation Kit

DNA	1 pmol of DNA ends
10X CutSmart Buffer	2 µl
Quick CIP	1 μΙ
Nuclease-free Water	То 20 µI
Incubation	37°C for 10 minutes
Heat Inactivation	80°C for 2 minutes

Note: Scale larger reaction volumes proportionally.

BLUNTING

- In some instances, the ends of the insert or vector require blunting
- · PCR with a proofreading polymerase will leave a predominantly blunt end
- T4 DNA Polymerase (NEB #M0203) or Klenow (NEB #M0210) will fill in a 5´ overhang and chew back a 3´ overhang
- The Quick Blunting Kit (NEB #E1201) is optimized to blunt and phosphorylate DNA ends for cloning in less than 30 minutes
- Analyze agarose gels with longwave UV (360 nM) to minimize UV exposure that may cause DNA damage

Blunting with the Quick Blunting Kit

•	
DNA	Up to 5 µg
10X Blunting Buffer	2.5 µl
dNTP Mix (1 mM)	2.5 µl
Blunt Enzyme Mix	1 μΙ
Nuclease-free Water	То 25 µI
Incubation	room temperature; 15 min for RE-digested DNA; 30 min for sheared/nebulized DNA or PCR products*
Heat Inactivation	70°C for 10 minutes

^{*} PCR-generated DNA must be purified before blunting using a purification kit (NEB #T1030), phenol extraction/ethanol precipitation, or gel extraction (NEB #T1020).

PHOSPHORYLATION

- For ligation to occur, at least one of the DNA ends (insert or vector) should contain a 5´ phosphate
- Primers are usually supplied non-phosphorylated; therefore, the PCR product will not contain a 5´ phosphate
- Digestion of DNA with a restriction enzyme will always produce a 5´ phosphate
- A DNA fragment can be phosphorylated by incubation with T4 Polynucleotide Kinase (NEB #M0201)

Phosphorylation With T4 PNK

DNA (20 mer)	1–2 μg
10X T4 PNK Buffer	5 μΙ
10 mM ATP	5 μl (1 mM final conc.)
T4 PNK	1 μl (10 units)
Nuclease-free Water	Το 50 μΙ
Incubation	37°C for 30 minutes

PURIFICATION OF VECTOR AND INSERT

 Purify the vector and insert by either running the DNA on an agarose gel and excising the appropriate bands or by using a spin column, such as Monarch DNA Gel Extraction Kit or PCR & DNA Cleanup Kit (NEB #T1020 or T1030)

- DNA can also be purified using β-Agarase I (NEB #M0392) with low melt agarose, or an appropriate spin column or resin
- Analyze agarose gels with longwave UV (360 nM) to minimize UV exposure that may cause DNA damage

LIGATION OF VECTOR AND INSERT

- Use a molar ratio of 1:3 vector to insert. Use NEBioCalculator to calculate molar ratios.
- If using T4 DNA Ligase (NEB #M0202) or the Quick Ligation Kit (NEB #M2200), thaw and resuspend the Ligase Buffer at room temp. If using Ligase Master Mixes, no thawing is necessary.
- The Quick Ligation Kit (NEB #M2200) is optimized for ligation of both sticky and blunt ends
- Instant Sticky-end Ligase Master Mix (NEB #M0370) is optimized for instant ligation of sticky/cohesive ends
- Blunt/TA Ligase Master Mix (NEB #M0367) is optimized for ligation of blunt or single base overhangs, which are the more challenging type of ends for T4 DNA Ligase
- Following ligation, chill on ice and transform
- DO NOT heat inactivate when using the Quick Ligation Buffer or Ligase Master Mixes, as this will inhibit transformation
- Electroligase (NEB #M0369) is optimized for ligation of both sticky and blunt ends and is compatible with electroporation (i.e., no cleanup step required)
- Improved Golden Gate Assembly can be achieved by selecting high fidelity overhangs [Potapov, V. et al. (2018) ACS Synth. Biol. 7(11), 2665–2674.

Ligation with the Quick Ligation Kit

Vector DNA (3 kb)	50 ng
Insert DNA (1 kb)	To 50 ng
2X Quick Ligation Buffer	10 μΙ
Quick T4 DNA Ligase	1 μΙ
Nuclease-free Water	20 μl (mix well)
Incubation	Room temperature for 5 minutes

Ligation with Instant Sticky-end Ligase Master Mix

Vector DNA (3 kb)	50 ng
Insert DNA (1 kb)	50 ng
Master Mix	5 μl
Nuclease-free Water	To 10 µl
Incubation	None

Ligation with Blunt/TA Ligase Master Mix

Vector DNA (3 kb)	50 ng
Insert DNA (1 kb)	50 ng
Master Mix	5 μl
Nuclease-free Water	To 10 µI
Incubation	Room temperature for 15 minutes

TRANSFORMATION

- To obtain tranformants in 8 hrs., use NEB Turbo Competent E. coli (NEB #C2984)
- If recombination is a concern, then use the recA⁻ strains NEB 5-alpha Competent E. coli (NEB #C2987), NEB-10 beta Competent E. coli (NEB #C3019) or NEB Stable Competent E. coli (NEB #C3040)
- NEB-10 beta Competent E. coli works well for constructs larger than 5 kb
- NEB Stable Competent E. coli (NEB #C3040) can be used for constructs with repetitive sequences such as lentiviral constructs
- If electroporation is required, use NEB 5-alpha (NEB #C2989) or NEB 10-beta (NEB #C3020) Electrocompetent E. coli
- Use pre-warmed selection plates
- · Perform several 10-fold serial dilutions in SOC for plating

Transformation with NEB 5-alpha Competent E. coli

DNA	1–5 μl containing 1 pg – 100 ng of plasmid DNA
Competent E. coli	50 μl
Incubation	On ice for 30 minutes
Heat Shock	Exactly 42°C for exactly 30 seconds
Incubation	On ice for 5 minutes Add 950 µl room temperature SOC 37°C for 60 minutes, with shaking

Troubleshooting Guide for Cloning

We strongly recommend running the following controls during transformations. These controls may assist in identifying which step(s) in the cloning workflow has failed.

- Transform 100 pg 1 ng of uncut vector to check cell viability, calculate transformation efficiency and verify the antibiotic resistance of the plasmid.
- Transform the cut vector to determine the amount of background due to undigested plasmid. The number of colonies in this control should be < 1% of the number of colonies in the uncut plasmid control transformation (from control #1).
- Transform a vector only ligation reaction. The ends of the vector should not be able to re-ligate because either they are incompatible (e.g., digested with two restriction enzymes that do not generate compatible ends) or the 5′ phosphate group has been removed in a dephosphorylation reaction (e.g., blunt ends treated with rSAP). This control transformation should yield the same number of colonies as control #2.
- Digest vector DNA with a single restriction enzyme, re-ligate and transform. The ends of the vector DNA should be compatible and easily joined during the ligation reaction, resulting in approximately the same number of colonies as control #1.

The cloning workflow often benefits from an accurate quantitation of the amount of DNAs that are being worked with. We recommend quantification of DNAs whenever possible.

PROBLEM	CAUSE	SOLUTION
	Cells are not viable	Transform an uncut plasmid (e.g., pUC19) and calculate the transformation efficiency of the competent cells. If the transformation efficiency is low (< 104) re-make the competent cells or consider using commercially available high efficiency competent cells.
	Incorrect antibiotic or antibiotic concentration	Confirm antibiotic and antibiotic concentration
	DNA fragment of interest is toxic to the cells	 Incubate plates at lower temperature (25–30°C). Transformation may need to be carried out using a strain that exerts tighter transcriptional control over the DNA fragment of interest (e.g., NEB 5-alpha F⁻ I^a Competent E. coli (NEB #C2992))
	If using chemically competent cells, the wrong heat-shock protocol was used	• Follow the manufacturer's specific transformation protocol (Note: going above the recommended temperature during the heat shock can result in competent cell death)
	If using electrocompetent cells, PEG is present in the ligation mix	Clean up DNA by drop dialysis prior to transformation with Monarch PCR & DNA Cleanup Kit (NEB #T1030) Try NEB's ElectroLigase (NEB #M0369)
	If using electrocompetent cells, arcing was	Clean up the DNA prior to the ligation step To the queste to get rid of any trapped air hubbles.
	observed or no voltage was registered	Tap the cuvette to get rid of any trapped air bubbles Be sure to follow the manufacturer's specified electroporation parameters
	Construct is too large	Select a competent cell strain that can be transformed efficiently with large DNA constructs (≥ 10 kb, we recommend trying NEB 10-beta Competent <i>E. coli</i> (NEB #C3019))
	J.	For very large constructs (> 10 kb), consider using electroporation
Few or no	Construct may be susceptible to recombination	Select a recA- strain such as NEB 5-alpha (NEB #C2987), NEB 10-beta (NEB #C3019) or NEB Stable (NEB #C3040) Competent E. coli
transformants	The insert comes directly from mammalian or plant DNA and contains methylated cytosines, which are degraded by many <i>E. coli</i> strains	• Use a strain that is deficient in McrA, McrBC and Mrr, such as NEB 10-beta Competent E. coli
	Too much ligation mixture was used	\bullet Use < 5 μl of the ligation reaction for the transformation
		Make sure that at least one fragment being ligated contains a 5´ phosphate moiety
	Inefficient ligation	Vary the molar ratio of vector to insert from 1:1 to 1:10. Use NEBiocalculator to calculate molar ratios
		• Purify the DNA to remove contaminants such as salt and EDTA with Monarch PCR & DNA Cleanup Kit (5 µg) (NEB #T1030)
		ATP will degrade after multiple freeze-thaws; repeat the ligation with fresh buffer
	3	Heat inactivate or remove the phosphatase prior to ligation
		Ligation of single base-pair overhangs (most difficult) may benefit from being carried out with Blunt/TA Master Mix (NEB #M0367), Quick Ligation Kit (NEB #M2200) or concentrated T4 DNA Ligase (NEB #M0202)
		• Test the activity of the ligase by carrying out a ligation control with Lambda-HindIII digested DNA (NEB #N0312)
		• Purify the DNA prior to phosphorylation with Monarch PCR & DNA Cleanup Kit (5 µg) (NEB #T1030). Excess salt, phosphate or ammonium ions may inhibit the kinase.
	Inefficient phosphorylation	• If the ends are blunt or 5´ recessed, heat the substrate/buffer mixture for 10 minutes at 70°C. Rapidly chill on ice before adding the ATP and enzyme, then incubate at 37°C.
		• ATP was not added. Supplement the reaction with 1 mM ATP, as it is required by T4 Polynucleotide Kinase (NEB #M0201)
		Alternatively, use 1X T4 DNA Ligase Buffer (contains 1 mM ATP) instead of the 1X T4 PNK Buffer

Troubleshooting Guide for Cloning (continued)

PROBLEM	CAUSE	SOLUTION		
		Heat inactivate or remove the restriction enzymes prior to blunting		
		• Clean up the PCR fragment prior to blunting with Monarch PCR & DNA Cleanup Kit (NEB #T1030)		
		Sonicated gDNA should be blunted for at least 30 minutes		
		• Do not use > 1 unit of enzyme/µg of DNA		
	Inefficient blunting	• Do not incubate for > 15 minutes		
	monoion blanting	• Do not incubate at temperatures > 12°C (for T4 DNA Polymerase, NEB #M0203) or > 24°C (for Klenow, NEB #M0210)		
		 Make sure to add a sufficient amount of dNTPs to the reaction (33 μM each dNTP for DNA Polymerase I, Large (Klenow) Fragment, NEB #M0210 and 100 μM each dNTP for T4 DNA Polymerase, NEB #M0203). 		
Few or no transformants		 When using Mung Bean Nuclease (NEB #M0250), incubate the reaction at room temperature. Do not use > 1 unit of enzyme/µg DNA or incubate the reaction > 30 minutes. 		
	Inefficient A-Tailing	Clean up the PCR prior to A-tailing. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030). High-fidelity enzymes will remove any non-templated nucleotides.		
		Check the methylation sensitivity of the enzyme(s) to determine if the enzyme is blocked by methylation of the recognition sequence		
	Restriction enzyme(s) didn't	Use the recommended buffer supplied with the restriction enzyme		
	cleave completely	Clean up the DNA to remove any contaminants that may inhibit the enzyme. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030).		
		When digesting a PCR fragment, make sure to have at least 6 nucleotides between the recognition site and the end of the DNA molecule		
	Antibiotic level used was too low	Increase the antibiotic level on plates to the recommended amount		
Colonies don't contain a plasmid	Alltiblotic level used was too low	Use fresh plates with fresh antibiotics		
oontain a piasinia	Satellite colonies were selected	Choose large, well-established colonies for analysis		
	Recombination of the plasmid has occurred	• Use a recA ⁻ strain such NEB 5-alpha, NEB 10-beta or NEB Stable Competent E. coli		
	Incorrect PCR amplicon was used	Optimize the PCR conditions		
	during cloning	Gel purify the correct PCR fragment. NEB recommends the Monarch DNA Gel Extraction Kit (NEB #T1020).		
Colonies contain the	Internal recognition site was present	Use NEBcutter to analyze insert sequence for presence of an internal recognition site		
wrong construct		• Incubate plates at lower temperature (25–30°C)		
	DNA fragment of interest is toxic to the cells	• Transformation may need to be carried out using a strain that exerts tighter transcriptional control of the DNA fragment of interest (e.g., NEB 5-alpha F´ Iq Competent E. coll)		
	Mutations are present in the sequence	Use a high-fidelity polymerase (e.g., Q5 High-Fidelity DNA Polymerase, NEB #M0491)		
	initiations are present in the sequence	Re-run sequencing reactions		
	Inefficient dephosphorylation	Heat inactivate or remove the restriction enzymes prior to dephosphorylation		
	Kinase is present/active	 Heat inactivate the kinase after the phosphorylation step. Active kinase will re-phosphorylate the dephosphorylated vector. 		
Too much background		Check the methylation sensitivity of the restriction enzyme(s) to be sure it is not inhibited by methylation of the recognition sequence		
	Restriction enzyme(s) didn't cleave completely	Use the recommended buffer supplied with the restriction enzyme		
	ologyo completely	Clean up the DNA to remove contaminants. (e.g., too much salt). NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030).		
	Antibiotic level is too low	Confirm the correct antibiotic concentration		
		Make sure at least one DNA fragment being ligated contains a 5´ phosphate		
	Inefficient ligation	Vary the molar ratios of vector to insert from 1:1 to 1:10. Use NEBioCalculator to calculate molar ratios.		
Ran the ligation on		Purify the DNA to remove contaminants such as salt and EDTA. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030).		
a gel and saw no		ATP will degrade after multiple freeze-thaws; repeat the ligation with fresh buffer		
ligated product		Heat inactivate or remove the phosphatase prior to ligation		
		Ligation of single base-pair overhangs (most difficult) may benefit from being carried out with Blunt/TA Master Mix, Duick Ligation Vit or consentrated TA DNA Ligace.		
		Quick Ligation Kit or concentrated T4 DNA Ligase • Test the activity of the ligase by carrying out a ligation control with Lambda-HindlII digested DNA		
The ligated DNA ran		- rest the activity of the flydse by carrying out a flydhoff conflict with Edificial-Afficial digested DIVA		
as a smear on an agarose gel	The ligase is bound to the substrate DNA	• Treat the ligation reaction with Proteinase K (NEB #P8107) prior to running on a gel		
	The restriction enzyme(s) is bound to the	Lower the number of units		
The digested DNA ran as a smear on an agarose gel	substrate DNA	Add SDS (0.1–0.5%) to the loading buffer to dissociate the enzyme from the DNA		
		Use fresh, clean running buffer		
	Nuclease contamination	Use a fresh agarose gel Class on the DNA NED consequents the Manager POR & DNA Classics (%) (NED (TAGGO))		
		Clean up the DNA. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030).		

PROBLEM	CAUSE	SOLUTION
		DNA isolated from a bacterial source may be blocked by Dam and Dcm methylation
		DNA isolated from eukaryotic source may be blocked by CpG methylation
	Cleavage is blocked by methylation	Check the methylation sensitivity of the enzyme(s) to determine if the enzyme is blocked by methylation of the recognition sequence
		• If the enzyme is inhibited by Dam or Dcm methylation, grow the plasmid in a dam-/dcm- strain (NEB #C2925)
		Enzymes that have low activity in salt-containing buffers (NEBuffer 3.1) may be salt sensitive, so clean up the DNA prior to digestion. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030).
	Salt inhibition	DNA purification procedures that use spin columns can result in high salt levels, which inhibit enzyme activity. Monarch kits (NEB #T1010, #T1020, #T1030) use columns that have been designed to minimize salt carry over into the eluted DNA, so using them can minimize this issue. To prevent this, DNA solution should be no more than 25% of total reaction volume.
	Inhibition by PCR components	• Clean up the PCR fragment prior to restriction digest. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030).
Incomplete restriction	Using the wrong buffer	Use the recommended buffer supplied with the restriction enzyme
enzyme digestion	Too few units of enzyme used	• Use at least 3–5 units of enzyme per µg of DNA
	Incubation time was too short	Increase the incubation time
	Digesting supercoiled DNA	Some enzymes have a lower activity on supercolled DNA. Increase the number of enzyme units in the reaction.
	Presence of slow sites	Some enzymes can exhibit slower cleavage towards specific sites. Increase the incubation time, 1-2 hours is typically sufficient.
	Two sites required	Some enzymes require the presence of two recognition sites to cut efficiently. For more information, visit the table "Restriction Enzymes Requiring Multi-sites" on neb.com.
	DNA is contaminated with an inhibitor	Assay substrate DNA in the presence of a control DNA. Control DNA will not cleave if there is an inhibitor present. Miniprep DNA is particularly susceptible to contaminants. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030).
		Clean DNA with a spin column, with Monarch PCR & DNA Cleanup Kit (NEB #T1030), resin or drop dialysis, or increase volume to dilute contaminant
	If larger bands than expected are seen in	Lower the number of units in the reaction
	the gel, this may indicate binding of the enzyme(s) to the substrate	Add SDS (0.1–0.5%) to the loading buffer to dissociate the enzyme from the substrate
		Use the recommended buffer supplied with the restriction enzyme
		Decrease the number of enzyme units in the reaction
	Star activity	Make sure the amount of enzyme added does not exceed 10% of the total reaction volume. This ensures that the total glycerol concentration does not exceed 5% v/v
		Decrease the incubation time. Using the minimum reaction time required for complete digestion will help prevent star activity.
Extra bands in the gel		• Try using a High-Fidelity (HF) restriction enzyme. HF enzymes have been engineered for reduced star activity.
Extra bands in the ger	Partial restriction enzyme digest	Enzymes that have low activity in salt-containing buffers (e.g., NEBuffer 3.1) may be salt sensitive. Make sure to clean up the DNA prior to digestion. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030).
		DNA purification procedures that use spin columns can result in high salt levels, which inhibit enzyme activity. Monarch kits (NEB #T1010, #T1020, #T1030) use columns that have been designed to minimize salt carry over into the eluted DNA, so using them can minimize this issue. To prevent this, DNA solution should be no more than 25% of total reaction volume
		• Clean-up the PCR fragment prior to restriction digest. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030).
		Use the recommended buffer supplied with the restriction enzyme
		• Use at least 3–5 units of enzyme per µg of DNA
		Digest the DNA for 1–2 hours
	Used the wrong primer sequence	Double check the primer sequence
	Incorrect annealing temperature	Use the NEB Tm calculator to determine the correct annealing temperature (www.neb.com/TmCalculator)
No PCR fragment	Incorrect extension temperature	• Each polymerase type has a different extension temperature requirement. Follow the manufacturer's recommendations.
amplified	Too few units of polymerase	Use the recommended number of polymerase units based on the reaction volume Each polymerase has a different primer concentration requirement. Make sure to follow the manufacturar's recommendations.
	Incorrect primer concentration Mg ²⁺ levels in the reaction are not optimal	 Each polymerase has a different primer concentration requirement. Make sure to follow the manufacturer's recommendations. Titrate the Mg²⁺ levels to optimize the amplification reaction. Follow the manufacturer's recommendations.
	Difficult template	With difficult templates, try different polymerases and/or buffer combinations
The PCR reaction is a smear on a gel and sare larger than expected it may indicate binding of the enzyme(s) to the DNA * Add SDS (0.1–0.5%) to the loading buffer to dissociate the enzyme from the DNA		
•	Annealing temperature is too low	Use the NEB Tm calculator to determine the annealing temperature of the primers
	Mg ²⁺ levels in the reaction are not optimal	• Titrate the Mg ²⁺ levels to optimize the amplification reaction. Make sure to follow the manufacturer's recommendations.
Extra bands in PCR reaction	Additional priming sites are present	Double check the primer sequence and confirm it does not bind elsewhere in the DNA template
. Off fouction	Formation of primer dimers	Primer sequence may not be optimal. Additional primers may need to be tested in the reaction.
	Incorrect polymerase choice	• Try different polymerases and/or buffer combinations

Optimization Tips for Your Cloning Reactions

New England Biolabs offers a wide selection of reagents for your cloning experiments. For more information, visit **ClonewithNEB.com**. The following tips can be used to help optimize each step in your cloning workflow. Tips for restriction enzyme digestion and amplification can be found on pages 290 and 337, respectively.

CDNA SYNTHESIS

Starting Material

- Intact RNA of high purity is essential for generating cDNA for cloning applications
- Total RNA or mRNA can be used in the reverse transcription reaction. Total RNA is generally sufficient for cDNA synthesis reactions. However, if desired, mRNA can be easily obtained using a PolyA Spin mRNA Isolation Kit (NEB #S1560) or Magnetic mRNA Isolation Kit (NEB #S1550).
- The amount of RNA required for cDNA cloning depends on the abundance of the transcript-of-interest. In general, 1 ng to 1 µg total RNA or 0.1–100 ng mRNA are recommended.

Product Selection

 Streamline your reaction setup by using the ProtoScript II First Strand cDNA Synthesis Kit (NEB #E6560). This kit combines ProtoScript II Reverse Transcriptase (NEB #M0360), a thermostable M-MuLV (RNase H-) Reverse Transcriptase, and recombinant RNase Inhibitor in an enzyme Master Mix, along with a separate Reaction Mix containing dNTPs. Additionally, the kit contains two optimized reverse transcription primer mixes.

Yield

- ProtoScript II Reverse Transcriptase is capable of generating cDNA of more than 10 kb up to 48°C. We recommend 42°C for routine reverse transcription.
- You can increase the yield of a long cDNA product by doubling the amount of enzyme and dNTPs

Additives

 For most RT-PCR reactions, RNase H treatment is not required. But for some difficult amplicons or sensitive assays, add 2 units of E. coli RNase H to the reaction and incubate at 37°C for 20 minutes

PHOSPHORYLATION

Enzyme

- T4 Polynucleotide Kinase (NEB #M0201) and T4 DNA Ligase (NEB #M0202) can be used together in the T4 DNA Ligase Buffer
- T4 Polynucleotide Kinase is inhibited by high levels of salt (50% inhibition by 150 mM NaCl), phosphate (50% inhibition by 7 mM phosphate) and ammonium ions (75% inhibited by 7 mM (NH₄)₂SO₄)
- If using T4 Polynucleotide Kinase and working with 5´-recessed ends, heat the reaction mixture for 10 min at 70°C, chill rapidly on ice before adding the ATP (or Ligase Buffer containing ATP) and enzyme, then incubate at 37°C

Additives

The addition of PEG 8000 (up to 5%) can improve results

DEPHOSPHORYLATION

Enzyme

- When dephosphorylating a fragment following a restriction enzyme digest, a DNA clean up step is required if the restriction enzyme(s) used is NOT heat inactivatable. We recommend the Monarch PCR & DNA Cleanup Kit (NEB #T1030).
- When working with the Quick Dephosphorylation Kit (NEB #M0508), rSAP (NEB #M0371) or AP (NEB #M0289), which are heat-inactivatable enzymes, a DNA clean-up step after dephosphorylation is not necessary prior to the

ligation step. However, when using CIP (NEB #M0290), a clean-up step (e.g., Monarch PCR & DNA Cleanup Kit, NEB #T1030) prior to ligation is necessary.

Additives

 AP requires the presence of Zn²⁺ in the reaction, so don't forget to supplement the reaction with 1X Antarctic Phosphatase Reaction Buffer when using other NEBuffers

BLUNTING/END REPAIR

Enzyme

- Make sure that you choose the correct enzyme to blunt your fragment.
 The Quick Blunting Kit (NEB #E1201), T4 DNA Polymerase (NEB #M0203) and DNA Polymerase I, Large (Klenow) Fragment (NEB #M0210) will fill 5' overhangs and degrade 3' overhangs. Mung Bean Nuclease (NEB #M0250) degrades 5' overhangs.
- T4 DNA Polymerase and DNA Polymerase I, Large (Klenow) Fragment are active in all NEBuffers. Please remember to add dNTPs.

Clean-up

- When trying to blunt a fragment after a restriction enzyme digestion, if the
 restriction enzyme(s) used are heat inactivable, then a clean-up step prior to
 blunting is not needed. Alternatively, if the restriction enzyme(s) used are not
 heat inactivable, a DNA clean-up step is recommended prior to blunting.
- When trying to blunt a fragment amplified by PCR, a DNA clean-up step is necessary prior to the blunting step to remove the nucleotides and polymerase

 When trying to dephosphorylate a fragment after the blunting step, you will need to add a DNA clean-up step (e.g., Monarch PCR & DNA Cleanup Kit, NEB #T1030) after the blunting and before the addition of the phosphatase

Temperature

 When trying to blunt a fragment with Mung Bean Nuclease, the recommended temperature of incubation is room temperature, since higher temperatures may cause sufficient breathing of the dsDNA ends that the enzyme may degrade some of the dsDNA sequence. The number of units to be used and time of incubation may be determined empirically to obtain best results.

Heat Inactivation

Mung Bean Nuclease reactions should not be heat inactivated. Although
Mung Bean Nuclease can be inactivated by heat, this is not recommended
because the DNA begins to "breathe" before the Mung Bean Nuclease is
inactivated and undesirable degradation occurs at breathing sections. Purify
DNA by phenol/chloroform extraction and ethanol precipitation or spin
column purification [e.g., Monarch PCR & DNA Cleanup Kit (NEB #T1030)].

A-TAILING

 If the fragment to be tailed has been amplified with a high-fidelity polymerase, the DNA needs to be purified prior to the tailing reaction. For this we recommend the Monarch PCR & DNA Cleanup Kit (NEB T1030). Otherwise, any high-fidelity polymerase present in the reaction will be able to remove any non-templated nucleotides added to the end of the fragments.

DNA LIGATION

Reaction Buffers

- T4 DNA Ligase Buffer (NEB #B0202) should be thawed on the bench or in the palm of your hand, and not at 37°C (to prevent breakdown of ATP)
- · Once thawed, T4 DNA Ligase Buffer should be placed on ice
- Ligations can be performed in any of the four standard restriction endonuclease NEBuffers or in T4 Polynucleotide Kinase Buffer (NEB #B0201) supplemented with 1 mM ATP
- When supplementing with ATP, use ribo-ATP (NEB #P0756). Deoxyribo-ATP will inhibit ligation.
- Before ligation, completely inactivate restriction enzyme by heat inactivation, spin column (e.g., Monarch PCR & DNA Cleanup Kit, NEB #T1030) or Phenol/EtOH purification

DNA

- Either heat inactivate (AP, SAP, Quick Dephosphorylation Kit) or remove phosphatase (CIP) before ligation
- Keep total DNA concentration between 1–10 μg/ml
- Vector:Insert molar ratios between 1:1 and 1:10 are optimal for single insertions.
 Use NEBioCalculator at NEBioCalculator.neb.com to calculate molar ratios.
- For cloning more than one insert, we recommend the NEBuilder® HiFi DNA Assembly Master Mix (NEB #E2621) or Cloning Kit (NEB #E5520)

 If you are unsure of your DNA concentration, perform multiple ligations with varying ratios

Ligase

- For cohesive-end ligations, standard T4 DNA Ligase. Instant Sticky-end Ligase Master Mix or the Quick Ligation Kit are recommended.
- For blunt and single-base overhangs the Blunt/TA Ligase Master Mix is recommended
- For ligations that are compatible with electroporation, Electroligase is recommended
- Standard T4 DNA Ligase can be heat inactivated at 65°C for 20 minutes
- Do not heat inactivate the Quick Ligation Kit or the ligase master mixes

Transformation

- Add between 1–5 μl of ligation mixture to competent cells for transformation
- · Extended ligation with PEG causes a drop off in transformation efficiency
- Electroporation is recommended for larger constructs (> 10,000 bp). Dialyze samples or use a spin column first if you have used the Quick Ligation Kit or ligase master mixes
- For ligations that are compatible with electroporation, Electroligase is recommended.

TRANSFORMATION

Thawing

- · Cells are best thawed on ice
- DNA should be added as soon as the last trace of ice in the tube disappears
- Cells can be thawed by hand, but warming above 0°C decreases efficiency

DNA

 Up to 10 µl of DNA from a ligation mix can be used with only a 2-fold loss of efficiency

Incubation & Heat Shock

- Incubate on ice for 30 minutes. Expect a 2-fold loss in transformation efficiency for every 10 minutes this step is shortened.
- Both temperature and time are specific to the transformation volume and vessel. Typically, 30 seconds at 42°C is recommended, except when using BL21 (NEB #C2530) which requires exactly 10 seconds.

Outgrowth

- Outgrowth at 37°C for 1 hour is best for cell recovery and for expression
 of antibiotic resistance. Expect a 2-fold loss in transformation efficiency for
 every 15 minutes this step is shortened.
- SOC and NEB 10-beta/Stable Outgrowth Medium give 2-fold higher transformation efficiency than LB medium
- Incubation with shaking or rotation results in 2-fold higher transformation efficiency

Plating

- Selection plates can be used warm or cold, wet or dry with no significant effects on transformation efficiency
- Warm, dry plates are easier to spread and allow for the most rapid colony formation

DNA Contaminants to Avoid

CONTAMINANT	REMOVAL METHOD	
Detergents	Ethanol precipitate	
Phenol	Extract with chloroform and ethanol precipitate	
Ethanol or Isopropanol	Dry pellet before resuspending	
PEG	Column purify (e.g., Monarch PCR & DNA Cleanup Kit) or phenol/chloroform extract and ethanol precipitate	

Troubleshooting Guide for DNA Cleanup & Plasmid Purification using Monarch® Kits

PROBLEM	PRODUCT	POSSIBLE CAUSE	SOLUTION	
No DNA purified	Monarch Plasmid Miniprep Kit	Buffers added incorrectly	Add buffers in the correct order so that the sample is bound, washed and eluted in the correct sequence Ensure ethanol was added to Plasmid Wash Buffer 2	
	(NEB #T1010)	Plasmid loss during culture growth	Ensure proper antibiotic and concentration was used to maintain selection during culture growth	
	Monarch DNA Gel Extraction Kit (NEB #T1020)		5 11 11 11 11 11 11 11 11 11 11 11 11 11	
	Monarch PCR & DNA Cleanup Kit (5 μg) (NEB #T1030)	Ethanol not added to wash buffer	Ensure the proper amount of ethanol was added to Monarch DNA Wash Buffer	
		Incomplete lysis	Pellet must be completely resuspended before addition of Plasmid Lysis Buffer (B2) — color should change from light to dark pink Avoid using too many cells; this can overload the column. If culture volume is larger than recommended, scale up buffers B1-B3.	
		Plasmid loss during culture growth	Ensure proper antibiotic and concentration was used to maintain selection during culture growth	
	Monarch Plasmid Miniprep Kit	Low-copy plasmid selected	Increase amount of cells processed and scale buffers accordingly	
	(NEB #T1010)	Lysis of cells during growth	Harvest culture during transition from logarithmic growth to stationary phase (-12-16 hours)	
		Incomplete neutralization	Invert tube several times until color changes to yellow	
		Incomplete elution	Deliver Elution Buffer directly to center of column Larger elution volumes and longer incubation times can increase yield For elution of plasmids > 10 kb, heat the DNA Elution Buffer to 50°C and extend incubation time to 5 minutes	
Low DNA yield		Buffers added incorrectly	Be sure that buffers have been reconstituted correctly and that reagents have been added in the correct order	
		Gel slice not fully dissolved	Undissolved agarose may clog the column and interfere with binding, Incubate in Monarch Gel Dissolving Buffer for proper time and temperature.	
	Monarch DNA Gel Extraction Kit (NEB #T1020)	Gel dissolved above 60°C	Dissolve gel slice in specified range (37-55°C). Higher temperatures can denature DNA	
	(NEC TITOLO)	Incomplete elution during preparation	Deliver Elution Buffer directly to center of column Larger elution volumes and longer incubation times can increase yield For elution of DNA > 10 kb, heat the Elution Buffer to 50°C and extend incubation time to 5 minutes Multiple rounds of elution can also be performed	
		Buffers added incorrectly	Be sure that buffers have been reconstituted correctly and that reagents have been added in the correct order	
	Monarch PCR & DNA Cleanup Kit (5 μg) (NEB #T1030)	Incomplete elution during preparation	Deliver Elution Buffer directly to center of column Larger elution volumes and longer incubation times can increase yield For elution of DNA > 10 kb, heat the Elution Buffer to 50°C and extend incubation time to 5 minutes Multiple rounds of elution can also be performed	
	Monarch Plasmid Miniprep Kit (NEB #T1010)	Plasmid degradation	Be cautious of strains with high levels of endogenous endonuclease (e.g., HB101 and JM 100 series)	
		Plasmid is denatured	Limit incubation with Plasmid Lysis Buffer (B2) to two minutes, as NaOH in the buffer can denature the plasmid	
Low DNA quality		gDNA contamination	Use careful inversion mixing after cell lysis to avoid shearing of host cell chromosomal DNA. Do not vortex.	
		RNA contamination	Incubate sample in neutralization buffer for the full 2 minutes. For cell culture volumes > 3 ml, increase the spin after neutralization to 5 minutes.	
		Improper storage	Elute DNA in DNA Elution Buffer or nuclease-free water, and store at -20°C. Do not store in solutions containing magnesium.	
		Ethanol has been carried over	Centrifuge final wash for 1 minute to ensure complete removal Ensure column tip does not come in contact with flow through	
	Monarch Plasmid Miniprep Kit	Excessive salt in sample	Use both plasmid wash buffers and do not skip wash steps	
Low DNA performance	(NEB #T1010)	Excessive carbohydrate has been carried over	Avoid strains with high amounts of endogenous carbohydrate (e.g., HB101 and JM 100 series). Be sure to follow protocol and include Plasmid Wash Buffer 1 step.	
		Gel slice not fully dissolved	Undissolved agarose may leach salts into the eluted DNA	
	Monarch DNA Gel Extraction Kit (NEB #T1020)	Ethanol has been carried over	Centrifuge final wash for 1 minute to ensure complete removal Ensure column tip does not come in contact with flow through	
		Trace amounts of salts have been carried over	Ensure column tip does not come in contact with new tube for elution	
	Monarch PCR & DNA Cleanup Kit (5 μg) (NEB #T1030)	Ethanol has been carried over	Centrifuge final wash for 1 minute to ensure complete removal Ensure column tip does not come in contact with flow through	
		Trace amounts of salts have been carried over	Ensure column tip does not come in contact with new tube	

Guidelines for Choosing Sample Input Amounts When Using the Monarch Genomic Purification Kit

Genomic DNA yield, purity and integrity vary immensely based on sample type, input amount and sample condition. Below, we have provided some empirical yield, purity, and DIN data from a wide variety of sample types, as well as guidance on the maximal input amounts for each of those samples when using the Monarch Genomic DNA Purification Kit. It is very important not to overload the column and the buffer system when extracting and purifying gDNA, as DNA yields, purity, integrity, and length may suffer.

SAMPLE TYPE	RECOMMENDED Input amount	TYPICAL YIELD (μg)	DIN	MAXIMUM INPUT AMOUNT
TISSUE*				
Tail (mouse)	10 mg	12–20	8.5-9.5	25 mg
Ear (mouse)	10 mg	18–21	8.5-9.5	10 mg
Liver (mouse and rat)	10 mg	15-30	8.5-9.5	15 mg
Kidney (mouse)	10 mg	10-25	8.5-9.5	10 mg
Spleen (mouse)	10 mg	30-70	8.5-9.5	10 mg
Heart (mouse)	10 mg	9–10	8.5-9.5	25 mg
Lung (mouse)	10 mg	14–20	8.5-9.5	15 mg
Brain (mouse and rat)	10 mg	4–10	8.5-9.5	12 mg
Muscle (mouse and rat)	10 mg	4–7	8.5-9.5	25 mg
Muscle (deer)	10 mg	5	8.5-9.5	25 mg
BL00D**				
Human (whole)	100 μΙ	2.5-4	8.5-9.5	100 μΙ
Mouse	100 μΙ	1–3	8.5–9.5	100 μΙ
Rabbit	100 μΙ	3–4	8.5-9.5	100 μΙ
Pig	100 μΙ	3.5-5	8.5-9.5	100 μΙ
Guinea pig	100 μΙ	3–8	8.5–9.5	100 μΙ
Cow	100 μΙ	2–3	8.5–9.5	100 μΙ
Horse	100 μΙ	4–7	8.5-9.5	100 μΙ
Dog	100 μΙ	2–4	8.5-9.5	100 μΙ
Chicken (nucleated)	10 μΙ	30-45	8.5-9.5	10 µl
CELLS				
HeLa	1 x 10 ⁶ cells	7–9	9.0-9.5	5 x 10 ⁶ cells
HEK293	1 x 10 ⁶ cells	7–9	9.0-9.5	5 x 10 ⁶ cells
NIH3T3	1 x 10 ⁶ cells	6–7.5	9.0-9.5	5 x 10 ⁶ cells
BACTERIA				
E. coli (Gram-negative)	2 x 10° cells	6–10	8.5-9.0	2 x 10° cells
Rhodobacter sp. (Gram-negative)	2 x 10° cells	6–10	8.5-9.0	2 x 10° cells
B. cereus (Gram-positive)	2 x 10° cells	6–9	8.5-9.0	2 x 10° cells
ARCHAEA				
T. kodakarensis	2 x 10° cells	3–5	8.5-9.0	2 x 10° cells
YEAST		<u> </u>		
S. cerevisiae	5 x 10 ⁷ cells	0.5-0.6	8.5-9.0	5 x 10 ⁷ cells
SALIVA/BUCCAL CELLS***				
Saliva (human)	200 µI	2–3	7.0-8.0	500 µІ
Buccal swab (human)	1 swab	5–7	6.0-7.0	1 swab

^{*} Tissue gDNA yields are shown for frozen tissue powder, frozen tissue pieces and RNAlater-stabilized tissue pieces. Though frozen tissue powder results in highly-intact gDNA, lower yields can be expected than when using frozen or RNAlater-stabilized tissue pieces. Residual nuclease activity in tissue pieces will cut the gDNA, resulting in a slightly smaller overall size; however, this gDNA is optimal for silica-based purification.

^{**} Human whole blood samples stabilized with various anticoagulants (e.g., EDTA, citrate and heparin) and various counter-ions were evaluated and results were comparable in all cases. Additionally, all indicated blood samples were tested both as fresh and frozen samples, yielding comparable results. Human samples were donated by healthy individuals; yields from unhealthy donors may differ.

^{***} Buccal swabs and saliva samples partially consist of dead cell material with degraded gDNA. Therefore, the purified gDNA from those samples will naturally have lower DIN values.

Troubleshooting Guide for Genomic DNA Purification using the Monarch Kit

PROBLEM	CAUSE	SOLUTION			
LOW YIELD	UNUUL				
Cells	Frozen cell pellet was thawed and/or resuspended too abruptly	• Thaw cell pellets slowly on ice and flick tube several times to release the pellet from bottom of tube. Use cold PBS, and resuspend gently by pipetting up and down 5–10 times until pellet is dissolved			
CONS	Cell Lysis Buffer was added concurrently with enzymes	Add Proteinase K and RNase A to sample and mix well before adding the Cell Lysis Buffer			
	Blood was thawed, allowing for DNase activity	Keep blood samples frozen and add Proteinase K, RNase A and Blood Lysis Buffer directly to the frozen samples			
Blood	Blood sample is too old	• Fresh (unfrozen) whole blood should not be older than 1 week. Older samples will show progressive DNA degradation and loss of yield.			
	Formation of hemoglobin precipitates	• Species with high hemoglobin content (e.g., guinea pig) may accumulate insoluble hemoglobin complexes that clog the membrane. Reduce Proteinase K lysis time from 5 to 3 minutes.			
	Tissue pieces are too large	• Cut starting material to the smallest possible pieces or grind with liquid nitrogen. In large tissue pieces, nucleases will destroy the DNA before the Proteinase K can lyse the tissue.			
	Membrane is clogged with tissue fibers	 Proteinase K digestion of fibrous tissues (e.g., muscle, heart, skin, ear clips), brain tissue and all RNAlater-stabilized tissues leads to the release of small indigestible protein fibers, which block the binding sites of the silica membrane. To remove fibers, centrifuge lysate at maximum speed for 3 minutes, as indicated in the protocol. For ear clips and brain tissue, use no more than 12–15 mg input material. 			
Tissue	Sample was not stored properly	 Samples stored for long periods of time at room temperature, 4°C or -20°C, will show degradation and loss of gDNA. Flash freeze tissue samples with liquid nitrogen or dry ice and store at -80°C. Alternatively, use stabilizing reagents to protect the gDNA. 			
	Genomic DNA was degraded (common in DNase-rich tissues)	 Organ tissues (e.g., pancreas, intestine, kidney, liver) contain significant amounts of nucleases. Store properly to prevent DNA degradation. Keep on ice during sample preparation. Refer to the protocol for the recommended amount of starting material and Proteinase K to use. 			
	Column is overloaded with DNA	 Some organ tissues (e.g., spleen, kidney, liver) are extremely rich in genomic DNA. Using inputs larger than recommended will result in the formation of tangled, long-fragment gDNA that cannot be eluted from the silica membrane. Reduce the amount of input material. 			
	Incorrect amount of Proteinase K added	• Most samples are digested with 10 µl Proteinase K, but for brain, kidney and ear clips, use 3 µl.			
DNA DEGRADA	ATION	Construction of the Lorentz of the Character of the Construction o			
	Tissue samples were not stored properly	 Samples stored for long periods of time at room temperature, 4°C or -20°C, will show degradation and loss of gDNA. Flash freeze tissue samples with liquid nitrogen or dry ice and store at -80°C. Alternatively, use stabilizing reagents to protect the gDNA. 			
Tissue	Tissue pieces are too large	Cut starting material to the smallest possible pieces or grind with liquid nitrogen. In large tissue pieces, nucleases will degrade the DNA before Proteinase K can lyse the tissue.			
	High DNase content of soft organ tissue	• Organ tissues (e.g., pancreas, intestine, kidney, liver) contain significant amounts of nucleases. Store properly to prevent DNA degradation. Keep on ice during sample preparation. Refer to the protocol for the recommended amount of starting material and Proteinase K to use.			
	Blood sample is too old	• Fresh (unfrozen) whole blood should not be older than 1 week. Older samples will show progressive DNA degradation and loss of yield.			
Blood	Blood was thawed, allowing for DNase activity	Keep frozen blood samples frozen and add enzymes and lysis buffer directly to the frozen samples			
SALT CONTAN	INATION				
	Guanidine Thiocyanate salt from the binding buffer was carried over into the eluate	 When transferring the lysate/binding buffer mix, avoid touching the upper column area with the pipet tip and always pipet carefully onto the silica membrane. Avoid transferring any foam that may have been present in the lysate; foam can enter into the cap area of the spin column. Close the caps gently to avoid splashing the mixture into the upper column area and move the samples with care in and out of the centrifuge. If salt contamination is a concern, invert the columns a few times (or vortex briefly) with gDNA Wash Buffer as indicated in the protocol. 			
PROTEIN CON	TAMINATION	(* * * * * * * * * * * * * * * * * * *			
	Incomplete digestion	 Cut samples to the smallest possible pieces. Incubate sample in the lysis buffer for an extra 30 minutes to 3 hours to degrade any remaining protein complexes. 			
Tissue	Membrane is clogged with tissue fibers	 Proteinase K digestion of fibrous tissues (e.g., muscle, heart, skin, ear clips), brain tissue and all RNAlater-stabilized tissues leads to the release of small, indigestible protein fibers, which block the binding sites of the silica membrane. To remove fibers, centrifuge the lysate at maximum speed for 3 minutes as indicated in the protocol. For ear clips and brain tissue, use no more than 12–15 mg input material. 			
Blood	High hemoglobin content	 Some blood samples (e.g., horse) are rich in hemoglobin, evidenced by their dark red color. Extend lysis time by 3-5 minutes for best purity results. 			
DIUUU	Formation of hemoglobin precipitates	Species with high hemoglobin content (e.g., guinea pig) may accumulate insoluble hemoglobin complexes that clog the membrane. Reduce Proteinase K lysis from 5 to 3 minutes.			
RNA CONTAM	INATION				
Tissue	Too much input material	• DNA-rich tissues (e.g., spleen, liver and kidney) will become very viscous during lysis and may inhibit RNase A activity. Do not use more than the recommended input amount.			
	Lysis time was insufficient	Extend lysis time by 30 minutes to 3 hours after the tissue piece has completely dissolved			
TISSUE DIGES	TION TAKES TOO LONG	Out the use along to the small set associated size as wind with limit at the second size of the state of the second size of the			
	Tissue pieces too large	Cut tissue pieces to the smallest possible size or grind with liquid nitrogen before starting lysis Vertex to release pieces from the table better and immediately offer addition Proteins & Kond Tiggue Lypis Puffer.			
	Tissue pieces are stuck to bottom of tube	Vortex to release pieces from the tube bottom, and immediately after adding Proteinase K and Tissue Lysis Buffer Lieu recommended input amount			
TISSUE LYSAT	Too much starting material E APPEARS TURBID	Use recommended input amount			
AGGGE ETGAL	Formation of indigestible fibers	• Proteinase K digestion of fibrous tissues (e.g., muscle, heart, skin, ear clips), brain tissue and all RNAlater-stabilized tissues leads to the release of small indigestible protein fibers, which block the binding sites of the silica membrane. To remove fibers, centrifuge lysate at maximum speed for 3 minutes, as indicated in the protocol. For ear clips and brain tissue, use no more than 12–15 mg input material.			
RATIO A ₂₆₀ /A ₂₃₀ > 2.5					
	Slight variations in EDTA concentration in eluates	• EDTA in elution buffer may complex with cations like Mg²+ and Ca²+ samples present in genomic DNA, which may lead to higher than usual A₂m/A₂m ratio. In some cases, this ratio exceeds a value of 3.0 and is consistent with highly pure samples. In these cases, the elevated value does not have any negative effect on downstream applications.			

Guidelines for Choosing Sample Input Amounts When Using the Monarch Total RNA Miniprep Kit

RNA yield, purity, and integrity vary immensely based on sample type, input amount and sample condition. Below, we have provided some empirical yield, purity, and RIN data from a wide variety of sample types, as well as guidance on the maximal input amounts for each of those samples when using the Monarch Total RNA Miniprep Kit. It is very important not to overload the column when extracting and purifying RNA, as yields, purity and integrity may suffer.

SAMPLE TYPE(1)		INPUT	AVERAGE YIELD (μg)	OBSERVED RIN	MAXIMUM Starting Material
CULTURED CELI	LS				
HeLa		1 x 10° cells	12–15	9-10	1 x 10 ⁷ cells
HEK 293		1 x 10° cells	12–14	9-10	1 x 10 ⁷ cells
NIH3T3		1 x 10° cells	8–12	9-10	1 x 10 ⁷ cells
MAMMALIAN BL	LOOD (2)				
Human	Fresh	200 μΙ	0.5-1.0	7-8	3 ml
	Frozen	200 μΙ	0.5–1.0	7-8	3 ml
	Stabilized	200 μΙ	0.5-1.0	7-8	3 ml
Rat	Frozen	100 μΙ	5.6	9	1 ml*
BLOOD CELLS					
PBMC (isolated fro	om 5 ml whole blood)	5 ml	3	7	1 x 10 ⁷ cells
TISSUE					
Rat liver	Frozen pulverized	10 mg	25	8–9	20 mg
	Stabilized solid	10 mg	50–60	8–9	20 mg
Rat spleen (stabilized solid with bead homogenizer)		10 mg	40–50	9	20 mg
Rat kidney (frozen pulverized)		10 mg	7–10	9	50 mg
Rat brain	Frozen pulverized	10 mg	2–3	8–9	50 mg
	Stabilized solid	10 mg	0.5–1.5	8–9	50 mg
	Stabilized solid with bead homogenizer	10 mg	5–8	8–9	50 mg
Rat muscle (frozen	pulverized)	10 mg	2–3	8–9	50 mg
Mouse muscle	Frozen pulverized	10 mg	3	8–9	50 mg
	Powder with bead homogenizer	10 mg	5	7–8	50 mg
	Stabilized solid with bead homogenizer	10 mg	8–10	9	50 mg
Mouse heart (stabilized solid w/bead homogenizer)		10 mg	5–6	8–9	50 mg
YEAST					
S. cerevisiae	Frozen with bead homogenizer	1 x 10 ⁷ cells	50	9–10**	5 x 10 ⁷ cells
	Fresh with Zymolyase®	1 x 10 ⁷ cells	60	9**	5 x 10 ⁷ cells
BACTERIA					
E. coli	Frozen	1 x 10° cells	5	10	1 x 10 ⁹ cells
	Frozen with bead homogenizer	1 x 10° cells	10	10	1 x 109 cells
	Frozen with lysozyme	1 x 10° cells	70	10	1 x 109 cells
B. cereus	Frozen with lysozyme	1 x 10 ⁸ cells	20–30	9	1 x 10 ⁹ cells
	Frozen with bead homogenizer	1 x 10 ⁸ cells	8	9–10	1 x 10° cells
PLANT					
Corn leaf (frozen pulverized with bead homogenizer)		100 mg	45	8	100 mg
Tomato leaf (frozen pulverized with bead homogenizer)		100 mg	30	8	100 mg

⁽¹⁾ RNA for other blood samples, including drosophila, zebrafish embryos/larvae, plasma, serum, saliva, buccal swabs and nucleated blood have been successfully purified with this kit; protocols are available in the product manual.

 $^{^{\}mbox{\scriptsize (2)}}$ A protocol for nucleated blood (e.g., birds, reptiles) is also available.

 $^{^{\}star}\,$ Mouse blood also has a maximum input of 1 ml.

^{**} S.cerevisiae total RNA was run on an Agilent® Nano 600 Chip using plant assay.

Troubleshooting Guide for Total RNA Extraction & Purification Using Monarch Kits

PROBLEM	CAUSE	SOLUTION
Clogged column	Insufficient sample disruption or homogenization	Increase time of sample digestion or homogenization Centrifuge sample after Proteinase K digestion or homogenization to pellet debris and use only supernatant for next steps Use larger volume of DNA/RNA Protection Reagent (NEB #T2011) and/or RNA Lysis Buffer (NEB #T2012) for sample disruption and homogenization. See sample-specific protocols in the product manual.
	Too much sample	Reduce amount of starting material to match kit specifications to ensure buffer amounts are sufficient and column is not overloaded. See Guidelines for Choosing Sample Input Amounts on page 357.
	Incomplete elution	After addition of Nuclease-free Water (NEB #B1500) to column matrix, incubate 5-10 min at room temperature and then centrifuge to elute Perform a second elution (note: this will dilute sample)
	Sample is degraded	Store input sample at -80°C prior to use Use Monarch DNA/RNA Protection Reagent (NEB #T2011) to maintain RNA integrity during storage
Low RNA yield	Insufficient disruption or homogenization	Increase time of sample digestion or homogenization Centrifuge sample after Proteinase K digestion or homogenization to pellet debris and use only supernatant for next steps Use larger volume of DNA/RNA Reagent (NEB #T2011) and/or RNA Lysis Buffer (NEB #T2012) for sample disruption and homogenization. See sample specific protocol in the product manual. For Proteinase K treated samples, doubling Proteinase K (from 5% to 10%) may lead to an increase in RNA yield
	Too much sample	Reduce amount of starting material to match kit specifications to ensure buffer amounts are sufficient and column is not overloaded. See Guidelines for Choosing Sample Input Amounts on page 357.
	Starting material not handled/stored properly	 Store input sample at -80°C prior to use. Degradation of RNA may occur if sample is not flash frozen or protected by a preservation reagent. Use Monarch DNA/RNA Protection Reagent (NEB #T2011) to maintain RNA integrity during storage.
RNA degradation	Deviation from the stated protocol may expose RNA to unwanted RNase activities	Refer to the General Guidelines for working with RNA in the product manual
	RNase contamination of eluted materials or kit buffers may have occurred	See General Guidelines for working with RNA in the product manual for advice on reducing risks of contamination
	Low A _{280/280} values indicate residual protein in the purified sample	Ensure the Proteinase K step was utilized for the recommended time. Ensure samples have no debris prior to addition of ethanol and loading onto RNA Purification Column.
Low OD ratios	Low A _{280/230} values indicate residual guanidine salts have been carried over during elution	Ensure wash steps are carried out prior to eluting sample. Use care to ensure the tip of the column does not contact the flow-through after the final wash. If unsure, please repeat centrifugation. When reusing collection tubes, blot rim of tube on a Kimwipe prior to reattachment to the column to remove any residual wash buffer.
	Genomic DNA not removed by column	Perform optional on-column DNase I treatment to remove unwanted gDNA from lysed sample Perform in-tube/off-column DNase I treatment to remove gDNA
DNA contamination	Too much sample	 Reduce amount of starting material to match kit specifications to ensure buffer amounts are sufficient and column is not overloaded. See Guidelines for Choosing Sample Input Amounts on page 357.
Low performance of RNA in downstream steps	Salt and/or ethanol carryover has occurred	Use care to ensure the tip of the RNA Purification Column does not contact the flow-through after the final wash. If unsure, please repeat centrifugation. Be sure to spin the RNA Purification Column for 2 minutes following the final wash with RNA Wash Buffer When reusing collection tubes, blot rim of tube on a Kimwipe prior to reattachment to the column to remove any residual wash buffer Add additional wash step and/or extend spin time for final wash
Unusual spectrophotometric	RNA concentration is too low for spectrophotometric analysis	 For more concentrated RNA, elute with 30 μl of nuclease-free water Increase amount of starting material (within kit specifications). See Guidelines for Choosing Sample Input Amounts on page 357.
readings	Silica fines in eluate	Re-spin eluted samples and pipet aliquot from the top of the liquid to ensure the A _{260/230} is unaffected by possible elution of silica particles

IPPENDIX

Genetic Markers

A *genotype* indicates the genetic state of the DNA in an organism. It is a theoretical construct describing a genetic situation that explains the observed properties (phenotype, see below) of a strain. *E. coli* genotypes list only genes that are defective (1). If a gene is not mentioned, then it is not known to be mutated*. ***. Prophages and plasmids that were present in the original K-12 strain (F, λ , e14, rac) are normally listed only if absent. However, for simplicity, we have not listed λ except when it is present, and we have listed F and its variants in all cases. Parentheses or brackets surround a prophage or plasmid when listed. Genes are given three-letter, lower-case, italicized names (e.g., *dam*) that are intended to be mnemonics suggesting the function of the gene (here, **DN**A adenine methylase). If the same function is affected by several genes, the different genes are distinguished with uppercase italic letters (e.g., *recA*, *recB*, *recC*, *recD* all affect **rec**ombination). Proper notation omits superscript + or – in a genotype, but these are sometimes used redundantly for clarity, as with F'lac-proA*B*. Deletion mutations are noted as Δ , followed by the names of deleted genes in parentheses, [e.g., $\Delta(lac$ -proA*B*. Deletion arabic numerals (e.g., *hsdR17*) and may be characterized as am=amber (UAG) mutation or t= inactive at high temperature, as appropriate. Some common alleles [e.g., $\Delta(lac$ -pro)X111] break the rules. If two strains' genotypes list a gene with the same allele number, they should carry exactly the same mutation.

The *phenotype* of a strain is an observable behavior, e.g., Lac⁻ fails to grow on lactose as a sole carbon source. Phenotypes are capitalized and in Roman type, and the letters are always followed by superscript + or - (or sometimes r, resistant, or s, sensitive). Although phenotypes do not, strictly speaking, belong in a genotype, they are sometimes included following the genotype designation when the former is not obvious from the latter [e.g., rpsL104 (Str')—gene name from ribosomal protein, small subunit, S12, confers resistance to streptomycin].

Some common genes of interest are described below and on the next page; a catalogue of genetically defined genes can be found in reference 2 and on the very useful internet site maintained by the *E. coli* Genetic Stock Center (CGSC) at Yale University http://cgsc.biology.yale.edu/. Additional information from CGSC can be obtained from curator Mary Berlyn by e-mail cgsc@yale.edu.

- * Most E. coli laboratory strains have been heavily mutagenized over forty years of study, and different lines may carry different, so far undiscovered, mutations that may or may not affect your situation. For this reason, it is sometimes useful to try more than one line or strain background in your experiments.
- ** E. coli B and its derivatives are naturally Lon- and Dcm-. We have listed this in brackets even though it is the wild type state for these strains.

dam

Endogenous adenine methylation at
GATC sequences is abolished. dam strains
have a high recombination frequency,
express DNA repair functions constitutively, and are poorly transformed by
Dam-modified plasmids. Used for making
DNA susceptible to cleavage by some
restriction enzymes (e.g., BcII).

dcm Endogenous cytosine methylation at CCWGG sequences is abolished. Used for making DNA susceptible to cleavage by some restriction enzymes (e.g., Avall).

dnaJ

One of several "chaperonins" is inactive.
This defect has been shown to stabilize certain mutant proteins expressed in
E. coli.

dut dUTPase activity is abolished. This mutation, in combination with ung, allows incorporation of uracil into DNA. Used for oligonucleotide mutagenesis.

endA Activity of nonspecific Endonuclease I is abolished. DNA preparations are thought to be of higher quality when prepared from endA strains.

e14 An excisable prophage-like element, present in K-12 but missing from many derivatives. e14 carries the *mcrA* gene among others, therefore e14- strains are McrA-.

F A low-copy number self-transmissible plasmid. F' factors carry portions of the *E. coli* chromosome, most notably the *lac* operon and *proAB* on F' *lac-proA*B**.

An iron uptake receptor is mutated. This mutation confers resistance to phage T1 (ferric hydroxamate uptake). Former name is *tonA*.

gal The ability to metabolize galactose is abolished.

glnV See supE.

fhuA

hsd\$

gyrA A point mutation in DNA gyrase, subunit A. This mutation confers resistance to the antibiotic nalidixic acid.

hflA This mutation results in high frequency lysogenization by λ .

hsdR, DNA that does not contain methylation

of certain sequences is recognized as foreign by EcoKl or EcoBl and restricted (degraded). These enzymes recognize different sequences and are encoded by different alleles of *hsdRMS*. *hsdR* mutations abolish restriction but not protective methylation (r-m¹), while *hsdS* mutations abolish both (r-m-). DNA made in the latter will be restricted when introduced into a wild-type strain.

References

(1) Demerec et al. (1966) Genetics, 54, 61-76.

(2) Berlyn, M.K.B. (1996). In F. C. Niedhardt et al. (Ed.), Escherichia coli and Salmonella: cellular and molecular biology, (2nd ed.), Vol. 2, (pp. 1715–1902). ASM Press.

(3) Raleigh, E.A. et al. (1991) *J. Bacteriol.*, 173, 2707–2709.

laciq The *lac* repressor is overproduced, turning off expression from P*lac* more completely.

lacZ β-galactosidase activity is abolished.

IacZ:: The phage T7 RNA polymerase **T7gene 1** (= gene 1) is inserted into the lacZ gene.

lacY Lactose permease activity is abolished.

 $\Delta(lac)$ = deletion; there are four common deletions involving lac:

 $\Delta(lacZ)$ M15 expresses a fragment that complements the lac α -fragment encoded by many vectors. These vectors will yield blue color on X-Gal only if the host carries Δ M15.

 Δ U169, Δ X111, and Δ X74 all delete the entire lac operon from the chromosome, in addition to varying amounts of flanking DNA. Δ X111 deletes proAB as well, so that the cell requires proline for growth on minimal medium, unless it also carries $F'lac\ proA^*B^*$.

Ion Activity of a protease responsible for degrading aberrant proteins is abolished. Some eukaryotic proteins are stabilized in Ion strains. E. coli B naturally lacks Lon.

Genetic Markers (continued)

IysY The lysozyme gene from the T7 bacteriophage is mutated. The mutation K128Y eliminates lysozyme activity, but the mutant protein still binds to and inhibits T7 RNA polymerase.

malB

The malB region encompasses the genes malEFG and malK lamB malM. Δ(malB) deletes most or all of this region and eliminates expression of Maltose Binding Protein (MalE).

ncrA,
ncrBC

A restriction system that requires methyl cytosine is abolished. DNA containing methylcytosine in some sequences is restricted by Mcr*. dcm-modified DNA is not restricted by Mcr*. Δ(mcrC-mrr) deletes six genes: mcrC-mcrB-hsdS-hsdM-hsdR-mrr, mcrA is lost with e14.

mrr A restriction system that requires cytosine or adenine methylation is abolished; however, dam⁻, dcm⁻ or EcoKI-modified DNA is not restricted by Mrr⁺. The methylcytosine-dependent activity is also known as McrF (3).

mtl The ability to metabolize the sugar alcohol mannitol is abolished.

ompT Activity of outer membrane protease (protease VII) is abolished.

phoA Activity of alkaline phosphatase is abolished.

prc See tsp.

recA Homologous recombination is abolished; particularly desirable when working with sequences containing direct repeats > 50 bp.

recB, Exonuclease and recombination activity of Exonuclease V is abolished. Homologous recombination is much reduced in recB recC strains that are not also sbcB or sbcA. Stability of inverted repeat sequences is enhanced in recB recC strains, especially if they are also sbcB sbcC. Plasmid replication may be aberrant.

recD Exonuclease activity of ExoV is abolished, but recombination activity is elevated. Inverted repeat sequences in λ can be propagated in recD strains. Plasmid replication is aberrant.

recF Plasmid-by-plasmid homologous recombination is abolished.

recJ Plasmid-by-plasmid homologous recombination is abolished.

relA1 Lacks ppGpp synthesis during the stringent response to amino acid starvation; activity of ATP:GTP 3´-pyrophosphotransferase (EC2.7.6.5) is abolished.

rfbD Lacks functional TDP-rhamnose synthetase, and thus does not synthesize the cell surface O-antigen.

rpoH (also known as htpR) Lack of this heat-shock transcription factor abolishes expression of some stress-induced protease activities in addition to lon. Some cloned proteins are more stable in rpoHam supCts strains at high temperature.

sbcB Exo I activity is abolished. Strains carrying recB recC and sbcB are usually also sbcC. These quadruple mutant strains are recombination-proficient and propagate inverted repeats in λ, but plasmid replication is aberrant.

sbcC Usually found with recB recC sbcB. However, strains carrying sbcC alone are recombination-proficient and stably propagate inverted repeats both in λ and in plasmids.

sulA Mutations in this gene allows cells to divide and recover from DNA damage in a lon mutant background (<u>suppressor</u> of Lon).

supC(ts) A thermosensitive tyrosine-inserting ochre (UAA) and amber (UAG) suppressor tRNA. Nonsense mutations in the same strain are suppressed only at low temperatures. Now called tyrT. **supE** A glutamine-inserting amber (UAG) suppressor tRNA; required for growth of some phage vectors. Now called **glnV**.

supF A tyrosine-inserting amber (UAG) suppressor tRNA; required for lytic growth of S7 or S100 λ phage, such as λgt11. Now called tyrT.

thi-1 The ability to synthesize thiamine is abolished (vitamin B1).

traD The self-transmissibility of the F factor is severely reduced.

tsp A periplasmic protease that may degrade secreted or cytoplasmically overexpressed proteins after lysis is abolished. Now called prc.

tsx Confers resistance to bacteriophage T6.

tyrT See supC, supF.

ung Uracil N-glycosylase activity is abolished.
Uracil incorporated into DNA is removed by Ung*, leaving baseless site. See dut.

xyl The ability to metabolize the sugar xylose is abolished.

(P1) The cell carries a P1 prophage. Such strains express the P1 restriction system.

(P2) The cell carries a P2 prophage. This allows selection against Red* Gam*λ (Spi⁻ selection).

(\$0) The cell carries the lambdoid prophage \$0. A defective \$0 prophage carrying the lac M15 deletion is present in some strains.

(Mu) Mu prophage; Mud means the phage is defective.

Enhancing Transformation Efficiency

Transformation efficiency is defined as the number of colony forming units (cfu) that would be produced by transforming 1 µg of plasmid into a given volume of competent cells. However, 1 µg of plasmid is rarely transformed. Instead, efficiency is routinely calculated by transforming 100 pg—1 ng of highly purified supercoiled plasmid under ideal conditions. Transformation Efficiency (TE) is calculated as: TE = Colonies/µg/Dilution. Efficiency calculations can be used to compare cells or ligations. Our recommended protocols and tips are presented here to help you achieve maximum results.

Recommended Protocols

High Efficiency Transformation Protocol

- 1. Thaw cells on ice for 10 minutes
- 2. Add 1 pg-100 ng of plasmid DNA (1-5 μ I) to cells and mix without vortexing
- 3. Place on ice for 30 minutes
- Heat shock at 42°C for 10–30 seconds or according to recommendations. For BL21, use exactly 10 seconds.
- 5. Place on ice for 5 minutes
- Add 950 µl of room temperature SOC or NEB 10-beta/Stable Outgrowth Medium
- 7. Place at 37°C for 60 minutes. Shake vigorously (250 rpm) or rotate.
- 8. Mix cells without vortexing and perform several 10-fold serial dilutions in SOC or NEB 10-beta/Stable Outgrowth Medium.
- Spread 50–100 µI of each dilution onto pre-warmed selection plates and incubate overnight at 37°C (30°C for SHuffle® strains) or according to recommendations

5 Minute Transformation Protocol

(10% efficiency compared to above protocol)

- 1. Thaw cells in your hand
- 2. Add 1 pg-100 ng of plasmid DNA (1-5 µl) to cells and mix without vortexing
- 3. Place on ice for 2 minutes
- 4. Heat shock at 42°C for 30 seconds or according to recommendations.
- 5. Place on ice for 2 minutes
- Add 950 µl of room temperature SOC or NEB 10-beta/Stable Outgrowth Medium. Immediately spread 50–100 µl onto a selection plate and incubate overnight at 37-42°C. (30°C for SHuffle strains) NOTE: Selection using antibiotics other than ampicillin may require some outgrowth prior to plating.

Transformation Tips

Thawing

- Cells are best thawed on ice
- DNA should be added as soon as the last trace of ice in the tube disappears
- Cells can be thawed by hand, but warming above 0°C decreases efficiency

Incubation of DNA with Cells on Ice

 Incubate on ice for 30 minutes. Expect a 2-fold loss in TE for every 10 minutes this step is shortened.

Heat Shock

 Both temperature and time are specific to the transformation volume and vessel. Typically, 30 seconds at 42°C is recommended.

Outgrowth

- Outgrowth at 37°C for 1 hour is best for cell recovery and for expression of antibiotic resistance. Expect a 2-fold loss in TE for every 15 minutes this step is shortened.
- SOC gives 2-fold higher TE than LB medium
- · Incubation with shaking or rotation results in 2-fold higher TE

Plating

- Selection plates can be used warm or cold, wet or dry with no significant effects on TE
- Warm, dry plates are easier to spread and allow for the most rapid colony formation

DNA

- DNA should be purified and resuspended in water or TE Buffer
- Up to 10 µl of DNA from a ligation mix can be used with only a 2-fold loss of efficiency
- Purification by either a spin column or phenol/chloroform extraction and ethanol precipitation is ideal
- The optimal amount of DNA is lower than commonly recognized. Using clean, supercoiled pUC19, the efficiency of transformation is highest in the 100 pg—1 ng range. However, the total colonies which can be obtained from a single transformation reaction increase up to about 100 ng.

DNA Contaminants to Avoid

CONTAMINANT	REMOVAL METHOD
Detergents	Ethanol precipitate
Phenol	Extract with chloroform and ethanol precipitate
Ethanol or Isopropanol	Dry pellet before resuspending
PEG	Column purify or phenol/chloroform extract and ethanol precipitate
DNA binding proteins (e.g., ligase)	Column purify or phenol/ chloroform extract and ethanol precipitate

Electroporation Tips

NEB Turbo (NEB #C2986), NEB 5-alpha (NEB #C2989) and NEB 10-beta (NEB #C3020) Competent *E. coli* Strains are available as electrocompetent cells. The following tips will help maximize transformation efficiencies.

- · Pre-chill electroporation cuvettes and microcentrifuge tubes on ice
- · Thaw cells on ice and suspended well by carefully flicking the tubes
- Once DNA is added, electroporation can be carried out immediately. It is not necessary to incubate DNA with cells. The maximum recommended volume of a DNA solution to be added is 2.5 µl. Addition of a large volume of DNA decreases transformation efficiency.
- DNA should be purified and suspended in water or TE. Transformation
 efficiency is > 10-fold lower for ligation mixtures than the control pUC19
 plasmid due to the presence of ligase and salts. If used directly, ligation
 reactions should be heat-inactivated at 65°C for 20 min and then diluted
 10-fold. For optimal results, spin columns are recommended for clean up of
 ligation reactions.
- Electroporation conditions vary with different cuvettes and electroporators.
 If you are using electroporators not specified in the protocol, you may need to optimize the electroporation conditions. Cuvettes with 1mm gap are recommended (e.g., BTX Model 610/613 and Bio-Rad #165-2089). Higher voltage is required for cuvettes with 2 mm gap.
- Arcing may occur due to high concentration of salts or air bubbles
- It is essential to add recovery medium to the cells immediately after electroporation. One minute delay can cause a 3-fold reduction in efficiency.
- Cold and dry selection plates lead to lower transformation efficiency. Prewarm plates at 37°C for 1 hour. Using 37°C pre-warmed recovery medium increases the efficiency by about 20%.
- Refreeze unused cells in a dry ice/ethanol bath for 5 min and then store at -80°C. Do not use liquid nitrogen. Additional freeze-thaw cycles result in lower transformation efficiency.



Protein Expression with T7 Express Strains

T7 Protein Expression

- Transform expression plasmid into a T7 expression strain. Plate out on antibiotic selection plates and incubate overnight at 37°C (24 hours at 30°C for SHuffle strains).
- 2. Resuspend a single colony in 10 ml liquid culture with antibiotic
- 3. Incubate at 37°C until OD₆₀₀ reaches 0.4–0.6
- Induce with 40 µl of a 100 mM stock of IPTG (final conc. = 0.4 mM) and induce for 2 hours at 37°C (4 hours at 30°C or 16°C overnight for SHuffle strains)
- Check expression by Coomassie stained protein gel, Western Blot or activity assay. Check expression in the total cell extract (soluble + insoluble) and the soluble fraction alone.
- For large scale, inoculate 1 L of liquid medium (with antibiotic) with a freshly grown colony or 10 ml of freshly grown culture. Incubate at 37°C (30°C for SHuffle strains) until OD₆₀₀ reaches 0.4–0.6. Add IPTG to 0.4 mM. Induce 2 hours at 37°C or 15°C overnight (4 hours at 30°C or 16°C overnight for SHuffle strains).

Troubleshooting Tips

No Colonies or No Growth in Liquid Culture

- Even though T7 expression is tightly regulated, there may be a low level of basal expression in the T7 Express host. If toxicity of the expressed protein is likely, transformation of the expression plasmid should be carried out in a more tightly controlled expression strain:
 - In P strains over-expression of the LacP repressor reduces basal expression of the T7 RNA polymerase
 - In IysY strains, mutant T7 lysozyme is produced which binds to T7 RNA polymerase, reducing basal expression of the target protein. Upon induction, newly made T7 RNA polymerase titrates out the lysozyme and results in expression of the target protein.
- Incubation at 30°C or room temperature may also alleviate toxicity issues
- Check antibiotic concentration (test with control plasmid)

No Protein Visible on Gel or No Activity

- Check for toxicity the cells may have eliminated or deleted elements in the expression plasmid. If this is the case, test F and/or IysY strains to reduce basal expression.
- Culture cells for protein induction. Just before induction, plate a sample on duplicate plates with and without antibiotic selection. If toxicity is an issue, significantly fewer colonies will be seen on plates containing antibiotic (indicating that the plasmid has been lost) compared to plates without antibiotic.

Induced Protein is Insoluble

T7 expression often leads to very high production of protein that can result in the target protein becoming insoluble. In this case:

- Induce at lower temperatures (12–15°C overnight)
- Reduce IPTG concentration to 0.01 mM 0.1 mM
- Induce for less time (as little as 15 minutes)
- Induce earlier in growth (OD₆₀₀ = 0.3 or 0.4)





DNA/RNA Input Guidelines for NGS Library Prep

DNA SAMPLE INPUT GUIDELINES

Integrity of DNA

• Start with as high quality DNA as possible. The quality of the input material directly affects the quality of the library. Absorbance measurements can be used as an indication of DNA purity. Ideally, the ratio of the absorbance at 260 nm to 280 nm should be between 1.8–2.0. However, measurements can be affected by the presence of RNA or small nucleic acid fragments. A DNA Integrity Number can be determined using the Agilent TapeStation® and qPCR-based methods can also provide a measurement of DNA integrity.

Quantitation of DNA

 It is important to quantify accurately the DNA sample prior to library construction. Fluorescence-based detection which utilizes dsDNA-specific dyes, such as the Qubit® from Life Technologies, is more accurate than UV spectrometer-based measurements, as the presence of RNA or other contaminants can result in overestimation of the amount of the DNA sample.

RNA SAMPLE INPUT GUIDELINES

Integrity of RNA

- It is important to start with high quality RNA. The use of degraded RNA can result in low yield or failure to generate libraries. We recommend determining RNA quality using the RNA Integrity Number (RIN) estimated by the Agilent® Bioanalyzer® or similar instrumentation. Ideally, the RNA sample should have a RIN value higher than 7, enabling use of poly(A) mRNA or rRNA depletion protocols. Degraded RNA with RIN values as low as 1-2 can be used if specific protocols are followed.
- RNA should be completely free of DNA. DNase digestion of the purified RNA with RNase-free DNase is recommended.

Quantitation of RNA

• It is important to quantify accurately the RNA sample prior to library construction. The concentration can be estimated with the Agilent Bioanalyzer on a pico or nano chip. Alternatively, RNA concentration can be determined by measuring the absorbance at 260 nm (A₂₆₀) in a spectrophotometer such as a NanoDrop®. However, free nucleotides or other organic compounds routinely used to extract RNA will also absorb UV light near 260 nm and will result in an over-estimation of the RNA concentration.

BEAD-BASED CLEAN-UPS AND SIZE SELECTION

Integrity of DNA

- · Be careful not to disturb the bead pellet when transferring material
- Be sure to vortex the beads just before use they should be a uniform suspension
- Do not over-dry the beads. This can make resuspension difficult and reduce yield.
- Bead-based clean-ups and size-selection are explained in animations and videos available on our website
- Use a magnet that is strong enough to separate the beads completely and quickly

INDICES

- Open only one index primer vial at a time, to minimize the risk of contamination
- When you are using a subset of the indices supplied in a kit, or using indices from more than one kit, it is important to optimize the combination of indices used, to ensure balanced sequencing reads.

We provide recommendations for NEBNext index combinations at **NEBNext.com**.

Labeling with SNAP-tag® Technology-Troubleshooting Guide

APPLICATION	PROBLEM	CAUSE	SOLUTION
	No labeling	Fusion protein not expressed	Verify transfection Check expression of fusion protein via Western blot or SDS-PAGE with Vista Green label
	Weak labeling	Poor expression and/or insufficient exposure of fusion protein to substrate	Increase substrate concentration Increase incubation time
Cellular		Rapid turnover of fusion protein	Analyze samples immediately or fix cells directly after labeling Label at lower temperature (4°C or 16°C)
Labeling	High background	Non-specific binding of substrates	Reduce substrate concentration and/or incubation time Allow final wash step to proceed for up to 2 hours Include fetal calf serum or BSA during labeling
	Signal strongly reduced	Instability of fusion protein	Fix cells Switch tag from N-terminus to C-terminus or vice versa
	after short time	Photobleaching	Add commercially available anti-fade reagent Reduce illumination time and/or intensity
	Precipitation	Insoluble fusion	Test from pH 5.0 to 10.0 Optimize salt concentration [50 to 250 mM] Add 0.05 to 0.1% Tween 20
Labeling in Solution	Weak or no labeling	Exhaustive labeling has not been achieved	Increase incubation time to 2 hrs at 25°C or 24 hrs at 4°C Reduce the volume of protein solution labeled Check expression of fusion protein via SDS-PAGE with Vista Green label
	Loss of activity	Instability of fusion protein	Reduce labeling time Decrease labeling temperature (4°C or 16°C)

Cellular Imaging & Analysis FAQs

- **Q.** How does SNAP-tag® labeling differ from using GFP fusion proteins?
- A. GFP and SNAP-tag are both valuable technologies used to visualize proteins in live cells. GFP is an intrinsically fluorescent protein derived from Aequorea victoria while SNAP-tag is derived from hAGT, a human DNA repair protein. In contrast to GFP, the fluorescence of SNAP-tag fusions can be readily turned on with the addition of a variety of fluorescent probes added directly to the culture media. Substituting different fluorophores or other functionalities (biotin, magnetic beads, blocking agents) requires no new cloning or expression, merely incubation of the appropriate substrate with cells, cell lysates or recombinant proteins.
- Q. What is the difference between SNAP- and CLIP-tag™?
- A. SNAP-tag and CLIP-tag are both derived from O⁶ -alkylguanine-DNA-alkyltransferase (hAGT). SNAP-tag recognizes O⁶-labeled benzylguanine substrates while CLIP-tag recognizes O²-labeled benzylcytosine substrates. Each tag transfers the label from the substrate to itself, resulting in specific covalent labeling. In creating the tags, hAGT has been engineered to no longer interact with DNA, but rather with derivatives of the free benzylguanine or benzylcytosine substrates. The tags exhibit no cross-reactivity with one another, enabling researchers to simultaneously label fusion proteins containing SNAP- and CLIP-tags with different fluorophores in live cells.

- Q. Can I clone my protein as a fusion to the N- or C-terminus of the tags?
- A. Yes. SNAP- and CLIP-tags can be fused to either the N- or C-terminus of a protein of interest. However, to label surface proteins on the outside of cells, the SNAP-tag or CLIP-tag must be cloned so that it is oriented to the extracellular surface of the plasma membrane. In this orientation, the tag is accessible to its fluorophore conjugated substrate.
- Q. Are the substrates toxic to cells?
- A. No toxicity has been noted by proliferation or viability assays when using up to 20 μM substrate for 2 hours. Most of the substrates can be incubated with cells for 24 hours up to a concentration of 20 μM without significant toxicity.
- Q. How stable is the labeled protein in mammalian cells?
- A. The stability of the tagged protein in the cell is dependent upon the stability of protein of interest. Labeled SNAP-tag fusion protein has been detected for up to 2 days in mammalian cells.

- Q. Are SNAP-tag substrates stable to fixation?
- A. Yes. SNAP-tag substrates are derived from organic fluorophores which are stable to fixation. Fluorescently-labeled SNAP-tag fusion proteins do not lose signal intensity in contrast to some GFP spectral variants. After labeling the SNAP-tag fusion proteins, the cells can be fixed with standard fixation methods such as para-formal-dehyde, ethanol, methanol, methanol/acetone etc. without loss of signal.
- **Q.** What conditions are recommended for SNAP-tag labeling in vitro?
- A. The SNAP-tag labeling reaction is tolerant of a wide range of buffers. The requirements of the fusion partner should dictate the buffer selected. The following buffer guidelines are recommended: pH between 5.0 and 10.0, monovalent salts (e.g. sodium chloride) between 50 mM and 250 mM and at least 1 mM DTT. Non-ionic detergents can be added to 0.5% v/v if required, but SDS and other ionic detergents should be avoided entirely because they inhibit the activity of the SNAP-tag. Metal chelating reagents (e.g., EDTA and EGTA) also inhibit SNAP-tag activity and should be avoided.

APPEN

Frequencies of Restriction Sites in Sequenced DNAs

The table below summarizes the frequencies with which restriction enzyme sites occur in eleven commonly used DNA molecules. Detailed restriction maps can be found on subsequent pages. The sites listed in these tables were identified by computer analyses of published sequences. Although we have tried to ensure their accuracy, the sites have not necessarily been confirmed by experimentation. When the same specificity is displayed by several enzymes, the site is listed by

the name of the enzyme that is available from New England Biolabs.

Other enzymes with the same specificity are listed in the table of isoschizomers on page 311–327. Enzymes not available from NEB are listed with an (x). If NEB offers an HF version of that enzyme, it is indicated by a red dot (•). Recognition sequences are written 5′ to 3′.

ENZYME	SITE	ADENO-2	LAMRDA	M13MP18	pBR322	pKLAC2	pMAL-P5X	pSNAP _e	pTXB1	pTYB21	pUC19	T7	
Aarl (x)	CACCTGC	9	12	0	0	0	0	1	0	0	0	5	١
Aatll	GACGTC	3	10	0	1	0	0	5	1	0	1	1	l
Accl	GTMKAC	17	9	1	2	5	1	2	5	3	1	33	١
Acc65I	GGTACC	8	2	1	0	0	0	2	0	1	1	5	l
Acil	CCGC	582	516	42	67	81	81	75	102	102	34	199	١
AcII	AACGTT	3	7	2	4	2	5	3	12	13	2	19	ı
Acul	CTGAAG	23	40	0	2	8	4	5	2	2	2	1	1
Afel	AGCGCT	13	2	1	4	0	2	0	1	1	0	0	ı
AfIII	CTTAAG	4	3	0	0	0	0	0	0	0	0	19	1
AfIIII	ACRYGT	25	20	3	1	4	2	5	3	4	1	23	ı
Agel (•)	ACCGGT	5	13	0	0	1	0	1	1	0	0	2	1
Ahdl	GACNNNNNGTC	9	9	0	1	1	2	1	2	2	1	14	ı
Alel	CACNNNNGTG	10	20	1	0	1	0	1	0	2	0	8	1
Alul	AGCT	158	143	27	17	38	28	27	30	31	16	140	ı
Alwl	GGATC	35	58	3	12	18	12	20	15	17	10	1	١
AlwNI	CAGNNNCTG	25	41	1	1	5	2	4	1	2	1	15	ı
Apal	GGGCCC	12	1	0	0	0	1	1	1	1	0	0	1
ApaLI	GTGCAC	7	4	0	3	3	6	4	4	4	3	1	ı
ApeKI	GCWGC	179	199	10	21	27	25	20	27	26	12	116	١
Apol (•)	RAATTY	29	58	11	0	19	5	3	5	6	1	13	l
Ascl	GGCGCGCC	2	2	0	0	0	0	1	0	0	0	0	١
Asel	ATTAAT	3	17	7	1	5	4	7	10	10	3	12	l
AsiSI	GCGATCGC	1	0	0	0	0	0	0	0	0	0	0	١
Aval	CYCGRG	40	8	2	1	2	1	2	3	1	1	4	l
Avall	GGWCC	73	35	1	8	7	9	6	6	7	2	54	١
AvrII	CCTAGG	2	2	0	0	0	0	1	0	0	0	3	l
BaeGI	GKGCMC	45	10	1	3	8	8	8	6	5	3	16	١
Bael	ACNNNNGTAYC	5	10	3	0	1	0	0	0	1	0	3	l
BamHI (•)	GGATCC	3	5	1	1	1	1	1	1	1	1	0	١
Banl	GGYRCC	57	25	7	9	7	4	8	8	7	4	33	l
Banll	GRGCYC	57	7	2	2	5	2	5	3	4	1	1	١
Bbsl (•)	GAAGAC	27	24	0	3	3	3	2	4	5	0	38	l
BbvCl	CCTCAGC	9	7	2	0	0	0	0	0	0	0	10	١
Bbvl	GCAGC	179	199	10	21	27	25	20	27	26	12	116	l
Bccl	CCATC	62	145	14	9	22	16	8	14	20	3	121	١
BceAl	ACGGC	80	115	7	3	13	11	8	13	12	2	47	l
Bcgl	CGANNNNNNTGC	10	28	0	3	6	4	1	4	6	1	19	١
BciVI	GTATCC	9	26	0	2	3	4	3	4	4	2	23	l
BcII (•)	TGATCA	5	8	0	0	2	2	1	1	2	0	1	١
BcoDI	GTCTC	60	37	5	3	11	8	4	8	9	4	95	l
Bfal	CTAG	54	13	5	5	19	3	14	8	8	4	60	١
BfuAl	ACCTGC	39	41	3	1	5	4	4	2	3	1	18	l
BfuCl (x)	GATC	87	116	6	22	35	23	31	24	27	15	6	١
Bgll	GCCNNNNNGGC	20	29	1	3	3	1	7	2	2	2	2	l
BgIII	AGATCT	11	6	1	0	1	1	2	0	1	0	1	١
Blpl	GCTNAGC	8	6	0	0	0	1	0	1	1	0	20	l
BmgBl	CACGTC	15	17	0	0	1	1	2	0	0	0	8	1
Bmrl	ACTGGG	22	4	1	5	2	5	6	11	11	2	6	I
Bmtl (•)	GCTAGC	4	1	0	1	4	0	1	1	1	0	1	1
Bpml	CTGGAG	32	25	2	4	3	4	1	5	6	1	23	I
Bpu10l	CCTNAGC	23	19	4	1	0	1	2	0	2	0	39	1
	CTTGAG	19	13	4		7			7				I
BpuEl		19	13	0	6 1	3	5	9	2	9	4 1	56	1
Bsal (•)	GGTCTC						2	2				29	I
BsaAl	YACGTR	22	14	5	1	4	0	3	2	4	0	35	1
BsaBI	GATNNNNATC	2	21	2	1	2	2	1	1	2	0	7	l
BsaHl	GRCGYC	44	40	1	6	6	5	8	12	8	3	8	1
BsaJI	CCNNGG	234	105	9	8	18	10	17	15	16	5	85	l
BsaWI	WCCGGW	28	81	6	5	8	7	5	8	6	3	32	

Frequencies of Restriction Sites (continued)

ENZYME	SITE	ADENO-2	LAMBDA	M13MP18	pBR322	pKLAC2	pMAL-P5X	pSNAP _F	pTXB1	pTYB21	pUC19	T7
BsaXI	ACNNNNNCTCC	29	19	4	0	3	1	1	2	3	1	12
BseRI	GAGGAG	63	19	1	0	2	0	3	0	0	0	13
BseYI	GCTGGG	31	32	3	2	3	4	5	4	4	1	29
Bsgl	GTGCAG	34	41	0	1	4	6	3	5	4	0	21
BsiEl	CGRYCG	50	22	3	7	11	8	6	9	7	5	17
BsiHKAI	GWGCWC	38	28	3	8	8	9	9	7	7	5	24
BsiWI (•)	CGTACG	4	1	0	0	0	1	0	1	0	0	0
BsII	CCNNNNNNNGG	216	176	17	20	18	16	26	31	27	6	90
Bsml	GAATGC	10	46	1	1	3	1	5	1	0	0	15
BsmBI	CGTCTC	21	14	1	1	2	2	0	2	2	2	16
BsmFI	GGGAC	59	38	2	4	4	1	5	4	4	0	46
BsoBI	CYCGRG	40	8	2	1	2	1	2	3	1	1	4
Bsp1286I	GDGCHC	105	38	5	10	16	11	15	11	10	5	40
BspCNI	CTCAG	75	80	24	7	10	10	9	20	23	5	142
BspDI	ATCGAT	2	15	2	1	2	0	0	0	0	0	3
BspEl	TCCGGA	8	24	0	1	1	2	0	1	1	0	0
BspHI	TCATGA	3	8	1	4	2	1	2	2	2	3	13
BspMI	ACCTGC	39	41	3	1	5	4	4	2	3	1	18
BspQI	GCTCTTC	7	10	0	1	2	1	3	1	1	1	4
BspUI(x)	GCSGC	232	181	7	21	25	18	27	22	23	7	40
Bsrl	ACTGG	86	110	19	18	23	26	19	32	30	11	118
BsrBI	CCGCTC	28	17	4	2	6	4	6	9	9	3	17
BsrDI	GCAATG	14	44	3	2	7	4	4	4	4	2	18
BsrFI	RCCGGY	40	61	1	7	9	2	6	11	7	1	3
BsrGI (•)	TGTACA	5	5	1	0	1	0	1	1	2	0	13
BssHII	GCGCGC	52	6	0	0	1	2	2	1	1	0	1
BssKI (x)	CCNGG	233	185	11	16	25	27	28	46	42	12	11
BssSI	CACGAG	11	8	0	3	5	3	4	2	4	3	31
BstAPI	GCANNNNNTGC	20	34	0	2	3	2	0	3	2	1	12
BstBI	TTCGAA	1	7	0	0	2	0	0	0	1	0	7
BstEII (•)	GGTNACC	10	13	0	0	1	1	0	1	1	0	1
BstNI	CCWGG	136	71	7	6	15	14	19	19	19	5	2
BstUI	CGCG	303	157	17	23	26	31	19	41	35	10	65
BstXI	CCANNNNNTGG	10	13	0	0	2	4	1	4	3	0	11
BstYI	RGATCY	22	21	2	8	12	9	11	10	12	7	2
BstZ17I	GTATAC	3	3	0	1	3	0	1	1	1	0	8
Bsu36I	CCTNAGG	7	2	1	0	1	1	1	0	0	0	30
Btgl	CCRYGG	82	46	2	2	4	3	6	1	3	0	26
BtgZl	GCGATG	23	45	4	3	4	6	6	4	3	0	24
Btsl	GCAGTG	22	34	1	2	7	5	5	4	4	3	20
BtsCI	GGATG	78	150	4	12	20	17	7	12	12	5	97
Cac8l	GCNNGC	285	238	28	31	33	32	41	49	45	14	104
Clal	ATCGAT	2	15	2	1	2	0	0	0	0	0	3
CspCl	CAANNNNGTGG	6	7	1	0	1	0	1	0	0	0	9
CviAII	CATG	183	181	14	26	38	23	21	23	29	11	148
CviKI-1	RGCY	680	692	103	73	131	86	119	112	116	45	562
CviQI	GTAC	83	113	19	3	15	7	14	6	10	3	168
Ddel	CTNAG GATC	97 87	104 116	30 6	8 22	17 35	11 23	11 31	20 24	26 27	6 15	282 6
Dpnl	GATC	87	116	6	22	35	23	31	24	27	15	
DpnII												6
Drall (•)	TTTAAA CACNNNGTG	12 10	13	5 1	3	5 2	1	6	3 1	3 1	3	9 16
(/	GACNNNNNGTC	6	10 3	1	0 2	5	2	3 2	4	4	2	
Drdl Eael	YGGCCR	70	39	3	6	10	5	15	4	5	3	11 2
	CGGCCG	19	2	0	1	4	1	2	2	2	0	0
Eagl Earl	CTCTTC	29	34	2	2	11	6	4	3	4	3	46
Ecil	GGCGGA	29 29	34	2	4	6	6	8	9	11	3	46
Eco53KI	GAGCTC	29 16	32	1	0	2	1	2	0	1	3 1	0
EcoNI	CCTNNNNNAGG	10	9	0	1	3	0	0	2	1	0	1
Eco0109I	RGGNCCY	44	3	0	4 7	1 7	2	5	1	2 5	1	22
EcoP15I	CAGCAG	50	72	4			10	6	6		3	36
EcoRI (•)	GAATTC	5	5	1	1	1	1	1	1	1	1	0
	GATATC	9	21	0	1	1	1	1	1	1	0	0
EcoRV (•) Esp3I	CGTCTC	21	14	1	1	2	2	0	2	2	2	16

ENZYME	SITE	ADENO-2	LAMBDA	M13MP18	pBR322	pKLAC2	pMAL-P5X	pSNAP _F	pTXB1	pTYB21	pUC19	T7	
Faul	CCCGC	147	90	10	10	14	17	11	28	28	5	24	
Fnu4HI	GCNGC	411	380	17	42	52	43	47	49	49	19	156	
Fokl	GGATG	78	150	4	12	20	17	7	12	12	5	97	
Fsel	GGCCGGCC	3	0	0	0	0	0	0	0	0	0	0	
Fspl	TGCGCA	17	15	1	4	3	2	2	1	1	2	7	
Haell	RGCGCY	76	48	6	11	6	9	3	7	7	3	26	
HaellI	GGCC	216	149	15	22	31	23	36	34	36	11	68	
Hgal	GACGC	87	102	7	11	10	12	7	20	18	4	70	
Hhal	GCGC	375	215	26	31	36	39	27	41	39	17	103	
HinP1I	GCGC	375	215	26	31	36	39	27	41	39	17	103	
Hincll	GTYRAC	25	35	1	2	9	7	4	7	6	1	61	
HindIII (•)	AAGCTT	12	6	1	1	1	1	4	0	1	1	0	
Hinfl	GANTC	72	148	26	10	31	9	11	16	20	6	218	
Hpal	GTTAAC	6	14	0	0	3	1	1	1	2	0	18	
Hpall	CCGG	171	328	18	26	32	25	24	50	40	13	58	
Hphl	GGTGA	99	168	18	12	15	19	14	21	21	7	102	
Hpy99I	CGWCG	61	102	8	9	14	9	13	18	14	5	29	
Hpy166II	GTNNAC	116	125	10	8	29	20	13	27	28	5	199	
Hpy188I	TCNGA	80	170	31	15	24	19	17	19	26	10	153	
Hpy188III	TCNNGA	103	185	28	19	32	22	25	27	29	13	173	
HpyAV	CCTTC	84	106	14	10	24	14	11	16	18	6	110	
HpyCH4III	ACNGT	122	187	31	14	25	20	15	18	17	8	174	
HpyCH4IV	ACGT	83	143	22	10	21	10	19	23	26	5	170	
HpyCH4V	TGCA	207	273	18	21	39	28	30	26	25	13	116	
Kasl	GGCGCC	20	1	1	4	1	1	1	1	1	1	2	
KpnI (•)	GGTACC	8	2	1	0	0	0	2	0	1	1	5	
Mbol	GATC	87	116	6	22	35	23	31	24	27	15	6	
Mboll	GAAGA	113	130	10	11	38	15	14	14	17	8	140	
Mfel (•)	CAATTG	4	8	0	0	2	1	2	1	1	0	8	
Mlul (•)	ACGCGT	5	7	0	0	0	1	2	2	1	0	1	
MluCl	AATT	87	189	62	8	43	22	19	31	44	7	79	
Mlyl	GAGTC	40	61	8	4	17	5	6	11	10	4	115	
Mmel	TCCRAC	25	18	3	4	8	3	5	4	4	2	33	
Mnll	CCTC	397	262	62	26	56	24	41	35	39	13	342	
Mscl	TGGCCA	17	18	1	1	0	1	2	0 32	1 41	0	2	
Msel MsII	TTAA CAYNNNNRTG	115 35	195 62	63	15 7	41 10	24 10	23 6	32 7	9	13 3	207 38	
	CMGCKG	95	75	4	6	14	11	8	10	11	6	35	
MspA1I	CCGG	171	328	18	26	32	25		50	40	13	58	
Mspl Mwol	GCNNNNNNGC	391	347	19	34	33	30	24 42	41	35	13	170	
Nael	GCCGGC	13	1	1	4	3	0	2	5	5	0	0	
Narl	GGCGCC	20	1	1	4	1	1	1	1	1	1	2	
Ncil	CCSGG	97	114	4	10	10	13	9	27	23	7	9	
Ncol (•)	CCATGG	20	4	0	0	10	1	3	0	1	0	1	
Ndel	CATATG	2	7	3	1	1	1	1	1	1	1	7	
NgoMIV	GCCGGC	13	1	1	4	3	0	2	5	5	0	0	
Nhel (•)	GCTAGC	4	1	0	1	4	0	1	1	1	0	1	
Nialli	CATG	183	181	14	26	38	23	21	23	29	11	148	
NIaIV	GGNNCC	178	82	18	24	22	14	20	22	24	11	99	
NmeAIII	GCCGAG	17	8	0	3	2	3	2	2	3	1	14	
Notl (•)	GCGGCCGC	7	0	0	0	1	1	1	1	1	0	0	
Nrul (•)	TCGCGA	5	5	0	1	1	0	1	1	0	0	3	
Nsil (•)	ATGCAT	9	14	0	0	6	0	1	0	0	0	8	
Nspl	RCATGY	41	32	6	4	9	3	5	5	6	3	24	
Nt.BstNBI	GAGTC	40	61	8	4	17	5	6	11	10	4	115	
Nt.CviPII	CCD	4148	4641	570	457	806	570	609	716	743	251	3575	
Pacl	TTAATTAA	1	0	1	0	0	0	1	0	0	0	1	
PaeR7I	CTCGAG	6	1	0	0	1	0	1	1	0	0	0	
Pcil	ACATGT	9	2	3	1	3	1	2	1	1	1	6	
PfIFI	GACNNNGTC	12	2	0	1	1	1	1	1	2	0	1	
PfIMI	CCANNNNTGG	18	14	0	2	3	1	5	2	3	0	8	
Phol (x)	GGCC	216	149	15	22	31	23	36	34	36	11	68	
(**)	GAGTC	40	61	8	4	17	5	6	11	10	4	115	
Plel	(3A(311)	411											
PleI PluTI							1	1	1				
PleI PluTI Pmel	GGCGCC GTTTAAAC	20	1 2	1 0	4 0	1 0				1 1	1 0	2	

Frequencies of Restriction Sites (continued)

ENZYME	SITE	ADENO-2	LAMBDA	M13MP18	pBR322	pKLAC2	pMAL-P5X	pSNAP _F	pTXB1	pTYB21	pUC19	T7
PpuMI	RGGWCCY	23	3	0	2	1	2	1	0	0	0	12
PshAl	GACNNNNGTC	2	7	0	1	1	0	0	1	2	0	6
Psil	TTATAA	4	12	2	0	2	1	1	1	1	0	5
PspGI	CCWGG	136	71	7	6	15	14	19	19	19	5	2
Psp0MI	GGGCCC	12	1	0	0	0	1	1	1	1	0	0
PspXI	VCTCGAGB	3	1	0	0	0	0	1	1	0	0	0
PstI	CTGCAG	30	28	1	1	3	1	4	1	1	1	0
Pvul (•)	CGATCG	7	3	1	1	3	2	1	1	1	2	0
Pvull (•)	CAGCTG	24	15	3	1	3	3	3	3	3	2	3
Rsal	GTAC	83	113	19	3	15	7	14	6	10	3	168
RsrII	CGGWCCG	2	5	0	0	0	1	1	0	0	0	1
Sacl (•)	GAGCTC	16	2	1	0	2	1	2	0	1	1	0
SacII	CCGCGG	33	4	0	0	2	0	1	1	1	0	0
Sall (•)	GTCGAC	3	2	1	1	1	1	1	4	1	1	0
Sapl	GCTCTTC	7	10	0	1	2	1	3	1	1	1	4
Sau3AI	GATC	87	116	6	22	35	23	31	24	27	15	6
Sau96I	GGNCC	164	74	4	15	14	20	21	26	28	6	79
Sbfl (•)	CCTGCAGG	3	5	1	0	1	1	1	0	1	1	0
Scal (•)	AGTACT	5	5	0	1	2	1	2	1	2	1	4
ScrFI	CCNGG	233	185	11	16	25	27	28	46	42	12	11
SexAl	ACCWGGT	9	5	0	0	3	0	0	0	0	0	0
SfaNI	GCATC	85	169	7	22	18	20	17	23	19	8	96
SfcI	CTRYAG	47	38	7	4	9	4	10	6	7	4	48
Sfil	GGCCNNNNNGGCC	3	0	0	0	0	0	1	0	0	0	1
Sfol	GGCGCC	20	1	1	4	1	1	1	1	1	1	2
SgrAl	CRCCGGYG	6	6	0	1	0	0	0	1	0	0	0
Smal	CCCGGG	12	3	1	0	1	0	1	0	0	1	0
SmII	CTYRAG	29	17	4	6	8	5	10	8	9	4	75
SnaBl	TACGTA	0	1	1	0	1	0	1	0	0	0	13
Spel (•)	ACTAGT	3	0	0	0	0	0	1	1	1	0	2
Sphl (•)	GCATGC	8	6	1	1	2	0	2	2	2	1	0
Sspl (•)	AATATT	5	20	6	1	6	2	1	3	5	1	6
Stul	AGGCCT	11	6	0	0	1	0	0	1	1	0	1
Styl (•)	CCWWGG	44	10	0	1	4	1	4	2	4	0	36
StyD4I	CCNGG	233	185	11	16	25	27	28	46	42	12	11
Swal	ATTTAAAT	1	0	1	0	0	0	1	1	1	0	1
Taql	TCGA	50	121	12	7	32	16	15	22	20	4	111
Tatl(x)	WGTACW	19	24	5	2	5	1	8	3	6	2	37
Tfil	GAWTC	32	87	18	6	14	4	5	5	10	2	103
Tsel	GCWGC	179	199	10	21	27	25	20	27	26	12	116
Tsp45I	GTSAC	73	81	9	9	9	7	5	12	11	4	108
TspMI	CCCGGG	12	3	1	0	1	0	1	0	0	1	0
TspRI	CASTG	83	119	9	11	22	14	16	16	14	10	94
Tth1111	GACNNNGTC	12	2	0	1	1	1	1	1	2	0	1
Xbal	TCTAGA	5	1	1	0	1	0	1	1	1	1	3
Xcml	CCANNNNNNNNTGG	14	12	0	0	1	3	0	3	4	0	8
Xhol	CTCGAG	6	1	0	0	1	0	1	1	0	0	0
Xmal	CCCGGG	12	3	1	0	1	0	1	0	0	1	0
XmnI	GAANNNNTTC	5	24	2	2	3	1	3	7	8	1	12
Zral	GACGTC	3	10	0	1	0	0	5	1	0	1	1
	a	•	10			0	,			0		'

APPENDIX

Lambda

48,502 base pairs GenBank Accession #: NC_001416 See page 118 for ordering information.


There are no restriction sites for the following enzymes: AsiSI, Fsel, I-Ceul, I-Scel, Notl, PI-Pspl, PI-Scel, Pacl, Sfil, Spel, Sffl(x), Swal (x) = enzyme not available from NEB

Lambda (λ) is a large, temperate *E. coli* bacteriophage with a linear, largely double-stranded DNA genome (1-5). At each end, the 5´ strand overhangs the 3´ strand by 12 bases. These single-stranded overhangs are complementary and anneal to form a cos site following entry into a host cell. Once annealed, the genome is a circular, completely double-stranded molecule which serves as a template for rolling-circle replication.

Many laboratory strains of lambda are derivatives of the strain λ $c1857ind^1$ Sam^7 , which contains four point mutations relative to the wild type strain. The ind^1 mutation in the c1 gene creates a new HindIII site at 37584 not present in the wild type. All lambda products sold by NEB are λ $c1857ind^1$ Sam^7 .

Numbering of the genome sequence begins at the first (5´-most) base of the left end (bottom of diagram below) and continues rightward from late genes *nu1* and *A* towards the early genes. The map below shows the positions of all known ORFs larger than 200 codons.

Enzymes with unique restriction sites are shown in **bold** type and enzymes with two restriction sites are shown in regular type. Location of sites of all NEB restriction enzymes for select plasmids can be found on the NEB website (choose Tools & Resources > DNA Sequences and Maps tool). Restriction site coordinates refer to the position of the 5´-most base on the top strand in each recognition sequence.



- (1) Echols, H. and Murialdo, H. (1978) *Microbiol. Rev.*, 42, 577–591.
- (2) Szybalski, E.H. and Szybalski, W. (1979) Gene, 7, 217-270
- (3) Daniels, D.L., de Wet, J.R. and Blattner, F.R. (1980) *J. Virol.*, 33, 390–400.
- (4) Sanger, F., Coulson, A.R., Hong, G.F., Hill, D.F. and Petersen, G.B. (1982) *J. Mol. Biol.*, 162, 729–773.
- (5) Daniels, D.L. et al. (1983). In R.W. Hendrix, J.W. Roberts, F.W. Stahl and R.A. Weisberg (Eds.), Lambda II: Appendix, New York: Cold Spring Harbor Press.

LacZ

BamHI

Accl Hincll Sall

LacZ

Accl Hincll Smal EcoRI Pstl HindIII mp9/pUC9

LacZ ECORI Saci Xmai Bamhi Xbai Sali Psti Hindili TCGTATGTTGTGGGAATTGTGGGGATAACAATTTCACACAGGAAACAGCTATGGAATTGGAGCTCGCGGGGATCGTTGAGACTGGGAAGTTGGCACTGGCGTTGTTTACAACGTCGTGAAAACCCTGGGAAAACCCTGGCG MetThrWetIleThrAsnSerSerSerProGlyspProLeuGluSerThrCysSerProAluSerJarayaValValLeuGlnArgArgArgAspTrpGluAsnProGly Smal Xmal BamHI Xbal mp10/pUC12

6164
TCGTATGTTGTGTGGGGATAACAATTTCACACAGGAAACAGCTATGACCATGATTACGCCAAGCTTGGGGTGCAGGTCGACTCTAGAGGATCCCCGGGCGAGCTCGAGCTTGACAGCTTGGGAAAACCCTGGCG
MetThrMetileTGTropserleuglucspprokrgAlaSerSerAsnSerLeudlucspprokrgAlaSerAsnSerLeudlucspprokrgAlaSerAsnSerLeudlucspangasptrogly Sacl EcoRI Smal Xbal BamHl Xmal Accl Hincl I Sall Pst Hind

mp11/pUC13

Smal Kpnl Xmal mp18/pUC18

LacZ 6164
TCGTATGTTGTGTGTGTGAGCGGATAACAATTTCACACAGGAAACAGCTATGACTACGCCAAGGTTGCATGCCTGCAGGTCGACTCTAGAGGATCCCGGGGTACCGGGGTACCGGAATTCACTGGCCGTCGTTTTACAACGTCGTGAATGATGGCTGGAAAACCCTGGCG Sacl EcoRI Smal BamHl Xmal Kpnl Xbal Accl Hincll Sall Pstl HindIII Sphl

370

M13mp18

GenBank Accession #: X02513 Revised sequence file available at www.neb.com. See page 118 for ordering information.

There are no restriction sites for the following enzymes: Aarl(x), Aatll, Acul, Aflll, Agel, Ahdl, Apal, Apall, Ascl, AsiSl, Avrll, Bbsl, Bcgl, BciVl, Bcll, Blpl, BmgBl, Bmtl, Bsal, Bsgl, BsiWl, BspEl, BspQl, BssHll, BssSl, BstAPl, BstBl, BstEll, BstXl, BstZ171, Eagl, EcoNI, EcoO109I, EcoRV, Fsel, FspAl(x), Hpal, I-Ceul, I-Scel, Mfel, Mlul, Ncol, Nhel, NmeAlll, Notl, Nrul, Nsil, PI-Pspl, PI-Scel, PaeR7I, PfiFI, PfiMI, Pmel, PmlI, PpuMI, PshAl, PspOMI, PspXI, RsrIl, SacIl, SanDl(x), Sapl, Scal, SexAl, Sfil, SgrAl, Spel, Srfl(x), Stul, Styl, Tth1111, Xcml, Xhol, Zral

(x) = enzyme not available from NEB

M13 is a filamentous *E. coli* bacteriophage specific for male (F factor-containing) cells. Its genome is a circular, single-stranded DNA molecule 6407 bases in length, and contains 10 genes. A double-stranded form (RF) arises as an intermediate during DNA replication.

The M13mp phage vectors, derived from M13, contain the $lacZ\alpha$ gene and differ from each other by the cloning sites embedded within it. The location of cloning sites inside this gene allows screening for insertions using α -complementation. The map of M13mp18, whose multiple cloning site (MCS) was later employed to construct the plasmid pUC19, is shown below; sequences of the MCS region from other M13mp vectors are shown on the previous page. M13mp19 is identical to M13mp18 except that the MCS region (6231-6288) is inverted.

The complete nucleotide sequences of M13mp18 and M13mp19 have recently been determined at New England Biolabs (1), resulting in several nucleotide changes relative to the previous sequence data (2,3).


Enzymes with unique restriction sites are shown in **bold** type. Location of sites of all NEB restriction enzymes for select plasmids can be found on the NEB website (choose Tools & Resources > DNA Sequences and Maps tool). Restriction site coordinates refer to the position of the 5′-most base on the top strand in each recognition sequence.

Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons.

M13 origin of replication arrows indicate the direction of synthesis of both the (+) and (-) strands.

Feature	Description	Coordinates
gene II	replication	6848-831 (cw)
gene X	replication	496-831
gene V	replication	843-1106
gene VII	minor coat protein	1108-1209
gene IX	minor coat protein	1206-1304
gene VIII	major coat protein	1301-1522
gene III	minor coat protein	1578-2852
gene VI	minor coat protein	2855-3193
gene I	phage assembly	3195-4241
gene XI (I*)	phage assembly	3915-4241
gene IV	phage assembly	4219-5499
ori	M13 origin (+) of replication	5487-5867
$lacZ\alpha$	for $lpha$ -complementation	6216-6722
MCS	multiple cloning site	6230-6286

(cw) = clockwise





- (1) Stewart, F.J. (2002) unpublished observations.
- (2) Messing, J. et al. (1977) Proc. Natl. Acad. Sci. USA, 74, 3652-3646.
- (3) Yanisch-Perron, C., Vieira, J. and Messing, J. (1985) Gene, 33, 103-119.

GenBank Accession #: J01749 See page 118 for ordering information.

There are no restriction sites for the following enzymes: Aarl(x), Acc65I, AfIII, AgeI, AleI, ApaI, ApuI, AscI, AsiSI, AvrII, BaeI, BbvCI, BcII, BgIII, BlpI, BmgBI, BsaXI, BseRI, BsiWI, BsrGI, BssHII, BstBI, BstEII, BstXI, Bsu36I, CspCI, DraIII, Eco53KI, FseI, HpaI, I-CeuI, I-SceI, KpnI, MfeI, MluI, NcoI, NotI, NsiI, PI-PspI, PI-SceI, PacI, PaeR7I, PmeI, PmII, PsiI, PspOMI, PspXI, RsrII, SacI, SacII, SanDI(x), SbfI, SexAI, SfiI, SmaI, SnaBI, SpeI, SrII(x), StuI, SwaI, TspMI, XbaI, XcmI, XhoI, XmaI

(x) = enzyme not available from NEB

pBR322 is an *E. coli* plasmid cloning vector containing the origin of replication from pMB1 (a plasmid in the ColE1 compatibility group; 1–3). The *rop* gene product, which regulates plasmid replication by stabilizing the interaction between RNAI and RNAII transcripts, maintains the copy number at about 20 per cell. However, pBR322 can be amplified with chloramphenicol or spectinomycin (4).


Enzymes with unique restriction sites are shown in **bold** type and enzymes with two restriction sites are shown in regular type. Location of sites of all NEB restriction enzymes for select plasmids can be found on the NEB website (choose Tools & Resources > DNA Sequences and Maps tool). Restriction site coordinates refer to the position of the 5′-most base on the top strand in each recognition sequence.

Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons.

Origin of replication coordinates include the region from the -35 promoter sequence of the RNAII transcript to the RNA/ DNA switch point. *bla* (Apⁿ) gene coordinates include the signal sequence.

Feature	Coordinates	Source
tet (TcR)	86-1276	pSC101
bla (ApR)	4153-3293	Tn3
rop	1915-2106	pMB1
origin	3122-2534	pMB1

ori = origin of replication Ap = ampicillin, Tc = tetracycline



- (1) Bolivar, F. et al. (1977) Gene, 2, 95-113.
- (2) Sutcliffe, J.G. (1979) Cold Spring Harb. Symp. Quant. Biol., 43, 77-90.
- (3) Watson, N. (1988) Gene, 70, 399-403.
- (4) Sambrook, J., Fritsch, E.F., and Maniatis, T. (1989). Molecular Cloning: A Laboratory Manual, (2nd ed.), Cold Spring Harbor, Cold Spring Harbor Laboratory Press.

GenBank Accession #: EU196354 See page 215 for ordering information.

There are no restriction sites for the following enzymes: Aarl(x), Aatll, Acc65l, Afel, Aflll, Apal, Ascl, AsiSI, AvrII, BbvCI, BIpI, Bpu10I, BsiWI, Fsel, FspAl(x), I-Ceul, I-Scel, Kpnl, Mlul, Mscl, Pacl. PI-Pspl. PI-Scel. Pmel. Pmll. PspOMI. PspXI, RsrII, SanDI(x), SfiI, SgrAI, SpeI, SrfI(x), Swal 7ral

(x) = enzyme not available from NEB

pKLAC2 is an expression vector capable both of replication in E. coli and stable integration into the genome of the yeast Kluyveromyces lactis (1). It is designed for high-level expression of recombinant protein in K. lactis using the K. lactis Protein Expression Kit (NEB #E1000). pKLAC2 contains a universal multiple cloning site (MCS) that is compatible with all NEB expression systems.

In E. coli, it replicates using the pMB1 origin of replication from pBR322 (although the rop gene is missing) and carries the bla (Ap^R) marker for selection with ampicillin. Upon transformation of K. lactis GG799 competent cells (NEB #C1001), SacII- or BstXI-linearized pKLAC2 integrates into the K. lactis chromosome at the LAC4 locus. Yeast transformants can be selected using the acetamidase selectable marker (amdS), which is expressed from the yeast ADH1 promoter. Acetamidase expressed from pKLAC2 permits transformed cells to utilize acetamide as a sole nitrogen source on defined medium (2).

The multiple cloning site (MCS) is positioned to allow translational fusion of the K. lactis lpha-mating factor secretion domain $(\alpha$ -MF) to the N-terminus of the recombinant target protein. This directs the fusion protein to the general secretory pathway, but the α -MF domain is cleaved off in the Golgi apparatus by the Kex protease, resulting in secretion of the recombinant protein alone.

Expression of the recombinant fusion protein is driven by the K. lactis LAC4 promoter, which has been modified to be transcriptionally silent in E. coli (1). This facilitates the cloning of proteins that are toxic to E. coli. This promoter is split such that when pKLAC2 is cleaved with SacII or BstXI, the recombinant protein and selectable marker are flanked by the two halves of the promoter. When these ends recombine with the LAC4 promoter in the K. lactis chromosome, the result is integration of the recombinant fusion protein (driven by the LAC4 promoter) and amdS upstream of the LAC4 gene (driven by a duplicate copy of the LAC4 promoter) (2).

Enzymes with unique restriction sites are shown in **bold** type and selected enzymes with two restriction sites are shown in regular type. Location of sites of all NEB restriction enzymes for select plasmids can be found on the NEB website (choose Tools

& Resources > DNA Sequences and Maps tool). Restriction site coordinates refer to the position of the 5´-most base on the top strand in each recognition sequence.


Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons. Components of coordinated regions are indented below the region itself.

pMB1 origin of replication coordinates include the region from the -35 promoter sequence of the RNAII transcript to the RNA/DNA switch point. Promoter and transcription terminator coordinates represent cloned regions and not necessarily the precise functional elements

Feature	Coordinates	Source
expression region:		
α-mating factor		
leader sequence	14-349	K. lactis
MCS	257-354	_
LAC4TT region	371-953	K. lactis
AdH1 promoter region	1010-1712	S. cerevisiae
amdS	1713-3359	A. nidulans
LAC4 promoter		
region (5´ end)	4068-4648	K. lactis
origin	5102-5690	pMB1
bla (ApR)	6721-5861	Tn3
LAC4 promoter		
region (3´ end)	7475-9107	K. lactis (modified)
MCS LAC4TT region AdH1 promoter region amdS LAC4 promoter region (5' end) origin bita (Ap") LAC4 promoter	257-354 371-953 1010-1712 1713-3359 4068-4648 5102-5690 6721-5861	K. lactis S. cerevisiae A. nidulans K. lactis pMB1 Tn3

ori = origin of replication Ap = ampicillin

TT = transcription terminator



BstXI 4319

381-392.

(1) Colussi, P.A. and Taron, C.H. (2005) Appl. Environ. Microbiol., 71, 7092-7098. (2) van Ooyen, A.J. et al. (2006) FEMS Yeast Res., 6,

Sequence file available at www.neb.com. See page 213 for ordering information.

Feature	Coordinates	Source
lacl ^q	81-1163	E. coli
P _{tac}	1406-1433	-
expression ORF	1528-2832	-
malE	1528-2703	E. coli
MCS	2709-2832	-
bla (ApR)	3162-4022	Tn3
origin	4110-4698	pMB1
rop	5068-5259	pMB1

There are no restriction sites for the following enzymes: Aarl(x), Aatll, Acc65I, AfIII, Agel, Alel, AscI, AsiSI, AvrII, Bael, BbvCI, Bmtl, BsaAI, BseRI, BspDI, BsrGI, BstBI, BstZ171, Clal, CspCI, DrallI, EcoNI, Fsel, I-Ceul, I-Scel, Kpnl, Nael, NgoMIV, Nhel, NruI, Nsil, PI-PspI, PI-Scel, PacI, PaeR7I, Pmel, PmII, PshAI, PspXI, SacII, SanDI(x), SexAI, SfiI, SgrAI, SmaI, SnaBI, Spel, SphI, Srfi(x), StuI, SwaI, TspMI, XbaI, XhoI, XmaI, ZraI

(x) = enzyme not available from NEB

pMAL-p5X is an *E. coli* plasmid cloning vector designed for recombinant protein expression and purification using the pMAL Protein Fusion and Purification System (NEB #E8200) (1–3). It contains the pMB1 origin of replication from pBR322 and is maintained at a similar copy number to pBR322.

The multiple cloning site (MCS) is positioned to allow translational fusion of the *E. coli* maltose binding protein (MBP, encoded by the *malE* gene) to the N-terminus of the cloned target protein. The pMAL-p5 and -c5 series of vectors differs from the -p4 and -c4 series in that they contain a universal multiple cloning site (MCS) that is compatible with other NEB expression systems and is followed by stop codons in all three reading frames. In addition, $lacZ\alpha$ and the M13 origin have been removed. In these vectors, MBP has been engineered for tighter binding to amylose. This allows easy purification of the fusion protein, and the MBP domain can be subsequently removed using Factor Xa protease (3).


Transcription of the gene fusion is controlled by the inducible "tac" promoter ($P_{\rm tac}$). Basal expression from $P_{\rm tac}$ is minimized by the binding of the Lac repressor, encoded by the lacla gene, to the lac operator immediately downstream of $P_{\rm tac}$. A portion of the rnB operon containing two terminators, derived from the vector pKK233-2, prevents transcription originating from $P_{\rm tac}$ from interfering with plasmid functions.

pMAL-c5-series vectors are identical to the pMAL-p5-series vectors above except for a deletion of the *malE* signal sequence (nt 1531-1605) (1).

Enzymes with unique restriction sites are shown in **bold** type, and enzymes with two restriction sites are shown in regular type. Location of sites of all NEB restriction enzymes for select plasmids can be found on the NEB website (choose Tools & Resources > DNA Sequences and Maps tool). Restriction site coordinates refer to the position of the 5′-most base on the top strand in each recognition sequence.

Open reading frame (ORF) coordinates are in the form "translational start — translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons.

The pMB1 origin of replication includes the region from the -35 promoter sequence of the RNAII transcript to the RNA/DNA switch point (labeled "ori") and the rop gene, which controls expression of the RNAII transcript. b/a (Ap^R) gene coordinates include the signal sequence.



- (1) Guan, C. et al. (1987) Gene, 67, 21-30.
- (2) Maina, C.V. et al. (1988) Gene, 74, 365-373.
- (3) Riggs, P.D. (1992). In F.M. Ausubel, et al. (Eds.), Current Prot. in Molecular Biol. New York: John Wiley & Sons, Inc.

pMiniT 2.0

Sequence available at www.neb.com See page 91 for more information.

Feature	Coordinates	Source
Constitutive promoter	1-214	pNK2138
SP6 promoter	479-496	SP6
Toxic minigene	541-549	_
Synthetic T7 promoter	619-602	T7
bla (ApR)	733-1593	Tn3
origin	1764-2352	pUC19

There are no restriction sites for the following enzymes: Absl(x), Acc65l, Accl, AfIII, Agel, Ajul(x), Alel, Alol(x), Apal, Arsl(x), Ascl, AsiSI, AvrII, Bael, BanII, Barl(x), Bbsl, BbvCI, BcII, BgIII, Blpl(x), BmgBl, Bmtl, Bpll(x), Bpu10l, Bsal, BsaAl, BsaBI, BseRI, BsgI, BsiWI, BsmFI, BsmI, BspDI, BspEI, BsrGI, BssHII, BstAPI, BstBI, BstEII, BstXI, BstZ17I, Bsu36I, Btgl, Clal, CspCI, DrallI, Eco53kl, EcoNI, EcoO109I, EcoRV, Fall(x), Fsel, FspAl(x), Hincll, Hindlll, Hpal, Kasl, Kfll(x), Kpnl, MauBl(x), Mfel, Mlul, Mrel(x), Mscl, Mtel(x), Nael, Narl, Ncol, NgoMIV, Nhel, Nsil, Pasl(x), PfIFI, PfIMI, Pfol(x), PluTI, PmII, PpuMI, PshAI, PsiI, PspOMI, Psrl(x), Pvull, Rsrll, Sacl, Sacll, Sall, SexAl, Sfil, Sfol, SgrAl, SgrDl(x), Smal, SnaBl, Spel, Sphl, Srfl(x), Stul, Styl, Swal, TspMI, Tth111I, Xbal, Xcml, Xmal


(x) = enzyme not available from NEB

pMiniT 2.0 is an *E. coli* plasmid cloning vector designed for cloning blunt-ended or single-base overhang PCR products, or amplicons, using the NEB PCR Cloning Kit (NEB #E1202, #E1203). The pMiniT2.0 also enables *in vitro* transcription using SP6 and T7 promoters. It is compatible with Golden Gate Assembly as the Bsal site has been removed from the Ampicillin resistance gene.


In *E. coli*, it replicates using the pMB1 origin of replication from pUC19 and carries the *bla* (ApR) marker for selection with ampicillin. pMiniT2.0 contains a toxic minigene that is under the control of a constitutive promoter. If the pMiniT 2.0 vector recircularizes without

an insert, the toxic minigene it will cause lethal inhibition of protein synthesis and no colony will result. If the pMiniT 2.0 Vector carries an insert, a colony will grow.

The map shown below displays the construct formed if no insert is present. Unique restriction sites are shown in **bold**. Additional restriction sites that can be used for subcloning are also shown. Expanded box below shows location of sequencing primers, restriction sites for subcloning or linearization for *in vitro* transcription, RNA Polymerase promoter sequences and placement of insertion site within the toxic minigene.



Features within Sequence Flanking the Toxic Minigene/Cloning Site:



Sequence file available at www.neb.com.

There are no restriction sites for the following enzymes: Aarl(x), Acc65I, AccI, AfeI, AfIII, AgeI, AleI, ApaI, AsiSI, AvaI, AvrII, BaeI, BbsI, BcII, BfuAI, BgIII, BIpI, BmgBI, BmtI, BsaAI, BsaBI, BsgI, BsiWI, BsmFI, BsmI, BsoBI, BspDI, BspEI, BspMI, BsrGI, BstBI, BstEII, BstXI, BstZ17I, Bsu36I, BtgI, BtgZI, ClaI, CspCI, DralII, EagI, EcoNI, EcoRV, FseI, FspAI(x), HincII, HpaI, I-CeuI, I-SceI, KpnI, MfeI, MluI, MscI, NaeI, NcoI, NgoMIV, NheI, NotI, NruI, NsiI, PI-PspI, PI-SceI, PaeR7I, PfIFI, PfIMI, PmII, PpuMI, PshAI, PsiI, PspOMI, PspXI, RsrII, SacII, SaII, SanDI(x), SexAI, SfiI, SgrAI, SmaI, SnaBI, SpeI, SphI, SrfI(x), StuI, StyI, SwaI, TspMI, Tth111I, XcmI, XhoI, XmaI

(x) = enzyme not available from NEB

pNEB206A is an *E. coli* plasmid vector designed for fast and efficient cloning of PCR products to be used in conjunction with USER Enzyme (NEB #M5505; 1). It is derived from pNEB193 containing the high-copy pUC19 origin of replication and $lacZ\alpha$ gene for screening of insertions at the cloning site using α -complementation (2).

The plasmid is supplied in a linearized form 2,706 bp in length (with bp 438-453 excised from the circular form), flanked by two noncomplementary 8-base 3′ overhangs at the intended cloning site. Amplification with deoxyuridine-containing primers and subsequent treatment (as defined in the protocol "Cloning with USER Enzyme" found on our website), results in PCR products with 5′ overhangs complementary to those in pNEB206A. These products can be directionally cloned into pNEB206A at high efficiency without the use of restriction enzymes or DNA ligase, forming recombinant circular molecules.

Enzymes with unique restriction sites are shown in **bold** type, and enzymes with two restriction sites are shown in regular type. Location of sites of all NEB restriction enzymes for select plasmids can be found on the NEB website (choose Tools &


Resources > DNA Sequences and Maps tool). Restriction site coordinates refer to the position of the 5′-most base on the top strand in each recognition sequence. Coordinates on the map and in the tables refer to the 2,722 bp circular plasmid prior to linearization and can be used to calculate relative distances.

Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons.

Origin of replication coordinates include the region from the -35 promoter sequence of the RNAII transcript to the RNA/ DNA switch point. *bla* (Apⁿ) gene coordinates include the signal sequence. Cloning site coordinates include those bases in the circular form that are single-stranded in or missing from the supplied linear form.

Feature	Coordinates	Source
lacZa.	505-146	-
cloning site	430-461	-
origin	1491-903	pUC19
bla (ApR)	2522-1662	Tn3

ori = origin of replication Ap = ampicillin



pNEB206A (linearized form) cloning site:



- (1) Bitinaite, J. and Vaiskunaite, R. (2003) unpublished observations.
- (2) Yanisch-Perron, C. et al. (1985) Gene, 33, 103-119.

pSNAP_f

Sequence file available at www.neb.com. See page 282 for ordering information.

There are no restriction sites for the following enzymes: Absl(x), Afel, AfIII, Ajul(x), Alfl(x), Alol(x), AsiSI, Bael, Barl(x), BbvCI, Blpl, Bpll(x), BsiWI, BsmBI, BspDI, BspEI, BstAPI, BstBI, BstEII, ClaI, EcoNI, Esp3I, FseI, FspAI(x), KfII(x), MauBI(x), Mrel(x), PasI(x), PfoI(x), PshAI, Psrl(x), SexAl, SgrAl, Srfl(x), Stul, Xcml (x) = enzyme not available from NEB

pSNAP, Vector is a mammalian expression plasmid intended for the cloning and stable or transient expression of SNAP-tag® protein fusions in mammalian cells. This plasmid encodes SNAP,, a SNAP-tag protein, which is expressed under control of the CMV promoter. SNAP, is an improved version of the SNAPtag which exhibits faster labeling kinetics. The SNAP-tag is a novel tool for protein research, allowing the specific, covalent attachment of virtually any molecule to a protein of interest. The SNAP-tag is a small protein based on human O6-alkylguanine-DNA-alkyltransferase (hAGT). SNAP-tag substrates are derivatives of benzyl purines and benzyl pyrimidines. In the labeling reaction, the substituted benzyl group of the substrate is covalently attached to the SNAP-tag. Use of this system involves two steps: sub-cloning and expression of the protein of interest as a SNAP-tag fusion, and labeling of the fusion with the SNAP-tag substrate of choice. Further details are provided with the SNAP-Cell Starter Kit (NEB #E9100) and SNAP-Surface Starter Kit (NEB #E9120).

Codon usage of the gene is optimized for expression in mammalian cells. pSNAP, contains two multiple cloning sites to allow cloning of the fusion partner as a fusion to the N- or C-terminus of the SNAP-tag. The expression vector has an Internal Ribosome Entry Site (IRES) and a neomycin resistance gene downstream of the SNAP-tag for the efficient selection of stable transfectants.

Enzymes with unique restriction sites are shown in **bold** type. Location of sites of all NEB restriction enzymes for select plasmid can be found on the NEB website (choose Tools &

Resources > DNA Sequences and Maps tool). Restriction site coordinates refer to the position of the 5'-most base on the top strand in each recognition sequence.

Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons. Component genes or regions of fusion ORFs are indented below the ORF itself.

pUC19 origin of replication coordinates include the region from the -35 promoter sequence of the RNAII transcript to the RNA/DNA switch point. bla (ApR) gene coordinates include the signal sequence.


Feature	Coordinates	Source
CMV promoter	251-818	-
expression region	915-1564	-
MCS1	915-965	-
SNAP,	969-1514	-
MCS2	1515-1564	-
IRES	1910-2500	ECMV
Neo ^R	2536-3339	Tn5
origin	4094-4682	pUC19
bla (ApR)	4853-5713	Tn3

ori = origin of replication

Ap = ampicillin

Neo = neomycin

IRES = internal ribosomal entry site



MCS₂

AscI

MCS₁

Xhol Pmel Pacl BamHI CCTGCA GGCGGATCCG CGTTTAAACT CGAGGTTAAT TAATGAGCGG CCGC GGACGT CCGCCTAGGC GCAAATTTGA GCTCCAATTA ATTACTCGCC GGCG... Sequence file available at www.neb.com. See page 214 for ordering information.

Feature	Coordinates	Source
bla (ApR)	140-1000	Tn3
M13 origin	1042-1555	M13
origin	1666-2254	pMB1
rop	2814-2623	pMB1
lacl	4453-3371	E. coli
T7 promoter	5637-5654	T7
expression ORF	5725-6558	-
MCS	5722-5775	-
Mxe GyrA intein	5776-6369	M. xenopi
CBD	6400-6558	B. circulans

ori = origin of replication Ap = ampicillin

There are no restriction sites for the following enzymes: Aarl(x), Acc651, AfIII, Alel, AscI, AsiSI, AvrII, Bael, BbvCI, BgIII, BmgBI, Bpu10I, BseRI, BspDI, BstBI, Bsu36I, ClaI, CspCI, Eco53KI, FseI, FspAI(x), HindIII, I-CeuI, I-SceI, KpnI, MscI, NcoI, NsiI, PI-PspI, PI-SceI, PacI, PmII, PpuMI, RsrII, SacI, SanDI(x), SbfI, SexAI, SfiI, SmaI, SnaBI, SrfI(x), TspMI, XmaI

(x) = enzyme not available from NEB

pTXB1 is an *E. coli* plasmid cloning vector designed for recombinant protein expression, purification, and ligation using the IMPACT™ Kit (NEB #E6901) (1,2). It contains the pMB1 origin of replication from pBR322 and is maintained at a similar copy number to pBR322; in addition, pTXB1 also contains an M13 origin of replication.

The multiple cloning site (MCS) is positioned to allow translational fusion of the *Mxe* GyrA intein tag to the C-terminus of the cloned target protein (2,3). The chitin binding domain (CBD) from *B. circulans*, fused to the C-terminus of the intein, facilitates purification of the intein-target protein precursor.


Transcription of the gene fusion is controlled by the inducible T7 promoter, requiring $E.\ coli$ strains containing integrated copies of the T7 RNA polymerase gene [e.g., C2566 or BL21(DE3)] for expression. Basal expression from the T7 promoter is minimized by the binding of the Lac repressor, encoded by the lacl gene, to the lac operator immediately downstream of the T7 promoter (4). Translation of the fusion utilizes the translation initiation signal (Shine Dalgarno sequence) from the strongly expressed T7 gene 10 protein (ϕ 10).


pTXB1 and pTXB3 are identical except for the MCS regions: pTXB1 contains an Ndel site, and pTXB3 an Ncol site, overlapping the initiating methionine codon of the intein fusion gene. The N-terminal cysteine residue ("Cys,") of the intein is shaded.

Enzymes with unique restriction sites are shown in **bold** type. Location of sites of all NEB restriction enzymes for select plasmids can be found on the NEB website (choose Tools & Resources > DNA Sequences and Maps tool). Restriction site coordinates refer to the position of the 5´-most base on the top strand in each recognition sequence.

Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons. Component genes or regions of fusion ORFs are indented below the ORF itself.

pMB1 origin of replication coordinates include the region from the -35 promoter sequence of the RNAII transcript to the RNA/DNA switch point. For the M13 origin, the arrow shows the direction of synthesis of the (+) strand, which gets packaged into phage particles. *bla* (Ap[®]) gene coordinates include the signal sequence.





Spel
ATCACGGGAGATGCACTAGTTGCCCTACCCGAGGGCGAGTCGGTACGCATCGCCGACATCGTGCCG...

- (1) Chong, S. et al. (1997) *Gene*, 192, 271–281
- (2) Evans, T.C., Benner and Xu, M.-Q. (1998) *Protein Sci.*, 7, 2256–2264.
- (3) Southworth, M.W. et al. (1999) *Biotech-niques*, 27, 110–120.
- (4) Dubendorff, J.W. and Studier, F.W. (1991)J. Mol. Biol., 219, 45–59.

pTYB21

Sequence file available at www.neb.com. See page 214 for ordering information.

Feature	Coordinates	Source
bla (ApR)	140-1000	Tn3
M13 origin	1042-1555	M13
origin	1666-2254	pMB1
rop	2814-2623	pMB1
lacl	4453-3371	E. coli
T7 promoter	5637-5654	T7
expression ORF	5725-7368	_
MCS	7301-7361	_
Sce VMA intein	5770-7299	S. cerevisiae
CBD	6595-6747	B. circulans

ori = origin of replication Ap = ampicillin

There are no restriction sites for the following enzymes: Aarl(x), Aatll, Aflll, Agel, Ascl, AsiSI, AvrII, BbvCI, BmgBI, BseRI, BsiWI, BsmI, BspDI, Bsu36I, Clal, CspCI, Fsel, FspAI(x), I-Ceul, I-Scel, Nrul, Nsil, PI-Pspl, PI-Scel, Pacl, PaeR7I, PpuMI, PspXI, RsrII, SanDI(x), SexAI, SfiI, SgrAI, Smal, SnaBl, Srfl(x), TspMl, Xhol, Xmal, Zral (x) = enzyme not available from NEB

pTYB21 is an E. coli plasmid cloning vector designed for recombinant protein expression and purification using the IMPACT™ Kit (NEB #E6901) (1,2). It contains the pMB1 origin of replication from pBR322 and is maintained at a similar copy number to pBR322; in addition, pTYB21 also contains an M13 origin of replication.

The multiple cloning site (MCS) is positioned to allow translational fusion of the Sce VMA intein tag to the N-terminus of the cloned target protein (2). The chitin binding domain (CBD) from B. circulans, facilitates purification of the intein-target protein precursor.

Transcription of the gene fusion is controlled by the inducible T7 promoter, requiring *E. coli* strains containing integrated copies of the T7 RNA polymerase gene [e.g., C2566 or BL21(DE3)] for expression. Basal expression from the T7 promoter is minimized by the binding of the Lac repressor, encoded by the lacl gene, to the lac operator immediately downstream of the T7 promoter (3). Translation of the fusion utilizes the translation initiation signal


(Shine Dalgarno sequence) from the strongly expressed T7 gene 10 protein (φ10).

pTYB21 contains a Sapl site which allows for cloning of a target gene without any extra amino acids. pTYB22 is identical to pTYB21 except for the MCS regions (see below). pTYB22 contains an Ndel site overlapping the initiating methionine codon of the intein fusion gene, pTYB21 differs from pTYB11 in that it contains a universal MCS that is compatible with all NEB expression systems.

Enzymes with unique restriction sites are shown in **bold** type. Location of sites of all NEB restriction enzymes for select plasmids can be found on the NEB website (Technical Reference > DNA Sequences and Maps). Restriction site coordinates refer to the position of the 5´-most base on the top strand in each recognition sequence.

Open reading frame (ORF) coordinates are in the form "translational start - translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons. Component genes or regions of fusion ORFs are indented below the ORF itself.

pMB1 origin of replication coordinates include the region from the -35 promoter sequence of the RNAII transcript to the RNA/DNA switch point. For the M13 origin, the arrow shows the direction of synthesis of the (+) strand, which gets packaged into phage particles. bla (ApR) gene coordinates include the signal sequence.





- (1) Chong et al. (1996) J. Biol. Chem., 271, 22159-22168
- (2) Chong et al. (1998) NAR, 26, 5109-5115.
- (3) Dubendorff, J.W. and Studier, F.W. (1991) J. Mol. Biol., 219, 45-59.

GenBank Accession #: L09137 See page 118 for ordering information.

Feature	Coordinates	Source
lacZ $lpha$	469-146	_
origin	1455-867	pMB1 (mutant)
bla (ApR)	2486-1626	Tn3

ori = origin of replication Ap = ampicillin

There are no restriction sites for the following enzymes: Aarl(x), Afel, Aflll, Agel, Alel, Apal, Ascl, AsiSl, Avril, Bael, Bbsl, BbvCl, Bcll, Bglll, Blpl, BmgBl, Bmtl, Bpu10l, BsaAl, BsaBl, BseRl, Bsgl, BsiWl, BsmFl, Bsml, BspDl, BspEl, BsrGl, BssHll, BstBl, BstEll, BstXl, Bst2171, Bsu36l, Btgl, BtgZl, Clal, CspCl, Dralli, Eagl, EcoNl, EcoRV, Fsel, FspAl(x), Hpal, I-Ceul, I-Scel, Mfel, Mlul, Mscl, Nael, Nool, NgoMIV, Nhel, Notl, Nrul, Nsil, PI-Pspl, PI-Scel, Pacl, PaeR7I, PfiFl, PfiMl, Pmel, Pmll, PpuMl, PshAl, Psil, PspOMl, PspXl, RsrIl, SacIl, SanDl(x), SexAl, Sfil, SgrAl, SnaBl, Spel, SrfI(x), Stul, Styl, Swal, Tth1111, Xcml, Xhol

(x) = enzyme not available from NEB

pUC19 is a small, high-copy number $\it E.~coli$ plasmid cloning vector containing portions of pBR322 and M13mp19 (1). It contains the pMB1 origin of replication from pBR322, but it lacks the $\it rop$ gene and carries a point mutation in the RNAII transcript (G 2975 in pBR322 to A 1308 in pUC19; 2). These changes together result in a temperature-dependent copy number of about 75 per cell at 37°C and > 200 per cell at 42°C (2,3). The multiple cloning site (MCS) is in frame with the $\it lacZ\alpha$ gene, allowing screening for insertions using $\it \alpha$ -complementation.

 $\,$ pUC18 is identical to pUC19 except that the MCS region (nt 397-454) is inverted.

pNEB193 is also identical to pUC19 except for the addition of several restriction endonuclease sites to the MCS. Its total length is 2713 bp.


Enzymes with unique restriction sites are shown in **bold** type, and enzymes with two restriction sites are shown in regular type. Location of sites of all NEB restriction enzymes for select plasmids can be found on the NEB website (choose Tools & Resources > DNA Sequences and Maps tool). Restriction site coordinates refer to the position of the 5´-most base on the top strand in each recognition sequence.

Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons.

Origin of replication coordinates include the region from the -35 promoter sequence of the RNAII transcript to the RNA/DNA switch point. *bla* (Ap[®]) gene coordinates include the signal sequence.

References

- (1) Yanisch-Perron, C., Vieira, J. and Messing, J. (1985) *Gene*, 33, 103–119.
- (2) Lin-Chao, S., Chen, W.-T. and Wong, T.-T. (1992) Mol. Microbiol., 6, 3385–3393
- (3) Miki, T. et al. (1987) Protein Eng., 1, 327-332.



pUC19 MCS



pNEB193 MCS

Ec agtgA	oRI ATT	Sa CGA		Kp CGG	nl	Smal CCG		Bssl scl GCG			mHI ATC			TAA	Xb GTC		Sal AGT			Pmel iTTT		Sbfl Ps CCT		Sp .GGC			ndIII GCT	TGG	cgt	aat	cat	ggt	cat
	40	0		4	10			420			43	0		4	40			450			46	0		4	70			480			490	0	Τ
s	N	S	S	Р	٧	R	Р	R	Α	Р	D	K	Ι	L	D	L	Т	S	Q	K	F	R	С	Α	Н	L	S	Р	T	Ι	М	Т	М
																						<u> </u>				– Ia	ιcΖ α	tran	slati	onal	star	t —	

T7

39,937 base pairs GenBank Accession #: NC_001604 Not currently available from NEB.

There are no restriction sites for the following enzymes: Afel, Apal, Ascl, AsiSI, BamHI, BsiWI, BspEI, Eagl, Eco53KI, EcoRI, EcoRV, Fsel, HindIII, I-Ceul, I-Scel, Nael, NgoMIV, Notl, PI-PspI, PI-Scel, PaeR7I, PspOMI, PspXI, PstI, PvuI, Sacl, SacII, SalI, SbfI, SexAI, SgrAI, Smal, SphI, SrfI(x), TspMI, XhoI, Xmal


(x) = enzyme not available from NEB

T7 is a lytic *E. coli* bacteriophage with a linear, double-stranded DNA genome containing 56 genes (1-4). Genes are classified as early or late based on the order of transcription in the infected host and their dependence on host or phage RNA polymerase.

Numbering of the sequence begins at the first (5'-most) base of the left end (bottom of the diagram below) and continues rightward (upward) in the direction of early to late genes. The map below shows the positions of all known ORFs larger than 200 codons.


Enzymes with unique restriction sites are shown in **bold** type, and enzymes with two restriction sites are shown in regular type. Location of sites of all NEB restriction enzymes for select plasmids can be found on the NEB website (choose Tools & Resources > DNA Sequences and Maps tool). Restriction site coordinates refer to the position of the 5´-most base on the top strand in each recognition sequence.

- Oakley, J.L. and Coleman, J.E. (1977) Proc. Natl. Acad. Sci. USA, 74, 4266–4270.
- (2) Stahl, S.J. and Zinn, K. (1981) J. Mol. Biol., 148, 481-485.
- (3) Dunn, J.J. and Studier, F.W. (1981) *J. Mol. Biol.*, 148, 303–330.
- (4) Dunn, J.J. and Studier, F.W. (1983) J. Mol. Biol., 166, 477–535



Amino Acid Structures

Each amino acid is accompanied by its three- and one-letter code, residue molecular weight (actual molecular weight minus water) and side-chain pK, where appropriate.



APPENDIX

DNA Base Pairs

The structures of the adenosine:thymidine and guanosine:cytidine base pairs are shown in the context of the ribose phosphodiester backbones. The numbering schemes of the ribose and nucleotide moieties are indicated. Arrows indicate the polarity of each strand from 5′ to 3′.

Nucleic Acid Data

Average weight of a DNA basepair (sodium salt) = 650 daltons

1.0 A_{260} unit ds DNA = 50 μ g/ml = 0.15 mM (in nucleotides)

1.0 A_{260}^{200} unit ss DNA = 33 μ g/ml = 0.10 mM (in nucleotides)

 $1.0 A_{260}^{200}$ unit ss RNA = $40 \mu g/ml = 0.11 mM$ (in nucleotides)

MW of a double-stranded DNA molecule = (# of base pairs) x (650 daltons/base pair)

Moles of ends of a double-stranded DNA molecule = 2 x (grams of DNA) / (MW in daltons)

Moles of ends generated by restriction endonuclease cleavage:

a) circular DNA molecule: 2 x (moles of DNA) x (number of sites)

b) linear DNA molecule: 2 x (moles of DNA) x (number of sites) + 2 x (moles of DNA)

 $1 \mu g$ of 1000 bp DNA = 1.52 pmol = 9.1 x 10^{11} molecules

 $1 \mu g$ of pUC18/19 DNA (2686 bp) = 0.57 pmol = 3.4×10^{11} molecules

1 μ g of pBR322 DNA (4361 bp) = 0.35 pmol = 2.1 x 10¹¹ molecules

1 μ g of M13mp18/19 DNA (7249 bp) = 0.21 pmol = 1.3 x 10¹¹ molecules

1 μg of λ DNA (48502 bp) = 0.03 pmol = 1.8 x 10¹⁰ molecules

1 pmol of 1000 bp DNA = $0.66 \mu g$

1 pmol of pUC18/19 DNA (2686 bp) = $1.77 \mu g$

1 pmol of pBR322 DNA (4361 bp) = 2.88 µg

1 pmol of M13mp18/19 DNA (7249 bp) = $4.78 \mu g$

1 pmol of λ DNA (48502 bp) = 32.01 μ g

1.0 kb DNA = coding capacity for 333 amino acids ≈ 37,000 dalton protein

10,000 dalton protein ≈ 270 bp DNA

50,000 dalton protein ≈ 1.35 kb DNA

Isotope Data

	Particle		
Isotope	Emitted	Half Life	
¹⁴ C	β	5,730 years	1 Ci = 1,000 mCi
3H	β	12.3 years	1 mCi = 1,000 μCi
125	γ	60 days	$1 \mu Ci = 2.2 \times 10^6 \text{ disintegrations/minute}$
³² P	β	14.3 days	1 Becquerel = 1 disintegration/second
³³ P	β	25 days	$1 \mu Ci = 3.7 \times 10^4 Becquerels$
³⁵ S	β	87.4 days	1 Becquerel = 2.7 x 10 ⁻⁵ μCi

Acids and Bases			Molecular	Specific	% by	Conc Reagent	
	Compound	Formula	Weight	Gravity	Weight	Molarity	
	Acetic acid, glacial	CH₃COOH	60.0	1.05	99.5	17.4	
	Formic acid	HC00H	46.0	1.20	90	23.4	
	Hydrochloric acid	HCI	36.5	1.18	36	11.6	
	Nitric acid	HNO ₃	63.0	1.42	71	16.0	
	Perchloric acid	HCIO ₄	100.5	1.67	70	11.6	
	Phosphoric acid	H ₃ PO ₄	98.0	1.70	85 06	18.1	
	Sulfuric acid	H ₂ SO ₄	98.1 35.0	1.84 0.90	96 28	18.0	
	Ammonium hydroxide Potassium hydroxide	NH₄OH KOH	56.1	1.52	20 50	14.8 13.5	
	Sodium hydroxide	NaOH	40.0	1.53	50	19.1	
	β-mercaptoethanol	HSCH ₂ CH ₂ OH		1.11	100	14.3	
Protein Data			Bacterial Cells: <i>E.</i>	coli or Salmonei	la typhimur	rium	
			Cell Data	per co		per liter at 10º cells per ml	
	The control of the co	C 4 Pt 0					
	Theoretical maximum yield		Wet Weight	9.5 x 10		0.95 g	
	(10 ⁹ cells /ml) if protein of		Dry Weight	2.8 x 10		0.28 g	
	0.1% of total protein: 3 2.0% of total protein: 3		Total Protein Volume	1.55 x 10 1.15 µm³ = 1		0.15 g	
	50.0% of total protein:		Protein Conc. in th				
Common Plasmid				Gene Product		Molecular Weight	
		Gene		# of Residues		(daltons)	
Gene Products		<i>tet</i> (pBR322)		401		43,267	
		<i>amp</i> (pBR322	hla)	286		45,267 31,515	
		kan (pACYC1		264		29,047	
		cam (pACYC1		219		25,663	
		lacZα (pUC19		107		12,232	
		lacZ	,	1,023		116,351	
				,			
 Nucleotide			Molecular	λ max	Ab	osorbance at λ max	
		Compound		· 			
Nucleotide Physical Properties		•	Weight	λ max (pH 7.0		M solution (pH 7.0)	
		Compound ATP CTP		λ max		M solution (pH 7.0) 15,400	
		ATP CTP	Weight 507.2 483.2	λ max (pH 7.0 259 271		M solution (pH 7.0) 15,400 9,000	
		ATP	Weight 507.2	λ max (pH 7.0 259		M solution (pH 7.0) 15,400 9,000 13,700	
		ATP CTP GTP	Weight 507.2 483.2 523.2	λ max (pH 7.0 259 271 253		M solution (pH 7.0) 15,400 9,000 13,700 10,000	
		ATP CTP GTP UTP	Weight 507.2 483.2 523.2 484.2	λ max (pH 7.0 259 271 253 262		M solution (pH 7.0) 15,400 9,000 13,700	
		ATP CTP GTP UTP dATP	Weight 507.2 483.2 523.2 484.2 491.2	λ max (pH 7.0 259 271 253 262 259		M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200	
		ATP CTP GTP UTP dATP dCTP	Weight 507.2 483.2 523.2 484.2 491.2 467.2	λ max (pH 7.0 259 271 253 262 259 271		M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300	
Physical Properties	pH of Tris Buffer	ATP CTP GTP UTP dATP dCTP dCTP dGTP dTTP	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ max (pH 7.0 259 271 253 262 259 271 253 267	1 1	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600	
Physical Properties pH vs Temperature	5°C 25°C	ATP CTP GTP UTP dATP dCTP dCTP dGTP dTTP	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ max (pH 7.0 259 271 253 262 259 271 253 267	e Gel	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600 Resolution	
Physical Properties	5°C 25°C 7.76 7.20	ATP CTP GTP UTP dATP dCTP dCTP dGTP dTTP (0.05 M) 37°C 6.91	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ max (pH 7.0 259 271 253 262 259 271 253 267	e Gel	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600 Resolution um Resolution	
Physical Properties pH vs Temperature	5°C 25°C 7.76 7.20 7.89 7.30	ATP CTP GTP UTP dATP dCTP dGTP dTTP (0.05 M) 37°C 6.91 7.02	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ max (pH 7.0 259 271 253 262 259 271 253 267	e Gel	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600 Resolution um Resolution near DNA (kb)	
Physical Properties pH vs Temperature	5°C 25°C 7.76 7.20 7.89 7.30 7.97 7.40	ATP CTP GTP UTP dATP dCTP dGTP dTTP	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ. max (pH 7.0 259 271 253 262 259 271 253 267 Agarose	e Gel	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600 Resolution um Resolution near DNA (kb) 30 to 1.0	
Physical Properties pH vs Temperature	5°C 25°C 7.76 7.20 7.89 7.30 7.97 7.40 8.07 7.50	ATP CTP GTP UTP dATP dCTP dGTP dTTP (0.05 M) 37°C 6.91 7.02 7.12 7.22	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ. max (pH 7.0 259 271 253 262 259 271 253 267 Agarose	e Gel Optimi	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600 Resolution mear DNA (kb) 30 to 1.0 12 to 0.8	
Physical Properties pH vs Temperature	5°C 25°C 7.76 7.20 7.89 7.30 7.97 7.40 8.07 7.50 8.18 7.60	ATP CTP GTP UTP dATP dCTP dGTP dTTP (0.05 M) 37°C 6.91 7.02 7.12 7.22 7.30	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ. max (pH 7.0 259 271 253 262 259 271 253 267 Agarose % Gel 0.5 0.7 1.0	e Gel Optimi	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600 Resolution um Resolution near DNA (kb) 30 to 1.0 12 to 0.8 10 to 0.5	
Physical Properties pH vs Temperature	5°C 25°C 7.76 7.20 7.89 7.30 7.97 7.40 8.07 7.50	ATP CTP GTP UTP dATP dCTP dGTP dTTP (0.05 M) 37°C 6.91 7.02 7.12 7.22	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ. max (pH 7.0 259 271 253 262 259 271 253 267 Agarose % Gel 0.5 0.7 1.0 1.2	e Gel Optimi	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600 Resolution um Resolution near DNA (kb) 30 to 1.0 12 to 0.8 10 to 0.5 7 to 0.4	
Physical Properties pH vs Temperature	5°C 25°C 7.76 7.20 7.89 7.30 7.97 7.40 8.07 7.50 8.18 7.60 8.26 7.70 8.37 7.80	ATP CTP GTP UTP dATP dCTP dGTP dTTP (0.05 M) 37°C 6.91 7.02 7.12 7.22 7.30 7.40 7.52	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ. max (pH 7.0 259 271 253 262 259 271 253 267 Agarose % Gel 0.5 0.7 1.0	e Gel Optimi	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600 Resolution um Resolution near DNA (kb) 30 to 1.0 12 to 0.8 10 to 0.5	
Physical Properties pH vs Temperature	5°C 25°C 7.76 7.20 7.89 7.30 7.97 7.40 8.07 7.50 8.18 7.60 8.26 7.70	ATP CTP GTP UTP dATP dCTP dGTP dTTP (0.05 M) 37°C 6.91 7.02 7.12 7.22 7.30 7.40	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ. max (pH 7.0 259 271 253 262 259 271 253 267 Agarose % Gel 0.5 0.7 1.0 1.2	e Gel Optimi	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600 Resolution um Resolution near DNA (kb) 30 to 1.0 12 to 0.8 10 to 0.5 7 to 0.4	
Physical Properties pH vs Temperature	5°C 25°C 7.76 7.20 7.89 7.30 7.97 7.40 8.07 7.50 8.18 7.60 8.26 7.70 8.37 7.80 8.48 7.90 8.58 8.00 8.68 8.10	ATP CTP GTP UTP dATP dCTP dGTP dTTP (0.05 M) 37°C 6.91 7.02 7.12 7.22 7.30 7.40 7.52 7.62 7.71 7.80	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ. max (pH 7.0 259 271 253 262 259 271 253 267 Agarose % Gel 0.5 0.7 1.0 1.2	e Gel Optimi	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600 Resolution um Resolution near DNA (kb) 30 to 1.0 12 to 0.8 10 to 0.5 7 to 0.4	
Physical Properties pH vs Temperature	5°C 25°C 7.76 7.20 7.89 7.30 7.97 7.40 8.07 7.50 8.18 7.60 8.26 7.70 8.37 7.80 8.48 7.90 8.58 8.00 8.68 8.10 8.78 8.20	ATP CTP GTP UTP dATP dCTP dGTP dTTP (0.05 M) 37°C 6.91 7.02 7.12 7.22 7.30 7.40 7.52 7.62 7.71 7.80 7.91	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ. max (pH 7.0 259 271 253 262 259 271 253 267 Agarose % Gel 0.5 0.7 1.0 1.2	e Gel Optimi	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600 Resolution um Resolution near DNA (kb) 30 to 1.0 12 to 0.8 10 to 0.5 7 to 0.4	
Physical Properties pH vs Temperature	5°C 25°C 7.76 7.20 7.89 7.30 7.97 7.40 8.07 7.50 8.18 7.60 8.26 7.70 8.37 7.80 8.48 7.90 8.58 8.00 8.68 8.10 8.78 8.20 8.88 8.30	ATP CTP GTP UTP dATP dCTP dGTP dTTP (0.05 M) 37°C 6.91 7.02 7.12 7.22 7.30 7.40 7.52 7.62 7.71 7.80 7.91 8.01	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ. max (pH 7.0 259 271 253 262 259 271 253 267 Agarose % Gel 0.5 0.7 1.0 1.2	e Gel Optimi	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600 Resolution um Resolution near DNA (kb) 30 to 1.0 12 to 0.8 10 to 0.5 7 to 0.4	
Physical Properties pH vs Temperature	5°C 25°C 7.76 7.20 7.89 7.30 7.97 7.40 8.07 7.50 8.18 7.60 8.26 7.70 8.37 7.80 8.48 7.90 8.58 8.00 8.68 8.10 8.78 8.20 8.88 8.30 8.98 8.40	ATP CTP GTP UTP dATP dCTP dGTP dTTP (0.05 M) 37°C 6.91 7.02 7.12 7.22 7.30 7.40 7.52 7.62 7.71 7.80 7.91 8.01 8.10	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ. max (pH 7.0 259 271 253 262 259 271 253 267 Agarose % Gel 0.5 0.7 1.0 1.2	e Gel Optimi	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600 Resolution um Resolution near DNA (kb) 30 to 1.0 12 to 0.8 10 to 0.5 7 to 0.4	
Physical Properties pH vs Temperature	5°C 25°C 7.76 7.20 7.89 7.30 7.97 7.40 8.07 7.50 8.18 7.60 8.26 7.70 8.37 7.80 8.48 7.90 8.58 8.00 8.68 8.10 8.78 8.20 8.88 8.30	ATP CTP GTP UTP dATP dCTP dGTP dTTP (0.05 M) 37°C 6.91 7.02 7.12 7.22 7.30 7.40 7.52 7.62 7.71 7.80 7.91 8.01	Weight 507.2 483.2 523.2 484.2 491.2 467.2 507.2	λ. max (pH 7.0 259 271 253 262 259 271 253 267 Agarose % Gel 0.5 0.7 1.0 1.2	e Gel Optimi	M solution (pH 7.0) 15,400 9,000 13,700 10,000 15,200 9,300 13,700 9,600 Resolution um Resolution near DNA (kb) 30 to 1.0 12 to 0.8 10 to 0.5 7 to 0.4	

Common Buffer Chart

The following chart lists some of the common buffers used in biology. The useful buffer range is the pK $_a$ \pm 0.5–1 pH unit. The buffering capacity decreases beyond this range.

COMMON NAME	pK _a AT 25°C	MOLECULAR WEIGHT	CHEMICAL FORMULA	CHEMICAL NAME
Phosphate	2.12	98.00	H ₃ PO ₄ -	-
Acetate	4.76	60.00	CH_3CO_2H or $C_2H_4O_2$	-
MES	6.15	195.20	$C_6H_{13}NO_4S$	2-(N-morpholino)ethanesulfonic acid
PIPES	6.76	302.40	$C_8H_{18}N_2O_6S_2$	piperazine-N,N´-bis(2-ethanesulfonic acid)
Imidazole	6.95	68.08	$C_3H_4N_2$	1,3-Diaza-2,4-cyclopentadiene
MOPS	7.20	209.30	$C_7 H_{15} NO_4 S$	3-(N-morpholino)propanesulfonic acid
Phosphate	7.21	97.00	$H_2PO_4^-$	-
TES	7.40	229.20	$C_6H_{15}NO_6S$	N-Tris(hydroxymethyl)methyl-2-aminoethanesulfonic acid
HEPES	7.48	238.30	$C_8 H_{18} N_2 O_4 S$	4-(2-hydroxyethyl)-1-piperazineethanesulfonic acid
Tricine	8.05	179.20	$C_6H_{13}NO_5$	N-(2-Hydroxy-1,1-bis(hydroxymethyl)ethyl)glycine
Tris	8.06	121.14	$C_4H_{11}NO_3$	Tris(hydroxymethyl)methylamine
Bicine	8.35	163.20	C ₆ H ₁₃ NO ₄	N,N-bis(2-hydroxyethyl)glycine
TAPS	8.43	243.30	$C_7H_{17}NO_6S$	N-Tris(hydroxymethyl)methyl-3-aminopropanesulfonic acid
Phosphate	12.67	96.00	HPO ₄ ²⁻	-



NEB expanded its manufacturing footprint by opening a facility in Rowley, MA for production of GMP-grade materials. Pictured here are several of the team members.

#	
1 kb DNA Ladder	167
1 kb Plus DNA Ladder	167
3´-Biotin-GTP	191
5´ Deadenylase	199
3´ Desthiobiotin-GTP	191
5´ DNA Adenylation Kit	199
9° N DNA Ligase	
50 bp DNA Ladder	167
100 bp DNA Ladder	167
5-hmC Identification	
5-hydroxymethyluridine DNA Kinase	99
5-mC Identification	262, 273
5-methyl-dCTP	271

A	
α -N-Acetylgalactosaminidase	243
α1-2 Fucosidase	
α1-2.3.4.6 Fucosidase	
α1-2,4,6 Fucosidase 0	245
α1-3.4 Fucosidase	
α1-3,6 Galactosidase	245
α1-3,4,6 Galactosidase	
α1-6 Mannosidase	
α1-2,3 Mannosidase	
α1-2,3,6 Mannosidase	
α2-3 Neuraminidase S	
α2-3,6,8 Neuraminidase	
α2-3,6,8,9 Neuraminidase A	
Aatll	
AbaSI	
Accl	
Acc65I	
Acil	
Acids and Bases	
AcII	
Activity Chart for Restriction Enzymes	
Activity of DNA Modifying Enzymes in CutSmart Buffer	
Acul	
Acyclonucleotide Set	
Adaptors, Sequencing	
Adenosine 5´-Triphosphate (ATP)	
Afel	
AfIII	
AfIIII	
Afu Uracil-DNA Glycosylase (UDG)	111
Agarase (β)	112
Agel	21
Agel-HF	21
Ahdl	21
AleI-v2	21
Alkaline Phosphatase, Calf Intestinal (CIP)	110
Alul	21
Alul Methyltransferase	114
Alwl	21
AlwNI	21
Amino Acid Structures	382
Amplification, DNA	60-81
Amylose Magnetic Beads	213
Amylose Resin	213
Amylose Resin High Flow	213
AMV Reverse Transcriptase	80, 191, 192
Antarctic Phosphatase	
Antarctic Thermolabile UDG	112
Antibodies	
Anti-CBD Monoclonal Antibody	
Anti-M13 pIII Monoclonal Antibody	260
Anti-MBP Magnetic Beads	
Anti-MBP Monoclonal Antibody	213

Anti-SNAP-tag Antibody, Polyclonal	281
Goat Anti-Mouse IgG Magnetic Beads	219
Goat Anti-Rabbit IgG Magnetic Beads	219
Goat Anti-Rat IgG Magnetic Beads	219
Anti-MBP Magnetic Beads	213
Apal	22
ApaLI	22
APE I Repair Protein	109
ApeKI	22
Apol	22
Apol-HF	22
Appendix	
Apyrase	101
ArtCorps	6
Ascl	22
Asel	22
AsiSI	22
AspN, Endoproteinase	253
ATP	79
Aval	23
Avall	23
Avoiding Ribonuclease Contamination	184
Avrll	23

В	
β-Agarase I	112
β1-3 Galactosidase	245
β1-3,4 Galactosidase	
β1-4 Galactosidase S	246
β- <i>N</i> -Acetylglucosaminidase S	243
β-N-Acetylhexosaminidase _t	
Bacteriodes Heparinase I, II, III	
Bael	
BaeGI	
BamHI	
BamHI-HF.	
BamHI Methyltransferase	
Banl	
Banll	
Basic Research at NEB	
Bbsl	
BbsI-HF	
Bbvl	
BbvCI BC-NH2	
Bccl Scale S	
BceAl	
Bcgl	
BciVI	
BcII	
BcII-HF	
BcoDI	
Bfal	
BfuAl	
BG-GLA-NHS	
BG-Maleimide	
BG-NH2	283
BG-PEG-NH2	283
BgII	25
BgIII	25
BioBrick Assembly Kit	91
Biotin Labels	282
Bisulfite Sequencing	141, 263
BL21 Competent E. coli	217, 224, 225, 229
BL21(DE3) Competent E.coli	217, 224, 225, 229
Blpl	
Blue Gel Loading Dye	57
Blue Loading Buffer Pack	171
Blue Prestained Protein Std., Broad range (11–250 kDa)	171
Blunt/TA Ligase Master Mix	96

Blunting, DNA		
Exonuclease T	,	
Klenow DNA Polymerase		
Quick Blunting Kit		
T4 DNA Polymerase		
T4 Polynucleotide Kinase BmgBl		
Bmrl		
Bmtl		
Bmtl-HF		
Bpml		
Bpu10l		
BpuEl		26
Bsal		27
Bsal-HFv2		27
BsaAl		27
BsaBl		
BsaHI		
BsaJI		
BsaWI		
BsaXI		
BseRI		
BseYI		
Bsgl		
BsiEl		
BsiHKAI		
BsiWI		
BSII BSII		
Bsml		
BsmAl		
BsmBl		
BsmFl .		
BsoBl		
Bsp1286I		
BspCNI		
BspDI		
BspEl		
BspHI		
BspMI		30
BspQI		30
Bsrl		30
BsrBl		30
BsrDI		
BsrFI-v2		
BsrGl		
BsrGI-HF		
BssHII		
BssSI-v2		
Bst 2.0, Bst 2.0 WarmStart DNA Polymerases		
Bst 3.0 DNA Polymerase BstAPI		
BStBI		
Bst DNA Polymerase, Full Length.		
Bst DNA Polymerase, Large Fragment		
BstEII		
BstFII-HF		
BstNI		
BstUl		
BstXI		
BstYI		32
BstZ17I-HF		32
Bsu DNA Polymerase, Large Fragment		
Bsu36l		32
Btgl		32
BtgZI		32
Btsl-v2		33
BtsIMutl		
BtsCl		33
Buffers		
Common Buffer Chart		
Diluent, Restriction Endonucleases	57,	293

Double Digestion, Restriction Enzymes	29
Reaction, Restriction Endonucleases	56, 293
Thermophilic Polymerases, Reaction Buffers	79
Tris Ruffer nH of as a Function of Temperature	38,

C	
Cac8l	33
Calf Intestinal Alkaline Phosphatase (CIP)	
cAMP-dependent Protein Kinase	257
Cap Analogs, RNA	190
Capping System, Vaccinia	
Cas9 Nuclease, S. pyogenes	,
Casein Kinase II	257
cDNA Synthesis	00 101 100
AMV Reverse Transcriptase M-MuLV Reverse Transcriptase	
Primers	
ProtoScript II Reverse Transcriptase	
ProtoScript II First Strand cDNA Synthesis Kit	
ProtoScript First Strand Synthesis Kit	
Selection Chart	
WarmStart RTx Reverse Transcriptase	80, 191, 193
cDNA, Primers	193
Cellular Analysis	
Cellular Imaging & Analysis	
Labeling with SNAP-tag Technology/Troubleshooting Guide	
ChIP-Seq, Library Preparation for Next Generation Sequencing	
Chitin Resin.	
Chitin Magnetic Beads	
Clal Cleavage Close to the End of DNA Fragments	
Cleavage of Supercoiled DNA.	
CLIP-Biotin	
CLIP-Cell 505	
CLIP-Cell Block	
CLIP-Cell TMR-Star	280
CLIP-Surface 488	280
CLIP-Surface 547	280
CLIP-Surface 647	280
CLIP-tag	
Overview	
Comparison to GFP	
Vectors	282
Cloning BioBrick Assembly Kit	0.1
Cloning Plasmids	
Competent Cells for Cloning	
Enhancing Transformation Efficiency	,
Getting Started with Molecular Cloning	
Gibson Assembly Cloning Kit	
NEB Golden Gate Assembly Mix	90
NEB PCR Cloning Kit	91
NEBuilder HiFi DNA Assembly Cloning Kit	
Optimization Tips	
Q5 Site Directed Mutagenesis Kit	
Traditional Cloning Quick Guide	
Troubleshooting Guide	
Workflow Comparison	80
One Tag DNA Polymerase	65
Tag DNA Polymerase	
Color Prestained Protein Standard, Broad Range (10–250 kDa)	
Combinatorial Libraries	
Combined Bisulfite Restriction Analysis, Substrates for	
Comet Assay	
Common Buffer Chart	
Compatible Cohesive Ends	332-333
Competent Cells	
Competent Cell Cloning and Protein Expression Strain Selection Chart	224
Comparent Call Clening Selection Chart	110

Competent Cells Strain Properties	224
Competent E. coli	
BL21	217, 224, 225, 229
BL21(DE3)	
dam-/ dcm-	224, 225, 228
Enhancing Transformation Efficiency	361
Lemo21(DE3)	217, 224, 225, 229
NEB 5-alpha (High Efficiency)	
NEB 5-alpha (Subcloning Efficiency)	
NEB 5-alpha Electrocompetent	
NEB 5-alpha F´/a (High Efficiency)	
NEB 10-beta (High Efficiency)	
NEB 10-beta Electrocompetent	
NEB Express (High Efficiency)	
NEB Express / (High Efficiency)	
NEB Stable Competent E. coli	
NEB Turbo (High Efficiency)	
NEB Turbo Electrocompetent	224, 225, 226
NiCo21(DE3)	
Protein Expression with T7 Express Strains	362
SHuffle Express	217, 224, 225, 232
SHuffle T7	217, 224, 225, 233
SHuffle T7 Express	217, 224, 225, 232
SHuffle T7 Express lysY	217, 224, 225, 233
T7 Express (High Efficiency).	217, 224, 225, 231
T7 Express /ysY (High Efficiency)	217, 224, 225, 231
T7 Express /ysY//q (High Efficiency)	217, 224, 225, 231
Competitor Cross Reference Tool	288
Control DNA	
CpG Methyltransferase (M.SssI)	114
Cre Recombinase	113
Creative Action Institute	6
Cross Index of Recognition Sequences.	304-305
CspCI	
Customized Solutions	4
CviAII	
CviKI-1	
CviQI	34

U	
5´ Deadenylase	199
3´ Desthiobiotin-GTP	191
5´ DNA Adenylation Kit	199
7-deaza-dGTP	
dam-/ dcm- Competent E. coli	217, 224, 225, 228
Dam, Dcm and CpG Methylation	
dam Methyltransferase	
dATP Solution	
Ddel	
Deep Vent DNA Polymerase	
Deep Vent (exo ⁻) DNA Polymerase	
Deoxynucleotide (dNTP) Solution Set Deoxynucleotide (dNTP) Solution Mix	
Diluent Buffers (Restriction Endonuclease)	
Distributors	,
DNA	IIISIUG DAGK GOVGI
Adenylation	199
Base Pair Diagrams	383
Cleanup, Monarch Products for	126, 127
Cloning Plasmids	118
Enrichment	149-151
Extraction	
Fragmentation	, ,
Input Guidelines for Library Prep	
Ladders	
Library Preparation	
Ligases	
Loading Dye	
Maps	
Markers	169

Methylation Analysis	
Modifying Enzymes	84-119
Phosphatases	99-100
Plasmid Miniprep, Monarch Products for	125
Polymerases	
Properties of DNA Repair Enzymes	107
Reagents for Library Preparation	134-161
Substrates	118
NA Markers and Ladders	
φX174 DNA-HaelII Digest	169
50 bp DNA Ladder	167
1 kb Plus DNA Ladder	167
100 bp DNA Ladder	167
1 kb DNA Ladder	167
Fast DNA Ladder	167
Lambda DNA-BstEII Digest	169
Lambda DNA-HindIII Digest	169
Lambda DNA-Mono Cut Mix	169
Lambda PFG Ladder	168
Low Molecular Weight DNA Ladder	167
Midrange PFG Ladder	168
pBR322 DNA-BstNI Digest	169
pBR322 DNA-Mspl Digest	169
PCR Marker	167
Quick-Load 1 kb Plus DNA Ladder	168
Quick-Load 1 kb DNA Ladder	168
Quick-Load 1 kb Extend DNA Ladder	168
Quick-Load 100 bp DNA Ladder	
Quick-Load Purple 1 kb Plus DNA Ladder	
Quick-Load Purple 1 kb DNA Ladder	
Quick-Load Purple 100 bp DNA Ladder	
Quick-Load Purple 50 bp DNA Ladder	
Quick-Load Purple Low Molecular Weight DNA Ladder	
Supercoiled DNA Ladder	
Tridye 1 kb Plus DNA Ladder	
Tridye 1 kb DNA Ladder	
Tridye 100 bp DNA Ladder	
NA Modifying Enzymes, Activity in CutSmart Buffer	
NA Modifying Enzymes and Cloning Technologies	
DNA Polymerase I	
DNA Polymerase I, Large (Klenow) Fragment	77
DNA Repair	
Glycosylases on Various Damaged Bases	
NEBNext FFPE DNA Repair Mix	
PreCR Repair Mix	
Properties of DNA Repair Enzymes	
Proteins	
NA Library Preparation Reagents for Next Generation Sequencing	
NA Sequences and Maps Tool	
Nase I (RNase-Free)	
Onmt1, Methyltransferase	
Oouble Digestion	
Oouble Digest Finder	
)pnl	
)pnll	
Oral	
OrallI-HF	
Ordl	
ISRNA L'adder	
ISRNA Ligase (T4 RNA Ligase 2)	197

t	
Eael	34
Eagl	35
Eagl-HF	35
Earl	35
Ecil	35
E. coli DNA Ligase	98
E. coli Ribosome.	216
E. coli RNA Polymerase, Core Enzyme	189

E. coli RNA Polymerase, Holoenzyme.	
E. coli Poly(A) Polymerase Eco53KI	
EcoNI	
Eco0109I	
EcoP15I	35
EcoRI	
EcoRI-HF	
EcoGII Methyltransferase EcoRI Methyltransferase	
EcoRV	
EcoRV-HF	36
EGCasel	
ElectroLigase	
EM-Seq	141
0-Glycosidase	242
Endo D.	239
Endo F2	
Endo F3	
Endo H, Endo H, Endo S	
PNGase A	
PNGase F (Glycerol Free, Recombinant)	
Rapid PNGase F (non-reducing format)	241
PNGase F, Remove-iT	
Rapid PNGase F Endonucleases	241
APE I	109
Common Applications	
DNase I (RNase-Free)	202
Endonuclease III (Nth)	
Endonuclease IV	
Endonuclease V Endonuclease VIII	
Homing	
Mung Bean Nuclease	
Nicking Endonucleases	
Properties of	
Restriction Endonucleases	
T7 Endonuclease I	
Thermostable FEN1	
Tma Endonuclease III	
Tth Endonuclease IV	
Endoproteinase AspN Endoproteinase GluC	
Endoproteinase LysC	
EnGen Lba Cas12a (Cpf1)	
EnGen Mutation Detection Kit	
EnGen Sau Cas9	
EnGen Spy Case Nickase	
EnGen Spy Cas9 Nickase	
EnGen Spy dCas9 (SNAP-tag)	
Enhancing Transformation Efficiency	
Enterokinase, light chain	
Environmental Commitment Enzymatic Methyl-Seq	
Enzyme Finder	
Epigenetics	
EpiMark 5-hmC and 5-mC Analysis Kit	
EpiMark Bisulfite Conversion Kit	
EpiMark Hot Start <i>Taq</i> DNA Polymerase EpiMark Methylated DNA Enrichment Kit	
EpiMark N6-Methyladenosine Enrichment Kit	203
EpiMark Nucleosome Assembly Kit	
E. coli DNA Polymerase I	
Esp3l	
ET SSB (Extreme Thermostable ssBinding Protein)	
Exo-CIP Rapid PCR Cleanup Kit	الا

Exoglycosidases	
α-N-Acetylgalactosaminidase	243
α1-2 Fucosidase	244
α1-2,3,4,6 Fucosidase	244
α1-2,4,6 Fucosidase 0	245
α1-3,4 Fucosidase	244
lpha1-3,6 Galactosidase	245
α1-3,4,6 Galactosidase	245
α1-6 Mannosidase	247
α1-2,3 Mannosidase	246
α1-2,3,6 Mannosidase	246
α2-3 Neuraminidase S.	247
α2-3,6,8 Neuraminidase	247
α2-3,6,8,9 Neuraminidase A.	247
β-N-Acetylglucosaminidase S	243
β-N-Acetylhexosaminidase,	244
β1-3 Galactosidase	245
β1-3,4 Galactosidase	246
β1-4 Galactosidase S	246
Exonucleases	
Common Applications	103
Exonuclease I (<i>E. coli</i>)	103
Exonuclease III (E. coli)	104
Exonuclease T	105, 202
Exonuclease V (RecBCD)	104
Exonuclease VII	104
Exonuclease VIII, truncated	104
Lambda Exonuclease	105
Nuclease BAL-31	106
Properties of	102
RecJ _f	106
T7 Exonuclease	107

F	
Factor Xa Protease	254
Fast DNA Ladder	167
Fatl	36
Faul	36
Fetuin	242
Fluorescent Substrates for Protein Labeling	
Fnu4HI	36
Fokl	37
Fpg	110
Freezer Program Locator, NEBNow	289
Freezer Program, On-site	4
Frequencies of Restriction Sites	
Fsel	37
Fspl	
FspEl	
Fucosidase	
Furin	

G	
Gel Extraction, Monarch Products for	126, 12
Gel Loading Dye, Blue (6X)	5
Gel Loading Dye, Orange (6X)	57
Gel Loading Dye, Purple (6X)	57, 16
Gel Loading Dye, Purple (6X), No SDS	57, 166
Gene Assembly	88-90
General Guidelines for PCR Optimization	33
General Molecular Biology Data	
Genetic Code	382
Genetic Markers	
Genomic DNA Substrates	
HeLa Genomic DNA	27-
Genome Editing	174-179
Gibson Assembly Cloning Kit	8
Gibson Assembly Master Mix	89

GIODAI BUSINESS Development	4
Globin & rRNA Depletion	148
GluC, Endoproteinase	253
Glycobiology and Protein Tools	236-257
Glucosylation (T4-BGT)	262
Glycosidases	
Endoglycosidases	238-242
Exoglycosidases	243-247
Glycosylases	
Afu Uracil-DNA Glycosylase (UDG)	111
Endonuclease III (Nth)	
Endonuclease IV	
Endonuclease V	110
Endonuclease VIII	110
Fpg	
hAAG	
T4 PDG	
Tma Endonuclease III	
Tth Endonuclease IV	
Uracil-DNA Glycosylase (UDG)	
Goat Anti-Mouse IgG Magnetic Beads	
Goat Anti-Rabbit IgG Magnetic Beads	
Goat Anti-Rat IgG Magnetic Beads	
Golden Gate Assembly Kit	
GpC Methyltransferase (M.CviPI)	114

П	
hAAG	111
Haell	37
HaeIII	37
HaellI Methyltransferase	115
Half Life of Radioisotopes	383
Heat Inactivation of Restriction Endonucleases	293-298
HeLa Genomic DNA	271
Helper Phage, M13K07	118
Hemo Klen <i>Taq</i>	67
Heparin Lyases	
Bacteriodes Heparinase I	249
Bacteriodes Heparinase II	
Bacteriodes Heparinase III	
HF Restriction Enzymes	
Hgal	
Hhal	
Hhal Methyltransferase	115
HiFi <i>Taq</i> DNA Ligase	97
High-Fidelity PCR	
Deep Vent DNA Polymerase	
One Taq DNA Polymerase	
One Taq Hot Start DNA Polymerase	
Q5 High-Fidelity DNA Polymerase	
Q5 Hot Start High-Fidelity DNA Polymerase	
Phusion High-Fidelity DNA Polymerase	
Phusion Hot Start Flex	
Vent DNA Polymerase	
High-Fidelity (HF) Restriction Enzymes	
HiScribe SP6 RNA Synthesis Kit	
HiScribe T7 ARCA mRNA Kit (with or without tailing)	
HiScribe T7 High Yield RNA Synthesis Kit	
HiScribe T7 Quick High Yield RNA Synthesis Kit	, ,
Hincll	
HindIII	
HindIII-HF	
Hinfl HinP1I	
Histones	
EpiMark Nucleosome Assembly Kit.	
H1º Human, Recombinant.	
H2A Human, Recombinant	
H2A/H2B Dimer	
H2B Human, Recombinant	
HZD HUHIAH, NEGUHIDHIAHL	200

Histone H3.1 Human, Recombinant	268
H3.1/H4 Tetramer	267
Histone H3.2 Human, Recombinant	269
Histone H3.3 Human, Recombinant	269
Histone H4 Human, Recombinant	269
Histone Methyltransferases	270
Hi-T7 RNA Polymerase	188
Homing Endonucleases	
I-Ceul	
I-Scel	
PI-Pspl	
PI-Scel.	
Hot Start PCR	
EpiMark Hot Start DNA Polymerase	
Hot Start <i>Tag</i> DNA Polymerase	
Hot Start Tag 2X Master Mix.	
LongAmp Hot Start DNA Polymerase	67
LongAmp Hot Start <i>Tag</i> 2X Master Mix	67
One Tag Hot Start DNA Polymerase	
One Tag Hot Start 2X Master Mixes	
Phusion Hot Start Flex DNA Polymerase	
Phusion Hot Start Flex 2X Master Mix	64
Q5 Hot Start High-Fidelity DNA Polymerase	
Q5 Hot Start High-Fidelity 2X Master Mix	
Hpal	
Hpall	
Hpall Methyltransferase.	116
Hphl	
Hpy99I	
Hpy166II	39
Hpy188I	39
Hpy188III	39
HpyAV	
HpyCH4III	
HpyCH4IV	39
HpyCH4V	39
hSMUG1	112
Human Alkyladenine Glycosylase (hAAG)	111
Human DNA (cytosine-5) MTase (Dnmt1)	
Human RNase Inhibitor	
Hydrophillic Streptavidin Magnetic Beads.	
Hydroxymethylcytosine, Identification	262-263

I-Ceul	55
Icon Descriptions	
IdeZ Protease (IgG-specific)	248
Immobilized Secondary Antibodies	219
Immunomagnetic Isolation	218
IMPACT Systems	
IMPACT Kit	212, 214
IMPACT Vectors	214
IMPACT Vectors Comparison	214
Inhibitors, RNA	202, 203
Inorganic Pyrophosphatase (E. coli)	189
Inorganic Pyrophosphatase, Thermostable	189
Inorganic Pyrophosphatase, (yeast)	189
Instant Sticky-end Ligase Master Mix	95
Inteins	214
Interactive Tools	
Internet Technical Support	288
In vitro Protein Synthesis, PURExpress	216
In vitro Transcription, Protein	216
I-Scel	
ISO 9001:2000, ISO 13485, ISO 14001	
IsoAmp II Universal tHDA Kit	
Isoschizomer List	
Isothermal Amplification	
Isotope Data	

K	
Kasl	39
Kinases	
cAMP-dependent Protein Kinase (PKA)	257
Casein Kinase II (CK2)	257
T4 Polynucleotide Kinase	98
T4 Polynucleotide Kinase (3´ phosphatase minus)	98
K. lactis Protein Expression	
K. lactis Competent Cells	215
K. lactis Protein Expression Kit	212, 215
K. lactis Vectors	215
Klenow Fragment	77
Klenow Fragment (3´→5´ exo⁻)	77
KLD Enzyme Mix	
Kpnl	40
Kpnl-HF	40

L	
Labeling with SNAP-tag/Troubleshooting	364
Ladders	
DNA	166-168, 170
RNA	170, 207
Lambda	
BstEII Digest	
DNA	
DNA (N ⁶ -methyladenine-free)	
Exonuclease (I, T, III, λ).	
HindIII Digest	
Map of	
Mono Cut Mix	
PFG Ladders	
Protein Phosphatase (Lambda PP)	
Lemo21(DE3) Competent E. coli	
Lesion Bypass Polymerase, <i>Sulfolobus</i> DNA Polymerase IV	
Library Preparation, DNA	
Adaptors and Primers	
DNA Enrichment, Microbiome Kit	
DNA Input Guidelines	
DNA Repair	
Fragmentation	
Illumina Compatible	
Ion Torrent Compatible	
Methylated – DNA	
Modules	
Target Enrichment	
Library Preparation, RNA	
Illumina Compatible	
Modules	
rRNA Depletion	
Library Quantitation.	
Licensing	4
Ligases 9°N DNA Ligase	00
E. coli DNA Ligase	
ElectroLigase	
HiFi <i>Tag</i> DNA Ligase	
NEBNext Ultra II Ligation Module	
RtcB Ligase	, , ,
SplintR Ligase	
T3 DNA Ligase	
T4 DNA Ligase	
T7 DNA Ligase	
T4 RNA Ligase 1	
T4 RNA Ligase 2	
T4 RNA Ligase 2, truncated	
T4 RNA Ligase 2, truncated K227Q	
T4 RNA Ligase 2, truncated KQ	
Tag DNA Ligase	
Thermostable 5´ App DNA/RNA	
momodable of hpp bringmin	130

Ligase Master Mixes	
Blunt/TA Ligase Master Mix	96
Instant Sticky-end Ligase Master Mix	95
Ligase Selection Chart, DNA	93, 94
Ligase Selection Chart, RNA	196
Ligation	
DNA	93-98
Protein	214
RNA	195-199
Ligation Kits	
Quick Ligation Kit	95
Loading Buffers	171
Loading Dye, Blue	57
Loading Dye, Orange	57
Loading Dye, Purple	
LongAmp Hot Start Taq 2X (Polymerase, Master Mix)	62, 67
LongAmp Taq (Polymerase, Master Mix, Kit)	62, 67
Long PCR	
LongAmp <i>Taq</i>	62, 67
LongAmp Hot Start <i>Taq</i>	62, 67
Phusion High-Fidelity DNA Polymerase	62, 64
Phusion Hot Start Flex DNA Polymerase	62, 64
Q5 High-Fidelity DNA Polymerase	62, 63
Q5 Hot Start High-Fidelity DNA Polymerase	62, 63
Low Range ssRNA Ladder	170, 207
Low Molecular Weight DNA Ladder	167
LpnPI	264
LunaScript RT SuperMix Kit	70
Luna Universal One-Step RT-qPCR Kit	71
Luna Universal Probe One-Step RT-qPCR Kit	71
Luna Universal Probe qPCR Master Mix	70
Luna Universal qPCR Master Mix.	70
Luna qPCR & RT-qPCR Optimization Tips	339-340
Luna qPCR & RT-qPCR Troubleshooting Guide	341-342

M	
5-methyl-dCTP	79
M13K07 Helper Phage	118
M13mp18 DNA	118
M13mp/pUC Series, Maps of	370
M13mp18, Map of	371
Magnetic Matrices (beads)	
Amylose	213
Anti-MBP	213
Chitin	219
Goat Anti-Mouse	219
Goat Anti-Rabbit	219
Goat Anti-Rat	219
Hydrophilic Streptavidin	218
Oligo (dT) ₂₅ Magnetic Beads	205
Protein A	218
Protein G	
SNAP-Capture	
Streptavidin	,
Magnetic mRNA Isolation Kit	
Magnetic Protein Purification	218, 219
Magnetic Separation Racks	
6-Tube Magnetic Separation Rack	
12-Tube Magnetic Separation Rack	
50 ml Magnetic Separation Rack	
96-Well Microtiter Plate Magnetic Separation Rack	219
Maps	
Lambda	
M13mp/pUC Series.	
M13mp18	
pBR322	
pKLAC2	
pMAL-p5X	
pMiniT 2.0	
pNEB206A	376

pSNAP _f	377
pTXB1	
pTYB21	
pUC19T7	
Markers and Ladders (DNA, RNA, Protein)	
Markers, Molecular Weight	
DNA	166-168, 170
RNA	,
Protein	
Mass Spectrometry Standards, Trypsin Digested BSA	252
Hot Start <i>Tag</i> 2X Master Mix	66
Long Amp Hot Start <i>Taq</i> 2X Master Mix	
LongAmp Taq 2X Master Mix	67
Luna Universal Probe qPCR Master Mix	
Luna Universal qPCR Master Mix	
NEBNext Ultra II Q5 Master Mix	
One Tag 2X Master Mix with GC Buffer	
One Taq 2X Master Mix with Standard Buffer	65
One Taq Hot Start 2X Master Mix with GC Buffer	
One Tag Hot Start 2X Master Mix with Standard Buffer	
One Taq Hot Start Quick-Load 2X Master Mix with GC Buffer One Taq Hot Start Quick-Load 2X Master Mix with Standard Buffer	
One Tag Quick-Load 2X Master Mix with GC Buffer	
One Tag Quick-Load 2X Master Mix with Standard Buffer	
Phusion High-Fidelity PCR Master Mix with GC Buffer	64
Phusion High-Fidelity PCR Master Mix with HF Buffer	
Phusion Hot Start Flex 2X Master Mix	
Q5 High-Fidelity 2X Master MixQ5 Hot Start High-Fidelity 2X Master Mix	
Quick-Load <i>Taq</i> 2X Master Mix	
Tag 2X Master Mix	
Taq 5X Master Mix	66
WarmStart Colorimetric LAMP Master Mix (DNA & RNA)	
Mbol	40
Mboll	40
Mboll	
Mboll McrBC Methyl-seq	265
McrBC	265 141
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG)	
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR	
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases	265 141 264-265 334-336 271
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR	265 141 264-265 334-336 271
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul	265 141 264-265 334-336 271 114
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.SssI) dam	265 141 264-265 334-336 271 114 114 114 115
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.Sssl) dam EcoGII	265 141 264-265 334-336 271 114 114 114 115 115
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.Sssl) dam EcoGII EcoRI	265 141 264-265 334-336 271 114 114 115 115 115
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.Sssl) dam EcoGII	265 141 264-265 334-336 271 114 114 115 115 115 1114
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.Sssl) dam EcoGII EcoRI GpC (M.CviPI)	265 141 264-265 334-336 271 114 114 114 115 115 115 115 1114
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.SssI) dam EcoGII EcoRI GpC (M.CviPI) HaelII Hhal	265 141 264-265 334-336 271 114 114 115 115 115 115 115 115 115 11
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.SssI) dam EcoCII EcoRI GpC (M.CviPI) HaeIII Hhal HpaII HpaII Human DNA (cytosine-5) MTase (Dnmt1)	265 141 264-265 334-336 271 114 114 115 115 115 115 115 115 115 11
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.SssI) dam EcoGII EcoRI GpC (M.CviPI) HaeIII Hhal HpaII Human DNA (cytosine-5) MTase (Dnmt1) MspI	265 141 264-265 334-336 271 114 114 115 115 115 115 116 270 116
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.SssI) dam EcoCII EcoRI GpC (M.CviPI) HaeIII Hhal HpaII HpaII Human DNA (cytosine-5) MTase (Dnmt1)	265 141 264-265 334-336 271 114 114 115 115 115 115 116 270 116 116
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.Sssl) dam EcoGII EcoRI GpC (M.CviPI) HaeIII Hhal HpaII Human DNA (cytosine-5) MTase (Dnmt1) MspI TaqI	265 141 264-265 334-336 271 114 114 115 115 115 115 116 116 270 116 40
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.Sssl) dam EcoGII EcoRI GpC (M.CviPI) HaeIII Hhal HpaII Human DNA (cytosine-5) MTase (Dnmt1) MspI TaqI Mfel Mfel-HF Microbiome Enrichment	265 141 264-265 334-336 271 114 114 115 115 115 115 116 116 40 40 152,159
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.SssI) dam EcoGII EcoRI GpC (M.CviPI) HaelII Hhal HpaII Human DNA (cytosine-5) MTase (Dnmt1) Mspl Taq1 Mfel-HF Microbiome Enrichment Micrococcal Nuclease	265 141 264-265 334-336 271 114 114 115 115 115 116 270 116 116 40 40 152,159 105
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.SssI) dam EcoGII EcoRI GpC (M.CviPI) HaelII Hhal HpaII Human DNA (cytosine-5) MTase (Dnmt1) MspI Taq1 Mfel Mfel Mfel Mfel Mfel Mfel Mferococcal Nuclease microRNA Marker	265 141 264-265 334-336 271 114 114 115 115 115 115 116 270 116 40 40 152,159 105 170,207
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.SssI) dam EcoGII EcoRI GpC (M.CviPI) HaelII Hhal HpaII Human DNA (cytosine-5) MTase (Dnmt1) Mspl Taq1 Mfel-HF Microbiome Enrichment Micrococcal Nuclease	265 141 264-265 334-336 271 114 114 115 115 115 115 116 270 116 40 40 152,159 105 170,207
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR. Methyltransferases Alul BamHI CpG (M.SssI) dam EcoGII EcoRI GpC (M.CviPI) HaeIII Hhal HpaII Human DNA (cytosine-5) MTase (Dnmt1) MspI TaqI Mfel Mfel-Mfel-HF Microbiome Enrichment Micrococcal Nuclease microRNA Marker Midrange PFG Marker	265 141 264-265 334-336 271 114 114 115 115 115 115 116 270 116 40 40 40 152,159 170,207 168 207
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.SssI) dam EcoGII EcoRI GpC (M.CviPI) HaeIII HhaI HhaI HpaII Human DNA (cytosine-5) MTase (Dnmt1) MspI TaqI Mfel-HF Microbiome Enrichment Micrococcal Nuclease microRNA Marker miRNA Cap 2'-O-Methyltransferase Mlul	265 141 264-265 334-336 271 114 114 115 115 115 115 116 270 116 40 40 152,159 170,207 168 207 190 40
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.SssI) dam EcoGII EcoRI GpC (M.CviPI) HaeIII HhaI HhaI HpaII Human DNA (cytosine-5) MTase (Dnmt1) Mspl Taq1 Mfel Mfel-HF Microbiome Enrichment Micrococcal Nuclease microRNA Marker miRNA Cap 2'-0-Methyltransferase Mlul Mlul-HF	265 141 264-265 334-336 271 114 114 115 115 115 115 116 116 270 116 40 40 152,159 105 170,207 190 40 40 40 40 40 40 40 40 40 40 40 40 40
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.SssI) dam EcoGII EcoRI GpC (M.CviPI) HaelII Hhal HpaII Human DNA (cytosine-5) MTase (Dnmt1) Mspl Taq1 Mfel-HF Microbiome Enrichment Micrococcal Nuclease microRNA Marker Midrange PFG Marker miRNA Cap 2'-O-Methyltransferase Miul Mlul-HF MiuCI	265 141 264-265 334-336 271 114 114 115 115 115 116 270 116 40 40 40 152,159 105 170,207 168 207 190 40 40 40 40 40 40 40 40 40 40 40 40 40
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.SssI) dam EcoGII EcoRI GpC (M.CviPI) HaelII Hhal HpaII Human DNA (cytosine-5) MTase (Dnmt1) MspI Taq1 Mfel Mfel Mfel Mfel-Mfers Micrococcal Nuclease microRNA Marker Midrange PFG Marker miRNA Cap 2'-O-Methyltransferase Mlul Miul-HF MluCI Miyl	265 141 264-265 334-336 271 114 114 115 115 115 115 116 270 116 40 40 152,159 105 170,207 168 207 190 40 40 40 40 41
McrBC Methyl-seq Methylation-Dependent Restriction Enzymes Methylation Sensitivity (dam, dcm, CpG) Methylation Specific PCR Methyltransferases Alul BamHI CpG (M.SssI) dam EcoGII EcoRI GpC (M.CviPI) HaelII Hhal HpaII Human DNA (cytosine-5) MTase (Dnmt1) Mspl Taq1 Mfel-HF Microbiome Enrichment Micrococcal Nuclease microRNA Marker Midrange PFG Marker miRNA Cap 2'-O-Methyltransferase Miul Mlul-HF MiuCI	265 141 264-265 334-336 271 114 114 115 115 115 115 116 270 116 116 40 40 152,159 170,207 168 207 190 40 40 40 41

Modifying Enzymes, Activity in Cutsmart Buffer		299
Monarch Collection Tubes II		
Monarch DNA Cleanup Binding Buffer		
Monarch DNA Cleanup Columns (5 μg)		
Monarch DNA Elution Buffer	126,	127
Monarch DNA Gel Extraction Kit		
Monarch DNA/RNA Protection Reagent		130
Monarch DNA Wash Buffer		
Monarch Gel Dissolving Buffer		126
Monarch gDNA Binding Buffer		128
Monarch gDNA Blood Lysis Buffer		128
Monarch gDNA Cell Lysis Buffer		
Monarch gDNA Elution Buffer		
Monarch gDNA Purification Columns		128
Monarch gDNA Removal Column		130
Monarch gDNA Tissue Lysis Buffer		128
Monarch gDNA Wash Buffer		
Monarch Microfuge Tube EcoRack		127
Monarch Nucleic Acid Purification Kits	.119, 122-	13
Monarch PCR & DNA Cleanup Kit (5 µg)		127
Monarch Plasmid Lysis Buffer (B2)		
Monarch Plasmid Miniprep Columns		
Monarch Plasmid Miniprep Kit		125
Monarch Plasmid Neutralization Buffer (B3)		125
Monarch Plasmid Resuspension Buffer (B1)		125
Monarch Plasmid Wash Buffer 1		125
Monarch Plasmid Wash Buffer 2		125
Monarch RNA Cleanup Binding Buffer		13
Monarch RNA Cleanup Columns (50 μg, 500 μg)		131
Monarch RNA Cleanup Kit (10 μg, 50 μg, 500 μg)		131
Monarch RNA Cleanup Wash Buffer		13
Monarch RNA Lysis Buffer		130
Monarch RNA Priming Buffer		130
Monarch RNA Purification Column		
Monarch Total RNA Miniprep Enzyme Pack		
Monarch Total RNA Miniprep Kit		
Monarch RNA Wash Buffer		
Monarch RNase A		
Monarch Troubleshooting Guide for DNA		
Monarch Troubleshooting Guide for gDNA.		
Monarch Troubleshooting Guide for RNA.		
Monarch gDNA Kit Sample Input Amounts		
Monarch RNA Kit Sample Input Amounts		
mRNA Isolation	,	
Mscl		
Msel		4
MsII		
Mspl		
Mspl Methyltransferase		
MspA11		
MspJI		264
Multiplex PCR		·
Hemo Klen Taq.		
Multiplex PCR 5X Master Mix		
One Tag DNA Polymerase		
Phusion High-Fidelity DNA Polymerase		
Q5 High-Fidelity DNA Polymerase		
Mung Bean Nuclease Murine RNase Inhibitor		
Mutation Detection		
Mwol	,	
WINVI		42

9°N DNA Ligase 98 N6-Methyladenosine Enrichment 204 Nael 42 Narl 42 Nb.BbvCl 53 Nb.Bsml 53 Nb.BsrDl 54	N	
V6-Methyladenosine Enrichment 204 Nael 42 Narl 42 Nb.BbvCl 53 Nb.Bsml 53	9°N DNA Ligase	98
Vlael 42 Vlarl 42 Vlb. BbvCl 53 Vlb. Bsml 53		
\text{Vb.BbvCl} 53 \text{Vb.Bsml} 53		
Vb.Bsml 53	Narl	42
	Nb.BbvCI	53
Nb.BsrDI	Nb.Bsml	53
	Nb.BsrDI	54

All, Des Ol			F.4	NEDNI I Colorio		000	, .
Nb.BssSI				NEBNext Selector			
Nb.Btsl				NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina			
Ncol				NEBNext Small RNA Library Prep Set for Illumina (Multiplex Compatible)			
Ncol-HF.				NEBNext Ultra Directional RNA Library Prep Kit for Illumina			
Ndel				NEBNext Ultra DNA Library Prep Kit for Illumina			
NEB 5-alpha (High Efficiency) Competent <i>E. coli</i>				NEBNext Ultra End Repair/dA-Tailing Module.			
NEB 5-alpha (Subcloning Efficiency) Competent <i>E. coli</i>				NEBNext Ultra Ligation Module			
NEB 5-alpha Electrocompetent <i>E. coli</i>				NEBNext Ultra RNA Library Prep Kit for Illumina			
NEB 5-alpha F [*] / _P (High Efficiency) Competent <i>E. coli</i>				NEBNext Ultra II Directional RNA Library Prep Kit for Illumina			
NEB 10-beta (High Efficiency) Competent <i>E. coli</i>				NEBNext Ultra II Directional RNA Library Prep	72, 170,	100,	200
NEB 10-beta Electrocompetent <i>E. coli</i>				with Sample Purification Beads	42 143 1	160	2Ω·
NEB 10-beta/Stable Outgrowth Medium				NEBNext Ultra II Directional RNA Second Strand Synthesis Module			
NEB Cloning Competent <i>E. coli</i> Sampler				NEBNext Ultra II DNA Library Prep Kit for Illumina.			
NEB Express (High Efficiency) Competent <i>E. coli</i> .				NEBNext Ultra II DNA Library Prep with Sample Purification Beads.			
NEB Express / (High Efficiency) Competent <i>E. coli</i>	,			NEBNext Ultra II End Repair/dA-Tailing Module			
NEB Golden Gate Assembly Kit (Bsal-HFv2)				NEBNext Ultra II FS DNA Library Prep Kit			
NEB Golden Gate Assembly Tool.				NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads			
NEB PCR Cloning Kit				NEBNext Ultra II FS DNA Module			
NEBpubs				NEBNext Ultra II Ligation Module			
				NEBNext Ultra II Non-directional RNA Second Strand Synthesis Module			
NEB Stable Competent <i>E. coli</i> NEB Tools				NEBNext Ultra II Non-directional Riva Second Strand Synthesis Module			
NEB Tube Opener				NEBNext Ultra II Q5 master mix. NEBNext Ultra II RNA First Strand Synthesis Module			
NEB Turbo Electrocompetent <i>E. coli</i> .				NEBNext Ultra II RNA First Strand Synthesis Module. NEBNext Ultra II RNA Library Prep Kit for Illumina			
•				NEBNext Ultra II RNA Library Prep Kit for Humina			
NEB Turbo (High Efficiency) Competent E. coli							
NEBaseChanger				NEBnow Locator NEBuffers for Restriction Enzymes			
NEBcloner							
NEBcutter				NEBuilder Assembly Tool			
NEBioCalculator				NEBuilder HiFi DNA Assembly Bundle for Large Fragments			
NEBNext Adaptors & Primers for Illumina				NEBuilder HiFi DNA Assembly Cloning Kit			
NEBNext Adaptor Dilution Buffer				NEBuilder HiFi DNA Assembly Master Mix			
NEBNext ChIP-Seq Library Prep Master Mix Set for Illumina				Neuraminidase			
NEBNext dA-Tailing Module				New England Biolabs Foundation			
NEBNext Direct Cancer Hot Spot Panel				N-Glycan Sequencing Kit			
NEBNext Custom Ready Panels				NgoMIV			
NEBNext Direct BRCA1/BRCA2 Panel				Nhel			
NEBNext dsDNA Fragmentase				Nhel-HF			4
NEBNext dsDNA Fragmentase Reaction Buffer v2				Nicking Endonucleases			
NEBNext DNA Library Prep Master Mix Set for Illumina				Nb.BbvCl			
NEBNext End Repair Module				Nb.Bsml			
NEBNext Enzymatic Methyl-seq Kit				Nb.BsrDI			
NEBNext Enzymatic Methyl-seq Conversion Module				Nb.BssSl			
NEBNext FFPE DNA Repair Mix				Nb.Btsl			
NEBNext Globin RNA & rRNA Depletion Kit				Nt.Alwl			
NEBNext Globin & rRNA Depletion Kit w/RNA Sample Purif. Beads				Nt.BbvCl			
NEBNext High-Fidelity 2X PCR Master Mix				Nt.BsmAl			
NEBNext Library Dilution Buffer				Nt.BspQl			
NEBNext Library Quant Kit for Illumina				Nt.BstNBI			
NEBNext Magnetic Separation Rack			.148	Nt.CviPII			5
NEBNext Microbiome DNA Enrichment Kit	1	52,	159	NiCo21(DE3) Competent E. coli.	17, 224, 2	225,	23
NEBNext Magnesium RNA Fragmentation Module				Ni-NTA Magnetic Beads			21
NEBNext Multiplex Oligos (Dual Index Primers Set 1)	146, 1	58,	160	NIaIII			4
NEBNext Multiplex Oligos (Dual Index Primers Set 2)				NIaIV			4
NEBNext Multiplex Oligos for Enzymatic Methyl-seq (Unique Dual Index Primer Pa				NmeAllI			
NEBNext Multiplex Oligos for Illumina (96 Index Primers)	,			Notl			
NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs)				NotI-HF			4
NEBNext Multiplex Oligos for Illumina (Index Primers Set 1, 2, 3, 4)				Nrul			
NEBNext Multiplex Oligos (Methylated Adaptor, Index Primers Set 1)				Nrul-HF			
NEBNext Multiplex Small RNA Library Prep Kit for Illumina (Index Primers 1-48)				Nsil			
VEBNext Multiplex Small RNA Library Prep Set for Illumina (Set 1)				Nsil-HF			
VEBNext Multiplex Small RNA Library Prep Set for Illumina (Set 2)				Nspl			
NEBNext Poly(A) mRNA Magnetic Isolation Module.				Nt.Alwl			
VEBNext Q5U Master Mix				Nt.BbvCl			
VEBNext Quick Ligation Module				Nt.BsmAl			
NEBNext Reagents for Next Generation Sequencing			. 55	Nt.BspQI			
Illumina Compatible	-156 1	58	160	Nt.BstNBI			
Ion Torrent.				Nt.CviPII			
Selection Charts				Nucleases			
NEBNext RNA First Strand Synthesis Module				BAL-31 Nuclease			10
NEBNext rRNA Depletion Kit (Human/Mouse/Rat)	, 10, I	50,	101	Cas9			
(with or without RNA Sample Purification Beads)	117 1	en	161	DNase I (RNase-Free)			
	147						111

Endonuclease VIII	
Exonuclease I (E. coli)	
Exonuclease III (E. coli)	
Exonucleases and Endonucleases, Properties of	
Exonuclease T	
Exonuclease V	
Exonuclease VII	
Exonuclease VIII, truncated	
Homing	
Lambda Exonuclease	
Micrococcal	
Mung Bean Nuclease Nicking Endonucleases	
Nuclease P1	
RecJ.	
Restriction Endonucleases	
T5 Exonuclease	
T7 Exonuclease	
Nuclease-free Water	
Nucleic Acid Data	
Nucleic Acid Purification	
Nucleic Acid Purification Troubleshooting Guide for DNA	
Nucleic Acid Purification Troubleshooting Guide for gDNA	
Nucleic Acid Purification Troubleshooting Guide for RNA	
Nucleosome Control DNA	
Nucleoside Digestion Mix.	
Vucleotides	
7-deaza-dGTP	79
ATP	79
Acyclonucleotide Set	79
Deoxynucleotide (dNTP) Solution Set	79
Deoxynucleotide (dNTP) Solution Mix	
5-methyl-dCTP	79
Physical Properties.	381
Ribonucleotide Solution Mix Ribonucleotide Solution Set	79
Ribonucleotide Solution Mix	79
Ribonucleotide Solution MixRibonucleotide Solution Set	
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Decan Genome Legacy Center	
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Decan Genome Legacy Center DEM	
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase	
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase Digos, NEBNext (Multiplex)	79 79 79 6 4 242 146
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center OEM O-Glycosidase Digos, NEBNext (Multiplex) Digo d(T) ₂₆ Cellulose Beads	
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Decan Genome Legacy Center DEM D-Glycosidase Digos, NEBNext (Multiplex) Digo d(T) 25 Cellulose Beads Digo d(T) 25 Magnetic Beads	
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Decan Genome Legacy Center DEM D-Glycosidase Diligos, NEBNext (Multiplex) Diligo d(T) 25 Cellulose Beads Diligo d(T) 25 Magnetic Beads Diligo d(T) 23 VN Dnline Tools	
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Decan Genome Legacy Center DEM D-Glycosidase Diligos, NEBNext (Multiplex) Diligo d(T) 25 Cellulose Beads Diligo d(T) 25 Magnetic Beads Diligo d(T) 23 VN Dnline Tools	
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Decan Genome Legacy Center DEM D-Glycosidase Diligos, NEBNext (Multiplex) Diligo d(T) 25 Cellulose Beads Diligo d(T) 25 Magnetic Beads Diligo d(T) 23 VN Dnline Tools Dne Taq DNA Polymerase (Enzyme, Master Mix, Kit)	79 79 79 66 4 242 146 205 205 205 288-289 62, 65
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Dean Genome Legacy Center DEM D-Glycosidase Diligos, NEBNext (Multiplex) Diligo d(T) 25 Cellulose Beads Diligo d(T) 25 Magnetic Beads Diligo d(T) 27 WN Dnline Tools Dne Taq DNA Polymerase (Enzyme, Master Mix, Kit) Dne Taq Hot Start DNA Polymerase (Enzyme, Master Mix)	79 79 79 66 4 242 146 205 205 205 205 225 265 62,65
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Dean Genome Legacy Center DEM D-Glycosidase Diligos, NEBNext (Multiplex) Diligo d(T) 25 Cellulose Beads Diligo d(T) 25 Magnetic Beads Diligo d(T) 23 VN Dnline Tools Dne Taq DNA Polymerase (Enzyme, Master Mix, Kit) Dne Taq Hot Start DNA Polymerase (Enzyme, Master Mix) Dne Taq One-Step RT-PCR Kit	79 79 79 66 4 242 146 205 205 193 288-289 62, 65 62, 65
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase Digos, NEBNext (Multiplex) Digo d(T) 25 Cellulose Beads Digo d(T) 25 Magnetic Beads Digo d(T) 27 VN One Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq One-Step RT-PCR Kit One Taq One-Step RT-PCR Kit	79 79 79 66 4 242 146 205 205 193 288-289 62, 65 62, 65
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase Oligos, NEBNext (Multiplex) Oligo d(T) 25 Cellulose Beads Oligo d(T) 25 Magnetic Beads Oligo d(T) 27 VN Online Tools One Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq One-Step RT-PCR Kit One Taq One-Step RT-PCR Kit	75 75 75 24 244 146 205 205 193 288-286 62, 65 62, 65 77
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase Oligos, NEBNext (Multiplex) Oligo d(T) 25 Cellulose Beads Oligo d(T) 25 Magnetic Beads Oligo d(T) 27 VN Online Tools One Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq One-Step RT-PCR Kit One Taq One-Step RT-PCR Kit One Taq RT-PCR Kit One Taq RT-PCR Kit Optimization Tips	75 75 75 24 244 245 205 205 288-285 62, 65 62, 65 77 72
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase Dligos, NEBNext (Multiplex) Digo d(T) 25 Cellulose Beads Dligo d(T) 25 Magnetic Beads Digo d(T) 23 VN Inline Tools One Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq Hot Start DNA Polymerase (Enzyme, Master Mix) Dne Taq One-Step RT-PCR Kit One Taq RT-PCR Kit One Taq RT-PCR Kit Optimization Tips Cloning.	75 75 75 75 75 75 76 77 78 78 79 79 79 70 70 70 70 70 70 70 70 70 70
Ribonucleotide Solution Mix. Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase Digos, NEBNext (Multiplex) Digo d(T) 25 Cellulose Beads Digo d(T) 25 Magnetic Beads Digo d(T) 23 VN Dnline Tools One Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq One-Step RT-PCR Kit Done Taq RT-PCR Kit Dotimization Tips Cloning Luna qPCR	79 79 79 79 79 79 242 144 205 205 208 288-289 62, 65 62, 65 72 72 73 352-353 339 340
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase Digos, NEBNext (Multiplex) Digos (T) 25 Cellulose Beads Digo d(T) 25 Cellulose Beads Digo d(T) 23 VN Dolline Tools Dine Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq Hot Start DNA Polymerase (Enzyme, Master Mix) Dine Taq RT-PCR Kit Diptimization Tips Cloning Luna qPCR Luna RT-qPCR	75 75 75 75 75 75 75 75 75 75
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase Oligos, NEBNext (Multiplex) Oligo d(T) 25 Cellulose Beads Oligo d(T) 25 Magnetic Beads Oligo d(T) 28 Magnetic Beads Oligo d(T) 28 VN Online Tools One Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq One-Step RT-PCR Kit Optimization Tips Cloning Luna qPCR Luna RT-qPCR PCR Restriction Endonucleases	79 79 79 79 79 79 79 79 79 79 79 79 79 7
Ribonucleotide Solution Mix Ribonucleotide Solution Set Ocean Genome Legacy Center DEM O-Glycosidase Digos, NEBNext (Multiplex) Digo d(T) 25 Cellulose Beads Digo d(T) 25 Magnetic Beads Digo d(T) 27 VN Diline Tools One Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq One-Step RT-PCR Kit One Taq RT-PCR Kit Optimization Tips Cloning Luna qPCR Luna RT-qPCR PCR Restriction Endonucleases Orange Gel Loading Dye	75 75 75 75 75 75 76 76 76 77 77 77 77 352-353 339 340 337 340 352-353
Ribonucleotide Solution Mix Ribonucleotide Solution Set Ocean Genome Legacy Center DEM D-Glycosidase Digos, NEBNext (Multiplex) Digo d(T) 25 Cellulose Beads Digo d(T) 25 Magnetic Beads Digo d(T) 27 WN Donline Tools Done Taq DNA Polymerase (Enzyme, Master Mix, Kit) Dne Taq Hot Start DNA Polymerase (Enzyme, Master Mix) Dne Taq One-Step RT-PCR Kit Dote Taq Cone-Step RT-PCR Kit Dote Taq Cone-Step RT-PCR Kit Dote Taq RT-PCR Kit Dote Taq RT-PCR Kit Dote Taq RT-PCR Luna RT-qPCR Luna RT-qPCR PCR Restriction Endonucleases Drange Gel Loading Dye	79 79 79 79 79 79 79 79 79 79 79 79 79 7
Ribonucleotide Solution Mix Ribonucleotide Solution Set Ocean Genome Legacy Center DEM O-Glycosidase Digos, NEBNext (Multiplex) Digo d(T) 25 Cellulose Beads Digo d(T) 25 Magnetic Beads Digo d(T) 27 VN Diline Tools One Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq One-Step RT-PCR Kit One Taq RT-PCR Kit Optimization Tips Cloning Luna qPCR Luna RT-qPCR PCR Restriction Endonucleases Orange Gel Loading Dye	75 75 75 75 75 75 76 76 76 77 77 77 77 352-353 339 340 337 340 352-353
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase Dligos, NEBNext (Multiplex) Dligo d(T) 25 Cellulose Beads Dligo d(T) 25 Magnetic Beads Dligo d(T) 27 WN Done Taq DNA Polymerase (Enzyme, Master Mix, Kit) Done Taq One-Step RT-PCR Kit Done Taq One-Step RT-PCR Kit Dit mization Tips Cloning Luna qPCR Luna RT-qPCR PCR Restriction Endonucleases Drange Gel Loading Dye Dredering Information	75 75 75 75 75 24 24 24 20 20 28 28 62, 66 72 72 352-353 334 344 337 290
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase Dligos, NEBNext (Multiplex) Dligo d(T) 25 Cellulose Beads Dligo d(T) 25 Cellulose Beads Dligo d(T) 27 Magnetic Beads Dligo d(T) 27 VN Dnline Tools One Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq DNA Polymerase (Enzyme, Master Mix) Dne Taq Pto Start DNA Polymerase (Enzyme, Master Mix) Dne Taq Tar-PCR Kit Dne Taq RT-PCR Kit Dptimization Tips Cloning Luna qPCR Luna RT-qPCR PCR Restriction Endonucleases Orange Gel Loading Dye Drdering Information	75 75 75 75 75 75 75 75 75 75
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase Diigos, NEBNext (Multiplex) Digo d(T) 25 Cellulose Beads Diigo d(T) 25 Magnetic Beads Diigo d(T) 23 VN Online Tools One Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq DNA Polymerase (Enzyme, Master Mix) Dne Taq Pone-Step RT-PCR Kit Dote Taq RT-PCR Kit Dotimization Tips Cloning Luna qPCR Luna RT-qPCR PCR Restriction Endonucleases Orange Gel Loading Dye Ordering Information Pacl	79 79 79 79 79 79 79 79 79 79 79 79 79 7
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase Digos, NEBNext (Multiplex) Digo d(T) 25 Cellulose Beads Digo d(T) 26 Magnetic Beads Digo d(T) 27 WN Online Tools One Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq Hot Start DNA Polymerase (Enzyme, Master Mix) One Taq One-Step RT-PCR Kit Optimization Tips Cloning Luna qPCR Luna RT-qPCR PCR Restriction Endonucleases Orange Gel Loading Dye Ordering Information Pacl PagRNA Binding Protein Pacl PagRNI	79 79 79 79 79 79 79 79 79 79 79 79 79 7
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase Digos, NEBNext (Multiplex) Digo d(T) 25 Cellulose Beads Digo d(T) 25 Magnetic Beads Dilgo d(T) 27 WA Online Tools One Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq Hot Start DNA Polymerase (Enzyme, Master Mix) One Taq One-Step RT-PCR Kit Optimization Tips Cloning Luna qPCR Luna RT-qPCR PCR Restriction Endonucleases Orange Gel Loading Dye Ordering Information Pacl PageR71 Partnering with NEB	79 79 79 79 79 79 79 79 79 79 79 79 79 7
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center DEM O-Glycosidase Digos, NEBNext (Multiplex) Digo d(T) 25 Cellulose Beads Digo d(T) 25 Magnetic Beads Dilgo d(T) 27 WA Online Tools One Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq Hot Start DNA Polymerase (Enzyme, Master Mix) One Taq One-Step RT-PCR Kit Optimization Tips Cloning Luna qPCR Luna RT-qPCR PCR Restriction Endonucleases Orange Gel Loading Dye Ordering Information Pacl PageR71 Partnering with NEB	79 79 79 79 79 79 79 79 79 79 79 79 79 7
Ribonucleotide Solution Mix Ribonucleotide Solution Set O Ocean Genome Legacy Center OEM O-Glycosidase Oligos, NEBNext (Multiplex) Oligo d(T) 25 Cellulose Beads Oligo d(T) 25 Magnetic Beads Oligo d(T) 25 Magnetic Beads Oligo d(T) 27 WN Online Tools One Taq DNA Polymerase (Enzyme, Master Mix, Kit) One Taq DNA Polymerase (Enzyme, Master Mix) One Taq One-Step RT-PCR Kit Optimization Tips Cloning Luna qPCR Luna RT-qPCR PCR Restriction Endonucleases Orange Gel Loading Dye Ordering Information P pp 19 siRNA Binding Protein Pacl Pacl Pacl Pacl Pacl Partnering with NEB pBR322	79 79 79 79 79 79 79 79 79 79 79 79 79 7

Mspl Digest	
PcilpCLIP, Vector	
PCR	202
Cleanup	81, 127
Cloning Kit	
Deep Vent DNA Polymerase	68
Deep Vent (exo-) DNA Polymerase.	68
Deoxynucleotides	
EpiMark Hot Start <i>Taq</i>	
Guidelines for Optimization	
Hemo Klen Taq	
Hot Start <i>Tag</i> DNA Polymerase	
LongAmp Hot Start <i>Tag</i> 2X Master Mix.	
LongAmp Hot Start <i>Taq</i> DNA Polymerase	
LongAmp <i>Taq</i> DNA Polymerase	
LongAmp <i>Taq</i> 2X Master Mix	
LongAmp <i>Taq</i> PCR Kit	
Marker	167
Multiplex PCR 5X Master Mix	
One Taq 2X Master Mix with Standard Buffer	
One Tag DNA Polymerase	
One Tag Hot Start 2X Master Mix with GC Buffer	
One Taq Hot Start 2X Master Mix with Standard Buffer One Taq Hot Start DNA Polymerase	
One <i>Tag</i> Hot Start Quick-Load 2X Master Mix with GC Buffer	
One Tag Hot Start Quick-Load 2X Master Mix with Standard Buffer	
One Tag Quick-Load 2X Master Mix with Standard Buffer	
PCR Buffers, Activity of Restriction Enzymes in	344-345
PCR Optimization, Guidelines	337
Phusion High-Fidelity DNA Polymerase	
Phusion High-Fidelity PCR Kit	
Phusion High-Fidelity PCR Master Mix with GC Buffer	
Phusion High-Fidelity PCR Master Mix with HF Buffer	
Phusion Hot Start Flex 2X Master Mix Phusion Hot Start Flex DNA Polymerase	
Polymerase Selection Chart	
PreCR Repair Mix	
Q5 High-Fidelity 2X Master Mix	
Q5 High-Fidelity DNA Polymerase	
Q5 High-Fidelity PCR Kit	63
Q5 Hot Start High-Fidelity 2X Master Mix	
Q5 Hot Start High-Fidelity DNA Polymerase	
qPCR	
Quick-Load <i>Taq</i> 2X Master Mix	
RT-qPCR Selection Chart	
Selection Chart	
Tag 5X Master Mix	
Taq DNA Polymerase with Standard Taq Buffer	
Taq DNA Polymerase with Standard Taq (Mg-free) Buffer	
Taq DNA Polymerase with ThermoPol Buffer	
Taq PCR Kit	
Troubleshooting Guide	338
Vent DNA Polymerase	
Vent (exo ⁻) DNA Polymerase	68
PCR from Blood	07
Hemo Klen <i>Tag.</i> PCR Selection Tool	
Per Selection 1001 Peptide Display Cloning System, Ph.D.	
Peptide Ligation	
Performance Chart for Restriction Enzymes	
PfIFI	
PfIMI	
Phage Display Libraries, Ph.D. Kits	
Ph.D. Peptide Display Cloning System	
Ph.D7 Phage Display Peptide Library Kit	
Ph.D12 Phage Display Peptide Library	
Ph.D12 Phage Display Peptide Library Kit	
Ph.DC7C Phage Display Peptide Library Kitphi29 DNA Polymerase	
אווע באווע ו חואיום באוום אווע האווע ה	

φΧ174	
HaellI Digest	
RF I, RF II DNA	
Virion DNA Phosphatases	118
Phosphatase, Antarctic	110
Phosphatase, Calf Intestinal (CIP)	
Phosphatase, Shrimp Alkaline (rSAP)	
Quick Dephosphorylation Kit	
Phosphatases, Protein	
Lambda Protein Phosphatase (λ-PPase)	
Phusion High-Fidelity DNA Polymerase (Enzyme, Master Mix, Kit)	
Phusion Hot Start Flex DNA Polymerase (Enzyme, Master Mix)	
pKLAC2, Map of	
PI-PspI PI-Scel	
Plasmid Miniprep, Monarch Products for	
Plel	
PluTI	
pMAL-c5X Vector	
pMAL-p5X, Map of	374
pMAL-p5X Vector	
pMAL Protein Fusion and Purification System	
pMiniT 2.0, Map of	
Pmel	
PmII	
p-Nitrophenylphosphate (PNPP)pNEB206A, Map of	
PNGase A	
PNGase F	
PNGase F (Glycerol Free, Recombinant, Rapid, non-reducing)	
PNGase F, Remove-iT	
Polbase	289
Poly(A) Polymerase, E. coli	188
polyA Spin mRNA Isolation Kit	
Poly(U) Polymerase	188
Polymerases, DNA Activity of Postriction Engagement in BCR Buffers	244 245
Activity of Restriction Enzymes in PCR Buffers	
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase.	73
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase	
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase.	
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase Bst DNA Polymerase (Large Fragment, Full Length)	
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase. Bst 3.0 DNA Polymerase. Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase. Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases.	
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase. Bst 3.0 DNA Polymerase. Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase. Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases. DNA Polymerase I (E. coli)	
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase (E. coli) DNA Polymerase I, Large (Klenow) Fragment	
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase. Bst 3.0 DNA Polymerase. Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase. Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases. DNA Polymerase I (E. colir) DNA Polymerase I, Large (Klenow) Fragment. Epimark Hot Start Taq DNA Polymerase.	73 73 73 73 78 68 77 77
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase. Bst 3.0 DNA Polymerase (Large Fragment, Full Length) Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase. Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases. DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase. Guidelines for Optimization.	73 73 73 73 78 68 77 77 68 68
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase. Bst 3.0 DNA Polymerase. Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase. Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases. DNA Polymerase I (E. coli') DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization. Hemo KlenTaq	73 73 73 73 78 68 77 77 77 68 337 62,67
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq	73 73 73 73 78 68 77 77 68 337 62, 67 62, 66
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Klenow Fragment (3'→5' exo⁻)	73 73 73 73 78 68 77 77 68 337 62,67 62,66
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq	73 73 73 73 73 78 68 68 77 77 68 337 62,67 62,66
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq Klenow Fragment (3´→5´ exo⁻) LongAmp Taq DNA Polymerases	73 73 73 73 73 78 78 68 68 337 62,67 62,66 77 62,66 62,67
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq Klenow Fragment (3'→5' exo⁻) LongAmp Taq DNA Polymerases One Taq DNA Polymerases	73 73 73 73 73 73 78 68 77 77 68 337 62,67 62,66 77 62,65 62,65 62,67 52 62,67
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I (E coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq Klenow Fragment (3'→5' exo⁻) LongAmp Taq DNA Polymerases One Taq DNA Polymerases phi29 DNA Polymerase Phusion High-Fidelity DNA Polymerases QS High-Fidelity DNA Polymerases	73 73 73 73 73 78 68 77 77 68 337 62,67 62,66 62,66 62,66 62,64 62,63
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase. Bst 3.0 DNA Polymerase. Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase. Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases. DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase. Guidelines for Optimization. Hemo KlenTaq. Klenow Fragment (3'→5' exo⁻) LongAmp Taq DNA Polymerases. One Taq DNA Polymerase. phi29 DNA Polymerase. phusion High-Fidelity DNA Polymerases. Q5 High-Fidelity DNA Polymerases. Reaction Buffers.	73 73 73 73 73 73 78 68 77 77 77 68 3337 62,67 62,66 77 62,67 62,65 62,64 62,63 79
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase. Bst 3.0 DNA Polymerase. Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase. Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases. DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization. Hemo KlenTaq. Klenow Fragment (3'→5' exo⁻) LongAmp Taq DNA Polymerases One Taq DNA Polymerases phi29 DNA Polymerase phusion High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Reaction Buffers. Selection Chart	73 73 73 73 73 78 68 77 77 77 68 3337 62,67 62,66 77 62,65 75 62,64 62,63 79 62 62
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase. Bst 3.0 DNA Polymerase. Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase. Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases. DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment. Epimark Hot Start Taq DNA Polymerase. Guidelines for Optimization. Hemo KlenTaq. Hot Start Taq. Klenow Fragment (3'→5' exo⁻) LongAmp Taq DNA Polymerases. One Taq DNA Polymerases. phi29 DNA Polymerase. Phusion High-Fidelity DNA Polymerases. Q5 High-Fidelity DNA Polymerases. Reaction Buffers. Selection Chart. Sulfolobus DNA Polymerase IV.	73 73 73 73 73 78 68 77 77 77 68 68 77 62,67 62,67 62,67 62,67 62,67 77 62,67 62,67 77 62,67 62,67 77 62,67 62,67 77 62,67 62,67 62,67 62,67 62,67 62,67 62,67 62,67 62,67 62,67 62,67 62,67
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq Klenow Fragment (3'→5' exo⁻) LongAmp Taq DNA Polymerases One Taq DNA Polymerase Phusion High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Reaction Buffers Selection Chart Sulfolobus DNA Polymerase IV. T4 DNA Polymerase	73 73 73 73 73 78 78 68 68 77 77 77 68 3337 62,67 62,66 77 62,65 62,65 62,64 62,63 75 62,64 76 76 76
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq Klenow Fragment (3'→5' exo⁻) LongAmp Taq DNA Polymerases One Taq DNA Polymerases phi29 DNA Polymerase Phusion High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Reaction Buffers Selection Chart Sulfolobus DNA Polymerase IV T4 DNA Polymerase T7 DNA Polymerase T7 DNA Polymerase	73 73 73 73 73 78 78 68 68 77 77 68 337 62,67 62,66 77 62,65 62,64 62,63 79 62,67 77 78
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq Klenow Fragment (3'→5' exo⁻) LongAmp Taq DNA Polymerases One Taq DNA Polymerase Phusion High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Reaction Buffers Selection Chart Sulfolobus DNA Polymerase IV. T4 DNA Polymerase	73 73 73 73 73 73 78 78 68 77 77 68 337 62,67 62,66 77 62,66 62,63 62,63 79 62,67 77 78
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq Klenow Fragment (3'→5' exo') LongAmp Taq DNA Polymerases One Taq DNA Polymerases phi29 DNA Polymerase Phusion High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Reaction Buffers Selection Chart Sulfolobus DNA Polymerase IV T4 DNA Polymerase T7 DNA Polymerase Taq DNA Polymerase Taq DNA Polymerase Taq DNA Polymerase	73 73 73 73 73 73 78 68 77 77 68 337 62,67 62,66 77 62,66 62,65 62,64 62,63 79 62 62 76 77 78 62 62 76 77 78
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerases (Le coli) DNA Polymerases I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq Klenow Fragment (3'→5' exo⁻) LongAmp Taq DNA Polymerases One Taq DNA Polymerases phi29 DNA Polymerase Phusion High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Reaction Buffers Selection Chart Sulfolobus DNA Polymerase IV T4 DNA Polymerase Taq DNA Polymerase Tag DNA Polymerase	73 73 73 73 73 73 78 68 77 77 68 337 62,67 62,66 62,66 62,66 62,63 79 62 76 77 78 62 76 77 78 62 76 77 78 78
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq Klenow Fragment (3'→5' exo⁻) LongAmp Taq DNA Polymerases One Taq DNA Polymerases Phusion High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Reaction Buffers Selection Chart Sulfolobus DNA Polymerase T7 DNA Polymerase Taq DNA Polymerase Terminal Transferase Therminator DNA Polymerase	73 73 73 73 73 73 78 68 68 77 77 62,67 62,66 62,66 62,63 79 62,63 79 62 76 77 78 62,67 77 78 62,67
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase Bst DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I (E coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Klenow Fragment (3'→5' exo⁻) LongAmp Taq DNA Polymerases One Taq DNA Polymerases phi29 DNA Polymerase Phusion High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Reaction Buffers Selection Chart Sulfolobus DNA Polymerase IV T4 DNA Polymerase T7 DNA Polymerase Taq DNA Polymerase Taq DNA Polymerase Tag DNA Polymerase	73 73 73 73 73 73 78 68 68 77 77 62,67 62,66 62,66 62,63 79 62,63 79 62 76 77 78 62,67 77 78 62,67
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq DNA Polymerases One Taq DNA Polymerases One Taq DNA Polymerases One Taq DNA Polymerases Phusion High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Reaction Buffers Selection Chart Sulfolobus DNA Polymerase T7 DNA Polymerase T7 DNA Polymerase Taq DNA Polymerase Terminal Transferase Therminator DNA Polymerase Troubleshooting guide Vent DNA Polymerase Polymerases, RNA E. coli Poly(A) Polymerase	73 73 73 73 73 73 78 68 68 77 77 68 3337 62,67 62,66 77 62,65 62,64 62,63 79 62,67 77 78 62,66 77 78 62,66 62,63
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq Klenow Fragment (3'→5' exo⁻) LongAmp Taq DNA Polymerases One Taq DNA Polymerases phi29 DNA Polymerase Phusion High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Reaction Buffers Selection Chart Sulfolobus DNA Polymerase IV T4 DNA Polymerase T7 DNA Polymerase Taq DNA Polymerase Terminal Transferase Therminator DNA Polymerase Troubleshooting guide Vent DNA Polymerase Polymerases, RNA E. coli Poly(A) Polymerase Poly(U) Polymerase	73 73 73 73 73 73 73 78 78 68 77 77 68 3337 62,67 62,66 77 62,66 62,63 79 62,64 62,63 79 62,67 77 78 62,66 77 78 62,66 78 76 62,66 78 8338 68
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I, Large (Klenow) Fragment Deep Vent DNA Polymerases DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq Klenow Fragment (3'→5' exo') LongAmp Taq DNA Polymerases One Taq DNA Polymerases One Taq DNA Polymerase Phusion High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Reaction Buffers Selection Chart Sulfolobus DNA Polymerase IV T4 DNA Polymerase T7 DNA Polymerase Taq DNA Polymerase Terminal Transferase Therminator DNA Polymerase Troubleshooting guide Vent DNA Polymerase Polymerases, RNA E coli Poly(A) Polymerase Polymerases Polymerases, RNA E coli Poly(A) Polymerase Polymerase SP6 RNA Polymerase	73 73 73 73 73 73 78 86 68 77 77 68 337 62,67 62,66 77 62,66 62,63 79 62,64 62,63 79 62,67 77 78 62,66 77 78 62,66 77 78 62,66 77 78 62,66 78 78 62,66 78 78 82,86
Activity of Restriction Enzymes in PCR Buffers Bst 2.0 DNA Polymerase Bst 3.0 DNA Polymerase (Large Fragment, Full Length) Bst 2.0 WarmStart DNA Polymerase Bsu DNA Polymerase Large Fragment Deep Vent DNA Polymerases DNA Polymerase I (E. coli) DNA Polymerase I, Large (Klenow) Fragment Epimark Hot Start Taq DNA Polymerase Guidelines for Optimization Hemo KlenTaq Hot Start Taq Klenow Fragment (3'→5' exo⁻) LongAmp Taq DNA Polymerases One Taq DNA Polymerases phi29 DNA Polymerase Phusion High-Fidelity DNA Polymerases Q5 High-Fidelity DNA Polymerases Reaction Buffers Selection Chart Sulfolobus DNA Polymerase IV T4 DNA Polymerase T7 DNA Polymerase Taq DNA Polymerase Terminal Transferase Therminator DNA Polymerase Troubleshooting guide Vent DNA Polymerase Polymerases, RNA E. coli Poly(A) Polymerase Poly(U) Polymerase	73 73 73 73 73 73 73 78 68 68 77 77 68 3337 62,67 62,66 77 62,66 62,65 62,64 62,63 79 62 62 76 62 63 64 62 63 68 78 76 68 78 76 68 78 76 88 88 88 88 88

Polymerases and Amplification Technologies	
PpuMI	45
PreCR Repair Mix	76
Primers for cDNA Synthesis	
Oligo d(T)18 mRNA Primer	193
Oligo d(T)23 VN	
Random Primer 6	
Random Primer 9	
Random Primer Mix	
Properties of Exonucleases and Endonucleases	
Properties of DNA Repair Enzymes	107
Properties of Protein Phosphatases	256
Proteases	248, 254, 255
Protein A Magnetic Beads	
Protein Data	
Protein Deglycosylation Mix II	
Protein Expression & Purification Selection Chart	
·	
Protein Expression & Purification Technologies	
Protein Expression with T7 Express Strains	
Protein Fusion and Purification System	
Protein G Magnetic Beads	218
Protein Kinases	257
Proteinase K	254
Protein Ligation	214
Protein Markers and Ladders	
Blue Prestained Protein Std., Broad range (10–200 kDa)	
Color Prestained Protein Std., Broad Range (10–250 kDa)	
Unstained Protein Std., Broad Range (10–200 kDa)	
,	
Protein Phosphatases	
Protein Standards	
Proteome Analysis	
ProtoScript II First Strand cDNA Synthesis Kit	
ProtoScript II Reverse Transcriptase	20 101 102
Trotoochpt in neverse manscriptase	00, 131, 132
ProtoScript First Strand cDNA Synthesis Kit. PshAl	80, 191, 194
ProtoScript First Strand cDNA Synthesis Kit. PshAl	80, 191, 194 45
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil	80, 191, 194 45 46
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP _f Vector	80, 191, 194 45 46 282
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP _f Vector. pSNAP _f Vector, Map of	80, 191, 194 45 46 282 377
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP _T Vector pSNAP _T Vector, Map of pSNAP-tag(T7)-2 Vector	80, 191, 194 45 46 282 377 282
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP _T Vector pSNAP _T Vector, Map of pSNAP-tag(T7)-2 Vector PspGI	
ProtoScript First Strand cDNA Synthesis Kit. PshAI PsiI pSNAP _t Vector. pSNAP _t Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI	
ProtoScript First Strand cDNA Synthesis Kit. PshAI PsiI pSNAP _t Vector pSNAP _t Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI	
ProtoScript First Strand cDNA Synthesis Kit. PshAI PsiI pSNAP _t Vector. pSNAP _t Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI	
ProtoScript First Strand cDNA Synthesis Kit. PshAI PsiI pSNAP _t Vector pSNAP _t Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI	
ProtoScript First Strand cDNA Synthesis Kit. PshAI PsiI pSNAP _t Vector pSNAP _t Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI PstI	
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP _f Vector pSNAP _f Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI PstI PstI-HF pTWIN1 Vector	
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP _f Vector. pSNAP _f Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI PstI PstI-HF pTWIN1 Vector, Map of	
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP _f Vector. pSNAP _f Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI PstI PstI-HF pTWIN1 Vector. pTXB1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector	
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP, Vector. pSNAP, Vector, Map of pSNAP-tag(T7)-2 Vector PspGl PspOMI PspXI PstI PstI-HF pTWIN1 Vector. pTXB1 Vector, Map of. pTXB1 Vector pTYB21 Vector	
ProtoScript First Strand cDNA Synthesis Kit. PshAI PsiI pSNAP ₁ Vector pSNAP ₁ Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI PstI PstI PstI-HF pTWIN1 Vector pTXB1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector pTYB21 Vector	80, 191, 194 45 46 282 377 282 46 46 46 46 46 41 214 378 214 214 379
ProtoScript First Strand cDNA Synthesis Kit. PshAI PsiI pSNAP _t Vector pSNAP _t Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI PstI PstI-HF pTWIN1 Vector pTXB1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector pTYB21 Vector, Map of pUC19 DNA	
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP _t Vector pSNAP _t Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI PstI PstI-HF pTWIN1 Vector pTXB1 Vector, Map of pTXB1 Vector pTXB1 Vector pTYB21 Vector pTYB21 Vector pTYB21 Vector pTYB21 Vector pTYB21 Vector pUC19 DNA pUC19, Map of	
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP _I Vector pSNAP _I Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI PstI PstI-HF pTWIN1 Vector pTXB1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector pTYB21 Vector Map of pUC19 DNA pUC19, Map of PURExpress In Vitro Protein Synthesis Kit	
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP _I Vector pSNAP _I Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspGI PspOMI PspXI PstI PstI-HF pTWIN1 Vector pTXB1 Vector, Map of pTXB1 Vector pTYB21 Vector Map of pUC19 DNA pUC19, Map of PURExpress In Vitro Protein Synthesis Kit	
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP ₁ Vector pSNAP ₂ Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspGI PspDMI PspXI Pstl Pstl-HF pTWIN1 Vector pTXB1 Vector, Map of pTXB1 Vector pTYB21 Vector pTYB21 Vector pTYB21 Vector pTYB21 Vector pTYB21 Vector Map of pUC19 DNA pUC19 MAp of PURExpress In Vitro Protein Synthesis Kit PURExpress Δ (aa, tRNA) Kit PURExpress Δ RF123 Kit	
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP _I Vector pSNAP _I Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspGI PspOMI PspXI PstI PstI-HF pTWIN1 Vector pTXB1 Vector, Map of pTXB1 Vector pTYB21 Vector Map of pUC19 DNA pUC19, Map of PURExpress In Vitro Protein Synthesis Kit	
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP ₁ Vector pSNAP ₂ Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspGI PspDMI PspXI Pstl Pstl-HF pTWIN1 Vector pTXB1 Vector, Map of pTXB1 Vector pTYB21 Vector pTYB21 Vector pTYB21 Vector pTYB21 Vector pTYB21 Vector Map of pUC19 DNA pUC19 MAp of PURExpress In Vitro Protein Synthesis Kit PURExpress Δ (aa, tRNA) Kit PURExpress Δ RF123 Kit	
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP ₁ Vector pSNAP ₂ Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI Pstl Pstl-HF pTWIN1 Vector pTXB1 Vector, Map of pTXB1 Vector pTYB21 Vector pTYB21 Vector pTYB21 Vector pTYB21 Vector Synthesis Kit PURExpress Δ (aa, tRNA) Kit PURExpress Δ Ribosome Kit.	80, 191, 194 45 46 282 377 282 46 46 46 46 214 378 214 214 379 118 380 212, 216 216 216 216
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP ₁ Vector. pSNAP ₁ Vector, Map of pSNAP-tag(T7)-2 Vector PspGl PspOMI PspXI PstI PstI-HF pTWIN1 Vector. pTXB1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector pTYB21 Vec	80, 191, 194 45 46 282 377 282 46 46 46 46 41 378 214 214 379 118 380 212, 216 216 216 216 216 57, 166
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP ₁ Vector pSNAP ₁ Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI PstI PstI-HF pTWIN1 Vector pTXB1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector pTYB21 Vector pUC19 DNA pUC19, Map of PURExpress In Vitro Protein Synthesis Kit PURExpress A RF123 Kit PURExpress A Ribosome Kit PURExpress Disulfide Bond Enhancer Purple Gel Loading Dye Pvul	80, 191, 194 45 46 282 377 282 46 46 46 46 46 214 378 214 214 217 379 118 380 212, 216 216 216 216 216 216 46
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP _I Vector pSNAP _I Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI PstI PstI-HF pTWIN1 Vector pTXB1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector pTYB21 Vector pTYB21 Vector pTYB21 Vector pTYB21 Vector, Map of pUC19 DNA pUC19, Map of PURExpress A (aa, tRNA) Kit PURExpress A RF123 Kit PURExpress Disulfide Bond Enhancer. Purple Gel Loading Dye Pvul Pvul-HF	80, 191, 194 45 46 282 377 282 46 46 46 46 214 378 214 214 379 118 380 212, 216 216 216 216 216 216 46 46 46 46 46 46 46 46 46 46 46 46 46
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP _I Vector pSNAP _I Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI PstI PstI-HF pTWIN1 Vector pTXB1 Vector, Map of pTXB1 Vector pTXB1 Vector pTYB21 Vector	
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP _I Vector pSNAP _I Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspJMI PstI PstI-HF pTWIN1 Vector pTXB1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector pTYB21 Vector pTYB21 Vector, Map of pUC19 DNA pUC19, Map of PURExpress A Rif123 Kit PURExpress A Rif123 Kit PURExpress Disulfide Bond Enhancer Purple Gel Loading Dye Pvul Pvul Pvul PvulI-HF	
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP ₁ Vector. pSNAP ₁ Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI PstI PstI-HF pTWIN1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector, pTYB21 Vector pTYB21 Vector, Map of pUC19 DNA pUC19 DNA pUC19, Map of PURExpress In Vitro Protein Synthesis Kit PURExpress A (aa, tRNA) Kit PURExpress A RF123 Kit PURExpress A Ribosome Kit PURExpress Disulfide Bond Enhancer. Purple Gel Loading Dye Pvul Pvul-HF Pvull-HF Pvull-HF Pyrophosphatase, Inorganic (yeast).	80, 191, 194 45 46 282 377 282 46 46 46 46 214 378 214 214 379 118 380 212, 216 216 216 46 47 47
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP ₁ Vector. pSNAP ₁ Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI PstI PstI-HF pTWIN1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector, Map of pUC19 DNA pUC19 DNA pUC19, Map of PURExpress A Vitro Protein Synthesis Kit PURExpress A RF123 Kit PURExpress A Ribosome Kit. PURExpress Disulfide Bond Enhancer Purple Gel Loading Dye Pvul Pvul-HF Pyrophosphatase, Inorganic (yeast). Pyrophosphatase, Inorganic (yeast). Pyrophosphatase, Inorganic (yeast).	80, 191, 194 45 46 282 377 282 46 46 46 46 214 378 214 214 379 118 380 212, 216 216 216 46 47 47 47
ProtoScript First Strand cDNA Synthesis Kit. PshAl Psil pSNAP ₁ Vector. pSNAP ₁ Vector, Map of pSNAP-tag(T7)-2 Vector PspGI PspOMI PspXI PstI PstI-HF pTWIN1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector, Map of pTXB1 Vector, pTYB21 Vector pTYB21 Vector, Map of pUC19 DNA pUC19 DNA pUC19, Map of PURExpress In Vitro Protein Synthesis Kit PURExpress A (aa, tRNA) Kit PURExpress A RF123 Kit PURExpress A Ribosome Kit PURExpress Disulfide Bond Enhancer. Purple Gel Loading Dye Pvul Pvul-HF Pvull-HF Pvull-HF Pyrophosphatase, Inorganic (yeast).	80, 191, 194 45 46 282 377 282 46 46 46 46 46 214 378 214 214 379 118 380 212, 216 216 216 46 47 47 189 189

Q	
Quality	2
Q5 High-Fidelity DNA Polymerases (Enzyme, Master Mix, Kit)	62, 63, 176
Q5 Hot Start High-Fidelity DNA Polymerases (Enzyme, Master Mix)	62, 63
Q5 Site Directed Mutagenesis Kit (with or without comp cells)	92, 176
Q5U Hot Start High-Fidelity DNA Polymerase	63
qPCR	
Luna Universal qPCR Master Mix	
Luna Universal Probe qPCR Master Mix	70
Optimization Tips	339
Troubleshooting Guide	341
Quick Blunting Kit	
Quick Dephosphorylation Kit	,
Quick Ligation Kit	
Quick-Load DNA Ladders	
Quick-Load Purple DNA Ladders	
Quick-Load <i>Taq</i> 2X Master Mix	66

R	
Random Primers	193
Random Primer Mix	193
Rapid PNGase F Antibody Standard	241
Rapid PNGase F, (non-reducing format)	
Reaction Buffers (Polymerases)	79
Reaction Buffers (Restriction Enzymes)	56
Reagents for DNA Library Preparation	134-161
Reagents for RNA Library Preparation	134-161
REBASE	289
RecA, RecA _r	117
RecJ _r	106
Recleavable Blunt Ends	331
Recleavable Filled-in 5´ Overhangs	330
Recognition Sequence	
Alphabetized List	309-310
Cross Index	304-309
Isoschizomer	311-327
Multiple	306-307
Palindromic Penta-nucleotide	307-309
Recombinant Albumin, Molecular Biology Grade	56
Recombinase, Cre	113
Reference Appendix	286-385
Remove-iT PNGase F.	241
Repair Proteins	76, 107-112
Restriction Endonucleases	
Activity Chart	293-298
Activity of DNA Modifying Enzymes in CutSmart Buffer	299
Activity of Enzymes at 37°C	299
Activity of Restriction Enzymes in PCR Buffers	334-345
Alphabetical Listing	16, 309-310
Cleavage of Supercoiled DNA	329
Cleavage Close to the End of DNA Fragments	
Compatible Cohesive Ends	332-333
Cross Index of Recognition Sequences	304-306
Diluent Buffers.	57, 293
Double Digestion	291
Extended Digests	328-329
Frequency of Sites in Sequenced DNA Molecules	365-368
Generating New Cleavage Sites	330
Heat Inactivation	
High Fidelity (HF)	301
Isoschizomer List	311-327
Methylation-Dependent	264-265
Methylation Sensitivity (Dam, Dcm, CpG)	
Multiple Recognition Sequences	
Nonpalindromic Recognition Sequences	
Optimizing Reactions	
Performance Chart	
Reaction Buffers	n

· ·	330
Recognition Sequences	
Reduced Star Activities of HF Enzymes	
Star Activity	
Survival in a Reaction	
Time-Saver Qualified	
Troubleshooting Guide	
Restriction Maps.	
Reverse Transcription	
AMV Reverse Transcriptase	80, 191, 192
M-MuLV Reverse Transcriptase	, ,
One Taq One-Step RT-PCR Kit	
One Tag RT-PCR Kit	
ProtoScript First Strand cDNA Synthesis Kit	
ProtoScript II First Strand cDNA Synthesis Kit ProtoScript II Reverse Transcriptase	
WarmStart RTx Reverse Transcriptase	
Ribonuclease Contamination	
RNA Depletion	
Ribonucleoside Vanadyl Complex	
Ribonucleotide Solution (Mix, Set)	
Ribosomal RNA Depletion	
Ribosome, E. coli	
RNA 5´ Pyrophosphohydrolase (RppH)	
RNA Cap Structure Analogs	
RNA Library Preparation for NGS	142-148, 160-161, 203
RNA Ligases Reported Activities and Applications	106
RtcB Ligase	
Selection Chart.	
T4 RNA Ligase 1	
T4 RNA Ligase 2	
T4 RNA Ligase 2, truncated	195-197
T4 RNA Ligase 2, truncated K227Q, KQ	
Thermostable 5´ App DNA/RNA Ligase	
RNA Loading Dye, 2X	207
RNA Markers and Ladders	
	470.007
dsRNA Ladder	,
dsRNA Ladder Low Range ssRNA Ladder	170, 207
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases	170, 207 170, 207 170, 207 170, 207
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase.	170, 207 170, 207 170, 207 170, 207 188 188
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases <i>E. coli</i> Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases <i>E. coli</i> Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase T7 RNA Polymerase RNA Reagents	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase T7 RNA Polymerase RNA Reagents RNA Reagents	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase. SP6 RNA Polymerase. T3 RNA Polymerase. T7 RNA Polymerase. RNA Reagents RNAse I, RNase III, Shortcut	170, 207 170, 207 170, 207 170, 207 188 188 188 188 188 188 207 200 201
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase. SP6 RNA Polymerase. T3 RNA Polymerase. T7 RNA Polymerase. RNA Reagents RNAse II, RNase III, Shortcut RNase Contamination, Avoiding	170, 207 170, 207 170, 207 170, 207 188 188 188 188 188 207 200 201 184
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNA Reagents RNAse I, RNase III, Shortcut RNase Contamination, Avoiding RNase H	170, 207 170, 207 170, 207 170, 207 188 188 188 188 188 200 200 201 184
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNA Reagents RNA Reagents RNAse II, RNase III, Shortcut RNase H RNase H	170, 207 170, 207 170, 207 170, 207 188 188 188 188 2188 200 201 184 200 200
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNA Reagents RNAse I, RNase III, Shortcut RNase Contamination, Avoiding RNase H	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNA Reagents RNAse I, RNase II, Shortcut RNase Contamination, Avoiding RNase H RNase HII RNase Inhibitor, Human Placenta	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNA Reagents RNAse I, RNase III, Shortcut RNase Contamination, Avoiding RNAse HII RNase IIIII RNAse IIII	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase T7 RNA Polymerase RNA Reagents RNAse III, Shortcut RNase Contamination, Avoiding RNase H RNase HIII RNAse IIII RNASE IIIIIIIIIIIIIIIIIIIIIIIIIIIIIIIIIIII	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase T7 RNA Polymerase RNA Reagents RNAse I, RNase III, Shortcut RNase Contamination, Avoiding RNase H RNase HIII RNase Inhibitor, Human Placenta RNase Inhibitor, Murine RNase T RNA Library Preparation RNA Polymerase RNA Polymerases RNA Polymerases RNA Polymerases RNA Polymerases RNA Polymerases RNA Polymerases RNA Synthesis	170, 207 170, 207 170, 207 170, 207 188 188 188 188 188 200 201 201 203 202 105, 202 142-161 188-189
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNA Reagents RNAse II, RNase III, Shortcut RNase Contamination, Avoiding RNase H RNase HII RNase Inhibitor, Human Placenta RNAse I Inhibitor, Murine RNASE T RNA Polymerases RNA Polymerase RNA Polymerase RNA Polymerase RNA Synthesis EnGen sgRNA Synthesis Kit	170, 207 170, 207 170, 207 170, 207 188 188 188 188 188 200 201 184 200 200 201 184 185 185 185 187 187 176, 179
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase 37 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNA Reagents RNAse II, RNase III, Shortcut RNase Contamination, Avoiding RNase HII RNase Inhibitor, Human Placenta RNase Inhibitor, Murine RNase T RNA Polymerases RNA Polymerase RNA Polymerase RNASH RNA	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNAse II, Shortcut RNase III, Shortcut RNase Contamination, Avoiding RNase HII RNase Inhibitor, Human Placenta RNase Inhibitor, Murine RNase T RNA Polymerase RNA Polymerase RNA Synthesis EnGen sgRNA Synthesis Kit HiScribe SP6 RNA Synthesis Kit HiScribe T7 ARCA mRNA Kit (with or without Tailing)	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNA Reagents RNase II, Shortcut RNase Contamination, Avoiding RNase H RNase HIII RNase Inhibitor, Human Placenta RNase Inhibitor, Murine RNase T RNA Library Preparation RNA Polymerases RNA Synthesis EnGen sgRNA Synthesis Kit HiScribe SP6 RNA Synthesis Kit HiScribe T7 ARCA mRNA Kit (with or without Tailing) HiScribe T7 Quick High Yield RNA Synthesis Kit	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNA Reagents RNase II, Shortcut RNase Contamination, Avoiding RNase H RNase HIII RNase Inhibitor, Human Placenta RNase Inhibitor, Murine RNase T RNA Library Preparation RNA Polymerases RNA Synthesis EnGen sgRNA Synthesis Kit HiScribe T7 ARCA mRNA Kit (with or without Tailing) HiScribe T7 Quick High Yield RNA Synthesis Kit HiScribe T7 Quick High Yield RNA Synthesis Kit HiScribe T7 High Yield RNA Synthesis Kit	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNA Reagents RNase II, Shortcut RNase Contamination, Avoiding RNase H RNase HIII RNase Inhibitor, Human Placenta RNase Inhibitor, Murine RNase T RNA Library Preparation RNA Polymerases RNA Synthesis EnGen sgRNA Synthesis Kit HiScribe SP6 RNA Synthesis Kit HiScribe T7 ARCA mRNA Kit (with or without Tailing) HiScribe T7 Quick High Yield RNA Synthesis Kit	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNA Reagents RNase II, Shortcut RNase Contamination, Avoiding RNase H RNase HII RNase Inhibitor, Human Placenta RNase Inhibitor, Murine RNase T RNA Polymerases RNA Polymerase RNA Library Preparation RNA Polymerases RNA Polymerases RNA Synthesis EnGen sgRNA Synthesis Kit HiScribe T7 ARCA mRNA Kit (with or without Tailing) HiScribe T7 Quick High Yield RNA Synthesis Kit HiScribe T7 High Yield RNA Synthesis Kit	170, 207 170, 207 170, 207 170, 207 188 188 188 188 188 188 200 200 201 184 200 203 202 105, 202 142-161 188-189 176, 179 185, 187 176, 185, 187 176, 185, 187
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNA Reagents RNAse III, Shortcut RNase III, Shortcut RNase Contamination, Avoiding RNase HI RNase Inhibitor, Human Placenta RNase Inhibitor, Murine RNase T RNA Library Preparation RNA Polymerases RNA Synthesis EnGen sgRNA Synthesis Kit HiScribe T7 ARCA mRNA Kit (with or without Tailing) HiScribe T7 Quick High Yield RNA Synthesis Kit Ribonucleotides Routine PCR LongAmp Taq DNA Polymerase One Taq DNA Polymerase	
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNA Reagents RNAse II, Shortcut RNase III, Shortcut RNase III, Shortcut RNase Contamination, Avoiding RNase H RNase Inhibitor, Human Placenta RNase Inhibitor, Murine RNase T RNA Library Preparation RNA Polymerases RNA Synthesis Engen sgRNA Synthesis Kit HiScribe T7 ARCA mRNA Kit (with or without Tailing) HiScribe T7 Unick High Yield RNA Synthesis Kit HiScribe T7 High Yield RNA Synthesis Kit HiScribe T7 High Yield RNA Synthesis Kit Ribonucleotides Routine PCR LongAmp Taq DNA Polymerase Phusion DNA Polymerase	170, 207 170, 207 170, 207 170, 207 188 188 188 188 188 188 188 200 201 184 200 203 202 105, 202 142-161 188-189 176, 179 185, 187 176, 185, 187 176, 185, 187 176, 185, 187 185 62, 67 62, 65
dsRNA Ladder Low Range ssRNA Ladder microRNA Marker ssRNA Ladder RNA Polymerases E. coli Poly(A) Polymerases Poly(U) Polymerase SP6 RNA Polymerase T3 RNA Polymerase T3 RNA Polymerase T7 RNA Polymerase RNA Reagents RNAse III, Shortcut RNase III, Shortcut RNase Contamination, Avoiding RNase HI RNase Inhibitor, Human Placenta RNase Inhibitor, Murine RNase T RNA Library Preparation RNA Polymerases RNA Synthesis EnGen sgRNA Synthesis Kit HiScribe T7 ARCA mRNA Kit (with or without Tailing) HiScribe T7 Quick High Yield RNA Synthesis Kit Ribonucleotides Routine PCR LongAmp Taq DNA Polymerase One Taq DNA Polymerase	170, 207 170, 207 170, 207 170, 207 188 188 188 188 188 188 188 200 201 184 200 203 202 105, 202 142-161 188-189 176, 179 185, 187 176, 185, 187 176, 185, 187 176, 185, 187 185 62, 67 62, 65 62, 64 62, 63

RppH	199
Rsal	47
RsrII	47
RtcB Ligase	195, 196, 198
RT-PCR	
One Taq One-Step RT-PCR Kit	72
One Taq RT-PCR Kit	72
RT-qPCR	
Luna Universal One-Step RT-qPCR Kit	71
Luna Universal Probe One-Step RT-qPCR Kit	71
Optimization Tips for Luna RT-qPCR	340
Troubleshooting Guide for Luna RT-qPCR	342

S 7-deaza-dGTP Sacl .. SacI-HF SacII. S-adenosylmethionine (SAM) 56 Sall. .47 Sall-HF. Sample Preparation Next Gen Sequencing.. .122-131 Nucleic Acid Purification . Target Enrichment DNA Enrichment... Sapl. Sau3Al Sau96I. Sbfl.. Sbfl-HF .48 Scal-HF ScrFI. SexAl SfaNI Sfcl. .49 Sfil. Sfol. .49 SgrAl. .49 sgRNA Synthesis. Shipping Information6 ShortCut RNase III Shrimp Alkaline Phosphatase (rSAP) .100 .217, 224, 225, 232 SHuffle Competent E. coli (Express, T7 Express). Single-stranded DNA Binding Protein..... siRNA Binding Protein, P19. Site-Directed Mutagenesis Kit, Q5 Smal SmII . ..50 SnaBl SNAP Analysis, Substrates for..... SNAP-Biotin... .282 SNAP-Capture Magnetic Beads .283 SNAP-Capture Pull-Down Resin .280 SNAP-Cell 430... SNAP-Cell 505-star .280 SNAP-Cell 647-SiR.. .280 SNAP-Cell Block. .281 .280 SNAP-Cell Oregon Green.. .279 SNAP-Cell Starter Kit... SNAP-Cell TMR-Star... SNAP-Surface 488.... .280 SNAP-Surface 549. .280 SNAP-Surface 594... .280 SNAP-Surface 649. .280 SNAP-Surface Alexa Fluor 488. .280 SNAP-Surface Alexa Fluor 546. .280 SNAP-Surface Alexa Fluor 647 SNAP-Surface Block... .281

SNAP-tag	
Comparison to GFP	279
Labeling with SNAP-tag/Troubleshooting Guide	364
Overview	278
SNAP-tag Purified Protein	281
Starter Kit	279
Vectors	282
SNAP-Vista Green	282
Sodium Orthovanadate (Vanadate).	256
SP6 RNA Polymerase	188
Spel	50
Spel-HF	50
Sphl	50
SphI-HF.	50
SplintR Ligase.	199
Srfl	50
Sspl	50
SspI-HF.	50
ssRNA Ladder	170, 207
ssRNA Ligase (T4 RNA Ligase 1)	195, 196, 197
SssI (CpG) Methyltransferase	114
Standards, Protein	171
Star Activity	300
Streptavidin Magnetic Beads	205
Stul	51
Styl	51
Styl-HF	51
StyD4I	
Sulfolobus DNA Polymerase IV (Lesion Bypass Polymerase)	
Supercoiled DNA, Cleavage of	
Supercoiled DNA Ladder	
Survival in a Reaction	328-329
Swal 64	
Synthetic Biology	
BioBrick Assembly Kit	
Gibson Assembly	
Golden Gate Assembly	
USER Enzyme	91

-	
T3 DNA Ligase	96
T3 RNA Polymerase	
T4 DNA Ligase	
T4 DNA Polymerase	
T4 Gene 32 Protein	
T4 PDG (T4 Endonuclease V)	111
T4 Phage β-glucosyltransferase (T4-BGT)	262
T4 Polynucleotide Kinases	98
T4 RNA Ligase 1 (ssRNA Ligase)	195, 196, 197
T4 RNA Ligase 2 (dsRNA Ligase)	195, 196, 197
T4 RNA Ligase 2, truncated	195, 196, 197
T4 RNA Ligase 2, truncated K227Q and KQ	195, 196, 198
T4 RNA Ligase Selection Chart	196
T7 DNA Ligase	96
T7 DNA Polymerase	78
T7 Endonuclease I	113, 178
T5 Exonuclease	107
T7 Exonuclease	107
T7 Express Competent E. coli, SHuffle	217, 224, 225, 232
T7 Express (High Efficiency), Competent E. coli	
T7 Express <i>lysY</i> Competent <i>E. coli</i> , SHuffle	
T7 Express, <i>lysY</i> (High Efficiency), Competent <i>E. coli</i>	
T7 Express, <i>lysY/l^q</i> (High Efficiency), Competent <i>E. coli</i>	
T7 Express Strains, Protein Expression with	
T7, Map of	
T7 RNA Polymerase	
Taq [~] l	
Taql Methyltransferase	
Taq 2X Master Mix	
Taq 2X Master Mix, Quick-Load	66

Tag 5X Master Mix	66
Taq DNA Ligase	97
Tag DNA Polymerase with Standard Tag Buffer	
Tag DNA Polymerase with Standard Tag (Mg-free) Buffer	66
Tag DNA Polymerase with ThermoPol Buffer	66
Tag PCR Kit	66
Tag Reaction Buffers	79
Target Enrichment	
Technical Support	286
TelN Protelomerase	113
Template Switching RT Enzyme Mix	
Terminal Transferase	
TEV Protease	255
Tfil	
Therminator DNA Polymerase	
Thermolabile Exonuclease I	
Thermolabile USER II Enzyme.	
Thermolabile Proteinase K.	
Thermophiles at 37°C, Activity of	
Thermophilic Helicase Dependent Amplification	
Thermophilic Polymerases, Diluent Buffers	
Thermophilic Polymerases, Reaction Buffers	
Thermostable 5´ App DNA/RNA Ligase	
Thermostable DNA Ligase (<i>Taq</i>)	, ,
Thermostable FEN1	
Thermostable Inorganic Pyrophosphatase	
Thermostable Ligase Reaction Temperature Calculator	
Thermostable RNase H.	
Time-Saver Qualified Restriction Enzymes	
Tm Calculator	
Tma Endonuclease III	
Topoisomerase	
Cre Recombinase	113
Topoisomerase I	
Transformation Efficiencies, Enhancing	
TriDye DNA Ladders	
Troubleshooting Guides	
Cloning	349-351
Genomic DNA Purification	
Luna qPCR	
Luna RT-qPCR	
DNA Cleanup & Plasmid Purification	
PCR	
Restriction Enzymes	
SNAP-tag Technology	
Total RNA Extraction & Purification	
Trypsin-digested BSA MS Standard (CAM-modified)	
Trypsin-ultra, Mass Spectrometry Grade	
Tsel	
Tsp45I	
TspMI.	
TspRI	
Tte UvrD Helicase	
Tth1111	
Tth Endonuclease IV	

U	
UDG	111, 112
Universal miRNA Cloning Linker	
Unstained Protein Standard, Broad Range (10-200 kDa)	
Uracil-DNA Glycosylase (UDG)	
Uracil-DNA Glycosylase (UDG), Afu	
USER Enzyme	
	•
V	
Vaccinia Capping System	190
(Vanadate) Sodium Orthovanadate	256
Vectors	
Cloning	118
IMPACT	
K. lactis	215
pMAL	
SNAP/CLIP	
Vent DNA Polymerase	
Vent (exo ⁻) DNA Polymerase	
Vista Label	
W	
WarmStart Colorimetric LAMP Master Mix (DNA & RNA)	7.1
WarmStart LAMP Kit (DNA & RNA)	
WarmStart RTx Reverse Transcriptase	
Water, Nuclease-free	
vvalei, ivuolease-iiee	
X	
Xbal	
Xcml	
Xhol	
Xmal	
Xmnl	
XRN-1	201
V	
Υ	
Yeast Medium Pack	215
Yeast Expression	215
yDcpS	
•	
7	
L	

Zral.

.53



Environmental stewardship is one of the founding principles of NEB. We feel it is important to protect and preserve our natural resources, including biodiversity. In addition to the content found in this catalog, we invite you to take a look at other NEB resources that address this topic.

Visit www.neb.com/EnvironmentalPhilosophy and download our Social and Environmental Responsibility Brochure.





Visit **www.neb.com/NEBTV** view our NEB TV episodes that discuss the conservation of biodiversity (Episode 25) and green laboratory practices (Episode 22).

Visit www.neb.com/NEBpodcast and listen to our podcast featuring a group of sustainable science experts that participated in a panel discussion during our Go Green Symposium (Episode 12).





Visit **LabConscious.com**, a community dedicated to promoting green lab initiatives to reduce waste, use green chemistry, conserve water and save energy. Join the discussion!



For 45 years, the New England Biolabs catalog has been a resource for scientists around the world. This catalog features a collection of mini-reviews that discuss conservation of biodiversity. As part of this effort, NEB will offer support to the organization below that protects biodiversity.



Reforestation and preservation activities:

Safe Habitat for Chimpanzees – Project 1169

Kibale National Park is one of the last remaining tropical forest blocks in Uganda. It harbours the greatest variety of primates found anywhere in East Africa. It is home to almost 1,500 chimpanzees, Uganda's largest population of this endangered species. The variety of further species, like forest elephants, wild cats, birds and plants represent the intact biodiversity.

Through reforestation and preservation activities, the project helps protect biodiversity and mitigate climate change on a total area of 10,000 hectares.

To learn more please visit: https://fpm.climatepartner.com/project/1169/en







Environmental Philosophy & the NEB Catalog

The NEB Catalog & Technical Reference is printed with sustainability in mind.

- The catalog is printed on Forest Stewardship Council (FSC) Certified, recycled paper.
- . The printing facility is also FSC Certified.
- The facility uses environmental-friendly inks and coatings

As with previous catalogs, we have achieved "ClimatePartner" certification by working with the ClimatePartner Company, a world-leading provider of carbon reduction solutions. This means that the unavoidable CO₂ generated by the printing and distribution of this catalog have been reduced to net zero through verified carbon offset projects. Offsets will be used to support Safe Habitat for Chimpanzees.

More information on carbon neutral production as well as details about the Biomass project Nr 1067 supported by this NEB catalog can be found at www.climatepartner.com.

Using 30% post-consumer fiber instead of virgin fiber for this catalog resulted in the following:

- 760 trees saved
- 3,000 lbs of solid waste not generated
- 160,000 gallons of water conserved
- 430,000 pounds of greenhouse gas emissions avoided
- 340,000,000 British Thermal Units (BTUs) of energy not consumed

Note: Statistics and calculations are based on the Environmental Paper Network Paper Calculator v3.2 and US-print production.





DISTRIBUTORS

ARGENTINA MIGLIORE LACLAUSTRA S.R.L. Tel. (11) 4372 9045 Email: info@migliorelaclaustra.com.ar www.migliorelaclaustra.com.ar

AUSTRALIA

NEW ENGLAND BIOLABS (AUSTRALIA) PTY See back cover

BAHRAIN

NEW ENGLAND BIOLABS Tel. +44(0)1462 420616 Email: exports.uk@neb.com

BELGIUM/LUXEMBURG/ THE NETHERLANDS

BIOKÉ

Tel. +31 (0)71 720 0220 Tel. (Be) Toll Free 0800 71640 Email: info@bioke.com www.bioke.com

BRAZIL

UNISCIENCE DO BRAZIL Tel. 11 3622 23 20 Email: info@unisciencce.com www.uniscience.com

VERITAS BIOTECNOLOGIA Tel. +55 11 2338 1016 Email: vendas@veritasbio.com.br

BULGARIA

Elta 90M Ltd. Tel. (+359) 2 983 96 49 Email: office@elta90.eu www.elta90.com

CAMEROON

INOABA BIOTEC WEST AFRICA LTD.
Tel. +234 803 403 5281-3
Email: info@inqababiotec.ng
www.inqababiotec.ng

CANADA

NEW ENGLAND BIOLABS, LTD. See back cover

CHILE

GENE X-PRESS Tel. +56 22 2454805 Email: info@genex.cl www.genexpress.cl

CHINA

NEW ENGLAND BIOLABS See back cover

COLOMBIA

SUMINISTROS CLINICOS ISLA SAS Tel. 571-6250580 Email: scliente@islasas.com www.islasas.com

CROATIA

DIAGNOSTICA SKALPELI d.o.o. Tel. +385 1 3778 484 Email: skalpeli@skalpeli.hr www.skalpeli.hr

CZECH REPUBLIC

BIOTECH a.s. Tel. +420 272 701 739 Email: info@ibiotech.cz www.ibiotech.cz

DENMARK BIONORDIKA DENMARK A/S Tel. +45 3956 2000 Email: info@bionordika.dk www.bionordika.dk

EGYPT SIGMA SCIENTIFIC SERVICES CO. Tel. +2 01006604079 Email: sigma.s@hotmail.com www.sigmaeg-co.com

FINLAND

BIONORDIKA OY Tel. +358 (0)207 410 270 Email: info@bionordika.fi www.bionordika.fi

FRANCE

NEW ENGLAND BIOLABS FRANCE See back cover

GHANA

INQABA BIOTEC WEST AFRICA LTD. Tel. +233 30 396 9327 Email: akwaaba@inqababiotec.ng www.inqababiotec.ng

GREECE

BIOLINE SCIENTIFIC Tel. 0030 210 5226547 Email: info@bioline.gr

GERMANY/AUSTRIA New England Biolabs GmbH See back cover

HONG KONG

BIOARROW TECHNOLOGY LTD. Tel. +852 2179-7168/2179-7193

HUNGARY

Email: info@kvalitex.hu www.kvalitex.hu

ICELAND

MEDOR EHF Tel. +354-412-7000 Email: medor@medor.is www.medor.is

INDIA

IMPERIAL LIFE SCIENCE (P) Tel. +91-124-4559800-99 info@imperialls.com www.imperialls.com

INDONESIA

PT. GENECRAFT LABS
Tel. +62-21-5890-3119
Email: sales@genecraftlabs.com
www.genecraftlabs.com

IRAQ

NEW ENGLAND BIOLABS Tel. +44(0)1462 420616 Email: exports.uk@neb.com

IRELAND, REPUBLIC OF

BRENNAN AND COMPANY Tel. +353 1 295 2501 Email: enquiries@brennanco.ie www.brennanco.ie

ISRAEL

ORNAT BIOCHEMICALS & LABORATORY EQUIPMENT LTD.
Tel. +972-8-9477077
Email: ornatbio@ornat.co.il
www.ornat.co.il

ITALY EUROCLONE S.P.A. Tel. +39 02 38.19.51 Email: info.bio@euroclone.it www.euroclonegroup.it

JAPAN

NEW ENGLAND BIOLABS JAPAN, INC. See back cover

JORDAN

BURJOUAN EST. FOR MEDICAL & LAB SUPPLIES Tel. +962-6-5350949 Email: info@burjouan.com www.burjouan.com

KENYA

INQABA BIOTEC EAST AFRICA LTD. Tel. +254 73 537 0693 Email: info@inqababiotec.co.ke www.inqababiotec.co.ke

KOREA, REPUBLIC OF KORAM BIOTECH CORP.

Tel. 82-2-556-0311 Email: koram@korambiotech.com www.korambiotech.com

Tel. +82 2 21057020 Email: info@philekorea.co.kr www.philekorea.co.kr

KUWAIT NEW ENGLAND BIOLABS Tel. +44(0)1462 420616 Email: exports.uk@neb.com

I FRANON

NEW ENGLAND BIOLABS Tel. +44(0)1462 420616 Email: exports.uk@neb.com

LIBYA

NEW ENGLAND BIOLABS Tel. +44(0)1462 420616 Email: exports.uk@neb.com

NEXT GENE SCIENTIFIC SDN BHD. Tel. +603-8051 5106 Email: sales@nextgene.my

MEXICO

QUIMICA VALANER S.A. DE C.V. Tel. (52) 5525 5725 Email: servicioaclientes@valaner.com www.valaner.com

NEW ZEALAND GENESEARCH PTY LTD Freephone: 0800 870 424 Tel. +64 9 887 0424 Email: orders@genesearch.co.nz www.genesearch.co.nz

NIGERIA

INQABA BIOTEC WEST AFRICA LTD. Tel. +234 805 882 7272 Email: info@inqababiotec.ng www.inqababiotec.ng

NORWAY

BIONORDIKA AS Tel. +47 23 03 58 00 Email: info@bionordika.no www.bionordika.no

OMAN NEW ENGLAND BIOLABS Tel. +44(0)1462 420616 Email: exports.uk@neb.com

PAKISTAN PAKISTAN HOSPITAL AND INDUSTRIAL LABORATORY CO. (PHILCO) Tel. 0092-42-37324276 Email: philcopakistan@yahoo.com www.philcopak.com

POI AND

Tel. +48 22 3359884 Email: labjot@labjot.com www.labjot.com

PORTUGAL

WERFEN Email: dac-pt@werfen.com www.pt.werfen.com

QATAR

QATAR SCIENTIFIC CO. WLL Tel.+974 4462 1165 Email: info@qsneoscience.com

ROMANIA

ELTA 90 MEDICAL RESEARCH SRL Tel. +40/21/2322694 Email: office@elta90mr.ro www.elta90mr.ro

RUSSIA and CIS

SKYGEN Tel. +7 (495) 215-02-22 Email: info@skygen.com www.skygen.com

SAUDI ARABIA (KSA)

NEW ENGLAND BIOLABS Tel. +44(0)1462 420616 Email: exports.uk@neb.com

SERBIA

ALFA GENETICS d.o.o.
Tel. + 381 11 208 4252
Email: office@alfagenetics.rs
www.alfagenetics.rs/

SENEGAL

INQABA BIOTEC WEST AFRICA LTD. Tel: +234 805 882 7272 Email: info@inqababiotec.ng www.inqababiotec.ng

SINGAPORE

NEW ENGLAND BIOLABS, PTE. LTD. See back cover

SLOVAK REPUBLIC

BioTech s.r.o. +421 2 54 77 44 88 info@ibiotech.sk http://www.ibiotech.sk

<u>SLO</u>VENIA

Mikro+Polo d.o.o Tel. +386 2 614 33 00 info@mikro-polo.si www.mikro-polo.si

SOUTH AFRICA

Inqaba Biotechnical Industries (Pty) Ltd. Tel. +27 12 343 5829 Email: info@inqababiotec.co.za www.inqababiotec.co.za

SPAIN

Tel. +34 902 20 30 90 Email: customerservice-es @werfen.com www.es.werfen.com

SWEDEN

BIONORDIKA SWEDEN AB Tel. +46 8 30 60 10 Email: info@bionordika.se www.bionordika.se

SWITZERLAND

Tel. +41 (0)61 486 80 80 Email: info@bioconcept.ch www.bioconcept.ch

TAIWAN
TAIGEN BIOSCIENCE CORP. Tel. (02) 28802913 Email: order@taigen.com www.taigen.com

TANZANIA

INQABA BIOTEC EAST AFRICA LTD Tel. +254 73 537 0693 Email: info@inqababiotec.co.ke www.inqababiotec.co.ke

THAILAND

Tel. (02) 412-5672 Email: info@theetrad.com www.theetrad.com

TURKEY EKA BIOLAB TEKNOLOJI Tel. +90/312 248 13 24 Email: info@ekabiolab.com www.ekabiolab.com

UNITED ARAB EMIRATES

NEW ENGLAND BIOLABS Tel. +44(0)1462 420616 Email: exports.uk@neb.com

UNITED KINGDOM
NEW ENGLAND BIOLABS See back cover

URUGUAY

BIKO SA Tel. +598 2204 27 43 Email: importaciones@biko.com.uy www.biko.com.uy

VENEZUELA

BIOPRODUCTS, INC. dba CORPOINTER Tel. 58-212-763-5141 Email: jmleal@bioproducts.net mchinea@bioproducts.net mlamanna@bioproducts.net www.bioproducts.net

VIFTNAM

BCE VIETNAM Tel. 84 4 3234 5666 Email: tamnt@bcevietnam.com www.bcevietnam.com.vn

YEMEN NEW ENGLAND BIOLABS Tel. +44(0)1462 420616 Email: exports.uk@neb.com

USA

New England Biolabs, Inc.
240 County Road
Ipswich, MA 01938-2723
Telephone (978) 927-5054
Toll Free (USA Orders) 1-800-632-5227
Toll Free (USA Tech) 1-800-632-7799
Fax (978) 921-1350
email: info@neb.com

AUSTRALIA

New England Biolabs (Australia) PTY Telephone: +61 (1800) 934218 email: info.au@neb.com www.nebiolabs.com.au

CANADA

New England Biolabs, Ltd.
Telephone (905) 665-4632
Toll Free 1-800-387-1095
Fax (905) 665-4630
Fax Toll Free 1-800-563-3789
email: info.ca@neb.com
www.neb.ca

CHINA

New England Biolabs (Beijing), Ltd. Telephone 010-82378265/82378266 Fax 010-82378262 email: info@neb-china.com www.neb-china.com

FRANCE

New England Biolabs France Telephone +33 (0) 1 69 90 87 15 Free Call 0800 100 632 Free Fax 0800 100 610 email: info.fr@neb.com www.neb-online.fr

GERMANY/AUSTRIA

New England Biolabs GmbH
Telephone +49 (0) 69/305 23140
Free Call 0800/246 5227 (Germany)
Free Call 00800/246 52277 (Austria)
Fax +49 (0) 69/305 23149
Free Fax 0800/246 5229 (Germany)
email: info.de@neb.com
www.neb-online.de

JAPAN

New England Biolabs Japan, Inc. Telephone +81 (0) 3 5669 6191 Fax +81 (0) 3 5669 6192 email: info.jp@neb.com www.nebj.jp

SINGAPORE

New England Biolabs, Pte. Ltd. Telephone +65 638 59623 Fax +65 638 59617 email: sales.sg@neb.com www.neb.sq

UNITED KINGDOM

New England Biolabs (UK), Ltd. Telephone (01462) 420616 Call Free 0800 318486 Fax (01462) 421057 Fax Free 0800 435682 email: info.uk@neb.com www.neb.uk.com

WANT TO EXPERIENCE THE NEB CATALOG IN AUGMENTED REALITY?



DOWNLOAD

the NEB Augmented Reality (AR) app for iPhone® or iPad® at the Apple® App Store or for Android™ on Google Play™.



EXPLORE

the NEB Catalog & Technical Reference and keep an eye out for the augmented reality icon.



DISCOVER

videos, tutorials and immersive experiences by scanning the page using the AR app.

